

INVERTA™ Implants Product Catalogue





Southern Implants is a leading provider of unique and innovative dental implant products with a focus on top-end professional users who want more choices. Southern's expertise in research, development and manufacturing of dental implants allows us to provide Innovative Treatment Solutions that will reduce treatment times and improve patient outcomes.

Striving for excellence and meeting customer needs, has led to our wide product range characterized by Unique and Innovative products which include;

- Multiple interfaces, to suit customer preference.
- INVERTA™ implant, featuring a body-shift design, engineered for primary stability and suitable for immediate loading.
- Co-Axis®, sub-crestal angle correcting implant, available in angulations of 12, 24 & 36° and various internal and external connections.
- MAX implant, specifically designed for immediate molar tooth replacement.
- The ZYGAN and ZYGEX implants for severely resorbed maxilla and craniofacial reconstruction.

Our product portfolio is in synchronized evolution with protocol improvements and technological advances.

 $My \, since re\, thanks\, to\, all\, specialists, dentists\, and\, technicians\, who\, put\, their\, trust\, in\, our\, company.$

Graham Blackbeard

Managing Director, Southern Implants

DEEP CONICAL - INVERTA™

Deep Conical INVERTA Ø3.5mm - Ø4.5mm Implants Ø4.0mm - Ø5.0mm Implants Surgical Information Prosthetics Information Site Preparation. Instrumentation.	-
CO-AXIS® Deep Conical INVERTA Ø3.5mm - Ø4.5mm 12° Co-Axis Implants Ø4.0mm - Ø5.0mm 12° Co-Axis Implants Surgical Information Prosthetics Information Site Preparation Instrumentation	Page 12 Page 12 Page 12 Page 13 Page 14 Page 15
EXTERNAL HEX - INVERTA	
External Hex INVERTA Ø3.5mm - Ø4.5mm Implants Ø4.0mm - Ø5.0mm Implants Surgical Information. Prosthetics Information Site Preparation Instrumentation.	Page 18 Page 18 Page 18 Page 19 Page 20 Page 21
CO-AXIS External Hex INVERTA Ø3.5mm - Ø4.5mm 12° Co-Axis Implants Ø4.0mm - Ø5.0mm 12° Co-Axis Implants Surgical Information. Prosthetics Information. Site Preparation. Instrumentation.	Page 24 Page 24 Page 25 Page 26 Page 27
Torque Table for Southern Screws.	Page 28
Driver Information	Page 29
Instrument Trays I-IV-EG.	Page 30
Implant dimensions and information Deep Conical INVERTA. Implant dimensions and information External Hex INVERTA.	Page 32 Page 33
Additional Instrument Information	Page 34
Explanation of Symbols.	Page 35



 ${\sf NOTE:} \ \ \, \textbf{--Images are for illustration purposes only and do not necessarily accurately represent the product.}$

- All dimensions in this catalogue are in mm, unless otherwise specified.
- Not all products are cleared for sale in all countries.

DEEP CONICAL INVERTATION

Straight

Ø3.5mm-Ø4.5mm

Ø4.0mm-Ø5.0mm

INVERTA[™] Deep Conical Ø3.5mm-Ø4.5mm



Implants are available in lengths of:

ITEM CODE	IMPLANT LENGTHS (in mm) L1
IV-DC35-4510	10.0
IV-DC35-4511	11.5
IV-DC35-4513	13.0
IV-DC35-4515	15.0

Coronal Ø D1	Maximum Ø D2
3.6mm	4.5mm

INVERTA Deep Conical Ø4.0mm-Ø5.0mm



Implants are available in lengths of:

ITEM CODE	IMPLANT LENGTHS (in mm) L1
IV-DC40-5010	10.0
IV-DC40-5011	11.5
IV-DC40-5013	13.0
IV-DC40-5015	15.0


Coronal Ø D1	Maximum Ø D2
4.0mm	5.0mm


NOTE: Implant dimensions and information - page 32


Surgical Components




Prosthetic Flowchart







Ø4.0mm-Ø5.0mm Tapered (IV-DC40-50)



NOTE: Site preparation sequence recommended by Southern Implants does not replace the judgement and experience of the surgeon.

Drill Information





1.22 Universal Hex Driver

PLACEMENT TOOL CONVERTERS I-WI-CST I-WI-SS For Handpiece inserts (Latch-type) featuring the WV&H hex I-HD-22U-M IV-TR-45xx IV-TR-50xx IV-TR-50xx

11/13/15

11/13/15

NOTE: - Refer to CAT-1191 for more information on wrench insert converters.
- Most instruments are available in various lengths




DEEP CONICAL INVERTA TM

Co-Axis®

Ø3.5mm-Ø4.5mm

Ø4.0mm-Ø5.0mm

INVERTA[™] Deep Conical Co-Axis[®] 12° Ø3.5mm-Ø4.5mm




Implants are available in lengths of:

ITEM CODE	IMPLANT LENGTHS (in mm)	
	L1	L2
IV-DC3512D-4511	11.5	10.86
IV-DC3512D-4513	13.0	12.36
IV-DC3512D-4515	15.0	14.36

Coronal Ø D1	Maximum Ø D2
3.6mm	4.5mm

INVERTA Deep Conical Co-Axis 12° Ø4.0mm-Ø5.0mm



Implants are available in lengths of:

ITEM CODE	IMPLANT LENGTHS (in mm)	
	L1	L2
IV-DC4012D-5011	11.5	10.86
IV-DC4012D-5013	13.0	12.36
IV-DC4012D-5015	15.0	14.36


Coronal Ø D1	Maximum Ø D2
4.0mm	5.0mm


NOTE: Implant dimensions and information - page 32

Surgical Components




Prosthetic Flowchart





Ø3.5mm-Ø4.5mm Tapered Co-Axis® (IV-DC3512d-45)

(illustrations are for 13mm implants)



*Final Tapered Drill Position

PLEASE NOTE:

Point 1


This corner of the drill is to be at bone level.

Point 2

This corner of the drill will be subcrestal.



Ø4.0mm-Ø5.0mm Tapered Co-Axis® (IV-DC4012d-50)



*Final Tapered Drill Position


PLEASE NOTE:

Point 1

This corner of the drill is to be at bone level.

Point 2

This corner of the drill will be subcrestal.



 $NOTE: Site \ preparation \ sequence \ recommended \ by \ Southern \ Implants \ does \ not \ replace \ the \ judgement \ and \ experience \ of \ the \ surgeon.$

Drill Information





PLACEMENT TOOL





I-WI-DC4-S/M





For Handpiece inserts (Latch-type) featuring the W&H hex

CONVERTERS

I-WI-SS



For SQUARE connection of fixture mounts and instruments



I-WI-SH

For HEX connection of fixture mounts

DRIVERS

TRY-IN / TAPS

I-HD-22U-M



1.22 Universal Hex Driver

IV-TR-45xx IV-TR-50xx







11/13/15

NOTE: - Refer to CAT-1191 for more information on wrench insert converters.

- Most instruments are available in various lengths



EXTERNAL HEX INVERTA TM

Straight

Ø3.5mm-Ø4.5mm

Ø4.0mm-Ø5.0mm

INVERTA[™] External Hex Ø3.5mm-Ø4.5mm



Implants are pre-mounted and available in lengths of:

ITEM CODE	IMPLANT LENGTHS (in mm) L1
IV-EX35-4510	10.0
IV-EX35-4511	11.5
IV-EX35-4513	13.0
IV-EX35-4515	15.0

Coronal Ø D1	Maximum Ø D2
3.5mm	4.5mm

INVERTA External Hex Ø4.0mm-Ø5.0mm



Implants are pre-mounted and available in lengths of:

ITEM CODE	IMPLANT LENGTHS (in mm) L1
IV-EX40-5010	10.0
IV-EX40-5011	11.5
IV-EX40-5013	13.0
IV-EX40-5015	15.0

Coronal Ø D1	Maximum Ø D2
4.0mm	5.0mm


NOTE: Implant dimensions and information - page 33

Surgical Components

COVER SCREW HEALING ABUTMENTS TBN WBN Ø3.6 Ø4.5 12/3/4/6/8 lengths 2/3/4/6 lengths


NOTE: TORQUE TO 20Ncm

Prosthetic Flowchart




Ø3.5mm-Ø4.5mm Tapered (IV-EX35-45)

(illustrations are for 13mm implants)



Ø4.0mm-Ø5.0mm Tapered (IV-EX40-50)



NOTE: Site preparation sequence recommended by Southern Implants does not replace the judgement and experience of the surgeon.

Drill Information





PLACEMENT TOOL CONVERTERS

I-CON-X

Fixture Mount Driver



I-WI-CST

For Handpiece inserts (Latch-type) featuring the W&H hex.



I-WI-SH

For HEX connection of fixture mounts



I-WI-SS

For SQUARE connection of fixture mounts and instruments

DRIVERS TRY-IN / TAPS



Cover Screw Driver 0.9 Hex







NOTE: - Refer to CAT-1191 for more information on wrench insert converters.
- Most instruments are available in various lengths



EXTERNAL HEX INVERTA TM

Co-Axis®

Ø3.5mm-Ø4.5mm

Ø4.0mm-Ø5.0mm

INVERTA[™] External Hex Co-Axis[®] 12° Ø3.5mm-Ø4.5mm




Implants are pre-mounted and available in lengths of:

ITEM CODE	IMPLANT LE	
IV-EX3512D-4511	11.5	10.78
IV-EX3512D-4513	13.0	12.28
IV-EX3512D-4515	15.0	14.28

Coronal Ø D1	Maximum Ø D2
3.5mm	4.5mm

INVERTA External Hex Co-Axis 12° Ø4.0mm-Ø5.0mm



Implants are pre-mounted and available in lengths of:

ITEM CODE	IMPLANT LE L1	NGTHS (in mm) L2
IV-EX4012D-5011	11.5	10.72
IV-EX4012D-5013	13.0	12.22
IV-EX4012D-5015	15.0	14.22

Coronal Ø D1	Maximum Ø D2
4.0mm	5.0mm


NOTE: Implant dimensions and information - page 33

Surgical Components




NOTE: TORQUE TO 20Ncm

Prosthetic Flowchart




Ø3.5mm-Ø4.5mm Tapered Co-Axis® (IV-EX3512d-45)


(illustrations are for 13mm implants)





Ø4.0mm-Ø5.0mm Tapered Co-Axis® (IV-EX4012d-50)





NOTE: Site preparation sequence recommended by Southern Implants does not replace the judgement and experience of the surgeon.

Drill Information





PLACEMENT TOOL CONVERTERS

I-CON-X

For Handpiece inserts (Latch-type) featuring the W&H hex.

I-WI-CST

I-WI-SH

For HEX For SQUARE connection of fixture mounts fixture mounts and instruments

I-WI-SS

DRIVERS TRY-IN / TAPS

Fixture Mount Driver



NOTE: - Refer to CAT-1191 for more information on wrench insert converters.
- Most instruments are available in various lengths

EXTERNAL HEX

2 SERIES Z Screws 1.22 Hex (M2)

TORQUE: 32-40Ncm Head Diameter: 2.70mm

TSHZ2

Laboratory Screws



BSH2*

3 SERIES Z Screws 1.22 Hex (M2)



TORQUE: 32-40Ncm **Head Diameter:** 2.40mm

TSHZ3

Note: Screw TORQUE with PEEK Prosthetics for < Ø4.0mm implant interfaces torque to 15Ncm ≥ Ø4.0mm implant interfaces torque to 20Ncm

BSH3

Laboratory Screw

DC (DEEP CONICAL)

1.22 Hex



DC4 TORQUE: 20Ncm **Head Diameter:** 2.35mm

TS-DC4-16

Laboratory Screw



Note: Screw TORQUE with PEEK Prosthetics: DC4 interface 20Ncm

1.22 Hex Passive Abutment Screws



TORQUE: 20Ncm Head Diameter: 2.40mm

PA-DC4-16B* (*anodized blue)

ANGLED COMPACT CONICAL ABUTMENT SCREWS

DC (Deep Conical)

1,22 Hex



TS-DC4-MC

TORQUE: 20Ncm **Head Diameter:** 2.0mm

External Hex

Series 9Z Screws (M2)

Unigrip



TSUZ9

GSUZ9

TORQUE: 20Ncm Head Diameter:

Digital Laboratory Analoge Screw

1.22 Hex



TORQUE: Finger tighten

LAD-S

NOTE: Screw supplied with all Digital Analogues

1 SERIES SCREW OPTIONS

Unigrip

1 22 Hex



TORQUE: 10-15Ncm **Head Diameter** 2.25mm

TSU1 GSU1 TSH₁ GSH1 Laboratory Screw

1.22 Hex



BSH1* *(Blackened and for laboratory use only)

Screw Head Connections

Hex



Unigrip



- NOTE:
- Due to design revisions screw tips may be flat or rounded.
- Always ensure that the correct screw is used for the relevant implant and component.
- Refer to CAT-1003 for alternative slotted 1 Series screws

Drivers



Abutment Drivers (for Compact Conical Abutments)



Cover Screw Drivers

0.9 HEX




Drivers for 1 Series Screws (for Compact Conical Prosthetics)


1.22 HEX



NOTE: Refer to CAT-1203 for alternative drivers

I-IV-EG for surgical placement of INVERTA Implants






The surgical kit has an intuitive layout to guide the surgeon through the drill sequence.


All instruments and tooling used during the procedure must be maintained in good condition, and cleaned and sterilized prior to to use. Please consult the Southern Implants Cleaning and Sterilization Procedure Guidelines (CAT-1039) for guidance concerning the maintenance of drills, instruments, and surgical trays.

											F 4 4		F
RANGE										N N	IMPLANI LENGIRS	ב ב ב	2
		CORONAL DIAMETER	MAXIMUM DIAMETER	PROSTHETICS	CONNECTION DIAMETER	THREAD PITCH	APEX DIAMETER	BUCCAL PLATFORM	PLATFORM ANGLE	10	11,5	13	15
IV-DC35-45	Ø3.5mm-Ø4.5mm	3.6	4.5		2.8	9.0	2.6	0.40	.0	>	>	>	>
IV-DC40-50	Ø4.0mm-Ø5.0mm	4.0	5.0		2.8	9.0	2.6	09'0	.0	>	>	>	>
IV-DC3512D-45	Ø3.5mm-Ø4.5mm	3.6	4.5		2.8	9'0	2.6	0.52	12°		>	^	>
IV-DC4012D-50	Ø4.0mm-Ø5.0mm	4.0	2'0		2.8	9'0	2.6	0.72	12°		>	^	^



Not all products are cleared for sale in all countries.

RANGE										IMPI	IMPLANT LENGTHS	LENG	THS
		CORONAL	MAXIMUM DIAMETER	Prosthetics	CONNECTION DIAMETER	THREAD PITCH	APEX DIAMETER	BUCCAL PLATFORM	PLATFORM ANGLE	10	11.5	13	15
IV-EX35-45	Ø3.5mm-Ø4.5mm	3.5	4.5		3,4	9'0	2.6	0.04	°0	>	7	>	7
IV-EX40-50	Ø4.0mm-Ø5.0mm	4.0	5.0		3.4	9'0	2.6	0.32	00	>	>	>	7
IV-EX3512D-45	Ø3.5mm-Ø4.5mm	3.5	4.5		3,4	9'0	2.6	0.12	12°		٨	٨	٨
IV-EX4012D-50	Ø4.0mm-Ø5.0mm	4.0	2.0		3,4	9'0	2.6	0.52	12°		٨	٨	٨



BONE MILLS

DC (Deep Conical)

CODES I-HBM-DC35 I-HBM-DC40

NOTE: \bullet Refer to CAT-1009 for more information on Bone Mills

For Polishing Protector Caps refer to CAT-1010

EXTERNAL HEX

HANDPIECE



I-BNM-45 I-BNM-55

CODES

Bone Mill Protector Cap

I-BNM-CAP

The following symbols are used on packaging labels and they indicate the following:

- 1 Manufacturer
- 2 Colour code indicating platform diameter
- 3 Implant image
- 4 Implant details and size
- 5 STERILE R Sterilization using Irradiation
 - Do not Resterilize
 - Consult instruction for use
 - ② Do not reuse
 - **C E** CE mark and notified body number
 - **Expiry date**
 - Sterile unless package is opened or damaged
- 6 2D Bar coding Contains the GTIN, Expiry Date and LOT Number
- 7 Patient sticker for documentation purposes (to be used by health care provider on patient file)
- 8 Prescription device

CAUTION: FEDERAL LAW RESTRICTS THE DEVICE TO SALE BY OR ON THE ORDER OF A LICENCED HEALTH CARE PROVIDER.









For more information, please contact your Southern Implants Representative or visit southernimplants.com



South Africa - Headquarters 1 Albert Road, Irene, RSA **T:** +27-12-667-1046 | **E:** info@southernimplants.com

Subsidiaries

Australia

Southern Implants Australia T: +61-2-8076-9337 E: info@southernimplants.com.au **United Kingdom**

Southern Implants UK
T: +44-20-8899-6845 / 6 / 7
E: info@southernimplants.co.uk

USA and Canada

Southern Implants North America Inc. T: +1-561-472-0990 E: customercare@southernimplants.com

Southern Implants Products are distributed world-wide, please visit southern implants.com for a list of Distributors.

©2019 Southern Implants (Pty) Ltd. All rights reserved. Co-Axis, Southern Implants, and Subcrestal Angle Correctionare registered trademarks and Body-Shift, INVERTA and PROVATA are trademarks of Southern Implants (Pty) Ltd. Not all implants are available in all markets.