

SIDIGITAL Product Catalogue





SIDigital is Southern Implants' comprehensive digital solution for CAD/CAM procedures.

Precision is one of the most important aspects of any dental restoration, and this becomes more significant when it comes to digital dentistry. Southern's digital solutions provide digital tools which accurately facilitate the digital workflow between dentist and laboratory for CAD/CAM procedures. This is achieved through the use of digital libraries and our unique abutment options, to cater for single unit and multiple unit restorations, all in the digital workflow.

Dental professionals wanting to work with SIDigital workflow procedures, will need to download and install the SIDigital libraries into either 3Shape, Exocad or Dental Wings software systems. These digital libraries link the intra-oral/desktop scan, model creator, design software and milling of the restoration.

My sincere thanks to all specialists, dentists and technicians who give continual feedback, suggestions and input. Our products are an interpretation of your needs.

Graham Blackbeard

Managing Director, Southern Implants

Introduction	Page 05
SIDigital Libraries and Abutment Selection	
	Page 10
scan flags	-
Lab Analogues	Page 16
ABUTMENT SELECTION	Faye 20
CIA Abutments	Page 26
CAB Abutments	Page 28
TIB Abutments	Page 30
Procedures for Lise	
CIA CAB and TIB Abutments using SIDigital libraries	Page 35
TIB Abutments used with Sirona	Page 37
Passive Abutments.	Page 38
Procedure for Lise	
Procedule for use	Page /1
1 assive Abuthent	i aye 4 i
Final Milling	Page 46
Southern Implants CAD/CAM workflow for Exocad 3Shane Dental Wings & Sirona	
External Hey	Page 50
M-Series (Internal Hex)	Page 62
PROVATA [™]	Page 66
DC (Deen Conical)	Page 72
TBI-NEX	Page 80
IT (Internal Octagon) Range	Page 90
Driver & Screw Information	Page 94
Instrument Information	Page 95
Polishing Protector Caps.	Page 96
Instructions for Use	Page 97
Explanation of Symbols	Page 99

Please note: • Images are for illustration purposes only and do not necessarily accurately represent the product.

All dimensions in this catalogue are in mm, unless otherwise specified.

· Not all products are cleared for sale in all countries.

INTRODUCTION

SIDIGITAL is Southern Implants' digital offering for use in CAD/CAM procedures.

Precision is the most important aspect in any dental restoration and this is more relevant with SIDigital. Southern's digital solutions provide digital tools which accurately facilitate the digital workflow between dentist and laboratory during CAD/CAM procedures. This is achieved with the use of digital libraries and Southern Implants' unique abutment options, which allows for single and multiple unit restorations all in a digital workflow.

Dental professionals working with digital workflow procedures are required to download and install SIDigital libraries into 3Shape, Excocad, and/or Dental Wings software systems. Digital libraries are the link between the intraoral / desktop scan, model creator, design software, and restoration milling.





SIDIGITAL LIBRARIES

SIDigital libraries are available for all dental professionals and are free to download. A digital workflow with digital libraries makes it easier to get consistent and accurate results when designing restorations. This is due to the digital file being the link to the different steps in the workflow, and the components in the library are exact replicas of the physical parts. By using the libraries, you are enabled to create complete dental restorations with an optimal fit.

- The libraries are available for 3Shape, Exocad, and/or Dental Wings.
- SIDigital libraries need to be imported into these systems.

Procedure to register online:

- Go to www.southernimplants.com
- Click on SIDidital solutions
- Scroll down and click "CAD/CAM library Access & Login"



- · Complete the registration form while taking care to provide the correct email address to receive updates.
- · Accept terms and conditions, and submit.

NB! Only register ONCE!

🗇 Southern Implants, Dental Ir 🗖 Southern Implants CA	D, × +	-	٥	\times
$\label{eq:constraint} \leftarrow \ \rightarrow \ \heartsuit \ \ \ \textcircled{\ \ \ } \ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \$	□ ☆ =	\mathbb{Z}	۵	
	Email Access Code:		L	igin 🍈
	~ S >			
	SOUTHERNIMPLANTS			
	CAD/CAM LIBRARIES			
	You need an access code to download our CAD/CAM libraries. You can request a code by Illing in the form below. The information will be reviewed by our expresentatives who will send you an access code once approved. If you already have an access code you can login on the top right of this page.			
	CUSTOMER SERVICE			
	sidigital (Seo Aversinglants.com			
	NOTE			
	Do not register more than once. This causes issues with data capture. Keep your usemane and password safe and use it for all further downloads.			
	PASSWORD REQUEST			
	Personal Information:			
	First Name: Las Name: Prode: Prod: Prod: Prode: Prode: Prode: Prode: Prode: Prode: Pro			
	In which country do you practice?			
	Country not listed? What year of practice do you have? Other? Which CADRAN eyderms do you use? Any system ratilisted? Which Can What year which you use?			
	DIGITIAL DATA LICENSE AGREEMENT			~

The dental professional will then be provided with their unique "password" or "access code" as login details (via e-mail).

Procedure to download libraries:

- To download the libraries visit www.southernimplants.com and follow the same steps as above to get to the registration page.
- · At the top right corner enter your email address and access code obtained from SIDigital and login.

🗇 Southern Implants, Dental Ir 🗧 Southern Implants CA	D ₂ × +	- ø ×
$\leftarrow \rightarrow$ O $ $ \triangle si-cadcam.appspot.com		
	TS S	Access Code:
	SOUTHERNIMPLANTS	
	CAD/CAM LIBRARIES	
	You need an access code to download our CADICAM Renaries. You can request a code by Bling in the fram below. The information will be reviewed by our representatives who will send you an access code once approved. If you already have an access code you can logit on the up offort of this gap.	
	CUSTOMER SERVICE	
	sidghai@exutheminglants.com	
	NOTE	
	Do not register more than once. This causes issues with data capture. Keep your usemame and password safe and use it for all further downloads.	
	PASSWORD REQUEST	
	Personal Information:	
	First Name: Last Name: Last Name: Practice / Company: Email: Email: Practice Practice Practice Date of Registration (Vyythmixit)	
	Practice Information:	
	In which country do you practice? Country not issen? Either type of practice do you have? Deve? Uther CADICAM systems do you use? Any systems mail issen? Either of our implant systems do you use?	
	DIGITIAL DATA LICENSE AGREEMENT	

Library installation instructions:



• Select either 3Shape, Exocad or Dental Wings.

NOTE: USA users must download the 3Shape FDA encrypted library



· SIDigital libraries are packaged as ".exe extension" type since 3Shape, Exocad, and Dental Wings all have different file types.

· Open Southern Implants and click on the library to download. (If a screen comes up prompting if the library should be run click "Run anyway".)

😼 Setup - SI 3Shape Library —	×
License Agreement Please read the following important information before continuing.	BOUTHERN MPLANTS
Please read the following License Agreement. You must accept the terms of this agreement before continuing with the installation.	
DIGITAL DATA LICENSE AGREEMENT	^
Between	
Southern Implants Pty (Ltd) 1 Albert Road Irene	~
 I accept the agreement I do not accept the agreement 	
Next >	Cancel

Accept the agreement, and click next.

Information Please read the following important information before continuing.	BOUTHER IMPLANT
When you are ready to continue with Setup, click Next.	
Changelog: 2018.06.13 R 1.2 -ADDITION OF INT-HEX WIDE CONNECTIONS -PASSIVE SCREW CHANNEL WIDENED 2018.05.07 R 1.1 -ADDITION OF NEW TIB's -TIB ANTI-ROTATION MADE TIGHTER 2018.03.22 R 1 -CIA NE MATCH STOCK 2017.12.06 - TIB & CIA BUG FIXED	~

A log of changes and updates made to the libraries appears, continue by clicking next.

SIDIGITAL LIBRARIES

👸 Setup - SI 3Shape Library - 🗆 🗙
Select Components Which components should be installed?
Select the components you want to install; dear the components you do not want to install. Click Next when you are ready to continue.
Full Installation
< Back Next > Cancel

• The next window will ask for full installation, continue by clicking next.

Setup is now ready to begin installing SI 3Shape Library on your computer.		157
Click Install to continue with the installation, or click Back if you want to review or change any settings. Setup type: Full Installation Selected components: Full Installation Full Installation V	hape Library on your computer.	BOUTHERN IMPLANTS
Setup type: Full Installation Selected components: Full Installation <	, or click Back if you want to revie	ew or
Selected components: Full Installation		^
Full Installation		
<		
<		
< >		
		~
		>
		hape Library on your computer. , or click Back if you want to revi

· Install the library.

· Click on Finish

L Windows (Ci)				
	PVSW			
	Southern Implants	18/10/2017 10:11	File folder	

- The library will download on the local disk "C:" drive, in a folder labelled "Southern Implants".
- When installing the library into the software of choice, download the library from this "Southern Implants" folder on the local disk "C:" drive.

Note:

- **Do not register more than once on the website.** If you have, please email <u>sidigital@southernimplants.com</u> for assistance.
- You can re-download the library as many times as you want with your access code. This code does not expire and can be used numerous times.

DIGITALLY OBTAINING THE IMPLANT/ANALOGUE POSITION

Intended use:

Southern Implants' scan flags are intended to be used during the digital impression procedure, in order to obtain the exact position and orientation of the respective dental implant or laboratory analogue, and transfer this position digitally during CAD/CAM scanning procedures. This aids the CAD/CAM software to align the CAD/CAM components to the implant in a digital format.

Description:

The scan flags are pre-manufactured and are available in different connection types. The Scan flags are made from Titanium and is supplied with a screw to secure the scan flag onto the dental implant or analogue.

EXTERNAL HEX

	Order Number	to be used with	Driver
	SFT-EX-30	Ø3mm External Hex implant (Piccolo)	
	SFT-EX-34	Ø3.25mm External Hex implant (IBN / IBNT)	
2	SFT-EX-40	Ø4mm External Hex implant (IB / IBT) & MAX-6	
	SFT-EX-50	Ø5mm External Hex implant & MAX-7	I-HD-22U-M (1.22/1.27 Universal)
	SFT-EX-60	Ø6mm External Hex implant & MAX-8	

TRI-NEX

	Order Number	to be used with	Driver
5 	SFT-EL-35	Ø3.5mm TRI-NEX implant	
st Lu l'	SFT-EL-43	Ø4.3mm TRI-NEX implant	
	SFT-EL-50	Ø5mm TRI-NEX implant & TRI-MAX-7	I-UGI-M (Unigrip)
	SFT-EL-60	Ø6mm TRI-NEX implant & TRI-MAX-8 & 9	

SCAN FLAGS

DCRA	NGE		
	Order Number	to be used with	Driver
	SFT-DC3	Ø3mm Deep Conical implant	
	SFT-DC4	Ø3.5mm & Ø4mm Deep Conical implant	
	SFT-DC5	Ø5mm Deep Conical implant	I-HD-M (1.22 hex)
	NAL HEX RAN	IGE (M-SERIES & PROVATA®)	
	Order Number	to be used with	Driver

	Ø3.7mm / Ø4.2mm / Ø5mm Internal Hex implants (M-series)	(TORCOLA)
SFT-M	Ø4mm / Ø5mm Provata	
	Ø6mm PROMAX implants	
SFT-Z	Ø7mm / Ø8mm / Ø9mm PROMAX implants	I-HD-27-M (1.27 hex)

IT (INTERNAL OCTAGON) RANGE

C	Order Number	to be used with	Driver
e la constante de la constante	SFT-IT	Ø4.8mm IT implants & MAXIT7	
	SFT-IT6	Ø6.5mm IT implants & MAXIT8 & 9	I-SCS-M (Torx)

ABUTMENT LEVEL (COMPACT CONICAL ABUTMENT)

Procedure for use:

The scan flag can be scanned directly on the implant with an intraoral scanner or indirectly on a model. Scan flags are used with digital libraries available for 3Shape, Exocad, and/or Dental Wings. A scan flag is attached to the implant or analogue before scanning to provide the exact orientation of the dental implant in the software.

- · that there is no damage.
- Attach the correct Scan flag to the dental implant, abutment or lab analogue. Check that a proper fit is obtained, and hand tighten the screw with the appropriate driver.
- The patient is scanned with the scan flag attached to an implant, using an intraoral scanner. The laboratory model is scanned with the scan flag attached to the laboratory analogue, using a desktop scanner.
- The scan flag is removed from the implant or analogue.
- The scan flag in the digital format is now matched and aligned with the corresponding scan flag from the library file.
- The software recognizes the orientation of the scan flag to the implant or analogue. This allows the software to know where to place the abutment for the design stage.

Note: Follow the instructions of the scanner used, and for scanning procedures.

Materials:

Scan flag: Titanium Screw: Titanium Grade 4 (anodized)

FOR 3D PRINTED OR MILLED MODELS

Intended use

Southern Implants' digital laboratory analogues are intended to be used as dental implant replicas that has to be inserted into a 3D printed or a milled model in order to duplicate the location, orientation, and restorative platform of the implant placed in the mouth during CAD/CAM procedures.

Description

Digital laboratory analogues are pre-manufactured components and are available in different connections. The items are to be fitted into a 3D printed or a milled model. Digital analogues are manufactured from Titanium, and are a two-piece component. The main component is anodized to match the prosthetic diameter of the implant being restored, and is supplied with a screw to secure the seat of the analogue in the model.

NOTE: Southern Implants digital analogues are for single use only. Re-use of the analogue can result in loss of accuracy.

EXTERNAL HEX

	Order Number	Scan flag	to be used with	Driver
, a	LAD-IP	SFT-EX-30	Ø3.0mm External Hex implant	
8	LAD-IBN	SFT-EX-34	Ø3.25mm External Hex implant	
	LAD-IB	SFT-EX-40	Ø4mm External Hex implant & MAX-6	
y å	LAD-BA	SFT-EX-50	Ø5mm External Hex implant & MAX-7	I-HD-M (1.22 hex)
	LAD-BBB	SFT-EX-60	Ø6mm External Hex implant & MAX-8	

TRI-NEX

	Order Number	Scan flag	to be used with	Driver
	LAD-L-35	SFT-EL-35	Ø3.5mm TRI-NEX implant	154283
	LAD-L-43	SFT-EL-43	Ø4.3mm TRI-NEX implant	
	LAD-L-50	SFT-EL-50	Ø5mm TRI-NEX implant & TRI-MAX-7	
a a	LAD-L-60	SFT-EL-60	Ø6mm TRI-NEX implant & TRI-MAX-8 & 9	۱۰۰تا طח-۱۱ (1.22 hex)

DC RANGE

 Order Number	Scan flag	to be used with	Driver
LAD-DC3	SFT-DC3	Ø3mm Deep Conical implant	18572-1
LAD-DC4	SFT-DC4	Ø3.5mm & Ø4mm Deep Conical implant	
LAD-DC5	SFT-DC5	Ø5mm Deep Conical implant	I-HD-M (1.22 hex)

INTERNAL HEX RANGE (M-SERIES & PROVATA®)

	Order Number	Scan flag	to be used with	Driver
			Ø3.7mm / Ø4.2mm / Ø5mm Internal Hex implants (M-series)	10002-1
	LAD-M	SFT-M	Ø4mm / Ø5mm Provata	U
			Ø6mm PROMAX implants	
Ц К Å	LAD-M-P45	SFT-M	Ø5mm Provata & Ø6mm PROMAX implants	I-HD-M (1.22 hex)
	LAD-Z	SFT-Z	Ø7mm / Ø8mm / Ø9mm PROMAX implants	

IT (INTERNAL OCTAGON) RANGE

	Order Number	Scan flag	to be used with	Driver
1	LAD-ITS	SFT-IT	Ø4.8mm IT implants & MAXIT7	
	LAD-IT6	SFT-IT6	Ø6.5mm IT implants & MAXIT8 & 9	I-HD-M (1.22 hex)

ABUTMENT LEVEL (COMPACT CONICAL ABUTMENT)

	Order Number	Scan flag	to be used with	Driver
	LAD-MC	SFT-MC-48	Ø4.8mm platform Compact conical abutments	
í de la comercia de l	LAD-MCW	SFT-MC-60	Ø6mm platform Compact conical abutments	I-HD-M

(1.22 hex)

Procedures:

- 1. The implant position is obtained digitally by an intraoral scan of the patient with the scan flag attached to the implant.
- 2. Remove the scan flag from the implant/s, and replace with healing abutment/s.
- 3. The scan is then imported into the design software.

Design procedures:

- 4. The scan flag in the digital form is now matched and aligned with the corresponding scan flag in the library.
- 5. Model design software is chosen and the digital laboratory analogue is positioned digitally in the model.
- 6. The software will guide the user through the steps to design and complete the model.

Milling Procedures:

- 7. The digital stl file of the model is sent to a 3d printer or a milling machine to print/mill.
- 8. The digital analogues are inserted from the top into the open channel of the 3d printed/milled model. A soft tissue model is always recommended.
- 9. The model is turned upside down, and the screw is inserted apically and finger tightened with a 1.22mm hex driver to the analogue. This will secure the analogue into the printed or milled model.

The digital analogue is inserted into the channel from the top of the model.

The model is turned upside down and a hex driver (I-HD-S/M/L) is used to secure the analogue to the model with the analogue screw.

Completed soft tissue model with apical screws inserted into the model.

Notes:

- The digital analogue must be used with the corresponding digital analogue in the libraries.
- The screw must not be torqued only finger tightened.

Material:

Digital lab analogue: Titanium Digital lab analogue screw: Stainless steel

The Concept

The CIA Abutment was originally designed for use during direct scanning procedures with CAD/CAM systems as a screw channel could not be added during the design procedures (thus the need for the PEEK pin). In modern systems, this function is readily available and the PEEK pin can be discarded. The CIA Abutment features a unique narrow platform with a sandblasted surface to aid in retention during cementation procedures.

EXTERNAL HEX

	Order Nu	ımber		to be used with	Scan flag	Prosthetic screw options	Retaining screw torque
Engaging		Non-Engaging				1.22 Hex	
	CIA-EX-34		CIA-NX-34	Ø3.25mm External Hex implant	SFT-EX-34		
	CIA-EX-40		CIA-NX-40	Ø4mm External Hex implant & MAX-6	SFT-EX-40	3 Series Titanium & Gold Alloy screw	32-40Nom
	CIA-EX-50		CIA-NX-50	Ø5mm External Hex implant & MAX-7	SFT-EX-50	(Refer to Catalogue for codes & Instrument detail)	32-40NGIII
	CIA-EX-60		CIA-NX-60	Ø6mm External Hex implant & MAX-8	SFT-EX-60		

TRI-NEX

	Order Nu	mber		to be used with	Scan flag	Prosthetic screw options	Retaining screw torque
Engaging		Non-Engaging				Unigrip	
	CIA-EL-35		CIA-NL-35	Ø3.5mm TRI-NEX implant	SFT-EL-35	TS-L-18 TS-L-18C GS-L-18	32Ncm
	CIA-EL-43		CIA-NL-43	Ø4.3mm TRI-NEX implant	SFT-EL-43		
	CIA-EL-50		CIA-NL-50	Ø5mm TRI-NEX implant & TRI-MAX-7	SFT-EL-50	TS-L-20 TS-L-20C GS-L-20	32-40Ncm
PT P	CIA-EL-60	-	CIA-NL-60	Ø6mm TRI-NEX implant & TRI-MAX-8 & 9	SFT-EL-60		

	Order Nu	umber		to be used with	Scan flag	Prosthetic screw options	Retaining screw torque
Engaging		Non-Engaging				1.27 Hex (Titanium) Torx (Gold)	
1950h		8101		Ø3.7mm / Ø4.2mm / Ø5mm Internal Hex implants (M-series)			
	CIA-EM		CIA-NM	Ø4mm / Ø5mm Provata	SFT-M		
				Ø6mm PROMAX implants		TS-Z-18 GS-Q-18	32Ncm
	CIA-EZ		CIA-NZ	Ø7mm / Ø8mm / Ø9mm PROMAX implants	SFT-Z		

INTERNAL HEX RANGE (M-SERIES & PROVATA®)

IT (INTERNAL OCTAGON) RANGE

	Order Nu	ımber		to be used with	Scan flag	Prosthetic screw options	Retaining screw torque
Engaging		Non-Engaging				Torx	
<u>I</u>	CIA-ITS		CIA-ITS-NE	Ø4.8mm IT implants & MAXIT7	SFT-IT	TSIT2	22 40Nom
4	CIA-ITS6		CIA-ITS6-NE	Ø6.5mm IT implants & MAXIT8 & 9	SFT-IT6	GSIT2	32-40NCM

Note: refer to product catalogue for available instruments and screw types.

Description

CIAAbutments are pre-manufactured components which are available in engaging and non-engaging versions for different connection systems. The CIAAbutment is made from Titanium and includes a PEEK pin, to aid in creating screw channels with certain CAD software during design procedures. These abutments are connected directly to an endosseous implant for use as an aid in prosthetic rehabilitation. Engaging versions have an anti-rotational "lobe" and are designated for single unit cases.

Material:

CIA Abutment: Titanium Titanium Scanning Pin: Polyetheretherketome (PEEK) Abutment screws: Titanium alloy (Ti-90%, AI-6%, V-4%) or Gold alloy (Au-61%; Ag-16.5%; Pt-13.5%; Cu-9%)

The Concept

The CAB Abutment was designed for Internal Hex (M-Series & Provata) and DC implant systems. These newer implant systems came in during the age of Digital dentistry. The concept of the CAB abutment is to be used with CAD/CAM systems where the screw access hole can be created in the design steps. This eliminates the need for a PEEK scanning pin as with the CIAAbutment.

DC RANGE

	Order Number	to be used with	Scan flag	Abutment screw options	Final screw torque
Engaging				1.22 Hex	
ŧ	CAB-DC3	Ø3mm Deep Conical implant	SFT-DC3	TS-DC3-14	15Ncm
ŧ	CAB-DC4	Ø3.5mm & Ø4mm Deep Conical implant	SFT-DC4	TS-DC4-16	20Ncm
	CAB-DC5	Ø5mm Deep Conical implant	SFT-DC5	TS-DC5-20	25-32Ncm

INTERNAL HEX RANGE (M-SERIES & PROVATA®)

Order Number				to be used with	Scan flag	Abutment screw options	Final screw torque
Engaging		Non-Engaging				1.27 Hex (Titanium) Torx (Gold)	
_				Ø3.7mm / Ø4.2mm / Ø5mm Internal Hex implants (M-series)		× 7	
Ļ	CAB-M	4	CAB-NM	Ø4mm / Ø5mm Provata	SET	TS-Z-18	
				Ø6mm PROMAX implants	3F 1-IM		32Ncm
Ļ	CAB-M-P45	5		Ø5mm Provata & Ø6mm PROMAX implants		GS-Q-18	
4	CAB-Z		Ø7mm /	Ø8mm / Ø9mm PROMAX implants	SFT-Z		

Description

CAB Abutments are pre-manufactured components which are available in engaging and non-engaging versions for the standard platform internal hex abutments. The CAB Abutment is made from Titanium and anodized yellow. The abutment has retention grooves to aid in retention during cementation procedures. These abutments are connected directly to an endosseous implant, aiding in prosthetic rehabilitation.

Material:

CAB Abutment:Titanium or Titanium alloy (Ti-6AL-4V)Abutment screws:Titanium alloy Ti-90%, Al-6%, V-4% or Gold alloy Au-61%, Ag-16.5%, Pt-13.5%, Cu -9%

The Concept

The TIB Abutment was designed for use with SIRONA (CEREC) systems. The TIB Abutment range is also used with digital libraries in 3Shape, Exocad, and/or Dental Wings. This makes the TIB Abutment the most extended CAD/CAM abutment solution at Southern Implants. Each abutment is available in 3 different collar heights as well as being available in engaging and non-engaging versions.

Note: Screws sold separately for all abutments. TIB Abutments are not available for Ø3 DC implants and Ø3mm External Hex implants.

EXTERNAL HEX

0.6mm	Order Number Collar Height	3mm	to be used with	Scan flag	3rd party Scan Body	Abutment screw	3rd party Software Selection
Engaging	1.000	omm			2003		
TIB-EX-34	TIB-EX-34-C1.5	TIB-EX-34-C3	(2) 25mm External Llaw			TSHZ3 (3 Series screws. Refer to CAT-2020 for codes & instrument	
Non-Engaging TIB-NX-34	TIB-NX-34-C1.5	TIB-NX-34-C3	implants	3F1-LX-34	L	information) Final torque	DO 3.4L
						32-40Ncm	
Engaging TIB-EX-40	TIB-EX-40-C1.5	TIB-EX-40-C3	Ø4mm External Hex implant & MAX-6	SFT-EX-40	L	TSHZ3 (3 Series screws. Refer to CAT-2020 for codes & instrument	BO 4.1L
Non-Engaging TIB-NX-40	TIB-NX-40-C1.5	TIB-NX-40-C3				information) Final torque 32-40Ncm	
Engaging TIB-EX-50	TIB-EX-50-C1.5	TIB-EX-50-C3	Ø5mm External Hex implant & MAX-7	SFT-EX-50	L	TSHZ3 (3 Series screws. Refer to CAT-2020 for codes & instrument information	BO 5.0L
TIB-NX-50	TIB-NX-50-C1.5	TIB-NX-50-C3				Final torque 32-40Ncm	
Engaging TIB-EX-60	TIB-EX-60-C1.5	TIB-EX-60-C3	Ø6mm External Hex implant & MAX-8	SFT-EX-60		TSHZ3 (3 Series screws. Refer to CAT-2020 for codes & instrument	
Non-Engaging TIB-NX-60	TIB-NX-60-C1.5	TIB-NX-60-C3			L	information) Final torque 32-40Ncm	NB RS 6.0L

TRI-NE	ΞX						
0.6mm	Order Number Collar Height 1.5mm	3mm	to be used with	Scan flag	3rd party Scan Body	Abutment screw	3rd party Software Selection
Engaging TIB-EL-35	TIB-EL-35-C1.5	TIB-EL-35-C3	Ø3.5mm TRI-NEX implant	SFT-EL-35	L	TS-L-18 TS-L-18C GS-L-18 Final torque 32Ncm	NB RS 3.5L
Engaging TIB-EL-43	TIB-EL-43-C1.5	TIB-EL-43-C3	Ø4.3mm TRI-NEX implant	SFT-EL-43	L	TS-L-20 TS-L-20C GS-L-20 Final torque 32-40Ncm	NB RS 4.3L
Engaging TIB-EL-50	TIB-EL-50-C1.5	TIB-EL-50-C3	Ø5mm TRI-NEX implant & TRI-MAX-7	SFT-EL-50	L	TS-L-20 TS-L-20C GS-L-20 Final torque 32-40Ncm	NB RS 5.0L
Engaging TIB-EL-60	TIB-EL-60-C1.5	TIB-EL-60-C3	Ø6mm TRI-NEX implant & TRI-MAX-8 & 9	SFT-EL-60	L	TS-L-20 TS-L-20C GS-L-20 Final torque 32-40Ncm	NB RS 6.0L

0.6mm	Order Number Collar Height 1.5mm	3mm	to be used with	Scan flag	3rd party Scan Body	Abutment screw	3rd party Software Selection
Engaging							
TIB-DC4	TIB-DC4-C1.5	TIB-DC4-C3	03.5 & 04.0 DC implants	SET DC4	I	TS-DC4-16	AT OS 3 5/4 01
Non-Engaging				3F1-DC4	L	Final taxaua	AT US 3.3/4.0L
TIB-NDC4	TIB-NDC4-C1.5	TIB-NDC4-C3				Final torque	
						20Ncm	
Engaging							
TIB-DC5	TIB-DC5-C1.5	TIB-DC5-C3	GE 0 DO invites (057.005		TS-DC5-20	
Non-Engaging			05.0 DC Impiants	SFT-DC5	L	Final torque	AT OS 4.5/5.0L
TIB-NDC5	TIB-NDC5-C1.5	TIB-NDC5-C3				25-32Ncm	

DC-RANGE (DEEP CONICAL)

INTERNAL HEX RANGE (M-SERIES & PROVATA®)

0.6mm	Order Number Collar Height 1.5mm	3mm	to be used with	Scan flag	3rd party Scan Body	Abutment screw	3rd party Software Selection
Engaging TIB-M	TIB-M-C1.5	TIB-M-C3	Ø3.7mm / Ø4.2mm / Ø5mm Internal Hex implants (M-series)			TS-Z-18	
TT.		1	Ø4mm/Ø5mm Provata	SFT-M	L	GS-Q-18	NB RS 4.3L
TIB-NM	TIB-NM-C1.5	TIB-NM-C3	Ø5mm Provata & Ø6mm PROMAX implants			Final torque 32Ncm	
Engaging TIB-M-P45	j		Ø5mm Provata & Ø6mm PROMAX implants	SFT-M	L	TS-Z-18 GS-Q-18 Final torque 32Ncm	NB RS 4.3L
Engaging TIB-Z	TIB-Z-C1.5	TIB-Z-C3	Ø7mm / Ø8mm / Ø9mm PROMAX implants	SFT-Z	L	TS-Z-18 GS-Q-18	B0 5.0L
Non-Engaging TIB-NZ	TIB-NZ-C1.5	TIB-NZ-C3				Final torque 32Ncm	

IT (INTERNAL OCTAGON)

Description

TIB Abutments are pre-manufactured components which are available in engaging and non-engaging versions for different connection systems. The TIB Abutment is made from Titanium and is anodized yellow. It is available in a 0.6mm, 1.5mm, and 3mm collar height. The abutment has retention grooves to aid in retention during cementation procedures. These abutments are connected directly to an endosseous implant for aiding in prosthetic rehabilitation.

Material:

TiB Abutment:Titanium or Titanium alloy (Ti-6AL-4V)Abutment screws:Titanium alloy (Ti-90%, Al-6%, V-4%) or Gold alloy (Au-61%; Ag-16.5%; Pt-13.5%; Cu-9%)

CIA, CAB, and TIB Abutments using SIDigital libraries

- · Download Southern Implants' digital library for 3Shape, Exocad, and/or Dental Wings by registering on www.southernimplants.com
- · Load the libraries into the CAD/CAM system by following their instructions.
- The implant position is obtained digitally either by an intraoral scan of the patient with scan flags attached to the endosseous implant, or a desktop scan of the dental model with the scan flags attached to the laboratory analogue.
- Remove the scan flag from the implants or from the model. (Replace it with healing abutments if used intra-orally).
- The scan is then imported into the design software.
- The scan flag in the digital form is now matched and aligned with the corresponding scan flag in the library.

- The software recognizes the position of the Scan flag to the implant/s or analogue/s as well as where to place the analogues/abutment for the design steps.
 - For a printed/milled model, model design software is used and the digital laboratory analogues are positioned digitally in the model.
- The software will guide the user through the steps to design and complete the model.
- The digital analogues are then inserted into the model.

The abutment/s is chosen from the library. The user will determine the desired tooth or teeth positions needed for the abutment/s, and the path of insertion will be determined. This will govern if a screw retained or cemented restoration is possible.

- The software will guide the user through the steps to complete the restoration. Selection of the relevant restorative material will then be required; the most common material to use is Zirconia.
- The restoration design is completed and the design file is sent to a milling unit or a production facility.
- The milled component/s is to be sintered to the specifications of the manufacturer if Zirconia is used. Once sintered, the abutment, the custom abutment or restoration is cemented to the prefabricated Titanium abutment by closing the screw channel temporarily to keep the screw channel free of cement.

Finishing procedures:

Polishing Protector Caps

These items are used to protect the interface and edges of various components, milled structures, and abutments when polishing. They are attached to the relevant abutments using the appropriate laboratory screw. Polishing then commences according to standard laboratory procedures. (Refer to page 96 for information available for Polishing Protector Caps).

Final Clinical procedures:

- 1. Clean and disinfect the restoration as applicable per the restorative material manufacturer's instructions.
- 2. a. For screw retained unit/s: Place and tighten the restoration screw. Verify the correct seating of the restoration using radiographic imaging. Tighten the restoration, using a torque wrench, to the specified torque value for the applicable prosthetic screw.
 - b. For cement retained unit/s: Place and tighten the restoration screw. Verify the correct seating of the substructure/custom post using radiographic imaging. Tighten the abutment with the substructure/custom post, using a torque wrench, to the specified torque value for the applicable prosthetic screw. Close the screw access holes and cement the final prosthesis in place.
TIB abutments used with Sirona

Scanning procedure:

- 1. Use the corresponding third party Scan body as indicated on pages 30 to 33 or the Worklfow from page 49.
- 2. Choose either the white or grey Scan body depending on the Cerec scanner (refer to Cerec manuals for more detail).
- 3. Align the guide grooves in the Scan body with the TIB Abutment and make sure that it seats flush without any gaps.





- 4. Select the correct scan body and corresponding TiBase in the software.
- 5. Take a scan with your scanner like normal. Ensure that the top of the scan body is completely captured.
- 6. Note: The diameter of the TIB Abutment must not be reduced. Shortening of the TIB is not permitted.
- 7. The software will guide the user to design the shape of the restoration.

Milling Procedures:

- 1. Mill the shape from meso block size L(see the" InLab 3D for Abutment User Manual").
- 2. Carefully follow the instructions and information in the operating instructions for meso blocks.
- 3. Send the design files to a milling unit or production facility accepting Sirona designs, or send the design to your InLab specialist.



The Concept:

The Passive Abutment concept allows one to achieve predictable passive fit of milled structures in a practical and repeatable way and thus eliminates the need for complex and intensive laboratory procedures usually undertaken to improve the fit. The Passive fit is achieved by luting a pre-machined Titanium interface component into the finished prosthesis, using the laboratory master model as the blueprint for fit. No additional clinical steps are required.

EXTERNAL HEX

Order Number		to be used with	Scan flag	Abutment screw	Final torque
Engaging SP16	Non-Engaging SP17	Ø3mm External Hex implants	SFT-EX-30	TS-P-16	25-32Ncm
100 A	(B)				
Engaging SBN16	Non-Engaging SBN-17-TT	Ø3.4mm External Hex implants	SET-EX-34		
A	20073				
Engaging	Non-Engaging			2 Series screws	
SB16	SB-17-TT	Ø4mm External Hex implant & MAX-6	SFT-EX-40	(Refer to CAT-2020	32-40Ncm
A	10013	·		information)	
Engaging	Non-Engaging				
SBA16	SBA-17-TT	Ø5mm External Hex implant & MAX-7	SFT-EX-50		
	LEA.				
Engaging	Non-Engaging				
SBBB16	SBBB-17-TT	Ø6mm External Hex implant & MAX-8	SFT-EX-60		
.0000.	_fma.				
TRI-NEX					
Order Number		to be used with	Scan flag	Abutment screw	Final torque

				30101	
Engaging PA-EL-35	Non-Engaging PA-NL-35	Ø3.5mm TRI-NEX implant	SFT-EL-35	PA-L-18 PA-L-18G	32Ncm
Engaging PA-EL-43	Non-Engaging PA-NL-43	Ø4.3mm TRI-NEX implant	SFT-EL-43		00.40N
Engaging PA-EL-50 -	Non-Engaging PA-NL-50	Ø5mm TRI-NEX implant & TRI-MAX-7	SFT-EL-50	PA-L-20 PA-L-20G	32-40Ncm
Engaging PA-EL-60 -	Non-Engaging PA-NL-60	Ø6mm TRI-NEX implant & TRI-MAX-8 & 9	SFT-EL-60		

DC-RANGE (DEEP CONICAL)

Order Number		to be used with	Scan flag	Abutment screw	Final torque
Engaging PA-DC3	Non-Engaging PA-NDC3	Ø3mm Deep conical implant	SFT-DC3	PA-DC3-14T	15Ncm
Engaging PA-DC4	Non-Engaging PA-NDC4	Ø3.5 & 4mm Deep conical implant	SFT-DC4	PA-DC4-16T	20Ncm
Engaging PA-DC5	Non-Engaging PA-NDC5	Ø5mm Deep conical implant	SFT-DC5	PA-DC5-20T	25-32Ncm

INTERNAL HEX RANGE (M-SERIES & PROVATA®)

Order Number		to be used with	Scan flag	Abutment screw	Final torque
Engaging	Non-Engaging	Ø3.7mm / Ø4.2mm / Ø5mm			
PA-EM-S	PA-NM-S	Internal Hex implants (M-series)			
1.5	a	Ø4mm/Ø5mm Provata	SFT-M		
_		Ø5mm Provata & Ø6mm PROMAX implants			
Engaging PA-EM-SP45 		Ø6mm PROMAX implants		PA-Z-18T PA-Q-18G	32Ncm
Engaging PA-EZ	Non-Engaging PA-NZ	Ø7mm /Ø8mm /Ø9mm PROMAX implants	SFT-Z	-	

IT (INTERNAL OCTAGON) RANGE

Order Number		to be used with	Scan flag	Abutment screw	Final torque
Engaging ITS-PA	Non-Engaging	Ø4.8mm IT implants & MAXIT7	SFT-IT		20 (0)
Engaging ITS6-PA	Non-Engaging ITS6-PA-NE	Ø6mm IT implants & MAXIT8 & 9	SFT-IT6	15-11-PA	32-40NCM

ABUTMENT LEVEL (COMPACT CONICAL ABUTMENT)

Order Number	to be used with	Scan flag	Abutment screw	Final torque
Engaging PA-MC-48	Ø4.8mm Comp Conical	SFT-MC-48	1 Series Titanium screws.	10.45Nem
Engaging PA-MC-60	Ø6mm Comp Conical	SFT-MC-60	for codes & instrument information)	10-15NCM

Description

Passive Abutments are pre-manufactured for aiding prosthetic rehabilitation. These abutments are connected directly to an endosseous implant, or connect the prosthesis to a compact conical abutment.



Detailed Description

The Passive Abutment consists of four components:

- 1. Plastic cylinder this component is not used during CAD/CAM procedures.
- 2. Titanium interfacial component/abutment –this pre-machined component forms the final interface between the milled restoration and the implant.
- 3. Luting screw –this small screw is used to clamp the interfacial component onto the laboratory analogue during the process of luting the casting onto the interfacial component. The component is not used clinically.
- 4. Prosthetic screw –this screw retains the completed prosthesis to the implant at final placement and provides a compressive force across the cement line.

NOTE: Only in the DC range is the prosthetic screw included, for all other systems the prosthetic screw is sold separately.



Materials:

Passive Abutment:Titanium grade,2,3,4 or 5Plastic Sleeve:Polyoxymethylene (POM)Luting screw:Medical grade PEEK (DES-3000)Abutment screws:Titanium alloy Ti-90%, Al-6%, V-4% or Gold Alloy Au-61%, Ag-16.5%, Pt-13.5%, Cu-9%

Overview of Use:

- The passive abutment is selected in the software, and the restoration is designed and milled.
- The milled structure may then undergo further laboratory processing e.g. ceramic firing, finishing and polishing before being assembled with the interfacial component.
- The Passive Abutment (Titanium interfacial component) is kept separate from the manufacturing of the milled structure and is, therefore, not subjected to degradation by heat-cycles or finishing procedures. The integrity of the machined part maintains its original condition.
- The finished milled structure is assembled with the interfacial ring by luting the laboratory model, before placement in the patient's mouth. For assembly, the Titanium interfacial component is clamped to the analogue on the master model by means of the luting screw. The luting screw ensures that the interfacial component is held in full contact with the analogue.
- The finished prosthesis is then luted to the clamped interfacial ring using resin cement. In this way the resin cement serves as a space filler between the milled structure and the interfacial ring. This compensates for any minor fitting and finishing discrepancies, and eliminated any misfit of the milled structure to the implant.
- At placement in the mouth, the prosthetic screw retains the prosthesis to the implant and maintains a compressive force over the cement line. The cement line is, therefore, not responsible for retention of the prosthesis. It is merely a space filler. The luting screw is discarded after the luting procedure.



Digital workflow by using digital libraries

Scanning procedure:

- 1. Download Southern Implants digital library for 3Shape, Exocad, and/or Dental Wings by registering on www.southernimplants.com.
- 2. Load the libraries into the CAD/CAM system by following their instructions.
- 3. The implant position is obtained digitally either by an intraoral scan of the patient with scan flags attached to the endosseous implant/compact conical abutments, or a desktop scan of the dental model with the scan flags attached to the laboratory analogue.
- 4. The scan is then imported into the design software.

Designing procedures:

- 5. The scan flag in the digital form is now matched and aligned with the corresponding scan flag in the library.
- 6. The software recognizes the position of the scan flag to the implant or analogue.
- 7. The Passive abutment/s is chosen from the library. The user will determine the desired tooth or teeth positions needed for the abutment/s, and the path of insertion will be determined. This will govern if a screw retained or cemented restoration is possible.
- 8. The software will guide the user through the steps to complete the restoration.
- 9. Selection of the relevant restorative material will then be required; the most common material to use is Zirconia.

Milling & sintering procedures:

10. The milled custom abutment or screw retained restoration is to be sintered to the specifications of the manufacturer if Zirconia was used. Once sintered the restoration is cemented to the prefabricated Titanium Passive Abutment by following the following steps:

Luting procedures for Passive Abutments:

- 1. It is highly recommended to use a model with removable soft issue mask. This will allow easy access to the analogues for further lab procedures, and will greatly ease latter assembly and processing procedures.
- 2. The luting procedure and concept stays the same when using passive abutments on compact conical level, multiple unit restorations and single unit restorations. For illustration purposes the following steps will show the steps applicable to all.





3. Finger-tighten, with the PEEK luting screw, the titanium interfacial component to the laboratory analogue. Using the corresponding driver as indicated on pages 38-39, and see page 76 for more driver information.

NOTE: Sandblasting of the Titanium interfacial component is not recommended as it rounds of the corners and the exact machined dimensions will be lost.



The milled structure can easily be fitted and removed from the interfacial abutment without the need to remove or replace the PEEK luting screw. This is due to the PEEK luting screw attaching the titanium interfacial abutment to the analogue, and not the milled structure to the analogue.

If it is needed to screw retain the milled structure to the lab analogue to undertake further laboratory procedures, one of the PEEK luting screws can be
removed and a laboratory screw can be used to attach the prosthesis to the master model. This will ensure that the prosthesis does not move during
articulation, opposite set up, or any other laboratory procedures.

- 5. It is important to ultrasonically clean or steam clean the following before the luting procedure starts:
- The titanium interfacial components
- The short luting screws
- The fitting surfaces of the prosthesis
- Also clean the analogues (Implant Replicas) in the model by brushing with soap and water or steam cleaning to remove any debris which may interfere with perfect seating of the interfacial components.
- 6. Clean the restoration in an ultrasonic unit for about 1 minute. Dry it with oil-free air. Apply Monobond Plus (Ivoclar Vivadent) to the cleaned surfaces of the titanium interfacial component and restoration using a brush or microbrush. Allow the Monobond Plus to react for 60 seconds and disperse with compressed air.



7. Lute the Titanium interfacial abutment to the milled structure by applying a self-cure resin cement or dual cure resin cement (e.g. Rely X by 3M) to the surface of all of the titanium rings. (NB: refrigeration of self-cure resin cements will usually lengthen working time for ease of use on multi-unit structures).



8. Limit the amount of resin cement being applied to the angle between the horizontal plane and vertical plane of the titanium ring. This will avoid excess cement extruding upwards through the screw hole in the milled structure and inadvertently locking the luting screw into the cement. Definitely avoid placing any cement in the area immediately around the head of the luting screw.

In the event that cement locks the luting screw in position, a diameter 1.5mm round burr is rotated into the hex of the screw. This usually separates the screw head from the shaft and frees the prosthesis. Take care not to damage the components. The prosthesis can then be removed.

9. Fit the prosthesis over the titanium interfacial component/abutment and settle the prosthesis firmly into place with finger pressure to extrude excess cement. Arch milled structures can be left seated under their own weight to allow cement to harden. Light cure for 60 seconds, un-mount the restoration and remove excess cement using a sharp blade, probe, or hand scaler instrument to make polishing easier. (Extruded cement breaks away easily in large pieces from the outer polished surfaces of the structure and titanium ring).



10. Once resin cement has hardened, remove all luting screws and then remove any prosthetic retaining screws so that the prosthesis can be lifted from the model. Attach polishing protectors to each of the fitting surfaces of the cemented titanium interfacial component. Polish the remaining cement line using a fine-edged lens-shaped rubber wheel. You will notice that the cement line is often not of constant thickness. This variation is indicative of the extent of misfit which existed and has now been corrected by the cement space of the Passive Abutment. (This effect is not as prominent in milled structures as it is in conventional cast structures).

FINISHING PROCEDURES:

Polishing Protector Caps

These items are used to protect the interface and edges of various components, milled structures, and abutments when polishing. They are attached to the relevant Passive Abutments, using the appropriate laboratory screw. Polishing then commences according to standard laboratory procedures. (Refer to CAT-1010 for information on Polishing Protector Caps).



All caps are made from a hardened stainless steel to ensure longer life.

1. Once polishing is completed, remove protector caps and replace the casting on the cleaned model analogues to inspect and verify the quality of fit obtained. (Resin cement is best cleaned from analogues using a brush with alcohol)

The fit would be expected to be excellent in all areas, but, in the unlikely event that a luting error has occurred, the offending titanium ring may be removed, cleaned and re-cemented to the prosthesis as required. A titanium ring can easily be removed by forcing a sharp blade into the cement line, or by punching out the ring using the shaft of a lab hand piece drill applied through the screw access hole (place the bridge rings down on a folded towel for padding and give the drill shaft a sharp tap).

VERY IMPORTANT: as this technique relies absolutely on the accuracy of the master model to achieve passive fit of the prosthesis, it is vital that accurate impression techniques be used and that the quality and condition of the model and analogues be maintained at all times.



Milling & sintering procedures:

- The restoration design is completed and the design file is sent to a milling unit or a production facility.
- The milled component/s is to be sintered to the specifications of the manufacturer if Zirconia is used. Once sintered, the abutment, custom abutment, or restoration is cemented to the prefabricated Titanium abutment by closing the screw channel temporarily to keep the screw channel free of cement.
- . The milled substructure milled is cemented onto the abutments and the final restoration constructed and cemented onto this.

FINISHING PROCEDURES:

Polishing Protector Caps

These items are used to protect the interface and edges of various components, milled structures, and abutments when polishing. They are attached to the relevant abutments, using the appropriate laboratory screw. Polishing then commences according to standard laboratory procedures. (Refer to CAT-1010 for information available for Polishing Protector Caps).



All caps are made from a hardened stainless steel to ensure longer life.

Final Clinical procedures:

- 1. Clean and disinfect the restoration as applicable per the restorative material manufacturer's instructions.
- 2. a. For screw retained unit/s: Place and tighten the restoration screw. Verify the correct seating of the restoration using radiographic imaging. Tighten the restoration, using a torque wrench, to the torque value specified for the applicable prosthetic screw.
 - b. For cement retained unit/s: Place and tighten the restoration screw. Verify the correct seating of the substructure/custom post using radiographic imaging. Tighten the substructure/custom post, using a torque wrench, to the torque value specified for the applicable prosthetic screw. Close the screw access holes and cement the final prosthesis in place.

WORKFLOWS

EXTERNAL HEX

Ø3.0 Prosthetic Platform



Compact Conical Abutment Main Selection SI PASSIVE COMPACT CONICAL **PASSIVE Abutment** (Compact Conical) SFT-MC-48 LAD-MC АРМС Secondary Selection PA-MC-48 NE 1 Series PA-MC-48 Non-Engaging Screws 2 28 1-0 ДП\ **В** 1/3/4/5.5

50

Indirect

Not available

Ø3.25 Prosthetic Platform

Retaining

Screws

Software

Selection











Ø3.25 Prosthetic Platform

	Implants	
CODE	LENGTHS	
IBN	8.5 / 10 / 11.5 / 13 / 15 / 18	
IBNT	8.5 / 10 / 11.5 / 13 / 15 / 18	
IBNT12d-	8.5 / 10 / 11.5 / 13 / 15 / 18	
MSc-IBNT12d-	8.5 / 10 / 11.5 / 13 / 15 / 18	
CODE	LENGTHS	
IBR12d-	8.5 / 10 / 11.5 / 13 / 15 / 18	
MSc-IBR12d-	8.5 / 10 / 11.5 / 13 / 15 / 18	
IBR24d-	8.5 / 10 / 11.5 / 13 / 15 / 18 / 20 / 22 / 24	
MSc-IBR24d-	8.5 / 10 / 11.5 / 13 / 15 / 18 / 20 / 22 / 24	
NVERTA®	,	

CODE	LENGTHS
IV-EX35-45	10 / 11.5 / 13 / 15
IV-EX3512D-45	11.5 / 13 / 15
CODE	LENGTHS
CODE	LENGTHS 10 / 11.5 / 13 / 15



Ø4.0 Prosthetic Platform







Ø4.0 Prosthetic Platform

	Implants	
CODE	LENGTHS	
IBS	7 / 8.5 / 10 / 11.5 / 13 / 15 / 18 / 20	
IBT	6 / 8.5 / 10 / 11.5 / 13 / 15	
IBPS-	8.5 / 10 / 11.5 / 13 / 15	
IBT12d-	8.5 / 10 / 11.5 / 13 / 15 / 18	
MSc-IBT12d-	8.5 / 10 / 11.5 / 13 / 15 / 18	
BAR12d-	8.5 / 10 / 11.5 / 13 / 15 / 18	
MSc-BAR12d-	8.5 / 10 / 11.5 / 13 / 15 / 18	
BAR24d-	8.5 / 10 / 11.5 / 13 / 15 / 18	
MSc-BAR24d-	8.5 / 10 / 11.5 / 13 / 15 / 18	
BAR36d-	8.5 / 10 / 11.5 / 13 / 15 / 18	
MSc-BAR36d-	8.5 / 10 / 11.5 / 13 / 15 / 18	
MAX-6-	6 / 7 / 9 / 11	
MSc-MAX-6-	6 / 7 / 9 / 11	



Ø5.0 Prosthetic Platform





Ø5.0 Prosthetic Platform

	Implants		
CODE	LENGTHS		
BA	6 / 7 / 8.5 / 10 / 11.5 / 13 / 15 / 18		
BAT	6 / 8.5 / 10 / 11.5 / 13 / 15 / 18		
MSc-BAT	6 / 8.5 / 10 / 11.5 / 13 / 15 / 18		ті
BAT12d-	8.5 / 10 / 11.5 / 13 / 15 / 18		т
MSc-BAT12d-	8.5 / 10 / 11.5 / 13 / 15 / 18	-	TII T
MAX-7-	7 / 9 / 11		
MSc-MAX-7-	7/9/11		



Ø6.0 Prosthetic Platform





Ø6.0 Prosthetic Platform

	Implants
CODE	LENGTHS
BBBS	7 / 8.5 / 10 / 11.5 / 13 / 15
BBBT	6 / 8.5 / 10 / 11.5 / 13 / 15 / 18
BBBT12d-	10 / 11.5 / 13 / 15 / 18
BBBT24d-	10 / 11.5 / 13 / 15
MAX-8 [™] -	7 / 9 / 11
MSc-MAX-8 [™] -	7 / 9 / 11



INTERNAL HEX

Standard Prosthetic Platform





Indirect



Standard Prosthetic Platform



CAD/CAM Workflow (Exocad, 3Shape, Dental Wings)

Platform Matched Prosthetics (for Ø5mm implants)

Direct

Implants	Scanning Flag	Digital Analogues	Prosthetic Components	Software Selection	Retaining Screws
			TIB Abutments	Main Selection SI TIB INT-HEX Secondary Selection TIB-M-P45 EN	TS-Z-18
CODE LENGTHS					
IM-T50 08 / 10 / 11 / 13 / 15			CAB Abutments	Main Selection	Gold BS-Z-18
	SFT-M	LAD-M-P45	CAB-M-P45 Engaging	Secondary Selection CAB-M-P45 EN	Brass *(Blackened and for laboratory use only)
					PA-M-18T
			PASSIVE Abutment	Main Selection SI PASSIVE INT-HEX	Titanium PA-Q-18G
			PA-EM-SP45 Engaging	Secondary Selection PA-EM-SP45 EN	Gold PA-M-18B



Standard Prosthetic Platform









Standard Prosthetic Platform

Implants		TIB Abutment		3rd Party Scanbody	Software Selection	Retaining Screws
						TS-Z-18
	LENGTHS	ТІВ-М	Engaging			Titanium
MSC-PRO4	08 / 10 / 11 / 13 / 15 / 18	TIB-M-C1.5 TIB-M-C3	Engaging Engaging			(1.27 Hex)
PRO5	08 / 10 / 11 / 13 / 15 / 18	TIB-NM TIB-NM-C1.5	Non-Engaging Non-Engaging			GS-Q-18
MSC-PRO5	08 / 10 / 11 / 13 / 15 / 18	TIB-NM-C3	Non-Engaging			
PRO12D4	08 / 10 / 11 / 13 / 15 / 18	-				
MSC-PRO12D4	08 / 10 / 11 / 13 / 15 / 18				NB RS 4.3L	(Quad)
PRO12D5	08 / 10 / 11 / 13 / 15 / 18					DC 7 49 *
MSC-PRO12D5	08 / 10 / 11 / 13 / 15 / 18					B3-2-10 "
PROMAX6	07 / 09 / 11					
Platform Matched Prosthetics (for Ø5mm and Ø6mm implants)		TIB-M-P45	Engaging			Brass (1.27 Hex) *(Blackened. For laboratory use only)
CODE	LENGTHS					
PRO5	08 / 10 / 11 / 13 / 15 / 18					
PROMAX6	07 / 09 / 11					
			L			

CAD/CAM Workflow (Exocad, 3Shape, Dental Wings)

Platform Matched Prosthetics (for Ø5mm and Ø6mm implants)

Direct

Implants		Scanning Flag	Digital Analogues	Prosthetic Components	Software Re Selection S	etaining Screws
				TIB Abutments	Main Selection SI TIB INT-HEX Secondary Selection TIB-M-P45 EN	TS-Z-18
CODE	LENGTHS				_	-
PRO5	08 / 10 / 11 / 13 / 15 / 18					Gold
MSC-PRO5	08 / 10 / 11 / 13 / 15 / 18			CAB Abutments	Main Selection SI CAB INT-HEX	BS-Z-18
PROMAX6	07 / 09 / 11	SFT-M	LAD-M-P45		Secondary Selection	
		- -		CAD-IM-P43 Engaging	P	Brass (Blackened and for laboratory use only) 'A-M-18T
				PASSIVE Abutment	Main Selection SI PASSIVE INT-HEX	Titanium
				-P-L PA-EM-SP45 Engaging	Secondary Selection PA-EM-SP45 EN	Gold
	_				_	Brass*

Wide Prosthetic Platform



CODE



Indirect



Wide Prosthetic Platform

Implants				
CODE	LENGTHS			
PROMAX7	07 / 09 / 11			

07 / 09 / 11

07 / 09 / 11

PROMAX8

PROMAX9

TIB-Z Engaging TIB-ZC-1.5 Engaging TIB-XC-3 Non-Engaging TIB-XC-1.5 Non-Engaging TIB-XC-1.5 Non-Engaging TIB-XC-1.5 Non-Engaging Discourse B0 5.0L BS-Z18- Image: State of the	TIB Abuti	ment	3rd Party Scanbody	Software Selection	Retaining Screws
	TIB-Z TIB-Z-C1.5 TIB-NZ TIB-NZ-C1.5 TIB-NZ-C3	Engaging Engaging Engaging Non-Engaging Non-Engaging Non-Engaging	L	— B0 5.0L	TS-Z-18 U Titanium (1.27 Hex) GS-Q-18 U Gold (Quad) BS-Z-18 * U D BS-Z-18 * U D Srass (1.27 Hex) "(Bhorder Srass) "(Bhorder) "Bootory" "(Bhorder) "Bootory" "(Bhorder) "(Dhorder) "(Bhorder) "(Dhorder)

DEEP CONICAL

Ø3.0 Prosthetic Platform






Not available

Ø3.5 Prosthetic Platform







Implants

Ø3.5 Prosthetic Platform

CODE	LENGTHS	
DCT35	08 / 09 / 11 / 13 / 15	
DCC35	08 / 09 / 11 / 13 / 15	
DCT40	06 / 08 / 09 / 11 / 13 / 15	
DCC40	06 / 08 / 09 / 11 / 13 / 15	
DCT40xx-12d	08 / 09 / 10 / 11 / 13 / 15	
DCC40xx-12d	08 / 09 / 10 / 11 / 13 / 15	
DCT50xx-12d	08 / 09 / 10 / 11 / 13 / 15	
DCC50xx-12d	08 / 09 / 10 / 11 / 13 / 15	

INVERTA®

CODE	LENGTHS
IV-DC35-45	10 / 11.5 / 13 / 15
IV-DC3512D-45	11 / 13 / 15
CODE	LENGTHS
IV-DC40-50	10 / 11.5 / 13 / 15
IV-DC4012D-50	11 / 13 / 15

Software Selection Retaining Screws **TIB Abutment** 3rd Party Scanbody TS-DC4-16 ų, 1000 TIB-DC4 TIB-DC4-C1.5 TIB-DC4-C3 Engaging Engaging Engaging Titanium TIB-NDC4 TIB-NDC4-C1.5 TIB-NDC4-C3 Non-Engaging Non-Engaging Non-Engaging AT OS 3.5/4.0L BS-DC4-16*

Ø5.0 Prosthetic Platform





CODE

DCT50

DCC50

Ø5.0 Prosthetic Platform



TRINEX

Ø3.5 Prosthetic Platform

TS-L-18

- Manada

Titanium

TS-L-18C

Coated GS-L-18

Gold BS-L-18*

Brass (Blackened

for laborate use only)

PA-L-18 ίΠ

and a

Titanium PA-L-18G ίΠ)

Concession of the local division of the loca

Gold

PA-L-18B*

Brass



Indirect



Ø3.5 Prosthetic Platform



Ø4.3 Prosthetic Platform

Screws TS-L-20

SHOW SHOW

Titanium

TS-L-20C

Coated

GS-L-20

A DATE

Gold BS-L-20

Brass *(Blackened for laboratory use only)

PA-L-20

あたる

Titanium

PA-L-20G

Gold

PA-L-20B

Brass









Ø4.3 Prosthetic Platform

CODE	LENGTHS
IA-LH-43-	8 / 10 / 11.5 / 13 / 16
IA-LHS-43-	8 / 10 / 11.5 / 13 / 15
IA50-12d	10 / 11.5 / 13 / 16

Implants



Ø5.0 Prosthetic Platform





Indirect



Ø5.0 Prosthetic Platform

	Implants	
CODE	LENGTHS	
IA-LH-50-	8 / 10 / 11.5 / 13 / 16	
IA-LHS-50-	8 / 10 / 11.5 / 13 / 15	
TRI-MAX7-	7 / 9 / 11	TIE



Ø6.0 Prosthetic Platform



Ø6.0 Prosthetic Platform

	Implants	
CODE	LENGTHS	
IA-LH-60-	8 / 10 / 11.5 / 13 / 16	
TRI-MAX8	7 / 9 / 11	
TRI-MAX9	7 / 9 / 11	



IT RANGE

Ø4.8 Prosthetic Platform



Ø4.8 Prosthetic Platform

	Implants	TIB Abutment
	6	
CODE	LENGTHS	
ITC3	06 / 08 / 10 / 12 / 14	
ITT4	08 / 10 / 12 / 14	
ITC4	06 / 08 / 10 / 12 / 14	TIB-ITS TIB-ITS-C1.5 TIB-ITS-C3
ITC5	06 / 08 / 10 / 12 / 14	TIB-ITSNE TIB-ITSNE-C1.5
ITT5	08 / 10 / 12 / 14	TIB-ITSNE-C3
ITST12d-4xxf	08 / 10 / 12 / 14	L
ITST12d-5xxf	08 / 10 / 12 / 14	-
MAXIT7-	7 / 9 / 11	



Ø6.5 Prosthetic Platform



Ø6.5 Prosthetic Platform

	Implants	
CODE	LENGTHS	
ITC6-5	06 / 08 / 10 / 12 / 14	
ITT6-5	08 / 10 / 12 / 14	
ITT6	08 / 10 / 12 / 14	
MAXIT [®] 8-	7 / 9 / 11	
MAXIT [®] 9-	7 / 9 / 11	-







ABUTMENT DRIVERS & SCREW INFORMATION



INSTRUMENT INFORMATION

Ratchet & Torque Attachment Wrench



I-FME-XS / M / L



Fixture Mount extension (hex top)



I-WI-CS / L

Converts a hand-piece/latch driven instrument with W&H hex to fit a Southern torque/ratchet wrench



For Handpiece (latch-type) inserts featuring the W&H hex





For SQUARE connection of fixture mounts and instruments

I-WI-SH



For HEX connection of fixture mounts

I-WI-SL



For Handpiece inserts (Latch-type) without the W&H hex.

IEX tion of

These items are used to protect the interface and edges of various components, castings, and abutments when polishing. They are attached to the relevant cast components, using the appropriate laboratory screw. Polishing then commences according to standard laboratory procedures.

- Short (S) Polishing Protector Caps are indicated when polishing multiple unit restorations.
- Long Polishing Protector Caps are indicated when polishing single unit restorations. All caps are made from a hardened stainless steel to ensure longer life.

	CODE	S	Ø implant
PPC-BN	PPC-IBN		Ø3.25mm
PPC-NPI	PPC-INPi		Ø3.25mm
PC-IB	PPC-IB		Ø4.0mm
PPC-BA	PPC-BA		Ø5.0mm
PPC-B3B	PPC-BBB		Ø6.0mm
PPC-MAX9	PPC-MAX9		Ø9.0mm

DC IMPLANT RANGE

	CODE	S	Ø implant
P - (4×3)	PPC-DC3		Ø3.0mm
P1/2-P214	PPC-DC4	0.045	Ø4.0mm
Provinces	PPC-DC5	033	Ø5.0mm

INTERNAL HEX (M-SERIES / PROVATA) IMPLANT RANGE

	CODE	S	Ø implant
FPC rd	РРС-М	20	Standard restorative interface
P 52	PPC-Z		Wide restorative interface

TRI-NEX IMPLANT RANGE

EXTERNAL HEX IMPLANT RANGE

	CODE	S	Ø implant
PPC-35	PPC-35	353	Ø3.5mm
PPC-43	PPC-43	438	Ø4.3mm
PFC 50	PPC-50	.5CS	Ø5.0mm
PPC-50	PPC-60	605	Ø6.0mm

IT IMPLANT RANGE

	CODE	S	Ø implant
PPC-ITS	PPC-ITS		Ø4.8mm restorative interface
PPC-IIS6	PPC-ITS6		Ø6.5mm restorative interface

STANDARD ABUTMENTS

PPC-AB-1	- replicates standard abutments with Ø4.50mm restorative interface
PPC-AWB-1	- replicates wide standard abutments with Ø5.50mm restorative interface

CONICAL AND COMPACT CONICAL ABUTMENTS

	PPC-MC-1	- replicates compact conical abutments
	PPC-MCW-1	- replicates wide compact conical abutments

CONTRAINDICATIONS

Do no use in:

- patients who are allergic to or have hypersensitivity to Titanium, Titanium Alloy or Titanium
- · patients who are medically unfit for implant procedures (e.g. uncontrolled diabetes and untreated infection in nearby bone)
- · cases where an adequate number of implants cannot be placed to achieve full functional support for a prosthesis
- · cases where there is insufficient blood supply to the implant site
- patients with insufficient mental health precluding patient co-operation
- patients who abuse drugs or alcohol
- · cases where a pre-operative screening exposes possible risks to the healing of the bone or soft-tissue
- patients who by nature of their condition, occupation or activity will be unable to keep the implant site clean
- patients where site specific contraindications exist.

WARNINGS

All Southern Implants products are intended to be used by appropriately trained and licensed professionals. For the safe and effective use of Osseointegrated Fixtures, it is strongly suggested that specialised training be undertaken, including hands-on training to learn proper technique and radiographic evaluations. THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING. Responsibility for proper patient selection, adequate training, experience in the placement of implants, and providing appropriate information for informed consent, rests with the practitioner. Improper technique can result in damage to anatomical structures during placement, implant failure and/or loss of supporting bone. Southern Implants will not accept liability for damage caused by improper implant treatment.

Southern Implants osseointegrated fixtures have only been validated for use with the corresponding Southern Implants abutments and accessories. Although care has been taken to create interfaces that are equivalent to similar products on the market, Southern Implants cannot guarantee outcomes obtained using components from other manufacturers. Southern Implants will not accept liability for damage caused by improper selection of incompatible abutments and accessories.

Failure to recognize actual drill lengths relative to radiographic measurements can result in permanent injury to nerves or other vital structures.

Do not re-use Implants, Cover screws, Temporary Abutments and Abutments. These are single-use products. DO NOT re-sterilize or autoclave these components. Re-using these components may result in damage on the surface or critical dimensions. This may result in performance and compatibility issues. The removal of proteins from metal (such as Titanium) is extremely difficult and re-use can lead to secondary infections. Sterility is assured unless the container or seal is damaged or opened.

Electro-surgery should not be attempted around metal implants, as they are conductive.

One hundred percent implant success cannot be guaranteed. Non-observance of the indicated limitations of use and working steps may result in failure. Minimizing trauma to the host tissue increases the potential for successful osseointegration.

CAUTIONS

New and experienced Implant users should do training before using a new system or attempt to do a new treatment method.

Take special care when treating patients who have local or systemic factors that could affect the healing of the bone and soft tissue. (I.e. poor oral hygiene, uncontrolled diabetes, are on steroid therapy, smokers, infection in the nearby bone and patients who had oro-facial radiotherapy.)

Thorough screening of prospective implant candidates must be performed including:

A comprehensive medical and dental history.

Visual and radiological inspection to determine adequate bone dimensions, anatomical landmarks, occlusal conditions, periodontal status, and adequacy of bone. Bruxism and unfavourable jaw relations must be taken into account.

Proper pre-operative planning with a good team approach between well trained surgeons, restorative dentists and lab technicians is essential for successful implant treatment.

Minimizing the trauma to the host tissue increases the potential for successful osseointegration.

Electro-surgery should not be attempted around metal implants, as they are conductive.

Re-use of single-use items may result in damage on the surface or to critical dimensions.

During procedure

Care must be taken that parts are not swallowed during any of the procedures, thus rubber-dam application is recommended when appropriate. Care must be taken to apply the correct tightening torque of abutments and abutment screws.

Post-procedure

Regular patient follow-up, and proper oral hygiene must be achieved are essential for favourable long-term results.

PACKAGING

- 1) Components are packed in a peel pouch. Labeling information is located on the bottom half of the pouch inside the packet.
- 2) Non-sterile components used in the laboratory are supplied clean but not sterile. These are: laboratory analogs, passive abutments, cast waxing sleeves, casting precision tools and gold abutments with plastic sleeves.

CLEANING

- Refer to CAT-1039
- Final restoration should be cleaned and disinfected, as per restorative material manufacturer's instructions, before intra oral use.
- Used surgical trays including grommets must be cleaned with suitable disinfectants.
- Multiple-part instruments must be disassembled prior to cleaning and sterilization.
- Internal debris/residue on instruments must be removed with a soft brush.
- Instruments should be inspected, cleaned separately and discarded if damaged.
- Best results are achieved if instruments are cleaned by material type.
- Instruments and trays can be cleaned and disinfected in a dedicated instrument washer or alternatively by hand, followed by an ultrasonic bath with a detergent
 appropriate for surgical instruments.
- · Instruments and trays must be rinsed and dried thoroughly.

STERILIZATION

- · Refer to CAT-1039
- Pre-vacuum sterilization method: Steam sterilize the component at 132°C (270°F) at 180-220 kPa for 4 minutes, or at 135°C (275°F) at 180-220 kPa for 3 minutes. Dry for at least 20 minutes in the chamber. Only an appropriate regulatory authority approved sterilizer and wrap/pouch for steam sterilization must be used.
- · It is the responsibility of the user to establish whether or not their sterilizer is appropriate regulatory authority approved to meet the recommended parameters.
- The product must be stored in a dry place in the original packaging at room temperature and not exposed to direct sunlight. Incorrect storage may influence device characteristics.

STORAGE

- · Maintain original packaging integrity in storage.
- Devices should be stored at room temperature and must not be exposed to direct sunlight.
- Incorrect storage may influence device characteristics.
- Packaging, after re-sterilization, should be dry before storage to avoid corrosion and degradation of cutting edges.

DISPOSAL

Disposal of the device and its packaging shall follow local regulations and environmental requirements, taking different contamination levels into account.

MAGNETIC RESONANCE (MR) SAFETY INFORMATION

This device has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artefact in the MR environment. The safety of this device in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

CAUTION: (USA ONLY)

United States Federal Law restricts this device to sale to, or on the order of, a licensed dentist or physician.

DISCLAIMER OF LIABILITY

This product is part of the Southern Implants' product range and should only be used with the associated original products, according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants' product range and take full responsibility for the correct indications and use of this product. Southern Implants does not assume liability for damage due to incorrect use. Please note that some Southern Implants' products may not be cleared or released for sale in all markets.

For more information, please visit your Southern Implants Representative or visit southernimplants.com

The following symbols are used on our packaging labels and they indicate the following:



2

For more information, please contact your Southern Implants Representative or visit southernimplants.com



Subsidiaries

Australia Southern Implants Australia T: +61-2-8076-9337 E: info@southernimplants.com.au United Kingdom Southern Implants UK T: +44-20-8899-6845 / 6 / 7 E: info@southernimplants.co.uk

South Africa - Headquarters 1 Albert Road, Irene, RSA T: +27-12-667-1046 | E: info@southernimplants.com

> USA and Canada Southern Implants North America Inc. T: +1-561-472-0990 E: customercare@southernimplants.com

Southern Implants Products are distributed world-wide, please visit southernimplants.com for a list of Distributors.

©2019 Southern Implants (Pty) Ltd. All rights reserved. Co-Axis, Southern Implants, and Subcrestal Angle Correctionare registered trademarks and Body-Shift, INVERTA and PROVATA are trademarks of Southern Implants (Pty) Ltd. Not all implants are available in all markets.