



SOUTHERNIMPLANTS®

Innovative Treatment Solutions


SIREAL

Guided Surgery System



Product Catalogue





Southern Implants® is a leading provider of unique and innovative dental implant products with a focus on top-end professional users who want more choices. Southern's expertise in research, development and manufacturing of dental implants allows us to provide Innovative Treatment Solutions that will reduce treatment times and improve patient outcomes.

Striving for excellence and meeting customer needs, has led to our wide product range characterised by Unique and Innovative products which include:

- Multiple interfaces, to suit customer preference.
- INVERTA® implant, featuring a Body-Shift™ design, engineered for primary stability and suitable for immediate loading.
- Co-Axis®, Subcrestal Angle Correction® implants, available in angulations of 12°, 24° and 36° and various internal and external connections.
- MAX implant, specifically designed for immediate molar tooth replacement.
- The ZYGAN®, ZYGEX and ZYGIN implants for severely resorbed maxilla and craniofacial reconstruction.
- The Machined Surface Coronally (MSC) dental implant surface treatment offers practitioners an innovative way to take advantage of the best characteristics of both smooth and moderately rough implant surfaces.

Our product portfolio is in synchronised evolution with protocol improvements and technological advances.

My sincere thanks to all specialists, dentists and technicians who put their trust in our company.


Graham Blackbeard
Managing Director, Southern Implants

CONTENTS

Introduction	Page 04
SIREAL GUIDED SURGERY	
SIREAL Standard	Page 05
Universal Guided Surgery Tool.....	Page 06
Offset Sleeves.....	Page 06
Assembly & use of the universal guided surgery tool.....	Page 07
Step 1: Setting the offset.....	Page 07
Step 2: Inserting an instrument.....	Page 07
Choosing the correct offset sleeve.....	Page 08
SIREAL Standard guide sleeves.....	Page 09
Placement Tools.....	Page 10
Cortical Perforator.....	Page 10
Surgical Tray.....	Page 11
Clinical Procedure.....	Page 12
Site Preparation Protocols.....	Page 14
Guided Workflows	
External Hex.....	Page 16
PROVATA®.....	Page 20
Internal Hex (M-series).....	Page 22
TRI-NEX®.....	Page 24
Deep Conical.....	Page 26
INVERTA® External Hex.....	Page 28
INVERTA® Deep Conical.....	Page 30
SIREAL Narrow	Page 32
Narrow Diameter Guided Surgery.....	Page 33
SIREAL Order Guide.....	Page 34
Explanation of Labeling Symbols.....	Page 35

For more information scan the below



or visit

SOUTHERNIMPLANTS.COM

INTRODUCTION

The **Southern Implants** guided surgery solution, **SIREAL** guide, provides a complete computer-assisted dental implant planning and placement solution for Southern Implants' tapered implants. This is achieved by virtual prosthesis and on-screen design of a surgical guide, enabling prosthetically driven implant placement.

Surgical guide types

The surgical guide type selection depends on the dental professional's preference, patient anatomy and the available planning software.

There are three types of surgical guides:



NOTE: all surgical guides are patient specific and consists of a 3D printed or milled acrylic guide and metal guide sleeves.

Treatment planning

Diagnostic and patient specific conditions influence the guided treatment plan. The type of restoration, provisionals, number of implants and imaging procedures must be taken into consideration during planning.

The following considerations should be reviewed during pre-planning:

- quantity, quality and health of both soft and hard tissues.
- occlusal analysis.
- oral hygiene assessment.
- the patient's vertical opening of the mouth needs to be sufficient to accommodate the instruments used during guided surgery.

CT scanning

Several imaging technologies are available to accurately scan data. The dental professional and/or radiologist, needs to follow the instructions of the imaging system used.

Warning: there may be distortion in the CT image data. These distortions could lead to fit and trajectory problems. It is recommended to validate the guide fit and trajectory by taking a CT scan of the patient wearing the guide before surgery. Open the CT scan image to review both the position and orientation of the guide sleeve. Measure guide sleeve distance and orientation in the CT scan and compare it to the offset/prolongation selected during the planning phase.

The dental professional must follow Southern Implants sleeve offsets and prolongations, failing to do so will result in patient injury. The guide manufacturer ensures compatibility with Southern Implants guided instruments by using SIREAL Guide sleeves, and instruments positioned according to offsets and prolongations described in this manual.

Verify the fit of the guide by seating it on the patient's jaw. It is recommended to validate the fit and sleeve position with a CT scan of the patient with the guide in-situ. If the guide was manufactured on a stone model, the inaccuracy of the model or poor image quality from the scan data may result in the guide not fitting. Should there be a variance, do not proceed, remake the guide.

After fixing the guide into place, proceed using SIREAL drills and instrumentation to prepare the osteotomy. The surgical protocol together with the surgical guide will govern which instruments are required to prepare each implant site.

Please note:

- images are for illustration purposes only and do not necessarily accurately represent the product.
- all dimensions in this catalogue are in mm, unless otherwise specified.
- not all products are cleared for sale in all countries.

SIREAL STANDARD



Implant planning



Sleeves

Ø6.2 (outer) / Ø5.2 (inner)



Ø5.1 offset sleeve

9 mm; 10.5 mm; 12 mm; 14 mm




Commission Guide



SIREAL GUIDED SURGERY: THE CONCEPT

UNIVERSAL GUIDED SURGERY TOOL

The universal guided surgery tool, I-DE-GS4, from Southern Implants is the solution to the SIREAL guided surgery system. This tool allows clinicians to utilise their standard Southern Implants drill kit, and convert the drills into guided surgery drills.



The universal guided surgery tool features a “twist lock” mechanism which locks the latch grip of the drills and placement tools into the I-DE-GS4 tool.

The W&H hex allows for handpieces with the W&H connection to engage the tool, allowing torque to be applied through the instrument.

Note: High torque can only be applied to instruments with a W&H hex to a maximum of their specific torque rating, and no higher than 70 Ncm. Instruments and drills without the W&H hex (universal tools) do not exceed 40 Ncm.

OFFSET SLEEVES



Four offset sleeves are available: **9 mm**, **10.5 mm**, **12 mm** and **14 mm**. This is to accommodate the patient’s vertical opening or adjacent teeth height that could interfere with the guide sleeve.

Offset is measured from the implant platform to the top of the guide sleeve.

SIREAL GUIDED SURGERY: THE CONCEPT

Assembly and use of the universal guided surgery tool:

Follow instructions illustrated below, with the handpiece latch pointing to the right.

Step 1: Setting the offset sleeve

The offset is the distance between the implant platform to the top of the surgical guide sleeve.


Select the offset sleeve that corresponds to one of the offset lengths (9 mm, 10.5 mm, 12 mm and 14 mm) during the implant planning and guide design phase.

Screw the offset sleeve in a **clockwise** direction onto the universal guided surgery tool, with the latch facing to the right as illustrated below.



Step 2: Inserting a drill or placement tool

The universal guided surgery tool is designed primarily for use with drills and placement tools, to allow for partially or fully guided surgeries.



Ensure that the twist lock on the I-DE-GS4 is rotated fully in the **clockwise** direction before inserting the drill or placement tool. Insert the drill/placement tool until it seats inside the tool (this might require rotation until the seat lines up with the latch).

Once the latch is seated, rotate the twist lock in an **anti-clockwise** direction until locked into position.

To release the drill or placement tool, rotate the twist lock 45° in a **clockwise** direction and separate the drill/placement tool from the I-DE-GS4.

SIREAL GUIDED SURGERY: THE CONCEPT

Choosing the correct offset and offset sleeve


The example cases below utilises an 11 mm length implant.

In a molar site, the vertical opening of the patient will be limited. When placing the implant in a posterior site, it's best to go for a shorter total drill length. Select the 9 mm offset sleeve (the top of the guide sleeve is 9 mm from the implant platform), and that is where the drill will stop.

An 11 mm implant + 9 mm offset = 20 mm total drill length (which is the total length from the top of the guide sleeve to the apex of the osteotomy).

In an anterior case where the patient has long dentition and you can't fit the sleeve in between the adjacent teeth, lift the sleeve offset to 14 mm above the planned implant platform.


The 11 mm implant + 14 mm offset = 25 mm total drill length.



Illustrating how to determine the correct offset sleeve to determine the maximum drill length

This planning and sleeve selection can be utilized for both the 10.5 mm and 12 mm.

SIREAL Offset is measured from implant platform to top of the sleeve.



Scan or visit southernimplants.com for online calculator.


SIREAL GUIDED SURGERY: THE CONCEPT

SIREAL STANDARD GUIDE SLEEVES

Closed sleeves

Indications for use:

- should be used with tapered implants.
- Available in 4 mm and 9 mm lengths.



Fixation Pins and Fixation Pin Sleeves

Indications for use:

- the Southern Implants® fixation pins are used to stabilise the surgical guide.
- a D-12T-M15 drill is used to drill through the fixation pin sleeve while guide is in situ. After drilling, insert the pin.
- it is recommended to use 3 pins for full arch guide. If a tooth supported guide requires additional stability, a minimum of 2 pins should be used.
- pin(s) (I-D12-GP) must not interfere with placement tool (I-DE-GS4) or drill trajectory.
- vertical opening and anatomical constraints of the patient must be considered when designing the guide with fixation pins.


C-Sleeves (open sleeves)

Indications for use:

- should be used with Co-Axis® tapered implants.
- should be positioned to allow access to the fixture mount screw. This will assist the user to remove the fixture mount when placing the implant fully guided.
- used in the posterior region where vertical opening is a challenge. The universal guided surgery tool and drill can be implemented from the side which allows additional space saving of the offset distance.
- allows irrigation at the osteotomy site while drilling.

NOTE:

- the lip on the guide sleeve, adds 0.25 mm, this does not need to be taken into consideration as most Southern Implants drills extend 1 mm longer.
- **always plan for at least 2 mm from nerves /anatomical structures.**



CLINICAL PROCEDURE

Step 1: insert the cortical perforator (D-GS-CP) through the guide sleeve to perforate the cortical bone.



WARNING: this tool fits directly into the handpiece and should not be used with the I-DE-GS4 tool.

Step 2: insert the D-3SPADE-1.8M into the I-DE-GS4 and proceed with pilot drilling.



Step 3 (optional): place the D-20T-M10 into the I-DE-GS4.



Step 4: place the D-DCT3011 into the I-DE-GS4.



CLINICAL PROCEDURE

Step 5: place the D-DCT3511 into the I-DE-GS4.





Step 6: place the D-DCT4011 into the I-DE-GS4.



Step 7: place the placement tool I-HDC4-GS into the I-DE-GS4.




Step 8: good primary stability will govern if immediate loading can be done or not.




NOTE: for implant torque, refer to surgical manual.

Ø3.0 mm Tapered (DCT30)


DCT3009



DCT3011




DCT3013




Ø3.5 mm Tapered (DCT35)


DCT3508




DCT3509




DCT3511



DCT3513




DCT3515




Ø4.0 mm Tapered (DCT40)


DCT4006




DCT4008




DCT4009




DCT4011



DCT4013




DCT4015




Ø5.0 mm Tapered (DCT50)


DCT5009




DCT5011



DCT5013



DCT5015



INSTRUMENTATION

PLACEMENT TOOLS

This is for placement of Southern Implants tapered implants, partially and fully guided.


Partially guided

The surgical guide is designed to prepare the osteotomy from pilot to final drill. Implant placement is not done through the surgical guide.

Fully guided

The surgical guide is designed to prepare the osteotomy from pilot to final drill, as well as placing the implant through the surgical guide.

TAPERED IMPLANTS SUPPORTED FOR FULLY GUIDED SURGERY



Insertion tools supplied with PEEK bit.

Important: the PEEK bits (I-PBIT-H16/ L18 /L20) should be replaced on a regular basis. Items sold separately. General wear and tear are to be expected with regular use.


TAPERED IMPLANTS SUPPORTED FOR PARTIAL GUIDED SURGERY




NOTE: these placement tools are specifically developed for the SIREAL Guided System. The standard tools should NOT be used for fully guided surgery.

Cortical Perforator

The cortical perforator, D-GS-CP or D-GS-CP4, is used to initiate the osteotomy by perforating the cortical plate at the planned implant position.



(for use with the standard 05.1 guide sleeves)



(for use with the narrow 04.1 guide sleeves)


NOTE: this tool fits directly into the handpiece and **should not** be used with the I-DE-GS4 tool.

SURGICAL TRAY

I-DC-EG Deep Conical Instrument Tray

(for demonstration purposes only)


Pilot Drills
 D-3Spade-1.8M
 D-3SPADE-CP
 D-16T



Ø2.0 mm Twist Drill
 D-20T-M10
 D-20T-M15
 D-DC20
 D-220C


Important:

- **D-20TM10:** 10 mm twist drills are indicated for implant lengths of 11.5 mm and shorter.
- **D-20TM15 or D-DC20:** 15 mm twist drills are indicated for 12 mm implants and longer.



Fixation Pin Drill
 D-12T-M15

Minimum of 2 required per surgery

Primary Drills (Ø3.0, Ø3.5, Ø4.0, Ø5.0)

Additional Drills (6.1, 8, 9, 11, 13, 15)

Placement Tools

Drivers

Cortical Perforators

D-GS-CP
 D-GS-CP4



Placement Tools

I-HDC3-GS
 I-HDC4-GS
 I-HDC5-GS
 Handpiece Insert



Tray Pins

I-GS-TP
 I-GS-TPG

Pins to insert into surgical trays to hold offset sleeves in place (I-GSS04/05-set).




Dedicated Tapered Drills

Ø3.0 D-DCT3009 D-DCT3011 D-DCT3013	Ø3.5 D-DCT3508 D-DCT3509 D-DCT3511 D-DCT3513 D-DCT3515	Ø4.0 D-DCT4006 D-DCT4008 D-DCT4009 D-DCT4011 D-DCT4013 D-DCT4015	Ø5.0 D-DCT5009 D-DCT5011 D-DCT5013 D-DCT5015
--	--	---	---




TISSUE CUTTERS (Optional)

SIREAL standard




Narrow



CAUTION: when drilling close to crucial anatomical landmarks, consider that the drill preparation site may be up to 1 mm deeper than the corresponding implant length.


TAPERED DRILLS:

- SIREAL universal guided surgery tool is only to be used with Southern Implants standard length tapered drills.
- do not use long shaft drills. It will drill deeper than the planned depth.
- long drills can be identified by an "L" in the product code. For example: D-40TP-11.5-L



INVERTA® DRILLS:

When placing INVERTA® implants with SIREAL, the guided INVERTA drills laser marked with "GS" must be utilised.



Implant Code	Initiate the osteotomy		Drill sequence per bone densities	
	Cortical Perforator	Initial Drill	Optional Drill for medium and dense cortical bone	Final Drill
IP8.5	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-30TP-8.5
IP10	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-30TP-10
IP11.5	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-30TP-11.5
IP13	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-30TP-13
IP15	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-30TP-15

Implant Code	Initiate the osteotomy		Drill sequence per bone densities		
	Cortical Perforator	Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill Final Drill for soft bone	Final Drill for medium and dense bone
IBNT8.5	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-30TP-8.5	D-33TP-8.5
IBNT10	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-30TP-10	D-33TP-10
IBNT11.5	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-30TP-11.5	D-33TP-11.5
IBNT13	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-30TP-13	D-33TP-13
IBNT15	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-30TP-15	D-33TP-15
IBNT18	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-30TP-18	D-33TP-18

Implant Code	Initiate the osteotomy		Drill sequence per bone densities			
	Cortical Perforator	Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill Final Drill for soft bone	Final Drill for medium and dense bone	Final Drill for medium and dense bone
IBT6	D-GS-CP4	D-3SPADE-GS	D-20FM10	D-30TP-8.5	D-40TP-6	D-40TP-6
IBT8.5	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-33TP-8.5	D-40TP-8.5	D-40TP-8.5
IBT10	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-33TP-10	D-40TP-10	D-40TP-10
IBT11.5	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-33TP-11.5	D-40TP-11.5	D-40TP-11.5
IBT13	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-33TP-13	D-40TP-13	D-40TP-13
IBT15	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-33TP-15	D-40TP-15	D-40TP-15
IBT18	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-33TP-18	D-40TP-18	D-40TP-18

Implant Code	Initiate the osteotomy		Drill sequence per bone densities			
	Cortical Perforator	Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill Final Drill for soft bone	3rd Drill Final Drill for soft bone	Final Drill for medium and dense bone
BAT6	D-GS-CP	D-3SPADE-GS	D-20FM10	D-40TP-6	D-50TP-6	D-50TP-6
BAT8.5	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-33TP-8.5	D-40TP-8.5	D-50TP-8.5
BAT10	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-33TP-10	D-40TP-10	D-50TP-10
BAT11.5	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-33TP-11.5	D-40TP-11.5	D-50TP-11.5
BAT13	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-33TP-13	D-40TP-13	D-50TP-13
BAT15	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-33TP-15	D-40TP-15	D-50TP-15
BAT18	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-33TP-18	D-40TP-18	D-50TP-18

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities				Final Drill NOT through the guide Final Drill for medium and dense bone
		Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill	3rd Drill	4th Drill		
BBBT6	D-GS-CP	D-3SPADE-GS	D-20FM10	D-33TP-8.5	D-40TP-6	D-50TP-6	D-60TP-6	
BBBT8.5	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-33TP-8.5	D-40TP-8.5	D-50TP-8.5	D-60TP-8.5	
BBBT10	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-33TP-10	D-40TP-10	D-50TP-10	D-60TP-10	
BBBT11.5	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-33TP-11.5	D-40TP-11.5	D-50TP-11.5	D-60TP-11.5	
BBBT13	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-33TP-13	D-40TP-13	D-50TP-13	D-60TP-13	
BBBT15	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-33TP-15	D-40TP-15	D-50TP-15	D-60TP-15	
BBBT18	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-33TP-18	D-40TP-18	D-50TP-18	D-60TP-18	

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities			
		Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill	Final Drill for soft bone	Final Drill for medium and dense bone	
IBNT12D-8.5	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-30TP-8.5	D-33TP-8.5	D-33TP-8.5	
IBNT12D-10	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-30TP-10	D-33TP-10	D-33TP-10	
IBNT12D-11.5	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-30TP-11.5	D-33TP-11.5	D-33TP-11.5	
IBNT12D-13	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-30TP-13	D-33TP-13	D-33TP-13	
IBNT12D-15	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-30TP-15	D-33TP-15	D-33TP-15	
IBNT12D-18	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-30TP-18	D-33TP-18	D-33TP-18	

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities			
		Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill	Final Drill for soft bone	Final Drill for medium and dense bone	
IBT12D-8.5	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-33TP-8.5	D-40TP-8.5	D-40TP-8.5	
IBT12D-10	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-33TP-10	D-40TP-10	D-40TP-10	
IBT12D-11.5	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-33TP-11.5	D-40TP-11.5	D-40TP-11.5	
IBT12D-13	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-33TP-13	D-40TP-13	D-40TP-13	
IBT12D-15	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-33TP-15	D-40TP-15	D-40TP-15	
IBT12D-18	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-33TP-18	D-40TP-18	D-40TP-18	

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities			
		Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill	Final Drill for soft bone	Final Drill for medium and dense bone	
IBR12D-8.5	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-33TP-8.5	D-40TP-8.5	D-40TP-8.5	
IBR12D-10	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-33TP-10	D-40TP-10	D-40TP-10	
IBR12D-11.5	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-33TP-11.5	D-40TP-11.5	D-40TP-11.5	
IBR12D-13	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-33TP-13	D-40TP-13	D-40TP-13	
IBR12D-15	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-33TP-15	D-40TP-15	D-40TP-15	
IBR12D-18	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-33TP-18	D-40TP-18	D-40TP-18	

Scan or visit
SOUTHERNIMPLANTS.COM
for online calculator.



Implant Code	Cortical Perforator	Drill sequence per bone densities			
		Initiate the osteotomy Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill Final Drill for soft bone	Final Drill for medium and dense bone
IBR24D-8.5	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-33TP-8.5	D-40TP-8.5
IBR24D-10	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-33TP-10	D-40TP-10
IBR24D-11.5	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-33TP-11.5	D-40TP-11.5
IBR24D-13	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-33TP-13	D-40TP-13
IBR24D-15	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-33TP-15	D-40TP-15
IBR24D-18	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-33TP-18	D-40TP-18

Implant Code	Cortical Perforator	Drill sequence per bone densities			
		Initiate the osteotomy Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill	3rd Drill Final Drill for soft bone
BAT12D-10	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-33TP-10	D-40TP-10
BAT12D-11.5	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-33TP-11.5	D-40TP-11.5
BAT12D-13	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-33TP-13	D-40TP-13
BAT12D-15	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-33TP-15	D-40TP-15
BAT12D-18	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-33TP-18	D-40TP-18

Implant Code	Cortical Perforator	Drill sequence per bone densities			
		Initiate the osteotomy Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill	3rd Drill Final Drill for soft bone
BAR12D-8.5	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-33TP-8.5	D-40TP-8.5
BAR12D-10	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-33TP-10	D-40TP-10
BAR12D-11.5	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-33TP-11.5	D-40TP-11.5
BAR12D-13	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-33TP-13	D-40TP-13
BAR12D-15	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-33TP-15	D-40TP-15
BAR12D-18	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-33TP-18	D-40TP-18

Implant Code	Cortical Perforator	Initiate the osteotomy				Drill sequence per bone densities		
		Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill	3rd Drill	Final Drill for soft bone	Final Drill for medium and dense bone	
BAR24D-8.5	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-33TP-8.5	D-40TP-8.5	D-50TP-8.5		
BAR24D-10	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-33TP-10	D-40TP-10	D-50TP-10		
BAR24D-11.5	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-33TP-11.5	D-40TP-11.5	D-50TP-11.5		
BAR24D-13	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-33TP-13	D-40TP-13	D-50TP-13		
BAR24D-15	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-33TP-15	D-40TP-15	D-50TP-15		
BAR24D-18	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-33TP-18	D-40TP-18	D-50TP-18		

Implant Code	Cortical Perforator	Initiate the osteotomy				Drill sequence per bone densities		
		Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill	3rd Drill	Final Drill for soft bone	Final Drill for medium and dense bone	
BAR36D-8.5	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-33TP-8.5	D-40TP-8.5	D-50TP-8.5		
BAR36D-10	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-33TP-10	D-40TP-10	D-50TP-10		
BAR36D-11.5	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-33TP-11.5	D-40TP-11.5	D-50TP-11.5		
BAR36D-13	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-33TP-13	D-40TP-13	D-50TP-13		
BAR36D-15	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-33TP-15	D-40TP-15	D-50TP-15		
BAR36D-18	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-33TP-18	D-40TP-18	D-50TP-18		

Implant Code	Cortical Perforator	Initiate the osteotomy				Drill sequence per bone densities				Final Drill NOT through the guide Final Drill for medium and dense bone
		Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill	3rd Drill	4th Drill	Final Drill for soft bone			
BBBT1 2D-10	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-33TP-10	D-40TP-10	D-50TP-10	D-60TP-10			
BBBT1 2D-11.5	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-33TP-11.5	D-40TP-11.5	D-50TP-11.5	D-60TP-11.5			
BBBT1 2D-13	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-33TP-13	D-40TP-13	D-50TP-13	D-60TP-13			
BBBT1 2D-15	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-33TP-15	D-40TP-15	D-50TP-15	D-60TP-15			
BBBT1 2D-18	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-33TP-18	D-40TP-18	D-50TP-18	D-60TP-18			

Implant Code	Cortical Perforator	Initiate the osteotomy				Drill sequence per bone densities				Final Drill NOT through the guide Final Drill for medium and dense bone
		Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill	3rd Drill	4th Drill	Final Drill for soft bone			
BBBT24D-10	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-33TP-10	D-40TP-10	D-50TP-10	D-60TP-10			
BBBT24D-11.5	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-33TP-11.5	D-40TP-11.5	D-50TP-11.5	D-60TP-11.5			
BBBT24D-13	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-33TP-13	D-40TP-13	D-50TP-13	D-60TP-13			
BBBT24D-15	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-33TP-15	D-40TP-15	D-50TP-15	D-60TP-15			
BBBT24D-18	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-33TP-18	D-40TP-18	D-50TP-18	D-60TP-18			

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities			Fully guided
		Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill Final Drill for soft bone	Final Drill for medium and dense bone	Final Drill for medium and dense bone	
PRO308	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-30TP-8.5	D-30TP-8.5	I-H-PRO3-GS	
PRO310	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-30TP-10	D-33TP-10	I-H-PRO3-GS	
PRO311	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-30TP-11.5	D-33TP-11.5	I-H-PRO3-GS	
PRO313	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-30TP-13	D-33TP-13	I-H-PRO3-GS	
PRO315	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-30TP-15	D-33TP-15	I-H-PRO3-GS	
PRO318	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-30TP-18	D-33TP-18	I-H-PRO3-GS	

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities			Fully guided
		Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill Final Drill for soft bone	Final Drill for medium and dense bone	Final Drill for medium and dense bone	
PRO408	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-33TP-8.5	D-40TP-8.5	I-H-PRO-GS	
PRO410	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-33TP-10	D-40TP-10	I-H-PRO-GS	
PRO411	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-33TP-11.5	D-40TP-11.5	I-H-PRO-GS	
PRO413	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-33TP-13	D-40TP-13	I-H-PRO-GS	
PRO415	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-33TP-15	D-40TP-15	I-H-PRO-GS	
PRO418	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-33TP-18	D-40TP-18	I-H-PRO-GS	

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities			Fully guided
		Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill	3rd Drill Final Drill for soft bone	Final Drill for medium and dense bone	
PRO508	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-33TP-8.5	D-40TP-8.5	I-H-PRO-GS	
PRO510	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-33TP-10	D-40TP-10	I-H-PRO-GS	
PRO511	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-33TP-11.5	D-40TP-11.5	I-H-PRO-GS	
PRO513	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-33TP-13	D-40TP-13	I-H-PRO-GS	
PRO515	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-33TP-15	D-40TP-15	I-H-PRO-GS	
PRO518	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-33TP-18	D-40TP-18	I-H-PRO-GS	

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities				Fully guided Implant placement
		Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill Final Drill for soft bone	Final Drill for medium and dense bone	Final Drill for medium and dense bone	Final Drill for medium and dense bone	
PRO12D408	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-33TP-8.5	D-40TP-8.5	D-40TP-8.5	I-H-PRO12D-GS	
PRO12D410	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-33TP-10	D-40TP-10	D-40TP-10	I-H-PRO12D-GS	
PRO12D411	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-33TP-11.5	D-40TP-11.5	D-40TP-11.5	I-H-PRO12D-GS	
PRO12D413	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-33TP-13	D-40TP-13	D-40TP-13	I-H-PRO12D-GS	
PRO12D415	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-33TP-15	D-40TP-15	D-40TP-15	I-H-PRO12D-GS	
PRO12D418	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-33TP-18	D-40TP-18	D-40TP-18	I-H-PRO12D-GS	

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities				Fully guided Implant placement
		Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill	3rd Drill Final Drill for soft bone	Final Drill for medium and dense bone	Final Drill for medium and dense bone	
PRO12D508	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-33TP-8.5	D-40TP-8.5	D-50TP-8.5	I-H-PRO12D-GS	
PRO12D510	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-33TP-10	D-40TP-10	D-50TP-10	I-H-PRO12D-GS	
PRO12D511	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-33TP-11.5	D-40TP-11.5	D-50TP-11.5	I-H-PRO12D-GS	
PRO12D513	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-33TP-13	D-40TP-13	D-50TP-13	I-H-PRO12D-GS	
PRO12D515	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-33TP-15	D-40TP-15	D-50TP-15	I-H-PRO12D-GS	
PRO12D518	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-33TP-18	D-40TP-18	D-50TP-18	I-H-PRO12D-GS	

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities			Fully guided
		Initial Drill	Optional Drill for medium and dense cortical bone.	Final Drill for medium and dense bone	Final Drill for medium and dense bone	Implant placement	
IM-T3708	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-MT3708	D-20FM10	I-H-PRO-GS	
IM-T3710	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-MT3710	D-20FM10	I-H-PRO-GS	
IM-T3711	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-MT3711	D-20FM10	I-H-PRO-GS	
IM-T3713	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-MT3713	D-20FM15	I-H-PRO-GS	
IM-T3715	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-MT3715	D-20FM15	I-H-PRO-GS	

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities			Fully guided
		Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill Final Drill for soft bone	Final Drill for medium and dense bone	Implant placement	
IM-T4208	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-MT3708	D-MT4208	I-H-PRO-GS	
IM-T4210	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-MT3710	D-MT4210	I-H-PRO-GS	
IM-T4211	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-MT3711	D-MT4211	I-H-PRO-GS	
IM-T4213	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-MT3713	D-MT4213	I-H-PRO-GS	
IM-T4215	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-MT3715	D-MT4215	I-H-PRO-GS	
IM-T4218	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-MT3718	D-MT4218	I-H-PRO-GS	

Implant Code	Cortical Perforator	Drill sequence per bone densities			Fully guided		
		Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill Final Drill for soft bone	3rd Drill Final Drill for soft bone	Final Drill for medium and dense bone	Implant placement
IM-T5008	D-GS-CP	D-3SPADE-1.8M	D-30FM10	D-MT3708	D-MT4208	D-MT5008	I-H-PRO-GS
IM-T5010	D-GS-CP	D-3SPADE-1.8M	D-30FM10	D-MT3710	D-MT4210	D-MT5010	I-H-PRO-GS
IM-T5011	D-GS-CP	D-3SPADE-1.8M	D-30FM10	D-MT3711	D-MT4211	D-MT5011	I-H-PRO-GS
IM-T5013	D-GS-CP	D-3SPADE-1.8M	D-30FM15	D-MT3713	D-MT4213	D-MT5013	I-H-PRO-GS
IM-T5015	D-GS-CP	D-3SPADE-1.8M	D-30FM15	D-MT3715	D-MT4215	D-MT5015	I-H-PRO-GS

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities			Fully guided Implant placement
		Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill Final Drill for soft bone	Final Drill for medium and dense bone	Final Drill for medium and dense bone	
IM-T4208-12d	D-GS-CP4	D-3SPADE-1.8M	D-20T-M10	D-30T-M10	D-MT3708	D-MT4208	I-H-PRO12D-GS
IM-T4210-12d	D-GS-CP4	D-3SPADE-1.8M	D-20T-M10	D-30T-M10	D-MT3710	D-MT4210	I-H-PRO12D-GS
IM-T4211-12d	D-GS-CP4	D-3SPADE-1.8M	D-20T-M10	D-30T-M10	D-MT3711	D-MT4211	I-H-PRO12D-GS
IM-T4213-12d	D-GS-CP4	D-3SPADE-1.8M	D-20T-M15	D-30T-M15	D-MT3713	D-MT4213	I-H-PRO12D-GS
IM-T4215-12d	D-GS-CP4	D-3SPADE-1.8M	D-20T-M15	D-30T-M15	D-MT3715	D-MT4215	I-H-PRO12D-GS
IM-T4218-12d	D-GS-CP4	D-3SPADE-1.8M	D-20T-M15	D-30T-M15	D-MT3718	D-MT4218	I-H-PRO12D-GS

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities		Fully guided
		Initial Drill	Optional Drill for medium and dense cortical bone.	Final Drill for medium and dense bone	Implant placement	
IA-LH-35-8	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-L-35-8	I-HLH-35GS	
IA-LH-35-10	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-L-35-10	I-HLH-35GS	
IA-LH-35-11.5	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-L-35-11.5	I-HLH-35GS	
IA-LH-35-13	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-L-35-13	I-HLH-35GS	
IA-LH-35-16	D-GS-CP4	D-3SPADE-1.8M	D-20FM15	D-L-35-16	I-HLH-35GS	

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities			Fully guided
		Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill Final Drill for soft bone	Final Drill for medium and dense bone	Implant placement	
IA-LH-43-8	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-L-35-8	D-L-43-8	I-HLH-43GS	
IA-LH-43-10	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-L-35-10	D-L-43-10	I-HLH-43GS	
IA-LH-43-11.5	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-L-35-11.5	D-L-43-11.5	I-HLH-43GS	
IA-LH-43-13	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-L-35-13	D-L-43-13	I-HLH-43GS	
IA-LH-43-16	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-L-35-16	D-L-43-16	I-HLH-43GS	

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities			Fully guided
		Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill	3rd Drill Final Drill for soft bone	Final Drill for medium and dense bone	Implant placement
IA-LH-50-8	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-L-35-8	D-L-43-8	D-L-50-8	I-HLH-50GS
IA-LH-50-10	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-L-35-10	D-L-43-10	D-L-50-10	I-HLH-50GS
IA-LH-50-11.5	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-L-35-11.5	D-L-43-11.5	D-L-50-11.5	I-HLH-50GS
IA-LH-50-13	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-L-35-13	D-L-43-13	D-L-50-13	I-HLH-50GS
IA-LH-50-16	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-L-35-16	D-L-43-16	D-L-50-16	I-HLH-50GS

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities			Final Drill NOT through the guide
		Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill	3rd Drill	4th Drill Final Drill for soft bone	Final Drill for medium and dense bone
IA-LH-60-8	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-L-35-8	D-L-43-8	D-L-50-8	D-L-60-8
IA-LH-60-10	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-L-35-10	D-L-43-10	D-L-50-10	D-L-60-10
IA-LH-60-11.5	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-L-35-11.5	D-L-43-11.5	D-L-50-11.5	D-L-60-11.5
IA-LH-60-13	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-L-35-13	D-L-43-13	D-L-50-13	D-L-60-13
IA-LH-60-16	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-L-35-16	D-L-43-16	D-L-50-16	D-L-60-16

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities			Fully guided Implant placement
		Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill Final Drill for soft bone	Final Drill for medium and dense bone	Final Drill for medium and dense bone	
IA43-12d-10	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-L-35-10	D-L-43-10	H-L43-12D-GS	
IA43-12d-11.5	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-L-35-11.5	D-L-43-11.5	H-L43-12D-GS	
IA43-12d-13	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-L-35-13	D-L-43-13	H-L43-12D-GS	
IA43-12d-16	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-L-35-16	D-L-43-16	H-L43-12D-GS	

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities			Fully guided Implant placement
		Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill	3rd Drill Final Drill for soft bone	Final Drill for medium and dense bone	
IA50-12d-10	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-L-35-10	D-L-43-10	D-L-50-10	H-L50-12D-GS
IA50-12d-11.5	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-L-35-11.5	D-L-43-11.5	D-L-50-11.5	H-L50-12D-GS
IA50-12d-13	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-L-35-13	D-L-43-13	D-L-50-13	H-L50-12D-GS
IA50-12d-16	D-GS-CP	D-3SPADE-1.8M	D-20FM15	D-L-35-16	D-L-43-16	D-L-50-16	H-L50-12D-GS

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities		Fully guided
		Initial Drill	Optional Drill for medium and dense cortical bone.	Final Drill	Final Drill	Implant placement
DCT3009	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-DCT3009	D-DCT3009	I-HDC3-GS
DCT3011	D-GS-CP4	D-3SPADE-1.8M	D-20FM10	D-DCT3011	D-DCT3011	I-HDC3-GS
DCT3013	D-GS-CP4	D-3SPADE-1.8M	D-DC20/D-20FM15	D-DC20/D-20FM15	D-DCT3013	I-HDC3-GS

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities			Fully guided
		Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill	Final Drill for soft bone	Final Drill for medium and dense bone	Implant placement
DCT3508	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-DCT3508	D-DCT3508	D-DCT3508	I-HDC4-GS
DCT3509	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-DCT3509	D-DCT3509	D-DCT3509	I-HDC4-GS
DCT3511	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-DCT3511	D-DCT3511	D-DCT3511	I-HDC4-GS
DCT3513	D-GS-CP	D-3SPADE-1.8M	D-DC20/D-20FM15	D-DCT3513	D-DCT3513	D-DCT3513	I-HDC4-GS
DCT3515	D-GS-CP	D-3SPADE-1.8M	D-DC20/D-20FM15	D-DCT3515	D-DCT3515	D-DCT3515	I-HDC4-GS

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities				Fully guided
		Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill	3rd Drill	Final Drill for soft bone	Final Drill for medium and dense bone	Implant placement
DCT4006	D-GS-CP	D-3SPADE-GS	D-16-T	D-DCT4006	D-DCT4006	D-DCT4006	D-DCT4006	I-HDC4-GS
DCT4008	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-DCT4008	D-DCT4008	D-DCT4008	D-DCT4008	I-HDC4-GS
DCT4009	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-DCT4009	D-DCT4009	D-DCT4009	D-DCT4009	I-HDC4-GS
DCT4011	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-DCT4011	D-DCT4011	D-DCT4011	D-DCT4011	I-HDC4-GS
DCT4013	D-GS-CP	D-3SPADE-1.8M	D-DC20/D-20FM15	D-DCT4013	D-DCT4013	D-DCT4013	D-DCT4013	I-HDC4-GS
DCT4015	D-GS-CP	D-3SPADE-1.8M	D-DC20/D-20FM15	D-DCT4015	D-DCT4015	D-DCT4015	D-DCT4015	I-HDC4-GS

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities				Fully guided	
		Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill	3rd Drill	4th Drill	Final Drill for soft bone	Final Drill for medium and dense bone	Implant placement
DCT5009	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-DCT5009	D-DCT5009	D-DCT5009	D-DCT5009	D-DCT5009	I-HDC5-GS
DCT5011	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-DCT5011	D-DCT5011	D-DCT5011	D-DCT5011	D-DCT5011	I-HDC5-GS
DCT5013	D-GS-CP	D-3SPADE-1.8M	D-DC20/D-20FM15	D-DCT5013	D-DCT5013	D-DCT5013	D-DCT5013	D-DCT5013	I-HDC5-GS
DCT5015	D-GS-CP	D-3SPADE-1.8M	D-DC20/D-20FM15	D-DCT5015	D-DCT5015	D-DCT5015	D-DCT5015	D-DCT5015	I-HDC5-GS

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities			
		Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill	3rd Drill Final Drill for soft bone	4th Drill Final Drill for soft bone	
DCT4008-12D	D-GS-CP	D-3SPADE-1.8M	D-20FM10		D-DCT3508	D-DCT4008	
DCT4009-12D	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-DC T3009	D-DCT3509	D-DCT4009	
DCT4011-12D	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-DC T3011	D-DCT3511	D-DCT4011	
DCT4013-12D	D-GS-CP	D-3SPADE-1.8M	D-DC20/D-20FM15	D-DC T3013	D-DCT3513	D-DCT4013	
DCT4015-12D	D-GS-CP	D-3SPADE-1.8M	D-DC20/D-20FM15		D-DCT3515	D-DCT4015	

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities			
		Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill	3rd Drill	4th Drill Final Drill for soft bone	Final Drill for medium and dense bone
DCT5009-12D	D-GS-CP	D-3SPADE-1.8M	D-20FM10		D-DCT3509	D-DCT4009	D-DCT5009
DCT5011-12D	D-GS-CP	D-3SPADE-1.8M	D-20FM10	D-DC T3009	D-DCT3511	D-DCT4011	D-DCT5011
DCT5013-12D	D-GS-CP	D-3SPADE-1.8M	D-DC20/D-20FM15	D-DC T3011	D-DCT3513	D-DCT4013	D-DCT5013
DCT5015-12D	D-GS-CP	D-3SPADE-1.8M	D-DC20/D-20FM15	D-DC T3013	D-DCT3515	D-DCT4015	D-DCT5015

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities	
		Initial Drill	Final Drill	Optional Drill for soft bone.	Final Drill for medium and dense bone
IVEX30-3710	D-GS-CP4	D-3SPADE-1.8M	D-IV3710GS		
IVEX30-3711	D-GS-CP4	D-3SPADE-1.8M	D-IV3711GS		
IVEX30-3713	D-GS-CP4	D-3SPADE-1.8M	D-IV3713GS		
IVEX30-3715	D-GS-CP4	D-3SPADE-1.8M	D-IV3715GS		

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities	
		Initial Drill	Optional Drill for soft bone.	Final Drill for medium and dense bone	Final Drill for medium and dense bone
IVEX35-4510	D-GS-CP	D-3SPADE-1.8M	D-IV3710GS	D-IV4510GS	
IVEX35-4511	D-GS-CP	D-3SPADE-1.8M	D-IV3711GS	D-IV4511GS	
IVEX35-4513	D-GS-CP	D-3SPADE-1.8M	D-IV3713GS	D-IV4513GS	
IVEX35-4515	D-GS-CP	D-3SPADE-1.8M	D-IV3715GS	D-IV4515GS	

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities		
		Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill for soft bone	Final Drill for medium and dense bone	Final Drill for medium and dense bone
IVEX40-5010	D-GS-CP	D-3SPADE-1.8M	D-IV3710GS	D-IV4510GS	D-IV5010GS	
IVEX40-5011	D-GS-CP	D-3SPADE-1.8M	D-IV3711GS	D-IV4511GS	D-IV5011GS	
IVEX40-5013	D-GS-CP	D-3SPADE-1.8M	D-IV3713GS	D-IV4513GS	D-IV5013GS	
IVEX40-5015	D-GS-CP	D-3SPADE-1.8M	D-IV3715GS	D-IV4515GS	D-IV5015GS	

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities			Final Drill NOT through the guide Final Drill for medium and dense bone
		Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill	3rd Drill for soft bone	Final Drill for soft bone	
IVEX52-6010	D-GS-CP	D-3SPADE-1.8M	D-IV3710GS	D-IV4510GS	D-IV5010GS	D-IV6010GS	
IVEX52-6011	D-GS-CP	D-3SPADE-1.8M	D-IV3711GS	D-IV4511GS	D-IV5011GS	D-IV6011GS	
IVEX52-6013	D-GS-CP	D-3SPADE-1.8M	D-IV3713GS	D-IV4513GS	D-IV5013GS	D-IV6013GS	
IVEX52-6015	D-GS-CP	D-3SPADE-1.8M	D-IV3715GS	D-IV4515GS	D-IV5015GS	D-IV6015GS	

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities	
		Initial Drill	Final Drill for medium and dense bone	Optional Drill for soft bone.	Final Drill for medium and dense bone
IV-EX3012D-3711	D-GS-CP4	D-3SPADE-1.8M	D-IV3711GS		
IV-EX3012D-3713	D-GS-CP4	D-3SPADE-1.8M	D-IV3713GS		
IV-EX3012D-3715	D-GS-CP4	D-3SPADE-1.8M	D-IV3715GS		

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities	
		Initial Drill	Optional Drill for soft bone.	Final Drill for medium and dense bone	Final Drill for medium and dense bone
IV-EX3512D-4510	D-GS-CP	D-3SPADE-1.8M	D-IV3710GS	D-IV4510GS	
IV-EX3512D-4511	D-GS-CP	D-3SPADE-1.8M	D-IV3711GS	D-IV4511GS	
IV-EX3512D-4513	D-GS-CP	D-3SPADE-1.8M	D-IV3713GS	D-IV4513GS	
IV-EX3512D-4515	D-GS-CP	D-3SPADE-1.8M	D-IV3715GS	D-IV4515GS	

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities		
		Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill Final Drill for soft bone	Final Drill for medium and dense bone	Final Drill for medium and dense bone
IV-EX4012D-5011	D-GS-CP	D-3SPADE-1.8M	D-IV3711GS	D-IV4511GS	D-IV5011GS	
IV-EX4012D-5013	D-GS-CP	D-3SPADE-1.8M	D-IV3713GS	D-IV4513GS	D-IV5013GS	
IV-EX4012D-5015	D-GS-CP	D-3SPADE-1.8M	D-IV3715GS	D-IV4515GS	D-IV5015GS	

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities			Final Drill NOT through the guide Final Drill for medium and dense bone
		Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill	3rd Drill Final Drill for soft bone	Final Drill for medium and dense bone	
IV-EX5212D-6011	D-GS-CP	D-3SPADE-1.8M	D-IV3711GS	D-IV4511GS	D-IV5011GS	D-IV6011GS	
IV-EX5212D-6013	D-GS-CP	D-3SPADE-1.8M	D-IV3713GS	D-IV4513GS	D-IV5013GS	D-IV6013GS	
IV-EX5212D-6015	D-GS-CP	D-3SPADE-1.8M	D-IV3715GS	D-IV4515GS	D-IV5015GS	D-IV6015GS	

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities		Fully guided	
		Initial Drill	Final Drill	Final Drill for medium and dense cortical bone.	Final Drill for medium and dense bone	Implant placement	Implant placement
IV-DC30-3710	D-GS-CP4	D-3SPADE-1.8M	D-IV3710GS				IHDC3-GS
IV-DC30-3711	D-GS-CP4	D-3SPADE-1.8M	D-IV3711GS				IHDC3-GS
IV-DC30-3713	D-GS-CP4	D-3SPADE-1.8M	D-IV3713GS				IHDC3-GS
IV-DC30-3715	D-GS-CP4	D-3SPADE-1.8M	D-IV3715GS				IHDC3-GS

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities		Fully guided	
		Initial Drill	2nd Drill for medium and dense cortical bone.	Final Drill for medium and dense bone	Final Drill for medium and dense bone	Implant placement	Implant placement
IV-DC35-4510	D-GS-CP	D-3SPADE-1.8M	D-IV4510GS	D-IV4510GS	D-IV4510GS		IHDC4-GS
IV-DC35-4511	D-GS-CP	D-3SPADE-1.8M	D-IV3711GS	D-IV4511GS	D-IV4511GS		IHDC4-GS
IV-DC35-4513	D-GS-CP	D-3SPADE-1.8M	D-IV3713GS	D-IV4513GS	D-IV4513GS		IHDC4-GS
IV-DC35-4515	D-GS-CP	D-3SPADE-1.8M	D-IV3715GS	D-IV4515GS	D-IV4515GS		IHDC4-GS

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities				Fully guided
		Initial Drill	2nd Drill for medium and dense cortical bone.	3rd Drill	Final Drill for soft bone	Final Drill for medium and dense bone	Final Drill for medium and dense bone	
IV-DC40-5010	D-GS-CP	D-3SPADE-1.8M	D-IV3710GS	D-IV4510GS	D-IV5010GS			IHDC4-GS
IV-DC40-5011	D-GS-CP	D-3SPADE-1.8M	D-IV3711GS	D-IV4511GS	D-IV5011GS			IHDC4-GS
IV-DC40-5013	D-GS-CP	D-3SPADE-1.8M	D-IV3713GS	D-IV4513GS	D-IV5013GS			IHDC4-GS
IV-DC40-5015	D-GS-CP	D-3SPADE-1.8M	D-IV3715GS	D-IV4515GS	D-IV5015GS			IHDC4-GS

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities				Final Drill NOT through the guide Final Drill for medium and dense bone
		Initial Drill	2nd Drill for medium and dense cortical bone.	3rd Drill	4th Drill	Final Drill for soft bone.	Final Drill for soft bone.	
IV-DC50-6010	D-GS-CP	D-3SPADE-1.8M	D-IV3710GS	D-IV4510GS	D-IV5010GS			D-IV6010GS
IV-DC50-6011	D-GS-CP	D-3SPADE-1.8M	D-IV3711GS	D-IV4511GS	D-IV5011GS			D-IV6011GS
IV-DC50-6013	D-GS-CP	D-3SPADE-1.8M	D-IV3713GS	D-IV4513GS	D-IV5013GS			D-IV6013GS
IV-DC50-6015	D-GS-CP	D-3SPADE-1.8M	D-IV3715GS	D-IV4515GS	D-IV5015GS			D-IV6015GS


Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities	
		Initial Drill	Optional Drill for medium and dense cortical bone.	Optional Drill for medium and dense cortical bone.	Final Drill for medium and dense bone
IV-DC3512D-4511	D-GS-CP	D-3SPADE-1.8M	D-IV3711GS	D-IV4511GS	
IV-DC3512D-4513	D-GS-CP	D-3SPADE-1.8M	D-IV3713GS	D-IV4513GS	
IV-DC3512D-4515	D-GS-CP	D-3SPADE-1.8M	D-IV3715GS	D-IV4515GS	

Implant Code	Cortical Perforator	Initiate the osteotomy		Drill sequence per bone densities		
		Initial Drill	Optional Drill for medium and dense cortical bone.	2nd Drill Final Drill for soft bone	Final Drill, for medium and dense bone	
IV-DC4012D-5011	D-GS-CP	D-3SPADE-1.8M	D-IV3711GS	D-IV4511GS	D-IV5011GS	
IV-DC4012D-5013	D-GS-CP	D-3SPADE-1.8M	D-IV3713GS	D-IV4513GS	D-IV5013GS	
IV-DC4012D-5015	D-GS-CP	D-3SPADE-1.8M	D-IV3715GS	D-IV4515GS	D-IV5015GS	

SIREAL NARROW



Implant planning



Sleeves

Ø5.1 (outer) / Ø4.2 (inner)



Ø4.1 offset sleeve

9 mm; 10.5 mm; 12 mm; 14 mm



Commission Guide



SIREAL NARROW

UNIVERSAL GUIDED SURGERY TOOL

The universal guided surgery tool has a diameter of 4.1 mm. The $\varnothing 4.1$ mm offset sleeves only fits directly on the threaded part of the I-DE-GS4. The shaft of the I-DE-GS4 guides through the narrow $\varnothing 4.2$ mm guide sleeves.



NARROW OFFSET SLEEVES




The universal guided surgery tool threads are used to engage the four varied Offset Sleeves for setting the depth of the tool. **NOTE:** the narrow offset sleeves do not have a cylinder attached as with the SIREAL standard offset sleeves. The shaft of the universal guided surgery tool acts as a $\varnothing 4.1$ mm cylinder, which fits through the narrow guide sleeves.



Four offset sleeves are available: **9 mm**, **10.5 mm**, **12 mm** and **14 mm**. This is to accommodate the patient's vertical opening or adjacent teeth height that could interfere with the guide sleeve.


NARROW GUIDE SLEEVES



UNIVERSAL GUIDED SURGERY TOOL




Ø5.1 STANDARD OFFSET SLEEVES



I-GSS05-9	Ø5.1 x 9 mm Offset Sleeve	(YELLOW)
I-GSS05-10.5	Ø5.1 x 10.5 mm Offset Sleeve	(GREEN)
I-GSS05-12	Ø5.1 x 12 mm Offset Sleeve	(BLUE)
I-GSS05-14	Ø5.1 x 14 mm Offset Sleeve	(PURPLE)

ADDITIONAL INFORMATION




- if driving on the latch, the recommended MAX torque is 40 Ncm.
- if driving on the W&H hex, the recommended MAX torque is 70Ncm Recommended to order 2.

The addition of the Ø5.1 Standard offset sleeve converts the cylinder diameter of the universal guided surgical tool from Ø4.1 mm to Ø5.1 mm.

UNIVERSAL GUIDED SURGERY TOOL

Fixed 10.5 mm offset



Features a built-in standard offset of 10.5 mm, eliminating the need for offset sleeves.

GUIDED SLEEVES

Ø6.2 (outer) / Ø5.2 (inner)



I-GS-6252-L4	4 mm	Closed Guide Sleeve
I-GS-6252-L9	9 mm	Closed Guide Sleeve
I-GSC-6252-L4	4 mm	Open Guide Sleeve
I-GSC-6252-L9	9 mm	Open Guide Sleeve

Why and when to use a C-guide sleeve?

- better irrigation at the osteotomy site.
- side entry: when the patient has limited vertical space (especially posteriorly).
- Co-Axis® compatibility
- fixture mount screw access.

Benefits of a 9 mm sleeve height

The SIREAL universal guided surgery tool will be guided by a longer guide sleeve length, which increases the guidance and stability.

NOTE: 9 mm guide sleeves can not be used with a 9 mm offset.


TISSUE CUTTERS



I-TC-GS4	Tissue Cutter Ø4.5 outer Ø4.1 inner
I-TC-GS5	Tissue Cutter Ø5.0 outer Ø4.7 inner

Can only be used with the standard offset sleeves and the standard guide sleeves.

CORTICAL PERFORATOR



D-GS-CP	Cortical Perforator Guide Surgery
---------	-----------------------------------

- compatible with Ø5.1 guided sleeves.
- recommended to be used as initial drill for 6 mm implants.
- do not use D-20TM10 to prevent preparing the osteotomy too deep.

UNIVERSAL GUIDED SURGERY TOOL



I-DE-GS4

Ø4.1 NARROW OFFSET SLEEVES



I-GSS04-9



I-GSS04-10.5



I-GSS04-12




I-GSS04-14


I-GSS04-9	Ø4.1 x 9 mm Offset Sleeve	(YELLOW)
I-GSS04-10.5	Ø4.1 x 10.5 mm Offset Sleeve	(GREEN)
I-GSS04-12	Ø4.1 x 12 mm Offset Sleeve	(BLUE)
I-GSS04-14	Ø4.1 x 14 mm Offset Sleeve	(PURPLE)

NARROW GUIDED SLEEVES

Ø5.1 (outer) / Ø4.2 (inner)




I-GS-5142-L4



I-GS-5142-L9



I-GSC-5142-L4



I-GSC-5142-L9


I-GS-5142-L4	4 mm	Closed Guide Sleeve
I-GS-5142-L9	9 mm	Closed Guide Sleeve
I-GSC-5142-L4	4 mm	Open Guide Sleeve
I-GSC-5142-L9	9 mm	Open Guide Sleeve

TISSUE CUTTER



I-TC-GS3	Tissue Cutter Ø3.9 outer Ø3.5 inner
----------	-------------------------------------

CORTICAL PERFORATOR




D-GS-CP4	Cortical Perforator Guide Surgery
----------	-----------------------------------

ADDITIONAL INFORMATION

- if driving on the latch, the recommended MAX torque is 40 Ncm.
- if driving on the W&H hex, the recommended MAX torque is 70Ncm Recommended to order 2.

The narrow offset sleeves utilise the Ø4.1 cylinder of the SIREAL tool to engage the NARROW sleeves.



Why and when to use a C-guide sleeve?

- better irrigation at the osteotomy site.
- side entry: when the patient has limited vertical space (especially posteriorly).
- Co-Axis® compatibility
- fixture mount screw access.

Benefits of a 9 mm guide sleeve height

The SIREAL universal guided surgery tool will be guided by a longer sleeve length, which increases the guidance and stability.

NOTE: 9 mm guide sleeves can not be used with a 9 mm offset.

Can only be used with the narrow offset sleeves and the narrow guide sleeves.

- compatible with Ø4.1 guided sleeves.
- recommended to be used as initial drill for 6 mm implants.
- do not use D-20TM10 to prevent preparing the osteotomy too deep.

GUIDE FIXATION PIN / SLEEVE / DRILL



I-D12-GS

Guide Fixation Pin



I-GS-2513-L6

Guide Fixation Pin Sleeve

I-D12-GS	Guide Fixation Pin
I-GS-2513-L6	Guide Fixation Pin Sleeve



D-12T-M15

Fixation Pin Drill

D-12T-M15	Twist Drill $\varnothing 1.2 \times 15 \text{ mm}$
-----------	--

INSTRUMENT TRAY INSERTS (to keep the Offset sleeves upright)




I-GS-TP

I-GS-TP	Guide Sleeve Pin GS Tray
---------	--------------------------

ADDITIONAL INFORMATION

Place at least 3 pins to stabilize the guide.



The guide pins, are inserted through the fixation guide sleeve, after drilling a hole with the D-12T-M15.

Due to the Narrow diameter of the $\varnothing 1.2$ twist drill:

- avoid lateral movement whilst drilling.
- have a spare available.








Compatible with silicone insert trays.

SIREAL ORDER GUIDE

INITIAL DRILLS

ADDITIONAL INFORMATION

	D-16-T Pilot drill - side cutting Ø1.6 mm	Drill is indicated when placing a 6 mm length implant. Do not use the D-20T-M10 initial drill.
	D-20T-M10 Initial drill - twist drill Ø2.0 mm	10 mm twist drill is indicated for implant lengths 11.5 mm and shorter. (Applicable to all implant ranges)
	D-20T-M15 Initial drill - twist drill Ø2.0 mm	15 mm twist drill is indicated for implant lengths 12 mm and longer. (Applicable to all implant ranges)
	D-DC20 Initial drill - twist drill Ø2.0 mm	D-DC20 drill is indicated for implant lengths 13 mm and longer. For shorter implants, use D-20T-M10.
	D-220C Initial drill - twist drill Ø2.0 mm	Indicated for IT implant lengths 12 mm and longer. For shorter implants, use D-20T-M10.

GUIDED DRILLS FOR INVERTA® GUIDED SURGERY

D-IV3708GS	Drill Taper IV GS Ø3.75 x 8 mm
D-IV3710GS	Drill Taper IV GS Ø3.75 x 10 mm
D-IV3711GS	Drill Taper IV GS Ø3.75 x 11.5 mm
D-IV3713GS	Drill Taper IV GS Ø3.75 x 13 mm
D-IV3715GS	Drill Taper IV GS Ø3.75 x 15 mm
D-IV4508GS	Drill Taper IV GS Ø4.5 x 8 mm
D-IV4510GS	Drill Taper IV GS Ø4.5 x 10 mm
D-IV4511GS	Drill Taper IV GS Ø4.5 x 11.5 mm
D-IV4513GS	Drill Taper IV GS Ø4.5 x 13 mm
D-IV4515GS	Drill Taper IV GS Ø4.5 x 15 mm
D-IV4518GS	Drill Taper IV GS Ø4.5 x 18 mm
D-IV5008GS	Drill Taper IV GS Ø5 x 8 mm
D-IV5010GS	Drill Taper IV GS Ø5 x 10 mm
D-IV5011GS	Drill Taper IV GS Ø5 x 11.5 mm
D-IV5013GS	Drill Taper IV GS Ø5 x 13 mm
D-IV5015GS	Drill Taper IV GS Ø5 x 15 mm
D-IV5018GS	Drill Taper IV GS Ø5 x 18 mm
D-IV6010GS	Drill Taper IV GS Ø6 x 10 mm
D-IV6011GS	Drill Taper IV GS Ø6 x 11.5 mm
D-IV6013GS	Drill Taper IV GS Ø6 x 13 mm
D-IV6015GS	Drill Taper IV GS Ø6 x 15 mm
D-IV6018GS	Drill Taper IV GS Ø6 x 18 mm



CAUTION:

- all INVERTA® tapered drills used with SIREAL must be part of the "GS" drill range.
- utilising the standard INVERTA® tapered drills will result in deeper site preparation than planned.

FULLY GUIDED IMPLANT PLACEMENT TOOLS

DC (Deep Conical)



IHDC3-GS IHDC4-GS IHDC5-GS

IHDC3-GS	Ø3.0 Placement Tool
IHDC4-GS	Ø4.0 Placement Tool
IHDC5-GS	Ø5.0 Placement Tool

Internal Hex and PROVATA®



IH-PRO-GS IH-PRO3-GS IH-PRO12D-GS

IH-PRO-GS	Placement Tool
IH-PRO3-GS	Placement Tool
IH-PRO12D-GS	Placement Tool, Co-Axis®

TRI-NEX®



IHLH-35GS IHLH-43GS IHLH-50GS IL43-12D-GS IL50-12D-GS

IHLH-35GS	Ø3.5 Placement Tool
IHLH-43GS	Ø4.3 Placement Tool
IHLH-50GS	Ø5.0 Placement Tool
IL43-12D-GS	Ø4.3 Placement Tool, Co-Axis®
IL50-12D-GS	Ø5.0 Placement Tool, Co-Axis®

ADDITIONAL INFORMATION



CAUTION:

- all fully guided implant placement tools used with the SIREAL universal guide tool must be part of the "GS" range.
- utilising the standard implant placement tools could result in deeper implant placement than planned.

PEEK bits



I-PBIT-H16



I-PBIT-L18

Insertion tools supplied with PEEK bit.

PEEK bits



I-PBIT-L18

for
Ø3.5 mm interface
instrumentation only.



I-PBIT-L20

for
Ø4.3 mm and Ø5.0 mm
interface
instrumentation.





















Insertion tool supplied with PEEK bit.





Important: the PEEK bits should be replaced on a regular basis.

Items sold separately.
General wear and tear are to be expected with regular use.

EXPLANATION OF SYMBOLS

The following symbols are used on packaging labels and they indicate the following:

- 1  Manufacturer
- 2  Colour code indicating platform diameter
- 3  Implant image
- 4  Implant details and size
- 5  Sterilisation using irradiation
-  European representative
-  Catalogue number
-  Batch code
-  Do not resterilise
-  Consult instruction for use
-  Do not reuse
-  CE mark and notified body number
-  Use by date
-  Date of manufacture
-  Do not use if package is damaged
-  Identifies the product as a medical device
-  MR Conditional / Magnetic Resonance Conditional
-  Single sterile barrier system
-  Double sterile barrier
- 6 **2D Bar coding**
Contains the GTIN, Use by date and LOT number
- 7 **Patient sticker for documentation purposes**
(to be used by health care provider on patient file)
- 8  Prescription device
CAUTION: FEDERAL LAW RESTRICTS THE DEVICE TO SALE BY OR ON THE ORDER OF A LICENCED HEALTH CARE PROVIDER.
- 9 **Product description**
(translated as per international standards)



For more information on Instructions for Use of our products, please scan the below,



or visit our website southernimplants.com/ifu

Images are for illustration purposes only and do not necessarily accurately represent the product.

For more information scan the below



to contact your Southern Implants Representative
or visit southernimplants.com



South Africa - Headquarters 1 Albert Road, Irene, RSA
T: +27-12-667-1046 | **E:** info@southernimplants.com

EC REP

Southern Implants Europe AB, Holmgatan 30, S-791 71, Falun, Sweden
E: ecrep@southernimplants.com

Subsidiaries

Australia

Southern Implants Australia
T: +61-2-8076-9337
E: info@southernimplants.com.au

Spain and Portugal

Southern Implants Ibérica
T: +34 935 053 507
E: info@southernimplants.es

United Kingdom and Ireland

Southern Implants UK
T: +44-20-8899-6845 / 6 / 7
E: info@southernimplants.co.uk

USA and Canada

Southern Implants North America Inc.
T: +1-561-472-0990
E: customercare@southernimplants.com

Southern Implants are distributed world-wide, please visit southernimplants.com for a list of Distributors.
For the latest and full assortment of Southern Implant® products, visit our current resources library at southernimplants.com.