**SOUTHERNIMPLANTS®**

Innovative Treatment Solutions

<b>English</b>	<b>INSTRUCTIONS FOR USE: Southern Implants® Reusable Instruments</b>
<b>Español</b>	<b>INSTRUCCIONES DE USO: Southern Implants® Instrumentos reutilizables</b>
<b>Italiano</b>	<b>ISTRUZIONI PER L'USO: Southern Implants® Strumenti riutilizzabili</b>
<b>Français</b>	<b>MODE D'EMPLOI : Southern Implants® Implants Instruments réutilisables</b>
<b>Deutsch</b>	<b>GEBRAUCHSANWEISUNG: Southern Implants® Wiederverwendbare Instrumente</b>
<b>Português</b>	<b>INSTRUÇÕES DE UTILIZAÇÃO: Southern Implants® Instrumentos Reutilizáveis</b>

ARCHIVED

**South Africa - Headquarters:** 1 Albert Road, Irene, 0062, RSA  
T: +27-12-667-1046 | E: [info@southernimplants.com](mailto:info@southernimplants.com)

EC REP

**Southern Implants Europe AB:** Holmgatan 30, S-791 71 Falun, Sweden  
T: +46 23 13300 | E: [ecrep@southernimplants.com](mailto:ecrep@southernimplants.com)

#### Subsidiaries

##### Australia

Southern Implants Australia  
T: +61-(0)-8-9466-2627  
E: [info@southernimplants.com.au](mailto:info@southernimplants.com.au)

##### Spain and Portugal

Southern Implants Iberica  
T: +34 935 053 507  
E: [info@southernimplants.es](mailto:info@southernimplants.es)

##### United Kingdom and Ireland

Southern Implants UK  
T: +44-20-8899-6845 / 6 / 7  
E: [info@southernimplants.co.uk](mailto:info@southernimplants.co.uk)

##### USA and Canada

Southern Implants North America Inc.  
T: +1-561-472-0990  
E: [customer care@southernimplants.com](mailto:customer care@southernimplants.com)

## Cleaning & Sterilization Procedure Guidelines

This document gives general guidance on how Southern Implants devices (components and instrumentation) that are suitable for cleaning and sterilization should be processed. These instructions are based on validated processes used to render a product free from viable microorganisms.

**NOTE:** The final responsibility for cleaning, disinfecting and sterilization techniques lies with the end-user (ISO17665).

### Do not re-use

Re-using single use devices may affect the integrity of the design, i.e. critical dimension may change - affecting fit or efficiency of the product.

### Do not re-sterilize

Changes in mechanical, physical or chemical characteristics may occur, affecting the design and performance of the item.

### Non-sterile

Indicates that the device is supplied clean but not sterile.

### Sterile

Indicates a device supplied sterilized by using irradiation.

### Do not use if package is opened or damaged

Medical devices should not be used if packaging is damaged or opened.

## Limitations on Re-usable items

A direct value for reusable instruments can not be given. Frequent processing may have minor effects on the instruments. The product life is normally determined by wear and damage during use, thus instruments if properly cared for and inspected after each use, can be reused many times. Maintain a checklist for these instruments recording the number of uses.

Prior to re-processing the device it should be thoroughly inspected and tested to determine its suitability for re-use.

**NOTE:** During use, handle drills and instruments with sterile tweezers to minimize contamination of the instrument tray and risk of damage to sterile surgical gloves.

## Containment

As soon as practically possible, remove all visible residue after use (bone, blood or tissue), by immersing the instrument in cold water. (Dried soil is difficult to remove).

## Pre-Cleaning

Disassemble instruments from hand pieces and all connecting parts from instruments i.e. ratchets, in order to clean soil from obstructed areas. Remove the PEEK bits from placement tools. Rinse with lukewarm water for 3 minutes, and remove hardened debris with a soft nylon brush. Avoid mechanical damage during cleaning.

## Manual Cleaning or Automated Cleaning

Prepare an ultrasonic bath with suitable detergent (i.e Steritech instrument cleaner -5% dilution), sonicate for 20 minutes (Alternative methods can be used if proven by the end user). Rinse with purified / sterile water.

**NOTE:** Always follow the instructions for use of the manufacturers of cleaning agents and disinfectants.

Load devices into a thermo-disinfector. Run the cleaning and disinfection cycle, followed by the drying cycle.

## Drying

Dry the instruments inside / outside with filtered compressed air or single use, lint free wipes. Pack the instruments as quickly as possible after removal into the storage container. If additional drying is necessary, dry in a clean location.

## Inspection

Do a visual inspection of the items to check for damages on instruments.

## Packaging

Use the correct packaging material as indicated for steam sterilization to ensure sterility is maintained. Double

packaging is recommended. Where appropriate, the cleaned, disinfected and checked devices can be assembled and placed in instrument trays as applicable. Instrument trays can be double wrapped or placed in sterilization bags.

## Sterilization

Southern Implants recommends the following procedure to sterilise the instruments prior to use/re-use:

Methods to sterilise the surgical instruments:

1. Pre-vacuum Sterilisation method: Steam sterilise the instruments at 132°C (270°F) at 180-220kPa for 4 minutes. Dry for at least 20 minutes in the chamber. Only an approved wrap or pouch for steam sterilisation must be used.
2. Pre-vacuum sterilisation method: Wrapped, steam sterilise at 135°C (275°F) for 3 minutes. Dry for 20 minutes in the chamber. Use a wrap or pouch that is cleared for the indicated steam sterilisation cycle.

**NOTE:** Users in the USA must ensure that the steriliser, wrap or pouch, and all steriliser accessories are cleared by the FDA, for the intended sterilisation cycle.

## Storage

Maintain packaging integrity to ensure sterility in storage. Packaging should be completely dry before storage to avoid corrosion and degradation of cutting edges.

## Intended use

### Osteotomy Guide

Are a non-sterile surgical instrument intended to be attached to a portion of bone (e.g. pelvis, femur, tibia, ulna) during an osteotomy procedure to outline the desired cut-line and mechanically support/guide the cutting device (e.g. scissoring cutting saw).

### Reusable Dental Implant Depth/Angle Gauge

A device designed to be used during the planning and post-osteotomy stages of a dental implantation procedure, prior to surgical implant placement/fixation, to measure/determine various lengths (e.g. drill depth of the osteotomy cavity, mucosal height, distance between implants/teeth for determining the mesiodistal position of the implant) and/or angles (e.g. direction of the perforation in relation to the antagonist arch), typically to help select an appropriate implant and/or final prosthetic abutment to be installed. It is typically graduated and may be a hand-held instrument or a stand-alone device placed intraorally; some types are intended for use during x-ray imaging.

### Dental osteotomes

A hand-held, dental, surgical instrument used during dental osteotomy for shaping and condensing bone, typically when performing a dental implant surgical procedure. It is held by the dentist who will use a surgical mallet or hammer to manually impart an impacting force to the proximal

end of the instrument. The device will have a straight or bayonet-shaped handle leading to a round or concave tip at the distal end, which may be marked with graduations to indicate depth differences. It is used for ridge expanding, bone condensing and sinus floor elevation procedures. This device is made of high-grade stainless steel.

**Basic Dental Retractors**

A hand-held, non-self-retaining dental instrument intended to be used to manually displace soft tissues of the oral cavity to improve their visualization and access, and to afford them protection during oral surgical procedures. It is typically a metallic device and includes a shaped/hooked blade with a handle/holding portion.

**Surgical drill guide - reusable**

A hand-held tubular surgical device designed to simultaneously place, angle, and guide a rotating surgical drill on hard tissue (e.g., bone), and to protect the surrounding soft tissue and dentition while the drill is in operation. It can have a handle at its proximal end for holding and manipulation, or it may be held by another device (e.g. orthopaedic jig, 3D-printed surgical guide). The device is typically dedicated to drill bit size, made of metal, and may have graduations along its length; it is not custom-made for a specific patient.

**Dental Torque Wrenches**

A hand-held dental tool intended for the manual application of rotation to tighten/loosen/remove a device (e.g. dental implant, spatial register (localizer/tracker/calibrator)) fitted into or close to the oral cavity, and/or to a surgical instrument, in association with a dental surgical procedure (e.g. dental guided surgery registration). It may have a built-in ratchet function, with forward and reverse direction, and a calibrated scale to which a recommended torque can be set; it is not a screwdriver. It is typically made of metallic materials and may have a standard socket connection to accommodate different tool bits/adaptors.

**Wrench Adaptors**

A component of a surgical wrench (typically ratchet type) intended to enable connection of the main body of the wrench to a variety of distal end pieces (bits). It is intended to be used during an invasive surgical procedure, and is not dedicated to dental use only. It may be referred to as a converter

**Intended User**

Maxillo-facial Surgeons, General Dentists, Orthodontists, Periodontist, Prosthodontists, and other appropriately trained and experienced implant users.

**Intended Environment**

The devices are intended to be used in a clinical environment such as an operating theater or a dentist consultation room.

**Intended Patient Population**

Patients that have lost one tooth or multiple teeth.

**Indications for Use**

Southern Implants (Pty) Ltd's manufactures and sells various reusable instruments that are intended to facilitate the insertion of restoration of its dental implants by clinicians.

**Contraindications**

- Patients under the age of 18
- Poor bone quality
- Blood disorders
- Infected implant site
- Vascular impairment
- Uncontrolled diabetes
- Drug or alcohol abuse
- Chronic high dose steroid therapy
- Anti-coagulant therapy
- Metabolic bone disease
- Radiotherapy treatment
- Sinus pathology

**Warnings**

THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING.

- For the safe and effective use of dental implants it is suggested that specialised training be undertaken, including hands-on training to learn proper technique, biomechanical requirements and radiographic evaluations.
- Responsibility for proper patient selection, adequate training, experience in the placement of implants, and providing appropriate information for informed consent rests with the practitioner. Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone.

**Cautions**

New and experienced Implant users should do training before using a new system or attempt to do a new treatment method. Take special care when treating patients who have local or systemic factors that could affect the healing of the bone and soft tissue (i.e. poor oral hygiene, uncontrolled diabetes, are on steroid therapy, smokers, infection in the nearby bone and patients who had oro-facial radiotherapy). Thorough screening of prospective implant candidates must be performed including:

- a comprehensive medical and dental history.
- visual and radiological inspection to determine adequate bone dimensions, anatomical landmarks, occlusal conditions and periodontal health.
- bruxism and unfavourable jaw relations must be taken into account.
- proper pre-operative planning with a good team approach between well trained surgeons, restorative dentists and lab technicians is essential for successful implant treatment.
- minimizing the trauma to the host tissue increases the potential for successful osseointegration.
- electro-surgery should not be attempted around metal implants, as they are conductive.

**Disposal**

Disposal of the device and its packaging; Follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

**Disclaimer of liability**

This product is part of the Southern Implants product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants product range and take full responsibility for the correct indications and use of this product. Southern Implants does not assume liability for damage due to incorrect use. Please note that some Southern Implants products may not be cleared or released for sale in all markets.

**Notice regarding serious incidents**

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and / or patient is established. The contact information for the manufacturer of this device to report a serious incident is as follows: [sicomplaints@southernimplants.com](mailto:sicomplaints@southernimplants.com)

**Basic UDI**

Product	Basic-UDI Number
Basic-UDI for Reusable Instruments	600954403876

**Related literature & catalogues**

- CAT-2004 - Tri-Nex® Implants Product Catalogue
- CAT-2005 - IT Implants Product Catalogue
- CAT-2010 - Osseointegrated Fixtures Catalogue
- CAT-2020 - External Hex Implants Product Catalogue
- CAT-2042 - Deep Conical Implants Product Catalogue
- CAT-2043 - Internal Hex Implants Product Catalogue
- CAT-2060 - PROVATA® Implants Product Catalogue
- CAT-2069 - INVERTA® Implants Product Catalogue
- CAT-2070 - Zygomatic Implants Product Catalogue

ARCHIVED

**Symbols and Warnings**

 Manufacturer: Southern Implants 1 Albert Rd, P.O Box 605 IRENE, 0062, South Africa. Tel: +27 12 667 1046	 2797	 Prescription device*	 Sterilization using Irradiation	 Non-sterile	 Caution	 Consult instruction for use	 Use by date (mm-yy)	 Do not reuse	 Do not re-sterilize	 Batch code	 Do not use if package is damaged	 Medical Device
* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.						Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.						
All rights reserved. Southern Implants®, the Southern Implants logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale.												

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