

## SIDIGITAL Library

### Exocad Installation instructions

## Download & Install the Southern Implants Exocad library in 3 Simple Steps

### STEP 1

LOG INTO  
SIDIGITAL PAGE

### STEP 2

DELETE OLD  
LIBRARIES  
(FOR EXISTING USERS ONLY)

### STEP 3

INSTALL NEW  
SIDIGITAL LIBRARY

## STEP 1


### STEP 1.1

#### Navigate to the CAD/CAM - SIDIGITAL page

Click [here](#) to go to the library registration page

Or

Visit [southernimplants.com](http://southernimplants.com) and navigate to the webpage by hovering over the CAD/CAM - SIDIGITAL tab on the menu bar and click on "CAD/CAM - SIDigital Library Registration & Access".



The screenshot shows the Southern Implants website. The top navigation bar includes links for Home, Treatment Solutions, Innovative Products, CAD/CAM - SIDigital, News & Events, Resource Library, and Contact Us. The CAD/CAM - SIDigital menu is open, showing options like SIDigital Catalog and Manual, CAD/CAM - SIDigital Library Registration & Access (highlighted), Scan Flag/Scan Body, CIA Scanning Abutments, Passive Abutment, TIB Abutments, Digital Laboratory Analogues, and SIREAL Guided Surgery Tool. The main content area features a 'Welcome to Southern Implants SIDigital' message, a 3shape logo, and an exocad logo. Below these are buttons for 'Installation Instructions' and 'Ins'. At the bottom, there are sections for 'New user:' with a 'Click here to register' button, and 'Existing and newly registered user log-in:' with a 'Click here to download' button. The footer contains contact information and social media links.


**STEP 1.2**  
**Create new account / log-in**

**FOR NEW USERS**

Create an account by clicking  
"Click here to register".

**FOR EXISTING USERS**

Navigate to the log in screen by clicking  
"Click here to download".



Complete and submit the form. A confirmation email with your login details will follow, check junk folder if not received within 48 hours. Navigate to the login page by clicking the link in the email or by navigating back to  
<https://southernimplants.com/cad-cam-download/>

From: sidigital@southernimplants.com  
Sent:  
To:  
Subject: SI Digital Library Login Details

Good day,

Thank you for registering to SIDigital. To download, click on this link [Access SIDigital](#)

Username:  
Password:

Kind regards,  
SI Digital Team

Enter your password and email address  
(remember to complete the CAPCHA verification)

Log in to download CAD/CAM libraries:

This content is password protected. To view it please enter your password below:

>> [Forgot password?](#)

Email:


Password:

☐ I'm not a robot



## STEP 1.3 Download Exocad Library


Click on the Exocad 'Download Latest Version' to initiate the download.

English

Home Treatment Solutions Innovative Products CAD/CAM - SIDigital News & Events Resource Library Contact Us


### Download SIDigital Libraries

Please choose your software: For any queries, contact: [sidigital@southernimplants.com](mailto:sidigital@southernimplants.com)




3 Shape  
USA Encrypted

Download latest version




3 Shape

Download latest version



Dental Wings

Download latest version



Exocad

Download latest version

3 Shape USA Encrypted Change Log

3 Shape Change Log

Dental Wings Change Log

Exocad Change Log

NOTE: A change-log is available of all revision changes to the library at the bottom of the screen.

## STEP 2


### DELETING YOUR PREVIOUS LIBRARY

(FOR EXISTING USERS ONLY)

#### STEP 2.1

##### Delete the Exocad library saved on your local computer


Navigate to your local hardrive (C:\) and delete the existing Southern Implants 'Exocad' folder.



For instances where the previous library was installed in a different location, it is essential that the user navigates to the correct directory and delete the library.

## Delete the Exocad abutment library saved on your local computer


Navigate to the Exocad abutment library folder (found in the Exocad installation folder at file directory DentalCADApp/library/implant) and delete all the Southern Implants abutment kits folders.



## STEP 2.3

## Delete the modelcreator kits saved on your local computer

Navigate to the Exocad abutment kit library folder (found in the Exocad installation folder at file directory DentalCADApp/library/modelcreator/implants) and delete all the Southern Implants modelcreator kits folders.




**The library is now deleted.**

## STEP 3

### STEP 3.1

#### Navigate to and unzip the Southern Implants library file


Navigate to your downloads folder on your computer and unzip the library file to a folder directory of your choice.



### STEP 3.2

#### Open the Exocad Implant folder



Open the unzipped folder and open the 'Implant' folder.



### STEP 3.3

#### Copy all the implant library folders to the Exocad Folder

Highlight all of the folders in the 'Implant' folder and copy them to the Exocad implant library folder (found in the Exocad installation folder at file directory DentalCADApp/library/implant).






### STEP 3.4

#### Open the Exocad's 'modelcreator' folder

Navigate back to the unzipped SI Exocad folder and open the 'modelcreator' folder.







### STEP 3.5

#### Copy all the model creator library folders to the Exocad Folder

Highlight all of the folders in the 'modelcreator' folder and copy them to the Exocad modelcreator implant library folder (found in the Exocad installation folder at file directory DentalCADApp/library/modelcreator/implants).




The Exocad library is now successfully installed.

## ADDITIONAL INFORMATION

### SIDIGITAL naming Scheme

Main software selection



Secondary software selection

