|           |                                                                                              |
|-----------|----------------------------------------------------------------------------------------------|
| English   | INSTRUCTIONS FOR USE: Southern Implants® Tapered Drills and Taps (reusable)                  |
| Español   | INSTRUCCIONES DE USO: Fresas cónicas y machos de roscar Southern Implants® (reutilizables)   |
| Italiano  | ISTRUZIONI PER L'USO: Southern Implants® Punte e Maschi Conici (riutilizzabili)              |
| Français  | MODE D'EMPLOI : Foret coniques et tarauds (réutilisables) Southern Implants®                 |
| Deutsch   | GEBRAUCHSANWEISUNGEN: Southern Implants® Kegelbohrer und Gewindeschneider (wiederverwendbar) |
| Português | INSTRUÇÕES DE UTILIZAÇÃO: Southern Implants® Brocas e machos cônicos (reutilizáveis)         |
| Ελληνικά  | ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ: Southern Implants® Κωνικά τρυπάνια και τάπες (επαναχρησιμοποιήσιμα)          |
| Svenska   | ANVISNINGAR FÖR ANVÄNDNING: Southern Implants® Koniska Borrar och Kranar (återanvändbara)    |

ARCHIVED

**South Africa - Headquarters:** 1 Albert Road, Irene, 0062, RSA  
T: +27-12-667-1046 | E: [info@southernimplants.com](mailto:info@southernimplants.com)

**Southern Implants Europe AB:** Holmgatan 30, S-791 71 Falun, Sweden  
T: +46 23 13300 | E: [ecrep@southernimplants.com](mailto:ecrep@southernimplants.com)

#### Subsidiaries

##### Australia

Southern Implants Australia  
T: +61-(0)-8-9466-2627  
E: [info@southernimplants.com.au](mailto:info@southernimplants.com.au)

##### Spain and Portugal

Southern Implants Iberica  
T: +34 935 053 507  
E: [info@southernimplants.es](mailto:info@southernimplants.es)

##### United Kingdom and Ireland

Southern Implants UK  
T: +44-20-8899-6845 / 6 / 7  
E: [info@southernimplants.co.uk](mailto:info@southernimplants.co.uk)

##### USA and Canada

Southern Implants North America Inc.  
T: +1-561-472-0990  
E: [customer care@southernimplants.com](mailto:customer care@southernimplants.com)

## Description

Southern Implants® tapered drills are intended to enlarge a channel made by a twist drill in order to prepare a conical osteotomy site for a tapered endosseous dental implant. This entails drilling a hole in the mandible or maxilla or making an existing hole larger in the conical dimension. The drills are for multiple use and supplied sterile. They are made from Stainless Steel and Titanium grade 4 or 5. The drills are available in diameters between Ø2.7 mm and Ø9 mm.

Southern Implants® taps are intended to be used to create a thread in the bone identical to that of the corresponding implant prior to placing the implant. The taps are for multiple use and supplied sterile. They are made from Titanium grade 4 or 5. The taps are available in diameters between Ø3 mm and Ø9 mm.

## Intended use

Tapered drills and taps are intended to aid in the treatment of partially or fully edentulous patients eligible for placement of one or more dental implants as a means of fixing a permanent or removable single crown, partial or full-arch dental prosthesis in the upper or lower jaw. The system constituents are classified as medical devices and are multiple use devices.

## Indications for use

The indications for use for the device are not dissimilar to that of dental implants. Southern Implants dental implants are intended for both one and two-stage surgical procedures in the following situations and with the following clinical protocols:

- replacing single and multiple missing teeth in the mandible and maxilla.
- immediate placement in extraction sites and in situations with a partially or completely healed alveolar ridge.
- immediate loading in all indications, except in single tooth situations on implants shorter than 8 mm or in soft bone (type IV) where implant stability may be difficult to obtain, and immediate loading may not be appropriate.

## Intended user

Maxillo-facial Surgeons, General Dentists, Orthodontists, Periodontists, Prosthodontists and other appropriately trained and experienced implant users.

## Intended environment

The devices are intended to be used in a clinical environment such as an operating theatre or a dentist consultation room.

## Intended patient population

The intended patient population for the device is not dissimilar to that of dental implant therapy. The intended patient population for implant therapy is partially or fully edentulous patients requiring prosthetic dental restoration in the upper or lower jaw. Restorations may comprise single teeth, partial or full bridges and are fixed restorations. Additionally, the intended population for some device variants includes those who have previously undergone dental implant therapy.

## Compatibility information

Southern Implants' implants tapered drills and taps are as described in Table A and B respectively. Details on site preparation sequences recommended by Southern Implants® can be found in the relevant product catalogues.

**Table A – Tapered Drills**

| Drill Code                         | Material                     | Coating (if any) | Number of Uses |
|------------------------------------|------------------------------|------------------|----------------|
| <b>EXTERNAL HEX &amp; PROVATA®</b> |                              |                  |                |
| D-30TP-XX                          | Titanium grade 5 (ASTM F136) | -                | Up to 10       |
| D-33TP-XX                          | Titanium grade 5 (ASTM F136) | -                | Up to 10       |
| D-40TP-XX                          | Stainless Steel (DIN 1.4197) | -                | Up to 10       |

|                                            |                              |                 |          |
|--------------------------------------------|------------------------------|-----------------|----------|
| D-50TP-XX                                  | Titanium grade 5 (ASTM F136) | -               | Up to 10 |
| D-60TP-XX                                  | Titanium grade 5 (ASTM F136) | -               | Up to 10 |
| <b>Guided</b>                              |                              |                 |          |
| D-33TP-XX-L                                | Titanium grade 5 (ASTM F136) | -               | Up to 10 |
| D-40TP-XX-L                                | Titanium grade 5 (ASTM F136) | -               | Up to 10 |
| D-50TP-XX-L                                | Titanium Alloy               | -               | Up to 10 |
| <b>Dense Bone Drills</b>                   |                              |                 |          |
| D-42TP-XX                                  | Stainless Steel (DIN 1.4197) | TiN             | Up to 10 |
| D-52TP-XX                                  | Stainless Steel (DIN 1.4197) | TiN             | Up to 10 |
| D-62TP-XX                                  | Stainless Steel (DIN 1.4197) | TiN             | Up to 10 |
| <b>TRI-NEX®</b>                            |                              |                 |          |
| D-L-35-XX                                  | Stainless Steel (DIN 1.4197) | -               | Up to 10 |
| D-L-43-XX                                  | Stainless Steel (DIN 1.4197) | -               | Up to 10 |
| D-L-50-XX                                  | Stainless Steel (DIN 1.4197) | -               | Up to 10 |
| D-L-60-XX                                  | Stainless Steel (DIN 1.4197) | -               | Up to 10 |
| <b>Soft Bone Drills</b>                    |                              |                 |          |
| DLS-35-XX                                  | Titanium grade 5 (ASTM F136) | -               | Up to 10 |
| DLS-43-XX                                  | Titanium grade 5 (ASTM F136) | -               | Up to 10 |
| DLS-50-XX                                  | Titanium grade 5 (ASTM F136) | -               | Up to 10 |
| DLS-60-XX                                  | Titanium grade 5 (ASTM F136) | -               | Up to 10 |
| <b>DEEP CONICAL (DC)</b>                   |                              |                 |          |
| D-DCT-30XX                                 | Titanium grade 5 (ASTM F136) | ALTiN           | Up to 10 |
| D-DCT-35XX                                 | Titanium grade 5 (ASTM F136) | ALTiN           | Up to 10 |
| D-DCT-40XX                                 | Titanium grade 5 (ASTM F136) | ALTiN           | Up to 10 |
| D-DCT-50XX                                 | Titanium grade 5 (ASTM F136) | ALTiN           | Up to 10 |
| <b>M-SERIES/INTERNAL HEX</b>               |                              |                 |          |
| D-MT37XX                                   | Titanium grade 5 (ASTM F136) | -               | Up to 10 |
| D-MT42XX                                   | Titanium grade 5 (ASTM F136) | -               | Up to 10 |
| D-MT50XX                                   | Titanium grade 5 (ASTM F136) | -               | Up to 10 |
| <b>INTERNAL OCTAGON (IT)</b>               |                              |                 |          |
| D-4XXT                                     | Stainless Steel (DIN 1.4197) | -               | Up to 10 |
| D-5XXT                                     | Stainless Steel (DIN 1.4197) | -               | Up to 10 |
| D-6XXT                                     | Stainless Steel (DIN 1.4197) | -               | Up to 10 |
| <b>INVERTA®</b>                            |                              |                 |          |
| D-IV2700                                   | Titanium grade 4 (ASTM F67)  | Anodized Silver | Up to 10 |
| D-IV37XX                                   | Titanium grade 4 (ASTM F67)  | Anodized Pink   | Up to 10 |
| D-IV45XX                                   | Titanium grade 5 (ASTM F136) | Anodized Blue   | Up to 10 |
| D-IV50XX                                   | Titanium grade 5 (ASTM F136) | Anodized Yellow | Up to 10 |
| D-IV60XX                                   | Titanium grade 5 (ASTM F136) | Anodized Silver | Up to 10 |
| <b>Guided</b>                              |                              |                 |          |
| D-IV37XXGS                                 | Titanium grade 4 (ASTM F67)  | Anodized Pink   | Up to 10 |
| D-IV45XXGS                                 | Titanium grade 5 (ASTM F136) | Anodized Blue   | Up to 10 |
| D-IV50XXGS                                 | Titanium grade 5 (ASTM F136) | Anodized Yellow | Up to 10 |
| D-IV60XXGS                                 | Titanium grade 5 (ASTM F136) | Anodized Silver | Up to 10 |
| <b>MAX (MAX, PROMAX, TRI-MAX®, MAXIT®)</b> |                              |                 |          |
| D-MAX6-XX                                  | Titanium grade 5 (ASTM F136) | -               | Up to 10 |
| D-70TP-XX                                  | Titanium grade 5 (ASTM F136) | -               | Up to 10 |
| D-70TP-XX-L                                | Titanium grade 5 (ASTM F136) | -               | Up to 10 |
| D-80TP-XX                                  | Titanium grade 5 (ASTM F136) | -               | Up to 10 |
| D-80TP-XX-L                                | Titanium grade 5 (ASTM F136) | -               | Up to 10 |
| D-90TP-XX                                  | Titanium grade 5 (ASTM F136) | -               | Up to 10 |
| D-90TP-XX-L                                | Titanium grade 5 (ASTM F136) | -               | Up to 10 |

| IE IMPLANTS           |                              |   |          |
|-----------------------|------------------------------|---|----------|
| D-40E-04F             | Titanium grade 5 (ASTM F136) | - | Up to 10 |
| SINGLE PLATFORM (SP1) |                              |   |          |
| D-SP30/L              | Titanium grade 5 (ASTM F136) | - | Up to 10 |
| D-SP35/L              | Titanium grade 5 (ASTM F136) | - | Up to 10 |
| D-SP40/L              | Titanium grade 5 (ASTM F136) | - | Up to 10 |
| D-SP45                | Titanium grade 5 (ASTM F136) | - | Up to 10 |
| D-SP50                | Titanium grade 5 (ASTM F136) | - | Up to 10 |
| Guided                |                              |   |          |
| D-SP30XX-GSN          | Titanium grade 5 (ASTM F136) | - | Up to 10 |
| D-SP30XX-GS           | Titanium grade 5 (ASTM F136) | - | Up to 10 |
| D-SP35XX-GSN          | Titanium grade 5 (ASTM F136) | - | Up to 10 |
| D-SP35XX-GS           | Titanium grade 5 (ASTM F136) | - | Up to 10 |
| D-SP40XX-GS           | Titanium grade 5 (ASTM F136) | - | Up to 10 |
| D-SP45XX-GS           | Titanium grade 5 (ASTM F136) | - | Up to 10 |
| D-SP50XX-GS           | Titanium grade 5 (ASTM F136) | - | Up to 10 |

Table B – Taps

| Drill Code              | Material                     | Coating (if any) | Number of Uses |
|-------------------------|------------------------------|------------------|----------------|
| EXTERNAL HEX & PROVATA® |                              |                  |                |
| D-TAP-IBN               | Titanium grade 5 (ASTM F136) | -                | Up to 10       |
| D-TAP-IBNT              | Titanium grade 5 (ASTM F136) | -                | Up to 10       |
| D-TAP-IBS               | Titanium grade 5 (ASTM F136) | -                | Up to 10       |
| D-TAP-IBT               | Titanium grade 5 (ASTM F136) | -                | Up to 10       |
| D-TAP-I4B               | Titanium grade 5 (ASTM F136) | -                | Up to 10       |
| D-TAP-BA                | Titanium grade 5 (ASTM F136) | -                | Up to 10       |
| D-TAP-BAT               | Titanium grade 5 (ASTM F136) | -                | Up to 10       |
| D-TAP-BBBS              | Titanium grade 5 (ASTM F136) | -                | Up to 10       |
| D-TAP-BBBT              | Titanium grade 5 (ASTM F136) | -                | Up to 10       |
| TRI-NEX®                |                              |                  |                |
| D-TAP-L-35              | Titanium grade 5 (ASTM F136) | -                | Up to 10       |
| D-TAP-L-43              | Titanium grade 5 (ASTM F136) | -                | Up to 10       |
| D-TAP-L-50              | Titanium grade 5 (ASTM F136) | -                | Up to 10       |
| D-TAP-L-60              | Titanium grade 5 (ASTM F136) | -                | Up to 10       |
| D-TAP-LS-35             | Titanium grade 5 (ASTM F136) | -                | Up to 10       |
| D-TAP-LS-43             | Titanium grade 5 (ASTM F136) | -                | Up to 10       |
| D-TAP-LS-50             | Titanium grade 5 (ASTM F136) | -                | Up to 10       |
| DEEP CONICAL (DC)       |                              |                  |                |
| D-TAP-DCC30             | Titanium grade 5 (ASTM F136) | -                | Up to 10       |
| D-TAP-DCC35             | Titanium grade 5 (ASTM F136) | -                | Up to 10       |
| D-TAP-DCC40             | Titanium grade 5 (ASTM F136) | -                | Up to 10       |
| D-TAP-DCC50             | Titanium grade 5 (ASTM F136) | -                | Up to 10       |
| INTERNAL OCTAGON (IT)   |                              |                  |                |
| D-TAP-ITT4              | Titanium grade 5 (ASTM F136) | -                | Up to 10       |
| D-TAP-ITT5              | Titanium grade 5 (ASTM F136) | -                | Up to 10       |
| D-TAP-ITT6              | Titanium grade 5 (ASTM F136) | -                | Up to 10       |
| D-TAP-ITC3              | Titanium grade 5 (ASTM F136) | -                | Up to 10       |
| D-TAP-ITC4              | Titanium grade 5 (ASTM F136) | -                | Up to 10       |
| D-TAP-ITC5              | Titanium grade 5 (ASTM F136) | -                | Up to 10       |
| INVERTA®                |                              |                  |                |
| IV-TR45XX               | Titanium Grade 4 (ASTM F67)  | Anodized Blue    | Up to 10       |
| IV-TR50XX               | Titanium Grade 4 (ASTM F67)  | Anodized Yellow  | Up to 10       |
| IV-TR60XX               | Titanium Grade 4 (ASTM F67)  | Anodized Silver  | Up to 10       |

| MAX (MAX, PROMAX, TRI-MAX®, MAXIT®) |                              |   |          |
|-------------------------------------|------------------------------|---|----------|
| D-TAP-MAX6-XX                       | Titanium grade 5 (ASTM F136) | - | Up to 10 |
| D-TAP-MAX7-XX                       | Titanium grade 5 (ASTM F136) | - | Up to 10 |
| D-TAP-MAX8-XX                       | Titanium grade 5 (ASTM F136) | - | Up to 10 |
| D-TAP-MAX9-XX                       | Titanium grade 5 (ASTM F136) | - | Up to 10 |
| D-TAP-MAX10-XX                      | Titanium grade 5 (ASTM F136) | - | Up to 10 |
| FIRST® ROHNER TECHNIQUE             |                              |   |          |
| D-TAP-35RT                          | Titanium grade 5 (ASTM F136) | - | Up to 10 |
| D-TAP-43RT                          | Titanium grade 5 (ASTM F136) | - | Up to 10 |

### Clinical benefits

The clinical benefits are not dissimilar to that of dental implant therapy. The clinical benefits of dental implant therapy include improved chewing function, speech, aesthetics and patient psychological wellbeing. Through this procedure patients can expect to have their missing teeth replaced and/or crowns restored.

The clinical benefits specific to drills include:

- Effective cutting of bone within acceptable time limits.
- Precise preparation of the implant site to the dimensions of the intended implant.
- Can remain effective after multiple uses to the extent stipulated by the manufacturer.

### Storage, cleaning and sterilisation

The component is supplied sterile (sterilised by gamma irradiation). Sterility is assured unless the container or seal is damaged or opened. If packaging is damaged do not use the product and contact your Southern representative or return to Southern Implants®. The devices must be stored in a dry place at room temperature and not exposed to direct sunlight. Incorrect storage may influence device characteristics.

### Limitations on reusable items

A direct value for reusable instruments cannot be given. Frequent processing may have minor effects on the instruments. The product life is normally determined by wear and damage during use, thus instruments if properly cared for and inspected after each use, can be reused many times. Maintain a checklist for these instruments recording the number of uses.

Prior to reprocessing the device, it should be thoroughly inspected and tested to determine its suitability for reuse.

**NOTE:** during use, handle drills and instruments with sterile tweezers to minimise contamination of the instrument tray and risk of damage to sterile surgical gloves.

### Containment

As soon as practically possible, remove all visible residue after use (bone, blood or tissue), by immersing the instrument in cold water (dried soil may be difficult to remove).

### Precleaning

Disassemble instruments from hand pieces and all connecting parts from instruments in order to clean soil from obstructed areas. Rinse with lukewarm water for 3 minutes and remove hardened debris with a soft nylon brush. Avoid mechanical damage during cleaning.

### Manual cleaning or automated cleaning

Prepare an ultrasonic bath with suitable detergent (i.e. Steritech instrument cleaner - 5% dilution), sonicate for 20 minutes (alternative methods can be used if proven by the end user). Rinse with purified/sterile water. Load devices into a thermo-disinfector. Run the cleaning and disinfection cycle, followed by the drying cycle.

**NOTE:** always follow the instructions for use of the manufacturers for cleaning agents and disinfectants.

**Drying**

Dry the drills with filtered compressed air or single use, lint free wipes. Pack the instruments as quickly as possible after removal into the storage container. If additional drying is necessary, dry in a clean location. Moisture on drills can cause corrosion and deterioration of the cutting edges.

**Inspection**

Do a visual inspection of the items to check for any damages.

**Packaging**

Use the correct packaging material as indicated for steam sterilisation to ensure sterility is maintained. Double packaging is recommended. Where appropriate, the cleaned, disinfected and checked devices can be assembled and placed in instrument trays as applicable. Instrument trays can be double wrapped or placed in sterilisation bags.

**Sterilisation**

Southern Implants® recommends one of the following procedures to sterilise the tapered drills and taps prior to use:

1. prevacuum sterilisation method: steam sterilise the abutments at 132°C (270°F) at 180 - 220 kPa for 4 minutes. Dry for at least 20 minutes in the chamber. Only an approved wrap or pouch for steam sterilisation must be used.
2. for users in the USA: prevacuum sterilisation method: wrapped, steam sterilise at 135°C (275°F) at 180 - 220 kPa for 3 minutes. Dry for 20 minutes in the chamber. Use a wrap or pouch that is cleared for the indicated steam sterilisation cycle.

**NOTE:** users in the USA must ensure that the steriliser, wrap or pouch, and all steriliser accessories are cleared by the FDA, for the intended sterilisation cycle.

**Storage**

Maintain packaging integrity to ensure sterility in storage. Packaging should be completely dry before storage to avoid corrosion and degradation of cutting edges.

**Contraindications**

The contraindications for this system are not dissimilar to those of dental implant therapy. Contra-indications to implant therapy include:

- patients medically unfit for oral surgical procedures.
- where an inadequate number of implants can be placed limiting the functional support of the prosthesis.
- patients under the age of 18 years.
- patients with poor bone quality.
- patients with blood disorders.
- presence of infection at the implant site.
- patients with vascular impairment.
- patients with uncontrolled diabetes.
- patients with drug or alcohol abuse dependencies.
- patients undergoing chronic high dose steroid therapy.
- patients undergoing anti-coagulant therapy.
- patients with metabolic bone disease.
- patients undergoing radiotherapy treatment.
- Patients with pure titanium, titanium alloy (Ti6Al4V), gold, palladium, or iridium allergies or hypersensitivities.

**Warnings and precautions**

THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING.

- For the safe and effective use of dental implants it is strongly suggested that specialised training be undertaken, including hands-on training to learn proper technique, biomechanical requirements and radiographic evaluations.

- Products must be secured against aspiration when handled intraorally. Aspiration of products may lead to infection or unplanned physical injury.

Responsibility for proper patient selection, adequate training, experience in the placement of implants and providing appropriate information for informed consent rests with the practitioner. Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone. Implant failure increases when implants are placed in irradiated bone as radiotherapy can result in progressive fibrosis of vessels and soft tissue, leading to diminished healing capacity.

It is important to be aware and avoid damage to vital structures such as nerves, veins and arteries. Injuries to vital anatomical structures may cause serious complications such as injury to the eye, nerve damage and excessive bleeding. It is essential to protect the infraorbital nerve. Failing to identify actual measurements relative to the radiographic data could lead to complications.

New and experienced implant users should do training before using a new system or attempting to do a new treatment method. Take special care when treating patients who have local or systemic factors that could affect the healing of the bone and soft tissue (i.e., poor oral hygiene, uncontrolled diabetes, are on steroid therapy, smokers, infection in the nearby bone and patients who had orofacial radiotherapy).

Thorough screening of prospective implant candidates must be performed including:

- a comprehensive medical and dental history.
- visual and radiological inspection to determine adequate bone dimensions, anatomical landmarks, occlusal conditions and periodontal health.
- bruxism and unfavourable jaw relations must be taken into account.
- proper preoperative planning with a good team approach between well trained surgeons, restorative dentists and lab technicians is essential for successful implant treatment.
- minimising the trauma to the host tissue increases the potential for successful osseointegration.
- electrosurgery should not be attempted around metal implants as they are conductive.

Should the device not operate as intended, it must be reported to the manufacturer of the device. The contact information for the manufacturer of this device to report a change in performance is: [sicomplaints@southernimplants.com](mailto:sicomplaints@southernimplants.com).

### Side effects

The side effects of the use of the system are not dissimilar to those of dental implant therapy. Possible side effects to implant therapy include:

- pain.
- swelling.
- phonetic difficulties.
- gingival inflammation.

Less common but more persistent symptoms include, but are not limited to:

- allergic reaction(s) to implant and/or abutment material.
- breakage of the implant and/or abutment.
- loosening of the abutment screw and/or retaining screw.
- infection requiring revision of the dental implant.
- nerve damage resulting in permanent weakness, numbness, or pain.
- histologic responses with possible macrophage and/or fibroblast involvement.
- fat emboli formation.
- loosening of the implant requiring revision surgery.
- perforation of the maxillary sinus.
- perforation of the labial and lingual plates.
- bone loss possibly resulting in revision or removal of the implant.

**Precaution: maintaining sterility protocol**

Implants are packaged as follows:

1. An outer package consisting of a rigid, clear box which acts as protection for the inner package.
2. The inner package consisting of a blister pack (clear plastic-formed blister base with a TYVEK “peel-back” lid).
3. Within the inner package, there is a hollow tube which contains one implant suspended from a titanium ring, this ensures the implant never touches the inside of the plastic tube.
4. Labelling information is located on the surface of the peel-back lid and on the outside of the rigid box.

Care must be taken to maintain the sterility of the implant by proper opening of the packaging and handling of the implant.

1. Open the implant package in the non-sterile field, with non-sterile gloves, tear the address label to open the box.
2. With non-sterile gloves, remove the inner blister pack. Do not place the plastic box or blister pack-lid onto the sterile field. The contents of this inner package are sterile.
3. The sealed blister is to be opened by an assistant (with nonsterile gloves): remove the TYVEK lid and drop or place the sterile tube onto the sterile field, open the tube cap and attach the implant placement tool onto the implant and carefully remove from the sterile tube. Do not touch the sterile implant.

Other sterile components are packed in a peel pouch or blister base with a “peel-back” lid. Labelling information is located on the bottom half of the pouch, inside the packet or on the surface of the peel-back lid. Sterility is assured unless the pouch is damaged or opened. Non-sterile components are supplied clean but not sterile in a peel pouch or blister base with peel-back lid. Labelling information is located on the bottom half of the pouch or on the surface of the peel-back lid.

**Notice regarding serious incidents**

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established.

The contact information for the manufacturer of this device to report a serious incident is as follows: [sicomplaints@southernimplants.com](mailto:sicomplaints@southernimplants.com).

**Materials**

Tapered Drills: Titanium Grade 4 (ASTM F67), Titanium grade 5 (ASTM F136) or Stainless Steel (DIN 1.4197)  
Taps: Titanium grade 5 (ASTM F136) or Titanium Grade 4 (ASTM F67)

**Disposal**

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

**Summary of Safety and Clinical Performance (SSCP)**

As required by the European Medical Device Regulation (MDR; EU2017/745), a Summary of Safety and Clinical Performance (SSCP) is available for perusal with regard to Southern Implants® product ranges.

The relevant SSCP can be accessed at <https://ec.europa.eu/tools/eudamed>.

**NOTE:** the above website will be available upon the launch of the European Database on Medical Devices (EUDAMED).

**Disclaimer of liability**

This product is part of the Southern Implants® product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has

to study the development of the Southern Implants® product range and take full responsibility for the correct indications and use of this product. Southern Implants® does not assume liability for damage due to incorrect use. Please note that some Southern Implants® products may not be cleared or released for sale in all markets.


### Basic UDI

| Product                                    | Basic-UDI Number |
|--------------------------------------------|------------------|
| Basic-UDI for Drills and Handpiece Devices | 6009544038759C   |
| Basic-UDI for Reusable Instruments         | 6009544038769E   |

### Related literature and catalogues

CAT-2004 - Tri-Nex® Implants Product Catalogue  
 CAT-2020 - External Hex Implants Product Catalogue  
 CAT-2042 - Deep Conical (DC) Implants Product Catalogue  
 CAT-2043 - Internal Hex (M-Series) Implants Product Catalogue  
 CAT-2060 - Internal Hex (PROVATA®) Implants Product Catalogue  
 CAT-2005 - IT (Internal Octagon) Implants Product Catalogue  
 CAT-2069 - INVERTA® Implants Product Catalogue  
 CAT-2068 - SIGuided Implants Product Catalogue  
 CAT-2010 - Osseointegrated Fixtures Product Catalogue  
 CAT-2091M - FIRST® Fibula Reconstruction System Technology  
 CAT-2093 - Single Platform Implants Product Catalogue

### Symbols and warnings

\* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

All rights reserved. Southern Implants®, the Southern Implants® logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants®. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale. It is the responsibility of the clinician to inspect the symbols that appear on the packaging of the product in use.

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