
English	INSTRUCTIONS FOR USE: Southern Implants® Tapered Drills and Taps
Español	INSTRUCCIONES DE USO: Fresas cónicas y machos Southern Implants®
Italiano	ISTRUZIONI PER L'USO: Southern Implants® Punte e Maschi Conici
Français	MODE D'EMPLOI : Forets et Tarauds coniques Southern Implants®
Deutsch	GEBRAUCHSANWEISUNGEN: Southern Implants® Kegelbohrer und Gewindeschneider
Português	INSTRUÇÕES DE UTILIZAÇÃO: Brocas e machos cónicos Southern Implants®
Ελληνικά	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ: Southern Implants® Κωνικά τρυπάνια και τάπες
Svenska	ANVISNINGAR FÖR ANVÄNDNING: Southern Implants® Koniska Borrar och Kranar
Română	INSTRUCȚIUNI DE UTILIZARE: Southern Implants® Burghie conice și tarozi

South Africa - Headquarters: 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

Subsidiaries

Australia

Southern Implants Australia
T: +61-(0)-8-9466-2627
E: info@southernimplants.com.au

Spain and Portugal

Southern Implants Iberica
T: +34 935 053 507
E: info@southernimplants.es

United Kingdom and Ireland

Southern Implants UK
T: +44-20-8899-6845 / 6 / 7
E: info@southernimplants.co.uk

USA and Canada

Southern Implants North America Inc.
T: +1-561-472-0990
E: customercare@southernimplants.com

Description

Southern Implants® Tapered Drills are intended to enlarge a channel made by a twist drill in order to prepare a conical osteotomy site for a tapered endosseous dental implant. This entails drilling a hole in the mandible or maxilla or making an existing hole larger in the conical dimension. The drills are available as multiple use and single-use versions, supplied either sterile or non-sterile. They are made from Stainless Steel and Titanium grade 4 or 5. The drills are available in diameters between Ø2.7 mm and Ø9 mm.

Southern Implants® Taps are intended to be used to create a thread in the bone identical to that of the corresponding implant prior to placing the implant. The taps are for multiple use and supplied sterile. They are made from Titanium grade 4 or 5. The taps are available in diameters between Ø3 mm and Ø9 mm.

Intended use

The Tapered Drills and Taps are intended to be used to prepare the osteotomy for implant placement.

Indications for use

Dental Implant therapy is indicated for patients with one or more missing teeth that wish to have the tooth/teeth restored with a dental implant and prosthetic tooth/teeth.

Intended user

The intended user for this system includes Maxillofacial Surgeons, General Dentists, Orthodontists, Periodontist, Prosthodontists and other appropriately trained and experienced implant users.

Intended environment

This device is intended to be used in a clinical environment such as an operating theatre or a dentist consultation room.

Intended patient population

The intended patient population for dental implant therapy includes partially or fully edentulous patients eligible or, otherwise, not contraindicated for implant placement that require prosthetic restoration or revision of existing restorations in the upper or lower jaw; where the planned restorations are to be fixed and comprise of single teeth, partial or full bridges.

Compatibility information

Southern Implants' implants tapered drills and taps are as described in Table A and B respectively. Details on site preparation sequences recommended by Southern Implants® can be found in the relevant product catalogues.

Table A – Tapered Drills

Drill Code	Material	Coating (if any)	Number of Uses	How Provided
EXTERNAL HEX & PROVATA®				
D-30TP-XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-33TP-XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-40TP-XX	Stainless Steel (DIN 1.4197)	-	Up to 10	Sterile
D-50TP-XX	Stainless Steel (DIN 1.4197)	-	Up to 10	Sterile
D-60TP-XX	Stainless Steel (DIN 1.4197)	-	Up to 10	Sterile
Guided				
D-33TP-XX-L	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-40TP-XX-L	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-50TP-XX-L	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
Dense Bone Drills				
D-42TP-XX	Stainless Steel (DIN 1.4197)	TiN	Up to 10	Sterile
D-52TP-XX	Stainless Steel (DIN 1.4197)	TiN	Up to 10	Sterile
D-62TP-XX	Stainless Steel (DIN 1.4197)	TiN	Up to 10	Sterile

TRI-NEX®				
D-L-35-XX	Stainless Steel (DIN 1.4197)	-	Up to 10	Sterile
D-L-43-XX	Stainless Steel (DIN 1.4197)	-	Up to 10	Sterile
D-L-50-XX	Stainless Steel (DIN 1.4197)	-	Up to 10	Sterile
D-L-60-XX	Stainless Steel (DIN 1.4197)	-	Up to 10	Sterile
Soft Bone Drills				
DLS-35-XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
DLS-43-XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
DLS-50-XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
DLS-60-XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
DEEP CONICAL (DC)				
D-DCT-30XX	Titanium grade 5 (ASTM F136)	ALTiN	Up to 10	Sterile
D-DCT-35XX	Titanium grade 5 (ASTM F136)	ALTiN	Up to 10	Sterile
D-DCT-40XX	Titanium grade 5 (ASTM F136)	ALTiN	Up to 10	Sterile
D-DCT-50XX	Titanium grade 5 (ASTM F136)	ALTiN	Up to 10	Sterile
M-SERIES/INTERNAL HEX				
D-MT37XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-MT42XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-MT50XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
INTERNAL OCTAGON (IT)				
D-4XXT	Stainless Steel (DIN 1.4197)	-	Up to 10	Sterile
D-5XXT	Stainless Steel (DIN 1.4197)	-	Up to 10	Sterile
D-6XXT	Stainless Steel (DIN 1.4197)	-	Up to 10	Sterile
INVERTA®				
D-IV2700	Titanium grade 4 (ASTM F67)	Anodized Silver	Up to 10	Sterile
D-IV37XX	Titanium grade 4 (ASTM F67)	Anodized Pink	Up to 10	Sterile
D-IV45XX	Titanium grade 5 (ASTM F136)	Anodized Blue	Up to 10	Sterile
D-IV50XX	Titanium grade 5 (ASTM F136)	Anodized Yellow	Up to 10	Sterile
D-IV60XX	Titanium grade 5 (ASTM F136)	Anodized Silver	Up to 10	Sterile
Guided				
D-IV37XXGS	Titanium grade 4 (ASTM F67)	Anodized Pink	Up to 10	Sterile
D-IV45XXGS	Titanium grade 5 (ASTM F136)	Anodized Blue	Up to 10	Sterile
D-IV50XXGS	Titanium grade 5 (ASTM F136)	Anodized Yellow	Up to 10	Sterile
D-IV60XXGS	Titanium grade 5 (ASTM F136)	Anodized Silver	Up to 10	Sterile
MAX (MAX, PROMAX, TRI-MAX®, MAXIT®)				
D-MAX6-XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-70TP-XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-70TP-XX-L	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-80TP-XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-80TP-XX-L	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-90TP-XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-90TP-XX-L	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
IE IMPLANTS				
D-40E-04F	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
SINGLE PLATFORM (SP1)				
D-SP30/L	Titanium grade 5 (ASTM F136)	ALTiN	Up to 10	Non-sterile
D-SP35/L	Titanium grade 5 (ASTM F136)	ALTiN	Up to 10	Non-sterile
D-SP40/L	Titanium grade 5 (ASTM F136)	ALTiN	Up to 10	Non-sterile
D-SP45/L	Titanium grade 5 (ASTM F136)	ALTiN	Up to 10	Non-sterile
D-SP50/L	Titanium grade 5 (ASTM F136)	ALTiN	Up to 10	Non-sterile
DS-SP30/L	Titanium grade 5 (ASTM F136)	Anodized Blue	1	Sterile
DS-SP35/L	Titanium grade 5 (ASTM F136)	Anodized Blue	1	Sterile
DS-SP40/L	Titanium grade 5 (ASTM F136)	Anodized Blue	1	Sterile
DS-SP45/L	Titanium grade 5 (ASTM F136)	Anodized Blue	1	Sterile

DS-SP50/L	Titanium grade 5 (ASTM F136)	Anodized Blue	1	Sterile
Guided				
D-SP30XX-GSN	Titanium grade 5 (ASTM F136)	ALTiN	Up to 10	Non-sterile
D-SP30XX-GS	Titanium grade 5 (ASTM F136)	ALTiN	Up to 10	Non-sterile
D-SP35XX-GSN	Titanium grade 5 (ASTM F136)	ALTiN	Up to 10	Non-sterile
D-SP35XX-GS	Titanium grade 5 (ASTM F136)	ALTiN	Up to 10	Non-sterile
D-SP40XX-GS	Titanium grade 5 (ASTM F136)	ALTiN	Up to 10	Non-sterile
D-SP45XX-GS	Titanium grade 5 (ASTM F136)	ALTiN	Up to 10	Non-sterile
D-SP50XX-GS	Titanium grade 5 (ASTM F136)	ALTiN	Up to 10	Non-sterile
DS-SP30XX-GSN	Titanium grade 5 (ASTM F136)	Anodized Blue	1	Sterile
DS-SP30XX-GS	Titanium grade 5 (ASTM F136)	Anodized Blue	1	Sterile
DS-SP35XX-GSN	Titanium grade 5 (ASTM F136)	Anodized Blue	1	Sterile
DS-SP35XX-GS	Titanium grade 5 (ASTM F136)	Anodized Blue	1	Sterile
DS-SP40XX-GS	Titanium grade 5 (ASTM F136)	Anodized Blue	1	Sterile
DS-SP45XX-GS	Titanium grade 5 (ASTM F136)	Anodized Blue	1	Sterile
DS-SP50XX-GS	Titanium grade 5 (ASTM F136)	Anodized Blue	1	Sterile

Table B – Taps

Drill Code	Material	Coating (if any)	Number of Uses	How Provided
EXTERNAL HEX & PROVATA®				
D-TAP-IBN	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-IBNT	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-IBS	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-IBT	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-I4B	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-BA	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-BAT	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-BBBS	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-BBBT	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
TRI-NEX®				
D-TAP-L-35	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-L-43	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-L-50	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-L-60	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-LS-35	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-LS-43	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-LS-50	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
DEEP CONICAL (DC)				
D-TAP-DCC30	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-DCC35	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-DCC40	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-DCC50	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
INTERNAL OCTAGON (IT)				
D-TAP-ITT4	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-ITT5	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-ITT6	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-ITC3	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-ITC4	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-ITC5	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
INVERTA®				
IV-TR45XX	Titanium Grade 4 (ASTM F67)	Anodized Blue	Up to 10	Sterile
IV-TR50XX	Titanium Grade 4 (ASTM F67)	Anodized Yellow	Up to 10	Sterile
IV-TR60XX	Titanium Grade 4 (ASTM F67)	Anodized Silver	Up to 10	Sterile

MAX (MAX, PROMAX, TRI-MAX®, MAXIT®)				
D-TAP-MAX6-XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-MAX7-XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-MAX8-XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-MAX9-XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-MAX10-XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
FIRST® ROHNER TECHNIQUE				
D-TAP-35RT	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-43RT	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile

Clinical benefits

The Tapered Drills and Taps are used as part of dental implant therapy, a procedure in which a patient's missing teeth and/or crowns are restored. As a result, the benefits of the Tapered Drills and Taps mirror those of dental implant therapy in general, including improved chewing function, speech, aesthetics, and patient psychological wellbeing.

The clinical benefits specific to drills include the effective cutting of bone within acceptable limits of time, precise preparation of the implant site in the dimensions of the intended implants and the ability to remain effective after multiple uses to the extent stipulated by the manufacturer.

The clinical benefit specific to the handpiece devices includes the successful achievement of its respective intended function.

Storage, cleaning and sterilisation

The component is supplied sterile (sterilised by gamma irradiation) and intended for single-use prior to the expiration date (see packaging label). Sterility is assured unless the container or seal is damaged or opened. If packaging is damaged do not use the product and contact your Southern representative or return to Southern Implants®. The devices must be stored in a dry place at room temperature and not exposed to direct sunlight. Incorrect storage may influence device characteristics. Do not reuse components indicated for single-use only. Reusing these components may result in:

- damage to the surface or critical dimensions, which may result in performance and compatibility degradation.
- adds the risk of cross-patient infection and contamination if single-use items are reused.

Southern Implants® does not accept any responsibility for complications associated with reused single-use components.

Limitations on reusable items

A direct value for reusable instruments cannot be given. Frequent processing may have minor effects on the instruments. The product life is normally determined by wear and damage during use, thus instruments if properly cared for and inspected after each use, can be reused many times. Maintain a checklist for these instruments recording the number of uses.

Prior to reprocessing the device it should be thoroughly inspected and tested to determine its suitability for reuse.

NOTE: during use, handle drills and instruments with sterile tweezers to minimise contamination of the instrument tray and risk of damage to sterile surgical gloves.

Containment

As soon as practically possible, remove all visible residue after use (bone, blood or tissue), by immersing the instrument in cold water (dried soil may be difficult to remove).

Precleaning

Disassemble instruments from hand pieces and all connecting parts from instruments in order to clean soil from obstructed areas. Remove the PEEK bits from placement tools. Rinse with lukewarm water for 3 minutes and remove hardened debris with a soft nylon brush. Avoid mechanical damage during cleaning.

Manual cleaning or automated cleaning

Prepare an ultrasonic bath with suitable detergent (i.e. Steritech instrument cleaner - 5% dilution), sonicate for 20 minutes (alternative methods can be used if proven by the end user). Rinse with purified/sterile water.

NOTE: always follow the instructions for use of the manufacturers for cleaning agents and disinfectants.

Load devices into a thermidisinfector. Run the cleaning and disinfection cycle, followed by the drying cycle.

Drying

Dry the instruments both inside and outside with filtered compressed air or single use, lint free wipes. Pack the instruments as quickly as possible after removal into the storage container. If additional drying is necessary, dry in a clean location.

Inspection

Do a visual inspection of the items to check for damages on instruments.

Packaging

Use the correct packaging material as indicated for steam sterilisation to ensure sterility is maintained. Double packaging is recommended. Where appropriate, the cleaned, disinfected and checked devices can be assembled and placed in instrument trays as applicable. Instrument trays can be double wrapped or placed in sterilisation bags.

Sterilisation

Southern Implants® recommends one of the following procedures to sterilise the device prior to use:

1. prevacuum sterilisation method: steam sterilise the device at 132°C (270°F) at 180 - 220 kPa for 4 minutes. Dry for at least 20 minutes in the chamber. Only an approved wrap or pouch for steam sterilisation must be used.
2. for users in the USA: prevacuum sterilisation method: wrapped, steam sterilise at 135°C (275°F) at 180 - 220 kPa for 3 minutes. Dry for 20 minutes in the chamber. Use a wrap or pouch that is cleared for the indicated steam sterilisation cycle.

NOTE: users in the USA must ensure that the steriliser, wrap or pouch, and all steriliser accessories are cleared by the FDA, for the intended sterilisation cycle.

Storage

Maintain packaging integrity to ensure sterility in storage. Packaging should be completely dry before storage to avoid corrosion and degradation of cutting edges.

Contraindications

Since the Tapered Drills and Taps are required or utilized as part of dental implant surgery, the complications for this device are not dissimilar to those of dental implant therapy as a whole. These contraindications include:

- patients medically unfit for oral surgical procedures.
- cases where an adequate number of implants cannot be placed to provide full functional support of the prosthesis.
- patients under the age of 18 years.
- patients with poor bone quality.
- patients with blood disorders.
- presence of infection at the implant site.
- patients with vascular impairment.
- patients with uncontrolled diabetes.
- patients with drug or alcohol abuse dependencies.
- patients undergoing radiotherapy treatment.
- patients undergoing chronic high dose steroid therapy.
- patients with metabolic bone disease.

- patients with pure titanium, titanium alloy (Ti6Al4V), gold, palladium, polyether ether ketone (PEEK) or iridium allergies or hypersensitivities.

Other than the above, there are no contraindications unique to this system.

Warnings and precautions

THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING

- For the safe and effective use of dental implants, new technologies/systems and the Tapered drills and Taps, it is strongly suggested that specialised training be undertaken. This training should include hands-on methods to gain competency on the proper technique, biomechanical requirements and radiographic evaluations required for the specific system.
- Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone.
- Use of the device with incompatible or non-corresponding devices can result in device failure or poor performance.
- When handled intraorally, devices must be adequately secured against aspiration. Aspiration of products may lead to infection or physical injury.
- Electrosurgery should not be attempted around metal implants as they are conductive.
- The use of non-sterile items can lead to secondary infections of the tissue or transfer infectious diseases.
- Failure to perform appropriate cleaning, re-sterilisation and storage procedures as per the Instructions for Use document can result in device damage or secondary infection.
- Exceeding the number of recommended uses for reusable devices can result in device damage, secondary infection or patient harm,
- Blunt drills may cause damage to the bone which could compromise osseointegration.
- For short implants, clinicians should closely monitor patients for any of the following conditions: peri implant bone loss, changes to implant's response to percussion, or radiographic changes in bone to implant contact along the implant's length. If the implant shows mobility or greater than 50% bone loss, the implant should be evaluated for possible removal. If the clinicians choose a short implant, then clinicians should consider a two-stage surgical approach, splinting a short implant to an additional implant, and placement of the widest possible fixture. Allow longer periods for osseointegration and avoid immediate loading.

Note: Training should be undertaken by both new and experienced implant users before using a new system or attempting to do a new treatment method.

It is important that thorough screening of prospective implant candidates is performed. This screening should include at least:

- a comprehensive medical and dental history.
- visual and radiological inspection to determine adequate bone dimensions, anatomical landmarks, occlusal conditions and periodontal health. Failure to accurately estimate lengths of drills relative to radiographic measurements can result in permanent injury to nerves or other vital structures.
- the evaluation and consideration of bruxism and unfavourable jaw relations.

Proper preoperative planning with a good team approach between well trained surgeons, restorative dentists and lab technicians is essential for successful implant treatment.

For successful implant treatment, it is important to:

- minimise the trauma to the host tissue as this increases the potential for successful osseointegration.
- be aware and avoid damage to vital anatomical structures such as nerves, veins and arteries. Injury to these structures may cause serious complications, including, injury to the eye, nerve damage and excessive bleeding. It is essential to protect the infraorbital nerve.
- identify actual measurements relative to the radiographic data as failing to do so could lead to complications.

Special care should be taken when treating patients with local or systemic risk factors that could affect the healing of the bone and soft tissue or otherwise increase the severity of the side effects, risk of complications and/or implant failure. These factors include:

- smoking.
- history of periodontal disease.
- poor oral hygiene.
- history of orofacial radiotherapy.**
- bruxism and unfavourable jaw relations.

**** The potential for implant failure and other complications increases when implants are placed in irradiated bone as radiotherapy can result in progressive fibrosis of vessels and soft tissue (i.e., osteoradionecrosis), leading to diminished healing capacity. Implant treatment of irradiated patients is dependent upon issues like the timing of implant placement in relation to the radiation therapy, anatomic site chosen for implant placement and radiation dosage at that site and consequent risk of osteoradionecrosis.**

The responsibility for proper patient selection, adequate training and experience in implant placement and the provision of the appropriate information required for informed consent rests with the practitioner.

Should the device not operate as intended, it must be reported to the manufacturer of the device. The contact information or the manufacturer of this device to report a change in performance is:

sicomplaints@southernimplants.com.

Side effects

Since the device is required or utilized as part of dental implant surgery, the expected side effects are not dissimilar to those of dental implant therapy as a whole. Common side effects of dental implant therapy include pain, local inflammation, phonetic difficulties and gingival inflammation. Other, less common side effects or complications to dental implant therapy include, but are not limited to: (1) wound dehiscence; (2) peri-implantitis; (3) transient weakness, numbness, and/or pain associated with mild nerve injury; (4) fat emboli formation; (5) marginal bone loss within acceptable limits; (6) allergic reaction(s) to implant and/or abutment material; (7) unspecified infection; (8) implant failure due to insufficient levels of osseointegration; (9) breakage of the implant and/or abutment; (10) loosening of the abutment screw and/or retaining screw; (11) complications requiring revision of the dental implant; (12) nerve damage resulting in permanent weakness, numbness, or pain; (13) loosening of the implant requiring revision surgery; (14) perforation of the maxillary sinus; (15) perforation of the labial and/or lingual plates and (16) bone loss possibly resulting in revision or removal of the implant.

Precaution: maintaining sterility protocol

Implants are packaged as follows:

1. An outer package consisting of a rigid, clear box which acts as protection for the inner package.
2. The inner package consisting of a blister pack (clear plastic-formed blister base with a TYVEK "peel-back" lid).
3. Within the inner package, there is a hollow tube which contains one implant suspended from a titanium ring, this ensures the implant never touches the inside of the plastic tube.
4. Labelling information is located on the surface of the peel-back lid and on the outside of the rigid box.

Care must be taken to maintain the sterility of the implant by proper opening of the packaging and handling of the implant.

1. Open the implant package in the non-sterile field, with non-sterile gloves, tear the address label to open the box.
2. With non-sterile gloves, remove the inner blister pack. Do not place the plastic box or blister pack-lid onto the sterile field. The contents of this inner package are sterile.
3. The sealed blister is to be opened by an assistant (with nonsterile gloves): remove the TYVEK lid and drop or place the sterile tube onto the sterile field, open the tube cap and attach the implant placement tool onto the implant and carefully remove from the sterile tube. Do not touch the sterile implant.

Other sterile components are packed in a peel pouch or blister base with a "peel-back" lid. Labelling information is located on the bottom half of the pouch, inside the packet or on the surface of the peel-back lid. Sterility is assured unless the pouch is damaged or opened. Non-sterile components are supplied clean but not sterile in a peel pouch or blister base with peel-back lid. Labelling information is located on the bottom half of the pouch or on the surface of the peel-back lid.

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established.

The contact information for the manufacturer of this device to report a serious incident is as follows:

sicomplaints@southernimplants.com.

Materials

Tapered Drills: Titanium Grade 4 (ASTM F67), Titanium grade 5 (ASTM F136)
or Stainless Steel (AISI 440A/DIN 1.4197)

Taps: Titanium grade 5 (ASTM F136) or Titanium Grade 4 (ASTM F67)

Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

Summary of Safety and Clinical Performance (SSCP)

As required by the European Medical Device Regulation (MDR; EU2017/745), a Summary of Safety and Clinical Performance (SSCP) is available for perusal with regard to Southern Implants® product ranges.

The relevant SSCP can be accessed at <https://ec.europa.eu/tools/eudamed>.

NOTE: the above website will be available upon the launch of the European Database on Medical Devices (EUDAMED).

Disclaimer of liability

This product is part of the Southern Implants® product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants® product range and take full responsibility for the correct indications and use of this product. Southern Implants® does not assume liability for damage due to incorrect use. Please note that some Southern Implants® products may not be cleared or released for sale in all markets.

Basic UDI

Product	Basic-UDI Number
Basic-UDI for Drills and Handpiece Devices	6009544038759C
Basic-UDI for Reusable Instruments	6009544038769E

Related literature and catalogues

CAT-2004 - Tri-Nex® Implants Product Catalogue
 CAT-2020 - External Hex Implants Product Catalogue
 CAT-2042 - Deep Conical (DC) Implants Product Catalogue
 CAT-2043 - Internal Hex (M-Series) Implants Product Catalogue
 CAT-2060 - Internal Hex (PROVATA®) Implants Product Catalogue
 CAT-2005 - IT (Internal Octagon) Implants Product Catalogue
 CAT-2069 - INVERTA® Implants Product Catalogue
 CAT-2068 - SIGuided Implants Product Catalogue
 CAT-2010 - Osseointegrated Fixtures Product Catalogue
 CAT-2091M - FIRST® Fibula Reconstruction System Technology
 CAT-2093 - Single Platform Implants Product Catalogue

Symbols and warnings

											
Manufacturer: Southern Implants® 1 Albert Rd, P.O Box 605 IRENE, 0062, South Africa. Tel: +27 12 667 1046	CE mark 2797	Prescription device*	Sterilised using irradiation	Non-sterile	Use by date (mm-yy)	Do not reuse	Do not resterilise	Catalogue number	Batch code	Medical device	Authorised representative in the European Community
											
Authorised representative for Switzerland	Date of manufacture	Magnetic Resonance conditional	Magnetic Resonance safe	Single sterile barrier system with protective packaging inside	Single sterile barrier system	Consult instruction for use	Caution	Keep away from sunlight	Do not use if package is damaged		

* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

All rights reserved. Southern Implants®, the Southern Implants® logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants®. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale. It is the responsibility of the clinician to inspect the symbols that appear on the packaging of the product in use.

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