English	INSTRUCTIONS FOR USE: Southern Implants® Tapered Drills and Taps
Español	INSTRUCCIONES DE USO: Fresas cónicas y machos Southern Implants®
Italiano	ISTRUZIONI PER L'USO: Southern Implants® Punte e Maschi Conici
Français	MODE D'EMPLOI : Forets et Tarauds coniques Southern Implants®
Deutsch	GEBRAUCHSANWEISUNG: Southern Implants® Kegelbohrer und Gewindeschneider
Português	INSTRUÇÕES DE UTILIZAÇÃO: Southern Implants® Brocas e machos cônicos
Ελληνικά	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ: Southern Implants® Κωνικά τρύπανα και τάπες
Svenska	BRUKSANVISNING: Southern Implants® Avsmalnande Borrar och Kranar
Română	INSTRUCȚIUNI DE UTILIZARE: Burghiuri și robinete conice Southern Implants®
Türkçe	KULLANMA TALİMATI: Southern Implants® Konik Frezler ve Kılavuzlar

South Africa - Headquarters: 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

Subsidiaries

Australia

Southern Implants Australia
T: +61-(0)-8-9466-2627
E: info@southernimplants.com.au

Spain and Portugal

Southern Implants Iberica
T: +34 935 053 507
E: info@southernimplants.es

United Kingdom and Ireland

Southern Implants UK
T: +44-20 8059 4490
E: info@southernimplants.co.uk

USA and Canada

Southern Implants North America Inc.
T: +1-561-472-0990
E: customercare@southernimplants.com

Description

Southern Implants® Tapered Drills are intended to enlarge a channel made by a twist drill in order to prepare a conical osteotomy site for a tapered endosseous dental implant. This entails drilling a hole in the mandible or maxilla or making an existing hole larger in the conical dimension. The drills are available as multiple use and single-use versions, supplied either sterile or non-sterile. They are made from Stainless Steel and Titanium grade 4 or 5. The drills are available in diameters between Ø2.7mm and Ø9mm.

Southern Implants® Taps are intended to be used to create a thread in the bone identical to that of the corresponding implant prior to placing the implant. The taps are for multiple use and supplied sterile. They are made from Titanium grade 4 or 5. The taps are available in diameters between Ø3mm and Ø9mm.

Intended use

The Southern Implants® Drills and Handpiece devices are intended to prepare the implant site/osteotomy and aid in the placement/removal of the endosseous implant and/or prosthetic device(s).

Specifically, the Southern Implants® Tapered Drills are intended to enlarge the initial channel created by a Twist Drill. The Southern Implants® Tap Drills are intended to cut threads into the bone, facilitate the insertion of the target implant(s) and/or reduce the torque required for implant placement.

Indications for use

The Southern Implants® Drills and Handpiece Devices are indicated for the functional and aesthetic treatment of patients who are eligible for implant-retained prostheses.

Intended user

Maxillofacial Surgeons, General Dentists, Orthodontists, Periodontists, Prosthodontists, and other appropriately trained and/or experienced medical professionals.

Intended environment

The devices are intended to be used in a clinical environment such as an operating theatre or a dentist consultation room.

Intended patient population

The Southern Implants® Drills and Handpiece devices are intended to be used in patients subject to implant treatment.

Compatibility information

Southern Implants' implants tapered drills and taps are as described in Table A and B respectively. Details on site preparation sequences recommended by Southern Implants® can be found in the relevant product catalogues.

Table A – Tapered Drills

Drill Code	Material	Coating (if any)	Number of Uses	How Provided
EXTERNAL HEX & PROVATA®				
D-30TP-XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-33TP-XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-40TP-XX	Stainless Steel (DIN 1.4197)	-	Up to 10	Sterile
D-50TP-XX	Stainless Steel (DIN 1.4197)	-	Up to 10	Sterile
D-60TP-XX	Stainless Steel (DIN 1.4197)	-	Up to 10	Sterile
Guided				
D-33TP-XX-L	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-40TP-XX-L	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-50TP-XX-L	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile

Dense Bone Drills				
D-42TP-XX	Stainless Steel (DIN 1.4197)	TiN	Up to 10	Sterile
D-52TP-XX	Stainless Steel (AISI 440A)	TiN	Up to 10	Sterile
D-62TP-XX	Stainless Steel (AISI 440A)	TiN	Up to 10	Sterile
TRI-NEX®				
D-L-35-XX	Stainless Steel (DIN 1.4197)	-	Up to 10	Sterile
D-L-43-XX	Stainless Steel (DIN 1.4197)	-	Up to 10	Sterile
D-L-50-XX	Stainless Steel (DIN 1.4197)	-	Up to 10	Sterile
D-L-60-XX	Stainless Steel (DIN 1.4197)	-	Up to 10	Sterile
Soft Bone Drills				
DLS-35-XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
DLS-43-XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
DLS-50-XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
DLS-60-XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
DEEP CONICAL (DC)				
D-DCT-30XX	Titanium grade 5 (ASTM F136)	ALTiN	Up to 10	Sterile
D-DCT-35XX	Titanium grade 5 (ASTM F136)	ALTiN	Up to 10	Sterile
D-DCT-40XX	Titanium grade 5 (ASTM F136)	ALTiN	Up to 10	Sterile
D-DCT-50XX	Titanium grade 5 (ASTM F136)	ALTiN	Up to 10	Sterile
M-SERIES/INTERNAL HEX				
D-MT37XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-MT42XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-MT50XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
INTERNAL OCTAGON (IT)				
D-4XXT	Stainless Steel (DIN 1.4197)	-	Up to 10	Sterile
D-5XXT	Stainless Steel (DIN 1.4197)	-	Up to 10	Sterile
D-6XXT	Stainless Steel (DIN 1.4197)	-	Up to 10	Sterile
INVERTA®				
D-IV2700	Titanium grade 4 (ASTM F67)	Anodized Silver	Up to 10	Sterile
D-IV37XX	Titanium grade 4 (ASTM F67)	Anodized Pink	Up to 10	Sterile
D-IV45XX	Titanium grade 5 (ASTM F136)	Anodized Blue	Up to 10	Sterile
D-IV50XX	Titanium grade 5 (ASTM F136)	Anodized Yellow	Up to 10	Sterile
D-IV60XX	Titanium grade 5 (ASTM F136)	Anodized Silver	Up to 10	Sterile
Guided				
D-IV37XXGS	Titanium grade 4 (ASTM F67)	Anodized Pink	Up to 10	Sterile
D-IV45XXGS	Titanium grade 5 (ASTM F136)	Anodized Blue	Up to 10	Sterile
D-IV50XXGS	Titanium grade 5 (ASTM F136)	Anodized Yellow	Up to 10	Sterile
D-IV60XXGS	Titanium grade 5 (ASTM F136)	Anodized Silver	Up to 10	Sterile
MAX (MAX, PROMAX, TRI-MAX®, MAXIT®)				
D-MAX6-XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-70TP-XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-70TP-XX-L	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-80TP-XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-80TP-XX-L	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile

D-90TP-XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-90TP-XX-L	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
IE IMPLANTS				
D-40E-04F	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
SINGLE PLATFORM (SP1)				
D-SP30/L	Titanium grade 5 (ASTM F136)	ALTiN	Up to 10	Non-sterile
D-SP35/L	Titanium grade 5 (ASTM F136)	ALTiN	Up to 10	Non-sterile
D-SP40/L	Titanium grade 5 (ASTM F136)	ALTiN	Up to 10	Non-sterile
D-SP45/L	Titanium grade 5 (ASTM F136)	ALTiN	Up to 10	Non-sterile
D-SP50/L	Titanium grade 5 (ASTM F136)	ALTiN	Up to 10	Non-sterile
DS-SP30/L	Titanium grade 5 (ASTM F136)	Anodized Blue	1	Sterile
DS-SP35/L	Titanium grade 5 (ASTM F136)	Anodized Blue	1	Sterile
DS-SP40/L	Titanium grade 5 (ASTM F136)	Anodized Blue	1	Sterile
DS-SP45/L	Titanium grade 5 (ASTM F136)	Anodized Blue	1	Sterile
DS-SP50/L	Titanium grade 5 (ASTM F136)	Anodized Blue	1	Sterile
Guided				
D-SP30XX-GSN	Titanium grade 5 (ASTM F136)	ALTiN	Up to 10	Non-sterile
D-SP30XX-GS	Titanium grade 5 (ASTM F136)	ALTiN	Up to 10	Non-sterile
D-SP35XX-GSN	Titanium grade 5 (ASTM F136)	ALTiN	Up to 10	Non-sterile
D-SP35XX-GS	Titanium grade 5 (ASTM F136)	ALTiN	Up to 10	Non-sterile
D-SP40XX-GS	Titanium grade 5 (ASTM F136)	ALTiN	Up to 10	Non-sterile
D-SP45XX-GS	Titanium grade 5 (ASTM F136)	ALTiN	Up to 10	Non-sterile
D-SP50XX-GS	Titanium grade 5 (ASTM F136)	ALTiN	Up to 10	Non-sterile
DS-SP30XX-GSN	Titanium grade 5 (ASTM F136)	Anodized Blue	1	Sterile
DS-SP30XX-GS	Titanium grade 5 (ASTM F136)	Anodized Blue	1	Sterile
DS-SP35XX-GSN	Titanium grade 5 (ASTM F136)	Anodized Blue	1	Sterile
DS-SP35XX-GS	Titanium grade 5 (ASTM F136)	Anodized Blue	1	Sterile
DS-SP40XX-GS	Titanium grade 5 (ASTM F136)	Anodized Blue	1	Sterile
DS-SP45XX-GS	Titanium grade 5 (ASTM F136)	Anodized Blue	1	Sterile
DS-SP50XX-GS	Titanium grade 5 (ASTM F136)	Anodized Blue	1	Sterile
BHORA GUIDED SURGERY				
D-30TP-BG	Titanium grade 5 (ASTM F136)	-	Up to 10	Non-sterile
D-33TP-BG	Titanium grade 5 (ASTM F136)	-	Up to 10	Non-sterile
D-35TP-BG	Titanium grade 5 (ASTM F136)	-	Up to 10	Non-sterile
D-40TP-BG	Titanium grade 5 (ASTM F136)	-	Up to 10	Non-sterile
D-45TP-BG	Titanium grade 5 (ASTM F136)	-	Up to 10	Non-sterile
D-50TP-BG	Titanium grade 5 (ASTM F136)	-	Up to 10	Non-sterile
D-30TP-BGN	Titanium grade 5 (ASTM F136)	-	Up to 10	Non-sterile
D-33TP-BGN	Titanium grade 5 (ASTM F136)	-	Up to 10	Non-sterile
D-35TP-BGN	Titanium grade 5 (ASTM F136)	-	Up to 10	Non-sterile
D-40TP-BGN	Titanium grade 5 (ASTM F136)	-	Up to 10	Non-sterile

Table B – Taps

Drill Code	Material	Coating (if any)	Number of Uses	How Provided
EXTERNAL HEX & PROVATA®				
D-TAP-IBN	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-IBNT	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-IBS	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-IBT	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-I4B	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-BA	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-BAT	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-BBBS	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-BBBT	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
TRI-NEX®				
D-TAP-L-35	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-L-43	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-L-50	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-L-60	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-LS-35	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-LS-43	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-LS-50	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
DEEP CONICAL (DC)				
D-TAP-DCC30	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-DCC35	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-DCC40	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-DCC50	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
INTERNAL OCTAGON (IT)				
D-TAP-ITT4	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-ITT5	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-ITT6	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-ITC3	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-ITC4	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-ITC5	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
MAX (MAX, PROMAX, TRI-MAX®, MAXIT®)				
D-TAP-MAX6-XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-MAX7-XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-MAX8-XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-MAX9-XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-MAX10-XX	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
FIRST® ROHNER TECHNIQUE				
D-TAP-35RT	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile
D-TAP-43RT	Titanium grade 5 (ASTM F136)	-	Up to 10	Sterile

Clinical performance

Southern Implants® Tapered Drills and Taps are designed to aid in implant site preparation by initiating and shaping a precisely profiled osteotomy for implant placement. Their specialised geometry, including appropriate lengths, diameters, and cutting features, ensures accurate and efficient site preparation.

Clinical benefits

Southern Implants® Tapered Drills and Taps play a vital role in implant-retained prosthetic treatment by enabling precise osteotomy preparation. While their clinical benefits are primarily tied to the overall success of the restorative treatment rather than their direct performance, their role in achieving accurate implant placement supports favourable outcomes.

The primary clinical benefits of this treatment approach include the effective restoration of the identified defect and enhanced patient satisfaction. Additionally, successful prosthetic rehabilitation can contribute to improved social interaction, psychosocial well-being, self-esteem, and confidence. Although these benefits are indirect and challenging to quantify, they are closely linked to positive treatment outcomes, such as optimal implant placement, functional restoration, and patient-reported satisfaction with both aesthetics and function.

Storage, cleaning and sterilisation

APPLICABLE FOR STERILE COMPONENTS:

The component is supplied sterile (sterilised by gamma irradiation) and may be intended for single-use prior to the expiration date (see packaging label). Sterility is assured unless the container or seal is damaged or opened. If packaging is damaged do not use the product and contact your Southern representative or return to Southern Implants®. The devices must be stored in a dry place at room temperature and not exposed to direct sunlight. Incorrect storage may influence device characteristics. Do not reuse components indicated for single-use only. Reusing these components may result in:

- damage to the surface or critical dimensions, which may result in performance and compatibility degradation.
- adds the risk of cross-patient infection and contamination if single-use items are reused.

Southern Implants® does not accept any responsibility for complications associated with reused single-use components.

APPLICABLE FOR REUSABLE COMPONENTS:

Limitations on reusable items

A direct value for reusable instruments cannot be given. Frequent processing may have minor effects on the instruments. The product life is normally determined by wear and damage during use, thus instruments if properly cared for and inspected after each use, can be reused many times. Maintain a checklist for these instruments recording the number of uses.

Prior to reprocessing the device, it should be thoroughly inspected and tested to determine its suitability for reuse.

NOTE: during use, handle drills and instruments with sterile tweezers to minimise contamination of the instrument tray and risk of damage to sterile surgical gloves.

Containment

As soon as practically possible, remove all visible residue after use (bone, blood or tissue), by immersing the instrument in cold water (dried soil may be difficult to remove).

Precleaning

Disassemble instruments from hand pieces and all connecting parts from instruments in order to clean soil from obstructed areas. Remove the PEEK bits from placement tools. Rinse with lukewarm water for 3 minutes and remove hardened debris with a soft nylon brush. Avoid mechanical damage during cleaning.

Manual cleaning or automated cleaning

Prepare an ultrasonic bath with suitable detergent (i.e. Steritech instrument cleaner - 5% dilution), sonicate for 20 minutes (alternative methods can be used if proven by the end user). Rinse with purified/sterile water.

NOTE: always follow the instructions for use of the manufacturers for cleaning agents and disinfectants.

Load devices into a thermidisinfector. Run the cleaning and disinfection cycle, followed by the drying cycle.

Drying

Dry the instruments both inside and outside with filtered compressed air or single use, lint free wipes. Pack the instruments as quickly as possible after removal into the storage container. If additional drying is necessary, dry in a clean location.

Inspection

Do a visual inspection of the items to check for damages on instruments.

Packaging

Use the correct packaging material as indicated for steam sterilisation to ensure sterility is maintained. Double packaging is recommended. Where appropriate, the cleaned, disinfected and checked devices can be assembled and placed in instrument trays as applicable. Instrument trays can be double wrapped or placed in sterilisation bags.

Sterilisation

Southern Implants® recommends one of the following procedures to sterilise the device prior to use:

1. prevacuum sterilisation method: steam sterilise the device at 132°C (270°F) at 180 - 220 kPa for 4 minutes. Dry for at least 20 minutes in the chamber. Only an approved wrap or pouch for steam sterilisation must be used.
2. for users in the USA: prevacuum sterilisation method: wrapped, steam sterilise at 135°C (275°F) at 180 - 220 kPa for 3 minutes. Dry for 20 minutes in the chamber. Use a wrap or pouch that is cleared for the indicated steam sterilisation cycle.

NOTE: users in the USA must ensure that the steriliser, wrap or pouch, and all steriliser accessories are cleared by the FDA, for the intended sterilisation cycle.

Storage

Maintain packaging integrity to ensure sterility in storage. Packaging should be completely dry before storage to avoid corrosion and degradation of cutting edges.

Contraindications

The contraindications of all device groups used as part of the specific treatment or procedure apply. Therefore, the contraindications of the systems/medical devices utilized as part of implant surgery/therapy should be noted and the relevant documents consulted.

Contraindications specific to the Southern Implants® Drills and Handpiece Devices:

- Allergies or hypersensitivity to chemical ingredients in the following materials: Titanium, Aluminium, Vanadium, Stainless Steel.

Warnings and precautions

IMPORTANT NOTICE: THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING

- To ensure the safe and effective use of dental implants, new technologies/systems and the Metal Abutment devices, it is strongly recommended that specialised training be undertaken. This training should include hands-on methods to gain competency on the proper technique, the system's biomechanical requirements and radiographic evaluations required for the specific system.

- Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone.
- Use of the device with incompatible or non-corresponding devices can result in the poor performance or failure of the device.
- When handling devices intraorally, it is imperative that they are adequately secured to prevent aspiration, as aspiration of products may lead to infection or physical injury.
- The use of non-sterile items can lead to secondary infections of the tissue or transfer infectious diseases.
- Failure to adhere to appropriate cleaning, re-sterilization, and storage procedures as outlined in the Instructions for Use can result in device damage, secondary infections, or patient harm.
- Exceeding the number of recommended uses for reusable devices can result in device damage, secondary infection or patient harm.
- The use of blunt drills may cause damage to the bone, potentially compromising osseointegration.

It is crucial to emphasize that training should be undertaken by both new and experienced implant users prior to utilizing a new system or attempting a new treatment method.

Patient Selection and Preoperative Planning

A comprehensive patient selection process and meticulous preoperative planning are essential for successful implant treatment. This process should involve consultation among a multidisciplinary team, including well-trained surgeons, restorative dentists, and laboratory technicians.

Patient screening should include, at a minimum, a thorough medical and dental history, as well as visual and radiological inspections to assess the presence of adequate bone dimensions, the positioning of anatomical landmarks, the presence of unfavourable occlusal conditions, and the periodontal health status of the patient.

For successful implant treatment, it is important to:

1. Minimize trauma to the host tissue, thereby increasing the potential for successful osseointegration.
2. Accurately identify measurements relative to radiographic data, as failure to do so may lead to complications.
3. Be vigilant in avoiding damage to vital anatomical structures, such as nerves, veins, and arteries. Injury to these structures may result in serious complications, including ocular injury, nerve damage, and excessive bleeding.

The responsibility for proper patient selection, adequate training and experience in implant placement, and the provision of appropriate information required for informed consent rests with the practitioner. By combining thorough screening of prospective implant candidates with a practitioner possessing a high level of competence in the use of the system, the potential for complications and severe side effects can be significantly reduced.

High Risk Patients

Special care should be exercised when treating patients with local or systemic risk factors that may adversely affect the healing of bone and soft tissue or otherwise increase the severity of side effects, the risk of complications, and/or the likelihood of implant failure. Such factors include:

- poor oral hygiene
- history of smoking/vaping/tobacco use
- history of periodontal disease
- history of orofacial radiotherapy**
- bruxism and unfavourable jaw relations
- use of chronic medications that may delay healing or increase the risk of complications including, but not limited to, chronic steroid therapy, anti-coagulant therapy, TNF- α blockers, bisphosphonate and cyclosporin

*** The potential for implant failure and other complications increases when implants are placed in irradiated bone, as radiotherapy can lead to progressive fibrosis of blood vessels and soft tissue (i.e., osteoradionecrosis), resulting in diminished healing capacity. Contributing factors to this increased risk include the timing of implant placement in relation to radiation therapy, the proximity of radiation exposure to the implant site, and the radiation dosage at that site.*

Should the device not operate as intended, it must be reported to the manufacturer of the device. The contact information or the manufacturer of this device to report a change in performance is: sicomplaints@southernimplants.com.

Side effects

The clinical outcome of treatment is influenced by multiple factors. The following side effects and residual risks are associated with this device group or the surgical procedures involved in its use. These may require further treatment, revision surgery, or additional consultations with the relevant medical professional.

- Allergic or hypersensitivity reaction(s)
- Anaesthesia, paraesthesia, hyperesthesia, and hypoesthesia (transient or permanent)
- Anatomical landmark damage
- Bruising
- Dental injury during surgery
- Fat emboli formation
- Fistula formation
- Gingival injury
- Infection (acute and/or chronic)
- Localized inflammation
- Loss or damage to adjacent dentition
- Mild thermal injury to adjacent structures (e.g. burn to the lip)
- Nasal bleeding (Zygomatic Implant procedures)
- Nerve damage
- Pain, discomfort or tenderness
- Perforation of the labial or lingual plates
- Peri-implantitis, peri-mucositis or otherwise poor peri-implant soft tissue health
- Sinus Perforation
- Soft tissue irritation
- Surgical side effects such as pain, inflammation, bruising and mild bleeding
- Thermal injury to bone resulting in insufficient levels of osseointegration
- Wound dehiscence or poor healing

Additionally, the normal side effects associated with anaesthesia should also be expected.

Note: these side effects and residual risks can vary in both severity and frequency.

Precaution: maintaining sterility protocol

APPLICABLE TO STERILE COMPONENTS

The devices are packed in a peel pouch or blister base with a “peel-back” lid. Labelling information is located on the bottom half of the pouch, inside the packet or on the surface of the peel-back lid. Sterility is assured unless the pouch is damaged or opened.

APPLICABLE TO NON-STERILE COMPONENTS

The devices are supplied clean but not sterile in a peel pouch or blister base with peel-back lid. Labelling information is located on the bottom half of the pouch or on the surface of the peel-back lid.

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established.

The contact information for the manufacturer of this device to report a serious incident is as follows:

sicomplaints@southernimplants.com.

Materials

Tapered Drills:	Commercially Pure Titanium Grade 4 according to ASTM F67 and ISO 5832-2, Titanium alloy Ti-6Al-4V (90% titanium, 6% aluminium, 4% vanadium) according to ASTM F136 and ISO 5832-3 or Stainless Steel according to AISI 440A/DIN 1.4197
Taps:	Titanium alloy Ti-6Al-4V (90% titanium, 6% aluminium, 4% vanadium) according to ASTM F136 and ISO 5832-3

Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

Summary of Safety and Clinical Performance (SSCP)

As required by the European Medical Device Regulation (MDR; EU2017/745), a Summary of Safety and Clinical Performance (SSCP) is available for perusal with regard to Southern Implants® product ranges.

The relevant SSCP can be accessed at <https://ec.europa.eu/tools/eudamed>.

NOTE: the above website will be available upon the launch of the European Database on Medical Devices (EUDAMED).

Disclaimer of liability

This product is part of the Southern Implants® product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants® product range and take full responsibility for the correct indications and use of this product. Southern Implants® does not assume liability for damage due to incorrect use. Please note that some Southern Implants® products may not be cleared or released for sale in all markets.

Basic UDI

Product	Basic-UDI Number
Basic-UDI for Drills and Handpiece Devices	6009544038759C

Related literature and catalogues

CAT-2004 - Tri-Nex® Implants Product Catalogue
 CAT-2020 - External Hex Implants Product Catalogue
 CAT-2042 - Deep Conical (DC) Implants Product Catalogue
 CAT-2043 - Internal Hex (M-Series) Implants Product Catalogue
 CAT-2060 - Internal Hex (PROVATA®) Implants Product Catalogue
 CAT-2005 - IT (Internal Octagon) Implants Product Catalogue
 CAT-2069 - Deep Conical INVERTA® Implants Product Catalogue
 CAT-2068 - SIGuided Implants Product Catalogue
 CAT-2010 - Osseointegrated Fixtures Product Catalogue
 CAT-2091M - FIRST® Fibula Reconstruction System Technology
 CAT-2093 - Single Platform Implants Product Catalogue
 CAT-2095 - External Hex INVERTA® Implants Product Catalogue

Symbols and warnings

											
Manufacturer: Southern Implants® 1 Albert Rd, P.O Box 605 IRENE, 0062, South Africa. Tel: +27 12 667 1046	CE mark 2797	Prescription device*	Sterilised using irradiation	Non-sterile	Use by date (mm-yy)	Do not reuse	Do not resterilise	Catalogue number	Batch code	Medical device	Authorised representative in the European Community
											
Authorised representative for Switzerland	Date of manufacture	Magnetic Resonance conditional	Magnetic Resonance safe	Single sterile barrier system with protective packaging inside	Single sterile barrier system	Consult instruction for use	Caution	Keep away from sunlight	Do not use if package is damaged		

* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

All rights reserved. Southern Implants®, the Southern Implants® logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants®. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale. It is the responsibility of the clinician to inspect the symbols that appear on the packaging of the product in use.

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