English	INSTRUCTIONS FOR USE: Southern Implants® Trephine Drills
Español	INSTRUCCIONES DE USO: Fresas de Trépano Southern Implants®
Italiano	ISTRUZIONI PER L'USO: Southern Implants® Punte Trephine
Français	NOTICE D'UTILISATION : Fraises creuses de Southern Implants®
Deutsch	GEBRAUCHSANWEISUNG: Southern Implants® Trepanbohrer
Português	INSTRUÇÕES DE UTILIZAÇÃO: Southern Implants® Brocas de trefina
Ελληνικά	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ: Τρύπανα διάτρησης της Southern Implants®
Svenska	BRUKSANVISNING: Southern Implants® Trephine borrar
Română	INSTRUCȚIUNI DE UTILIZARE: Burghiu Southern Implants® Trephine
Türkçe	KULLANMA TALİMATI: Southern Implants® Trefin Frezler

South Africa - Headquarters: 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

Subsidiaries

Australia

Southern Implants Australia
T: +61-(0)-8-9466-2627
E: info@southernimplants.com.au

Spain and Portugal

Southern Implants Iberica
T: +34 935 053 507
E: info@southernimplants.es

United Kingdom and Ireland

Southern Implants UK
T: +44-20 8059 4490
E: info@southernimplants.co.uk

USA and Canada

Southern Implants North America Inc.
T: +1-561-472-0990
E: customercare@southernimplants.com

Description

Southern Implants trephine drills are described in table A. These devices attach to a handpiece of an implant motor unit, and have a latch compatible to ISO 1797. This is in order to connect the trephine drill to the handpiece of an implant motor unit. Trephine drills have laser marked lines at different levels in millimeters to indicate the depth of cutting. The laser marking on the shaft of the drill is to indicate the internal diameter of the device. The trephine drills have cutting edges on the end of the drill with a hollow tube in the middle of the drill. When the drill is rotated by a dental handpiece, the cutting edges remove a thin annulus of bone around an implant, separating the core of bone from the surrounding bone. This device is supplied sterile and for single patient use.

Intended use

The Southern Implants® Drills and Handpiece devices are intended to prepare the implant site/osteotomy and aid in the placement/removal of the endosseous implant and/or prosthetic device(s).

Specifically, the Southern Implants® Trephine Drills are intended to harvest autogenous bone and/or remove osseointegrated implants.

Indications for use

The Southern Implants® Drills and Handpiece Devices are indicated for the functional and aesthetic treatment of patients who are eligible for implant-retained prostheses.

Intended user

Maxillofacial Surgeons, General Dentists, Orthodontists, Periodontists, Prosthodontists, and other appropriately trained and/or experienced medical professionals.

Intended environment

The devices are intended to be used in a clinical environment such as an operating theatre or a dentist consultation room.

Intended patient population

The Southern Implants® Drills and Handpiece devices are intended to be used in patients subject to implant treatment.

Compatibility information

Southern Implants' implants trephine drills are as described in Table A. Details on site preparation sequences recommended by Southern Implants® can be found in the relevant product catalogues.

Table A – Trephine Drills

D-TRE-3	D-TRE-4	D-TRE-5	D-TRE-7
D-TRE-3 15 13 7 10 Ø3.6mm (internal) Ø4.3mm (outer)	D-TRE4 Ø4.1mm (internal) Ø4.8mm (outer)	D-TRE5 Ø5.1mm (internal) Ø5.8mm (outer)	D-TRE7 Ø7.0mm (internal) Ø7.7mm (outer)

Clinical performance

Southern Implants® Trepine Drills are designed to aid in implant site preparation by initiating and shaping a precisely profiled osteotomy for implant placement. Their specialised geometry, including appropriate lengths, diameters, and cutting features, ensures accurate and efficient site preparation.

Clinical benefits

Southern Implants® Trepine Drills play a vital role in implant-retained prosthetic treatment by enabling precise osteotomy preparation. While their clinical benefits are primarily tied to the overall success of the restorative treatment rather than their direct performance, their role in achieving accurate implant placement supports favourable outcomes.

The primary clinical benefits of this treatment approach include the effective restoration of the identified defect and enhanced patient satisfaction. Additionally, successful prosthetic rehabilitation can contribute to improved social interaction, psychosocial well-being, self-esteem, and confidence. Although these benefits are indirect and challenging to quantify, they are closely linked to positive treatment outcomes, such as optimal implant placement, functional restoration, and patient-reported satisfaction with both aesthetics and function.

Storage, cleaning and sterilisation

APPLICABLE FOR STERILE COMPONENTS:

The component is supplied sterile (sterilised by gamma irradiation) and may be intended for single-use prior to the expiration date (see packaging label). Sterility is assured unless the container or seal is damaged or opened. If packaging is damaged do not use the product and contact your Southern representative or return to Southern Implants®. The devices must be stored in a dry place at room temperature and not exposed to direct sunlight. Incorrect storage may influence device characteristics. Do not reuse components indicated for single-use only. Reusing these components may result in:

- damage to the surface or critical dimensions, which may result in performance and compatibility degradation.
- adds the risk of cross-patient infection and contamination if single-use items are reused.

Southern Implants® does not accept any responsibility for complications associated with reused single-use components.

Contraindications

The contraindications of all device groups used as part of the specific treatment or procedure apply. Therefore, the contraindications of the systems/medical devices utilized as part of implant surgery/therapy should be noted and the relevant documents consulted.

Contraindications specific to the Southern Implants® Drills and Handpiece Devices:

- Allergies or hypersensitivity to chemical ingredients in the following materials: Titanium, Aluminium, Vanadium, Stainless Steel.

Warnings and precautions

IMPORTANT NOTICE: THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING

- To ensure the safe and effective use of dental implants, new technologies/systems and the Metal Abutment devices, it is strongly recommended that specialised training be undertaken. This training should include hands-on methods to gain competency on the proper technique, the system's biomechanical requirements and radiographic evaluations required for the specific system.
- Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone.
- Use of the device with incompatible or non-corresponding devices can result in the poor performance or failure of the device.
- When handling devices intraorally, it is imperative that they are adequately secured to prevent aspiration, as aspiration of products may lead to infection or physical injury.
- The use of non-sterile items can lead to secondary infections of the tissue or transfer infectious diseases.

- Failure to adhere to appropriate cleaning, re-sterilization, and storage procedures as outlined in the Instructions for Use can result in device damage, secondary infections, or patient harm.
- Exceeding the number of recommended uses for reusable devices can result in device damage, secondary infection or patient harm.
- The use of blunt drills may cause damage to the bone, potentially compromising osseointegration.

It is crucial to emphasize that training should be undertaken by both new and experienced implant users prior to utilizing a new system or attempting a new treatment method.

Patient Selection and Preoperative Planning

A comprehensive patient selection process and meticulous preoperative planning are essential for successful implant treatment. This process should involve consultation among a multidisciplinary team, including well-trained surgeons, restorative dentists, and laboratory technicians.

Patient screening should include, at a minimum, a thorough medical and dental history, as well as visual and radiological inspections to assess the presence of adequate bone dimensions, the positioning of anatomical landmarks, the presence of unfavourable occlusal conditions, and the periodontal health status of the patient.

For successful implant treatment, it is important to:

1. Minimize trauma to the host tissue, thereby increasing the potential for successful osseointegration.
2. Accurately identify measurements relative to radiographic data, as failure to do so may lead to complications.
3. Be vigilant in avoiding damage to vital anatomical structures, such as nerves, veins, and arteries. Injury to these structures may result in serious complications, including ocular injury, nerve damage, and excessive bleeding.

The responsibility for proper patient selection, adequate training and experience in implant placement, and the provision of appropriate information required for informed consent rests with the practitioner. By combining thorough screening of prospective implant candidates with a practitioner possessing a high level of competence in the use of the system, the potential for complications and severe side effects can be significantly reduced.

High Risk Patients

Special care should be exercised when treating patients with local or systemic risk factors that may adversely affect the healing of bone and soft tissue or otherwise increase the severity of side effects, the risk of complications, and/or the likelihood of implant failure. Such factors include:

- poor oral hygiene
- history of smoking/vaping/tobacco use
- history of periodontal disease
- history of orofacial radiotherapy**
- bruxism and unfavourable jaw relations
- use of chronic medications that may delay healing or increase the risk of complications including, but not limited to, chronic steroid therapy, anti-coagulant therapy, TNF- α blockers, bisphosphonate and cyclosporin

*** The potential for implant failure and other complications increases when implants are placed in irradiated bone, as radiotherapy can lead to progressive fibrosis of blood vessels and soft tissue (i.e., osteoradionecrosis), resulting in diminished healing capacity. Contributing factors to this increased risk include the timing of implant placement in relation to radiation therapy, the proximity of radiation exposure to the implant site, and the radiation dosage at that site.*

Should the device not operate as intended, it must be reported to the manufacturer of the device. The contact information or the manufacturer of this device to report a change in performance is: sicomplaints@southernimplants.com.

Side effects

The clinical outcome of treatment is influenced by multiple factors. The following side effects and residual risks are associated with this device group or the surgical procedures involved in its use. These may require further treatment, revision surgery, or additional consultations with the relevant medical professional.

- Allergic or hypersensitivity reaction(s)
- Anaesthesia, paraesthesia, hyperesthesia, and hypoesthesia (transient or permanent)
- Anatomical landmark damage
- Bruising
- Dental injury during surgery
- Fat emboli formation
- Fistula formation
- Gingival injury
- Infection (acute and/or chronic)
- Localized inflammation
- Loss or damage to adjacent dentition
- Mild thermal injury to adjacent structures (e.g. burn to the lip)
- Nasal bleeding (Zygomatic Implant procedures)
- Nerve damage
- Pain, discomfort or tenderness
- Perforation of the labial or lingual plates
- Peri-implantitis, peri-mucositis or otherwise poor peri-implant soft tissue health
- Sinus Perforation
- Soft tissue irritation
- Surgical side effects such as pain, inflammation, bruising and mild bleeding
- Thermal injury to bone resulting in insufficient levels of osseointegration
- Wound dehiscence or poor healing

Additionally, the normal side effects associated with anaesthesia should also be expected. Note: these side effects and residual risks can vary in both severity and frequency.

Precaution: maintaining sterility protocol

APPLICABLE TO STERILE COMPONENTS

The devices are packed in a peel pouch or blister base with a “peel-back” lid. Labelling information is located on the bottom half of the pouch, inside the packet or on the surface of the peel-back lid. Sterility is assured unless the pouch is damaged or opened.

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established.

The contact information for the manufacturer of this device to report a serious incident is as follows:

sicomplaints@southernimplants.com.

Materials

Trephine Drills: Titanium alloy Ti-6Al-4V (90% titanium, 6% aluminium, 4% vanadium) according to ASTM F136 and ISO 5832-3

Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

Summary of Safety and Clinical Performance (SSCP)

As required by the European Medical Device Regulation (MDR; EU2017/745), a Summary of Safety and Clinical Performance (SSCP) is available for perusal with regard to Southern Implants® product ranges.

The relevant SSCP can be accessed at <https://ec.europa.eu/tools/eudamed>.

NOTE: the above website will be available upon the launch of the European Database on Medical Devices (EUDAMED).

Disclaimer of liability

This product is part of the Southern Implants® product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants® product range and take full responsibility for the correct indications and use of this product. Southern Implants® does not assume liability for damage due to incorrect use. Please note that some Southern Implants® products may not be cleared or released for sale in all markets.


Basic UDI

Product	Basic-UDI Number
Basic-UDI for Drills and Handpiece Devices	6009544038759C

Related literature and catalogues

CAT-2004 - Tri-Nex Implants Product Catalogue
 CAT-2005 - IT Implants Product Catalogue
 CAT-2010 - Osseointegrated Fixtures Catalogue
 CAT-2020 - External Hex Implants Product Catalogue
 CAT-2042 - Deep Conical Implants Product Catalogue
 CAT-2043 - Internal Hex Implants Product Catalogue
 CAT-2060 - PROVATA® Implants Product Catalogue
 CAT-2069 - INVERTA® Implants Product Catalogue
 CAT-2070 - Zygomatic Implants Product Catalogue

Symbols and warnings

* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

All rights reserved. Southern Implants®, the Southern Implants® logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants®. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale. It is the responsibility of the clinician to inspect the symbols that appear on the packaging of the product in use.

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