English	INSTRUCTIONS FOR USE: Southern Implants® Passive abutments
Español	INSTRUCCIONES DE USO: Pilares pasivos Southern Implants®
Italiano	ISTRUZIONI PER L'USO: Southern Implants® Monconi passivi
Français	MODE D'EMPLOI : Piliers passifs Southern Implants®
Deutsch	GEBRAUCHSANWEISUNG: Southern Implants® Passive Abutments
Português	INSTRUÇÕES DE UTILIZAÇÃO: Pilares passivos Southern Implants®

ARCHIVED

South Africa - Headquarters: 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

Subsidiaries

Australia

Southern Implants Australia
T: +61-(0)-8-9466-2627
E: info@southernimplants.com.au

Spain and Portugal

Southern Implants Iberica
T: +34 935 053 507
E: info@southernimplants.es

United Kingdom and Ireland

Southern Implants UK
T: +44-20-8899-6845 / 6 / 7
E: info@southernimplants.co.uk

USA and Canada

Southern Implants North America Inc.
T: +1-561-472-0990
E: customer@southernimplants.com


Intended use

Southern Implants® dental abutments are intended to be used in the Maxilla or Mandible for supporting a prosthesis on endosseous implants in order to restore chewing function for the patient.

Description

The Southern Implants prosthetic range consist of abutments or cylinders for prosthetic rehabilitation on the Southern range of implants. These are pre-manufactured and can be attached direct or indirect to the implant as an aid in prosthetic rehabilitation. These instructions apply to Passive abutments with burn out sleeves. These components are supplied non-sterile.

Fig.1 - Passive abutment components.

Indications for use

Southern Implants dental implants are intended for both one and two-stage surgical procedures in the following situations and with the following clinical protocols:

- replacing single and multiple missing teeth in the mandible and maxilla.
- immediate placement in extraction sites and in situations with a partially or completely healed alveolar ridge.
- immediate loading in all indications, except in single tooth situations on implants shorter than 8mm or in soft bone (type IV) where implant stability may be difficult to obtain and immediate loading may not be appropriate.

Intended user

Dental Technicians, Maxillo-facial Surgeons, General Dentists, Orthodontists, Periodontists, Prosthodontists and other appropriately trained and experienced implant users.

Intended environment

The Passive abutments are intended to be used in a dental laboratory for making of the restoration and in a clinical environment such as an operating theatre or a dentist consultation room for final clinical procedures.

Intended patient population

This device is used in the dental restoration of partially or fully edentulous patients in the upper or lower jaw. Restorations may comprise of single teeth, partial or full bridges and may be fixed or removable.

Compatibility information

Southern Implants' implants should be restored with Southern components. In the Southern Implants range there are 5 implant

connections. The implant code and connection type can be identified by specific abbreviations in the product codes. Range identifiers are summarised in Table A.

Table A

Implant connection type	Compatible device
External Hex (EX)	Parts labelled SP16, SBN16, SB16, SBA16, SBBB16 and SMAX9h for engaging items
	Parts labelled SP17, SBN17, SB17, SBA17, SBBB17 and SMAX9nh for non-engaging items
TRI-NEX (EL) (Lobe)	Parts labelled PA-EL-(Ø) for engaging items
	Parts labelled PA-NL-(Ø) for non-engaging items
Deep Conical (DC)	Parts labelled PA-DC-(Ø) for engaging items
	Parts labelled PA-NDC-(Ø) for non-engaging items
Internal Hex (M-Series and PROVATA®) (M)(Z)	Parts labelled PA-3EM-S (used with PROVATA® Ø3.3 mm platform) for engaging items
	Parts labelled PA-3NM-S (used with PROVATA® Ø3.3 mm platform) for non-engaging items
	Parts labelled PA-EM-S (used with M-Series Ø3.75, 4.2 and 5.0 mm and PROVATA® Ø4.0, 5.0 and 6.0 mm platforms) for engaging items
	Parts labelled PA-NM-S (used with M-Series Ø3.75, 4.2 and 5.0 mm and PROVATA® Ø4.0, 5.0 and 6.0 mm platforms) for non-engaging items
	Parts labelled PA-EZ (used with PROVATA® Ø7.0, 8.0 and 9.0 mm platforms) for engaging items
	Parts labelled PA-NZ (used with PROVATA® Ø7.0, 8.0 and 9.0 mm platforms) for non-engaging items
IT (ITS) (ITS6) - Internal Octagon	Parts labelled ITS-PA (used with Ø4.8 mm platforms) for engaging item
	Parts labelled ITS-PA-ne (used with Ø4.8 mm platforms) for non-engaging items
	Parts labelled ITS6-PA (used with Ø6.5 mm platforms) for engaging items
	Parts labelled ITS6-PA-ne (used with Ø6.5 mm platforms) for non-engaging items
Abutment level (MC)	Parts labelled PA-MC-48 (used with Ø4.8 mm abutment platforms) for non-engaging items
	Parts labelled PA-MC-60 (used with Ø6.0 mm abutment platforms) for non-engaging items

Surgical procedures

1. Use a model with removable soft tissue mask, allowing easy access to the analogues for lab procedures.
2. Finger-tighten the titanium interfacial component and plastic sleeve to the analogue with the laboratory screw. Use the corresponding driver as indicated in Table D. Do not over-tighten to avoid distortion of the plastic. The waxing sleeve can be cut back or added to as needed. The wax-up is completed and sprued on the model.

Investing and casting

1. The retaining screw must be removed to allow the wax-up with plastic cylinders to be lifted from the model, leaving behind the loose titanium interfacial component.
2. Standard procedures are used for investing and casting. An appropriate casting alloy must be chosen, depending on whether a ceramic veneered bridge or cast bar is being manufactured. Alloys that are commonly used are: Degunorm, Argipal, Begopal 300, Begocer-G, Pors-on 4, Degudent G etc.
3. For complete burn-out: the plastic cylinder requires an oven temperature of 820°C for at least 45 minutes.
4. Devest ultrasonically as opposed to blasting with sand or glass beads. This helps preserve the sharp edges and fitting surfaces of the casting.

Refining the screw seat

1. The screw seat is the internal ledge in the casting where the head of the screw will seat. The cast surface of the screw seat will be rough due to the casting process and must therefore be refined using special hand-held reamers as indicated in Table B. The correct diameter of reamer must be chosen. This is an

important step to ensure proper seating and tightening of the final prosthetic screw.

NOTE: do not sandblast the titanium interfacial component.

Table B

Reamer diameter	Reamer Code	Use with	Screw head diameter
2.2	LT18-2.2	TS-DC3-MC	2.00
		TS-DC4-MC	
		TS-DC3-14	2.10
		TS-P-16	2.15
		TSS4	2.20
		GSH4	
		PA-DC3-14T	
		TS-Z-MC-18	
		TS-Z-MC-16	
		PA-M-16T	
2.3	LT18-2.3	TSH1	2.25
		TSU1	
		TSS1	
		GSH1	
		GSU1	
		GSS1	
		TU-MCL-18C	
		GU-MCL-18	
		TS-Z-18	2.30
		GS-Q-18	
		TS-DC4-16	
		TS-Z-16	
2.4	LT18-2.4	TS-DC5-20	2.35
		TSHZ3	2.40
		TSUZ3	
		TSSZ3	
		TSU3	
		TSS3	
		GSUZ3	
		GSQZ3	
		GSSZ3	
		GSS3	
		GSU3	
		GSQ3	
		TS-IT-TIB	
		PA-DC4-16T	
2.6	LT18-2.6	TU-MCL-20C	2.45
		GU-MCL-20	2.50
		TSUZ9	
		TSU9	
		GSUZ9	
		GSU9	
		TS-L-18	
		TS-L-18C	
		GS-L-18	
		TS-L-20	
		TS-L-20C	
		GS-L-20	
		PA-M-18T	2.60
		PA-Q-18G	

2.8	LT18-2.8	TSHZ2	2.70
		TSUZ2	
		TSU2	
		TSSZ2	
		TSS2	
		GSUZ2	
		GSQZ2	
		GSSZ2	
		GSIT2	
		GSS2	
		GSU2	
		GSQ2	
		TSIT2	
		TS-IT-PA	
		PA-L-18	
2.9	LT18-2.9	PA-L-18G	2.80
		PA-L-20	
		PA-L-20G	
		PA-DC5-20T	

Luting procedures for Passive Abutments

The luting procedure and concept stays the same when using passive abutments on compact conical level, multiple unit restorations and single unit restorations. The following steps are applicable to all.

- It is important to ultrasonically clean or steam clean the following before the luting procedure starts:
 - the titanium interfacial components.
 - the luting screws.
 - the fitting surfaces of the prosthesis.
 - clean the analogues in the model by brushing with soap and water or steam clean to remove any debris which may interfere with perfect seating of the interfacial components.
- Finger-tighten (with the PEEK luting screw) the titanium interfacial component to the laboratory analogue. Use the corresponding driver as indicated in Table C.
- The sleeve can easily be fitted and removed from the titanium interfacial component without the need to remove or replace the PEEK luting screw. This is due to the PEEK luting screw attaching the titanium interfacial component to the analogue and not the structure to the analogue. The PEEK luting screw can be removed and a laboratory screw can be used to attach the prosthesis to the master model. This will ensure that the prosthesis does not move during articulation, opposite set up, or any other laboratory procedures.
- Clean the restoration in an ultrasonic unit for about 1 minute. Dry with oil-free air. Apply a bonding agent like Monobond Plus (Ivoclar Vivadent) to the cleaned surfaces of the titanium interfacial component and restoration using a small brush. Allow the Monobond Plus to react for 60 seconds and disperse with compressed air.
- Lute the titanium interfacial abutment to the structure by applying a self-cure resin cement or dual cure resin cement (e.g. Rely X by 3M) to the surface of all of the titanium interfacial components.

NOTE: refrigeration of self-cure resin cements will usually lengthen working time for ease of use on multi-unit structures. Limit the amount of resin cement being applied to the angle between the horizontal plane and vertical plane of the titanium interfacial component. This will avoid excess cement extruding upwards through the screw hole in the structure and inadvertently locking the luting screw into the cement. Definitely avoid placing any cement in the area immediately around the head of the luting screw.

NOTE: in the event that cement locks the luting screw in position, a Ø1.5 mm round burr is rotated into the hex of the screw. This separates the screw head from the shaft and frees the prosthesis. Take care not to damage the components. The prosthesis can

then be removed.

6. Fit the prosthesis over the titanium interfacial component and settle the prosthesis firmly into place with finger pressure to extrude excess cement. Arch structures can be left seated under their own weight to allow cement to harden.
7. Light cure for 60 seconds.
8. Un-mount the restoration and remove excess cement using a sharp blade, probe or hand scaler instrument to make polishing easier.
9. Once resin cement has hardened, remove all luting screws and then remove any prosthetic retaining screws so that the prosthesis can be lifted from the model.
10. Attach polishing protectors to each of the fitting surfaces of the cemented titanium interfacial component. Polish the remaining cement line using a fine-edged lens-shaped rubber wheel. The cement line is often not of constant thickness. This variation is indicative of the extent of misfit which existed and has now been corrected by the cement space of the Passive Abutment.

Clinical procedures

The clinician receives the restoration from the laboratory.

1. Remove the healing abutment or temporary restoration.
2. Clean, disinfect and sterilise the restoration.
3. Insert the restoration into the patient's mouth.
4. Position the restoration on the implant making sure that the retentive elements of the implant/abutment connections are properly aligned.

Table C - Luting Screw

Driver type	External Hex	DC	Tri-Nex	Internal Hex	IT	Compact Conical screw
0.90 mm		DC3 only				✓
1.22 mm / 1.27 mm Universal driver	✓	✓		✓	✓	
1.22 mm hex driver	✓	✓			✓	
1.27 mm hex driver				✓		
Unigrip driver			✓			
Quad driver						
Blade driver						
Torx driver						

5. Fix the abutment to the implant with the correct screw using the applicable driver (Table D), torque the screw down to the value indicated in Table E.

Table D - Retaining Screw

Driver type	External Hex	DC	Tri-Nex	Internal Hex	IT	Compact Conical screw
1.22 mm / 1.27 mm Universal driver	✓	✓		✓		✓
1.22 mm hex driver	✓	✓				✓
1.27 mm hex driver	✓		✓	✓		✓
Unigrip driver	✓					✓
Quad driver	✓			Gold screws only		✓
Blade driver	✓					✓
Torx driver	✓				✓	✓

Table E

Direct to Implant	Torque
Ext-Hex	
Ø3.0 mm	25-32 Ncm
Ø3.25, 4.0, 5.0, 6.0, 7.0, 8.0 and 9.0 mm	32-40 Ncm
Tri-Nex	
Ø3.0mm	30 Ncm
Ø4.3, 5.0, 6.0, 7.0, 8.0 and 9.0 mm	32-40 Ncm
DC	
Ø3.0 mm	15 Ncm
Ø3.5 and 4.0 mm	20 Ncm
Ø5.0 mm	25-32 Ncm
Internal Hex (M-Series & PROVATA®)	
Ø3.75, 4.2 and 5.0 mm M-Series	32 Ncm
Ø3.3, 4.0, 5.0, 6.0, 7.0, 8.0 and 9.0 mm PROVATA® Implant	32 Ncm
IT Octagon	
Ø3.3, 4.1, 4.9, 4.0, 5.0, 6.0, 7.0, 8.0 and 9.0 mm IT Implant	32-40 Ncm
Abutment level	All Passive Abutments on Compact Conical Abutments: 10-15 Ncm

6. Verify the correct seating of the restoration using radiographic image.
7. Do not exceed the recommended torque value as this may result in failure of the screw, abutment or implant. Do not torque less than the recommended value, this may result in loosening of the abutment.
8. Close the screw access hole.
9. Cement the final prosthesis.

Clinical benefits

Through this procedure patients can expect to have their missing teeth replaced and/or crowns restored

Before surgery

All components, instruments and tooling used during the clinical or laboratory procedure must be maintained in good condition and care must be taken that instrumentation does not damage implants or other components.

During surgery

Take care that parts are not swallowed or aspirated during any of the procedures and apply the correct tightening torque to abutments and abutment screws.

Post-surgery

Regular patient follow-up and proper oral hygiene must be achieved to ensure favourable long-term results.

Storage, cleaning and sterilisation

The implants, healing screws and healing abutments are supplied sterile (sterilised by gamma irradiation) and intended for single-use prior to the expiration date (see packaging label). Sterility is assured unless the container or seal is damaged or opened. If packaging is damaged do not use the product and contact your Southern representative or return to Southern Implants. Do not reuse implants, cover screws, temporary abutments and abutments. Re-using these components may result in:

- damage on the surface or critical dimensions, which may result in performance and compatibility degradation.
- adds the risk of cross-patient infection and contamination if single-use items are reused.

Southern Implants does not accept any responsibility for complications associated with reused components.

Cleaning and disinfection

An implant restoration is a single or multiple-tooth implant crown, bridge or substructure, attached to a Southern Implants abutment or multiple abutments. Before intraoral use the final restoration needs to be cleaned and disinfected, as per restorative material manufacturer's instructions.

Sterilisation

Southern Implants recommends the following procedure to sterilise the restoration prior to use:

Methods to sterilise the restoration and abutment screw

1. Pre-vacuum sterilisation method: steam sterilise the abutments at 132°C (270°F) at 180-220kPa for 4 minutes. Dry for at least 20 minutes in the chamber. Only an approved wrap or pouch for steam sterilisation must be used.
2. Pre-vacuum sterilisation method: wrapped, steam sterilise at 135°C (275°F) for 3 minutes. Dry for 20 minutes in the chamber. Use a wrap or pouch that is cleared for the indicated steam sterilisation cycle.

NOTE: users in the USA must ensure that the steriliser, wrap or pouch, and all steriliser accessories are cleared by the FDA, for the intended sterilisation cycle.

Contraindications

Do not use in patients:

- who are medically unfit for dental implant procedures.
- where adequate numbers of implants could not be placed to achieve full functional support of the prosthesis.
- who are allergic or have hypersensitivity to pure titanium or titanium alloy (Ti-6Al-4V), gold, palladium, platinum, iridium or polyetheretherketone (PEEK).
- who are under the age of 18, have poor bone quality, blood disorders, infected implant site, vascular impairment, uncontrolled diabetes, drug or alcohol abuse, chronic high dose steroid therapy, anti-coagulant therapy, metabolic bone disease, radiotherapy treatment or sinus pathology.

Warnings

THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING.

- For the safe and effective use of dental implants it is suggested that specialised training be undertaken, including hands-on training to learn proper technique, biomechanical requirements and radiographic evaluations.
- Responsibility for proper patient selection, adequate training, experience in the placement of implants and providing appropriate information for informed consent rests with the practitioner. Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone.

Cautions

New and experienced implant users should do training before using a new system or attempting to do a new treatment method. Take special care when treating patients who have local or systemic factors that could affect the healing of the bone and soft tissue (i.e. poor oral hygiene, uncontrolled diabetes, are on steroid therapy, smokers, infection in the nearby bone and patients who had oro-facial radiotherapy). Thorough screening of prospective implant candidates must be performed including:

- a comprehensive medical and dental history.
- visual and radiological inspection to determine adequate bone dimensions, anatomical landmarks, occlusal conditions and periodontal health.
- bruxism and unfavourable jaw relations must be taken into account.
- proper pre-operative planning with a good team approach between well trained surgeons, restorative dentists and lab technicians is essential for successful implant treatment.
- minimizing the trauma to the host tissue increases the potential for successful osseointegration.
- electro-surgery should not be attempted around metal implants, as they are conductive.

Single use devices

Do not reuse devices indicated for single use (use the device prior to the expiration date).

Side effects

Potential side effects and temporary symptoms: pain, swelling, phonetic difficulties, gingival inflammation. More persistent symptoms: the risks and complications with implants include, but are not limited to: (1)

allergic reaction(s) to implant and/or abutment material; (2) breakage of the implant and/or abutment; (3) loosening of the abutment screw and/or retaining screw; (4) infection requiring revision of the dental implant; (5) nerve damage that could cause permanent weakness, numbness, or pain; (6) histologic responses possibly involving macrophages and/or fibroblasts; (7) formation of fat emboli; (8) loosening of the implant requiring revision surgery; (9) perforation of the maxillary sinus; (10) perforation of the labial and lingual plates and (11) bone loss possibly resulting in revision or removal.

Breakage

Implant and abutment fractures can occur when applied loads exceed the tensile or compressive strength of the material. Potential overloading conditions may result from: deficiencies in implant numbers, lengths and/or diameters to adequately support a restoration, excessive cantilever length, incomplete abutment seating, abutment angles greater than 30 degrees, occlusal interferences causing excessive lateral forces, patient parafunction (e.g. bruxing, clenching), loss or changes in dentition or functionality, inadequate prosthesis fit and physical trauma. Additional treatment may be necessary when any of the above conditions are present to reduce the possibility of hardware complications or failure.

Performance requirements

It is the responsibility of the clinician to instruct the patient on all appropriate contraindications, side effects and precautions as well as the need to seek the services of a trained dental professional if there are any changes in the performance of the implant (e.g. looseness of the prosthesis, infection or exudate around the implant, pain or any other unusual symptoms that the patient has not been told to expect).

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established. The contact information for the manufacturer of this device to report a serious incident is as follows: sicomplaints@southernimplants.com

Materials

Passive Abutment:	Titanium grade 2, 3, 4 or 5
Plastic Sleeve:	Polyoxymethylene (POM)
Luting screw:	Medical grade PEEK (polyetheretherketone)
Abutment screws:	Titanium alloy Ti-90%, Al-6%, V-4% Gold Alloy Au-61%, Ag-16.5%, Pt-13.5%, Cu-9%

Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

Disclaimer of liability

This product is part of the Southern Implants product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants product range and take full responsibility for the correct indications and use of this product. Southern Implants does not assume liability for damage due to incorrect use. Please note that some Southern Implants products may not be cleared or released for sale in all markets.

Healing

The healing time required for osseointegration depends on the individual and treatment protocol. It is the responsibility of the practitioner to decide when the implant can be restored. Good primary stability will govern if immediate loading can be done.

Implant care and maintenance

Potential implant patients should establish an adequate oral hygiene regime prior to implant therapy. Proper post-operative, oral hygiene and implant maintenance instructions must be discussed with the patient, as this will determine the longevity and health of the implants. The patient should maintain regular prophylaxis and evaluation appointments.

MR conditional

Non-clinical testing has demonstrated that the Southern Implants dental implants, metallic abutments and prosthetic screws are MR conditional. A patient with these devices can be safely scanned in a MR system meeting the following conditions:

- Static magnetic field of 1.5 Tesla and 3.0 Tesla only.
- Maximum spatial gradient magnetic field of 4000 Gauss/cm (40 T/m).
- Maximum MR system reported, head specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode) or whole body averaged specific absorption rate (wbSAR) of 1 W/kg.

Under the scan conditions defined above, the Southern Implants dental implants, abutments and prosthetic screws are expected to produce a maximum temperature rise of 5.8°C after 15 minutes of continuous scanning. In the non-clinical testing, the image artifact caused by the device extends approximately 20 mm from the Southern Implants dental implants, abutments and prosthetic screws, when imaged with a gradient echo pulse sequence and a 3.0 Tesla MRI system. Removable restorations should be taken out prior to scanning, as is done for watches, jewellery etc. Should there be no MR symbol on the product label, please note that this device has not been evaluated for safety and compatibility in the MR environment. This device has not been tested for heating or migration in the MR environment.

Basic UDI

Product	Basic-UDI Number
Basic-UDI for Metal Abutments	600954403872

Related literature and catalogues

CAT-2004 - Tri-Nex Implants Product Catalogue
 CAT-2005 - IT Implants Product Catalogue
 CAT-2020 - External Hex Implants Product Catalogue
 CAT-2042 - Deep Conical Implants Product Catalogue
 CAT-2043 - Internal Hex Implants Product Catalogue
 CAT-2060 - PROVATA® Implants Product Catalogue
 CAT-2069 - INVERTA® Implants Product Catalogue
 CAT-2070 - Zygomax Implants Product Catalogue

Symbols and Warnings

 Manufacturer: Southern Implants 1 Albert Rd, P.O Box 605 IRENE, 0062, South Africa. Tel: +27 12 667 1046	           	Prescription device* Non-sterile Consult instruction for use Use by date (mm-yy) Do not reuse Batch code Do not use if package is damaged Medical Device Authorize Representative in the European Community Catalog number Date of Manufacture
* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist. Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law. All rights reserved. Southern Implants®, the Southern Implants logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale.		

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