**SOUTHERNIMPLANTS®**

Innovative Treatment Solutions

|                  |   |
|------------------|---|
| <b>English</b>   | <b>INSTRUCTIONS FOR USE: Southern Implants® Temporary Titanium Abutments</b>          |
| <b>Español</b>   | <b>INSTRUCCIONES DE USO: Pilares temporales de titanio Southern Implants®</b>         |
| <b>Italiano</b>  | <b>ISTRUZIONI PER L'USO: Southern Implants® Monconi provvisori in titanio</b>         |
| <b>Français</b>  | <b>MODE D'EMPLOI : Southern Implants® Piliers temporaires en titane</b>               |
| <b>Deutsch</b>   | <b>GEBRAUCHSANWEISUNG: Temporäre Titan-Abutments von Southern Implants®</b>           |
| <b>Português</b> | <b>INSTRUÇÕES DE UTILIZAÇÃO: Pilares provisórios de titânio da Southern Implants®</b> |

ARCHIVED

**South Africa - Headquarters:** 1 Albert Road, Irene, 0062, RSA  
T: +27-12-667-1046 | E: [info@southernimplants.com](mailto:info@southernimplants.com)

EC REP

**Southern Implants Europe AB:** Holmgatan 30, S-791 71 Falun, Sweden  
T: +46 23 13300 | E: [ecrep@southernimplants.com](mailto:ecrep@southernimplants.com)

#### Subsidiaries

##### Australia

Southern Implants Australia  
T: +61-(0)-8-9466-2627  
E: [info@southernimplants.com.au](mailto:info@southernimplants.com.au)

##### Spain and Portugal

Southern Implants Iberica  
T: +34 935 053 507  
E: [info@southernimplants.es](mailto:info@southernimplants.es)

##### United Kingdom and Ireland

Southern Implants UK  
T: +44-20-8899-6845 / 6 / 7  
E: [info@southernimplants.co.uk](mailto:info@southernimplants.co.uk)

##### USA and Canada

Southern Implants North America Inc.  
T: +1-561-472-0990  
E: [customer care@southernimplants.com](mailto:customer care@southernimplants.com)

**Intended use**

Southern Implants® dental implant abutments are intended to be used in the Maxilla or Mandible for supporting a prosthesis on endosseous implants in order to restore chewing function for the patient.

**Description**

These are pre-manufactured dental implant abutments that can either be connected direct to an endosseous implant or connect the prosthesis to a compact conical abutment for use as an aid in temporary prosthetic rehabilitation. Refer to individual product catalogues for product characteristics and compatible accessories. The Titanium Abutments are provided sterile, however, it will no longer be sterile after modification.

**Indications for use**

Southern Implants dental implants are intended for both one and two-stage surgical procedures in the following situations and with the following clinical protocols:

- replacing single and multiple missing teeth in the mandible and maxilla.
- immediate placement in extraction sites and in situations with a partially or completely healed alveolar ridge.
- immediate loading in all indications, except in single tooth situations on implants shorter than 8 mm or in soft bone (type IV) where implant stability may be difficult to obtain and immediate loading may not be appropriate.

**Intended user**

Dental Technicians, Maxillo-facial Surgeons, General Dentists, Orthodontists, Periodontists, Prosthodontists and other appropriately trained and experienced implant users.

**Intended environment**

The devices are intended to be used in a clinical environment such as an operating theater or a dentist consultation room and in a dental laboratory for making of the restoration.

**Intended patient population**

This device is used in the dental restoration of partially or fully edentulous patients in the upper or lower jaw. Restorations may comprise of partial or full bridges, multi-unit cases and may be fixed or removable.

**Compatibility information**

Southern Implants' implants should be restored with Southern Implants components. In the Southern Implants range there are 5 implant connections. The implant code and connection type can be identified by specific abbreviations in the product codes. Range identifiers are summarised in table A.

**Table A**

| Implant connection type | Compatible device  |
|-------------------------|--|
| External Hex (TC)       | Parts labelled TCP1h, TCBN1h/5h, TCB1h/5h, TCBA1h/5h, TCBBB1h/5h and TCMAx9-1h for engaging items                    |
|                         | Parts labelled TCP1nh, TCBN1nh/5nh, TCB1nh/5nh, TC9, TCBA1nh/5nh, TCBBB1nh/5nh and TCMAx9-1nh for non-engaging items |
| Tri-Nex (EL) (Lobe)     | Parts labelled TC-EL-(Ø) for engaging items  |
|                         | Parts labelled TC-NL-(Ø) for non-engaging items  |
| Deep Conical (DC)       | Parts labelled TC-DC-(Ø) for engaging items  |
|                         | Parts labelled TC-NDC-(Ø) for non-engaging items   |
| Internal Hex (M)        | Parts labelled TC-M (used with Ø3.75, 4.20 and 5.00 mm platforms) for engaging items                                 |
|                         | Parts labelled TC-NM (used with Ø3.75, 4.20 and 5.00 mm platforms) for non-engaging items                            |

|                                |  |
|--------------------------------|--|
| Internal Hex (PROVATA®) (M)(Z) | Parts labelled TC-3M (used with Ø3.3 mm platform) for engaging items                       |
|                                | Parts labelled TC-3NM (used with Ø3.3 mm platform) for engaging items                      |
|                                | Parts labelled TC-M (used with Ø4.0, 5.0 and 6.0 mm platforms) for engaging items          |
|                                | Parts labelled TC-NM (used with Ø4.0, 5.0 and 6.0 mm platforms) for non-engaging items     |
|                                | Parts labelled TC-EZ-(*) (used with Ø7.0, 8.0 and 9.0 mm platforms) for engaging items     |
|                                | Parts labelled TC-NZ-(*) (used with Ø7.0, 8.0 and 9.0 mm platforms) for non-engaging items |
| IT (ITS) (ITS6) - Octagon      | Parts labelled ITS-TC1 (used with Ø4.8 mm platforms) for engaging items                    |
|                                | Parts labelled ITS-TC1ne (used with Ø4.8 mm platforms) for non-engaging items              |
|                                | Parts labelled ITS6-TC1 (used with Ø6.5 mm platforms), for engaging items                  |
|                                | Parts labelled ITS6-TC1ne (used with Ø6.5 mm platforms) for non-engaging items             |
| Abutment level                 | Parts labelled TMC1/5 and TMCSL (used with Ø4.8 mm abutment platforms) non-engaging        |
|                                | Parts labelled TMCW1/5 (used with Ø6.0 mm abutment platforms) non-engaging                 |

(\*) is indicative of various lengths available.

**Chairside procedure** (making a temporary restoration)

**NOTE:** modification of Titanium abutments can be done with a carbide burr or disk. It is recommended to do this extraorally and with copious irrigation during cutting.

1. Determine the cuff height that will be most suitable for the restoration. Titanium abutments are available in 1 mm and 5mm collar heights. DC (Deep Conical) and IT (octagon) are only available in 1 mm collar heights. PROVATA® (Internal Hex implant Promax ranges) are available in 1 mm and 3 mm collar heights.
2. Connect the abutment to the implant and modify the abutment to the correct occlusal height, with no additional angular correction. Modification of the abutment must be done with copious amounts of irrigation intra-orally (extraoral trimming of the abutment is the preferred recommendation). For single-unit use, do not reduce the post below a minimum height of 4 mm.
3. With a 5 mm Titanium abutment the collar can be trimmed to follow the contours of the soft tissue.
4. Close the screw channel hole in a way that will ensure the prosthetic screw can be retrieved.
5. Make a temporary restoration by using a pre-formed stent and suitable temporary material.
6. Unscrew the temporary prosthesis.
7. Make final adjustments.
8. Clean and disinfect the restoration as applicable per the restorative material manufacturer's instructions.
9. Attach the Temporary Abutment to the endosseous implant or compact conical abutment with the compatible prosthetic screw. Tighten the restoration using a manual torque wrench to the torque value specified in Table C.
10. Close screw access hole.
11. Cement final prosthesis if applicable.

**Laboratory procedures** (making a temporary restoration)

1. The laboratory receives the impression either implant level or abutment level.
2. The corresponding laboratory analogue is connected to the impression coping, fabricate a working model with removable gingival mask or soft tissue material.
3. The same steps as for clinical procedures will apply.

**Clinical procedures** (placing the temporary restoration)

The clinician receives the restoration from the laboratory.

4. Remove the healing abutment.
5. Clean, disinfect and sterilise the restoration as described.
6. Insert the restoration into the patient's mouth.
7. Position the restoration on the implant/abutment making sure that the retentive elements of the implant/abutment connections are

properly aligned.

**Table B**

| Driver type                        | External Hex | DC | Tri-Nex | Internal Hex     | IT | Compact Conical screw |
|------------------------------------|--------------|----|---------|------------------|----|-----------------------|
| 1.22 mm / 1.27 mm Universal driver | ✓            | ✓  |         | ✓                |    | ✓                     |
| 1.22 mm hex driver                 | ✓            | ✓  |         |                  |    | ✓                     |
| 1.27 mm hex driver                 |              |    |         | ✓                |    |                       |
| Unigrip driver                     | ✓            |    | ✓       |                  |    | ✓                     |
| Quad driver                        | ✓            |    |         | Gold screws only |    |                       |
| Blade driver                       | ✓            |    |         |                  |    | ✓                     |
| Torx driver                        |              |    |         |                  | ✓  |                       |

- Fix the abutment to the implant/abutment with the correct screw using applicable driver (Table B). Torque the screw down to the value indicated in Table C.

**Table C**

| Direct to Implant   | Torque   |
|---|--|
| <b>Ext-Hex</b>  |  |
| Ø3.0 mm   | 25-32 Ncm  |
| Ø3.25, 4.0, 5.0, 6.0, 7.0, 8.0 and 9.0 mm                     | 32-40 Ncm  |
| <b>Tri-Nex</b>  |  |
| Ø3.5mm  | 32 Ncm   |
| Ø4.3, 5.0, 6.0, 7.0, 8.0 and 9.0 mm                           | 32-40 Ncm  |
| <b>DC</b>   |  |
| Ø3.0 mm   | 15 Ncm   |
| Ø3.5, 4.0 mm  | 20 Ncm   |
| Ø5.0 mm   | 25-32 Ncm  |
| <b>Internal Hex (M-Series and PROVATA®)</b>                   |  |
| Ø3.75, 4.2, 5.0 mm M-Series                                   | 32 Ncm   |
| Ø3.3, 4.0, 5.0, 6.0, 7.0, 8.0 and 9.0 mm PROVATA® Implant     | 32 Ncm   |
| <b>IT Octagon</b>   |  |
| Ø3.3, 4.1, 4.9, 4.0, 5.0, 6.0, 7.0, 8.0 and 9.0 mm IT Implant | 32-40 Ncm  |
| <b>Abutment level</b>   | All Titanium cylinders on Compact Conical Abutments: 10-15 Ncm |

- Verify the correct seating of the restoration using radiographic image.
- Do not exceed the recommended torque value as this may result in failure of the screw, abutment or implant. Do not torque less than the recommended value, this may result in loosening of the abutment that can lead to abutment or implant failure.
- Close the screw access hole.
- Cement the temporary prosthesis if applicable.

#### Clinical benefits

Through this procedure patients can expect to have their missing teeth replaced and/or crowns restored.

#### Before surgery

All components, instruments and tooling used during the clinical or laboratory procedure must be maintained in good condition and care must be taken that instrumentation does not damage implants or other components.

#### During surgery

Care must be taken that parts are not swallowed during any of the procedures, a rubber-dam application is recommended when appropriate. Care must be taken to apply the correct tightening torque of abutments

and abutment screws.

#### Post-surgery

Regular patient follow-up and proper oral hygiene must be achieved to ensure favourable long-term results.

#### Storage, cleaning and sterilisation

The implants, cover screws and healing abutments are supplied sterile (sterilised by gamma irradiation) and intended for single-use prior to the expiration date (see packaging label). Sterility is assured unless the container or seal is damaged or opened. If packaging is damaged do not use the product and contact your Southern representative or return to Southern Implants. Do not reuse implants, cover screws, temporary abutments and abutments. Re-using these components may result in:

- damage on the surface or critical dimensions, which may result in performance and compatibility degradation.
- adds the risk of cross-patient infection and contamination if single-use items are reused.

Southern Implants does not accept any responsibility for complications associated with reused components.

#### Cleaning and disinfection

An implant restoration is a single or multiple-tooth implant crown, bridge or substructure, attached to a Southern Implants abutment or multiple abutments. Before intraoral use the final restoration needs to be cleaned and disinfected, as per restorative material manufacturer's instructions.

#### Sterilisation

Southern Implants recommends the following procedure to sterilise the restoration prior to use:

Methods to sterilise the restoration and abutment screw

- Pre-vacuum sterilisation method: steam sterilise the abutments at 132°C (270°F) at 180-220kPa for 4 minutes. Dry for at least 20 minutes in the chamber. Only an approved wrap or pouch for steam sterilisation must be used.
- Pre-vacuum sterilisation method: wrapped, steam sterilise at 135°C (275°F) for 3 minutes. Dry for 20 minutes in the chamber. Use a wrap or pouch that is cleared for the indicated steam sterilisation cycle.

**NOTE:** users in the USA must ensure that the steriliser, wrap or pouch and all steriliser accessories are cleared by the FDA, for the intended sterilisation cycle.

#### Contraindications

Do not use in patients:

- who are medically unfit for dental implant procedures.
- where adequate numbers of implants could not be placed to achieve full functional support of the prosthesis.
- who are allergic or have hypersensitivity to pure titanium or titanium alloy (Ti-6Al-4V), gold, palladium, platinum or iridium.
- who are under the age of 18, have poor bone quality, blood disorders, infected implant site, vascular impairment, uncontrolled diabetes, drug or alcohol abuse, chronic high dose steroid therapy, anti-coagulant therapy, metabolic bone disease, radiotherapy treatment and sinus pathology.

#### Warnings

THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING.

- For the safe and effective use of dental implants it is suggested that specialised training be undertaken, including hands-on training to learn proper technique, biomechanical requirements and radiographic evaluations.
- Responsibility for proper patient selection, adequate training, experience in the placement of implants and providing appropriate information for informed consent rests with the practitioner. Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone.

#### Cautions

New and experienced implant users should do training before using a new system or attempting to do a new treatment method. Take special care when treating patients who have local or systemic factors that could

affect the healing of the bone and soft tissue (i.e. poor oral hygiene, uncontrolled diabetes, are on steroid therapy, smokers, infection in the nearby bone and patients who had oro-facial radiotherapy). Thorough screening of prospective implant candidates must be performed including:

- a comprehensive medical and dental history.
- visual and radiological inspection to determine adequate bone dimensions, anatomical landmarks, occlusal conditions and periodontal health.
- bruxism and unfavourable jaw relations must be taken into account.
- proper pre-operative planning with a good team approach between well trained surgeons, restorative dentists and lab technicians is essential for successful implant treatment.
- minimizing the trauma to the host tissue increases the potential for successful osseointegration.
- electro-surgery should not be attempted around metal implants, as they are conductive.

#### Single use devices

Do not reuse devices indicated for single use. (Use the device prior to the expiration date).

#### Side effects

Potential side effects and temporary symptoms: pain, swelling, phonetic difficulties, gingival inflammation. More persistent symptoms: the risks and complications with implants include, but are not limited to: (1) allergic reaction(s) to implant and/or abutment material; (2) breakage of the implant and/or abutment; (3) loosening of the abutment screw and/or retaining screw; (4) infection requiring revision of the dental implant; (5) nerve damage that could cause permanent weakness, numbness, or pain; (6) histologic responses possibly involving macrophages and/or fibroblasts; (7) formation of fat emboli; (8) loosening of the implant requiring revision surgery; (9) perforation of the maxillary sinus; (10) perforation of the labial and lingual plates and (11) bone loss possibly resulting in revision or removal.

#### Breakage

Implant and abutment fractures can occur when applied loads exceed the tensile or compressive strength of the material. Potential overloading conditions may result from: deficiencies in implant numbers, lengths and/or diameters to adequately support a restoration, excessive cantilever length, incomplete abutment seating, abutment angles greater than 30 degrees, occlusal interferences causing excessive lateral forces, patient parafunction (e.g. bruxing, clenching), loss or changes in dentition or functionality, inadequate prosthesis fit and physical trauma. Additional treatment may be necessary when any of the above conditions are present to reduce the possibility of hardware complications or failure.

#### Changes in performance

It is the responsibility of the clinician to instruct the patient on all appropriate contraindications, side effects and precautions as well as the need to seek the services of a trained dental professional if there are any changes in the performance of the implant (e.g. looseness of the prosthesis, infection or exudate around the implant, pain or any other unusual symptoms that the patient has not been told to expect).

#### Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established. The contact information for the manufacturer of this device to report a serious incident is as follows: [sicomplaints@southernimplants.com](mailto:sicomplaints@southernimplants.com)

#### Materials

Titanium abutments: Titanium grade 2, 3, 4, or 5.

Abutment screws: Titanium alloy Ti-90%, Al-6%, V-4%

Gold Alloy Au-61%, Ag-16.5%, Pt-13.5%, Cu -9%

**NOTE:** DC and Internal Hex (PROVATA®) Titanium abutments are anodized gold in colour.

#### Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into

account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

#### Disclaimer of liability

This product is part of the Southern Implants product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants product range and take full responsibility for the correct indications and use of this product. Southern Implants does not assume liability for damage due to incorrect use. Please note that some Southern Implants products may not be cleared or released for sale in all markets.

#### Healing

The healing time required for osseointegration depends on the individual and treatment protocol. It is the responsibility of the practitioner to decide when the implant can be restored. Good primary stability will govern if immediate loading can be done.

#### Implant care and maintenance

Potential implant patients should establish an adequate oral hygiene regime prior to Implant therapy. Proper post operative oral hygiene and implant maintenance instructions must be discussed with the patient, as this will determine the longevity and health of the Implants. The patient should maintain regular prophylaxis and evaluation appointments.

#### MR conditional

Non-clinical testing has demonstrated that the Southern Implants dental implants, metallic abutments and prosthetic screws are MR conditional. A patient with these devices can be safely scanned in a MR system meeting the following conditions:

- Static magnetic field of 1.5 Tesla and 3.0 Tesla only.
- Maximum spatial gradient magnetic field of 4000 Gauss/cm (40 T/m).
- Maximum MR system reported, head specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode) or whole body averaged specific absorption rate (wbSAR) of 1 W/kg.

Under the scan conditions defined above, the Southern Implants dental implants, abutments and prosthetic screws are expected to produce a maximum temperature rise of 5.8°C after 15 minutes of continuous scanning. In the non-clinical testing, the image artifact caused by the device extends approximately 20 mm from the Southern Implants dental implants, abutments and prosthetic screws, when imaged with a gradient echo pulse sequence and a 3.0 Tesla MRI system. Removable restorations should be taken out prior to scanning, as is done for watches, jewellery etc. Should there be no MR symbol on the product label, please note that this device has not been evaluated for safety and compatibility in the MR environment. This device has not been tested for heating or migration in the MR environment.

#### Basic UDI

| Product                       | Basic-UDI Number |
|-------------------------------|------------------|
| Basic-UDI for Metal Abutments | 600954403872     |

#### Related literature and catalogues

CAT-2004 - Tri-Nex Implants Product Catalogue  
 CAT-2005 - IT Product Catalogue  
 CAT-2020 - External Hex Implants Product Catalogue  
 CAT-2042 - Deep Conical Implants Product Catalogue  
 CAT-2043 - Internal Hex Implants Product Catalogue  
 CAT-2060 - PROVATA® Implants Product Catalogue  
 CAT-2069 - INVERTA® Implants Product Catalogue  
 CAT-2070 - Zygomatic Implants Product Catalogue  
 CAT-2010 - Osseointegrated Fixtures Product Catalogue

ARCHIVED

**Symbols and Warnings**

|   |   |   |   |   |   |   |   |  |   |   |   |   |
|---|---|---|---|---|---|---|---|--|---|---|---|---|
| <br><b>Manufacturer: Southern Implants</b><br>1 Albert Rd, P.O Box 605 IRENE, 0062,<br>South Africa.<br>Tel: +27 12 667 1046 |  |  |  |  |  |  |  |  |  |  |  |  |
|   | Prescription device*  | Sterilization using Irradiation   | Consult instruction for use   | Use by date (mm-yy)   | Do not reuse  | Batch code  | Do not use if package is damaged  | Medical Device   | Authorize Representative in the European Community                                    | Catalog number  | Date of Manufacture   |   |

\* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

All rights reserved. Southern Implants®, the Southern Implants logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale.

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