English	INSTRUCTIONS FOR USE: Southern Implants® Cylindrical Step Twist Drills and Taps
Español	INSTRUCCIONES DE USO: Machos de roscar y brocas helicoidales escalonadas cilíndricas de Southern Implants®
Italiano	ISTRUZIONI PER L'USO: Southern Implants® Punte elicoidali a gradini cilindrici e maschiatore
Français	MODE D'EMPLOI : Forets étagés et tarauds hélicoïdaux cylindriques Southern Implants®
Deutsch	GEBRAUCHSANWEISUNGEN: Southern Implants® Zylindrische Stufenbohrer und Gewindeschneider
Português	INSTRUÇÕES DE UTILIZAÇÃO: Southern Implants® Brocas cilíndricas helicoidais em etapas e machos

South Africa - Headquarters: 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

Subsidiaries

Australia

Southern Implants Australia
T: +61-(0)-8-9466-2627
E: info@southernimplants.com.au

Spain and Portugal

Southern Implants Iberica
T: +34 935 053 507
E: info@southernimplants.es

United Kingdom and Ireland

Southern Implants UK
T: +44-20-8899-6845 / 6 / 7
E: info@southernimplants.co.uk

USA and Canada

Southern Implants North America Inc.
T: +1-561-472-0990
E: customer care@southernimplants.com

EC REP

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

CH REP

MedEnvoy Switzerland: Gotthardstrasse 28, 6302 Zug, Switzerland

Intended use

Twist drills are used to drill a cylindrical osteotomy for implant placement. Twist drills are manufactured from stainless steel or titanium grade 5. The twist drills are medical devices. The D-ZYG-CH-29 and D-ZYG-CH-29S are multiple use drills intended for multiple use on multiple patients. All the other twist drills are intended for single use and multiple use on a single patient.

The bone taps are intended to be used to prepare the osteotomy for implant placement. The bone taps are implant length specific. The taps are medical devices and are intended for multiple uses on multiple patients.

Description

Southern Implants cylindrical twist drills and cylindrical bone taps, are as described in Table 1 and Table 2, respectively. These drills are made of stainless steel, Titanium (Grade 4) or Titanium Alloy (Grade 5) and some

are AlTiN coated. The drills/taps have a latch dimension compatible to ISO 1797. This is in order to connect the drill/tap to the handpiece of an implant motor unit. Drill sizes are identified by laser markings on the shaft of the drill. The drills have laser markings to indicate the drill depth.

Cylindrical twist drills and taps indicated in Table 1, are provided sterile:

- twist drills indicated, with 1 Number of uses is for single use only.
- taps indicated with 10 Number of uses, can be used up to 10 times or when the cutting efficiency deteriorates.

Indications for use

Dental implant therapy is indicated for patients with one or more missing teeth that wish to have the tooth/teeth restored with a dental implant and a prosthetic tooth/teeth.

Table 1 - Twist Drills

Drill code	Material	Coating (if any)	Description of product	Number of uses
External Hex, Internal Hex (M-series & Provata), TRI-NEX				
D-20T-MXX	Stainless Steel	-	Twist Drill Ø2.0	1
D-25T-MXX	Stainless Steel	-	Twist Drill Ø2.5	1
D-29T-MXX	Stainless Steel	-	Twist Drill Ø2.85	1
D-30T-MXX	Stainless Steel	-	Twist Drill Ø3.0	1
D-33T-MXX	Stainless Steel	-	Twist Drill Ø3.25	1
D-35T-MXX	Stainless Steel	-	Twist Drill Ø3.5	1
D-40T-MXX	Stainless Steel	-	Twist Drill Ø4.0	1
D-43T-MXX	Stainless Steel	-	Twist Drill Ø4.3	1
D-46T-MXX	Stainless Steel	-	Twist Drill Ø4.6	1
D-50T-MXX	Stainless Steel	-	Twist Drill Ø5.0	1
D-53T-MXX	Stainless Steel	-	Twist Drill Ø5.3	1
D-56T-MXX	Stainless Steel	-	Twist Drill Ø5.6	1
IT (Internal Octagon)				
D-220C	Stainless Steel	-	Twist Drill Ø2.2	1
D-275C	Stainless Steel	-	Twist Drill Ø2.75	1
D-350C	Stainless Steel	-	Twist Drill Ø3.5	1
D-430C	Stainless Steel	-	Twist Drill Ø4.3	1
D-220C-L	Titanium	-	Twist Drill Ø2.2 Longer shaft drill	1
D-275C-L	Titanium	-	Twist Drill Ø2.75 Longer shaft drill	1
D-350C-L	Titanium	-	Twist Drill Ø3.5 Longer shaft drill	1
D-430C-L	Titanium	-	Twist Drill Ø4.3 Longer shaft drill	1
DC (Deep Conical)				
D-DC20	Stainless Steel	AlTiN	Twist Drill DC Ø2.2	1
D-DC25	Stainless Steel	AlTiN	Twist Drill DC Ø2.5	1

D-DC27	Stainless Steel	AlTiN	Twist Drill DC Ø2.7	1
D-DC29	Stainless Steel	AlTiN	Twist Drill DC Ø2.85	1
D-DC32	Stainless Steel	AlTiN	Twist Drill DC Ø3.2	1
D-DC34	Stainless Steel	AlTiN	Twist Drill DC Ø3.35	1
D-DC37	Stainless Steel	AlTiN	Twist Drill DC Ø3.7	1
D-DC39	Stainless Steel	AlTiN	Twist Drill DC Ø3.85	1
D-DC42	Stainless Steel	AlTiN	Twist Drill DC Ø4.2	1
D-DC47	Stainless Steel	AlTiN	Twist Drill DC Ø4.7	1
D-DC49	Stainless Steel	AlTiN	Twist Drill DC Ø4.8	1
Zygomatic, Oncology, Zygan & Zygex implants				
D-ZYG-27	Titanium	Anodised yellow	Twist Drill Zygomatic Ø2.7	1
D-ZYG-27S	Titanium	Anodised yellow	Twist Drill Zygomatic Ø2.7 (Short)	1
D-ZYG-27ST-GSM	Titanium	Anodised yellow	Twist Drill Zygomatic Ø2.7 (Medium)	1
D-ZYG-27ST-GSL	Titanium	Anodised yellow	Twist Drill Zygomatic Ø2.7 (Long)	1
D-ZYG-29	Titanium	Anodised yellow	Twist Drill Zygomatic Ø2.9	1
D-ZYG-29S	Titanium	Anodised yellow	Twist Drill Zygomatic Ø2.9 (Short)	1
D-ZYG-35	Titanium	Anodised yellow	Twist Drill Zygomatic Ø3.5	1
D-ZYG-35S	Titanium	Anodised yellow	Twist Drill Zygomatic Ø3.5 (Short)	1
Guided Surgery				
D-20T-GS-20	Titanium	-	Twist Drill Guided Surgery Ø1.95	1
D-20ST-GS-20	Titanium	-	Twist Drill Guided Surgery Ø1.95 Adjustable	1
D-20T-GS-23	Titanium	-	Twist Drill Guided Surgery Ø1.95	1
D-20T-GS-28	Titanium	-	Twist Drill Guided Surgery Ø1.95	1
D-20ST-GS-30	Titanium	-	Twist Drill Guided Surgery Ø1.95 Adjustable	1
D-28T-GS-20	Titanium	-	Twist Drill Guided Surgery Ø2.75	1
D-28T-GS-23	Titanium	-	Twist Drill Guided Surgery Ø2.75	1
D-28T-GS-28	Titanium	-	Twist Drill Guided Surgery Ø2.75	1

FIRST				
D-20T-32RT	Titanium	-	Twist Drill FIRST technique Ø2.0	1
D-29T-32RT	Titanium	-	Twist Drill FIRST technique Ø2.85	1
D-33T-32RT	Titanium	-	Twist Drill FIRST technique Ø3.25	1
D-35T-32RT	Titanium	-	Twist Drill FIRST technique Ø3.5	1
D-40T-32RT	Titanium	-	Twist Drill FIRST technique Ø4.0	1
EXTRA Oral				
D-20E-03F	Titanium	-	Twist Drill IE implant Ø2.0 / 3mm	1
D-20E-04F	Titanium	-	Twist Drill IE implant Ø2.0 / 4mm	1
D-20E-06F	Titanium	-	Twist Drill IE implant Ø2.0 / 6mm	1
D-30E-03F	Titanium	-	Twist Drill IE implant Ø3.0 / 3mm	1
D-30E-04F	Titanium	-	Twist Drill IE implant Ø3.0 / 4mm	1
D-30E-06F	Titanium	-	Twist Drill IE implant Ø3.0 / 6mm	1

Table 2 - Bone taps

Drill code	Material	Coating (if any)	Description of product	Number of uses
External Hex				
D-TAP-IBN	Titanium	-	Tap for Hard Bone IBN	Up to 10
D-TAP-IBS	Titanium	-	Tap for Hard Bone IBS	Up to 10
D-TAP-I4B	Titanium	-	Tap for Hard Bone I4B	Up to 10
D-TAP-BA	Titanium	-	Tap for Hard Bone BA	Up to 10
Tri-Nex				
D-TAP-LS-35	Titanium	-	Tap for Hard Bone Ø3.5	Up to 10
D-TAP-LS-43	Titanium	-	Tap for Hard Bone Ø4.3	Up to 10
D-TAP-LS-50	Titanium	-	Tap for Hard Bone Ø5.0	Up to 10
DC (Deep Conical)				
D-TAP-DCC30	Titanium	-	Tap for Hard Bone Ø3.0	Up to 10
D-TAP-DCC35	Titanium	-	Tap for Hard Bone Ø3.5	Up to 10
D-TAP-DCC40	Titanium	-	Tap for Hard Bone Ø4.0	Up to 10
D-TAP-DCC50	Titanium	-	Tap for Hard Bone Ø5.0	Up to 10
IT (Internal Octagon)				
D-TAP-ITC3	Titanium	-	Tap for Hard Bone Ø3.3	Up to 10
D-TAP-ITC4	Titanium	-	Tap for Hard Bone Ø4.1	Up to 10
D-TAP-ITC5	Titanium	-	Tap for Hard Bone Ø4.9	Up to 10
FIRST				
D-TAP-35RT	Titanium	-	Tap for Hard Bone Ø3.5 FIRST Technique	Up to 10
D-TAP-43RT	Titanium	-	Tap for Hard Bone Ø4.3 FIRST Technique	Up to 10

Intended user

Maxillo-facial Surgeons, General Dentists, Orthodontists, Periodontists, Prosthodontists and other appropriately trained and experienced implant users.

Intended environment

The cylindrical twist drills and cylindrical bone taps are intended to be used in a clinical environment such as an operating theater or a dentist consultation room.

Intended patient population

This device is used in the dental restoration of partially or fully edentulous patients in the upper or lower jaw. Restorations may comprise of single teeth, partial or full bridges and may be fixed or removable.

Clinical procedure

A proper clinical and radiological evaluation must be done to determine the bone dimensions and bone quality. Ensure that all instruments and drills are in a good condition.

Clinical benefits


Implant therapy results in restoration of chewing function, speech and aesthetics in patients with missing teeth.

Surgical procedure

Southern Implants provides the user with different drill options, for placement of cylindrical implants, refer to each individual product catalogue for drill protocols for different bone quality. Ensure that all instruments and drills are in a good condition. Blunt drills may cause damage to the bone which could compromise osseointegration.

The drill sizes are identified by different laser markings on the shaft of the drill. The drills have different laser markings on the body of the drill which corresponds to the length of implant being placed. Refer to individual product catalogues for more information.

(Fig. 1)


1. Drill in the planned direction to the full depth of implant length being placed as indicated on the markings on the drill. Cylindrical drills extend up to 1mm longer than the implant, when seated. Allow for this additional length when drilling near vital anatomical structures.

NOTE: connect the drill/tap latch to the handpiece. If the latch is not inserted fully into the handpiece the torque is applied to the latch, resulting in possible twisting of the latch or damage to the handpiece. Consult the instructions for use of the handpiece to ensure proper engagement of the latch.

2. Drill at sufficient speed (between 1000-1500rpm for Cylindrical Drills), under constant irrigation with sterile saline.
3. Use an up-and-down motion with the handpiece, without stopping the motor. This will allow the irrigation to flush away bone debris on the drill.
4. Insert the Direction Indicator (I-DI) after using the 2mm Twist drill, to verify the alignment with adjacent implants or teeth. A radiograph can be taken at this stage to verify depth and direction. If the drilling direction is incorrect, start a new direction with the initial drill.
5. Gradually enlarge the osteotomy in a stepwise approach to the desired diameter and depth, the depth can be determined by a depth gauge:
 - I-DG-20 for External Hex, Internal Hex, IT and TRI-Nex.
 - I-DG-24 for External Hex, Internal Hex, IT and TRI-Nex.
 - I-ZYG-DG-1 for Zygomatic Implants.
 - I-DG-DC for DC implants.

Figure 2: (DC implant used for illustration purposes)

6. During surgery the clinician will be able to assess the bone quality and should use dense bone protocols when necessary, to prepare the site. This is to avoid the implant getting stuck before it is properly seated in the osteotomy. It is recommended to undersize the osteotomy in soft bone, and use a bone tap in dense bone.
7. Preparing the site further should involve: ensure the drill goes to full depth and/or use of the optional bone tap to pre-tap the site. Tap at low speed (25rpm) and after tapping to full depth, switch the handpiece to reverse mode to remove the tap.

Some cylindrical twist drills and cylindrical bone taps can be re-used (refer to Table 1 and 2). It is recommended to maintain a log of these drills, recording the number of uses. Prior to re-processing these components, it should be thoroughly inspected and tested to determine its suitability for re-use.

NOTE:

- Southern Implants taps feature a W&H hex on the shaft, to achieve higher torque use a converter (I-WI-C-S) over the shaft of the tap. The converter will engage the W&H hex on the tap and convert the tap to be used with a Southern Implants torque wrench. This will avoid the latch from getting stuck in the handpiece.
- avoid lateral pressure (bending) on the drills during drilling procedures.
- lateral pressure to the drill can cause drill fracture.
- verify the drill is securely locked into the handpiece before drilling procedure starts.

Before surgery

All components, instruments and tooling used during the clinical or laboratory procedure must be maintained in good condition and care must be taken that instrumentation does not damage implants or other components.

During surgery

Care must be taken that parts are not swallowed during any of the procedures, a rubber-dam application is recommended when appropriate. All components, instruments and tooling used during the clinical or laboratory procedure must be maintained in good condition and care must be taken that instrumentation does not damage implants or other components or equipment.

Post-surgery


Regular patient follow-up, and proper oral hygiene must be achieved to ensure favourable long-term results.

Drill extension


When a drill extension is used (I-DE-K / I-DE-G), care must be taken to ensure that the latch is fully engaged to prevent distortion. See *Figure 3 and 4* below:

- drill extensions must NOT be used with Ø6mm and larger drills, use longer shaft drills instead.
- drill extensions must NOT be used with bone taps.

(Fig. 3)

(Fig. 3)

The orientations indicated in *Figure 3* and *Figure 4* ensure that the catch feature of the drill extension (circled in blue) slots into the latch groove of the drill (circled in green). This prevents the drill from sliding out of the drill extension.

NOTE:

- do not apply more than 40-45Ncm to any latch type drill/instrument, this could cause damage to the handpiece and latch of the instrument.
- blunt drills cause excessive torque and result in damage to the handpiece or drill latch.

Storage, cleaning and sterilisation

These devices are supplied sterile (sterilised by gamma irradiation). Sterility is assured unless the container or seal is damaged or opened. If packaging is damaged do not use the product and contact your Southern representative or return to Southern Implants. The devices must be stored in a dry place at room temperature and not exposed to direct sunlight. Incorrect storage may influence device characteristics.

Sterilisation

Southern Implants recommends the following procedure to sterilise the instruments prior to use/re-use:

Methods to sterilise these devices:

1. pre-vacuum sterilisation method: steam sterilise the instruments at 132°C (270°F) at 180-220kPa for 4 minutes. Dry for at least 20 minutes in the chamber. Only an approved wrap or pouch for steam sterilisation must be used.
2. pre-vacuum sterilisation method: wrapped, steam sterilise at 135°C (275°F) for 3 minutes. Dry for 20 minutes in the chamber. Use a wrap or pouch that is cleared for the indicated steam sterilisation cycle.

NOTE: users in the USA must ensure that the steriliser, wrap or pouch and all steriliser accessories are cleared by the FDA, for the intended sterilisation cycle.

Contraindications

Do not use in patients:

- who are medically unfit for dental implant procedures.
- where adequate numbers of implants could not be placed to achieve full functional support of the prosthesis.
- who are allergic or have hypersensitivity to pure titanium or titanium alloy (Ti-6Al-4V), gold, palladium, platinum or iridium.
- who are under the age of 18, have poor bone quality, blood disorders, infected implant site, vascular impairment, uncontrolled diabetes, drug or alcohol abuse, chronic high dose steroid therapy, anti-coagulant therapy, metabolic bone disease, radiotherapy treatment and sinus pathology.

Warnings

- The use of non-sterile items and reuse of single-use devices can lead to secondary infections of the tissue or transfer infectious diseases.
- Failure to perform appropriate cleaning, re-sterilisation and storage procedures as per the Instructions for Use document can result in device damage or secondary infection.
- For the safe and effective use of dental implants, it is recommended that specialised training be undertaken.

- Failure to accurately estimate lengths of drills relative to radiographic measurements can result in permanent injury to nerves or other vital structures.
- Use of the device with incompatible or non-corresponding devices can result in device failure or poor performance.

Cautions

New and experienced implant users should do training before using a new system or attempting to do a new treatment method. Take special care when treating patients who have local or systemic factors that could affect the healing of the bone and soft tissue. (i.e. poor oral hygiene, uncontrolled diabetes, are on steroid therapy, smokers, infection in the nearby bone and patients who had oro-facial radiotherapy).

Thorough screening of prospective implant candidates must be performed including:

- a comprehensive medical and dental history.
- visual and radiological inspection to determine adequate bone dimensions, anatomical landmarks, occlusal conditions and periodontal health.
- bruxism and unfavourable jaw relations must be taken into account.
- proper pre-operative planning with a good team approach between well-trained surgeons, restorative dentists and lab technicians is essential for successful implant treatment.
- minimizing the trauma to the host tissue increases the potential for successful osseointegration.
- electro-surgery should not be attempted around metal implants, as they are conductive.

Single use devices

Do not reuse devices indicated for single use. (Use the device prior to the expiration date).

Do not reuse implants, single use drills, cover screws, temporary abutments and abutments. Reusing these components may result in:

- damage on the surface or critical dimensions, which may result in performance and compatibility degradation.
- adds the risk of cross-patient infection and contamination if single-use items are reused.

Southern Implants does not accept any responsibility for complications associated with reused components.

Reusable devices

Prior to re-using this device, it needs to be inspected, if there are signs of visible corrosion, deformed or twisted connections, dull cutting edges, wear and damage, this device shall be disposed of. After inspection, and if re-use seems fit, the devices are cleaned and sterilised. If re-use seems fit:

- containment: as soon as practically possible, remove all visible residue after use (bone, blood or tissue), by immersing the instrument in cold water (Dried soil is difficult to remove).
- pre-cleaning: rinse with lukewarm water for 3 minutes and remove hardened debris with a soft nylon brush. Avoid mechanical damage during cleaning.
- manual cleaning or automated cleaning: prepare an ultrasonic bath with suitable detergent, sonicate for 20 minutes (alternative methods can be used if proven by the end user). Rinse with purified/sterile water. Load devices into a thermo-disinfector. Run the cleaning and disinfection cycle, followed by the drying cycle.

NOTE: always follow the instructions for use of the manufacturers of cleaning agents and disinfectants.

- Drying: dry the instruments with filtered compressed air or single use lint-free wipes. Pack the instruments as quickly as possible after removal into the storage container. If additional drying is necessary, dry in a clean location. Moisture on these devices can cause corrosion and deterioration of the cutting edges.
- Inspection: do a visual inspection of the items to check for any damage.
- Packaging: use the correct packaging material as indicated for steam sterilisation to ensure sterility is maintained. Double packaging is recommended.

Side effects

Potential side effects and temporary symptoms: pain, swelling, phonetic difficulties, gingival inflammation. More persistent symptoms: the risks and complications with implants include, but are not limited to: (1) allergic reaction(s) to implant and/or abutment material; (2) breakage of the implant and/or abutment; (3) loosening of the abutment screw and/or retaining screw; (4) infection requiring revision of the dental implant; (5) nerve damage that could cause permanent weakness, numbness, or pain; (6) histologic responses possibly involving macrophages and/or fibroblasts; (7) formation of fat emboli; (8) loosening of the implant requiring revision surgery; (9) perforation of the maxillary sinus; (10) perforation of the labial and lingual plates; and (11) bone loss possibly resulting in revision or removal.

Breakage

Breakage or deformation can result when applied loads exceed the strength of the material. Overloading may result from: excessive pressure exerted on the handpiece; excessive lateral forces due to non-concentric loading; excessive torque; incorrect fit of the latch attachment into the handpiece; use outside of manufacturer recommendations; use of the device with incompatible hardware. In the above cases, additional treatment may be necessary to reduce the likelihood of complications.

Changes in performance

It is the responsibility of the clinician to ensure all warnings, precautions and recommended techniques are adhered to. In case of unexpected or unusual changes in performance, the manufacturer should be contacted immediately.

Performance requirements

1. Can be used in the upper or lower jaw.
2. Cuts bone effectively.
3. Prepares implant site to precise dimensions of intended implant.
4. Reusable drill devices can be used for 10 osteotomy sites.
5. Does not cause adverse biologic events as a result of toxicity or particle release.

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established. The contact information for the manufacturer of this device to report a serious incident is as follows: sicomplaints@southernimplants.com

Materials

Drills/Taps:	Stainless Steel, Titanium (Grade 4) or Titanium Alloy (Ti-6AL-4V)
Drill Coating:	None, or Titanium Nitride (TiN), or Aluminium Titanium Nitride (AlTiN)

Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

Disclaimer of liability

This product is part of the Southern Implants product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants product range and take full responsibility for the correct indications and use of this product. Southern Implants does not assume liability for damage due to incorrect use. Please note that some Southern Implants products may not be cleared or released for sale in all markets.

Basic UDI

Product	Basic-UDI Number
Basic-UDI for Drills and Hand Piece Devices	600954403875
Basic-UDI for Reusable Instruments	600954403876

Related Literature and Catalogues

CAT-2004 - Tri-Nex Implants Product Catalogue

CAT-2005 - IT Implants Product Catalogue

CAT-2020 - External Hex Implants Product Catalogue

CAT-2042 - Deep Conical Implants Product Catalogue

CAT-2043 - Internal Hex Implants Product Catalogue

CAT-2060 - PROVATA® Implants Product Catalogue

CAT-2069 - INVERTA® Implants Product Catalogue

CAT-2070 - Zygomatic Implants Product Catalogue

ARCHIVED

Symbols and Warnings

 Manufacturer: Southern Implants 1 Albert Rd, P.O Box 605 IRENE, 0062, South Africa. Tel: +27 12 667 1046													
Prescription device* Sterilization using Irradiation Consult instruction for use Use by date (mm-yy) Do not reuse Batch code Do not use if package is damaged Medical Device Date of Manufacture Authorize Representative in the European Community Catalog number Authorized Representative for Switzerland	* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist. Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law. All rights reserved. Southern Implants®, the Southern Implants logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale.												

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