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
South Africa - Headquarters: 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

Subsidiaries

Australia

Southern Implants Australia
T: +61-(0)-8-9466-2627
E: info@southernimplants.com.au

Spain and Portugal

Southern Implants Iberica
T: +34 935 053 507
E: info@southernimplants.es

United Kingdom and Ireland

Southern Implants UK
T: +44-20-8899-6845 / 6 / 7
E: info@southernimplants.co.uk

USA and Canada

Southern Implants North America Inc.
T: +1-561-472-0990
E: customercare@southernimplants.com

Description

These drills are intended to enlarge the channel made by the round burr drill in order to prepare a cylindrical osteotomy site for a cylindrical endosseous dental implant. This entails drilling a hole in the mandible or maxilla or making an existing hole larger in diameter. The majority of these drills are for single use and supplied sterile. Twist drills are made from titanium grade 5 or stainless steel. The drills are available in diameters between Ø2 mm and Ø5.6 mm.

Intended use

Twist drills are intended to aid in the treatment of partially or fully edentulous patients eligible for placement of one or more dental implants as a means of fixing a permanent or removable single crown, partial or full-arch dental prosthesis in the upper or lower jaw. All of these twist drills are classified as medical devices and are intended for single use on a single patient.

Indications for use

The indications for use for the device are not dissimilar to that of dental implants. Southern Implants dental implants are intended for both one and two-stage surgical procedures in the following situations and with the following clinical protocols:

- replacing single and multiple missing teeth in the mandible and maxilla.
- immediate placement in extraction sites and in situations with a partially or completely healed alveolar ridge.
- immediate loading in all indications, except in single tooth situations on implants shorter than 8 mm or in soft bone (type IV) where implant stability may be difficult to obtain, and immediate loading may not be appropriate.

Intended user

Maxillo-facial Surgeons, General Dentists, Orthodontists, Periodontists, Prosthodontists and other appropriately trained and experienced implant users.

Intended environment

The devices are intended to be used in a clinical environment such as an operating theatre or a dentist consultation room.

Intended patient population

The intended patient population for the device is not dissimilar to that of dental implant therapy. The intended patient population for implant therapy is partially or fully edentulous patients requiring prosthetic dental restoration in the upper or lower jaw. Restorations may comprise single teeth, partial or full bridges and are fixed restorations. Additionally, the intended population for some device variants includes those who have previously undergone dental implant therapy.

Compatibility information

Southern Implants' implants twist drills are as described in Table A. Details on site preparation sequences recommended by Southern Implants® can be found in the relevant product catalogues.

Table A – Twist Drills

Drill Code	Material	Coating (if any)	Number of Uses
EXTERNAL HEX, INTERNAL HEX (M-SERIES & PROVATA®), TRI-NEX®			
D-12T-M15	Stainless Steel (DIN 1.4197)	-	1
D-20T-MXX	Stainless Steel (DIN 1.4197)	-	1
D-25T-MXX	Stainless Steel (DIN 1.4197)	-	1
D-29T-MXX	Stainless Steel (DIN 1.4197)	-	1
D-30T-MXX	Stainless Steel (DIN 1.4197)	-	1
D-33T-MXX	Stainless Steel (DIN 1.4197)	-	1
D-35T-MXX	Stainless Steel (DIN 1.4197)	-	1

D-40T-MXX	Stainless Steel (DIN 1.4197)	-	1
D-43T-MXX	Stainless Steel (DIN 1.4197)	-	1
D-46T-MXX	Stainless Steel (DIN 1.4197)	-	1
D-50T-MXX	Stainless Steel (DIN 1.4197)	-	1
D-53T-MXX	Stainless Steel (DIN 1.4197)	-	1
D-56T-MXX	Stainless Steel (DIN 1.4197)	-	1
INTERNAL OCTAGON (IT)			
D-220C	Stainless Steel (DIN 1.4197)	-	1
D-275C	Stainless Steel (DIN 1.4197)	-	1
D-350C	Stainless Steel (DIN 1.4197)	-	1
D-430C	Stainless Steel (DIN 1.4197)	-	1
D-220C-L	Stainless Steel (DIN 1.4197)	-	1
D-275C-L	Stainless Steel (DIN 1.4197)	-	1
D-350C-L	Stainless Steel (DIN 1.4197)	-	1
D-430C-L	Stainless Steel (DIN 1.4197)	-	1
DEEP CONICAL (DC)			
D-DC20	Stainless Steel (DIN 1.4197)	ALTiN	1
D-DC25	Stainless Steel (DIN 1.4197)	ALTiN	1
D-DC27	Stainless Steel (DIN 1.4197)	ALTiN	1
D-DC29	Stainless Steel (DIN 1.4197)	ALTiN	1
D-DC32	Stainless Steel (DIN 1.4197)	ALTiN	1
D-DC34	Stainless Steel (DIN 1.4197)	ALTiN	1
D-DC37	Stainless Steel (DIN 1.4197)	ALTiN	1
D-DC39	Stainless Steel (DIN 1.4197)	ALTiN	1
D-DC42	Stainless Steel (DIN 1.4197)	ALTiN	1
D-DC47	Stainless Steel (DIN 1.4197)	ALTiN	1
D-DC49	Stainless Steel (DIN 1.4197)	ALTiN	1
ZYGOMATIC			
D-ZYG-27	Titanium grade 5 (ASTM F136)	Anodized Yellow	1
D-ZYG-27S	Titanium grade 5 (ASTM F136)	Anodized Yellow	1
D-ZYG-27ST-GSM	Titanium grade 5 (ASTM F136)	Anodized Yellow	1
D-ZYG-27ST-GSL	Titanium grade 5 (ASTM F136)	Anodized Yellow	1
D-ZYG-29	Titanium grade 5 (ASTM F136)	Anodized Yellow	1
D-ZYG-29S	Titanium grade 5 (ASTM F136)	Anodized Yellow	1
D-ZYG-35	Titanium grade 5 (ASTM F136)	Anodized Yellow	1
D-ZYG-35S	Titanium grade 5 (ASTM F136)	Anodized Yellow	1
GUIDED SURGERY			
D-20T-GS-20	Titanium grade 5 (ASTM F136)	-	1
D-20T-GS-23	Titanium grade 5 (ASTM F136)	-	1
D-20T-GS-28	Titanium grade 5 (ASTM F136)	-	1
D-28T-GS-20	Titanium grade 5 (ASTM F136)	-	1
D-28T-GS-23	Titanium grade 5 (ASTM F136)	-	1
D-28T-GS-28	Titanium grade 5 (ASTM F136)	-	1
FIRST® ROHNER TECHNIQUE			
D-20T-32RT	Titanium grade 5 (ASTM F136)	-	1
D-29T-32RT	Titanium grade 5 (ASTM F136)	-	1
D-33T-32RT	Titanium grade 5 (ASTM F136)	-	1
D-35T-32RT	Titanium grade 5 (ASTM F136)	-	1
D-40T-32RT	Titanium grade 5 (ASTM F136)	-	1
IE IMPLANTS			
D-20E-03F	Titanium grade 5 (ASTM F136)	-	1
D-20E-04F	Titanium grade 5 (ASTM F136)	-	1
D-20E-06F	Titanium grade 5 (ASTM F136)	-	1
D-30E-03F	Titanium grade 5 (ASTM F136)	-	1
D-30E-04F	Titanium grade 5 (ASTM F136)	-	1

D-30E-06F	Titanium grade 5 (ASTM F136)	-	1
CORTICAL PERFORATORS			
D-GS-CP	Titanium grade 5 (ASTM F136)	-	1
D-GS-CP4	Titanium grade 5 (ASTM F136)	-	1

Clinical benefits

The clinical benefits are not dissimilar to that of dental implant therapy. The clinical benefits of dental implant therapy include improved chewing function, speech, aesthetics and patient psychological wellbeing. Through this procedure patients can expect to have their missing teeth replaced and/or crowns restored.

The clinical benefits specific to drills include:

- Effective cutting of bone within acceptable time limits.
- Precise preparation of the implant site to the dimensions of the intended implant.
- Can remain effective after multiple uses to the extent stipulated by the manufacturer.

Storage, cleaning and sterilisation

The component is supplied sterile (sterilised by gamma irradiation) and intended for single-use prior to the expiration date (see packaging label). Sterility is assured unless the container or seal is damaged or opened. If packaging is damaged do not use the product and contact your Southern representative or return to Southern Implants®. The devices must be stored in a dry place at room temperature and not exposed to direct sunlight. Incorrect storage may influence device characteristics. Do not reuse components indicated for single-use only. Reusing these components may result in:

- damage to the surface or critical dimensions, which may result in performance and compatibility degradation.
- adds the risk of cross-patient infection and contamination if single-use items are reused.

Southern Implants® does not accept any responsibility for complications associated with reused single-use components.

Contraindications

The contraindications for this system are not dissimilar to those of dental implant therapy. Contra-indications to implant therapy include:

- patients medically unfit for oral surgical procedures.
- where an inadequate number of implants can be placed limiting the functional support of the prosthesis.
- patients under the age of 18 years.
- patients with poor bone quality.
- patients with blood disorders.
- presence of infection at the implant site.
- patients with vascular impairment.
- patients with uncontrolled diabetes.
- patients with drug or alcohol abuse dependencies.
- patients undergoing chronic high dose steroid therapy.
- patients undergoing anti-coagulant therapy.
- patients with metabolic bone disease.
- patients undergoing radiotherapy treatment.
- Patients with pure titanium, titanium alloy (Ti6Al4V), gold, palladium, or iridium allergies or hypersensitivities.

Warnings and precautions

THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING.

- For the safe and effective use of dental implants it is strongly suggested that specialised training be undertaken, including hands-on training to learn proper technique, biomechanical requirements and radiographic evaluations.

- Products must be secured against aspiration when handled intraorally. Aspiration of products may lead to infection or unplanned physical injury.

Responsibility for proper patient selection, adequate training, experience in the placement of implants and providing appropriate information for informed consent rests with the practitioner. Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone. Implant failure increases when implants are placed in irradiated bone as radiotherapy can result in progressive fibrosis of vessels and soft tissue, leading to diminished healing capacity.

It is important to be aware and avoid damage to vital structures such as nerves, veins and arteries. Injuries to vital anatomical structures may cause serious complications such as injury to the eye, nerve damage and excessive bleeding. It is essential to protect the infraorbital nerve. Failing to identify actual measurements relative to the radiographic data could lead to complications.

New and experienced implant users should do training before using a new system or attempting to do a new treatment method. Take special care when treating patients who have local or systemic factors that could affect the healing of the bone and soft tissue (i.e., poor oral hygiene, uncontrolled diabetes, are on steroid therapy, smokers, infection in the nearby bone and patients who had orofacial radiotherapy).

Thorough screening of prospective implant candidates must be performed including:

- a comprehensive medical and dental history.
- visual and radiological inspection to determine adequate bone dimensions, anatomical landmarks, occlusal conditions and periodontal health.
- bruxism and unfavourable jaw relations must be taken into account.
- proper preoperative planning with a good team approach between well trained surgeons, restorative dentists and lab technicians is essential for successful implant treatment.
- minimising the trauma to the host tissue increases the potential for successful osseointegration.
- electrosurgery should not be attempted around metal implants as they are conductive.

Should the device not operate as intended, it must be reported to the manufacturer of the device. The contact information for the manufacturer of this device to report a change in performance is: sicomplaints@southernimplants.com.

Side effects

The side effects of the use of the system are not dissimilar to those of dental implant therapy. Possible side effects to implant therapy include:

- pain.
- swelling.
- phonetic difficulties.
- gingival inflammation.

Less common but more persistent symptoms include, but are not limited to:

- allergic reaction(s) to implant and/or abutment material.
- breakage of the implant and/or abutment.
- loosening of the abutment screw and/or retaining screw.
- infection requiring revision of the dental implant.
- nerve damage resulting in permanent weakness, numbness, or pain.
- histologic responses with possible macrophage and/or fibroblast involvement.
- fat emboli formation.
- loosening of the implant requiring revision surgery.
- perforation of the maxillary sinus.
- perforation of the labial and lingual plates.
- bone loss possibly resulting in revision or removal of the implant.

Precaution: maintaining sterility protocol

Implants are packaged as follows:

1. An outer package consisting of a rigid, clear box which acts as protection for the inner package.
2. The inner package consisting of a blister pack (clear plastic-formed blister base with a TYVEK “peel-back” lid).
3. Within the inner package, there is a hollow tube which contains one implant suspended from a titanium ring, this ensures the implant never touches the inside of the plastic tube.
4. Labelling information is located on the surface of the peel-back lid and on the outside of the rigid box.

Care must be taken to maintain the sterility of the implant by proper opening of the packaging and handling of the implant.

1. Open the implant package in the non-sterile field, with non-sterile gloves, tear the address label to open the box.
2. With non-sterile gloves, remove the inner blister pack. Do not place the plastic box or blister pack-lid onto the sterile field. The contents of this inner package are sterile.
3. The sealed blister is to be opened by an assistant (with nonsterile gloves): remove the TYVEK lid and drop or place the sterile tube onto the sterile field, open the tube cap and attach the implant placement tool onto the implant and carefully remove from the sterile tube. Do not touch the sterile implant.

Other sterile components are packed in a peel pouch or blister base with a “peel-back” lid. Labelling information is located on the bottom half of the pouch, inside the packet or on the surface of the peel-back lid. Sterility is assured unless the pouch is damaged or opened. Non-sterile components are supplied clean but not sterile in a peel pouch or blister base with peel-back lid. Labelling information is located on the bottom half of the pouch or on the surface of the peel-back lid.

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established.

The contact information for the manufacturer of this device to report a serious incident is as follows: sicomplaints@southernimplants.com.

Materials

Material type Titanium grade 5 (ASTM F136) or Stainless Steel (DIN 1.4197)

Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

Summary of Safety and Clinical Performance (SSCP)

As required by the European Medical Device Regulation (MDR; EU2017/745), a Summary of Safety and Clinical Performance (SSCP) is available for perusal with regard to Southern Implants® product ranges.

The relevant SSCP can be accessed at <https://ec.europa.eu/tools/eudamed>.

NOTE: the above website will be available upon the launch of the European Database on Medical Devices (EUDAMED).

Disclaimer of liability

This product is part of the Southern Implants® product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants® product range and take full responsibility for the correct indications and use of this product. Southern Implants® does not assume liability for damage due to incorrect use. Please note that some Southern Implants® products may not be cleared or released for sale in all markets.

Basic UDI

Product	Basic-UDI Number
Basic-UDI for Drills and Handpiece Devices	6009544038759C

Related literature and catalogues

CAT-2004 - Tri-Nex® Implants Product Catalogue
 CAT-2020 - External Hex Implants Product Catalogue
 CAT-2042 - Deep Conical (DC) Implants Product Catalogue
 CAT-2043 - Internal Hex (M-Series) Implants Product Catalogue
 CAT-2060 - Internal Hex (PROVATA®) Implants Product Catalogue
 CAT-2005 - IT (Internal Octagon) Implants Product Catalogue
 CAT-2070 - Zygomatic Implants Product Catalogue
 CAT-2069 - INVERTA® Implants Product Catalogue
 CAT-2068 - SiGuided Implants Product Catalogue
 CAT-2010 - Osseointegrated Fixtures Product Catalogue

Symbols and warnings

											
Manufacturer: Southern Implants® 1 Albert Rd, P.O Box 605 IRENE, 0062, South Africa. Tel: +27 12 667 1046	CE mark 2797	Prescription device*	Sterilised using irradiation	Non-sterile	Use by date (mm-yy)	Do not reuse	Do not resterilise	Catalogue number	Batch code	Medical device	Authorised representative in the European Community
											
Authorised representative for Switzerland	Date of manufacture	Magnetic Resonance conditional	Magnetic Resonance safe	Single sterile barrier system with protective packaging inside	Single sterile barrier system	Consult instruction for use	Caution	Keep away from sunlight	Do not use if package is damaged		

* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

All rights reserved. Southern Implants®, the Southern Implants® logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants®. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale. It is the responsibility of the clinician to inspect the symbols that appear on the packaging of the product in use.

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D-25T-MXX	Acier inoxydable (DIN 1.4197)	-	1
D-29T-MXX	Acier inoxydable (DIN 1.4197)	-	1
D-30T-MXX	Acier inoxydable (DIN 1.4197)	-	1
D-33T-MXX	Acier inoxydable (DIN 1.4197)	-	1
D-35T-MXX	Acier inoxydable (DIN 1.4197)	-	1
D-40T-MXX	Acier inoxydable (DIN 1.4197)	-	1
D-43T-MXX	Acier inoxydable (DIN 1.4197)	-	1
D-46T-MXX	Acier inoxydable (DIN 1.4197)	-	1
D-50T-MXX	Acier inoxydable (DIN 1.4197)	-	1
D-53T-MXX	Acier inoxydable (DIN 1.4197)	-	1
D-56T-MXX	Acier inoxydable (DIN 1.4197)	-	1
OCTOGONE INTERNE (IT)			
D-220C	Acier inoxydable (DIN 1.4197)	-	1
D-275C	Acier inoxydable (DIN 1.4197)	-	1
D-350C	Acier inoxydable (DIN 1.4197)	-	1
D-430C	Acier inoxydable (DIN 1.4197)	-	1
D-220C-L	Acier inoxydable (DIN 1.4197)	-	1
D-275C-L	Acier inoxydable (DIN 1.4197)	-	1
D-350C-L	Acier inoxydable (DIN 1.4197)	-	1
D-430C-L	Acier inoxydable (DIN 1.4197)	-	1
CONIQUE PROFOND (DC)			
D-DC20	Acier inoxydable (DIN 1.4197)	ALTiN	1
D-DC25	Acier inoxydable (DIN 1.4197)	ALTiN	1
D-DC27	Acier inoxydable (DIN 1.4197)	ALTiN	1
D-DC29	Acier inoxydable (DIN 1.4197)	ALTiN	1
D-DC32	Acier inoxydable (DIN 1.4197)	ALTiN	1
D-DC34	Acier inoxydable (DIN 1.4197)	ALTiN	1
D-DC37	Acier inoxydable (DIN 1.4197)	ALTiN	1
D-DC39	Acier inoxydable (DIN 1.4197)	ALTiN	1
D-DC42	Acier inoxydable (DIN 1.4197)	ALTiN	1
D-DC47	Acier inoxydable (DIN 1.4197)	ALTiN	1
D-DC49	Acier inoxydable (DIN 1.4197)	ALTiN	1
ZYGOMATIQUE			
D-ZYG-27	Titane grade 5 (ASTM F136)	Jaune anodisé	1
D-ZYG-27S	Titane grade 5 (ASTM F136)	Jaune anodisé	1
D-ZYG-27ST-GSM	Titane grade 5 (ASTM F136)	Jaune anodisé	1
D-ZYG-27ST-GSL	Titane grade 5 (ASTM F136)	Jaune anodisé	1
D-ZYG-29	Titane grade 5 (ASTM F136)	Jaune anodisé	1
D-ZYG-29S	Titane grade 5 (ASTM F136)	Jaune anodisé	1
D-ZYG-35	Titane grade 5 (ASTM F136)	Jaune anodisé	1
D-ZYG-35S	Titane grade 5 (ASTM F136)	Jaune anodisé	1
CHIRURGIE GUIDÉE			
D-20T-GS-20	Titane grade 5 (ASTM F136)	-	1
D-20T-GS-23	Titane grade 5 (ASTM F136)	-	1
D-20T-GS-28	Titane grade 5 (ASTM F136)	-	1
D-28T-GS-20	Titane grade 5 (ASTM F136)	-	1
D-28T-GS-23	Titane grade 5 (ASTM F136)	-	1
D-28T-GS-28	Titane grade 5 (ASTM F136)	-	1
TECHNIQUEFIRST® ROHNER			
D-20T-32RT	Titane grade 5 (ASTM F136)	-	1
D-29T-32RT	Titane grade 5 (ASTM F136)	-	1
D-33T-32RT	Titane grade 5 (ASTM F136)	-	1
D-35T-32RT	Titane grade 5 (ASTM F136)	-	1
D-40T-32RT	Titane grade 5 (ASTM F136)	-	1

IMPLANTS IE			
D-20E-03F	Titane grade 5 (ASTM F136)	-	1
D-20E-04F	Titane grade 5 (ASTM F136)	-	1
D-20E-06F	Titane grade 5 (ASTM F136)	-	1
D-30E-03F	Titane grade 5 (ASTM F136)	-	1
D-30E-04F	Titane grade 5 (ASTM F136)	-	1
D-30E-06F	Titane grade 5 (ASTM F136)	-	1
PERFORATEURS CORTICAUX			
D-GS-CP	Titane grade 5 (ASTM F136)	-	1
D-GS-CP4	Titane grade 5 (ASTM F136)	-	1

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D-30T-MXX	Edelstahl (DIN 1.4197)	–	1
D-33T-MXX	Edelstahl (DIN 1.4197)	–	1
D-35T-MXX	Edelstahl (DIN 1.4197)	–	1
D-40T-MXX	Edelstahl (DIN 1.4197)	–	1
D-43T-MXX	Edelstahl (DIN 1.4197)	–	1
D-46T-MXX	Edelstahl (DIN 1.4197)	–	1
D-50T-MXX	Edelstahl (DIN 1.4197)	–	1
D-53T-MXX	Edelstahl (DIN 1.4197)	–	1
D-56T-MXX	Edelstahl (DIN 1.4197)	–	1
Innenachtkant (IT)			
D-220C	Edelstahl (DIN 1.4197)	–	1
D-275C	Edelstahl (DIN 1.4197)	–	1
D-350C	Edelstahl (DIN 1.4197)	–	1
D-430C	Edelstahl (DIN 1.4197)	–	1
D-220C-L	Edelstahl (DIN 1.4197)	–	1
D-275C-L	Edelstahl (DIN 1.4197)	–	1
D-350C-L	Edelstahl (DIN 1.4197)	–	1
D-430C-L	Edelstahl (DIN 1.4197)	–	1
Tief konisch (DC)			
D-DC20	Edelstahl (DIN 1.4197)	ALTiN	1
D-DC25	Edelstahl (DIN 1.4197)	ALTiN	1
D-DC27	Edelstahl (DIN 1.4197)	ALTiN	1
D-DC29	Edelstahl (DIN 1.4197)	ALTiN	1
D-DC32	Edelstahl (DIN 1.4197)	ALTiN	1
D-DC34	Edelstahl (DIN 1.4197)	ALTiN	1
D-DC37	Edelstahl (DIN 1.4197)	ALTiN	1
D-DC39	Edelstahl (DIN 1.4197)	ALTiN	1
D-DC42	Edelstahl (DIN 1.4197)	ALTiN	1
D-DC47	Edelstahl (DIN 1.4197)	ALTiN	1
D-DC49	Edelstahl (DIN 1.4197)	ALTiN	1
Zygomatic			
D-ZYG-27	Titan Grad 5 (ASTM F136)	Gelb eloxiert	1
D-ZYG-27S	Titan Grad 5 (ASTM F136)	Gelb eloxiert	1
D-ZYG-27ST-GSM	Titan Grad 5 (ASTM F136)	Gelb eloxiert	1
D-ZYG-27ST-GSL	Titan Grad 5 (ASTM F136)	Gelb eloxiert	1
D-ZYG-29	Titan Grad 5 (ASTM F136)	Gelb eloxiert	1
D-ZYG-29S	Titan Grad 5 (ASTM F136)	Gelb eloxiert	1
D-ZYG-35	Titan Grad 5 (ASTM F136)	Gelb eloxiert	1
D-ZYG-35S	Titan Grad 5 (ASTM F136)	Gelb eloxiert	1
Geführte Chirurgie			
D-20T-GS-20	Titan Grad 5 (ASTM F136)	–	1
D-20T-GS-23	Titan Grad 5 (ASTM F136)	–	1
D-20T-GS-28	Titan Grad 5 (ASTM F136)	–	1
D-28T-GS-20	Titan Grad 5 (ASTM F136)	–	1
D-28T-GS-23	Titan Grad 5 (ASTM F136)	–	1
D-28T-GS-28	Titan Grad 5 (ASTM F136)	–	1
FIRST® ROHNER TECHNIK			
D-20T-32RT	Titan Grad 5 (ASTM F136)	–	1
D-29T-32RT	Titan Grad 5 (ASTM F136)	–	1
D-33T-32RT	Titan Grad 5 (ASTM F136)	–	1
D-35T-32RT	Titan Grad 5 (ASTM F136)	–	1
D-40T-32RT	Titan Grad 5 (ASTM F136)	–	1
IE Implantate			
D-20E-03F	Titan Grad 5 (ASTM F136)	–	1
D-20E-04F	Titan Grad 5 (ASTM F136)	–	1
D-20E-06F	Titan Grad 5 (ASTM F136)	–	1

D-30E-03F	Titan Grad 5 (ASTM F136)	–	1
D-30E-04F	Titan Grad 5 (ASTM F136)	–	1
D-30E-06F	Titan Grad 5 (ASTM F136)	–	1
Kortikale Perforatoren			
D-GS-CP	Titan Grad 5 (ASTM F136)	–	1
D-GS-CP4	Titan Grad 5 (ASTM F136)	–	1

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