SOUTHERNIMPLANTS®

Innovative Treatment Solutions

English	INSTRUCTIONS FOR USE: Drill Extension
Español	INSTRUCCIONES DE USO: Extensión de la fresa
Italiano	ISTRUZIONI PER L'USO: Prolunga per punte
Français	MODE D'EMPLOI: Prolongateur de foret
Deutsch	GEBRAUCHSANWEISUNGEN: Bohrer Verlängerung
Português	INSTRUÇÕES DE UTILIZAÇÃO: Extensão da broca
Ελληνικά	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ: Επέκταση τρυπανιού
Svenska	ANVISNINGAR FÖR ANVÄNDNING: Borrörlängare

ARCHIVED

South Africa - Headquarters: 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

Subsidiaries

Australia

Southern Implants Australia
T: +61-(0)-8-9466-2627
E: info@southernimplants.com.au

Spain and Portugal

Southern Implants Iberica
T: +34 935 053 507
E: info@southernimplants.es

United Kingdom and Ireland

Southern Implants UK
T: +44-20-8899-6845 / 6 / 7
E: info@southernimplants.co.uk

USA and Canada

Southern Implants North America Inc.
T: +1-561-472-0990
E: customer@southernimplants.com

EC REP

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

CH REP

MedEnvoy Switzerland: Gotthardstrasse 28, 6302 Zug, Switzerland

Intended use

The drill extension is used to increase the length of a latch type drill, providing sufficient clearance for the dental handpiece head from the surrounding anatomical structures.

Description

The drill extensions are made of hardened stainless steel. The drill extensions have a latch dimension compatible to ISO 1797. This is in order to connect the drill extension to the handpiece of an implant motor unit. The drill extensions have an internal and external W&H hex connection with dimensions compatible to ISO 1797. This in order to transmit higher torque values when necessary. The drill extensions are provided non-sterile and are multiple use instruments.

Table 1 - Drill extensions

Item code	Material	Description of product
I-DE-G	Stainless Steel	Drill Extension - Ganter (W&H compatible)
I-DE-GS4*	Stainless Steel	Drill Extension - Guided Surgery Tool
I-DE-K	Stainless Steel	Drill Extension - Komet
I-DE-MN	Stainless Steel	Drill Extension - Narrow variation
VXTR	Stainless Steel	Drill Extension - Versah variation (W&H compatible)

Indications for use

Dental implant therapy is indicated for patients with one or more missing teeth that wish to have the tooth/teeth restored with a dental implant and a prosthetic tooth/teeth.

Intended user

Maxillofacial Surgeons, General Dentists, Orthodontists, Periodontists, Prosthodontists and other appropriately trained and experienced implant users.

Intended environment

The cylindrical twist drills and cylindrical bone taps are intended to be used in a clinical environment such as an operating theater or a dentist consultation room.

Intended patient population

This device is used in the dental restoration of partially or fully edentulous patients in the upper or lower jaw. Restorations may comprise of single teeth, partial or full bridges and may be fixed or removable.

Clinical procedure

A proper clinical and radiological evaluation must be done to determine the bone dimensions and bone quality. Ensure that all instruments and drills are in a good condition.

Clinical benefits

Implant therapy results in restoration of chewing function, speech and aesthetics in patients with missing teeth.

Surgical procedure

The latch type drill/instrument should be placed into the drill extension and needs to be fully engaged. The flat side of the latch of the latch type drill/instrument needs to be aligned with the latch flat of the drill extension and the drill extension clip needs to be inside the groove of the drill/instrument's latch. See Figure 1 for correct alignment. If the latch type drill/instrument is not fully engaged, it can cause distortion of the drill/instrument's latch and the drill extensions' connection.

It is important that the:

- drill extensions must NOT be used with Ø6 mm and larger drills, use longer shaft drills instead.
- drill extensions must NOT be used with bone taps.


Figure 1: latch type instrument and drill extension showing the drill and extension alignment.

The orientations indicated in (Figure 1) ensure that the catch feature of the drill extension (circled in blue) slots into the latch groove of the drill (circled in green). This prevents the drill from sliding out of the drill extension.

NOTE:

- do not apply more than 40-45 Ncm to any latch type drill/instrument, this could cause damage to the handpiece and latch of the instrument.
- blunt drills cause excessive torque and result in damage to the handpiece or drill latch.

Before surgery

All components, instruments and tooling used during the clinical or laboratory procedure must be maintained in good condition and care must be taken that instrumentation does not damage implants or other components.

During surgery

Care must be taken that parts are not swallowed during any of the procedures, a rubber-dam application is recommended when appropriate. All components, instruments and tooling used during the clinical or laboratory procedure must be maintained in good condition and care must be taken that instrumentation does not damage implants or other components or equipment.

Post-surgery

Regular patient follow-up, and proper oral hygiene must be achieved to ensure favourable long-term results.

Storage, cleaning and sterilisation

These devices are supplied non-sterile. If packaging is damaged do not use the product and contact your Southern representative or return to Southern Implants. The devices must be stored in a dry place at room temperature and not exposed to direct sunlight. Incorrect storage may influence device characteristics.

Sterilisation

Southern Implants recommends the following procedure to sterilise the instruments prior to use/re-use:

Methods to sterilise these devices:

1. prevacuum sterilisation method: steam sterilise the instruments at 132°C (270°F) at 180-220 kPa for 4 minutes. Dry for at least 20 minutes in the chamber. Only an approved wrap or pouch for steam sterilisation must be used.
2. prevacuum sterilisation method: wrapped, steam sterilise at 135°C (275°F) for 3 minutes. Dry for 20 minutes in the chamber. Use a wrap or pouch that is cleared for the indicated steam sterilisation cycle.

NOTE: users in the USA must ensure that the steriliser, wrap or pouch and all steriliser accessories are cleared by the FDA, for the intended sterilisation cycle.

Contraindications

Do not use in patients:

- who are medically unfit for dental implant procedures.
- where adequate numbers of implants could not be placed to achieve full functional support of the prosthesis.
- who are allergic or have hypersensitivity to pure titanium or titanium alloy (Ti-6Al-4V), gold, palladium, platinum or iridium.
- who are under the age of 18, have poor bone quality, blood

disorders, infected implant site, vascular impairment, uncontrolled diabetes, drug or alcohol abuse, chronic high dose steroid therapy, anti-coagulant therapy, metabolic bone disease, radiotherapy treatment and sinus pathology.

Warnings

- The use of non-sterile items can lead to secondary infections of the tissue or transfer infectious diseases.
- Failure to perform appropriate cleaning, resterilisation and storage procedures as per the Instructions for Use document can result in device damage or secondary infection.
- For the safe and effective use of dental implants, it is recommended that specialised training be undertaken.
- Failure to accurately estimate lengths of drills relative to radiographic measurements can result in permanent injury to nerves or other vital structures.
- Use of the device with incompatible or non-corresponding devices can result in device failure or poor performance.

Cautions

New and experienced implant users should do training before using a new system or attempting to do a new treatment method. Take special care when treating patients who have local or systemic factors that could affect the healing of the bone and soft tissue (i.e. poor oral hygiene, uncontrolled diabetes, are on steroid therapy, smokers, infection in the nearby bone and patients who had orofacial radiotherapy).

Thorough screening of prospective implant candidates must be performed including:

- a comprehensive medical and dental history.
- visual and radiological inspection to determine adequate bone dimensions, anatomical landmarks, occlusal conditions and periodontal health.
- bruxism and unfavourable jaw relations must be taken into account.
- proper preoperative planning with a good team approach between well-trained surgeons, restorative dentists and lab technicians is essential for successful implant treatment.
- minimising the trauma to the host tissue increases the potential for successful osseointegration.
- electrosurgery should not be attempted around metal implants, as they are conductive.

Reusable devices

Prior to reusing this device, it needs to be inspected, if there are signs of visible corrosion, deformed or twisted connections, dull cutting edges, wear and damage, this device shall be disposed of. After inspection, and if re-use seems fit, the devices are cleaned and sterilised. If reuse seems fit:

- containment: as soon as practically possible, remove all visible residue after use (bone, blood or tissue), by immersing the instrument in cold water (dried soil is difficult to remove).
- pre-cleaning: rinse with lukewarm water for 3 minutes and remove hardened debris with a soft nylon brush. Avoid mechanical damage during cleaning.
- manual cleaning or automated cleaning: prepare an ultrasonic bath with suitable detergent, sonicate for 20 minutes (alternative methods can be used if proven by the end user). Rinse with purified/sterile water. Load devices into a thermo-disinfector. Run the cleaning and disinfection cycle, followed by the drying cycle.

NOTE: always follow the instructions for use of the manufacturers of cleaning agents and disinfectants.

- Drying: dry the instruments with filtered compressed air or single use lint-free wipes. Pack the instruments as quickly as possible after removal into the storage container. If additional drying is necessary, dry in a clean location. Moisture on these devices can cause corrosion and deterioration of the cutting edges.
- Inspection: do a visual inspection of the items to check for any damage.
- Packaging: use the correct packaging material as indicated for steam sterilisation to ensure sterility is maintained. Double packaging is recommended.

Side effects

Potential side effects and temporary symptoms: pain, swelling, phonetic difficulties, gingival inflammation. More persistent symptoms: the risks and complications with implants include, but are not limited to: (1) allergic reaction(s) to implant and/or abutment material; (2) breakage of the implant and/or abutment; (3) loosening of the abutment screw and/or retaining screw; (4) infection requiring revision of the dental implant; (5) nerve damage that could cause permanent weakness, numbness, or pain; (6) histologic responses possibly involving macrophages and/or fibroblasts; (7) formation of fat emboli; (8) loosening of the implant requiring revision surgery; (9) perforation of the maxillary sinus; (10) perforation of the labial and lingual plates; and (11) bone loss possibly resulting in revision or removal.

Breakage

Breakage or deformation can result when applied loads exceed the strength of the material. Overloading may result from: excessive pressure exerted on the handpiece; excessive lateral forces due to non-concentric loading; excessive torque; incorrect fit of the latch attachment into the handpiece; use outside of manufacturer recommendations; use of the device with incompatible hardware. In the above cases, additional treatment may be necessary to reduce the likelihood of complications.

Changes in performance

It is the responsibility of the clinician to ensure all warnings, precautions and recommended techniques are adhered to. In case of unexpected or unusual changes in performance, the manufacturer should be contacted immediately.

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established. The contact information for the manufacturer of this device to report a serious incident is as follows: sicomplaints@southernimplants.com

Materials

Drill Extension: Stainless Steel

Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

Disclaimer of liability

This product is part of the Southern Implants product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants product range and take full responsibility for the correct indications and use of this product. Southern Implants does not assume liability for damage due to incorrect use. Please note that some Southern Implants products may not be cleared or released for sale in all markets.

Basic UDI

Product	Basic-UDI Number
Basic-UDI for Drills and Hand Piece Devices	600954403875
Basic-UDI for Reusable Instruments	600954403876

Related Literature and Catalogues

CAT-2004 - Tri-Nex Implants Product Catalogue
 CAT-2005 - IT Implants Product Catalogue
 CAT-2020 - External Hex Implants Product Catalogue
 CAT-2042 - Deep Conical Implants Product Catalogue
 CAT-2043 - Internal Hex Implants Product Catalogue
 CAT-2060 - PROVATA® Implants Product Catalogue
 CAT-2069 - INVERTA® Implants Product Catalogue
 CAT-2070 - Zygomatic Implants Product Catalogue

ARCHIVED

Symbols and Warnings

 Manufacturer: Southern Implants 1 Albert Rd, P.O Box 605 IRENE, 0062, South Africa. Tel: +27 12 667 1046	 2797	 Prescription device*	 Non-sterile	 Consult instruction for use	 Batch code	 Do not use if package is damaged	 Medical Device	 Date of Manufacture	 Authorize Representative in the European Community	 Catalog number	 Authorized Representative for Switzerland
* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.											
All rights reserved. Southern Implants®, the Southern Implants logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale.											

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

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