
English	INSTRUCTIONS FOR USE: Southern Implants® Gold abutments
Español	INSTRUCCIONES DE USO: Pilares de oro Southern Implants®
Italiano	ISTRUZIONI PER L'USO: Southern Implants® Monconi in oro
Français	MODE D'EMPLOI : Piliers en or Southern Implants®
Deutsch	GEBRAUCHSANWEISUNG: Southern Implants® Gold-Abutments
Português	INSTRUÇÕES DE UTILIZAÇÃO: Pilares de ouro da Southern Implants®

South Africa - Headquarters: 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

Subsidiaries

Australia

Southern Implants Australia
T: +61-(0)-8-9466-2627
E: info@southernimplants.com.au

Spain and Portugal

Southern Implants Iberica
T: +34 935 053 507
E: info@southernimplants.es

United Kingdom and Ireland

Southern Implants UK
T: +44-20-8899-6845 / 6 / 7
E: info@southernimplants.co.uk

USA and Canada

Southern Implants North America Inc.
T: +1-561-472-0990
E: customer care@southernimplants.com

Intended use

Gold abutments are intended to be connected to an implanted dental implant to serve as support for a dental prosthesis.

Description

The gold abutments are pre-manufactured dental abutments that can either be connected direct to an endosseous implant or connect the prosthesis to a compact conical abutment for use as an aid in prosthetic rehabilitation. The gold abutment includes a plastic burn-out sleeve to aid in the laboratory procedure during wax-up procedures. Refer to individual product catalogues for product characteristics and compatible accessories. Non-engaging versions are indicated for multi-implant cases. Engaging versions are indicated for single implant cases and multi-unit custom abutment cases.

Indications for use

Southern Implants' dental implant systems are intended for surgical placement in the upper or lower jaw to provide a means for prosthetic attachment of crowns, bridges or overdentures utilizing delayed or immediate loading. Southern Implants' dental implants are intended for immediate function when good primary stability with appropriate occlusal loading is achieved.

Intended user

Dental Technicians, Maxillo-facial Surgeons, General Dentists, Orthodontists, Periodontists, Prosthodontists and other appropriately trained and experienced implant users.

Intended environment

The devices are intended to be used in a clinical environment such as an operating theater or a dentist consultation room.

Intended patient population

Patients that have lost one tooth or multiple teeth.

Compatibility information

Southern Implants' implants should be restored with Southern Implants' components. In the Southern Implants range of gold abutments, there are 5 implant connections and abutment-level options. The implant code and connection type can be identified by specific abbreviations in the product codes. Range identifiers are summarised in Table A.

Table A

Implant connection type	Compatible device
Internal Hex and PROVATA® (M/Z)	Parts labelled GC-EM or GC-EZ for engaging items
	Parts labelled GC-NM or GC-NZ for non-engaging items
External Hex (EX)	Parts labelled GC-EX-(Ø) for engaging items
	Parts labelled GC-NX-(Ø) for non-engaging items
Tri-Nex® (EL) (Lobe)	Parts labelled GC-EL-(Ø) for engaging items
	Parts labelled GC-NL-(Ø) for non-engaging items
Deep Conical (DC)	Parts labelled GC-DC-(Ø) for engaging items
	Parts labelled GC-NDC-(Ø) for non-engaging items
Internal Octagon (ITS) (ITS6)	Parts labelled ITS-GC1 (used with Ø4.8 mm platforms) for engaging items
	Parts labelled ITS-GC1ne (used with Ø4.8 mm platforms) for non-engaging items
	Parts labelled ITS6-GC1 (used with Ø6.5 mm platforms) for engaging items
	Parts labelled ITS6-GC1ne (used with Ø6.5 mm platforms) for non-engaging items
Abutment level	Parts labelled GMC1 (used with Ø4.8 mm abutment platforms) for non-engaging items
	Parts labelled GMCW1 (used with Ø6.0 mm abutment platforms) for non-engaging items

Laboratory procedures

This abutment has a machined gold base which fits directly to the implant or compact conical abutment, together with a waxing sleeve, to facilitate the wax up process.

NOTE: the gold abutment is threaded and the plastic waxing sleeve screws onto the abutment, ensure the two components are screwed together at all times.

1. Place the abutment onto the implant analogue, making sure that the retentive elements of the analogue-abutment connection are properly aligned. Fix the abutment by hand tightening the lab-screw to the analogue.
2. Trim the plastic sleeve to the correct occlusal height, with no additional angular correction and out of occlusion (retaining a minimal height of 4 mm for single units). Fabricate a cement-retained or screw-retained restoration. Add wax to achieve the desired substructure.

Investment

1. Clean the abutment before investing, do not cover the abutment margins in wax, do not use wetting agents.
2. When investing, avoid the use of a tension reliever as it leaves a residue, which can result in metal creep.
3. Sprue the restoration, making sure the screw channel is vertical to the investment ring.
4. The use of investment material for rapid heating methods (speed investment methods) is not recommended.
5. A phosphate bonded investment must be used.

WARNING: cast onto cylinder using precious metal.

Casting procedure

Gold abutment casting specifications:

1. The gold (Ceramicor®) abutment is not a porcelain-bonding alloy. The alloy is designed to have oxidation during casting-on procedures. Direct porcelain bonding could result in the porcelain turning green or chipping off due to a poor bond.
2. The Co-efficient of Expansion of Southern Implants Gold abutments is: 25-500°C:11.9x10⁻⁶ K⁻¹ and 25-600°C:12.2x10⁻⁶ K⁻¹
3. The melting temperature of the Gold abutment is approximately 1475°C.
4. The recommended cast on temperature is ± 920°C.
5. Casting at higher temperatures could result in distortion of the abutment.
6. Casting at a lower temperature of 650°C may result in a miscast.
7. The furnace should run up straight to the casting temperature to reduce investment expansion. Too much expansion could result in overflow of metal onto the fitting surface of the abutments.

Recommended casting alloys:

Contact your restorative material supplier for a compatible material to Southern Implants Gold abutments.

Commonly used metals for custom posts (non-ceramic bonding) are: Stabilor, Procast Y45-Argen, Argeno 1.

Commonly used metals for custom posts (ceramic bonding) are: Argeno 1, Degudent U, Argedent 3, Degunorm – Gold gate system, Bond-on 4, and Degudent G.

After casting procedure

1. After casting, the removal of the investment material is most critical. It must be removed ultrasonically. Do not sandblast or blast the fitting surface with glass beads. If this surface is blasted in any way, the precision-machined gold part may not fit the implant as intended. If the abutments are to be sandblasted cover it with a thin layer of wax or use a polishing protector cap.
2. Work off the restoration as normal. Take care not to damage the fitting surface.
3. Before sending the restoration and model to the clinician, clean the model thoroughly with soap, water and a toothbrush or steam clean. Then clean the abutment crown and its screw, preferably in the ultrasonic bath. Assemble the crown to the model and send it to the clinician.
4. Occlusion: implant borne crowns are usually kept out of occlusion by about one-tenth of a millimetre.

Clinical procedures

The clinician receives the restoration from the laboratory.

1. Remove the healing abutment or temporary restoration.
2. Clean, disinfect and sterilise the restoration as described above.
3. Insert the restoration into the patient's mouth.
4. Position the restoration on the implant making sure that the retentive elements of the implant/abutment connections are properly aligned.
5. Fix the abutment to the implant/abutment with the correct screw using applicable driver (Table B). Torque the screw down to the value indicated in Table C.

Table B

Driver type	External Hex	DC	Tri-Nex	Internal Hex	IT	Compact Conical screw
1.22 mm / 1.27 mm Universal driver	✓	✓		✓		✓
1.22 mm hex driver	✓	✓				✓
1.27 mm hex driver				✓		
Unigrip driver	✓		✓			✓
Quad driver	✓			Gold screws only		
Blade driver	✓					✓
Torx driver					✓	

Table C

Direct to Implant	Torque
External Hex	
Ø3.0 mm	25-32 Ncm
Ø3.25, 4.0, 5.0, 6.0, 7.0, 8.0 and 9.0 mm	32-40 Ncm
Tri-Nex®	
Ø3.5 mm	32 Ncm
Ø4.3, 5.0, 6.0, 7.0, 8.0 and 9.0 mm	32-40 Ncm
Deep Conical	
Ø3.0 mm	15 Ncm
Ø3.5 and 4.0 mm	20 Ncm
Ø5.0 mm	25-32 Ncm
Internal Hex (M-Series and PROVATA®)	
Ø3.75, 4.2 and 5.0 mm M-Series	32 Ncm
Ø4.0, 5.0, 6.0, 7.0, 8.0 and 9.0 mm PROVATA® Implant	32 Ncm
IT Octagon	
Ø3.3, 4.1, 4.9, 4.0, 5.0, 6.0, 7.0, 8.0 and 9.0 mm	32-40 Ncm
Abutment level	All gold cylinders on Compact Conical Abutments: 10-15 Ncm

6. Verify the correct seating of the restoration using radiographic imaging.
7. Do not exceed the recommended torque value as this may result in failure of the screw, abutment or implant. Do not torque less than the recommended value, this may result in loosening of the abutment that can lead to abutment or implant failure.
8. Close the screw access hole.
9. Cement the temporary prosthesis if applicable.

Clinical benefits

Patients can expect to have their missing teeth replaced and/or crowns restored.

Before surgery

All components, instruments and tooling used during the clinical or laboratory procedure must be maintained in good condition and care must be taken that instrumentation does not damage implants or other components.

During surgery

Care must be taken that parts are not swallowed during any of the procedures, a rubber-dam application is recommended when appropriate. Care must be taken to apply the correct tightening torque of abutments and abutment screws.

Post-surgery

Regular patient follow-up and proper oral hygiene must be achieved to ensure favourable long-term results.

Storage, cleaning and sterilisation

The implants, cover screws and healing abutments are supplied sterile (sterilised by gamma irradiation) and intended for single-use prior to the expiration date (see packaging label). Sterility is assured unless the container or seal is damaged or opened. If packaging is damaged, do not use the product and contact your Southern representative or return to Southern Implants. Do not reuse implants, cover screws, temporary abutments or abutments. Re-using these components may result in:

- damage to the surface or critical dimensions, which may result in performance and compatibility degradation.
- adds the risk of cross-patient infection and contamination if single-use items are reused.

Southern Implants does not accept any responsibility for complications associated with reused components.

Cleaning and disinfection

An implant restoration is a single or multiple-tooth implant crown, bridge or substructure, attached to a Southern Implants abutment or multiple abutments. Before intraoral use the final restoration needs to be cleaned and disinfected, as per restorative material manufacturer's instructions.

Sterilisation

Southern Implants recommends the following procedure to sterilise the restoration prior to use:

Methods to sterilise the restoration and abutment screw:

1. Pre-vacuum sterilisation method: steam sterilise the abutments at 132°C (270°F) at 180-220 kPa for 4 minutes. Dry for at least 20 minutes in the chamber. Only an approved wrap or pouch for steam sterilisation must be used.
2. Pre-vacuum sterilisation method: wrapped, steam sterilise at 135°C (275°F) for 3 minutes. Dry for 20 minutes in the chamber. Use a wrap or pouch that is cleared for the indicated steam sterilisation cycle.

NOTE: users in the USA must ensure that the steriliser, wrap or pouch and all steriliser accessories are cleared by the FDA, for the intended sterilisation cycle.

Contraindications

Do not use in patients:

- who are medically unfit for dental implant procedures.
- where adequate numbers of implants could not be placed to achieve full functional support of the prosthesis.
- who are allergic or have hypersensitivity to pure titanium or titanium alloy (Ti-6Al-4V), gold, palladium, platinum or iridium.
- who are under the age of 18, have poor bone quality, blood disorders, infected implant site, vascular impairment, uncontrolled diabetes, drug or alcohol abuse, chronic high dose steroid therapy, anti-coagulant therapy, metabolic bone disease, radiotherapy treatment or sinus pathology.

Warnings

THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING.

- For the safe and effective use of dental implants it is suggested that specialised training be undertaken, including hands-on training to learn proper technique, biomechanical requirements and radiographic evaluations.
- Responsibility for proper patient selection, adequate training, experience in the placement of implants and providing appropriate information for informed consent rests with the practitioner. Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone.

Cautions

New and experienced implant users should do training before using a new system or attempting to do a new treatment method. Take special care when treating patients who have local or systemic factors that could affect the healing of the bone and soft tissue (i.e. poor oral hygiene, uncontrolled diabetes, are on steroid therapy, smokers, infection in the nearby bone and patients who had oro-facial radiotherapy). Thorough screening of prospective implant candidates must be performed including:

- a comprehensive medical and dental history.
- visual and radiological inspection to determine adequate bone dimensions, anatomical landmarks, occlusal conditions and periodontal health.
- bruxism and unfavourable jaw relations must be taken into account.
- proper pre-operative planning with a good team approach between well trained surgeons, restorative dentists and lab technicians is essential for successful implant treatment.
- minimizing the trauma to the host tissue increases the potential for successful osseointegration.
- electro-surgery should not be attempted around metal implants, as they are conductive.

Single use devices

Do not reuse devices indicated for single use (use the device prior to the expiration date).

Side effects

Potential side effects and temporary symptoms: pain, swelling, phonetic difficulties, gingival inflammation. More persistent symptoms: the risks and complications with implants include, but are not limited to: (1) allergic reaction(s) to implant and/or abutment material; (2) breakage of the implant and/or abutment; (3) loosening of the abutment screw and/or retaining screw; (4) infection requiring revision of the dental implant; (5) nerve damage that could cause permanent weakness, numbness, or pain; (6) histologic responses possibly involving macrophages and/or fibroblasts; (7) formation of fat emboli; (8) loosening of the implant requiring revision surgery; (9) perforation of the maxillary sinus; (10) perforation of the labial and lingual plates; and (11) bone loss possibly resulting in revision or removal.

Breakage

Implant and abutment fractures can occur when applied loads exceed the tensile or compressive strength of the material. Potential overloading conditions may result from: deficiencies in implant numbers, lengths and/or diameters to adequately support a restoration, excessive cantilever length, incomplete abutment seating, abutment angles greater than 30°, occlusal interferences causing excessive lateral forces, patient parafunction (e.g. bruxing, clenching), loss or changes in dentition or functionality, inadequate prosthesis fit and physical trauma. Additional treatment may be necessary when any of the above conditions are present to reduce the possibility of hardware complications or failure.

Changes in performance

It is the responsibility of the clinician to instruct the patient on all appropriate contraindications, side effects and precautions as well as the need to seek the services of a trained dental professional if there are any changes in the performance of the implant (e.g. looseness of the prosthesis, infection or exudate around the implant, pain or any other unusual symptoms that the patient has not been told to expect).

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established. The contact information for the manufacturer of this device to report a serious incident is as follows: sicomplaints@southernimplants.com

Materials

Gold abutment:	Gold alloy (Au-60%, Pd-20%, Pt-19%, Ir-1%)
Plastic sleeve:	Polyoxymethylene (POM)
Abutment screws:	Titanium alloy (Ti-90%, Al-6%, V-4%) Gold Alloy (Au-61%, Ag-16.5%, Pt-13.5%, Cu-9%)

Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

Disclaimer of liability

This product is part of the Southern Implants product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants product range and take full responsibility for the correct indications and use of this product. Southern Implants does not assume liability for damage due to incorrect use. Please note that some Southern Implants products may not be cleared or released for sale in all markets.

Healing

The healing time required for osseointegration depends on the individual and treatment protocol. It is the responsibility of the practitioner to decide when the implant can be restored. Good primary stability will govern if immediate loading can be done.

Implant care and maintenance

Potential implant patients should establish an adequate oral hygiene regime prior to implant therapy. Proper post-operative oral hygiene and implant maintenance instructions must be discussed with the patient as this will determine the longevity and health of the implants. The patient should maintain regular prophylaxis and evaluation appointments.

MR conditional

Non-clinical testing has demonstrated that the Southern Implants dental implants, metallic abutments and prosthetic screws are MR conditional. A patient with these devices can be safely scanned in a MR system meeting the following conditions:

- Static magnetic field of 1.5 Tesla and 3.0 Tesla only.
- Maximum spatial gradient magnetic field of 4000 Gauss/cm (40 T/m).
- Maximum MR system reported, head specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode) or whole body averaged specific absorption rate (wbSAR) of 1 W/kg.

Under the scan conditions defined above, the Southern Implants dental implants, abutments and prosthetic screws are expected to produce a maximum temperature rise of 5.8°C after 15 minutes of continuous scanning. In the non-clinical testing, the image artifact caused by the device extends approximately 20 mm from the Southern Implants dental implants, abutments and prosthetic screws, when imaged with a gradient echo pulse sequence and a 3.0 Tesla MRI system. Removable restorations should be taken out prior to scanning, as is done for watches, jewellery etc. Should there be no MR symbol on the product label, please note that this device has not been evaluated for safety and compatibility in the MR environment. This device has not been tested for heating or migration in the MR environment.

Basic UDI

Product	Basic-UDI Number
Basic-UDI for Metal Abutments	600954403872

Related literature and catalogues

CAT-2004 - Tri-Nex® Implants Product Catalogue
 CAT-2005 - IT Product Catalogue
 CAT-2020 - External Hex Implants Product Catalogue
 CAT-2042 - Deep Conical Implants Product Catalogue
 CAT-2043 - Internal Hex Implants Product Catalogue
 CAT-2060 - PROVATA® Implants Product Catalogue
 CAT-2069 - INVERTA® Implants Product Catalogue
 CAT-2070 - Zygomatic Implants Product Catalogue

Symbols and Warnings

 Manufacturer: Southern Implants 1 Albert Rd, P.O Box 605 IRENE, 0062, South Africa. Tel: +27 12 667 1046	  2797 Prescription device*	 Consult instruction for use	 Do not reuse	 Batch code	 Do not use if package is damaged	 Medical Device	 Date of Manufacture	 Authorize Representative in the European Community	 Catalog number
* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.					Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.				
All rights reserved. Southern Implants®, the Southern Implants logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale.									

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