English	INSTRUCTIONS FOR USE: Southern Implants® Impression Copings
Español	INSTRUCCIONES DE USO: Cofias de impresión Southern Implants®
Italiano	ISTRUZIONI PER L'USO: Southern Implants® Cappette per Impronta
Français	NOTICE D'UTILISATION : Transferts d'empreintes de Southern Implants®
Deutsch	GEBRAUCHSANWEISUNG: Southern Implants® Abdruckpfosten
Português	INSTRUÇÕES DE UTILIZAÇÃO: Southern Implants® Coifas de impressão
Ελληνικά	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ: Southern Implants® Επικαλύψεις αποτύπωσης
Svenska	BRUKSANVISNING: Southern Implants® Avtryckskåpan
Română	INSTRUCȚIUNI DE UTILIZARE: Southern Implants® Pini de amprentare
Türkçe	KULLANMA TALİMATI: Southern Implants® İmpresyon Kopingleri

South Africa - Headquarters: 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

Subsidiaries

Australia

Southern Implants Australia
T: +61-(0)-8-9466-2627
E: info@southernimplants.com.au

Spain and Portugal

Southern Implants Iberica
T: +34 935 053 507
E: info@southernimplants.es

United Kingdom and Ireland

Southern Implants UK
T: +44-20 8059 4490
E: info@southernimplants.co.uk

USA and Canada

Southern Implants North America Inc.
T: +1-561-472-0990
E: customercare@southernimplants.com

Description

Southern Implants® Impression copings are premanufactured components that are connected directly to an endosseous implant and are used in single tooth or multiple unit reconstructions to take an impression of the patient's oral cavity. Impression copings come in a transfer (closed tray) and pickup (open tray) type configuration. These devices are supplied sterile and for single patient use.

Intended use

The Southern Implants® Impression Copings and Scan Flags are intended to facilitate the transfer of the orientation and position of placed endosseous implants or abutments to a lab model.

Indications for use

The Southern Implants® Impression Copings and Scan Flags are indicated for the prosthetic planning, design and manufacture.

Specifically, the Southern Implants® Impression Copings are indicated for manual impression procedures.

Intended user

Dental Technicians, Maxillofacial Surgeons, General Dentists, Orthodontists, Periodontists, Prosthodontists, and other appropriately trained and/or experienced medical professionals.

Intended environment

The devices are intended to be used in a dental laboratory as part of the restoration design and manufacture as well as in a clinical environment such as an operating theatre or a dentist consultation room.

Intended patient population

The Southern Implants® Impression Copings and Scan Flags are intended to be used in patients subject to implant-retained prosthetic restoration procedures.

Compatibility information

Southern Implants' implants should be restored with Southern Implants' components. In the Southern Implants' range there are 8 implant and/or abutment connections. The implant code and connection type can be identified by specific abbreviations in the product codes. Range identifiers are summarised in Table A.

Table A - Compatibility

Implant connection type	Compatible device
External Hex (EX)	Parts labelled CP (used with Ø3.00 mm platform) Parts labelled CBNU, CBNU-W and CBN70 (used with Ø3.25 mm platform) Parts labelled CBV, CBU, CBU-W, CB70 and CB75 (used with Ø4.00 mm platform) Parts labelled CBAV, CBAU, CBAU-M, CBAU-W and CBA75 (used with Ø5.00 mm platform) Parts labelled CBBBU, CBBB-W and CBBB75 (used with Ø6.00 mm platform) Parts labelled CMAX9 and CMAX9L (used with Ø7.00 mm platform)
TRI-NEX® (EL) (Lobe)	Parts labelled IC-L-(Ø), IC-L-(Ø)-W, ICT-L-(Ø) and ICT-L-(Ø)-W
Deep Conical (DC)	Parts labelled IC-DC(Ø), IC-DC(Ø)-W, IC-DCR(Ø), ICT-DC(Ø) and ICT-DC(Ø)-W
Internal Hex (M)	Parts labelled ICT-M, ICT-MW and IC-M (used with Ø3.75, 4.20 and 5.00 mm implants) Parts labelled ICT-MW-P45 and IC-MW-P45 (used with Ø5.00 mm implants)
Internal Hex PROVATA® (M/ Z)	Parts labelled SFT-PRO3 (used with Ø3.30 mm implants) Parts labelled ICT-M, ICT-MW and IC-M (used with Ø4.00, 5.00 and 6.00 mm implants) Parts labelled ICT-MW-P45 and IC-MW-P45 (used with Ø5.00 and 6.00 mm implants) Parts labelled ICT-Z, ICT-ZW and IC-Z (used with Ø7.00, 8.00 and 9.00 mm implants)
IT Internal Octagon (ITS/ ITS6)	Parts labelled ITS-CA5, ITS-TRC5 and ITSTRCNE (used with Ø4.80 mm platforms) Parts labelled ITS6-CA5, ITS6-TRC5 and ITS6-TRCNE (used with Ø6.50 mm platforms)
Single Platform (SP)	Parts labelled IC-SP(Ø) and ICT-SP(Ø) (used with Ø3.50, 4.00 and 5.00 mm implants) Parts labelled IC-SP(Ø)-PM and ICT-SP(Ø)-PM (used with Ø5.00 mm implants)
Compact Conical Abutments	Parts labelled CMC1, CMC2 and CMC-ZG-2 (used with Ø4.80 mm compact conical platform) Parts labelled CMCW1 and CMCW2 (used with Ø6.00 mm compact conical platform)

Clinical performance

The Impression Copings and Scan Flags are intended to aid in the prosthetic procedures associated with implant therapy. Consequently, the clinical performance of these devices is primarily defined by the success of the prosthetic outcome. This performance can be quantitatively assessed through prosthetic survival/success, which serves as an important indicator of their effectiveness and success in clinical practice.

Clinical benefits

The Southern Implants® Impression Copings and Scan Flags play a crucial role in ensuring the successful outcome of prosthetic procedures, though they do not provide a direct clinical benefit. However, the clinical benefits of the associated treatment, target implant, and/or prosthetic device(s) should be acknowledged, with relevant documentation reviewed as necessary.

Storage, cleaning and sterilisation

The component is supplied sterile (sterilised by gamma irradiation) and intended for single-use prior to the expiration date (see packaging label). Sterility is assured unless the container or seal is damaged or opened. If packaging is damaged do not use the product and contact your Southern representative or return to Southern Implants®. The devices must be stored in a dry place at room temperature and not exposed to direct sunlight. Incorrect storage may influence device characteristics. Do not reuse components indicated for single-use only. Reusing these components may result in:

- damage to the surface or critical dimensions, which may result in performance and compatibility degradation.
- adds the risk of cross-patient infection and contamination if single-use items are reused.

Southern Implants® does not accept any responsibility for complications associated with reused single-use components.

Contraindications

The contraindications of the specific target implant and prosthetic devices apply. Therefore, the contraindications of the systems/medical devices utilized as part of implant surgery/therapy should be noted and the relevant documents consulted.

Contraindications specific to the Southern Implants® Impression Copings and Scan Flags:

- Allergies or hypersensitivity to chemical ingredients in the following materials: Titanium, Aluminium, Vanadium, Acetal/polyoxymethylene and Stainless Steel.

Other than the above, there are no contraindications unique to this system.

Warnings and precautions

IMPORTANT NOTICE: THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING

- To ensure the safe and effective use of dental implants, new technologies/systems and the subject devices, it is strongly recommended that specialised training be undertaken. This training should include hands-on methods to gain competency on the proper technique, the system's biomechanical requirements and radiographic evaluations required for the specific system.
- Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone.
- Use of the device with incompatible or non-corresponding devices can result in the poor performance or failure of the device.
- When handling devices intraorally, it is imperative that they are adequately secured to prevent aspiration, as aspiration of products may lead to infection or physical injury.
- The use of non-sterile items can lead to secondary infections of the tissue or transfer infectious diseases.
- Failure to adhere to appropriate cleaning, re-sterilization, and storage procedures as outlined in the Instructions for Use can result in device damage, secondary infections, or patient harm.
- Exceeding the number of recommended uses for reusable devices can result in device damage, secondary infection or patient harm.
- The use of blunt drills may cause damage to the bone, potentially compromising osseointegration.

It is crucial to emphasize that training should be undertaken by both new and experienced implant users prior to utilizing a new system or attempting a new treatment method.

Patient Selection and Preoperative Planning

A comprehensive patient selection process and meticulous preoperative planning are essential for successful implant treatment. This process should involve consultation among a multidisciplinary team, including well-trained surgeons, restorative dentists, and laboratory technicians.

Patient screening should include, at a minimum, a thorough medical and dental history, as well as visual and radiological inspections to assess the presence of adequate bone dimensions, the positioning of anatomical landmarks, the presence of unfavourable occlusal conditions, and the periodontal health status of the patient.

For successful implant treatment, it is important to:

1. Minimize trauma to the host tissue, thereby increasing the potential for successful osseointegration.
2. Accurately identify measurements relative to radiographic data, as failure to do so may lead to complications.
3. Be vigilant in avoiding damage to vital anatomical structures, such as nerves, veins, and arteries. Injury to these structures may result in serious complications, including ocular injury, nerve damage, and excessive bleeding.

The responsibility for proper patient selection, adequate training and experience in implant placement, and the provision of appropriate information required for informed consent rests with the practitioner. By combining thorough screening of prospective implant candidates with a practitioner possessing a high level of competence in the use of the system, the potential for complications and severe side effects can be significantly reduced.

High Risk Patients

Special care should be exercised when treating patients with local or systemic risk factors that may adversely affect the healing of bone and soft tissue or otherwise increase the severity of side effects, the risk of complications, and/or the likelihood of implant failure. Such factors include:

- poor oral hygiene
- history of smoking/vaping/tobacco use
- history of periodontal disease
- history of orofacial radiotherapy**
- bruxism and unfavourable jaw relations
- use of chronic medications that may delay healing or increase the risk of complications including, but not limited to, chronic steroid therapy, anti-coagulant therapy, TNF- α blockers, bisphosphonate and cyclosporin

*** The potential for implant failure and other complications increases when implants are placed in irradiated bone, as radiotherapy can lead to progressive fibrosis of blood vessels and soft tissue (i.e., osteoradionecrosis), resulting in diminished healing capacity. Contributing factors to this increased risk include the timing of implant placement in relation to radiation therapy, the proximity of radiation exposure to the implant site, and the radiation dosage at that site.*

Should the device not operate as intended, it must be reported to the manufacturer of the device. The contact information or the manufacturer of this device to report a change in performance is:

sicomplaints@southernimplants.com.

Side effects

The clinical outcome of treatment is influenced by multiple factors. The following side effects and residual risks are associated with this device group or the surgical procedures involved in its use. These may require further treatment, revision surgery, or additional consultations with the relevant medical professional.

- Allergic or hypersensitivity reaction(s)
- Bruising
- Dental injury during surgery

- Gingival injury
- Infection (acute and/or chronic)
- Localized inflammation
- Pain, discomfort or tenderness
- Soft tissue irritation
- Surgical side effects such as pain, inflammation, bruising and mild bleeding
- Prosthetic misfit or improper connection
- Sub-optimal aesthetic result

Additionally, the normal side effects associated with anaesthesia should also be expected. Note: these side effects and residual risks can vary in both severity and frequency.

Precaution: maintaining sterility protocol

The devices are packed in a peel pouch or blister base with a “peel-back” lid. Labelling information is located on the bottom half of the pouch, inside the packet or on the surface of the peel-back lid. Sterility is assured unless the pouch or blister is damaged or opened.

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established.

The contact information for the manufacturer of this device to report a serious incident is as follows:

sicomplaints@southernimplants.com.

Materials

Impression Copings & Impression Coping Screws:

Commercially Pure Titanium Grade 3 or 4 according to ASTM F67 and ISO 5832-2, Titanium alloy Ti-6Al-4V (90% titanium, 6% aluminium, 4% vanadium) according to ASTM F136 and ISO 5832-3, Stainless Steel or Acetal (POM)

Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

Summary of Safety and Clinical Performance (SSCP)

As required by the European Medical Device Regulation (MDR; EU2017/745), a Summary of Safety and Clinical Performance (SSCP) is available for perusal with regard to Southern Implants® product ranges.

The relevant SSCP can be accessed at <https://ec.europa.eu/tools/eudamed>.

NOTE: the above website will be available upon the launch of the European Database on Medical Devices (EUDAMED).

Disclaimer of liability

This product is part of the Southern Implants® product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants® product range and take full responsibility for the correct indications and use of this product. Southern Implants® does not assume liability for damage due to incorrect use. Please note that some Southern Implants® products may not be cleared or released for sale in all markets.


Basic UDI

Product	Basic-UDI Number
Basic-UDI for Impression Copings	6009544038789J

Related literature and catalogues

CAT-2004 - TRI-NEX® Implants Product Catalogue
 CAT-2005 - IT Implants Product Catalogue
 CAT-2010 - Osseointegrated Fixtures Catalogue
 CAT-2020 - External Hex Implants Product Catalogue
 CAT-2042 - Deep Conical Implants Product Catalogue
 CAT-2043 - Internal Hex Implants Product Catalogue
 CAT-2060 - PROVATA® Implants Product Catalogue
 CAT-2069 - Deep Conical INVERTA® Implants Product Catalogue
 CAT-2070 - Zygomatic Implants Product Catalogue
 CAT-2093 - Single Platform Implants Product Catalogue
 CAT-2095 - External Hex INVERTA® Implants Product Catalogue

Symbols and warnings

* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

All rights reserved. Southern Implants®, the Southern Implants® logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants®. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale. It is the responsibility of the clinician to inspect the symbols that appear on the packaging of the product in use.

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