English	INSTRUCTIONS FOR USE: Southern Implants® Overdenture Abutments
Español	INSTRUCCIONES DE USO: Pilares de sobredentadura Southern Implants®
Italiano	ISTRUZIONI PER L'USO: Southern Implants® Monconi per Overdenture
Français	NOTICE D'UTILISATION : Piliers de prothèse Southern Implants®
Deutsch	GEBRAUCHSANWEISUNG: Southern Implants® Deckprothesen-Abutments
Português	INSTRUÇÕES DE UTILIZAÇÃO: Southern Implants® Pilares de sobredentadura
Ελληνικά	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ: Southern Implants® Κολοβώματα επένθετης οδοντοστοιχίας
Svenska	BRUKSANVISNING: Southern Implants® Täckprotes Distanser
Română	INSTRUCȚIUNI DE UTILIZARE: Southern Implants® Bonturi pentru suprastructuri dentare
Türkçe	KULLANMA TALİMATI: Southern Implants® Üst Protez Abutmentleri

Southern Implants (Pty.) Ltd. 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

EC REP

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

Subsidiaries

Australia
Southern Implants Australia
T: +61-(0)-8-9466-2627
E: info@southernimplants.com.au

Spain and Portugal
Southern Implants Iberica
T: +34 935 053 507
E: info@southernimplants.es

United Kingdom and Ireland
Southern Implants UK
T: +44-20-805-94490
E: info@southernimplants.co.uk

USA and Canada
Southern Implants North America Inc.
T: +1-561-472-0990
E: customercare@southernimplants.com

Description

Overdenture Abutments are premanufactured abutments that are connected directly to an endosseous implant and is used in multiple unit reconstructions to retain a tissue-supported overdenture in patients with extensive bone or soft tissue loss. The coronal part of the abutment has a retentive feature that allows an overdenture to clip onto the abutment so that it can be inserted and removed by the patient. These components are supplied non-sterile.

Intended use

The Southern Implants® Metal Abutments are prosthetic components intended for connection to endosseous implant systems to aid in prosthetic rehabilitation and implant procedures.

Specifically, the Southern Implants® Overdenture Abutments are intended to provide a supportive base for a removable dental prosthesis on endosseous implants.

Indications for use

The Southern Implants® Metal Abutments are indicated for use as components of an endosseous implant system to support prosthetic rehabilitation in implant treatment procedures.

Specifically, the Southern Implants® Overdenture Abutments are indicated for connection to Southern Implants® fixtures as part of the functional and aesthetic restoration of the upper or lower jaw using a fixed, removable prosthesis.

Intended user

Dental technicians, Maxillofacial Surgeons, General Dentists, Orthodontists, Periodontists, Prosthodontists, and other appropriately trained and/or experienced medical professionals.

Intended environment

The devices are intended to be used in a dental laboratory as part of the restoration design and manufacture as well as in a clinical environment such as an operating theatre or a dentist consultation room.

Intended patient population

Southern Implants® Metal Abutments are intended to be used in patients indicated and eligible, or otherwise not contraindicated, for implant-retained prosthetic restorations.

Compatibility information

Southern Implants' implants should be restored with Southern components. In the Southern Implants range there are 7 implant connections, the implant code and connection type, which can be identified by specific abbreviations in the product codes. Range identifiers are summarised in Table A.

Table A

Implant connection type	Compatible device
External Hex (EX)	Parts labelled OAB(*), OAB-H(*), OABN-H(*), OBA(*), OBZ(*), EQ-IP- (*), EQ-IBN- (*) and EQ-IB- (*)
Tri-Nex® (EL) (Lobe)	Parts labelled Z-EQ-131TN35(*), Z-EQ-131TN43(*), Z-EQ-131TN5(*), LOB-35-(*), LOB-43-(*), LOB-50-(*), and LOB-60-(*)
Deep Conical (DC)	Parts labelled EQ-DC(ø)-(L)
Internal Hex (M)	Parts labelled EQ-M-(*) (used with ø3.75, 4.2 and 5.0mm platforms)
Internal Hex PROVATA® (M) (Z)	Parts labelled EQ-3M-(*) (used with Ø3.3 mm platform)
	Parts labelled EQ-M-(*) (used with ø4.0, 5.0 and 6.0mm platforms)
IT (ITS) (ITS6) - Internal Octagon	Parts labelled TOD(*), TOD6-(*), and Z-EQ-131IT48(*) (used with ø4.8mm platforms)
Single Platform (SP)	Parts labelled EQ-SP-(*) (used with Ø3.5, 3.75, 4.0, 4.5, 5.0 and 6.0 mm implants)

NOTE: (*) is indicative of various lengths available.

NOTE: TRI-NEX® (Lobe) and IT (Octagon) implant equator abutments are supplied with a smart box housing and retention caps. For External Hex, Internal Hex (M-Series and Provata®) and DC (Deep Conical) and Single Platform (SP1) ranges the housing and retentive caps needs to be ordered separately.

Clinical performance

The clinical performance of the Southern Implants® Overdenture Abutments is primarily defined by their ability to maintain stability and retention of the prosthetic restoration over time. This stability is achieved through the secure fixation of the restoration to the abutment, which can be accomplished via cementation or casting, as well as the mechanical fastening of the abutment to the implant through screws.

Clinical benefits

The clinical benefits associated with the Southern Implants® Overdenture Abutments are intrinsically linked to the overall success of implant-supported prosthetic treatments and encompass both direct and indirect clinical benefits. The direct clinical benefits of the Metal Abutments arise specifically from their use in the treatment process. These benefits include the enhancement of peri-implant soft tissue health and the aesthetic outcomes of the final restoration. The clinical benefits experienced by patients as a result of the implant-supported prosthetic treatments in which these abutments are utilized represent the indirect clinical benefits of the Southern Implants® Overdenture Abutments. While these benefits are substantial, they are not directly attributable to the devices themselves; rather, they reflect the contributions of the abutments to the overall treatment system. Additionally, successful treatment outcomes are anticipated to result in significant enhancements in patients' psychosocial well-being, self-esteem, and overall quality of life. Overall, the Southern Implants® Overdenture Abutments not only facilitate the mechanical aspects of dental restorations but also contribute to the holistic improvement of patient outcomes and experiences.

Storage, cleaning and sterilisation

The component is supplied sterile (sterilised by gamma irradiation) and is intended for single use prior to the expiration date (see packaging label). Sterility is assured unless the container or seal is damaged or opened. If packaging is damaged do not use the product and contact your Southern representative or return to Southern Implants®. The devices must be stored in a dry place at room temperature and not exposed to direct sunlight. Incorrect storage may influence device characteristics.

Do not reuse components indicated for single-use only. Reusing these components may result in:

- damage to the surface or critical dimensions, which may result in performance and compatibility degradation.
- adds the risk of cross-patient infection and contamination if single-use items are reused.

Southern Implants® does not accept any responsibility for complications associated with reused single-use components.

Contraindications

The contraindications of all device groups used as part of the specific treatment or procedure apply. Therefore, the contraindications of the systems/medical devices utilized as part of implant surgery/therapy should be noted and the relevant documents consulted.

The contraindications specific to the Southern Implants® Overdenture Abutments include:

- patients who are medically unfit for implant or prosthetic related procedures
- patients that present with an allergy or hypersensitivity to the following materials: Titanium, Aluminium, Vanadium, Cobalt Chrome, Cobalt Chrome Molybdenum, Gold, Platinum, Palladium, Iridium, Silver and/or Copper

Warnings and precautions

IMPORTANT NOTICE: THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING

- To ensure the safe and effective use of dental implants, new technologies/systems and the Metal Abutment devices, it is strongly recommended that specialised training be undertaken. This training should include hands-on methods to gain competency on the proper technique, the system's biomechanical requirements and radiographic evaluations required for the specific system.
- Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone.
- Use of the device with incompatible or non-corresponding devices can result in the poor performance or failure of the device.
- When handling devices intraorally, it is imperative that they are adequately secured to prevent aspiration, as aspiration of products may lead to infection or physical injury.
- The use of non-sterile items can lead to secondary infections of the tissue or transfer infectious diseases.
- Failure to adhere to appropriate cleaning, re-sterilization, and storage procedures as outlined in the Instructions for Use can result in device damage, secondary infections, or patient harm.
- Exceeding the number of recommended uses for reusable devices can result in device damage, secondary infection or patient harm.
- The use of blunt drills may cause damage to the bone, potentially compromising osseointegration.

It is crucial to emphasize that training should be undertaken by both new and experienced implant users prior to utilizing a new system or attempting a new treatment method.

Patient Selection and Preoperative Planning

A comprehensive patient selection process and meticulous preoperative planning are essential for successful implant treatment. This process should involve consultation among a multidisciplinary team, including well-trained surgeons, restorative dentists, and laboratory technicians.

Patient screening should include, at a minimum, a thorough medical and dental history, as well as visual and radiological inspections to assess the presence of adequate bone dimensions, the positioning of anatomical landmarks, the presence of unfavourable occlusal conditions, and the periodontal health status of the patient.

For successful implant treatment, it is important to:

1. Minimize trauma to the host tissue, thereby increasing the potential for successful osseointegration.
2. Accurately identify measurements relative to radiographic data, as failure to do so may lead to complications.
3. Be vigilant in avoiding damage to vital anatomical structures, such as nerves, veins, and arteries. Injury to these structures may result in serious complications, including ocular injury, nerve damage, and excessive bleeding.

The responsibility for proper patient selection, adequate training and experience in implant placement, and the provision of appropriate information required for informed consent rests with the practitioner. By combining thorough screening of prospective implant candidates with a practitioner possessing a high level of competence in the use of the system, the potential for complications and severe side effects can be significantly reduced.

High Risk Patients

Special care should be exercised when treating patients with local or systemic risk factors that may adversely affect the healing of bone and soft tissue or otherwise increase the severity of side effects, the risk of complications, and/or the likelihood of implant failure. Such factors include:

- poor oral hygiene
- history of smoking/vaping/tobacco use
- history of periodontal disease
- history of orofacial radiotherapy**
- bruxism and unfavourable jaw relations

- use of chronic medications that may delay healing or increase the risk of complications including, but not limited to, chronic steroid therapy, anti-coagulant therapy, TNF- α blockers, bisphosphonate and cyclosporin

*** The potential for implant failure and other complications increases when implants are placed in irradiated bone, as radiotherapy can lead to progressive fibrosis of blood vessels and soft tissue (i.e., osteoradionecrosis), resulting in diminished healing capacity. Contributing factors to this increased risk include the timing of implant placement in relation to radiation therapy, the proximity of radiation exposure to the implant site, and the radiation dosage at that site.*

Should the device not operate as intended, it must be reported to the manufacturer of the device. The contact information or the manufacturer of this device to report a change in performance is: sicomplaints@southernimplants.com.

Side effects

The clinical outcome of treatment is influenced by multiple factors. The following side effects and residual risks are associated with the device group or the surgical procedures involved in its use. They may necessitate further treatment, revision surgery or additional visits to the relevant medical professional. Furthermore, these side effects and residual risks can vary in both severity and frequency.

Therefore, the full list of known potential undesirable side-effects and residual risks that have been identified in relation to the dental abutment device group include:

- Abutment fracture
- Abutment screw fracture or distortion
- Allergic reaction(s) to the abutment material
- Anaesthesia, paraesthesia, hyperesthesia, and hypoesthesia (transient or permanent)
- Bleeding
- Bleeding on probing
- Bruising
- Complications requiring revision surgery
- Difficulty in retrieving the abutment
- Gingival inflammation
- Gingival recession
- Hyperplastic tissue response
- Implant failure due to insufficient levels of osseointegration
- Improper implant positioning resulting in prosthetic compromise
- Infection (acute and/or chronic)
- Localized inflammation
- Loosening of the abutment screw and/or retaining screw
- Loss or damage to adjacent teeth
- Marginal bone loss
- Micromovements and implant instability
- Misfit or improper connection at the implant-abutment interface
- Overloading of the abutment/implant
- Pain, tenderness or discomfort
- Peri-implantitis, peri-mucositis or otherwise poor peri-implant soft tissue health
- Periodontal inflammation
- Phonetic difficulties
- Prosthetic failure
- Soft tissue irritation
- Sub-optimal aesthetic result
- Wound dehiscence or poor healing

Precaution: maintaining sterility protocol

The Overdenture Abutments are packed in a peel pouch or blister base with a “peel-back” lid. Labelling information is located on the bottom half of the pouch, inside the packet or on the surface of the peel-back lid. Sterility is assured unless the pouch or blister is damaged or opened.

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established.

The contact information for the manufacturer of this device to report a serious incident is as follows:

sicomplaints@southernimplants.com.

Materials

Material Type Commercially Pure Titanium (Grade 4) and Titanium Alloy (Ti-6Al-4V), Cobalt Chrome

Commercially pure titanium (grade 4) according to ASTM F67 and ISO 5832-2

Chemical Components	Residuals* (Fe, O, C, N, H) in total	Titanium (Ti)
Composition, % (mass/mass)	<1.05	Balance

*Where, Fe = Iron, O = Oxygen, C= Carbon, N = Nitrogen, H = Hydrogen

Titanium alloy (Ti-6Al-4V) according to ASTM F136 and ISO 5832-3

Chemical Components	Aluminium (Al)	Vanadium (V)	Residuals* (Fe, O, C, N, H) in total	Titanium (Ti)
Composition, % (mass/mass)	5.50 – 6.75	3.50 – 4.50	<0.55	Balance

*Where, Fe = Iron, O = Oxygen, C= Carbon, N = Nitrogen, H = Hydrogen

Wrought Cobalt Chromium Molybdenum Alloy (Co-28Cr-6Mo) according to ASTM F1537

Chemical Components	Cobalt (Co)	Chromium (Cr)	Molybdenum (Mo)	Residuals* (Ni, Fe, Mn, Si, C, N, P, S) in total
Composition, % (mass/mass)	Balance	27.0 – 30.0	5.00 – 7.00	<4.20

*Where, Fe = Iron, Ni = Nickel, C= Carbon, N = Nitrogen, Mn = Manganese, Si = Silicon, P = Phosphorous, S = Sulphur

Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

Magnetic Resonance (MR) Safety

Nonclinical testing has demonstrated that the Southern Implants® Dental/Craniomaxillofacial implants, metallic abutments and prosthetic screws are MR conditional.

MRI Safety Information

MR Conditional

A patient with the Dental/Craniomaxillofacial implants, metallic abutments and prosthetic screws from Southern Implants® may be safely scanned under the following conditions. Failure to follow these conditions may result in injury to the patient.

Device Name	Dental/Craniomaxillofacial implants, metallic abutments and prosthetic screws
Static Magnetic Field Strength (B₀) [T]	1.5T or 3T
Maximum Spatial Field Gradient [T/m] [gauss/cm]	30 T/m (3000 gauss/cm)
Type of RF Excitation	Circularly polarized (CP) (i.e., quadrature driven)
RF Transmit Coil Type	Integrated Whole Body Transmit Coil
RF Receive Coil	Any receive-only RF coil may be used
Operating Mode	Normal Operating Mode
RF Conditions	Normal Operating mode: (Head SAR of 3.2 W/kg, 2 W/kg whole body)
Scan Duration	If the closest part of implant <25cm away from the isocenter in feet-head direction: a continuous scan time of 15 minutes will require a cooling delay of at least 5 minutes. If the closest part of implant ≥ 25cm away from the isocenter in feet-head direction: This is limited to 1 hour of scan duration.
MR Image Artifact	The presence of this implant may produce an image artifact.

Additional information:

1. Local transmit coil: Local transmit coil safety was not evaluated. Local transmit coil is only allowed if the implant is ≥ 25cm away from isocenter and for a wbSAR of up to 2W/kg.
2. Artifacts: MR Image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the implant. Tests following guidelines in ASTM F2119-24 showed artifact widths of ≤20 mm from the surface of the implant for the gradient echo sequence and ≤10 mm from the surface of the implant for the spin echo sequence.

Removable restorations should be taken out prior to scanning, as is done for watches, jewellery etc.

Should there be no MR symbol on the product label, please note that this device has not been evaluated for safety and compatibility in the MR environment.

Summary of Safety and Clinical Performance (SSCP)

As required by the European Medical Device Regulation (MDR; EU2017/745), a Summary of Safety and Clinical Performance (SSCP) is available for perusal with regard to Southern Implants® product ranges.

The relevant SSCP can be accessed at <https://ec.europa.eu/tools/eudamed>.

NOTE: the above website will be available upon the launch of the European Database on Medical Devices (EUDAMED). In the interim, email clinical@southernimplants.com to request a Summary of Safety and Clinical Performance (SSCP) for these devices.

Disclaimer of liability

This product is part of the Southern Implants® product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants® product range and take full responsibility for the correct indications and use of this product. Southern Implants® does not assume liability for damage due to incorrect use. Please note that some Southern Implants® products may not be cleared or released for sale in all markets.

Basic UDI

Product	Basic-UDI Number
Basic-UDI for Metal Abutments	60095440387296
Basic-UDI for MA Overdenture Abut Titanium	60095440501682
Basic-UDI for MA Overdenture Abut Chrome C	6009544050157Y

Related literature and catalogues

CAT-2004 - Tri-Nex® Implants Product Catalogue

CAT-2005 - IT Implants Product Catalogue

CAT-2010 - Osseointegrated Fixtures Catalogue

CAT-2020 - External Hex Implants Product Catalogue

CAT-2042 - Deep Conical Implants Product Catalogue

CAT-2043 - Internal Hex Implants Product Catalogue

CAT-2060 - PROVATA® Implants Product Catalogue

CAT-2069 - INVERTA® Implants Product Catalogue

CAT-2070 - Zygomatic Implants Product Catalogue

CAT-2093 - Single Platform Implants Product Catalogue

Symbols and warnings

Manufacturer

CE mark

Prescription device*

Sterilised using irradiation

Non-sterile

Use by date (mm-yy)

Do not reuse

Do not resterilize

Catalogue number

Batch code

Medical device

Authorised representative in the European Community

Authorised representative for Switzerland

Date of manufacture

Magnetic Resonance conditional

Magnetic Resonance safe

Single sterile barrier system with protective packaging inside

Single sterile barrier system

Consult instruction for use

Caution

Keep away from sunlight

Do not use if package is damaged

* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

All rights reserved. Southern Implants®, the Southern Implants® logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants®. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale. It is the responsibility of the clinician to inspect the symbols that appear on the packaging of the product in use.

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