
English	INSTRUCTIONS FOR USE: Southern Implants® Screw Remover Drills and Guides
Español	INSTRUCCIONES DE USO: Fresas y guías de extracción de tornillos de Southern Implants®
Italiano	ISTRUZIONI PER L'USO: Southern Implants® Punte e Guide Rimuovi-Vite
Français	MODE D'EMPLOI : Forets extracteurs de vis et guides Southern Implants®
Deutsch	GEBRAUCHSANWEISUNGEN: Southern Implants® Bohrer und Führungen zur Schraubenentfernung
Português	INSTRUÇÕES DE UTILIZAÇÃO: Southern Implants® Brocas e guias para remoção de parafusos
Ελληνικά	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ: Southern Implants® Τρύπανα και οδηγοί αφαίρεσης βιδών
Svenska	ANVISNINGAR FÖR ANVÄNDNING: Southern Implants® Skrub borttagare borrar och guider

ARCHIVED

South Africa - Headquarters: 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

Subsidiaries

Australia

Southern Implants Australia
T: +61-(0)-8-9466-2627
E: info@southernimplants.com.au

Spain and Portugal

Southern Implants Iberica
T: +34 935 053 507
E: info@southernimplants.es

United Kingdom and Ireland

Southern Implants UK
T: +44-20-8899-6845 / 6 / 7
E: info@southernimplants.co.uk

USA and Canada

Southern Implants North America Inc.
T: +1-561-472-0990
E: customercare@southernimplants.com

Description

Screw remover drills are intended to be used to remove broken screws from an implanted dental implant. The screw removers are intended for single use on a single patient. The screw remover has cutting edges or cutting threads that cut with a counter clockwise rotation when rotated by a dental handpiece. These cutting edges or threads transfer the counter clockwise torque to the broken screw and unscrew the thread of the broken screw. The screw removers are available in lengths between 7.0 mm and 8.97 mm.

Screw remover guides are a hand-held manual surgical tool designed to remove implanted screws that have a damaged screw recess by self-tapping into the defective recess/head and providing sufficient grip to rotate and remove the screw in the normal way. Devices typically have a conical tip that is placed into the implant connection to provide a channel through which the screw remover is guided during use. The device is also intended to protect the implant internal connection during the screw removal procedure.

Intended use

Drills and handpieces are intended to be used to prepare the osteotomy for implant placement.

Indications for use

Dental implant therapy is indicated for patients with one or more missing teeth that wish to have the tooth/teeth restored with a dental implant and prosthetic tooth/teeth.

Intended user

Maxillo-facial Surgeons, General Dentists, Orthodontists, Periodontists, Prosthodontists and other appropriately trained and experienced implant users.

Intended environment

The devices are intended to be used in a clinical environment such as an operating theatre or a dentist consultation room.

Intended patient population

The intended patient population for the device is not dissimilar to that of dental implant therapy. The intended patient population for implant therapy is partially or fully edentulous patients requiring prosthetic dental restoration in the upper or lower jaw. Restorations may comprise single teeth, partial or full bridges and are fixed restorations.

Compatibility information

Southern Implants screw remover drills and guides are described in Table A and Table B. The screw remover drills attach to a handpiece of an implant motor unit and have an ISO 1797 compatible latch. This is to connect the drill to the handpiece of an implant motor unit. The screw remover has a laser marking on the shaft of the drill, to indicate the product code. The I-SR-6 has a 1.15 mm cutting tip, the I-SR-3 screw remover has cutting threads that cut with a counter clockwise rotation when rotated by a dental handpiece. Drill guides are anodized to the prosthetic diameter of the implant as described in Table B.

Table A – Screw Remover Drills

Drill Code	Material	Coating (if any)	Number of Uses
I-SR-6	Stainless Steel (DIN 1.4197)	-	1
I-SR-3	Stainless Steel (DIN 1.4197)	-	1

Table B – Screw Remover Guides

Item Code	Material	Anodizing (if any)	Implant Interface	Number of Uses
EXTERNAL HEX				
I-SRG-MS-C-IP	Titanium Grade 4 (ASTM F67)	Pink	Ø3 Interface Implants	Up to 10
I-SRG-EXT-IBN	Titanium Grade 4 (ASTM F67)	Blue	Ø3.25 Interface Implants	Up to 10
I-SRG-EXT-HEX	Titanium Grade 4 (ASTM F67)	Silver	Other Ex Hex Interface Implants	Up to 10
TRI-NEX®				
I-SRG-L-35	Titanium Grade 4 (ASTM F67)	Pink/ Purple	Ø3.5 Interface Implants	Up to 10

I-SRG-L-43	Titanium Grade 4 (ASTM F67)	Yellow	Ø4.3 Interface Implants	Up to10
I-SRG-L-50	Titanium Grade 4 (ASTM F67)	Blue	Ø5.0 Interface Implants	Up to10
I-SRG-L-60	Titanium Grade 4 (ASTM F67)	Green	Ø6.0 Interface Implants	Up to10
DEEP CONICAL (DC)				
I-SRG-DC3	Titanium Grade 4 (ASTM F67)	Yellow	DC3 Interface Implants	Up to10
I-SRG-DC4	Titanium Grade 4 (ASTM F67)	Blue	DC4 Interface Implants	Up to10
I-SRG-DC5	Titanium Grade 4 (ASTM F67)	Pink/ Purple	DC5 Interface Implants	Up to10
INTERNAL HEX (M-SERIES & PROVATA®)				
I-SRG-M	Titanium Grade 4 (ASTM F67)	Yellow	Standard Internal Hex Interface Implants	Up to10
I-SRG-Z	Titanium Grade 4 (ASTM F67)	Yellow	Wide Internal Hex Interface Implants	Up to10
IT (INTERNAL OCTAGON)				
I-SRG-IT	Titanium Grade 4 (ASTM F67)	Silver	All IT Implants	Up to10

Surgical procedures for screw removal

Step 1

The use of a fine probe to gently rotate the broken screw in an anti-clockwise direction, may often turn the screw out, to the point that it may be gripped with a haemostat and rotated out fully. The screw is usually quite loose once it has fractured, and should be freely able to rotate.

Step 2

Should step 1 not be successful and the screw is not loose, the use of an ultrasonic scaler in the same manner as the probe, should be applied. Often, the ultrasound assists in the rotation/freeing up of the screw. Once the screw has been exposed to ultrasound, Step 1 should be applied a second time.

Step 3

If the ultrasound is not successful on its own. The correct drill guide must be selected for the corresponding implant. Place the guide on top of the implant interface, it must be ensured that the drill guide is correctly seated on the implant interface. Each drill guide has a square connection on the opposing end, allowing them to be stabilized using a ratchet or torque wrench. The screw removal drills rotate in a counter-clockwise direction. The first drill, (I-SR-6), cuts a pilot hole in the centre of the screw to be removed. It is possible that the drilling process of the I-SR-6 may thread the screw out of the implant, if the drill grips the screw while cutting the pilot hole.

Step 4

If step 3 is not successful, the I-SR-3 will be used. Once the pilot hole has been cut by the I-SR-6, the I-SR-3 is designed to wedge into and engage the internal walls of the pilot hole, and rotate the fractured screw out in a counter-clockwise direction. Should this not be successful, a small rose head bur may also be tried to rotate the screw out by hand. After successful screw removal, it is important to ensure that the internal thread of the implant has not been damaged. This can be done using a healing abutment or impression coping pin. Should there be significant resistance when checking the internal thread for damage, it may be necessary to clean the internal thread of the implant by carefully using a thread tap.

NOTE:

- The screw removal instrumentation may only be used with a high torque – low speed motor unit, in reverse set at 600 – 800 rpm for the I-SR-6 with copious irrigation, and 10 - 20 rpm for the I-SR-3. Do not overheat the implant via friction this can lead to failure of the implant.
- Gently insert the drill into the latch-grip handpiece, ensuring that it is fully engaged. If the latch is not properly engaged in the handpiece, the torque applied when activating the motor unit may result in the possible distortion of the latch or damage to the handpiece, due to the incorrect distribution of forces within the latch-grip mechanism.
- Consult the instructions for use of the relevant handpiece manufacturer, to ensure proper engagement of the latch.
- Do not apply more than 40-45Ncm to any latch type instrument. This could cause damage to the handpiece gearing mechanism, and/or latch of the instrument.

Clinical benefits

The clinical benefits of the device are not dissimilar to that of dental implant therapy. These include improved chewing function, speech, aesthetics and patient psychological wellbeing. Through this procedure patients can expect to have their missing teeth replaced and/or crowns restored.

The clinical benefits specific to drills include:

- Effective cutting of bone within acceptable time limits.
- Precise preparation of the implant site to the dimensions of the intended implant.
- Can remain effective after multiple uses to the extent stipulated by the manufacturer.

The clinical benefits of handpiece devices include:

- Successful achievement of respective intended function (e.g., implant insertion, drill extension, screw removal)

Storage, cleaning and sterilisation

For screw remover drills:

The component is supplied sterile (sterilised by gamma irradiation) and intended for single-use prior to the expiration date (see packaging label). Sterility is assured unless the container or seal is damaged or opened. If packaging is damaged do not use the product and contact your Southern representative or return to Southern Implants®. The devices must be stored in a dry place at room temperature and not exposed to direct sunlight. Incorrect storage may influence device characteristics. Do not reuse components indicated for single-use only. Reusing these components may result in:

- damage to the surface or critical dimensions, which may result in performance and compatibility degradation.
- adds the risk of cross-patient infection and contamination if single-use items are reused.

Southern Implants® does not accept any responsibility for complications associated with reused single-use components.

For screw remover guides:

Limitations on reusable items

A direct value for reusable instruments cannot be given. Frequent processing may have minor effects on the instruments. The product life is normally determined by wear and damage during use, thus instruments if properly cared for and inspected after each use, can be reused many times. Maintain a checklist for these instruments recording the number of uses.

Prior to reprocessing the device, it should be thoroughly inspected and tested to determine its suitability for reuse.

NOTE: during use, handle drills and instruments with sterile tweezers to minimise contamination of the instrument tray and risk of damage to sterile surgical gloves.

Containment

As soon as practically possible, remove all visible residue after use (bone, blood or tissue), by immersing the instrument in cold water (dried soil may be difficult to remove).

Precleaning

Disassemble instruments from hand pieces and all connecting parts from instruments in order to clean soil from obstructed areas. Remove the PEEK bits from placement tools. Rinse with lukewarm water for 3 minutes and remove hardened debris with a soft nylon brush. Avoid mechanical damage during cleaning.

Manual cleaning or automated cleaning

Prepare an ultrasonic bath with suitable detergent (i.e., Steritech instrument cleaner - 5% dilution), sonicate for 20 minutes (alternative methods can be used if proven by the end user). Rinse with purified/sterile water.

NOTE: always follow the instructions for use of the manufacturers for cleaning agents and disinfectants.

Load devices into a thermidisinfector. Run the cleaning and disinfection cycle, followed by the drying cycle.

Drying

Dry the instruments both inside and outside with filtered compressed air or single use, lint free wipes. Pack the instruments as quickly as possible after removal into the storage container. If additional drying is necessary, dry in a clean location.

Inspection

Do a visual inspection of the items to check for damages on instruments.

Packaging

Use the correct packaging material as indicated for steam sterilisation to ensure sterility is maintained. Double packaging is recommended. Where appropriate, the cleaned, disinfected and checked devices can be assembled and placed in instrument trays as applicable. Instrument trays can be double wrapped or placed in sterilisation bags.

Sterilisation

Southern Implants® recommends one of the following procedures to sterilise the restoration prior to use:

1. prevacuum sterilisation method: steam sterilise the abutments at 132°C (270°F) at 180 - 220 kPa for 4 minutes. Dry for at least 20 minutes in the chamber. Only an approved wrap or pouch for steam sterilisation must be used.
2. for users in the USA: prevacuum sterilisation method: wrapped, steam sterilise at 135°C (275°F) at 180 - 220 kPa for 3 minutes. Dry for 20 minutes in the chamber. Use a wrap or pouch that is cleared for the indicated steam sterilisation cycle.

NOTE: users in the USA must ensure that the steriliser, wrap or pouch, and all steriliser accessories are cleared by the FDA, for the intended sterilisation cycle.

Storage

Maintain packaging integrity to ensure sterility in storage. Packaging should be completely dry before storage to avoid corrosion and degradation of cutting edges.

Contraindications

Since the device is required or utilized as part of dental implant surgery, the contraindications for the device are not dissimilar to dental implant therapy as a whole.

The contraindications for dental implant therapy include:

- who are medically unfit for dental implant procedures
- where adequate numbers of implants could not be placed to achieve full functional support of the prosthesis,
- who are allergic or have hypersensitivity to pure titanium or titanium alloy (Ti-6Al-4V), gold, palladium, platinum, iridium or stainless steel.
- who are under the age of 18, have poor bone quality, blood disorders, infected implant site, vascular impairment, uncontrolled diabetes, drug or alcohol abuse, chronic high dose steroid therapy, anti-coagulant therapy, metabolic bone disease, radiotherapy treatment.

Warnings and precautions

THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING

- For the safe and effective use of dental implants, new technologies/systems and the device, it is strongly suggested that specialised training be undertaken. This training should include hands-on methods to gain competency on the proper technique, biomechanical requirements and radiographic evaluations required for the specific system.
- Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone.
- Use of the device with incompatible or non-corresponding devices can result in device failure or poor performance.
- When handled intraorally, devices must be adequately secured against aspiration. Aspiration of products may lead to infection or physical injury.
- Electrosurgery should not be attempted around metal implants as they are conductive.
- The use of non-sterile items can lead to secondary infections of the tissue or transfer infectious diseases.
- Failure to perform appropriate cleaning, resterilisation and storage procedures as per the Instructions for Use document can result in device damage or secondary infection.
- Blunt drills may cause damage to the bone which could compromise osseointegration.
- For short implants, clinicians should closely monitor patients for any of the following conditions: peri implant bone loss, changes to implant's response to percussion, or radiographic changes in bone to implant contact along the implant's length. If the implant shows mobility or greater than 50% bone loss, the implant should be evaluated for possible removal. If the clinicians choose a short implant, then clinicians should consider a two-stage surgical approach, splinting a short implant to an additional implant, and placement of the widest possible fixture. Allow longer periods for osseointegration and avoid immediate loading.

NOTE: Training should be undertaken by both new and experienced implant users before using a new system or attempting to do a new treatment method.

It is important that thorough screening of prospective implant candidates is performed. This screening should include:

- a comprehensive medical and dental history.
- visual and radiological inspection to determine adequate bone dimensions, anatomical landmarks, occlusal conditions and periodontal health. Failure to accurately estimate lengths of drills relative to radiographic measurements can result in permanent injury to nerves or other vital structures.
- the evaluation and consideration of bruxism and unfavourable jaw relations.

Proper preoperative planning with a good team approach between well trained surgeons, restorative dentists and lab technicians is essential for successful implant treatment.

For successful implant treatment, it is important to:

- minimise the trauma to the host tissue as this increases the potential for successful osseointegration.
- be aware and avoid damage to vital anatomical structures such as nerves, veins and arteries. Injury to these structures may cause serious complications, including, injury to the eye, nerve damage and excessive bleeding. It is essential to protect the infraorbital nerve.
- identify actual measurements relative to the radiographic data as failing to do so could lead to complications.

Special care should be taken when treating patients with local or systemic risk factors that could affect the healing of the bone and soft tissue or otherwise increase the severity of the side effects, risk of complications and/or implant failure. These factors include:

- smoking.
- history of periodontal disease.
- poor oral hygiene.
- history of orofacial radiotherapy.**
- bruxism and unfavourable jaw relations.

** The potential for implant failure and other complications increases when implants are placed in irradiated bone as radiotherapy can result in progressive fibrosis of vessels and soft tissue, leading to diminished healing capacity.

The responsibility for proper patient selection, adequate training and experience in implant placement and the provision of the appropriate information required for informed consent rests with the practitioner.

Should the device not operate as intended, it must be reported to the manufacturer of the device. The contact information or the manufacturer of this device to report a change in performance is: sicomplaints@southernimplants.com

Side effects

Since the device is required or utilized as part of dental implant surgery, the side effects of the device are not dissimilar to those of dental implant therapy as a whole. Common side effects of dental implant therapy include pain, inflammation, phonetic difficulties and gingival inflammation. Other, less common side effects or complications to dental implant therapy include, but are not limited to: (1) wound dehiscence; (2) peri-implantitis; (3) transient weakness, numbness, and/or pain associated with mild nerve injury; (4) fat emboli formation; (5) marginal bone loss within acceptable limits; (6) allergic reaction(s) to implant and/or abutment material; (7) unspecified infection; (8) implant failure due to insufficient levels of osseointegration; (9) breakage of the implant and/or abutment; (10) loosening of the abutment screw and/or retaining screw; (11) complications requiring revision of the dental implant; (12) nerve damage resulting in permanent weakness, numbness, or pain; (13) loosening of the implant requiring revision surgery; (14) perforation of the maxillary sinus; (15) perforation of the labial and/or lingual plates and (16) bone loss possibly resulting in revision or removal of the implant.

Precaution: maintaining sterility protocol

Implants are packaged as follows:

1. An outer package consisting of a rigid, clear box which acts as protection for the inner package.
2. The inner package consisting of a blister pack (clear plastic-formed blister base with a TYVEK "peel-back" lid).
3. Within the inner package, there is a hollow tube which contains one implant suspended from a titanium ring, this ensures the implant never touches the inside of the plastic tube.
4. Labelling information is located on the surface of the peel-back lid and on the outside of the rigid box.

Care must be taken to maintain the sterility of the implant by proper opening of the packaging and handling of the implant.

1. Open the implant package in the non-sterile field, with non-sterile gloves, tear the address label to open the box.
2. With non-sterile gloves, remove the inner blister pack. Do not place the plastic box or blister pack-lid onto the sterile field. The contents of this inner package are sterile.
3. The sealed blister is to be opened by an assistant (with nonsterile gloves): remove the TYVEK lid and drop or place the sterile tube onto the sterile field, open the tube cap and attach the implant placement tool onto the implant and carefully remove from the sterile tube. Do not touch the sterile implant.

Other sterile components are packed in a peel pouch or blister base with a "peel-back" lid. Labelling information is located on the bottom half of the pouch, inside the packet or on the surface of the peel-back lid. Sterility is assured unless the pouch is damaged or opened. Non-sterile components are supplied clean but not sterile in a peel pouch or blister base with peel-back lid. Labelling information is located on the bottom half of the pouch or on the surface of the peel-back lid.

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established.

The contact information for the manufacturer of this device to report a serious incident is as follows: sicomplaints@southernimplants.com.

Materials

Screw remover drills	Stainless Steel (DIN 1.4197)
Screw remover guides	Titanium Grade 4 (ASTM F67)

Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

Summary of Safety and Clinical Performance (SSCP)

As required by the European Medical Device Regulation (MDR; EU2017/745), a Summary of Safety and Clinical Performance (SSCP) is available for perusal with regard to Southern Implants® product ranges.

The relevant SSCP can be accessed at <https://ec.europa.eu/tools/eudamed>.

NOTE: the above website will be available upon the launch of the European Database on Medical Devices (EUDAMED).

Disclaimer of liability

This product is part of the Southern Implants® product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants® product range and take full responsibility for the correct indications and use of this product. Southern Implants® does not assume liability for damage due to incorrect use. Please note that some Southern Implants® products may not be cleared or released for sale in all markets.

Basic UDI

Product	Basic-UDI Number
Basic-UDI for Drills and Handpiece Devices	6009544038759C
Basic-UDI for Reusable Instruments	6009544038769E

Related literature and catalogues

CAT-2004 - Tri-Nex® Implants Product Catalogue
CAT-2020 - External Hex Implants Product Catalogue
CAT-2042 - Deep Conical (DC) Implants Product Catalogue
CAT-2043 - Internal Hex (M-Series) Implants Product Catalogue
CAT-2060 - Internal Hex (PROVATA®) Implants Product Catalogue
CAT-2005 - IT (Internal Octagon) Implants Product Catalogue
CAT-2070 - Zygomatic Implants Product Catalogue
CAT-2069 - INVERTA® Implants Product Catalogue
CAT-2092 - Soft Bone Implants Product Catalogue
CAT-2010 - Osseointegrated Fixtures Product Catalogue

Symbols and warnings

											
Manufacturer: Southern Implants® 1 Albert Rd, P.O Box 605 IRENE, 0062, South Africa. Tel: +27 12 667 1046	CE mark 2797	Prescription device*	Sterilised using irradiation	Non-sterile	Use by date (mm-yy)	Do not reuse	Do not resterilise	Catalogue number	Batch code	Medical device	Authorised representative in the European Community
											
Authorised representative for Switzerland	Date of manufacture	Magnetic Resonance conditional	Magnetic Resonance safe	Single sterile barrier system with protective packaging inside	Single sterile barrier system	Consult instruction for use	Caution	Keep away from sunlight	Do not use if package is damaged		

* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

All rights reserved. Southern Implants®, the Southern Implants® logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants®. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale. It is the responsibility of the clinician to inspect the symbols that appear on the packaging of the product in use.

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