
English	INSTRUCTIONS FOR USE: Southern Implants® PEEK Temporary Abutments
Español	INSTRUCCIONES DE USO: Pilares temporales PEEK Southern Implants®
Italiano	ISTRUZIONI PER L'USO: Southern Implants® Monconi Provvisori in PEEK
Français	MODE D'EMPLOI : Piliers provisoires PEEK Southern Implants®
Deutsch	GEBRAUCHSANWEISUNGEN: Southern Implants® Provisorische PEEK-Abutments
Português	INSTRUÇÕES DE UTILIZAÇÃO: Pilares temporários PEEK da Southern Implants®
Ελληνικά	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ: Southern Implants® PEEK Προσωρινά κολοβώματα
Svenska	ANVISNINGAR FÖR ANVÄNDNING: Southern Implants® PEEK Tillfälliga distanser

ARCHIVED

South Africa - Headquarters: 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

Subsidiaries

Australia

Southern Implants Australia
T: +61-(0)-8-9466-2627
E: info@southernimplants.com.au

Spain and Portugal

Southern Implants Iberica
T: +34 935 053 507
E: info@southernimplants.es

United Kingdom and Ireland

Southern Implants UK
T: +44-20-8899-6845 / 6 / 7
E: info@southernimplants.co.uk

USA and Canada

Southern Implants North America Inc.
T: +1-561-472-0990
E: customercare@southernimplants.com

Description

The Southern Implants® PEEK Temporary Abutments are premanufactured and are available in a variety of connections, engaging and non-engaging, to fit the implant systems manufactured by Southern Implants. They are used as an aid in manufacturing a prosthesis for prosthetic rehabilitation. They can either be used for direct connection to an endosseous implant or they can be used for connecting the prosthesis to a compact conical abutment. The PEEK abutments are provided sterile; however, it will no longer be sterile after modification.

Intended use

The Southern Implants® PEEK Temporary Abutment family is intended to aid in the treatment of partially or fully edentulous patients eligible for placement of one or more dental implants as a means of fixing a permanent or removable single crown, partial or full-arch dental prosthesis in the upper or lower jaw. These devices are premanufactured prosthetic components directly connected to endosseous dental implants of the mandible and/or maxilla and are intended to be used to protect the implant's internal threads and implant during the healing phase and to shape the soft tissue in preparation for the final prosthesis. The PEEK temporary abutments are attached to the implants by a retention screw made from titanium alloy. The system constituents are classified as medical devices and are intended for single use on a single patient.

Indications for use

The Southern Implants® PEEK Abutments are premanufactured prosthetic components directly connected to endosseous dental implants and are intended for provisional use up to 180 days as an aid in prosthetic rehabilitation.

Intended user

The intended user for this system includes Dental Technicians, Maxillofacial Surgeons, General Dentists, Orthodontists, Periodontists, Prosthodontists and other appropriately trained and experienced implant users.

Intended environment

This system is intended to be used in a dental laboratory for the manufacture of the restoration and in a clinical environment such as an operating theatre or a dentist consultation room.

Intended patient population

The intended patient population for the device is not dissimilar to that of dental implant therapy. The intended patient population for implant therapy is partially or fully edentulous patients requiring prosthetic dental restoration in the upper or lower jaw. Restorations may comprise single teeth, partial or full bridges and are fixed restorations.

Compatibility information

Southern Implants' implants should be restored with Southern Implants' components. In the Southern Implants' range there are 5 implant connections. The implant code and connection type can be identified by specific abbreviations in the product codes. Range identifiers are summarised in Table A.

Table A

Implant connection type	Compatible device
External Hex (EX)	Parts labelled PKIP2H, PKBN2H, PKB2H, PKBA2H, PKBBB2H, PKMAX9-2H, PKR-EX-30, PKR-TOP1/2-EX-30, PKR-EX-34, PKR-TOP1/2/3/4-EX-34, PKR-EX-40, PKR-TOP1/2/3/4-EX-40, PKR-EX-50, PKR-EX-60 and PKR-EX-70 for engaging items.
	Parts labelled PKIP2NH, PKBN2NH, PKB2NH, PKBA2NH, PKBBB2NH, PKMAX9-2NH, PKR-NX-30, PKR-NX-34, PKR-NX-40, PKR-NX-50, PKR-NX-60 and PKR-NX-70 for non-engaging items.
Tri-Nex (EL) (Lobe)	Parts labelled PKC-EL-(Ø) and PKR-EL-(Ø) for engaging items.
	Parts labelled PKC-NL-(Ø) and PKR-NL-(Ø) for non-engaging items.
Deep Conical (DC)	Parts labelled PKC-DC(Ø), PKR-DC(Ø), PKR-TOP-1/2-DC3, PKR-TOP1/2/3/4-DC4 and PKR-TOP1/2/3/4-DC5 for engaging items.
	Parts labelled PKC-NDC(Ø) and PKR-NDC(Ø) for non-engaging items.
Internal Hex (M)	Parts labelled PKC-M-2 and PKR-M (used with Ø3.75, Ø4.2, Ø5.0 mm platforms) for engaging items.
	Parts labelled PKC-NM-2 and PKR-NM (used with Ø3.75, Ø4.2 and Ø5.0 mm platforms) for non-

	engaging items.
Internal Hex PROVATA® (3M) (M) (Z)	Parts labelled PKC-3M-2 and PKR-3M (used with Ø3.3 mm platforms) for engaging items.
	Parts labelled PKC-3NM-2 and PKR-3NM (used with Ø3.3 mm platforms) for non-engaging items.
	Parts labelled PKC-M-2 and PKR-M (used with Ø4.0, Ø5.0 and Ø6.0 mm platforms) for engaging items.
	Parts labelled PKC-NM-2 and PKR-NM (used with Ø4.0, Ø5.0 and Ø6.0 mm platforms) for non-engaging items.
	Parts labelled PKC-Z-2 and PKR-Z (used with Ø7.0, Ø8.0 and Ø9.0 mm platforms) for engaging items.
	Parts labelled PKC-NZ-2 and PKR-NZ (used with Ø7.0, Ø8.0 and Ø9.0 mm platforms) for non-engaging items.
Internal Octagon (ITS) (IT6)	Parts labelled ITS-PKC1 and PKR-ITS (used with Ø4.8 mm platforms) for engaging items.
	Parts labelled ITS6-PKC1 and PKR-IT6 (used with Ø6.5 mm platforms) for engaging items.
Abutment Level	Parts labelled PKC-MC and PKR-MC-48 (used with Ø4.8 mm abutment platforms) for non-engaging items.
	Parts labelled PKC-MCW and PKR-MC-60 (used with Ø6.0 mm abutment platforms) for non-engaging items.

Table B summarises the recommended screw torque values of the Southern Screws when used with the PEEK abutments.

Table B

Range	Screw type	Torque
Abutment Level	1 Series (TSS1, TSU1, TSH1)	15 Ncm
External Hex (EX)	3 Series (TSSZ3, TSUZ3, TSHZ3)	< Ø4.0 mm implant interfaces: 15 Ncm ≥ Ø4.0 mm implant interfaces: 20 Ncm
	TS-P-16	15 Ncm
Tri-Nex (EL)	TS-L-18	15 Ncm
	TS-L-20	20 Ncm
Deep Conical (DC)	TS-DC3-14	15 Ncm
	TS-DC4-16	20 Ncm
	TS-DC5-20	20 Ncm
Internal Hex (M-Series & PROVATA®)	TS-Z-16	15 Ncm
	TS-Z-18	15 Ncm
Internal Octagon (IT)	TSIT2	20 Ncm

Clinical benefits

Through this procedure patients can expect to have their missing teeth replaced and/or crowns restored.

Storage, cleaning and sterilisation

The component is supplied sterile (sterilised by gamma irradiation) and intended for single-use prior to the expiration date (see packaging label). Sterility is assured unless the container or seal is damaged or opened. If packaging is damaged do not use the product and contact your Southern representative or return to Southern Implants®. The devices must be stored in a dry place at room temperature and not exposed to direct sunlight. Incorrect storage may influence device characteristics. Do not reuse components indicated for single-use only. Reusing these components may result in:

- damage to the surface or critical dimensions, which may result in performance and compatibility degradation.
- adds the risk of cross-patient infection and contamination if single-use items are reused.

Southern Implants® does not accept any responsibility for complications associated with reused single-use components.

Sterilisation

Southern Implants® recommends one of the following procedures to sterilise the restoration prior to use:

1. Pre-vacuum sterilisation method: steam sterilise the abutments at 132°C (270°F) at 180 - 220 kPa for 4 minutes. Dry for at least 20 minutes in the chamber. Only an approved wrap or pouch for steam sterilisation must be used.
2. Pre-vacuum sterilisation method: wrapped, steam sterilise at 135°C (275°F) at 180 - 220 kPa for 3 minutes. Dry for 20 minutes in the chamber. Use a wrap or pouch that is cleared for the indicated steam sterilisation cycle.

NOTE: users in the USA must ensure that the steriliser, wrap or pouch, and all steriliser accessories are cleared by the FDA, for the intended sterilisation cycle.

Contraindications

Do not use in:

- patients medically unfit for oral surgical procedures.
- cases where an inadequate number of implants can be placed to provide the required functional support for the prosthesis.
- minors (patients under the age of 18 years).
- patients with poor bone quality.
- patients with blood disorders.
- presence of infection at the implant site.
- patients with vascular impairment.
- patients with uncontrolled diabetes.
- patients with drug or alcohol abuse dependencies.
- patients undergoing chronic high dose steroid therapy.
- patients undergoing anti-coagulant therapy.
- patients with metabolic bone disease.
- patients undergoing radiotherapy treatment.
- patients with pure titanium, titanium alloy (Ti6Al4V), gold, palladium, iridium or polyetheretherketone (PEEK) allergies or hypersensitivities.

Warnings and precautions

THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING.

- For the safe and effective use of dental implants it is strongly suggested that specialised training be undertaken, including hands-on training to learn proper technique, biomechanical requirements and radiographic evaluations.
- Products must be secured against aspiration when handled intraorally. Aspiration of products may lead to infection or unplanned physical injury.

Responsibility for proper patient selection, adequate training, experience in the placement of implants and providing appropriate information for informed consent rests with the practitioner. Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone. Implant failure increases when implants are placed in irradiated bone as radiotherapy can result in progressive fibrosis of vessels and soft tissue, leading to diminished healing capacity.

It is important to be aware and avoid damage to vital structures such as nerves, veins and arteries. Injuries to vital anatomical structures may cause serious complications such as injury to the eye, nerve damage and excessive bleeding. It is essential to protect the infraorbital nerve. Failing to identify actual measurements relative to the radiographic data could lead to complications.

New and experienced implant users should do training before using a new system or attempting to do a new treatment method. Take special care when treating patients who have local or systemic factors that could affect the healing of the bone and soft tissue (i.e., poor oral hygiene, uncontrolled diabetes, are on steroid therapy, smokers, infection in the nearby bone and patients who had orofacial radiotherapy).

Thorough screening of prospective implant candidates must be performed including:

- a comprehensive medical and dental history.
- visual and radiological inspection to determine adequate bone dimensions, anatomical landmarks, occlusal conditions and periodontal health.
- bruxism and unfavourable jaw relations must be taken into account.
- proper preoperative planning with a good team approach between well trained surgeons, restorative dentists and lab technicians is essential for successful implant treatment.
- minimising the trauma to the host tissue increases the potential for successful osseointegration.
- electrosurgery should not be attempted around metal implants as they are conductive.

Should the device not operate as intended, it must be reported to the manufacturer of the device. The contact information for the manufacturer of this device to report a change in performance is: sicomplaints@southernimplants.com.

Side effects

Potential side effects and temporary symptoms: pain, swelling, phonetic difficulties and gingival inflammation. Less common but more persistent symptoms include, but are not limited to: (1) allergic reaction(s) to implant and/or abutment material; (2) breakage of the implant and/or abutment; (3) loosening of the abutment screw and/or retaining screw; (4) infection requiring revision of the dental implant; (5) nerve damage that could cause permanent weakness, numbness, or pain; (6) histologic responses possibly involving macrophages and/or fibroblasts; (7) formation of fat emboli; (8) loosening of the implant requiring revision surgery; and (9) bone loss possibly resulting in revision or removal.

Precaution: maintaining sterility protocol

Implants are packaged as follows:

1. An outer package consisting of a rigid, clear box which acts as protection for the inner package.
2. The inner package consisting of a blister pack (clear plastic-formed blister base with a TYVEK “peel-back” lid).
3. Within the inner package, there is a hollow tube which contains one implant suspended from a titanium ring, this ensures the implant never touches the inside of the plastic tube.
4. Labelling information is located on the surface of the peel-back lid and on the outside of the rigid box.

Care must be taken to maintain the sterility of the implant by proper opening of the packaging and handling of the implant.

1. Open the implant package in the non-sterile field, with non-sterile gloves, tear the address label to open the box.
2. With non-sterile gloves, remove the inner blister pack. Do not place the plastic box or blister pack-lid onto the sterile field. The contents of this inner package are sterile.
3. The sealed blister is to be opened by an assistant (with nonsterile gloves): remove the TYVEK lid and drop or place the sterile tube onto the sterile field, open the tube cap and attach the implant placement tool onto the implant and carefully remove from the sterile tube. Do not touch the sterile implant.

Other sterile components are packed in a peel pouch or blister base with a “peel-back” lid. Labelling information is located on the bottom half of the pouch, inside the packet or on the surface of the peel-back lid. Sterility is assured unless the pouch is damaged or opened. Non-sterile components are supplied clean but not sterile in a peel pouch or blister base with peel-back lid. Labelling information is located on the bottom half of the pouch or on the surface of the peel-back lid.

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established.

The contact information for the manufacturer of this device to report a serious incident is as follows:
sicomplaints@southernimplants.com.

Materials

Abutment Material Type	White radiopaque medical grade PEEK
Abutment Screws Material Type	Grade 5 Titanium Alloy (ASTM F136)

Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

MR safety

MR safe: The Southern Implants® PEEK abutments (excluding the prosthetic screws) are fabricated from a material that is not affected by exposure to MRI energy and are MR safe. However, restorations fabricated from this material may be assembled with metal implants, abutments and screws that can be affected by MRI energy.

Summary of Safety and Clinical Performance (SSCP)

As required by the European Medical Device Regulation (MDR; EU2017/745), a Summary of Safety and Clinical Performance (SSCP) is available for perusal with regard to Southern Implants® product ranges.

The relevant SSCP can be accessed at <https://ec.europa.eu/tools/eudamed>.

NOTE: the above website will be available upon the launch of the European Database on Medical Devices (EUDAMED).

Disclaimer of liability

This product is part of the Southern Implants® product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants® product range and take full responsibility for the correct indications and use of this product. Southern Implants® does not assume liability for damage due to incorrect use. Please note that some Southern Implants® products may not be cleared or released for sale in all markets.

Basic UDI

Product	Basic-UDI Number
Basic-UDI for PEEK Abutments	6009544038749A

Related literature and catalogues

CAT-2004 - Tri-Nex Implants Product Catalogue
 CAT-2005 - IT Implants Product Catalogue
 CAT-2020 - External Hex Implants Product Catalogue
 CAT-2042 - Deep Conical Implants Product Catalogue
 CAT-2043 - Internal Hex Implants Product Catalogue
 CAT-2060 - PROVATA® Implants Product Catalogue
 CAT-2070 - Zygomatic Implants Product Catalogue
 CAT-2069 - INVERTA® Implants Product Catalogue
 CAT-2092 - Soft Bone Implants Product Catalogue

Symbols and warnings

											
Manufacturer: Southern Implants® 1 Albert Rd, P.O Box 605 IRENE, 0062, South Africa. Tel: +27 12 667 1046	CE mark 2797	Prescription device*	Sterilised using irradiation	Non-sterile	Use by date (mm-yy)	Do not reuse	Do not resterilise	Catalogue number	Batch code	Medical device	Authorised representative in the European Community
											
Authorised representative for Switzerland	Date of manufacture	Magnetic Resonance conditional	Magnetic Resonance safe	Single sterile barrier system with protective packaging inside	Single sterile barrier system	Consult instruction for use	Caution	Keep away from sunlight	Do not use if package is damaged		

* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

All rights reserved. Southern Implants®, the Southern Implants® logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants®. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale. It is the responsibility of the clinician to inspect the symbols that appear on the packaging of the product in use.

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