
English	INSTRUCTIONS FOR USE: Southern Implants® PEEK Temporary Abutments
Español	INSTRUCCIONES DE USO: Southern Implants® PEEK Temporary Abutments
Italiano	ISTRUZIONI PER L'USO: Southern Implants® Monconi Provvisori in PEEK
Français	MODE D'EMPLOI : Piliers provisoires PEEK de Southern Implants®
Deutsch	GEBRAUCHSANWEISUNG: Southern Implants® PEEK-Temporärabutments
Português	INSTRUÇÕES DE UTILIZAÇÃO: Pilares temporários PEEK Southern Implants®
Ελληνικά	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ: Southern Implants® PEEK Προσωρινά κολοβώματα
Svenska	BRUKSANVISNING: Southern Implants® PEEK Tillfälliga Abutment
Română	INSTRUCȚIUNI DE UTILIZARE: Ebutmente temporare Southern Implants® PEEK
Türkçe	KULLANMA TALİMATI: Southern Implants® PEEK Geçici Abutmentleri

South Africa - Headquarters: 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

Subsidiaries

Australia

Southern Implants Australia
T: +61-(0)-8-9466-2627
E: info@southernimplants.com.au

Spain and Portugal

Southern Implants Iberica
T: +34 935 053 507
E: info@southernimplants.es

United Kingdom and Ireland

Southern Implants UK
T: +44-20-8899-6845 / 6 / 7
E: info@southernimplants.co.uk

USA and Canada

Southern Implants North America Inc.
T: +1-561-472-0990
E: customercare@southernimplants.com

Description

The Southern Implants® PEEK Temporary Abutments are premanufactured and are available in a variety of connections, engaging and non-engaging, to fit the implant systems manufactured by Southern Implants. They are used as a temporary aid to seal an implant system connection, allow for soft tissue healing/profiling and/or provide a temporary supportive base for provisional dental restorations. They can either be used for direct connection to an endosseous implant or they can be used for connecting the prosthesis to a compact conical abutment. The PEEK abutments are provided sterile; however, it will no longer be sterile after modification.

Intended use

The Southern Implants® PEEK Temporary Abutments are intended to be temporarily connected to an implant system in order to seal the connection, allow for soft tissue healing/profiling and/or provide a temporary supportive base for provisional dental restorations.

Indications for use

The Southern Implants® PEEK Temporary Abutments are indicated for provisional use (up to 180 days) in the functional and aesthetic restoration of patients eligible for implant-retained prostheses.

Intended user

Dental Technicians, Maxillofacial Surgeons, General Dentists, Orthodontists, Periodontists, Prosthodontists, and other appropriately trained and/or experienced medical professionals.

Intended environment

The devices are intended to be used in a dental laboratory as part of the restoration design and manufacture as well as in a clinical environment such as an operating theatre or a dentist consultation room.

Intended patient population

The device is intended to be used in partially or fully edentulous patients eligible, or otherwise not contraindicated, for implant-retained prosthetic dental restorations.

Compatibility information

Southern Implants' implants should be restored with Southern Implants' components. In the Southern Implants' range there are 7 implant connections. The implant code and connection type can be identified by specific abbreviations in the product codes. Range identifiers are summarised in Table A.

Table A

Implant connection type	Compatible Platform	Engaging/ Non-Engaging	Compatible device
External Hex (EX)	Ø3.0 mm	Engaging	Parts labelled PKIP2H, PKC-TOP1/2-EX-30, and HA-PKC1/2-EX-30
		Non-Engaging	Parts labelled PKIP2NH
	Ø3.25 mm	Engaging	Parts labelled PKBN2H, PKC-TOP1/2/3/4-EX-34, and HA-PKC1/2/3/4-EX-34
		Non-Engaging	Parts labelled PKBN2NH
	Ø4 mm	Engaging	Parts labelled PKB2H, PKC-TOP1/2/3/4-EX-40, and HA-PKC1/2/3/4-EX-40
		Non-Engaging	Parts labelled PKB2NH
	Ø5 mm	Engaging	Parts labelled PKBA2H
		Non-Engaging	Parts labelled PKBA2NH
Tri-Nex (EL) (Lobe)	Ø3.5 mm	Engaging	Parts labelled PKC-EL-35
		Non-Engaging	Parts labelled PKC-NL-35
	Ø4.3 mm	Engaging	Parts labelled PKC-EL-43
		Non-Engaging	

	Ø5 mm	Non-Engaging	Parts labelled PKC-NL-43
		Engaging	Parts labelled PKC-EL-50
		Non-Engaging	Parts labelled PKC-NL-50
	Ø6 mm	Engaging	Parts labelled PKC-EL-60
		Non-Engaging	Parts labelled PKC-NL-60
Deep Conical (DC)	Ø3.0 mm	Engaging	Parts labelled PKC-DC3, PKC-TOP-1/2-DC3, and HA-PKC1/2-DC3
		Non-Engaging	Parts labelled PKC-NDC3
	Ø4.0 mm	Engaging	Parts labelled PKC-DC4, PKC-TOP1/2/3/4-DC4, and HA-PKC1/2/3/4-DC4
		Non-Engaging	Parts labelled PKC-NDC4
	Ø5.0 mm	Engaging	Parts labelled PKC-DC5, PKC-TOP1/2/3/4-DC5, and HA-PKC1/2/3/4-DC5
		Non-Engaging	Parts labelled PKC-NDC5
Internal Hex (M)	Ø3.75, Ø4.2, and Ø5.0 mm	Engaging	Parts labelled PKC-M-2
	Ø3.75, Ø4.2, and Ø5.0 mm	Non-Engaging	Parts labelled PKC-NM-2
Internal Hex PROVATA® (3M) (M) (Z)	Ø3.3 mm	Engaging	Parts labelled PKC-3M-2
		Non-Engaging	Parts labelled PKC-3NM-2
	Ø4.0, Ø5.0 and Ø6.0 mm	Engaging	Parts labelled PKC-M-2
		Non-Engaging	Parts labelled PKC-NM-2
	Ø7.0, Ø8.0 and Ø9.0 mm	Engaging	Parts labelled PKC-Z-2
		Non-Engaging	Parts labelled PKC-NZ-2
Internal Octagon (ITS) (IT6)	Ø4.8 mm	Engaging	Parts labelled ITS-PKC1
	Ø6.5 mm	Engaging	Parts labelled ITS6-PKC1
Abutment Level	Ø4.8 mm	Non-Engaging	Parts labelled PKC-MC
	Ø6.0 mm	Non-Engaging	Parts labelled PKC-MCW

Table B summarises the recommended screw torque values of the Southern Screws when used with the PEEK abutments.

Table B

Range	Screw type	Torque
Abutment Level	1 Series (TSS1, TSU1, TSH1)	15 Ncm
External Hex (EX)	3 Series (TSSZ3, TSUZ3, TSHZ3)	< Ø4.0 mm implant interfaces: 15 Ncm ≥ Ø4.0 mm implant interfaces: 20 Ncm
	TS-P-16	15 Ncm
Tri-Nex (EL)	TS-L-18	15 Ncm
	TS-L-20	20 Ncm
Deep Conical (DC)	TS-DC3-14	15 Ncm
	TS-DC4-16	20 Ncm
	TS-DC5-20	20 Ncm
Internal Hex (M-Series & PROVATA®)	TS-Z-16	15 Ncm
	TS-Z-18	15 Ncm
Internal Octagon (IT)	TSIT2	20 Ncm

Clinical performance

The PEEK Abutments are intended to seal the implant connection and create a stable environment that fosters the integration of the implant into the surrounding bone and the healing of the soft tissue. Consequently, the clinical performance of the PEEK Abutments is primarily defined by their effectiveness in facilitating osseointegration and promoting optimal soft tissue healing. This performance can be quantitatively assessed through rates of successful

osseointegration, the health of the peri-implant soft tissue, the incidence of soft tissue irritation, and the emergence profile established by the healing abutments.

Clinical benefits

Due to their role in the healing phase of dental implant therapy, the Southern Implants® PEEK Abutments provide clinical benefits associated with the enhancement of peri-implant soft tissue health and the improvement of aesthetic outcomes for the final restoration. Additionally, due to their intrinsic link with the overall success of the treatment, these abutments indirectly afford the benefits associated with the overall treatment, reflecting their contribution to the overall treatment system. These indirect benefits can be assessed through treatment success rates, patient satisfaction metrics and improvements in quality of life.

Storage, cleaning and sterilisation

The component is supplied sterile (sterilised by gamma irradiation) and intended for single-use prior to the expiration date (see packaging label). Sterility is assured unless the container or seal is damaged or opened. If packaging is damaged do not use the product and contact your Southern representative or return to Southern Implants®. The devices must be stored in a dry place at room temperature and not exposed to direct sunlight. Incorrect storage may influence device characteristics. Do not reuse components indicated for single-use only. Reusing these components may result in:

- damage to the surface or critical dimensions, which may result in performance and compatibility degradation.
- adds the risk of cross-patient infection and contamination if single-use items are reused.

Southern Implants® does not accept any responsibility for complications associated with reused single-use components.

Sterilisation

Southern Implants® recommends one of the following procedures to sterilise the restorations components prior to use:

1. prevacuum sterilisation method: steam sterilise the abutments at 132°C (270°F) at 180 - 220 kPa for 4 minutes. Dry for at least 20 minutes in the chamber. Only an approved wrap or pouch for steam sterilisation must be used.
2. for users in the USA: prevacuum sterilisation method: wrapped, steam sterilise at 135°C (275°F) at 180 - 220 kPa for 3 minutes. Dry for 20 minutes in the chamber. Use a wrap or pouch that is cleared for the indicated steam sterilisation cycle.

NOTE: users in the USA must ensure that the steriliser, wrap or pouch, and all steriliser accessories are cleared by the FDA, for the intended sterilisation cycle.

Contraindications

The contraindications of all device groups used as part of the specific treatment or procedure apply. Therefore, the contraindications of the systems/medical devices utilized as part of implant surgery/therapy should be noted and the relevant documents consulted.

The contraindications specific to this device group include:

- Allergies or hypersensitivity to chemical ingredients in the following materials: Polyether Ether Ketone (PEEK), Titanium, Aluminium, Vanadium.

Warnings and precautions

IMPORTANT NOTICE: THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING

- To ensure the safe and effective use of dental implants, new technologies/systems and the Metal Abutment devices, it is strongly recommended that specialised training be undertaken. This training should include hands-on methods to gain competency on the proper technique, the system's biomechanical requirements and radiographic evaluations required for the specific system.

- Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone.
- Use of the device with incompatible or non-corresponding devices can result in the poor performance or failure of the device.
- When handling devices intraorally, it is imperative that they are adequately secured to prevent aspiration, as aspiration of products may lead to infection or physical injury.
- The use of non-sterile items can lead to secondary infections of the tissue or transfer infectious diseases.
- Failure to adhere to appropriate cleaning, re-sterilization, and storage procedures as outlined in the Instructions for Use can result in device damage, secondary infections, or patient harm.
- Exceeding the number of recommended uses for reusable devices can result in device damage, secondary infection or patient harm.
- The use of blunt drills may cause damage to the bone, potentially compromising osseointegration.

It is crucial to emphasize that training should be undertaken by both new and experienced implant users prior to utilizing a new system or attempting a new treatment method.

Patient Selection and Preoperative Planning

A comprehensive patient selection process and meticulous preoperative planning are essential for successful implant treatment. This process should involve consultation among a multidisciplinary team, including well-trained surgeons, restorative dentists, and laboratory technicians.

Patient screening should include, at a minimum, a thorough medical and dental history, as well as visual and radiological inspections to assess the presence of adequate bone dimensions, the positioning of anatomical landmarks, the presence of unfavourable occlusal conditions, and the periodontal health status of the patient.

For successful implant treatment, it is important to:

1. Minimize trauma to the host tissue, thereby increasing the potential for successful osseointegration.
2. Accurately identify measurements relative to radiographic data, as failure to do so may lead to complications.
3. Be vigilant in avoiding damage to vital anatomical structures, such as nerves, veins, and arteries. Injury to these structures may result in serious complications, including ocular injury, nerve damage, and excessive bleeding.

The responsibility for proper patient selection, adequate training and experience in implant placement, and the provision of appropriate information required for informed consent rests with the practitioner. By combining thorough screening of prospective implant candidates with a practitioner possessing a high level of competence in the use of the system, the potential for complications and severe side effects can be significantly reduced.

High Risk Patients

Special care should be exercised when treating patients with local or systemic risk factors that may adversely affect the healing of bone and soft tissue or otherwise increase the severity of side effects, the risk of complications, and/or the likelihood of implant failure. Such factors include:

- poor oral hygiene
- history of smoking/vaping/tobacco use
- history of periodontal disease
- history of orofacial radiotherapy**
- bruxism and unfavourable jaw relations
- use of chronic medications that may delay healing or increase the risk of complications including, but not limited to, chronic steroid therapy, anti-coagulant therapy, TNF- α blockers, bisphosphonate and cyclosporin

*** The potential for implant failure and other complications increases when implants are placed in irradiated bone, as radiotherapy can lead to progressive fibrosis of blood vessels and soft tissue (i.e., osteoradionecrosis), resulting in diminished healing capacity. Contributing factors to this increased risk include the timing of implant placement in relation to radiation therapy, the proximity of radiation exposure to the implant site, and the radiation dosage at that site.*

Should the device not operate as intended, it must be reported to the manufacturer of the device. The contact information or the manufacturer of this device to report a change in performance is: sicomplaints@southernimplants.com.

Side effects

The clinical outcome of treatment is influenced by various factors. The following side effects and residual risks are associated with the device group and may necessitate further treatment, revision surgery or additional visits to the relevant medical professional's office. Furthermore, these side effects and residual risks can occur with varying possible severities and frequencies.

Therefore, the full list of known potential undesirable side-effects and residual risks that have been identified in relation to the PEEK Abutments include:

- Abutment fracture
- Allergic reaction(s) to the abutment material
- Bleeding
- Bruising
- Complications requiring revision surgery
- Discomfort
- Gingival inflammation
- Gingival recession
- Hematoma or bruising
- Hyperplastic tissue response
- Infection (acute and/or chronic)
- Localized inflammation
- Misfit or improper connection at the implant-abutment interface
- Overloading of the abutment/implant
- Pain, tenderness or discomfort
- Peri-implantitis, peri-mucositis or otherwise poor peri-implant soft tissue health
- Periodontal inflammation
- Phonetic difficulties
- Soft tissue irritation
- Wound dehiscence or poor healing

Additionally, the normal side effects associated with anaesthesia should also be expected.

Precaution: maintaining sterility protocol

The PEEK Abutments are packed in a peel pouch or blister base with a "peel-back" lid. Labelling information is located on the bottom half of the pouch, inside the packet or on the surface of the peel-back lid. Sterility is assured unless the pouch or blister is damaged or opened.

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established.

The contact information for the manufacturer of this device to report a serious incident is as follows:
sicomplaints@southernimplants.com.

Materials

Abutment Material Type	White Medical Grade polyether ether ketone (PEEK)
Abutment Screws Material Type	Titanium alloy (Ti-6Al-4V) according to ASTM F136 and ISO 5832-3

Titanium alloy (Ti-6Al-4V) according to ASTM F136 and ISO 5832-3

Chemical Components	Aluminium (Al)	Vanadium (V)	Residuals* (Fe, O, C, N, H) in total	Titanium (Ti)
Composition, % (mass/mass)	5.50 – 6.75	3.50 – 4.50	<0.55	Balance

*Where, Fe = Iron, O = Oxygen, C= Carbon, N = Nitrogen, H = Hydrogen

Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

Magnetic Resonance (MR) Safety

MR safe: the Southern Implants® PEEK abutments (excluding the prosthetic screws) are fabricated from a material that is not affected by exposure to MRI energy and are MR safe. However, restorations fabricated from this material may be assembled with metal implants, abutments and screws that can be affected by MRI energy.

Summary of Safety and Clinical Performance (SSCP)

As required by the European Medical Device Regulation (MDR; EU2017/745), a Summary of Safety and Clinical Performance (SSCP) is available for perusal with regard to Southern Implants® product ranges.

The relevant SSCP can be accessed at <https://ec.europa.eu/tools/eudamed>.

NOTE: the above website will be available upon the launch of the European Database on Medical Devices (EUDAMED).

Disclaimer of liability

This product is part of the Southern Implants® product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants® product range and take full responsibility for the correct indications and use of this product. Southern Implants® does not assume liability for damage due to incorrect use. Please note that some Southern Implants® products may not be cleared or released for sale in all markets.

Basic UDI

Product	Basic-UDI Number
Basic-UDI for PEEK Abutments	6009544038749A
Basic-UDI for Prosthetic Screws Titanium	60095440501886

Related literature and catalogues

CAT-2004 - Tri-Nex Implants Product Catalogue
 CAT-2005 - IT Implants Product Catalogue
 CAT-2020 - External Hex Implants Product Catalogue
 CAT-2042 - Deep Conical Implants Product Catalogue
 CAT-2043 - Internal Hex Implants Product Catalogue
 CAT-2060 - PROVATA® Implants Product Catalogue
 CAT-2069 - INVERTA® Implants Product Catalogue
 CAT-2070 - Zygomatic Implants Product Catalogue
 CAT-2092 - Soft Bone Implants Product Catalogue
 CAT-2093 - Single Platform Implants Product Catalogue
 CAT-2095 - External Hex INVERTA Implants Product Catalogue
 CAT-2096 - External Hex PTERYGOID Implants Product Catalogue

Symbols and warnings

											
Manufacturer: Southern Implants® 1 Albert Rd, P.O Box 605 IRENE, 0062, South Africa. Tel: +27 12 667 1046	CE mark 2797	Prescription device*	Sterilised using irradiation	Non-sterile	Use by date (mm-yy)	Do not reuse	Do not resterilise	Catalogue number	Batch code	Medical device	Authorised representative in the European Community
											
Authorised representative for Switzerland	Date of manufacture	Magnetic Resonance conditional	Magnetic Resonance safe	Single sterile barrier system with protective packaging inside	Single sterile barrier system	Consult instruction for use	Caution	Keep away from sunlight	Do not use if package is damaged		

* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

All rights reserved. Southern Implants®, the Southern Implants® logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants®. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale. It is the responsibility of the clinician to inspect the symbols that appear on the packaging of the product in use.

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