English	INSTRUCTIONS FOR USE: Southern Implants® CAB abutments
Español	INSTRUCCIONES DE USO: Pilares CAB Southern Implants®
Italiano	ISTRUZIONI PER L'USO: Southern Implants® Monconi CAB
Français	MODE D'EMPLOI : Piliers CAB Southern Implants®
Deutsch	GEBRAUCHSANWEISUNG: Southern Implants® CAB-Abutments
Português	INSTRUÇÕES DE UTILIZAÇÃO: Southern Implants® Pilares CAB
Română	INSTRUCȚIUNI DE UTILIZARE: Bonturi Southern Implants® CAB

South Africa - Headquarters: 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

Subsidiaries

Australia

Southern Implants Australia
T: +61-(0)-8-9466-2627
E: info@southernimplants.com.au

Spain and Portugal

Southern Implants Iberica
T: +34 935 053 507
E: info@southernimplants.es

United Kingdom and Ireland

Southern Implants UK
T: +44-20-8899-6845 / 6 / 7
E: info@southernimplants.co.uk

USA and Canada

Southern Implants North America Inc.
T: +1-561-472-0990
E: customercare@southernimplants.com

Intended use

Custom Abutment Base (CAB) abutments are intended to be connected to an implanted dental implant to serve as support for a dental prosthesis.

Description

Custom Abutment Base (CAB) abutments are pre-manufactured components that are connected direct to an endosseous implant for use as an aid in prosthetic rehabilitation.

The CAB abutments are available in the DC, Internal Hex and PROVATA® product ranges. The CAB abutments are provided sterile.

NOTE: screws are sold separately for all abutments. CAB abutments are anodised yellow.

Indications for use

Southern Implants' dental implant systems are intended for surgical placement in the upper or lower jaw to provide a means for prosthetic attachment of crowns, bridges or overdentures utilizing delayed or immediate loading. Southern Implants' dental implants are intended for immediate function when good primary stability with appropriate occlusal loading is achieved.

Intended user

Dental Technicians, Maxillo-facial Surgeons, General Dentists, Orthodontists, Periodontists, Prosthodontists and other appropriately trained and experienced implant users.

Intended environment

The devices are intended to be used in a clinical environment such as an operating theater or a dentist consultation room.

Intended patient population

Patients that have lost one tooth or multiple teeth.

Compatibility information

The CAB abutments are available for Deep Conical, Internal Hex and PROVATA® as indicated in Table A.

Table A

Implant connection type	Compatible device
Deep Conical (DC)	Parts labelled CAB-DC(Ø) for engaging items
Internal Hex (M)	Parts labelled CAB-M (used with Ø3.75, 4.20 and 5.00 mm platforms) for engaging items
	Parts labelled CAB-NM (used with Ø3.75, 4.20 and 5.00 mm platforms) for non-engaging items
PROVATA® (M) (Z)	Parts labelled CAB-M (used with Ø4.0, 5.0 and 6.0 mm platforms) for engaging items
	Parts labelled CAB-NM (used with Ø4.0, 5.0 and 6.0 mm platforms) for non-engaging items
	Parts labelled CAB-M-P45 (used with Ø5.0 and 6.0 mm platform-matched) for engaging items
	Parts labelled CAB-Z (used with Ø7.0, 8.0 and 9.0 mm platforms) for engaging items

Clinical benefits

Patients can expect to have their missing teeth replaced and/or crowns restored.

Before surgery

All components, instruments and tooling used during the clinical or laboratory procedure must be maintained in good condition and care must be taken that instrumentation does not damage implants or other components.

During surgery

Care must be taken that parts are not swallowed during any of the procedures, a rubber-dam application is recommended when appropriate. Care must be taken to apply the correct tightening torque of abutments and abutment screws.

Post-surgery

Regular patient follow-up and proper oral hygiene must be achieved to ensure favourable long-term results.

Digital workflow by using SIDIGITAL libraries**Scanning procedure**

1. Download Southern Implants digital library for 3Shape, Dental Wings and Exocad by registering on www.southernimplants.com.
2. Load the libraries into the CAD/CAM system by following their instructions.
3. The implant position is obtained digitally by an intraoral scan of the patient with the scan flag attached to the endosseous implant or a desktop scan of the dental model with the scan flag attached to the laboratory analogue.
4. Remove the scan flag from the implants or from the model. Replace healing abutments.
5. The scan is then imported into the design software.

Designing procedures

1. The scan flag in the digital form is now matched and aligned with the corresponding scan flag in the library.
2. The software recognises the position of the scan flag to the implant or analogue.
3. The CAB abutment is chosen from the library.
4. The software will guide the user through the steps to complete the restoration.
5. Selection of the restorative material (the most common material to use is Zirconia).

Milling and sintering procedures

1. Follow the instructions for use of the CAD/CAM system and milling material being used.
2. The milled restoration is cemented to the prefabricated titanium CAB abutment. Close the screw channel prior to cementing to keep the screw channel free of cement.

Digital workflow by scanning directly without using digital libraries**Scanning procedures**

1. Attach the CAB abutment to the lab analogue on the master model and screw down with the appropriate laboratory screw and driver.
2. Temporarily cover the screw hole with wax.
3. The CAB abutment is anodised gold and may not need scanning powder/spray to properly scan.
4. Scan using usual scan routines and design procedures.

Designing and milling procedures

1. Follow the instructions for use of the CAD/CAM system and milling material being used.
2. The milled restoration is cemented to the prefabricated titanium CAB abutment. Close the screw channel temporarily to keep the screw channel free of cement.

After casting procedure

Before sending the restoration and model to the clinician, clean the model thoroughly with soap, water and a toothbrush or steam clean. Then clean the abutment crown and its screw, preferably in the ultrasonic bath. Assemble the crown to the model and send it to the clinician.

Clinical procedures (placing restoration)

The clinician receives the restoration from the laboratory.

1. Remove the healing abutment or temporary restoration.
2. Clean, disinfect and sterilise the restoration as described.
3. Insert the restoration into the patient's mouth.
4. Position the restoration onto the implant, making sure that the retentive elements of the implant/abutment connections are properly aligned.
5. Fix the abutment to the implant with the correct screw using applicable driver (Table B). Torque the screw down to the value indicated in Table C.

Table B

Driver type	DC	Internal Hex & PROVATA®
1.22 mm/1.27 mm Universal driver	✓	✓
1.22 mm hex driver	✓	

1.27 mm hex driver		✓
Quad driver		Gold screws only

Table C

Direct to Implant	Torque
DC	
Ø3.0 mm	15 Ncm
Ø3.5 and 4.0 mm	20 Ncm
Ø5.0 mm	25-32 Ncm
Internal Hex (M-Series & PROVATA®)	
Ø3.75, 4.2 and 5.0 mm M-Series Implant	32 Ncm
Ø4.0, 5.0, 6.0, 7.0, 8.0 and 9.0 mm PROVATA® Implant	32 Ncm

- Verify the correct seating of the restoration using radiographic image.
- Do not exceed the recommended torque value as this may result in failure of the screw, abutment or implant. Do not torque less than the recommended value, this may result in loosening of the abutment that can lead to abutment or implant failure.
- Close the screw access hole.

Storage, cleaning and sterilisation

The implants, cover screws and healing abutments are supplied sterile (sterilised by gamma irradiation) and intended for single-use prior to the expiration date (see packaging label). Sterility is assured unless the container or seal is damaged or opened. If packaging is damaged, do not use the product and contact your Southern representative or return to Southern Implants. Do not reuse implants, cover screws, temporary abutments or abutments. Re-using these components may result in:

- damage to the surface or critical dimensions, which may result in performance and compatibility degradation.
- adds the risk of cross-patient infection and contamination if single-use items are reused.

Southern Implants does not accept any responsibility for complications associated with reused components.

Cleaning and disinfection

An implant restoration is a single or multiple-tooth implant crown, bridge or substructure, attached to a Southern Implants abutment or multiple abutments.

Before intraoral use the final restoration needs to be cleaned and disinfected, as per restorative material manufacturer's instructions.

Sterilisation

Southern Implants recommends the following procedure to sterilise the restoration prior to use:

Methods to sterilise the restoration and abutment screw:

- Pre-vacuum sterilisation method: steam sterilise the abutments at 132°C (270°F) at 180-220 kPa for 4 minutes. Dry for at least 20 minutes in the chamber. Only an approved wrap or pouch for steam sterilisation must be used.
- Pre-vacuum sterilisation method: wrapped, steam sterilise at 135°C (275°F) for 3 minutes. Dry for 20 minutes in the chamber. Use a wrap or pouch that is cleared for the indicated steam sterilisation cycle.

NOTE: Users in the USA must ensure that the steriliser, wrap or pouch and all steriliser accessories are cleared by the FDA, for the intended sterilisation cycle.

Contraindications

Do not use in patients:

- who are medically unfit for dental implant procedures.
- where adequate numbers of implants could not be placed to achieve full functional support of the prosthesis.
- who are allergic or have hypersensitivity to pure titanium or titanium alloy (Ti-6Al-4V), gold, palladium, platinum or iridium.
- who are under the age of 18, have poor bone quality, blood

disorders, infected implant site, vascular impairment, uncontrolled diabetes, drug or alcohol abuse, chronic high dose steroid therapy, anti-coagulant therapy, metabolic bone disease, radiotherapy treatment or sinus pathology.

Warnings

THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING.

- For the safe and effective use of dental implants it is suggested that specialised training be undertaken, including hands-on training to learn proper technique, biomechanical requirements and radiographic evaluations.
- Responsibility for proper patient selection, adequate training, experience in the placement of implants and providing appropriate information for informed consent rests with the practitioner. Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone.

Cautions

New and experienced implant users should do training before using a new system or attempting to do a new treatment method. Take special care when treating patients who have local or systemic factors that could affect the healing of the bone and soft tissue (i.e. poor oral hygiene, uncontrolled diabetes, are on steroid therapy, smokers, infection in the nearby bone and patients who had oro-facial radiotherapy). Thorough screening of prospective implant candidates must be performed including:

- a comprehensive medical and dental history.
- visual and radiological inspection to determine adequate bone dimensions, anatomical landmarks, occlusal conditions and periodontal health.
- bruxism and unfavourable jaw relations must be taken into account.
- proper pre-operative planning with a good team approach between well trained surgeons, restorative dentists and lab technicians is essential for successful implant treatment.
- minimizing the trauma to the host tissue increases the potential for successful osseointegration.
- electro-surgery should not be attempted around metal implants, as they are conductive.

Single use devices

Do not reuse devices indicated for single use (use the device prior to the expiration date).

Side effects

Potential side effects and temporary symptoms: pain, swelling, phonetic difficulties, gingival inflammation. More persistent symptoms: the risks and complications with implants include, but are not limited to: (1) allergic reaction(s) to implant and/or abutment material; (2) breakage of the implant and/or abutment; (3) loosening of the abutment screw and/or retaining screw; (4) infection requiring revision of the dental implant; (5) nerve damage that could cause permanent weakness, numbness, or pain; (6) histologic responses possibly involving macrophages and/or fibroblasts; (7) formation of fat emboli; (8) loosening of the implant requiring revision surgery; (9) perforation of the maxillary sinus; (10) perforation of the labial and lingual plates; and (11) bone loss possibly resulting in revision or removal.

Breakage

Implant and abutment fractures can occur when applied loads exceed the tensile or compressive strength of the material. Potential overloading conditions may result from: deficiencies in implant numbers, lengths and/or diameters to adequately support a restoration, excessive cantilever length, incomplete abutment seating, abutment angles greater than 30 degrees, occlusal interferences causing excessive lateral forces, patient parafunction (e.g. bruxing, clenching), loss or changes in dentition or functionality, inadequate prosthesis fit and physical trauma. Additional treatment may be necessary when any of the above conditions are present to reduce the possibility of hardware complications or failure.

Changes in performance

It is the responsibility of the clinician to instruct the patient on all appropriate contraindications, side effects and precautions as well as the

need to seek the services of a trained dental professional if there are any changes in the performance of the implant (e.g. looseness of the prosthesis, infection or exudate around the implant, pain or any other unusual symptoms that the patient has not been told to expect).

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established. The contact information for the manufacturer of this device to report a serious incident is as follows: sicomplaints@southernimplants.com

Materials

CAB abutment: Titanium (grade 4, ASTM F67)
Screws: Titanium alloy (Ti-90%, Al-6%, V-4%)
Gold alloy (Au-61%, Ag-16.5%, Pt-13.5%, Cu -9%)

Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

Disclaimer of liability

This product is part of the Southern Implants product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants product range and take full responsibility for the correct indications and use of this product. Southern Implants does not assume liability for damage due to incorrect use. Please note that some Southern Implants products may not be cleared or released for sale in all markets.

Healing

The healing time required for osseointegration depends on the individual and treatment protocol. It is the responsibility of the practitioner to decide when the implant can be restored. Good primary stability will govern if immediate loading can be done.

Implant care and maintenance

Potential implant patients should establish an adequate oral hygiene regime prior to implant therapy. Proper post-operative oral hygiene and implant maintenance instructions must be discussed with the patient as this will determine the longevity and health of the implants. The patient should maintain regular prophylaxis and evaluation appointments.

MR conditional

Non-clinical testing has demonstrated that the Southern Implants dental implants, metallic abutments and prosthetic screws are MR conditional. A patient with these devices can be safely scanned in a MR system meeting the following conditions:

- Static magnetic field of 1.5 Tesla and 3.0 Tesla only.
- Maximum spatial gradient magnetic field of 4000 Gauss/cm (40 T/m).
- Maximum MR system reported, head specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode) or whole body averaged specific absorption rate (wbSAR) of 1 W/kg.

Under the scan conditions defined above, the Southern Implants dental implants, abutments and prosthetic screws are expected to produce a maximum temperature rise of 5.8°C after 15 minutes of continuous scanning. In the non-clinical testing, the image artifact caused by the device extends approximately 20 mm from the Southern Implants dental implants, abutments and prosthetic screws, when imaged with a gradient echo pulse sequence and a 3.0 Tesla MRI system. Removable restorations should be taken out prior to scanning, as is done for watches, jewellery etc. Should there be no MR symbol on the product label, please note that this device has not been evaluated for safety and compatibility in the MR environment. This device has not been tested for heating or migration in the MR environment.

Basic UDI

Product	Basic-UDI Number
Basic-UDI for Metal Abutments	600954403872

Related literature and catalogues

CAT-2042 - Deep Conical Implants Product Catalogue
CAT-2043 - Internal Hex Implants Product Catalogue
CAT-2060 - PROVATA® Implants Product Catalogue

Symbols and Warnings

 Manufacturer: Southern Implants 1 Albert Rd, P.O Box 605 IRENE, 0062,	 2797	 Prescription device*	 Sterilization using instruction Device Irradiation	 Consult (mm-yy) for use	 Use by date	 Do not reuse	 Batch if package is damaged	 Medic	 Authorize Representative in the European Community	 Catalog Number	 Date of Manufacture
---	---	---	---	--	--	---	--	---	---	---	--

* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.

All rights reserved. Southern Implants®, the Southern Implants logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

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