
English	INSTRUCTIONS FOR USE: Southern Implants® Single Platform (SP1) implants
Español	INSTRUCCIONES DE USO: Implantes de Plataforma única (SP1) Southern Implants®
Italiano	ISTRUZIONI PER L'USO: Southern Implants® Impianti per Piattaforma Singola (SP1)
Français	MODE D'EMPLOI : Implants Southern Implants® Single Platform (SP1)
Deutsch	GEBRAUCHSANWEISUNGEN: Southern Implants® Single Platform (SP1)-Implantate
Português	INSTRUÇÕES DE UTILIZAÇÃO: Southern Implants® Implantes de plataforma única (SP1)
Ελληνικά	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ: Εμφυτεύματα Ενιαίας Πλατφόρμας (SP1) της Southern Implants®
Svenska	ANVISNINGAR FÖR ANVÄNDNING: Southern Implants® Enkel Plattform (SP1) Implantat
Română	INSTRUCȚIUNI DE UTILIZARE: Southern Implants® Implanturi Platformă Unică (SP1)
Türkçe	KULLANMA TALİMATI: Southern Implants® Tek Platform (SP1) implantları

South Africa - Headquarters: 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

Subsidiaries

Australia

Southern Implants Australia
T: +61-(0)-8-9466-2627
E: info@southernimplants.com.au

Spain and Portugal

Southern Implants Iberica
T: +34 935 053 507
E: info@southernimplants.es

United Kingdom and Ireland

Southern Implants UK
T: +44-20-8899-6845 / 6 / 7
E: info@southernimplants.co.uk

USA and Canada

Southern Implants North America Inc.
T: +1-561-472-0990
E: customercare@southernimplants.com

Description

The Single Platform (SP1) implant is a self-tapping tapered implant made of commercially pure Grade 4 Titanium (UTS $\geq$ 900 MPa). All implants are surface-roughened using Southern Implants' proven surface. The surface has an average Sa value of 1.4 microns. Tapered implants facilitate good stability for cases involving immediate loading and/or soft bone. Cover screws and healing abutments are supplied separately.

Intended use

Southern Implants® Dental Implants are intended to be surgically placed in the upper or lower jaw to provide a support structure and point of attachment for prosthetic devices.

Indications for use

The Southern Implants® Dental Implants are indicated for the functional and aesthetic restoration of the upper or lower jaw in cases of partial or full edentulousness, congenital tooth loss, failing/problematic natural dentition, or as part of the revision/replacement of previous dental restorations.

The Southern Implants® Dental Implants are suitable for immediate loading in all indications where sufficient primary stability is achieved with appropriate occlusal loading, except in single tooth restorations on implants shorter than 8 mm or in soft bone (type IV) where implant stability may be difficult to obtain and/or immediate loading deemed inappropriate.

Intended user

Maxillofacial Surgeons, General Dentists, Orthodontists, Periodontists, Prosthodontists, and other appropriately trained and/or experienced medical professionals.

Intended environment

The devices are intended to be used in a clinical environment such as an operating theatre or a dentist consultation room.

Intended patient population

The intended patient population for dental implant therapy includes partially or fully edentulous patients eligible or, otherwise, not contraindicated for implant placement.

Compatibility information

Southern Implants' implants should be restored with Southern Implants' components. Use components that correspond to the connection type and prosthetic platform when restoring the Single Platform (SP1) implants. For further information see the Single Platform (SP1) Product Catalogue (CAT-2093).

Table A – Single Platform (SP1) implant range

Item Code	Implant Length Codes	Implant Lengths (mm)	Implant Diameter (mm)
MSC-SP35	08/ 10/ 11/ 13/ 16/ 18/ 20/ 22/ 24	8.0/ 10.0/ 11.5/ 13.0/ 16.0/ 18.0/ 20.0/ 22.0/ 24.0	Ø3.5
MSC-SP37	08/ 10/ 11/ 13/ 16/ 18/ 20/ 22/ 24	8.0/ 10.0/ 11.5/ 13.0/ 16.0/ 18.0/ 20.0/ 22.0/ 24.0	Ø3.75
MSC-SP40	08/ 10/ 11/ 13/ 16/ 18/ 20/ 22/ 24	8.0/ 10.0/ 11.5/ 13.0/ 16.0/ 18.0/ 20.0/ 22.0/ 24.0	Ø4.0
PT-SP40	10/ 11/ 13/ 15/ 18/ 20/ 22/ 24	10.0/ 11.5/ 13.0/ 15.0/ 18.0/ 20.0/ 22.0/ 24.0	Ø4.0
MSC-SP45	08/ 10/ 11/ 13/ 16/ 18	8.0/ 10.0/ 11.5/ 13.0/ 16.0/ 18.0	Ø4.5
MSC-SP50	08/ 10/ 11/ 13/ 16/ 18	8.0/ 10.0/ 11.5/ 13.0/ 16.0/ 18.0	Ø5.0
MSC-SP60	08/ 10/ 11/ 13/ 16	8.0/ 10.0/ 11.5/ 13.0/ 16.0	Ø6.0

Clinical performance

The Single Platform (SP1) Implants are intended to provide a support structure and point of attachment for prosthetic devices in the functional and aesthetic treatment of the partially or fully edentulous upper or lower jaw. The rigid and secure retention of the prosthetic devices are accomplished through the stability of the implant itself. This stability is achieved through the implant's surgical implantation into the bone and subsequent successful osseointegration. Consequently, the clinical performance of the Single Platform (SP1) Implants is primarily defined by the degree of osseointegration achieved.

Clinical benefits

The clinical benefits associated with the prosthetic restoration of edentulism are significant, as this treatment can greatly enhance the functional and aesthetic outcomes of the patient. Therefore, the clinical benefit of the Single Platform (SP1) Implants encompasses the successful restoration of the edentulous state and the patient's satisfaction with the outcome. Furthermore, since this treatment addresses the functional and/or aesthetic impact of partial or full edentulism, it is anticipated that a successfully restored patient will experience an improvement in masticatory function and improved psychosocial wellbeing.

Storage, cleaning and sterilisation

The component is supplied sterile (sterilised by gamma irradiation) and intended for single-use prior to the expiration date (see packaging label). Sterility is assured unless the container or seal is damaged or opened. If packaging is damaged do not use the product and contact your Southern representative or return to Southern Implants®. The devices must be stored in a dry place at room temperature and not exposed to direct sunlight. Incorrect storage may influence device characteristics. Do not reuse components indicated for single-use only. Reusing these components may result in:

- damage to the surface or critical dimensions, which may result in performance and compatibility degradation.
- adds the risk of cross-patient infection and contamination if single-use items are reused.

Southern Implants® does not accept any responsibility for complications associated with reused single-use components.

Contraindications

The contraindications of all device groups used as part of the specific treatment or procedure apply. Therefore, the contraindications of the systems/medical devices utilized as part of implant surgery/therapy should be noted and the relevant documents consulted.

The contraindications specific to this device group include:

- Patients medically unfit for oral surgical procedures,
- cases where an adequate number of implants cannot be placed to provide full functional support of the prosthesis,
- patients with inadequate bone quality/quantity,
- patients with uncontrolled bleeding disorders,
- presence of infection at the implant site,
- patients with severe vascular impairment,
- patients with uncontrolled endocrine disorders,
- patients undergoing high dose steroid therapy,
- patients with metabolic bone disease,
- patients with incomplete mandibular or maxillary growth,
- patients with a weakened immune system or inadequate wound healing capacity
- patients with allergies or hypersensitivity to the material used (titanium).

Warnings and precautions

IMPORTANT NOTICE: THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING

- To ensure the safe and effective use of dental implants, new technologies/systems and the Metal Abutment devices, it is strongly recommended that specialised training be undertaken. This training should include hands-on methods to gain competency on the proper technique, the system's biomechanical requirements and radiographic evaluations required for the specific system.
- Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone.
- Use of the device with incompatible or non-corresponding devices can result in the poor performance or failure of the device.
- When handling devices intraorally, it is imperative that they are adequately secured to prevent aspiration, as aspiration of products may lead to infection or physical injury.
- The use of non-sterile items can lead to secondary infections of the tissue or transfer infectious diseases.

- Failure to adhere to appropriate cleaning, re-sterilization, and storage procedures as outlined in the Instructions for Use can result in device damage, secondary infections, or patient harm.
- Exceeding the number of recommended uses for reusable devices can result in device damage, secondary infection or patient harm.
- The use of blunt drills may cause damage to the bone, potentially compromising osseointegration.

It is crucial to emphasize that training should be undertaken by both new and experienced implant users prior to utilizing a new system or attempting a new treatment method.

Patient Selection and Preoperative Planning

A comprehensive patient selection process and meticulous preoperative planning are essential for successful implant treatment. This process should involve consultation among a multidisciplinary team, including well-trained surgeons, restorative dentists, and laboratory technicians.

Patient screening should include, at a minimum, a thorough medical and dental history, as well as visual and radiological inspections to assess the presence of adequate bone dimensions, the positioning of anatomical landmarks, the presence of unfavourable occlusal conditions, and the periodontal health status of the patient.

For successful implant treatment, it is important to:

1. Minimize trauma to the host tissue, thereby increasing the potential for successful osseointegration.
2. Accurately identify measurements relative to radiographic data, as failure to do so may lead to complications.
3. Be vigilant in avoiding damage to vital anatomical structures, such as nerves, veins, and arteries. Injury to these structures may result in serious complications, including ocular injury, nerve damage, and excessive bleeding.

The responsibility for proper patient selection, adequate training and experience in implant placement, and the provision of appropriate information required for informed consent rests with the practitioner. By combining thorough screening of prospective implant candidates with a practitioner possessing a high level of competence in the use of the system, the potential for complications and severe side effects can be significantly reduced.

High Risk Patients

Special care should be exercised when treating patients with local or systemic risk factors that may adversely affect the healing of bone and soft tissue or otherwise increase the severity of side effects, the risk of complications, and/or the likelihood of implant failure. Such factors include:

- poor oral hygiene
- history of smoking/vaping/tobacco use
- history of periodontal disease
- history of orofacial radiotherapy**
- bruxism and unfavourable jaw relations
- use of chronic medications that may delay healing or increase the risk of complications including, but not limited to, chronic steroid therapy, anti-coagulant therapy, TNF- α blockers, bisphosphonate and cyclosporin

*** The potential for implant failure and other complications increases when implants are placed in irradiated bone, as radiotherapy can lead to progressive fibrosis of blood vessels and soft tissue (i.e., osteoradionecrosis), resulting in diminished healing capacity. Contributing factors to this increased risk include the timing of implant placement in relation to radiation therapy, the proximity of radiation exposure to the implant site, and the radiation dosage at that site.*

Should the device not operate as intended, it must be reported to the manufacturer of the device. The contact information or the manufacturer of this device to report a change in performance is: sicomplaints@southernimplants.com.

Side effects

The clinical outcome of treatment is influenced by various factors. The following side effects and residual risks are associated with the device group and may necessitate further treatment, revision surgery or additional visits to the

relevant medical professional's office. Furthermore, these side effects and residual risks can occur with varying possible severities and frequencies.

The common side effects and residual risks that are associated with the device include:

- Pain or tenderness
- Localized inflammation
- Soft tissue irritation
- Discomfort
- Hematoma or bruising
- Bleeding

Additionally, the normal side effects associated with anaesthesia should also be expected.

The less common and potentially more severe side effects and residual risks that are associated with the device include:

- Allergic reaction(s) to the implant material
- Anaesthesia, paraesthesia, hyperesthesia, and hypoesthesia (transient or permanent)
- Anatomical landmark damage
- Bleeding on probing
- Bruising
- Complications requiring revision surgery
- Damage to adjacent dentition
- Dental injury during surgery
- Hyperplastic soft tissue response
- Implant failure
- Implant fracture
- Improper implant positioning resulting in prosthetic compromise
- Infection (acute and/or chronic)
- Insufficient levels of osseointegration resulting in instability, mobility and/or failure
- Localized inflammation
- Loosening of the abutment screw and/or retaining screw
- Loss or damage to adjacent teeth
- Marginal bone loss
- Micromovements and implant instability
- Misfit or improper connection at the implant-abutment interface
- Overloading of the abutment/implant
- Pain or discomfort
- Peri-implantitis, peri-mucositis or otherwise poor peri-implant soft tissue health
- Periodontal inflammation
- Phonetic difficulties
- Prosthetic failure
- Soft tissue irritation
- Sub-optimal aesthetic result
- Surgical side effects such as pain, inflammation, bruising and mild bleeding
- Wound dehiscence or poor healing

Note: these side effects and residual risks can vary in both severity and frequency.

Precaution: maintaining sterility protocol

Implants are packaged as follows:

1. An outer package consisting of a rigid, clear box which acts as protection for the inner package.
2. The inner package consisting of a blister pack (clear plastic-formed blister base with a TYVEK "peel-back" lid).
3. Within the inner package, there is a hollow tube which contains one implant suspended from a titanium ring, this ensures the implant never touches the inside of the plastic tube.

- Labelling information is located on the surface of the peel-back lid and on the outside of the rigid box.

Care must be taken to maintain the sterility of the implant by proper opening of the packaging and handling of the implant.

- Open the implant package in the non-sterile field, with non-sterile gloves, tear the address label to open the box.
- With non-sterile gloves, remove the inner blister pack. Do not place the plastic box or blister pack-lid onto the sterile field. The contents of this inner package are sterile.
- The sealed blister is to be opened by an assistant (with nonsterile gloves): remove the TYVEK lid and drop or place the sterile tube onto the sterile field, open the tube cap and attach the implant placement tool onto the implant and carefully remove from the sterile tube. Do not touch the sterile implant.

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established.

The contact information for the manufacturer of this device to report a serious incident is as follows:
sicomplaints@southernimplants.com.

Materials

Material type Commercially pure titanium (grade 4, ASTM F67 and ISO5832-2, UTS $\geq$ 900 MPa)

Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

Magnetic Resonance (MR) Safety

Nonclinical testing has demonstrated that the Southern Implants® dental implants, metallic abutments and prosthetic screws are MR conditional.

MRI Safety Information	
 MR Conditional	
A patient with the dental implants, metallic abutments and prosthetic screws from Southern Implants® may be safely scanned under the following conditions. Failure to follow these conditions may result in injury to the patient.	
Device Name	dental implants, metallic abutments and prosthetic screws
Static Magnetic Field Strength (B₀) [T]	1.5T or 3T
Maximum Spatial Field Gradient [T/m] [gauss/cm]	30 T/m (3000 gauss/cm)
Type of RF Excitation	Circularly polarized (CP) (i.e., quadrature driven)
RF Transmit Coil Type	Integrated Whole Body Transmit Coil
RF Receive Coil	Any receive-only RF coil may be used
Operating Mode	Normal Operating Mode (See SAR, landmark restrictions below)
RF Conditions	Normal Operating mode for all landmarks: (Head SAR of 3.2 W/kg for head landmark, 2 W/kg whole body, and appropriate partial body SAR for other landmarks)

Scan Duration	For imaging landmarks above the diaphragm, a continuous scan time of 15 minutes will require a cooling delay of at least 5 minutes. No scantime restriction for landmarks below the diaphragm <i>Or</i> If the closest part of implant <25 cm away from the isocenter in feet-head direction: a continuous scan time of 15 minutes will require a cooling delay of at least 5 minutes. If the closest part of implant ≥ 25 cm away from the isocenter in feet-head direction: There is no limit on scan duration.
MR Image Artifact	The presence of this implant may produce an image artifact.

Additional information:

1. Local transmit coil: Local transmit coil safety was not evaluated. Local transmit coil is only allowed if the implant is ≥ 25 cm away from isocenter and for a wbSAR of up to 2W/kg.
2. Artifacts: MR Image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the implant. Tests following guidelines in ASTM F2119-24 showed artifact widths of ≤20 mm from the surface of the implant for the gradient echo sequence and ≤10 mm from the surface of the implant for the spin echo sequence.

Removable restorations should be taken out prior to scanning, as is done for watches, jewellery etc.

Should there be no MR symbol on the product label, please note that this device has not been evaluated for safety and compatibility in the MR environment.

Summary of Safety and Clinical Performance (SSCP)

As required by the European Medical Device Regulation (MDR; EU2017/745), a Summary of Safety and Clinical Performance (SSCP) is available for perusal with regard to Southern Implants® product ranges.

The relevant SSCP can be accessed at <https://ec.europa.eu/tools/eudamed>.

NOTE: the above website will be available upon the launch of the European Database on Medical Devices (EUDAMED).

Disclaimer of liability

This product is part of the Southern Implants® product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants® product range and take full responsibility for the correct indications and use of this product. Southern Implants® does not assume liability for damage due to incorrect use. Please note that some Southern Implants® products may not be cleared or released for sale in all markets.

Basic UDI

Product	Basic-UDI Number
Basic-UDI for General Dental Implants	6009544038699H

Related literature and catalogues

CAT-2093 - Single Platform (SP1) Implants Product Catalogue

Symbols and warnings

											
Manufacturer: Southern Implants® 1 Albert Rd, P.O Box 605 IRENE, 0062, South Africa. Tel: +27 12 667 1046	CE mark 2797	Prescription device*	Sterilised using irradiation	Non-sterile	Use by date (mm-yy)	Do not reuse	Do not resterilise	Catalogue number	Batch code	Medical device	Authorised representative in the European Community
											
Authorised representative for Switzerland	Date of manufacture	Magnetic Resonance conditional	Magnetic Resonance safe	Single sterile barrier system with protective packaging inside	Single sterile barrier system	Consult instruction for use	Caution	Keep away from sunlight	Do not use if package is damaged		

* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

All rights reserved. Southern Implants®, the Southern Implants® logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants®. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale. It is the responsibility of the clinician to inspect the symbols that appear on the packaging of the product in use.

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