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South Africa - Headquarters: 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

Subsidiaries

United Kingdom and Ireland
Southern Implants UK
T: +44-20 8059 4490
E: info@southernimplants.co.uk

USA and Canada
Southern Implants North America Inc.
T: +1-561-472-0990
E: customercare@southernimplants.com

EC REP

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

CH REP

MedEnvoy Switzerland: Gotthardstrasse 28, 6302 Zug, Switzerland

Description

The ZAGA® Screw Retained Abutment for Zygomatic implants is a multi-unit abutment that is connected directly to an endosseous implant and is used in multiple unit reconstructions when it is desirable to raise the prosthetic interface to a more apical position than that of an implant head for a screw retained restoration. This abutment is compatible with the Straumann® ZAGA™ Zygomatic Implant system for which it is designed. This abutment is available in a straight and 17° angled configuration. The ZAGA® Screw Retained Abutments for Zygomatic implants are only indicated for multi-unit cases. These abutments are provided sterile.

Intended use

The Southern Implants® Metal Abutments are prosthetic components intended for connection to endosseous implant systems to aid in prosthetic rehabilitation and implant procedures.

The ZAGA® Screw Retained Abutments for Zygomatic Implants are intended to provide a supportive base for connected prosthetic devices on endosseous zygomatic implants.

Indications for use

The Southern Implants® Metal Abutments are indicated for use as components of an endosseous implant system to support prosthetic rehabilitation in implant treatment procedures.

The ZAGA® Screw Retained Abutments for Zygomatic Implants are indicated for connection to implanted Straumann® ZAGA™ Zygomatic Implants in order to raise the prosthetic platform to a more coronal position as part of the functional and aesthetic restoration of the upper jaw using a fixed, **permanent** prosthesis.

Intended user

Dental technicians, Maxillofacial Surgeons, General Dentists, Orthodontists, Periodontists, Prosthodontists, and other appropriately trained and/or experienced medical professionals.

Intended environment

Straumann® ZAGA® Screw Retained Abutments are intended to be used in a dental laboratory as part of the restoration design and manufacture as well as in a clinical environment such as an operating theatre or a dentist consultation room.

Intended patient population

Southern Implants® Metal Abutments are intended to be used in patients indicated and eligible, or otherwise not contraindicated, for implant-retained prosthetic restorations.

Compatibility information

Use only original Southern Implants® components to restore Straumann® ZAGA™ Zygomatic ranges. Use components that correspond to the connection type, and prosthetic platform when restoring Zygomatic implants.

Table A

Straumann® Zygomatic ZAGA™ Implants				
Straumann® Zygomatic Implant ZAGA™ Flat	Straumann® Zygomatic Implant ZAGA™ Round	Cover Screw and Driver	Abutment and Driver	Prosthetic Screw and Driver
CH-ZC-30.0	CH-ZT-30.0	CH-CS (cover screw) I-CS-HD (driver)	CH-SRA-xx* (straight Screw-Retained Abutment)	I-HD-M (driver)
CH-ZC-32.5	CH-ZT-32.5			
CH-ZC-35.0	CH-ZT-35.0			
CH-ZC-37.5	CH-ZT-37.5			
CH-ZC-40.0	CH-ZT-40.0		046.401	
CH-ZC-42.5	CH-ZT-42.5		046.411 (driver)	
CH-ZC-45.0	CH-ZT-45.0			

CH-ZC-47.5	CH-ZT-47.5			
CH-ZC-50.0	CH-ZT-50.0			
CH-ZC-52.5	CH-ZT-52.5			
CH-ZC-55.0	CH-ZT-55.0			
CH-ZC-57.5	CH-ZT-57.5			
CH-ZC-60.0	CH-ZT-60.0			

*xx denotes the abutment collar height

Table B

Straumann® Zygomatic ZAGA™ Implants			
Straumann® Zygomatic Implant ZAGA™ Flat	Cover Screw and Driver	Abutment and Driver	Prosthetic Screw and Driver
CH-ZF-30.0	CH-CS (cover screw) I-CS-HD (driver)	CH-SRA-xx* (straight Screw-Retained Abutment) CH-SRA17-xx* (angled Screw-Retained Abutment) 046.401 046.411 (driver)	I-HD-M (driver)
CH-ZF-32.5			
CH-ZF-35.0			
CH-ZF-37.5			
CH-ZF-40.0			
CH-ZF-42.5			
CH-ZF-45.0			
CH-ZF-47.5			
CH-ZF-50.0			
CH-ZF-52.5			
CH-ZF-55.0			
CH-ZF-57.5			
CH-ZF-60.0			

*xx denotes the abutment collar height

First clinical procedure

1. Select and connect the appropriate abutment to the implant, using the dedicated abutment driver.
2. Place the abutment on the implant and tighten to the recommended torque (see Table C) with the appropriate driver and torque wrench (see Table A and Table B).
CAUTION: never exceed the recommended torque as overtightening may lead to screw fractures.
3. Verify the correct seating of the abutments using radiographic imaging.
4. Connect impression copings (Traditional: 025.0012, 025.0014; Digital: 025.0001) to the ZAGA® screw-retained abutment.
5. Take an open or closed tray impression and remove/transfer the impression copings to the impression. Alternatively, perform a digital scan of the impression copings, thereafter removing the digital impression copings.
6. Connect the temporary restoration directly to the ZAGA® screw-retained abutment.

Table C

Abutment Recommended Torque Values	
Abutment	Recommended Torque (Ncm)
CH-SRA-xx* (straight Screw-Retained Abutment)	35
CH-SRA17-xx* (angled Screw-Retained Abutment)	20

*xx denotes the abutment collar height

Second clinical procedure

1. Attach the prosthesis to the ZAGA® screw-retained abutment: place and tighten the restoration screw.
2. Verify the correct seating of the restoration using radiographic imaging.
3. Tighten the restoration using a manual torque wrench to 10-15 Ncm.
4. Close the screw access hole.
5. Screw-retain or cement final prosthesis, if applicable.

Clinical performance

The clinical performance of the Straumann® ZAGA® Screw Retained Abutments is primarily defined by their ability to maintain stability and retention of the prosthetic restoration over time. This stability is achieved through secure fixation of the restoration to the abutment, which can be accomplished via cementation or casting, as well as the mechanical fastening of the abutment to the implant through screws.

Clinical benefits

The clinical benefits associated with the Straumann® ZAGA® Screw Retained Abutments are intrinsically linked to the overall success of implant-supported prosthetic treatments and encompass both direct and indirect clinical benefits. The direct clinical benefits of the abutments arise specifically from their use in the treatment process. These benefits include the enhancement of peri-implant soft tissue health and the aesthetic outcomes of the final restoration. The clinical benefits experienced by patients as a result of the implant-supported prosthetic treatments in which these abutments are utilized represent the indirect clinical benefits of the Straumann® ZAGA® Screw Retained Abutments. While these benefits are substantial, they are not directly attributable to the devices themselves; rather, they reflect the contributions of the abutments to the overall treatment system. Additionally, successful treatment outcomes are anticipated to result in significant enhancements in patients' psychosocial well-being, self-esteem, and overall quality of life. Overall, the Straumann® ZAGA® Screw Retained Abutments not only facilitate the mechanical aspects of dental restorations but also contribute to the holistic improvement of patient outcomes and experiences.

Before surgery

All components, instruments and tooling used during the clinical or laboratory procedure must be maintained in good condition and care must be taken that instrumentation does not damage implants or other components.

During surgery

Care must be taken that parts are not swallowed during any of the procedures, a rubber-dam application is recommended when appropriate. Care must be taken to apply the correct tightening torque of abutments and abutment screws.

Post-surgery

Regular patient follow-up and proper oral hygiene must be achieved to ensure favourable long-term results.

Storage, cleaning and sterilisation

The component is supplied sterile (sterilised by gamma irradiation) and intended for single-use prior to the expiration date (see packaging label). Sterility is assured unless the container or seal is damaged or opened. If packaging is damaged do not use the product and contact your Straumann representative or return to Southern Implants®. The devices must be stored in a dry place at room temperature and not exposed to direct sunlight. Incorrect storage may influence device characteristics.

Do not reuse components indicated for single-use only. Reusing these components may result in:

- damage to the surface or critical dimensions, which may result in performance and compatibility degradation.
- adds the risk of cross-patient infection and contamination if single-use items are reused.

Southern Implants® does not accept any responsibility for complications associated with reused single-use components.

Contraindications

The contraindications of all device groups used as part of the specific treatment or procedure apply. Therefore, the contraindications of the systems/medical devices utilized as part of implant surgery/therapy should be noted and the relevant documents consulted.

The contraindications specific to this device group include:

- patients who are medically unfit for implant or prosthetic related procedures
- patients that present with an allergy or hypersensitivity to the following materials: Titanium, Aluminium, Vanadium, Cobalt Chrome, Cobalt Chrome Molybdenum, Gold, Platinum, Palladium, Iridium, Silver and/or Copper.

Warnings and precautions

IMPORTANT NOTICE: THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING

- To ensure the safe and effective use of dental implants, new technologies/systems and the Metal Abutment devices, it is strongly recommended that specialised training be undertaken. This training should include hands-on methods to gain competency on the proper technique, the system's biomechanical requirements and radiographic evaluations required for the specific system.
- Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone.
- Use of the device with incompatible or non-corresponding devices can result in the poor performance or failure of the device.
- When handling devices intraorally, it is imperative that they are adequately secured to prevent aspiration, as aspiration of products may lead to infection or physical injury.
- The use of non-sterile items can lead to secondary infections of the tissue or transfer infectious diseases.
- Failure to adhere to appropriate cleaning, re-sterilization, and storage procedures as outlined in the Instructions for Use can result in device damage, secondary infections, or patient harm.
- Exceeding the number of recommended uses for reusable devices can result in device damage, secondary infection or patient harm.
- The use of blunt drills may cause damage to the bone, potentially compromising osseointegration.

It is crucial to emphasize that training should be undertaken by both new and experienced implant users prior to utilizing a new system or attempting a new treatment method.

Patient Selection and Preoperative Planning

A comprehensive patient selection process and meticulous preoperative planning are essential for successful implant treatment. This process should involve consultation among a multidisciplinary team, including well-trained surgeons, restorative dentists, and laboratory technicians.

Patient screening should include, at a minimum, a thorough medical and dental history, as well as visual and radiological inspections to assess the presence of adequate bone dimensions, the positioning of anatomical landmarks, the presence of unfavourable occlusal conditions, and the periodontal health status of the patient.

For successful implant treatment, it is important to:

1. Minimize trauma to the host tissue, thereby increasing the potential for successful osseointegration.
2. Accurately identify measurements relative to radiographic data, as failure to do so may lead to complications.
3. Be vigilant in avoiding damage to vital anatomical structures, such as nerves, veins, and arteries. Injury to these structures may result in serious complications, including ocular injury, nerve damage, and excessive bleeding.

The responsibility for proper patient selection, adequate training and experience in implant placement, and the provision of appropriate information required for informed consent rests with the practitioner. By combining thorough

screening of prospective implant candidates with a practitioner possessing a high level of competence in the use of the system, the potential for complications and severe side effects can be significantly reduced.

High Risk Patients

Special care should be exercised when treating patients with local or systemic risk factors that may adversely affect the healing of bone and soft tissue or otherwise increase the severity of side effects, the risk of complications, and/or the likelihood of implant failure. Such factors include:

- poor oral hygiene
- history of smoking/vaping/tobacco use
- history of periodontal disease
- history of orofacial radiotherapy**
- bruxism and unfavourable jaw relations
- use of chronic medications that may delay healing or increase the risk of complications including, but not limited to, chronic steroid therapy, anti-coagulant therapy, TNF- α blockers, bisphosphonate and cyclosporin

*** The potential for implant failure and other complications increases when implants are placed in irradiated bone, as radiotherapy can lead to progressive fibrosis of blood vessels and soft tissue (i.e., osteoradionecrosis), resulting in diminished healing capacity. Contributing factors to this increased risk include the timing of implant placement in relation to radiation therapy, the proximity of radiation exposure to the implant site, and the radiation dosage at that site.*

Should the device not operate as intended, it must be reported to the manufacturer of the device. The contact information or the manufacturer of this device to report a change in performance is:

sicomplaints@southernimplants.com.

Side effects

The clinical outcome of treatment is influenced by various factors and there are additional potential side effects and residual risks associated with the device group. These may necessitate further treatment, revision surgery or additional visits to the relevant medical professional's office. Furthermore, these side effects and residual risks can vary in both severity and frequency.

Therefore, the full list of known potential undesirable side-effects and residual risks that have been identified in relation to the dental abutment device group include:

- Abutment fracture
- Abutment screw fracture or distortion
- Allergic reaction(s) to the abutment material
- Anaesthesia, paraesthesia, hyperesthesia, and hypoesthesia (transient or permanent)
- Bleeding
- Bleeding on probing
- Bruising
- Complications requiring revision surgery
- Difficulty in retrieving the abutment
- Gingival inflammation
- Gingival recession
- Hyperplastic tissue response
- Implant failure due to insufficient levels of osseointegration
- Improper implant positioning resulting in prosthetic compromise
- Infection (acute and/or chronic)
- Localized inflammation
- Loosening of the abutment screw and/or retaining screw
- Loss or damage to adjacent teeth

- Marginal bone loss
- Micromovements and implant instability
- Misfit or improper connection at the implant-abutment interface
- Overloading of the abutment/implant
- Pain, tenderness or discomfort
- Peri-implantitis, peri-mucositis or otherwise poor peri-implant soft tissue health
- Periodontal inflammation
- Phonetic difficulties
- Prosthetic failure
- Soft tissue irritation
- Sub-optimal aesthetic result
- Wound dehiscence or poor healing

Precaution: maintaining sterility protocol

The Straumann® ZAGA® Screw Retained Abutments are packed in a peel pouch or blister base with a “peel-back” lid. Labelling information is located on the bottom half of the pouch, inside the packet or on the surface of the peel-back lid. Sterility is assured unless the pouch or blister is damaged or opened.

Breakage

Implant and abutment fractures can occur when applied loads exceed the tensile or compressive strength of the material. Potential overloading conditions may result from: deficiencies in implant numbers, lengths and/or diameters to adequately support a restoration, excessive cantilever length, incomplete abutment seating, abutment angles greater than 30 degrees, occlusal interferences causing excessive lateral forces, patient parafunction (e.g. bruxing, clenching), loss or changes in dentition or functionality, inadequate prosthesis fit and physical trauma. Additional treatment may be necessary when any of the above conditions are present to reduce the possibility of hardware complications or failure.

Changes in performance

It is the responsibility of the clinician to instruct the patient on all appropriate contraindications, side effects and precautions as well as the need to seek the services of a trained dental professional if there are any changes in the performance of the implant (e.g. looseness of the prosthesis, infection or exudate around the implant, pain or any other unusual symptoms that the patient has not been told to expect).

Healing

The healing time required for osseointegration depends on the individual and treatment protocol. It is the responsibility of the practitioner to decide when the implant can be restored. Good primary stability will govern if immediate loading can be done.

Implant care and maintenance

Potential implant patients should establish an adequate oral hygiene regime prior to implant therapy. Proper post-operative, oral hygiene and implant maintenance instructions must be discussed with the patient, as this will determine the longevity and health of the implants. The patient should maintain regular prophylaxis and evaluation appointments.

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established.

The contact information for the manufacturer of this device to report a serious incident is as follows:

sicomplaints@southernimplants.com.

Materials

Material Type Commercially pure titanium (grade 4) according to ASTM F67 and ISO 5832-2 (UTS $\geq$ 900MPa)

Chemical Components	Residuals* (Fe, O, C, N, H) in total	Titanium (Ti)
Composition, % (mass/mass)	<1.05	Balance

Titanium alloy (Ti-6Al-4V) according to ASTM F136 and ISO 5832-3 (UTS $\geq$ 1015MPa)

Chemical Components	Aluminium (Al)	Vanadium (V)	Residuals* (Fe, O, C, N, H) in total	Titanium (Ti)
Composition, % (mass/mass)	5.50 – 6.75	3.50 – 4.50	<0.55	Balance

*Where, Fe = Iron, O = Oxygen, C= Carbon, N = Nitrogen, H = Hydrogen

Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

Magnetic Resonance (MR) Safety

Nonclinical testing has demonstrated that the Southern Implants® dental implants, metallic abutments and prosthetic screws are MR conditional.

MRI Safety Information

MR Conditional

A patient with the dental implants, metallic abutments and prosthetic screws from Southern Implants® may be safely scanned under the following conditions. Failure to follow these conditions may result in injury to the patient.

Device Name	dental implants, metallic abutments and prosthetic screws
Static Magnetic Field Strength (B₀) [T]	1.5T or 3T
Maximum Spatial Field Gradient [T/m] [gauss/cm]	30 T/m (3000 gauss/cm)
Type of RF Excitation	Circularly polarized (CP) (i.e., quadrature driven)
RF Transmit Coil Type	Integrated Whole Body Transmit Coil
RF Receive Coil	Any receive-only RF coil may be used
Operating Mode	Normal Operating Mode
RF Conditions	Normal Operating mode: (Head SAR of 3.2 W/kg, 2 W/kg whole body)
Scan Duration	If the closest part of implant <25cm away from the isocenter in feet-head direction: a continuous scan time of 15 minutes will require a cooling delay of at least 5 minutes. If the closest part of implant $\geq$ 25cm away from the isocenter in feet-head direction: This is limited to 1 hour of scan duration.
MR Image Artifact	The presence of this implant may produce an image artifact.

Additional information:

1. Local transmit coil: Local transmit coil safety was not evaluated. Local transmit coil is only allowed if the implant is $\geq 25\text{cm}$ away from isocenter and for a wbSAR of up to 2W/kg .
2. Artifacts: MR Image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the implant. Tests following guidelines in ASTM F2119-24 showed artifact widths of $\leq 20\text{ mm}$ from the surface of the implant for the gradient echo sequence and $\leq 10\text{ mm}$ from the surface of the implant for the spin echo sequence.

Removable restorations should be taken out prior to scanning, as is done for watches, jewellery etc.

Should there be no MR symbol on the product label, please note that this device has not been evaluated for safety and compatibility in the MR environment.

Summary of Safety and Clinical Performance (SSCP)

As required by the European Medical Device Regulation (MDR; EU2017/745), a Summary of Safety and Clinical Performance (SSCP) is available for perusal with regard to Southern Implants® product ranges.

The relevant SSCP can be accessed at <https://ec.europa.eu/tools/eudamed>.

NOTE: the above website will be available upon the launch of the European Database on Medical Devices (EUDAMED).

Compatible Components

The following components are compatible with the ZAGA® screw-retained abutments:

Table D

Product Number	Component Description
024.4323-04	Protective cap for Screw Retained Abutment, Ø4.6 mm, H 5.1mm, PEEK/TAN
024.4324-04	Protective cap for Screw Retained Abutment, Ø4.6 mm, H 6.6mm, PEEK/TAN
024.4325-04	Protective cap for Screw Retained Abutment, Ø4.6mm, H 8.1mm, PEEK/TAN
024.0020-04	Protective cap for Screw Retained Abutment, Ø4.6mm, H4.5mm, PEEK/TAN
024.0024	Temporary Coping for Screw Retained Abutments, Ø4.6mm, H 11.5mm, TAN
023.4751	Coping for Screw Retained Abutment, Bridge, Ø4.6mm, H 11mm, Ti
023.4754	Coping for Screw Retained Abutment, Bridge, Ø4.6mm, Ceramicor®/POM
023.4758	Burn-out Coping for Screw Retained Abutments, Ø4.6mm, H 11mm, POM
023.0028	Variobase® for Bridge/Bar Cylindrical Coping for Screw Retained Abutments, Ø4.6mm, H4mm, TAN
023.0032	Burn-out Coping for Variobase®, Ø4.6mm, H 11.1mm, POM
023.4752	Coping for Screw Retained Abutments, Ø4.6mm, H 5.5mm, Ti
023.4755	Coping for Screw Retained Abutments, Ø4.6mm, H5.5mm, Ceramicor®
023.4763	Occlusal screw for Titanium, Gold, Burn-out and Variobase® Copings, L 3.7mm, TAN

Additional information for these compatible components can be retrieved from <http://ifu.straumann.com/> and searching for the product number.

Straumann® Titanium, gold, burn-out, Variobase® and temporary copings are premanufactured and are available in non-engaging versions. These copings are to be fitted indirectly onto the Straumann® Zygomatic Implants via the screw-retained abutments for Zygomatic implants.

The titanium, gold, burn-out, Variobase® and temporary copings (023.4751, 023.4752, 023.4754, 023.4755, 023.4758, 023.0032, 023.0028 and 024.0024) are used in the prosthetic restoration of dental implants. The intended use of these copings is to serve as a base for the construction of bridge and bar constructs. Depending upon the clinician's

choice of prosthesis, the dental laboratory would design the final restoration according to the occlusal geometry of the abutment.

The titanium, gold, burn-out, Variobase® and temporary copings are available in a straight (0°) configuration only. Final copings may be placed into occlusion for implants with sufficient primary stability or for implants that are fully osseointegrated. The final bar or bridge restoration is screw-retained to the abutment(s).

Consult the individual product instructions for use documentation for the restorative procedures (visit <http://ifu.straumann.com/>).

NOTE: some of the compatible components for the Screw Retained Abutments may not be cleared or released for sale in all markets.

Disclaimer of liability

This product is part of the Southern Implants® product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants® product range and take full responsibility for the correct indications and use of this product. Southern Implants® does not assume liability for damage due to incorrect use. Please note that some Southern Implants® products may not be cleared or released for sale in all markets.


Basic UDI

Product	Basic-UDI Number
Basic-UDI for Metal Abutments	60095440387296
Basic-UDI for ZAGA® Screw Retained Abutments for Zygomatic Implants	6009544050137U

Related literature and catalogues

CAT-8049-STR - Straumann® ZAGA™ Zygomatic Implants
 CAT-8048-STR - Straumann® ZAGA™ Drills & Handpiece Devices
 CAT-8051-STR - Straumann® ZAGA™ Cover Screws
 CAT-8080-STR - Straumann® ZAGA™ Reusable Instruments
 CAT-8082-STR - Straumann® ZAGA™ Instrument Trays
 CAT-8083-STR - Straumann® ZAGA™ General Use Instruments

Symbols and warnings

* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

All rights reserved. Southern Implants®, the Southern Implants® logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants®. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale. It is the responsibility of the clinician to inspect the symbols that appear on the packaging of the product in use.

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