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South Africa - Headquarters: 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

Subsidiaries

United Kingdom and Ireland
Southern Implants UK
T: +44-20 8059 4490
E: info@southernimplants.co.uk

USA and Canada
Southern Implants North America Inc.
T: +1-561-472-0990
E: customercare@southernimplants.com

EC REP

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

CH REP

MedEnvoy Switzerland: Gotthardstrasse 28, 6302 Zug, Switzerland

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Description

The Drills and Handpiece devices feature either an ISO 1797 or W&H compliant latch dimension, allowing compatibility and connection to industry-standard handpieces of implant motor units. This design enables the application of torque and the execution of their respective functions through rotational motion for cutting, tightening screws, or placing implants, as applicable. The drills may be utilized with irrigation at the discretion of the operator.

Twist Drills

The cylindrical twist drills for the Straumann® Zygomatic Implant System are made of stainless steel or Grade 5 Titanium alloy. The drills have a latch dimension compatible to ISO 1797. This is in order to connect the drill to the handpiece of an implant motor unit. Drill sizes are identified by laser markings on the shaft of the drill. The drills have laser markings to indicate the drill depth. Twist drills indicated with 10 number of uses, can be used up to 10 times or when the cutting efficiency deteriorates. These drills are provided sterile.

Round Burr Drills

Round burr drills are intended to place an initial channel into bone before using the initial twist drill. The drills are for single use and supplied sterile. They are made from stainless steel.

Canal Mills

The ZAGA Side Cut Burr is a drill dedicated to the Straumann® Zygomatic Implant System. This drill has a side-cutting action as opposed to a cutting tip to aid in this function. This drill is made from Grade 5 Titanium alloy and is anodised yellow. The ZAGA Side Cut Burr is a single-use drill. The drills are provided sterile and for single patient use only.

Insertion Tools

The implant insertion tools are intended to hold a dental implant (e.g., orthodontic implant screw) at its distal end during the application of the implant to the patient's oral cavity. It is typically a small metal shaft designed to be mounted onto the distal end of a manual screwdriver or a handpiece latch (ISO 1797 or W&H) whereby rotation is applied to introduce the implant into a surgical cavity. The insertion tools are provided non-sterile.

Screw Remover Drills

Screw remover drills are intended to be used to remove broken screws from an implanted dental implant. The screw removers are intended for multiple use and supplied non-sterile. The screw remover has cutting edges or cutting threads that cut with a counter clockwise rotation when rotated by a dental handpiece. These cutting edges or threads transfer the counter clockwise torque to the broken screw and unscrew the thread of the broken screw.

Intended use

The Southern Implants® Drills and Handpiece devices are intended to prepare the implant site/osteotomy and aid in the placement/removal of the endosseous implant and/or prosthetic device(s).

Specifically, The Southern Implants® Twist Drills are intended to initiate an osteotomy and/or enlarge the channel created by a preceding drill.

The Southern Implants® Round Burr Drills are intended to initiate the osteotomy and/or create the initial channel in the bone.

The Southern Implants® Canal Mills are intended to remove excess bone from the alveolar ridge and enable canal creation.

The Southern Implants® Insertion Tool is intended to transfer torque from a handpiece device to a dental implant.

The Southern Implants® Screw Remover is intended to apply or transmit counterclockwise torque from a dental handpiece to a prosthetic screw.

Indications for use

The Southern Implants® Drills and Handpiece Devices are indicated for the functional and aesthetic treatment of patients who are eligible for implant-retained prostheses.

Specifically, the Southern Implants® Screw Remover is indicated in the removal of broken or defective prosthetic screws from an implanted endosseous implant.

Intended user

Maxillofacial Surgeons, General Dentists, Orthodontists, Periodontists, Prosthodontists, and other appropriately trained and/or experienced medical professionals.

Intended environment

The devices are intended to be used in a clinical environment such as an operating theatre or a dentist consultation room.

Intended patient population

The Southern Implants® Drills and Handpiece devices are intended to be used in patients subject to implant treatment.

Compatibility information

The devices are as described in Table A. Details on site preparation sequences recommended by Southern Implants® can be found in the relevant product catalogues.

Table A

Drill Code	Material	Description of Product	Coating (if any)	Number of Uses	Recommended Drill Speed	How Provided
TWIST DRILLS						
D-ZYG-27	Titanium grade 5 (ASTM F136)	Twist Drill Zygomatic Ø2.7	Anodized Yellow	1	1000 – 1500 rpm	Sterile
D-ZYG-27S	Titanium grade 5 (ASTM F136)	Twist Drill Zygomatic Ø2.7, Short	Anodized Yellow	1	1000 – 1500 rpm	Sterile
D-ZYG-29	Titanium grade 5 (ASTM F136)	Twist Drill Zygomatic Ø2.9	Anodized Yellow	1	1000 – 1500 rpm	Sterile
D-ZYG-29S	Titanium grade 5 (ASTM F136)	Twist Drill Zygomatic Ø2.9, Short	Anodized Yellow	1	1000 – 1500 rpm	Sterile
D-ZYG-CH-27	Stainless Steel (DIN 1.4197)	Twist Drill Zygomatic Ø2.7	-	Up to 10	1000 – 1500 rpm	Sterile
D-ZYG-CH-27S	Stainless Steel (DIN 1.4197)	Twist Drill Zygomatic Ø2.7, Short	-	Up to 10	1000 – 1500 rpm	Sterile
D-ZYG-CH-29	Stainless Steel (DIN 1.4197)	Twist Drill Zygomatic Ø2.9	-	Up to 10	1000 – 1500 rpm	Sterile
D-ZYG-CH-29S	Stainless Steel (DIN 1.4197)	Twist Drill Zygomatic Ø2.9, Short	-	Up to 10	1000 – 1500 rpm	Sterile
D-35T-M15	Stainless Steel (DIN 1.4197)	Twist Drill Ø3.5x15	-	1	1000 – 1500 rpm	Sterile
ROUND BURR DRILLS						
D-ZYG-RB	Stainless Steel (DIN 1.4197)	Round Burr Zygomatic Ø2.9	-	1	1000 – 1500 rpm	Sterile
CANAL MILLS						
CH-D-CM	Titanium grade 5 (ASTM F136)	Side Cut Burr Ø2.8-Ø3.9	Anodised Yellow	1	800 rpm	Sterile
INSERTION TOOLS (HANDPIECE)						
I-CON-X	Stainless Steel (DIN 1.4197)	W&H Hexed Latch Handpiece	-	Indefinite (provided the device is not damaged)	15 – 20 rpm	Non-Sterile

SCREW REMOVAL DRILLS						
I-SR-6	Stainless Steel (DIN 1.4197)	Screw Remover Drill Ø1.15	-	Up to 10	600 – 800 rpm (in reverse)	Non-Sterile

Surgical procedures for Twist Drills

Ensure that all instruments and drills are in a good condition. The drill sizes are identified by different laser markings on the shaft of the drill. The drills have different laser markings on the body of the drill which corresponds to the length of implant being placed. Refer to individual product catalogues for more information.

1. Drill in the planned direction to the full depth of the implant length being placed as indicated on the markings on the drill. Cylindrical drills extend up to 1mm longer than the implant, when seated. Allow for this additional length when drilling near vital anatomical structures.

NOTE: connect the drill latch to the handpiece. If the latch is not inserted fully into the handpiece, the torque is applied to the latch, resulting in possible twisting of the latch or damage to the handpiece. Consult the instructions for use of the handpiece to ensure proper engagement of the latch.

2. Drill at sufficient speed (between 1000-1500rpm for cylindrical drills), under constant irrigation with sterile saline.
3. Use an up-and-down motion with the handpiece, without stopping the motor. This will allow the irrigation to flush away bone debris on the drill.
4. Insert the direction indicator (ZYG-TR-55-35/ 45/ 52) to verify the alignment with adjacent implants or teeth. A radiograph can be taken at this stage to verify depth and direction.
5. The depth can be determined by a depth gauge: I-ZYG-DG-1 or CH-I-DG.
6. During surgery the clinician will be able to assess the bone quality and should adjust drill protocols, when necessary, to prepare the site. This is to avoid the implant getting stuck before it is properly seated in the osteotomy. It is recommended to undersize the osteotomy in soft bone.
7. Preparing the site further should involve: ensuring the drill goes to full depth.

NOTE:

- avoid lateral pressure (bending) on the drills during drilling procedures.
- lateral pressure to the drill can cause drill fracture.
- verify the drill is securely locked into the handpiece before the drilling procedure starts.
- For the twist drills which can be re-used, it is recommended to maintain a log of these drills, recording the number of uses. Prior to re-processing these components, it should be thoroughly inspected and tested to determine its suitability for re-use.
- do not apply more than 40-45Ncm of torque to any latch type drill/instrument, this could cause damage to the handpiece and latch of the instrument.
- blunt drills cause excessive torque and result in damage to the handpiece or drill latch.

Surgical procedures for Canal Mills

Case Planning

A proper clinical and radiological evaluation must be done to determine the bone dimensions and bone quality.

Procedure

1. Verify the drill is securely locked into the handpiece before drilling procedure starts.
2. Drill at sufficient speed (800rpm), with constant irrigation with sterile saline.
3. Use an up-and-down motion with the hand-piece, without stopping the motor. This will allow the irrigation to flush away bone debris on the drill.
4. Gradually enlarge the osteotomy in a stepwise approach to the desired diameter and depth.
5. During surgery the clinician will be able to assess the bone quality and should use dense bone protocols, when necessary, to prepare the site. This is to avoid the implant getting stuck before it is properly seated in the osteotomy.
6. Avoid excessive lateral pressure (bending) on the drills during drilling procedures. Excessive lateral pressure to the drill can cause drill fracture.

NOTE:

- connect the drill latch to the handpiece. If the latch is not inserted fully into the handpiece, the torque is applied to the latch, resulting in possible twisting of the latch or damage to the handpiece. Consult the instructions for use of the handpiece to ensure proper engagement of the latch.
- do not apply more than 40-45Ncm to any latch type drill/ instrument, this could cause damage to the handpiece and latch of the instrument.
- blunt drills cause excessive torque and result in damage to the handpiece or drill latch.
- drill extension section is removed.

Surgical procedures for Screw Removal Drills**Step 1**

The use of a fine probe to gently rotate the broken screw in an anti-clockwise direction, may often turn the screw out, to the point that it may be gripped with a haemostat and rotated out fully. The screw is usually quite loose once it has fractured, and should be freely able to rotate.

Step 2

Should step 1 not be successful and the screw is not loose, the use of an ultrasonic scaler in the same manner as the probe, should be applied. Often, the ultrasound assists in the rotation/freeing up of the screw. Once the screw has been exposed to ultrasound, Step 1 should be applied a second time.

Step 3

If the ultrasound is not successful on its own. The correct drill guide must be selected for the corresponding implant. Place the guide on top of the implant interface, it must be ensured that the drill guide is correctly seated on the implant interface. Each drill guide has a square connection on the opposing end, allowing them to be stabilized using a ratchet or torque wrench. The screw removal drills rotate in a counter-clockwise direction. The first drill, (I-SR-6), cuts a pilot hole in the centre of the screw to be removed. It is possible that the drilling process of the I-SR-6 may thread the screw out of the implant, if the drill grips the screw while cutting the pilot hole.

NOTE:

- The screw removal instrumentation may only be used with a high torque – low speed motor unit, in reverse set at 600 – 800 rpm for the I-SR-6 with copious irrigation. Do not overheat the implant via friction this can lead to failure of the implant.
- Gently insert the drill into the latch-grip handpiece, ensuring that it is fully engaged. If the latch is not properly engaged in the handpiece, the torque applied when activating the motor unit may result in the possible distortion of the latch or damage to the handpiece, due to the incorrect distribution of forces within the latch-grip mechanism.
- Consult the instructions for use of the relevant handpiece manufacturer, to ensure proper engagement of the latch.
- Do not apply more than 40-45 Ncm to any latch type instrument. This could cause damage to the handpiece gearing mechanism, and/or latch of the instrument.

Clinical performance

Southern Implants® drills are designed to aid in implant site preparation by initiating and shaping a precisely profiled osteotomy for implant placement. Their specialised geometry, including appropriate lengths, diameters, and cutting features, ensures accurate and efficient site preparation.

The Southern Implants® handpiece devices are designed to assist in the placement or removal of endosseous implants and prosthetic components. The Southern Implants® handpiece devices are designed to enable secure attachment and controlled torque transfer from the handpiece to the implant or abutment during placement.

Clinical benefits

Southern Implants® drills play a vital role in implant-retained prosthetic treatment by enabling precise osteotomy preparation. Southern Implants® handpiece devices play a vital role in implant-retained prosthetic treatment by enabling secure attachment and controlled torque transfer to the implant or abutment during placement.

While their clinical benefits are primarily tied to the overall success of the restorative treatment rather than their direct performance, their role in achieving accurate implant placement supports favourable outcomes.

The primary clinical benefits of this treatment approach include the effective restoration of the identified defect and enhanced patient satisfaction. Additionally, successful prosthetic rehabilitation can contribute to improved social interaction, psychosocial well-being, self-esteem, and confidence. Although these benefits are indirect and challenging to quantify, they are closely linked to positive treatment outcomes, such as optimal implant placement, functional restoration, and patient-reported satisfaction with both aesthetics and function.

Storage, cleaning and sterilisation

Applicable For Sterile Components:

The component is supplied sterile (sterilised by gamma irradiation) and may be intended for single-use prior to the expiration date (see packaging label). Sterility is assured unless the container or seal is damaged or opened. If packaging is damaged do not use the product and contact your Southern representative or return to Southern Implants®. The devices must be stored in a dry place at room temperature and not exposed to direct sunlight. Incorrect storage may influence device characteristics. Do not reuse components indicated for single-use only. Reusing these components may result in:

- damage to the surface or critical dimensions, which may result in performance and compatibility degradation.
- adds the risk of cross-patient infection and contamination if single-use items are reused.

Southern Implants® does not accept any responsibility for complications associated with reused single-use components.

Applicable For Reusable Components:

Limitations on reusable items

A direct value for reusable instruments cannot be given. Frequent processing may have minor effects on the instruments. The product life is normally determined by wear and damage during use, thus instruments if properly cared for and inspected after each use, can be reused many times. Maintain a checklist for these instruments recording the number of uses.

Prior to reprocessing the device, it should be thoroughly inspected and tested to determine its suitability for reuse.

NOTE: during use, handle drills and instruments with sterile tweezers to minimise contamination of the instrument tray and risk of damage to sterile surgical gloves.

Containment

As soon as practically possible, remove all visible residue after use (bone, blood or tissue), by immersing the instrument in cold water (dried soil may be difficult to remove).

Precleaning

Disassemble instruments from hand pieces and all connecting parts from instruments in order to clean soil from obstructed areas. Remove the PEEK bits from placement tools. Rinse with lukewarm water for 3 minutes and remove hardened debris with a soft nylon brush. Avoid mechanical damage during cleaning.

Manual cleaning or automated cleaning

Prepare an ultrasonic bath with suitable detergent (i.e. Steritech instrument cleaner - 5% dilution), sonicate for 20 minutes (alternative methods can be used if proven by the end user). Rinse with purified/sterile water.

NOTE: always follow the instructions for use of the manufacturers for cleaning agents and disinfectants.

Load devices into a thermodisinfector. Run the cleaning and disinfection cycle, followed by the drying cycle.

Drying

Dry the instruments both inside and outside with filtered compressed air or single use, lint free wipes. Pack the instruments as quickly as possible after removal into the storage container. If additional drying is necessary, dry in a clean location.

Inspection

Do a visual inspection of the items to check for damages on instruments.

Packaging

Use the correct packaging material as indicated for steam sterilisation to ensure sterility is maintained. Double packaging is recommended. Where appropriate, the cleaned, disinfected and checked devices can be assembled and placed in instrument trays as applicable. Instrument trays can be double wrapped or placed in sterilisation bags.

Sterilisation

Southern Implants® recommends one of the following procedures to sterilise the device prior to use:

1. prevacuum sterilisation method: steam sterilise the device at 132°C (270°F) at 180 - 220 kPa for 4 minutes. Dry for at least 20 minutes in the chamber. Only an approved wrap or pouch for steam sterilisation must be used.
2. for users in the USA: prevacuum sterilisation method: wrapped, steam sterilise at 135°C (275°F) at 180 - 220 kPa for 3 minutes. Dry for 20 minutes in the chamber. Use a wrap or pouch that is cleared for the indicated steam sterilisation cycle.

NOTE: users in the USA must ensure that the steriliser, wrap or pouch, and all steriliser accessories are cleared by the FDA, for the intended sterilisation cycle.

Storage

Maintain packaging integrity to ensure sterility in storage. Packaging should be completely dry before storage to avoid corrosion and degradation of cutting edges.

Contraindications

The contraindications of all device groups used as part of the specific treatment or procedure apply. Therefore, the contraindications of the systems/medical devices utilized as part of implant surgery/therapy should be noted and the relevant documents consulted.

Contraindications specific to the Southern Implants® Drills and Handpiece Devices:

- Allergies or hypersensitivity to chemical ingredients in the following materials: Titanium, Aluminium, Vanadium, Stainless Steel.

Warnings and precautions

IMPORTANT NOTICE: THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING

- To ensure the safe and effective use of dental implants, new technologies/systems and the subject devices, it is strongly recommended that specialised training be undertaken. This training should include hands-on methods to gain competency on the proper technique, the system's biomechanical requirements and radiographic evaluations required for the specific system.
- Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone.
- Use of the device with incompatible or non-corresponding devices can result in the poor performance or failure of the device.
- When handling devices intraorally, it is imperative that they are adequately secured to prevent aspiration, as aspiration of products may lead to infection or physical injury.
- The use of non-sterile items can lead to secondary infections of the tissue or transfer infectious diseases.

- Failure to adhere to appropriate cleaning, re-sterilization, and storage procedures as outlined in the Instructions for Use can result in device damage, secondary infections, or patient harm.
- Exceeding the number of recommended uses for reusable devices can result in device damage, secondary infection or patient harm.
- The use of blunt drills may cause damage to the bone, potentially compromising osseointegration.

It is crucial to emphasize that training should be undertaken by both new and experienced implant users prior to utilizing a new system or attempting a new treatment method.

Patient Selection and Preoperative Planning

A comprehensive patient selection process and meticulous preoperative planning are essential for successful implant treatment. This process should involve consultation among a multidisciplinary team, including well-trained surgeons, restorative dentists, and laboratory technicians.

Patient screening should include, at a minimum, a thorough medical and dental history, as well as visual and radiological inspections to assess the presence of adequate bone dimensions, the positioning of anatomical landmarks, the presence of unfavourable occlusal conditions, and the periodontal health status of the patient.

For successful implant treatment, it is important to:

1. Minimize trauma to the host tissue, thereby increasing the potential for successful osseointegration.
2. Accurately identify measurements relative to radiographic data, as failure to do so may lead to complications.
3. Be vigilant in avoiding damage to vital anatomical structures, such as nerves, veins, and arteries. Injury to these structures may result in serious complications, including ocular injury, nerve damage, and excessive bleeding.

The responsibility for proper patient selection, adequate training and experience in implant placement, and the provision of appropriate information required for informed consent rests with the practitioner. By combining thorough screening of prospective implant candidates with a practitioner possessing a high level of competence in the use of the system, the potential for complications and severe side effects can be significantly reduced.

High Risk Patients

Special care should be exercised when treating patients with local or systemic risk factors that may adversely affect the healing of bone and soft tissue or otherwise increase the severity of side effects, the risk of complications, and/or the likelihood of implant failure. Such factors include:

- poor oral hygiene
- history of smoking/vaping/tobacco use
- history of periodontal disease
- history of orofacial radiotherapy**
- bruxism and unfavourable jaw relations
- use of chronic medications that may delay healing or increase the risk of complications including, but not limited to, chronic steroid therapy, anti-coagulant therapy, TNF- α blockers, bisphosphonate and cyclosporin

*** The potential for implant failure and other complications increases when implants are placed in irradiated bone, as radiotherapy can lead to progressive fibrosis of blood vessels and soft tissue (i.e., osteoradionecrosis), resulting in diminished healing capacity. Contributing factors to this increased risk include the timing of implant placement in relation to radiation therapy, the proximity of radiation exposure to the implant site, and the radiation dosage at that site.*

Should the device not operate as intended, it must be reported to the manufacturer of the device. The contact information or the manufacturer of this device to report a change in performance is: sicomplaints@southernimplants.com.

Side effects

The clinical outcome of treatment is influenced by multiple factors. The following side effects and residual risks are associated with this device group or the surgical procedures involved in its use. These may require further treatment, revision surgery, or additional consultations with the relevant medical professional.

- Allergic or hypersensitivity reaction(s)
- Anaesthesia, paraesthesia, hyperesthesia, and hypoesthesia (transient or permanent)
- Anatomical landmark damage
- Bruising
- Dental injury during surgery
- Fat emboli formation
- Fistula formation
- Gingival injury
- Infection (acute and/or chronic)
- Localized inflammation
- Loss or damage to adjacent dentition
- Mild thermal injury to adjacent structures (e.g. burn to the lip)
- Nasal bleeding (Zygomatic Implant procedures)
- Nerve damage
- Pain, discomfort or tenderness
- Perforation of the labial or lingual plates
- Peri-implantitis, peri-mucositis or otherwise poor peri-implant soft tissue health
- Sinus Perforation
- Soft tissue irritation
- Surgical side effects such as pain, inflammation, bruising and mild bleeding
- Thermal injury to bone resulting in insufficient levels of osseointegration
- Wound dehiscence or poor healing

Additionally, the normal side effects associated with anaesthesia should also be expected.

Note: these side effects and residual risks can vary in both severity and frequency.

Precaution: maintaining sterility protocol

Applicable To Sterile Components

The devices are packed in a peel pouch or blister base with a “peel-back” lid. Labelling information is located on the bottom half of the pouch, inside the packet or on the surface of the peel-back lid. Sterility is assured unless the pouch or blister is damaged or opened.

Applicable To Non-Sterile Components

The devices are supplied clean but not sterile in a peel pouch or blister base with peel-back lid. Labelling information is located on the bottom half of the pouch or on the surface of the peel-back lid.

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established.

The contact information for the manufacturer of this device to report a serious incident is as follows:

sicomplaints@southernimplants.com.

Materials

Material type: Titanium alloy Ti-6Al-4V (90% titanium, 6% aluminium, 4% vanadium) according to ASTM F136 and ISO 5832-3 or Stainless Steel according to DIN 1.4197/DIN 1.4542 and ASTM F899-23

Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

Summary of Safety and Clinical Performance (SSCP)

As required by the European Medical Device Regulation (MDR; EU2017/745), a Summary of Safety and Clinical Performance (SSCP) is available for perusal with regard to Southern Implants® product ranges.

The relevant SSCP can be accessed at <https://ec.europa.eu/tools/eudamed>.

NOTE: the above website will be available upon the launch of the European Database on Medical Devices (EUDAMED).

Disclaimer of liability

This product is part of the Southern Implants® product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants® product range and take full responsibility for the correct indications and use of this product. Southern Implants® does not assume liability for damage due to incorrect use. Please note that some Southern Implants® products may not be cleared or released for sale in all markets.

Basic UDI

Product	Basic-UDI Number
Basic-UDI for Drills and Handpiece Devices	6009544038759C

Related literature and catalogues

CAT-8049-STR - IFU - Straumann® ZAGA™ Zygomatic Implants

CAT-8047-STR - IFU - Straumann® ZAGA™ Screw-Retained Abutments

CAT-8051-STR - IFU - Straumann® ZAGA™ Cover Screws

CAT-8082-STR - IFU - Straumann® ZAGA™ Instrument Trays

CAT-8080-STR - IFU - Straumann® ZAGA™ Reusable Instruments

Symbols and warnings

Manufacturer:
Southern
Implants®
1 Albert Rd, P.O
Box 605 IRENE,
0062,
South Africa.
Tel: +27 12 667
1046

CE mark

Prescription
device*

Sterilised using
irradiation

Non-sterile

Use by date
(mm-yy)

Do not
reuse

Do not
resterilize

Catalogue
number

Batch code

Medical
device

Authorised
representative in
the European
Community

Authorised
representative for
Switzerland

Date of
manufacture

Magnetic
Resonance
conditional

Magnetic
Resonance
safe

Single sterile barrier system
with protective packaging
inside

Single sterile
barrier system

Consult
instruction for
use

Caution

Keep away
from
sunlight

Do not use if
package is
damaged

* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

All rights reserved. Southern Implants®, the Southern Implants® logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants®. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale. It is the responsibility of the clinician to inspect the symbols that appear on the packaging of the product in use.

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