English	INSTRUCTIONS FOR USE: Straumann® ZAGA™ Cover Screws
Español	INSTRUCCIONES DE USO: Tornillos de cobertura ZAGA™ Straumann®
Italiano	ISTRUZIONI PER L'USO: Straumann® ZAGA™ Viti di Copertura
Français	MODE D'EMPLOI : Vis de couverture Straumann® ZAGA™
Deutsch	GEBRAUCHSANWEISUNGEN: Straumann® ZAGA™ Abdeckschrauben
Português	INSTRUÇÕES DE UTILIZAÇÃO: Parafusos de cobertura Straumann® ZAGA™
Ελληνικά	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ: Βίδες κάλυψης ZAGA™ της Straumann®
Svenska	ANVISNINGAR FÖR ANVÄNDNING: Straumann® ZAGA™ Täckskravar
Brazilian	INSTRUÇÕES DE USO: Straumann® ZAGA™ Parafusos da tampa
Hrvatski	UPUTE ZA UPOTREBU: Straumann® ZAGA™ vijci
Dansk	BRUGSANVISNING: Straumann® ZAGA™ dækskruer
Suomi	KÄYTTÖOHJEET: Straumann® ZAGA™ -suojaruuvit
Nederlands (NL)	GEBRUIKSAANWIJZING: Straumann® ZAGA™ afdekschroeven
Magyar	HASZNÁLATI ÚTMUTATÓ: Straumann® ZAGA™ fedőcsavarok
Lietuviškai	NAUDOJIMO NURODYMAI: Straumann® ZAGA™ dengiantieji varžtai
Polski	INSTRUKCJA UŻYTKOWANIA: Straumann® ZAGA™ śruby zamykające
Română	INSTRUCȚIUNI DE UTILIZARE: Straumann® ZAGA™ Șuruburi de acoperire
Slovenščina	NAVODILA ZA UPORABO: Zaporni vijaki Straumann® ZAGA™
Українська	ІНСТРУКЦІЇ З ЕКСПЛУАТАЦІЇ: Заглушки Straumann® ZAGA™
Русский	ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ: Заглушки Straumann® ZAGA™
Български	ИНСТРУКЦИИ ЗА УПОТРЕБА: Покривни винтове Straumann® ZAGA™
Česky	NÁVOD K POUŽITÍ: Krycí šrouby Straumann® ZAGA™
Eesti	KASUTUSJUHEND: Straumann® ZAGA™ kattekruvid
Íslenska	NOTKUNARLEIÐBEININGAR: Straumann® ZAGA™ hlífarskrúfur
Latviski	LIETOŠANAS INSTRUKCIJA: Straumann® ZAGA™ vāka skrūves
Norsk	BRUKSANVISNINGER: Straumann® ZAGA™ dekskruer
Slovenčina	NÁVOD NA POUŽITIE: Krycie skrutky Straumann® ZAGA™
Türkçe	KULLANMA TALİMATI: Straumann® ZAGA™ Kapak Vidaları
Српски	УПУТСТВО ЗА УПОТРЕБУ: Straumann® ZAGA™ завртњи поклопца

South Africa - Headquarters: 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

Subsidiaries

United Kingdom and Ireland
Southern Implants UK
T: +44-20 8059 4490
E: info@southernimplants.co.uk

USA and Canada
Southern Implants North America Inc.
T: +1-561-472-0990
E: customercare@southernimplants.com

EC REP

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

CH REP

MedEnvoy Switzerland: Gotthardstrasse 28, 6302 Zug, Switzerland

Description

Cover screws are pre-manufactured temporary dental implant components used during the initial healing phase following implant insertion. They connect to the implant and seal the connection of the implant. The cover screw has a low profile so that soft tissue can be sutured closed over it. Cover screws have built in screws for retention to the implant and are made of commercially pure titanium. All Southern Implants cover screws are provided sterile.

Intended use

The Southern Implants® Metal Abutments are prosthetic components intended for connection to endosseous implant systems to aid in prosthetic rehabilitation and implant procedures.

Specifically, the Straumann® ZAGA™ Cover Screws are intended to be temporarily connected to endosseous implants following their surgical placement in order to seal the implant connection and allow for osseointegration and soft tissue healing prior to functional loading.

Indications for use

The Southern Implants® Metal Abutments are indicated for use as components of an endosseous implant system to support prosthetic rehabilitation in implant treatment procedures.

The Straumann® ZAGA™ Cover Screws are indicated for 2-stage implant surgeries where immediate loading of the endosseous implant(s) is deemed inappropriate due to the limited primary stability of the implant, patient specific risk factors, treatment specific factors and/or the discretion of the user.

Intended user

Maxillofacial Surgeons, General Dentists, Orthodontists, Periodontists, Prosthodontists, and other appropriately trained and/or experienced medical professionals.

Intended environment

The devices are intended to be used in a clinical environment such as an operating theatre or a dentist consultation room.

Intended patient population

Southern Implants® Metal Abutments are intended to be used in patients indicated and eligible, or otherwise not contraindicated, for implant-retained prosthetic restorations.

Compatibility information

Use only original Southern Implants® components to restore Straumann® ZAGA™ Zygomatic ranges. Use components that correspond to the connection type, and prosthetic platform when restoring Zygomatic implants.

Table A

Straumann® Zygomatic ZAGA™ Implants				
Straumann® Zygomatic Implant ZAGA™ Flat	Straumann® Zygomatic Implant ZAGA™ Round	Cover Screw and Driver	Abutment and Driver	Prosthetic Screw and Driver
CH-ZC-30.0	CH-ZT-30.0	CH-CS (cover screw) I-CS-HD (driver)	CH-SRA-xx* (straight Screw-Retained Abutment) 046.401 046.411 (driver)	I-HD-M (driver)
CH-ZC-32.5	CH-ZT-32.5			
CH-ZC-35.0	CH-ZT-35.0			
CH-ZC-37.5	CH-ZT-37.5			
CH-ZC-40.0	CH-ZT-40.0			
CH-ZC-42.5	CH-ZT-42.5			
CH-ZC-45.0	CH-ZT-45.0			
CH-ZC-47.5	CH-ZT-47.5			
CH-ZC-50.0	CH-ZT-50.0			
CH-ZC-52.5	CH-ZT-52.5			
CH-ZC-55.0	CH-ZT-55.0			
CH-ZC-57.5	CH-ZT-57.5			
CH-ZC-60.0	CH-ZT-60.0			

*xx denotes the abutment collar height

Table B

Straumann® Zygomatic ZAGA™ Implants			
Straumann® Zygomatic Implant ZAGA™ Flat	Cover Screw and Driver	Abutment and Driver	Prosthetic Screw and Driver
CH-ZF-30.0	CH-CS (cover screw) I-CS-HD (driver)	CH-SRA-xx* (straight Screw-Retained Abutment) CH-SRA17-xx* (angled Screw-Retained Abutment) 046.401 046.411 (driver)	I-HD-M (driver)
CH-ZF-32.5			
CH-ZF-35.0			
CH-ZF-37.5			
CH-ZF-40.0			
CH-ZF-42.5			
CH-ZF-45.0			
CH-ZF-47.5			
CH-ZF-50.0			
CH-ZF-52.5			
CH-ZF-55.0			
CH-ZF-57.5			
CH-ZF-60.0			

*xx denotes the abutment collar height

Surgical Procedures

1. Select the appropriate cover screw.
2. Connect the cover screw to the implant and tighten the cover screws, with the applicable driver (Table A and Table B).
3. Torque the cover screw down to the value indicated in Table C.
4. Reposition the flap margins together and suture closed.

Table C

Abutment Recommended Torque Values	
Abutment	Recommended Torque
CH-CS	Finger tighten

Before surgery

All components, instruments and tooling used during the clinical or laboratory procedure must be maintained in good condition and care must be taken that instrumentation does not damage implants or other components.

During surgery

Care must be taken that parts are not swallowed during any of the procedures, a rubber-dam application is recommended when appropriate. Care must be taken to apply the correct tightening torque of abutments and abutment screws.

Post-surgery

Regular patient follow-up and proper oral hygiene must be achieved to ensure favourable long-term results.

Clinical performance

The clinical performance of the Straumann® ZAGA™ Cover Screws is primarily assessed through their effectiveness in facilitating osseointegration and promoting optimal soft tissue healing. These temporary abutments are designed to be securely connected to endosseous implants immediately following surgical placement where their primary function is to seal the implant connection and creating a stable environment that fosters the integration of the implant into the surrounding bone and the healing of the soft tissue.

Clinical benefits

The clinical benefits associated with the Straumann® ZAGA™ Cover Screws are intrinsically linked to the overall success of implant-supported prosthetic treatments and encompass both direct and indirect clinical benefits. The direct clinical benefits of the cover screws arise specifically from their use in the treatment process. These benefits include the

enhancement of peri-implant soft tissue health and the aesthetic outcomes of the final restoration. The clinical benefits experienced by patients as a result of the implant-supported prosthetic treatments in which these abutments are utilized represent the indirect clinical benefits of the Straumann® ZAGA™ Cover Screws. While these benefits are substantial, they are not directly attributable to the devices themselves; rather, they reflect the contributions of the abutments to the overall treatment system. Additionally, successful treatment outcomes are anticipated to result in significant enhancements in patients' psychosocial well-being, self-esteem, and overall quality of life. Overall, the Straumann® ZAGA™ Cover Screws not only facilitate the mechanical aspects of dental restorations but also contribute to the holistic improvement of patient outcomes and experiences

Storage, cleaning and sterilisation

The component is supplied sterile (sterilised by gamma irradiation) and intended for single-use prior to the expiration date (see packaging label). Sterility is assured unless the container or seal is damaged or opened. If packaging is damaged do not use the product and contact your Straumann representative or return to Southern Implants®. The devices must be stored in a dry place at room temperature and not exposed to direct sunlight. Incorrect storage may influence device characteristics.

Do not reuse components indicated for single-use only. Reusing these components may result in:

- damage to the surface or critical dimensions, which may result in performance and compatibility degradation.
- adds the risk of cross-patient infection and contamination if single-use items are reused.

Southern Implants® does not accept any responsibility for complications associated with reused single-use components.

Contraindications

The contraindications of all device groups used as part of the specific treatment or procedure apply. Therefore, the contraindications of the systems/medical devices utilized as part of implant surgery/therapy should be noted and the relevant documents consulted.

The contraindications specific to this device group include:

- patients who are medically unfit for implant or prosthetic related procedures
- patients that present with an allergy or hypersensitivity to the following materials: Titanium, Aluminium, Vanadium, Cobalt Chrome, Cobalt Chrome Molybdenum, Gold, Platinum, Palladium, Iridium, Silver and/or Copper.

Warnings and precautions

IMPORTANT NOTICE: THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING

- To ensure the safe and effective use of dental implants, new technologies/systems and the Metal Abutment devices, it is strongly recommended that specialised training be undertaken. This training should include hands-on methods to gain competency on the proper technique, the system's biomechanical requirements and radiographic evaluations required for the specific system.
- Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone.
- Use of the device with incompatible or non-corresponding devices can result in the poor performance or failure of the device.
- When handling devices intraorally, it is imperative that they are adequately secured to prevent aspiration, as aspiration of products may lead to infection or physical injury.
- The use of non-sterile items can lead to secondary infections of the tissue or transfer infectious diseases.
- Failure to adhere to appropriate cleaning, re-sterilization, and storage procedures as outlined in the Instructions for Use can result in device damage, secondary infections, or patient harm.
- Exceeding the number of recommended uses for reusable devices can result in device damage, secondary infection or patient harm.
- The use of blunt drills may cause damage to the bone, potentially compromising osseointegration.

It is crucial to emphasize that training should be undertaken by both new and experienced implant users prior to utilizing a new system or attempting a new treatment method.

Patient Selection and Preoperative Planning

A comprehensive patient selection process and meticulous preoperative planning are essential for successful implant treatment. This process should involve consultation among a multidisciplinary team, including well-trained surgeons, restorative dentists, and laboratory technicians.

Patient screening should include, at a minimum, a thorough medical and dental history, as well as visual and radiological inspections to assess the presence of adequate bone dimensions, the positioning of anatomical landmarks, the presence of unfavourable occlusal conditions, and the periodontal health status of the patient.

For successful implant treatment, it is important to:

1. Minimize trauma to the host tissue, thereby increasing the potential for successful osseointegration.
2. Accurately identify measurements relative to radiographic data, as failure to do so may lead to complications.
3. Be vigilant in avoiding damage to vital anatomical structures, such as nerves, veins, and arteries. Injury to these structures may result in serious complications, including ocular injury, nerve damage, and excessive bleeding.

The responsibility for proper patient selection, adequate training and experience in implant placement, and the provision of appropriate information required for informed consent rests with the practitioner. By combining thorough screening of prospective implant candidates with a practitioner possessing a high level of competence in the use of the system, the potential for complications and severe side effects can be significantly reduced.

High Risk Patients

Special care should be exercised when treating patients with local or systemic risk factors that may adversely affect the healing of bone and soft tissue or otherwise increase the severity of side effects, the risk of complications, and/or the likelihood of implant failure. Such factors include:

- poor oral hygiene
- history of smoking/vaping/tobacco use
- history of periodontal disease
- history of orofacial radiotherapy**
- bruxism and unfavourable jaw relations
- use of chronic medications that may delay healing or increase the risk of complications including, but not limited to, chronic steroid therapy, anti-coagulant therapy, TNF- α blockers, bisphosphonate and cyclosporin

*** The potential for implant failure and other complications increases when implants are placed in irradiated bone, as radiotherapy can lead to progressive fibrosis of blood vessels and soft tissue (i.e., osteoradionecrosis), resulting in diminished healing capacity. Contributing factors to this increased risk include the timing of implant placement in relation to radiation therapy, the proximity of radiation exposure to the implant site, and the radiation dosage at that site.*

Should the device not operate as intended, it must be reported to the manufacturer of the device. The contact information or the manufacturer of this device to report a change in performance is: sicomplaints@southernimplants.com.

Side effects

The clinical outcome of treatment is influenced by various factors and there are additional potential side effects and residual risks associated with the device group. These may necessitate further treatment, revision surgery or additional visits to the relevant medical professional's office. Furthermore, these side effects and residual risks can vary in both severity and frequency.

Therefore, the full list of known potential undesirable side-effects and residual risks that have been identified in relation to the dental abutment device group include:

- Abutment fracture
- Abutment screw fracture or distortion

- Allergic reaction(s) to the abutment material
- Anaesthesia, paraesthesia, hyperesthesia, and hypoesthesia (transient or permanent)
- Bleeding
- Bleeding on probing
- Bruising
- Complications requiring revision surgery
- Difficulty in retrieving the abutment
- Gingival inflammation
- Gingival recession
- Hyperplastic tissue response
- Implant failure due to insufficient levels of osseointegration
- Improper implant positioning resulting in prosthetic compromise
- Infection (acute and/or chronic)
- Localized inflammation
- Loosening of the abutment screw and/or retaining screw
- Loss or damage to adjacent teeth
- Marginal bone loss
- Micromovements and implant instability
- Misfit or improper connection at the implant-abutment interface
- Overloading of the abutment/implant
- Pain, tenderness or discomfort
- Peri-implantitis, peri-mucositis or otherwise poor peri-implant soft tissue health
- Periodontal inflammation
- Phonetic difficulties
- Prosthetic failure
- Soft tissue irritation
- Sub-optimal aesthetic result
- Wound dehiscence or poor healing

Precaution: maintaining sterility protocol

The Straumann® ZAGA™ Cover Screws are packed in a peel pouch or blister base with a “peel-back” lid. Labelling information is located on the bottom half of the pouch, inside the packet or on the surface of the peel-back lid. Sterility is assured unless the pouch or blister is damaged or opened.

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established.

The contact information for the manufacturer of this device to report a serious incident is as follows:

sicomplaints@southernimplants.com.

Materials

Material type Commercially pure titanium (grade 4) according to ASTM F67 and ISO 5832-2 (UTS ≥ 900MPa)

Chemical Components		Residuals* (Fe, O, C, N, H) in total	Titanium (Ti)
Composition, % (mass/mass)		<1.05	Balance

*Where, Fe = Iron, O = Oxygen, C= Carbon, N = Nitrogen, H = Hydrogen

Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

Magnetic Resonance (MR) Safety

Nonclinical testing has demonstrated that the Southern Implants® dental implants, metallic abutments and prosthetic screws are MR conditional.

MRI Safety Information	
 MR Conditional	
A patient with the dental implants, metallic abutments and prosthetic screws from Southern Implants® may be safely scanned under the following conditions. Failure to follow these conditions may result in injury to the patient.	
Device Name	dental implants, metallic abutments and prosthetic screws
Static Magnetic Field Strength (B₀) [T]	1.5T or 3T
Maximum Spatial Field Gradient [T/m] [gauss/cm]	30 T/m (3000 gauss/cm)
Type of RF Excitation	Circularly polarized (CP) (i.e., quadrature driven)
RF Transmit Coil Type	Integrated Whole Body Transmit Coil
RF Receive Coil	Any receive-only RF coil may be used
Operating Mode	Normal Operating Mode
RF Conditions	Normal Operating mode: (Head SAR of 3.2 W/kg, 2 W/kg whole body)
Scan Duration	If the closest part of implant <25cm away from the isocenter in feet-head direction: a continuous scan time of 15 minutes will require a cooling delay of at least 5 minutes. If the closest part of implant ≥ 25cm away from the isocenter in feet-head direction: This is limited to 1 hour of scan duration.
MR Image Artifact	The presence of this implant may produce an image artifact.

Additional information:

1. Local transmit coil: Local transmit coil safety was not evaluated. Local transmit coil is only allowed if the implant is ≥ 25cm away from isocenter and for a wbSAR of up to 2W/kg.
2. Artifacts: MR Image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the implant. Tests following guidelines in ASTM F2119-24 showed artifact widths of ≤20 mm from the surface of the implant for the gradient echo sequence and ≤10 mm from the surface of the implant for the spin echo sequence.

Removable restorations should be taken out prior to scanning, as is done for watches, jewellery etc.

Should there be no MR symbol on the product label, please note that this device has not been evaluated for safety and compatibility in the MR environment.

Summary of Safety and Clinical Performance (SSCP)

As required by the European Medical Device Regulation (MDR; EU2017/745), a Summary of Safety and Clinical Performance (SSCP) is available for perusal with regard to Southern Implants® product ranges.

The relevant SSCP can be accessed at <https://ec.europa.eu/tools/eudamed>.

NOTE: the above website will be available upon the launch of the European Database on Medical Devices (EUDAMED).

Disclaimer of liability

This product is part of the Southern Implants® product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants® product range and take full responsibility for the correct indications and use of this product. Southern Implants® does not assume liability for damage due to incorrect use. Please note that some Southern Implants® products may not be cleared or released for sale in all markets.


Basic UDI

Product	Basic-UDI Number
Basic-UDI for Metal Abutments	60095440387296
Basic-UDI for Cover Screws	60095440500883

Related literature and catalogues

CAT-8049-STR - Straumann® ZAGA™ Zygomatic Implants
 CAT-8048-STR - Straumann® ZAGA™ Drills & Handpiece Devices
 CAT-8047-STR - Straumann® ZAGA™ Screw Retained Abutments
 CAT-8080-STR - Straumann® ZAGA™ Reusable Instruments
 CAT-8082-STR - Straumann® ZAGA™ Instrument Trays
 CAT-8083-STR - Straumann® ZAGA™ General Use Instruments

Symbols and warnings

* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

All rights reserved. Southern Implants®, the Southern Implants® logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants®. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale. It is the responsibility of the clinician to inspect the symbols that appear on the packaging of the product in use.

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