
|                  |                                                                             |
|------------------|-----------------------------------------------------------------------------|
| <b>English</b>   | <b>INSTRUCTIONS FOR USE: Southern Implants® Healing abutments</b>           |
| <b>Español</b>   | <b>INSTRUCCIONES DE USO: Pilares de Cicatrización Southern Implants®</b>    |
| <b>Italiano</b>  | <b>ISTRUZIONI PER L'USO: Southern Implants® Monconi di guarigione</b>       |
| <b>Français</b>  | <b>MODE D'EMPLOI : Piliers de cicatrisation Southern Implants®</b>          |
| <b>Deutsch</b>   | <b>GEBRAUCHSANWEISUNGEN: Southern Implants® Gingivaformer</b>               |
| <b>Português</b> | <b>INSTRUÇÕES DE UTILIZAÇÃO: Pilares de cicatrização Southern Implants®</b> |
| <b>Ελληνικά</b>  | <b>ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ: Κολοβώματα επούλωσης της Southern Implants®</b>          |
| <b>Svenska</b>   | <b>ANVISNINGAR FÖR ANVÄNDNING: Southern Implants® Läkdistanser</b>          |
| <b>Română</b>    | <b>INSTRUCȚIUNI DE UTILIZARE: Southern Implants® Bonturi de vindecare</b>   |

**South Africa - Headquarters:** 1 Albert Road, Irene, 0062, RSA  
T: +27-12-667-1046 | E: [info@southernimplants.com](mailto:info@southernimplants.com)

**Southern Implants Europe AB:** Holmgatan 30, S-791 71 Falun, Sweden  
T: +46 23 13300 | E: [ecrep@southernimplants.com](mailto:ecrep@southernimplants.com)

#### Subsidiaries

##### Australia

Southern Implants Australia  
T: +61-(0)-8-9466-2627  
E: [info@southernimplants.com.au](mailto:info@southernimplants.com.au)

##### Spain and Portugal

Southern Implants Iberica  
T: +34 935 053 507  
E: [info@southernimplants.es](mailto:info@southernimplants.es)

##### United Kingdom and Ireland

Southern Implants UK  
T: +44-20-8899-6845 / 6 / 7  
E: [info@southernimplants.co.uk](mailto:info@southernimplants.co.uk)

##### USA and Canada

Southern Implants North America Inc.  
T: +1-561-472-0990  
E: [customercare@southernimplants.com](mailto:customercare@southernimplants.com)

## Description

Southern Implants® Healing abutments are premanufactured prosthetic components intended to be temporarily connected to endosseous dental implants or Compact Conical Abutments during the initial healing phase of dental implant surgery. These abutments are available in various diameters and lengths to enable suitable soft tissue emergence profile for the final prosthesis.

## Intended use

The Healing Abutments are intended to be used as a temporary component attached to endosseous implants following surgery to allow soft tissue healing.

## Indications for use

The Healing Abutments are indicated for use in cases of partial or full edentulousness, congenital tooth loss, or failing/problematic natural dentition, when a patient is dissatisfied with existing removable dentures, or as part of the revision/replacement of previous dental restorations.

## Intended user

The intended user for this system includes Dental Technicians, Maxillo-facial Surgeons, General Dentists, Orthodontists, Periodontists, Prosthodontists, and other appropriately trained and experienced implant users.

## Intended environment

This device is intended to be used in a dental laboratory as part of the restoration design and manufacture as well as in a clinical environment such as an operating theatre or a dentist consultation room.

## Intended patient population

The intended patient population for dental implant therapy includes partially or fully edentulous patients eligible or, otherwise, not contraindicated for implant placement that require prosthetic restoration or revision of existing restorations in the upper or lower jaw; where the planned restorations are to be fixed and comprise of single teeth, partial or full bridges.

## Compatibility information

Southern Implants' implants should be restored with Southern Implants' components. In the Southern Implants' range there are 8 implant / abutment connections. The implant code and connection type can be identified by specific abbreviations in the product codes. Range identifiers are summarised in Table A.

**Table A - Compatible**

| Implant connection type          | Compatible device                                                                                                                                                  |
|----------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| External Hex (EX)                | Parts labelled TPW*, TBN*, WBN*, TV4B*, TB*, WB*, T4B*, T5B*, TBA*, XBA*, WBA*, TV5BA*, T5BA*, T6BA*, T7BA*, TBBB*, WBBB*, T6BBB, T7BBB*, TB9MAX*, TMAX9* and TPN* |
| TRI-NEX® (EL) (Lobe)             | Parts labelled HA-L-(Ø)* and HA-L-(Ø)W*                                                                                                                            |
| Deep Conical (DC)                | Parts labelled HA-DC(Ø)*, HA-DCR(Ø)*, HA-DC(Ø)-W* and HA-DC(Ø)-N*                                                                                                  |
| Internal Hex (M)                 | Parts labelled HA-M-37*, HA-M-45* and HA-M-55* (used with Ø3.75, 4.2 and 5.0 mm platforms)                                                                         |
|                                  | Parts labelled HA-M-P45* (used with Ø5.0 and 6.0 mm during platform matching procedures)                                                                           |
| Internal Hex PROVATA® (3M/ M/ Z) | Parts labelled HA-3M-35* and HA-3M-45* (used with Ø3.3 mm platform)                                                                                                |
|                                  | Parts labelled HA-M-37*, HA-M-45* and HA-M-55* (used with Ø4.0, 5.0 and 6.0 mm platforms)                                                                          |
|                                  | Parts labelled HA-M-P45* (used with Ø5.0 and Ø6.0 mm platforms during platform matching procedures)                                                                |
|                                  | Parts labelled HA-Z6* and HA-Z8* (used with Ø7.0, 8.0 and 9.0 mm platforms)                                                                                        |
| Internal Octagon IT (ITS/ ITS6)  | Parts labelled TT* (used with Ø4.8 mm platforms)                                                                                                                   |
|                                  | Parts labelled TT6* (used with Ø6.5 mm platforms)                                                                                                                  |
| Single Platform (SP)             | Parts labelled HA-SP45* and HA-SP50* (used with Ø3.5, 4.0 and 5.0 mm platforms)                                                                                    |
|                                  | Parts labelled HA-SP55-PM* and HA-SP65-PM* (used with Ø5.0 mm platforms during platform matching procedures)                                                       |
| Compact Conical Abutments        | Parts labelled HMC* and HMCT7* (used with Ø4.8 mm abutment platforms)                                                                                              |
|                                  | Parts labelled HMCW* and HMCTW9* (used with Ø6.0 mm abutment platforms)                                                                                            |

(\*) is indicative of various collar heights available

### Clinical benefits

The Healing Abutments are used as part of dental implant therapy, a procedure in which a patient's missing teeth and/or crowns are restored. As a result, the benefits of the Healing Abutments mirror those of dental implant therapy in general, including improved chewing function, speech, aesthetics, and patient psychological wellbeing.

### Storage, cleaning and sterilisation

The component is supplied sterile (sterilised by gamma irradiation) and intended for single-use prior to the expiration date (see packaging label). Sterility is assured unless the container or seal is damaged or opened. If packaging is damaged do not use the product and contact your Southern representative or return to Southern Implants®. The devices must be stored in a dry place at room temperature and not exposed to direct sunlight. Incorrect storage may influence device characteristics. Do not reuse components indicated for single-use only. Reusing these components may result in:

- damage to the surface or critical dimensions, which may result in performance and compatibility degradation.
- adds the risk of cross-patient infection and contamination if single-use items are reused.

Southern Implants® does not accept any responsibility for complications associated with reused single-use components.

### Contraindications

Since the Healing Abutments are required or utilized as part of dental implant surgery, the contraindications are not dissimilar to those of dental implant therapy as a whole. These contraindications include:

- patients medically unfit for oral surgical procedures
- cases where an adequate number of implants cannot be placed to provide full functional support of the prosthesis.
- patients under the age of 18 years.
- patients with poor bone quality.
- patients with blood disorders.
- presence of infection at the implant site.
- patients with vascular impairment.
- patients with uncontrolled diabetes.
- patients undergoing chronic high dose steroid therapy.
- patients with metabolic bone disease.
- patients with pure titanium, titanium alloy (Ti6Al4V), gold, palladium, polyether ether ketone (PEEK) or iridium allergies or hypersensitivities.

Other than the above, there are no contraindications unique to this system.

### Warnings and precautions

THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING

- For the safe and effective use of dental implants, new technologies/systems and the Healing Abutments, it is strongly suggested that specialised training be undertaken. This training should include hands-on methods to gain competency on the proper technique, biomechanical requirements and radiographic evaluations required for the specific system.
- Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone.
- Use of the device with incompatible or non-corresponding devices can result in device failure or poor performance.
- When handled intraorally, devices must be adequately secured against aspiration. Aspiration of products may lead to infection or physical injury.
- Electrosurgery should not be attempted around metal implants as they are conductive.
- The use of non-sterile items can lead to secondary infections of the tissue or transfer infectious diseases.
- Failure to perform appropriate cleaning, re-sterilisation and storage procedures as per the Instructions for Use document can result in device damage or secondary infection.

- Exceeding the number of recommended uses for reusable devices can result in device damage, secondary infection or patient harm,
- Blunt drills may cause damage to the bone which could compromise osseointegration.
- For short implants, clinicians should closely monitor patients for any of the following conditions: peri implant bone loss, changes to implant's response to percussion, or radiographic changes in bone to implant contact along the implant's length. If the implant shows mobility or greater than 50% bone loss, the implant should be evaluated for possible removal. If the clinicians choose a short implant, then clinicians should consider a two-stage surgical approach, splinting a short implant to an additional implant, and placement of the widest possible fixture. Allow longer periods for osseointegration and avoid immediate loading.

**Note:** Training should be undertaken by both new and experienced implant users before using a new system or attempting to do a new treatment method.

It is important that thorough screening of prospective implant candidates is performed. This screening should include at least:

- a comprehensive medical and dental history.
- visual and radiological inspection to determine adequate bone dimensions, anatomical landmarks, occlusal conditions and periodontal health. Failure to accurately estimate lengths of drills relative to radiographic measurements can result in permanent injury to nerves or other vital structures.
- the evaluation and consideration of bruxism and unfavourable jaw relations.

Proper preoperative planning with a good team approach between well trained surgeons, restorative dentists and lab technicians is essential for successful implant treatment.

For successful implant treatment, it is important to:

- minimise the trauma to the host tissue as this increases the potential for successful osseointegration.
- be aware and avoid damage to vital anatomical structures such as nerves, veins and arteries. Injury to these structures may cause serious complications, including, injury to the eye, nerve damage and excessive bleeding. It is essential to protect the infraorbital nerve.
- identify actual measurements relative to the radiographic data as failing to do so could lead to complications.

Special care should be taken when treating patients with local or systemic risk factors that could affect the healing of the bone and soft tissue or otherwise increase the severity of the side effects, risk of complications and/or implant failure. These factors include:

- smoking.
- history of periodontal disease.
- poor oral hygiene.
- history of orofacial radiotherapy.\*\*
- bruxism and unfavourable jaw relations.

*\*\* The potential for implant failure and other complications increases when implants are placed in irradiated bone as radiotherapy can result in progressive fibrosis of vessels and soft tissue (i.e., osteoradionecrosis), leading to diminished healing capacity. Implant treatment of irradiated patients is dependent upon issues like the timing of implant placement in relation to the radiation therapy, anatomic site chosen for implant placement and radiation dosage at that site and consequent risk of osteoradionecrosis.*

The responsibility for proper patient selection, adequate training and experience in implant placement and the provision of the appropriate information required for informed consent rests with the practitioner.

Should the device not operate as intended, it must be reported to the manufacturer of the device. The contact information or the manufacturer of this device to report a change in performance is:

[sicomplaints@southernimplants.com](mailto:sicomplaints@southernimplants.com).

### Side effects

Since the device is required or utilized as part of dental implant surgery, the expected side effects are not dissimilar to those of dental implant therapy as a whole. Common side effects of dental implant therapy include pain, local

inflammation, phonetic difficulties and gingival inflammation. Other, less common side effects or complications to dental implant therapy include, but are not limited to: (1) wound dehiscence; (2) peri-implantitis; (3) transient weakness, numbness, and/or pain associated with mild nerve injury; (4) fat emboli formation; (5) marginal bone loss within acceptable limits; (6) allergic reaction(s) to implant and/or abutment material; (7) unspecified infection; (8) implant failure due to insufficient levels of osseointegration; (9) breakage of the implant and/or abutment; (10) loosening of the abutment screw and/or retaining screw; (11) complications requiring revision of the dental implant; (12) nerve damage resulting in permanent weakness, numbness, or pain; (13) loosening of the implant requiring revision surgery; (14) perforation of the maxillary sinus; (15) perforation of the labial and/or lingual plates and (16) bone loss possibly resulting in revision or removal of the implant.

### **Precaution: maintaining sterility protocol**

Implants are packaged as follows:

1. An outer package consisting of a rigid, clear box which acts as protection for the inner package.
2. The inner package consisting of a blister pack (clear plastic-formed blister base with a TYVEK “peel-back” lid).
3. Within the inner package, there is a hollow tube which contains one implant suspended from a titanium ring, this ensures the implant never touches the inside of the plastic tube.
4. Labelling information is located on the surface of the peel-back lid and on the outside of the rigid box.

Care must be taken to maintain the sterility of the implant by proper opening of the packaging and handling of the implant.

1. Open the implant package in the non-sterile field, with non-sterile gloves, tear the address label to open the box.
2. With non-sterile gloves, remove the inner blister pack. Do not place the plastic box or blister pack-lid onto the sterile field. The contents of this inner package are sterile.
3. The sealed blister is to be opened by an assistant (with nonsterile gloves): remove the TYVEK lid and drop or place the sterile tube onto the sterile field, open the tube cap and attach the implant placement tool onto the implant and carefully remove from the sterile tube. Do not touch the sterile implant.

Other sterile components are packed in a peel pouch or blister base with a “peel-back” lid. Labelling information is located on the bottom half of the pouch, inside the packet or on the surface of the peel-back lid. Sterility is assured unless the pouch is damaged or opened. Non-sterile components are supplied clean but not sterile in a peel pouch or blister base with peel-back lid. Labelling information is located on the bottom half of the pouch or on the surface of the peel-back lid.

### **Notice regarding serious incidents**

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established.

The contact information for the manufacturer of this device to report a serious incident is as follows:  
[sicomplaints@southernimplants.com](mailto:sicomplaints@southernimplants.com).

### **Materials**

Material type                      Commercially Pure Titanium (Grade 4)

### **Disposal**

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

### **MR safety**

Nonclinical testing has demonstrated that the Southern Implants® dental implants, metallic abutments and prosthetic screws are MR conditional.

A patient with these devices can be safely scanned in a MR system meeting the following conditions:

- static magnetic field of 1.5 Tesla and 3.0 Tesla only.
- maximum spatial gradient magnetic field of 3000 Gauss/cm (30 T/m).
- maximum MR system reported SAR corresponding to Normal Operating mode for all landmarks (Head SAR of 3.2 W/kg for head landmark, 2 W/kg whole body, and appropriate partial body SAR for other landmarks). For imaging landmarks above the thorax, a continuous scan time of 15 minutes will require a cooling delay of at least 5 minutes.
- in the non-clinical testing, the image artifact caused by the device extends approximately 20 mm from the Southern Implants' dental implants, abutments and prosthetic screws, when imaged with a gradient echo pulse sequence and a 3.0 Tesla MRI system.

Removable restorations should be taken out prior to scanning, as is done for watches, jewellery etc.

Should there be no MR symbol on the product label, please note that this device has not been evaluated for safety and compatibility in the MR environment.

### Summary of Safety and Clinical Performance (SSCP)

As required by the European Medical Device Regulation (MDR; EU2017/745), a Summary of Safety and Clinical Performance (SSCP) is available for perusal with regard to Southern Implants® product ranges.

The relevant SSCP can be accessed at <https://ec.europa.eu/tools/eudamed>.

**NOTE:** the above website will be available upon the launch of the European Database on Medical Devices (EUDAMED).

### Disclaimer of liability

This product is part of the Southern Implants® product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants® product range and take full responsibility for the correct indications and use of this product. Southern Implants® does not assume liability for damage due to incorrect use. Please note that some Southern Implants® products may not be cleared or released for sale in all markets.

### Basic UDI

| Product                       | Basic-UDI Number |
|-------------------------------|------------------|
| Basic-UDI for Metal Abutments | 60095440387296   |

### Related literature and catalogues

CAT-2004 - TRI-NEX® Implants Product Catalogue  
 CAT-2005 - IT Implants Product Catalogue  
 CAT-2020 - External Hex Implants Product Catalogue  
 CAT-2042 - Deep Conical Implants Product Catalogue  
 CAT-2043 - Internal Hex Implants Product Catalogue  
 CAT-2060 - PROVATA® Implants Product Catalogue  
 CAT-2069 - INVERTA® Implants Product Catalogue  
 CAT-2070 - Zygomatic Implants Product Catalogue  
 CAT-2093 - Single Platform (SP1) Implants Product Catalogue

## Symbols and warnings

|                                                                                                                                |                                                                                     |                                                                                     |                                                                                     |                                                                                     |                                                                                      |                                                                                       |                                                                                       |                                                                                       |                                                                                       |                                                                                       |                                                                                       |
|--------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
|                                             |  |  |  |  |   |    |   |  |  |  |  |
| Manufacturer:<br>Southern Implants®<br>1 Albert Rd, P.O<br>Box 605 IRENE,<br>0062,<br>South Africa.<br>Tel: +27 12 667<br>1046 | CE mark<br>2797                                                                     | Prescription<br>device*                                                             | Sterilised<br>using<br>irradiation                                                  | Non-sterile                                                                         | Use by date<br>(mm-yy)                                                               | Do not reuse                                                                          | Do not<br>resterilise                                                                 | Catalogue<br>number                                                                   | Batch code                                                                            | Medical<br>device                                                                     | Authorised<br>representative in the<br>European<br>Community                          |
|                                             |  |  |  |  |  |  |  |  |  |                                                                                       |                                                                                       |
| Authorised<br>representative for<br>Switzerland                                                                                | Date of<br>manufacture                                                              | Magnetic<br>Resonance<br>conditional                                                | Magnetic<br>Resonance<br>safe                                                       | Single sterile barrier system<br>with protective packaging<br>inside                | Single sterile<br>barrier<br>system                                                  | Consult<br>instruction<br>for use                                                     | Caution                                                                               | Keep away<br>from sunlight                                                            | Do not use if<br>package is<br>damaged                                                |                                                                                       |                                                                                       |

\* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

All rights reserved. Southern Implants®, the Southern Implants® logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants®. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale. It is the responsibility of the clinician to inspect the symbols that appear on the packaging of the product in use.

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