SOUTHERNIMPLANTS®

Innovative Treatment Solutions

English	INSTRUCTIONS FOR USE: Straumann® ZAGA™ Cylindrical Twist Drills
Español	INSTRUCCIONES DE USO: Fresas helicoidales cilíndricas Straumann® ZAGA™
Italiano	ISTRUZIONI PER L'USO: Straumann® ZAGA™ Punta elicoidali cilindriche
Français	INSTRUCTIONS D'UTILISATION : Straumann® ZAGA™ Forets cylindriques hélicoïdaux
Deutsch	GEBRAUCHSANWEISUNG: Straumann® ZAGA™ Zylinderspiralbohrer
Português	INSTRUÇÕES DE UTILIZAÇÃO: Straumann® ZAGA™ Brocas Helicoidais Cilíndricas

ARCHIVED

South Africa - Headquarters: 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

Subsidiaries

Australia

Southern Implants Australia
T: +61-(0)-8-9466-2627
E: info@southernimplants.com.au

Spain and Portugal

Southern Implants Iberica
T: +34 935 053 507
E: info@southernimplants.es

United Kingdom and Ireland

Southern Implants UK
T: +44-20-8899-6845 / 6 / 7
E: info@southernimplants.co.uk

USA and Canada

Southern Implants North America Inc.
T: +1-561-472-0990
E: customer care@southernimplants.com

Intended use

The ZAGA™ cylindrical twist drills are intended to be used to prepare the osteotomy for implant placement. Cylindrical twist drills are designed with laser markings on the body, corresponding to implant length to assist with drilling depth.

Intended user

Maxillo-facial Surgeons, General Dentists, Orthodontists, Periodontist, Prosthodontists and other appropriately trained and experienced implant users.

Intended environment

The cylindrical twist drills are intended to be used in a clinical environment such as an operating theater or a dentist consultation room.

Intended patient population

This device is used in the dental restoration of partially or fully edentulous patients in the upper or lower jaw. Restorations may comprise single teeth, partial or full bridges, and may be fixed or removable.

Description

The ZAGA™ cylindrical twist drills are dedicated to the Straumann® Zygomatic Implant System and described in Table 1. These drills are made of stainless steel. The drills have a latch dimension compatible to ISO 1797. This is in order to connect the drill to the handpiece of an implant motor unit. Drill sizes are identified by laser markings on the shaft of the drill. The drills have laser markings to indicate the drill depth. Twist drills indicated with 10 Number of uses, can be used up to 10 times or when the cutting efficiency deteriorates. These drills are provided sterile.

Table 1

ZAGA Implants				
Drill code	Material	Coating (if any)	Description of product	Number of uses
D-ZYG-CH-29	Stainless Steel	-	Twist Drill Zygomatic ø2.9	Up to 10
D-ZYG-CH-29S	Stainless Steel	-	Twist Drill Zygomatic ø2.9 (Short)	Up to 10

Indications for use of our implant system

The Straumann® Zygomatic Implants are intended to be implanted in the upper jaw arch to provide support for fixed dental prostheses in patients with partially or fully edentulous maxillae. All implants are appropriate for immediate loading when good primary stability is achieved and with appropriate occlusal loading. This implant system is not intended, nor should it be used, in conjunction with an angled abutment. These implants are not intended for single unit loading.

Indications for use of our cylindrical twist drills

ZAGA™ cylindrical twist drills are indicated for a step-wise drilling approach, when preparing an osteotomy, for Straumann® Zygomatic Implants, in soft, medium or dense bone, by following the drill protocols as recommended in the product catalogues.

Contraindications

Do not use in patients:

- who are medically unfit for dental implant procedures
- where adequate numbers of implants could not be placed to achieve full functional support of the prosthesis,
- who are allergic or have hypersensitivity to pure titanium or titanium alloy (Ti-6Al-4V), gold, palladium, platinum or iridium.
- who are under the age of 18, have poor bone quality, blood disorders, infected implant site, vascular impairment, uncontrolled diabetes, drug or alcohol abuse, chronic high dose steroid therapy, anti-coagulant therapy, metabolic bone disease, radiotherapy treatment.

Warnings

THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING.

- For the safe and effective use of dental implants it is suggested that specialised training be undertaken, including hands-on training to learn proper technique, biomechanical requirements

and radiographic evaluations.

- Responsibility for proper patient selection, adequate training, experience in the placement of implants, and providing appropriate information for informed consent rests with the practitioner. Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone.
- The use of non-sterile items can lead to secondary infections of the tissue or transfer infectious diseases.
- Blunt drills may cause damage to the bone which could compromise osseointegration.

Cautions

New and experienced Implant users should do training before using a new system or attempting to do a new treatment method. Take special care when treating patients who have local or systemic factors that could affect the healing of the bone and soft tissue (i.e. poor oral hygiene, uncontrolled diabetes, are on steroid therapy, smokers, infection in the nearby bone and patients who had oro-facial radiotherapy).

Thorough screening of prospective implant candidates must be performed including:

- a comprehensive medical and dental history.
- visual and radiological inspection to determine adequate bone dimensions, anatomical landmarks, occlusal conditions and periodontal health.
- bruxism and unfavourable jaw relations must be taken into account.
- proper pre-operative planning with a good team approach between well-trained surgeons, restorative dentists and lab technicians is essential for successful implant treatment.
- minimizing the trauma to the host tissue increases the potential for successful osseointegration.
- electro-surgery should not be attempted around metal implants as they are conductive.

During surgery

Care must be taken that parts are not swallowed during any of the procedures, a rubber-dam application is recommended when appropriate. All components, instruments and tooling used during the clinical or laboratory procedure must be maintained in good condition and care must be taken that instrumentation does not damage implants or other components or equipment.

Post-surgery

Regular patient follow-up, and proper oral hygiene must be achieved to ensure favourable long-term results.

Storage, cleaning & sterilisation

These devices are supplied sterile (sterilised by gamma irradiation). Sterility is assured unless the container or seal is damaged or opened. If packaging is damaged do not use the product and contact your Straumann® representative/ or return to Southern Implants. The devices must be stored in a dry place at room temperature and not exposed to direct sunlight. Incorrect storage may influence device characteristics.

Single use devices:

Do not reuse implants, single use drills, cover screws, temporary abutments or abutments. Reusing these components may result in:

- damage to the surface or critical dimensions, which may result in performance and compatibility degradation.
- adds the risk of cross-patient infection and contamination if single-use items are reused.

Southern Implants does not accept any responsibility for complications associated with reused components.

Reusable devices:

Prior to re-using this device, it needs to be inspected, if there are signs of visible corrosion, deformed or twisted connections, dull cutting edges, wear and damage, this device shall be disposed of. After inspection, and if re-use seems fit, the devices are cleaned and sterilised. If re-use seems fit:

- containment: As soon as practically possible, remove all visible residue after use (bone, blood or tissue), by immersing the instrument in cold water (Dried soil is difficult to remove).
- pre-Cleaning: Rinse with lukewarm water for 3 minutes,

and remove hardened debris with a soft nylon brush. Avoid mechanical damage during cleaning.

- manual cleaning or automated cleaning: prepare an ultrasonic bath with suitable detergent, sonicate for 20 minutes (alternative methods can be used if proven by the end user). Rinse with purified / sterile water. Load devices into a thermo-disinfector. Run the cleaning and disinfection cycle, followed by the drying cycle.

NOTE: Always follow the instructions for use of the manufacturers of cleaning agents and disinfectants.

- drying: dry the instruments with filtered compressed air or single use lint-free wipes. Pack the instruments as quickly as possible after removal into the storage container. If additional drying is necessary, dry in a clean location. Moisture on these devices can cause corrosion and deterioration of the cutting edges.
- inspection: Do a visual inspection of the items to check for any damage.
- packaging: Use the correct packaging material as indicated for steam sterilisation to ensure sterility is maintained. Double packaging is recommended.

Sterilisation

Southern Implants recommends the following procedure to sterilise the instruments prior to use/re-use:

Methods to sterilise these devices:

1. Pre-vacuum Sterilisation method: Steam sterilise the instruments at 132°C (270°F) at 180-220kPa for 4 minutes. Dry for at least 20 minutes in the chamber. Only an approved wrap or pouch for steam sterilisation must be used.
2. Pre-vacuum sterilisation method: Wrapped, steam sterilise at 135°C (275°F) for 3 minutes. Dry for 20 minutes in the chamber. Use a wrap or pouch that is cleared for the indicated steam sterilisation cycle.

NOTE: Users in the USA must ensure that the steriliser, wrap or pouch, and all steriliser accessories are cleared by the FDA, for the intended sterilisation cycle.

Clinical procedure

A proper clinical and radiological evaluation must be done to determine the bone dimensions and bone quality.

Surgical procedure

Ensure that all instruments and drills are in a good condition. Blunt drills may cause damage to the bone which could compromise osseointegration.

The drill sizes are identified by different laser markings on the shaft of the drill. The drills have different laser markings on the body of the drill which corresponds to the length of implant being placed (Figure 1). Refer to individual product catalogues for more information.

Figure 1

D-ZYG-CH-29

D-ZYG-CH-29S


1. Drill in the planned direction to the full depth of the implant length being placed as indicated on the markings on the drill. Cylindrical drills extend up to 1mm longer than the implant, when seated. Allow for this additional length when drilling near vital anatomical structures.

NOTE: Connect the drill latch to the handpiece. If the latch is not inserted fully into the handpiece the torque is applied to the latch, resulting in possible twisting of the latch or damage to the handpiece. Consult the instructions for use of the handpiece to ensure proper engagement of the latch.

2. Drill at sufficient speed (between 1000-1500rpm for Cylindrical Drills), under constant irrigation with sterile saline.
3. Use an up-and-down motion with the handpiece, without stopping the motor. This will allow the irrigation to flush away bone debris on the drill.
4. Insert the Direction Indicator (ZYG-TR-55-35/45/52/52.5) to verify the alignment with adjacent implants or teeth. A radiograph can be taken at this stage to verify depth and direction.
5. The depth can be determined by a depth gauge (Figure 2):
6.
 - I-ZYG-DG-1
 - CH-I-DG

Figure 2

CH-I-DG

7. 6. During surgery the clinician will be able to assess the bone quality and should adjust drill protocols when necessary, to prepare the site. This is to avoid the implant getting stuck before it is properly seated in the osteotomy. It is recommended to undersize the osteotomy in soft bone.
8. 7. Preparing the site further should involve: ensuring the drill goes to full depth.

NOTE:

- avoid lateral pressure (bending) on the drills during drilling procedures.
- lateral pressure to the drill can cause drill fracture.
- verify the drill is securely locked into the handpiece before the drilling procedure starts.

These twist drills can be re-used. It is recommended to maintain a log of these drills, recording the number of uses. Prior to re-processing these components, it should be thoroughly inspected and tested to determine its suitability for re-use.

NOTE:

- do not apply more than 40-45Ncm to any latch type drill/instrument, this could cause damage to the handpiece and latch of the instrument.
- blunt drills cause excessive torque and result in damage to the

handpiece or drill latch.

Materials

Drills: Stainless Steel 1.4197

Clinical benefits

Patients can expect to have their missing teeth replaced and/ or crowns restored.

Healing

The healing time required for osseointegration depends on the individual and treatment protocol. It is the responsibility of the practitioner to decide when the implant can be restored. Good primary stability will govern if immediate loading can be done.

Implant care and maintenance

Potential implant patients should establish an adequate oral hygiene regime prior to Implant therapy. Proper post-operative, oral hygiene and implant maintenance instructions must be discussed with the patient, as this will determine the longevity and health of the Implants. The patient should maintain regular prophylaxis and evaluation appointments.

Side effects

Potential Side Effects and Temporary symptoms: Pain, swelling, phonetic difficulties, gingival inflammation. More persistent symptoms: The risks and complications with implants include, but are not limited to: (1) allergic reaction(s) to implant and/or abutment material; (2) breakage of the implant and/or abutment; (3) loosening of the abutment screw and/ or retaining screw; (4) infection requiring revision of the dental implant; (5) nerve damage that could cause permanent weakness, numbness, or pain; (6) histologic responses possibly involving macrophages and/ or fibroblasts; (7) formation of fat emboli; (8) loosening of the implant requiring revision surgery; (9) perforation of the maxillary sinus; (10) perforation of the labial and lingual plates; and (11) bone loss possibly resulting in revision or removal.

Breakage

Implant and abutment fractures can occur when applied loads exceed the tensile or compressive strength of the material. Potential overloading conditions may result from; deficiencies in implant numbers, lengths and/ or diameters to adequately support a restoration, excessive cantilever length, incomplete abutment seating, abutment angles greater than 30 degrees, occlusal interferences causing excessive lateral forces, patient parafunction (e.g., bruxing, clenching), loss or changes in dentition or functionality, inadequate prosthesis fit, and physical trauma. Additional treatment may be necessary when any of the above conditions are present to reduce the possibility of hardware complications or failure.

Changes in performance

It is the responsibility of the clinician to instruct the patient on all appropriate contraindications, side effects, and precautions as well as the need to seek the services of a trained dental professional if there are any changes in the performance of the implant (e.g., looseness of the prosthesis, infection or exudate around the implant, pain, or any other unusual symptoms that the patient has not been told to expect).

Disposal

Disposal of the device and its packaging; Follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

Disclaimer of liability

This product is part of the Southern Implants product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants product range and take full responsibility for the correct indications and use of this product. Southern Implants does not assume liability for damage due to incorrect use. Please note that some Southern Implants products may not be cleared or released for sale in all markets.

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and / or patient is established. The contact information for the manufacturer of this device to report a serious incident is as follows: sicomplaints@southernimplants.com

Basic UDI

Product	Basic-UDI Number
Basic-UDI for Drills and Hand Piece Devices	600954403875
Basic-UDI for Reusable Instruments	600954403876

Related literature & catalogues

CAT-8049 - Straumann® ZAGA™ Zygomatic Implants

Symbols and Warnings

 Manufacturer: Southern Implants 1 Albert Rd, P.O Box 605 IRENE, 0062, South Africa. Tel: +27 12 667 1046	 2797	 Prescription device*	 Sterilization using Irradiation	 Non-sterile	 Caution	 Consult instruction for use	 Use by date (mm-yy)	 Do not reuse	 Do not re-sterilize	 Batch code	 Do not use if package is damaged	 Medical Device
* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.						Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.						
All rights reserved. Southern Implants®, the Southern Implants logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale												

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