English	INSTRUCTIONS FOR USE: Southern Implants® Implant Zygomatic Implants
Español	INSTRUCCIONES DE USO: Implantes Cigomáticos Southern Implants®
Italiano	ISTRUZIONI PER L'USO: Southern Implants® Impianti Zigomatici Impianti
Français	MODE D'EMPLOI : Implants Southern Implants® Implants Zygomaticques
Deutsch	GEBRAUCHSANWEISUNGEN: Southern Implants® Implant Zygomatic-Implantate
Português	INSTRUÇÕES DE UTILIZAÇÃO: Southern Implants® Implantes zigomáticos
Ελληνικά	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ: Southern Implants® Zygomatic Implants
Svenska	BRUKSANVISNING: Southern Implants® Zygomatic Implantat
Română	INSTRUCȚIUNI DE UTILIZARE: Implanturi Zigomatice Southern Implants®

South Africa - Headquarters: 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

EC REP

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

Subsidiaries

Australia
Southern Implants Australia
T: +61-(0)-8-9466-2627
E: info@southernimplants.com.au

Spain and Portugal
Southern Implants Iberica
T: +34 935 053 507
E: info@southernimplants.es

United Kingdom and Ireland
Southern Implants UK
T: +44-20-8899-6845 / 6 / 7
E: info@southernimplants.co.uk

USA and Canada
Southern Implants North America Inc.
T: +1-561-472-0990
E: customer care@southernimplants.com

Description

The Southern Implants Zygomatic range includes the standard Zygomatic, the ZYGAN, the ZYGIN, the Oncology and the ZYGEX Implants. The implants are up to 60mm long to enable anchorage in the zygoma and have a 55° head angle. They are made from biocompatible, commercially pure Titanium grade 4 and are available in a range of lengths to be used with a range of prosthetic components. All standard Zygomatic implants in this range present the MSc feature, a 6mm coronal area of specific roughness ($Sa = 0.6 \pm 0.2$).

Intended use

The zygomatic implants are intended to treat partially or fully edentulous patients with severely resorbed or absent maxillae for whom conventional implants are not an option as a means of fixing a permanent or removable dental or maxillofacial prosthesis.

The system constituents are classified as medical devices and are intended for single use on a single patient.

Indications for use

The Southern Implants Zygomatic standard implant, the ZYGAN (narrow apex), ZYGIN (narrow apex and coronal implant head), the Oncology and the ZYGEX (narrow apex) implants are intended to be implanted in the upper jaw to provide support for fixed or removable dental prostheses in patients with partially or fully edentulous maxilla. All implants are appropriate for immediate loading when good primary stability is achieved and with appropriate occlusal loading.

Intended user

Maxillo-facial Surgeons, General Dentists, Orthodontists, Periodontists, Prosthodontists, and other appropriately trained and experienced implant users.

Intended environment

The devices are intended to be used in a clinical environment such as an operating theatre or a dentist consultation room.

Intended patient population

The intended patient population for the device is not dissimilar to that of dental implant therapy.

The intended patient population for implant therapy is partially or fully edentulous patients requiring prosthetic dental restoration in the upper or lower jaw. Restorations may comprise single teeth, partial or full bridges and are fixed restorations.

Compatibility information

Use only original Southern Implants components to restore Southern Implants Zygomatic ranges. Use components that correspond to the connection type, and prosthetic platform when restoring Zygomatic implants. For further information please see Zygomatic Product Catalogue CAT-2070.

NOTE: angled Compact Conical abutments are not indicated to be used with Southern Implants Zygomatic implant ranges.

Table A - Compatible

Item code	Implant length codes (in mm)
ONC-55-	27.5N/ 32.5N/ 37.5N/ 42.5N/ 47.5N
ZYGEX-	30/ 32.5/ 35/ 37.5/ 40/ 42.5/ 45/ 47.5/ 50/ 52.2/ 55/ 57.5/ 60
ZYG-55-	35N/ 37.5N/ 40N/ 42.5N/ 45N/ 47.5N/ 50N/ 52.5N/ 55N/ 60N
ZYGAN-	30/ 32.5/ 35/ 37.5/ 40/ 42.5/ 45/ 47.5/ 50/ 52.2/ 55/ 57.5/ 60
ZYGIN-	30/ 32.5/ 35/ 37.5/ 40/ 42.5/ 45/ 47.5/ 50/ 52.2/ 55/ 57.5/ 60

Applicable Implants	Cover Screw and Driver	Abutment and Driver	Prosthetic Screw Driver
All Zygomatic Implants with the exception of the ZYGIN implant.	SCU2 (Cover screw) I-CS-HD (driver)	AMCZ (Screw retained abutment), I-HAD (Driver)	1 Series screw (prosthetic screw), I-HD-M (driver)
ZYGIN	CS-ZYG (Cover screw) I-CS-HD (Driver)	MC-ZYG (Screw retained abutment), I-HAD (Driver)	1 Series screw (prosthetic screw), I-HD-M (driver)

Clinical benefits

Clinical benefits of zygomatic implant therapy are that patients can expect to have their missing teeth replaced and/or crowns restored.

Storage, cleaning and sterilisation

The component is supplied sterile (sterilised by gamma irradiation) and intended for single-use prior to the expiration date (see packaging label). Sterility is assured unless the container or seal is damaged or opened. If packaging is damaged do not use the product and contact your Southern representative or return to Southern Implants®. The devices must be stored in a dry place at room temperature and not exposed to direct sunlight. Incorrect storage may influence device characteristics. Do not reuse components indicated for single-use only. Reusing these components may result in:

- damage to the surface or critical dimensions, which may result in performance and compatibility degradation.
- adds the risk of cross-patient infection and contamination if single-use items are reused.

Southern Implants® does not accept any responsibility for complications associated with reused single-use components.

Contraindications

Contra-indications include:

- patients medically unfit for oral surgical procedures.
- cases where an inadequate number of implants can be placed to provide the required functional support for the prosthesis.
- minors (patients under the age of 18 years).
- patients with poor bone quality.
- patients with blood disorders.
- presence of infection at the implant site.
- patients with vascular impairment.
- patients with uncontrolled diabetes.
- patients with drug or alcohol abuse dependencies.
- patients undergoing chronic high dose steroid therapy.
- patients undergoing anti-coagulant therapy.
- patients with metabolic bone disease.
- patients undergoing radiotherapy treatment.*
- patients with pure titanium, titanium alloy (Ti6Al4V), gold, palladium, or iridium allergies or hypersensitivities.
- sinus pathology.
- inadequate bone volume for zygomatic implants.

These contraindications are not dissimilar to those of general implant therapy. Proceeding with implant therapy in contraindicated cases will increase the incidence and/or risk of hazardous situations.

** is not a contraindication for the Oncology Zygomatic Implant.*

Warnings and precautions

THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING.

- For the safe and effective use of dental implants it is strongly suggested that specialised training be undertaken, including hands-on training to learn proper technique, biomechanical requirements and radiographic evaluations.
- Products must be secured against aspiration when handled intraorally. Aspiration of products may lead to infection or unplanned physical injury.

Responsibility for proper patient selection, adequate training, experience in the placement of implants and providing appropriate information for informed consent rests with the practitioner. Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone. Implant failure increases when implants are placed in irradiated bone as radiotherapy can result in progressive fibrosis of vessels and soft tissue, leading to diminished healing capacity.

Additionally, use of Zygomatic implants in bone tissue which has been irradiated as part of cancer therapy may result in the following:

- delayed or failed osseointegration of implants due to reduced bone vascularity, clinically expressed as osteoradionecrosis.
- tissue dehiscence and osteoradionecrosis.
- implant failure and loss.
- implant treatment of irradiated patients is dependent upon issues like the timing of implant placement in relation to the radiation therapy, anatomic site chosen for implant placement and radiation dosage at that site and consequent risk of osteoradionecrosis.

It is important to be aware and avoid damage to vital structures such as nerves, veins and arteries. Injuries to vital anatomical structures may cause serious complications such as injury to the eye, nerve damage and excessive bleeding. It is essential to protect the infraorbital nerve. Failing to identify actual measurements relative to the radiographic data could lead to complications.

New and experienced implant users should do training before using a new system or attempting to do a new treatment method. Take special care when treating patients who have local or systemic factors that could affect the healing of the bone and soft tissue (i.e. poor oral hygiene, uncontrolled diabetes, are on steroid therapy, smokers, infection in the nearby bone and patients who had orofacial radiotherapy).

Thorough screening of prospective implant candidates must be performed including:

- a comprehensive medical and dental history.
- visual and radiological inspection to determine adequate bone dimensions, anatomical landmarks, occlusal conditions and periodontal health.
- bruxism and unfavourable jaw relations must be taken into account.
- proper preoperative planning with a good team approach between well trained surgeons, restorative dentists and lab technicians is essential for successful implant treatment.
- minimising the trauma to the host tissue increases the potential for successful osseointegration.
- electrosurgery should not be attempted around metal implants as they are conductive.

Should the device not operate as intended, it must be reported to the manufacturer of the device. The contact information for the manufacturer of this device to report a change in performance is:

sicomplaints@southernimplants.com.

Side effects

For Zygomatic Implants risk matrix the following side effects have been identified:

- Pain
- Swelling
- phonetic difficulties
- gingival inflammation

Less common but more persistent symptoms include, but are not limited to:

- allergic reaction(s) to implant and/or abutment material
- breakage of the implant and/or abutment
- loosening of the abutment screw and/or retaining screw
- infection requiring revision of the dental implant
- nerve damage resulting in permanent weakness, numbness, or pain
- histologic responses with possible macrophage and/or fibroblast involvement
- fat emboli formation
- loosening of the implant requiring revision surgery
- perforation of the maxillary sinus
- perforation of the labial and lingual plates
- bone loss possibly resulting in revision or removal of the implant

Other than the above, there are no side-effects unique to this system.

Precaution: maintaining sterility protocol

1. Care must be taken to maintain the sterility of the implant by proper opening of the packaging and handling of the implant.
2. The outer rigid plastic box and the outside of the inner plastic blister tray-lid are not sterile; do not touch the outside with sterile gloves and do not place the plastic box or blister tray-lid onto the sterile field.
3. Inside the plastic box, the sealed inner plastic blister and peel back TYVEK lid is sterile only on the inside. The sealed blister is to be opened by an assistant (with non sterile gloves). Remove the TYVEK lid and do not touch the sterile implant.
4. Follow the instructions illustrated in (Figure 1 to Figure 4) to remove the sterile implant, maintaining sterility and to attach the fixture mount and implant to the handpiece.
5. Maintain the sterility of the implant, after opening the tray and removing the implant, until placement in the surgical site.

NOTE: white gloves and background represent non-sterile items. Blue gloves and background represent sterile items.

Figure 1: to open implant package in the non-sterile field, with non-sterile gloves, tear the tamper-proof label to open the box. With nonsterile gloves

Figure 2: the sealed blister must be opened by an assistant (non-sterile gloves). Peel back the TYVEK lid and present the open blister to the surgeon.

Figure 3: without touching the outside of the blister, the surgeon removes the implant holder. Take care to not touch the implant.


Figure 4: engage the handpiece placement tool (I-CON-X) onto the fixture mount and separate the implant from the titanium clip on the implant holder

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established.

The contact information for the manufacturer of this device to report a serious incident is as follows:
sicomplaints@southernimplants.com.

Materials

Material type	Commercially pure titanium grade 4 (ASTM F67 and ISO 5832-2, UTS ≥ 900MPa)
---------------	---

Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

MR safety

Nonclinical testing has demonstrated that the Southern Implants® dental implants, metallic abutments and prosthetic screws are MR conditional.

A patient with these devices can be safely scanned in a MR system meeting the following conditions:

- static magnetic field of 1.5 Tesla and 3.0 Tesla only.
- maximum spatial gradient magnetic field of 3000 Gauss/cm (30 T/m).
- maximum MR system reported SAR corresponding to Normal Operating mode for all landmarks (Head SAR of 3.2 W/kg for head landmark, 2 W/kg whole body, and appropriate partial body SAR for other landmarks). For imaging landmarks above the thorax, a continuous scan time of 15 minutes will require a cooling delay of at least 5 minutes.
- in the non-clinical testing, the image artifact caused by the device extends approximately 20 mm from the Southern Implants' dental implants, abutments and prosthetic screws, when imaged with a gradient echo pulse sequence and a 3.0 Tesla MRI system.

Removable restorations should be taken out prior to scanning, as is done for watches, jewellery etc.

Should there be no MR symbol on the product label, please note that this device has not been evaluated for safety and compatibility in the MR environment.

Summary of Safety and Clinical Performance (SSCP)

As required by the European Medical Device Regulation (MDR; EU2017/745), a Summary of Safety and Clinical Performance (SSCP) is available for perusal with regard to Southern Implants® product ranges.

The relevant SSCP can be accessed at <https://ec.europa.eu/tools/eudamed>.

NOTE: the above website will be available upon the launch of the European Database on Medical Devices (EUDAMED).

Disclaimer of liability

This product is part of the Southern Implants® product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants® product range and take full responsibility for the correct indications and use of this product. Southern Implants® does not assume liability for damage due to incorrect use. Please note that some Southern Implants® products may not be cleared or released for sale in all markets.

Basic UDI

Product	Basic-UDI Number
Basic-UDI for Zygomatic Implants	60095440387194

Related literature and catalogues

CAT-2070- Zygomatic Implants Product Catalogue

Symbols and warnings

											
Manufacturer: Southern Implants® 1 Albert Rd, P.O Box 605 IRENE, 0062, South Africa. Tel: +27 12 667 1046	CE mark 2797	Prescription device*	Sterilised using irradiation	Non-sterile	Use by date (mm-yy)	Do not reuse	Do not resterilise	Catalogue number	Batch code	Medical device	Authorised representative in the European Community
											
Authorised representative for Switzerland	Date of manufacture	Magnetic Resonance conditional	Magnetic Resonance safe	Single sterile barrier system with protective packaging inside	Single sterile barrier system	Consult instruction for use	Caution	Keep away from sunlight	Do not use if package is damaged		

* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

All rights reserved. Southern Implants®, the Southern Implants® logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants®. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale. It is the responsibility of the clinician to inspect the symbols that appear on the packaging of the product in use.

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