English	INSTRUCTIONS FOR USE: Southern Implants® Zygomatic Implants
Español	INSTRUCCIONES DE USO: Implantes Zigomáticos Southern Implants®
Italiano	ISTRUZIONI PER L'USO: Southern Implants® Impianti Zigomatici
Français	MODE D'EMPLOI : Implants zygomatiques de Southern Implants®
Deutsch	GEBRAUCHSANWEISUNG: Southern Implants® Zygomatic-Implantate
Português	INSTRUÇÕES DE UTILIZAÇÃO: Southern Implants® Implantes Zigomáticos
Ελληνικά	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ: Ζυγωματικά εμφυτεύματα της Southern Implants®
Svenska	ANVISNINGAR FÖR ANVÄNDNING: Southern Implants® Zygomatic Implantat
Română	INSTRUCȚIUNI DE UTILIZARE: Implanturi zigom atice Southern Implants®
Türkçe	KULLANMA TALİMATI: Southern Implants® Zigomatik İmplantları

South Africa - Headquarters: 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

Subsidiaries

Australia

Southern Implants Australia
T: +61-(0)-8-9466-2627
E: info@southernimplants.com.au

Spain and Portugal

Southern Implants Iberica
T: +34 935 053 507
E: info@southernimplants.es

United Kingdom and Ireland

Southern Implants UK
T: +44-20 8059 4490
E: info@southernimplants.co.uk

USA and Canada

Southern Implants North America Inc.
T: +1-561-472-0990
E: customer care@southernimplants.com

Description

The Southern Implants Zygomatic range includes the standard Zygomatic, the ZYGAN, the ZYGIN, the ZYGON, the Wide ZYGIN, the Oncology and the ZYGEX Implants. The implants are up to 60mm long to enable anchorage in the zygoma and have a 55° head angle. They are made from biocompatible, commercially pure Titanium grade 4 and are available in a range of lengths to be used with a range of prosthetic components. All standard Zygomatic implants in this range present the MSc feature, a 6mm coronal area of specific roughness ($Sa = 0.6 \pm 0.2$).

Intended use

The Southern Implants® Zygomatic Implants are intended to be surgically placed in the zygoma in order to provide a support structure and point of attachment for prosthetic devices.

Indications for use

The Southern Implants® Zygomatic Implants are indicated for the functional and aesthetic restoration of the partially or fully edentulous upper jaw using multiple-unit prosthetic devices in cases where the patient presents with a severely resorbed or absent maxilla or for whom conventional implants are not an option.

Unless stated otherwise in the specific indications, these devices can be placed with immediate function if sufficient primary stability is achieved with appropriate occlusal loading.

Intended user

Maxillofacial Surgeons, General Dentists, Orthodontists, Periodontists, Prosthodontists, and other appropriately trained and/or experienced medical professionals.

Intended environment

The devices are intended to be used in a clinical environment such as an operating theatre or a dentist consultation room.

Intended patient population

The intended patient population for zygomatic dental implant therapy includes partially or fully edentulous patients with severely atrophic maxillae for whom conventional implants are unsuited where these patients are eligible or, otherwise, not contraindicated for zygomatic implant placement.

Compatibility information

Use only original Southern Implants components to restore Southern Implants Zygomatic ranges. Use components that correspond to the connection type, and prosthetic platform when restoring Zygomatic implants. For further information please see Zygomatic Product Catalogue CAT-2070.

NOTE: angled Compact Conical abutments are not indicated to be used with Southern Implants Zygomatic implant ranges.

Table A - Compatible

Item code	Implant length codes (in mm)
ONC-55-	27.5N/ 32.5N/ 37.5N/ 42.5N/ 47.5N
ZYGEX-	30/ 32.5/ 35/ 37.5/ 40/ 42.5/ 45/ 47.5/ 50/ 52.2/ 55/ 57.5/ 60
ZYG-55-	35N/ 37.5N/ 40N/ 42.5N/ 45N/ 47.5N/ 50N/ 52.5N/ 55N/ 60N
ZYGAN-	30/ 32.5/ 35/ 37.5/ 40/ 42.5/ 45/ 47.5/ 50/ 52.2/ 55/ 57.5/ 60
ZYGIN-	30/ 32.5/ 35/ 37.5/ 40/ 42.5/ 45/ 47.5/ 50/ 52.2/ 55/ 57.5/ 60
ZYGIN-W-	30/ 32.5/ 35/ 37.5/ 40/ 42.5/ 45/ 47.5/ 50/ 52.2/ 55/ 57.5/ 60
ZYGON-	30/ 32.5/ 35/ 37.5/ 40/ 42.5/ 45/ 47.5/ 50/ 52.2/ 55/ 57.5/ 60

Applicable Implants	Cover Screw and Driver	Abutment and Driver	Prosthetic Screw Driver
ONC-55, ZYGEX, ZYG-55, ZYGAN	SCU2 (Cover screw) I-CS-HD (driver)	AMCZ (Screw retained abutment) I-HAD (Driver)	1 Series screw (prosthetic screw) I-HD-M (driver)

ZYGIN, ZYGIN-W, ZYGON	CS-ZYG (Cover screw) I-CS-HD (Driver)	MC-ZYG (Screw retained abutment) I-HAD (Driver)	1 Series screw (prosthetic screw) I-HD-M (driver)
--------------------------	--	--	--

Clinical performance

The Zygomatic Implants are intended to provide a rigid retention method for fixed or removable prostheses in the restorative treatment of the partially or fully edentulous upper jaw that presents as severely resorbed or absent. The rigid and secure retention of the prosthetic restoration is accomplished through the stability of the implant itself. This stability is achieved through the fixture's surgical implantation into the bone and subsequent successful osseointegration. Consequently, the clinical performance of the Zygomatic Implants is primarily defined by the degree of osseointegration achieved. This performance can be quantitatively assessed through implant survival, which serves as an important indicator of the implant's effectiveness and long-term success in clinical practice.

Clinical benefits

The clinical benefits associated with the prosthetic restoration of partially or fully edentulous maxillae are significant, as this treatment can greatly enhance the functional and aesthetic outcomes of the patient. Therefore, the clinical benefit of the Zygomatic Implants encompasses the successful restoration of the edentulous maxillae and the patient's satisfaction with the outcome. Furthermore, since this treatment addresses the functional impact of partial or full edentulism of the maxilla, it is anticipated that a successfully restored patient will experience an improvement in masticatory function.

Storage, cleaning and sterilisation

The component is supplied sterile (sterilised by gamma irradiation) and intended for single-use prior to the expiration date (see packaging label). Sterility is assured unless the container or seal is damaged or opened. If packaging is damaged do not use the product and contact your Southern representative or return to Southern Implants®. The devices must be stored in a dry place at room temperature and not exposed to direct sunlight. Incorrect storage may influence device characteristics.

Do not reuse components indicated for single-use only. Reusing these components may result in:

- damage to the surface or critical dimensions, which may result in performance and compatibility degradation.
- adds the risk of cross-patient infection and contamination if single-use items are reused.

Southern Implants® does not accept any responsibility for complications associated with reused single-use components.

Contraindications

The contraindications of all device groups used as part of the specific treatment or procedure apply. Therefore, the contraindications of the systems/medical devices utilised as part of implant surgery/therapy should be noted and the relevant documents consulted.

The contraindications specific to this device group include:

- patients medically unfit for oral surgical procedures
- cases where an adequate number of implants cannot be placed to provide full functional support of the prosthesis
- patients with inadequate bone quality/quantity
- patients with uncontrolled bleeding disorders
- presence of infection at the implant site
- patients with severe vascular impairment
- patients with uncontrolled endocrine disorders
- patients undergoing high dose steroid therapy
- patients with metabolic bone disease
- patients with incomplete mandibular or maxillary growth
- patients with a weakened immune system or inadequate wound healing capacity

- patients with allergies or hypersensitivity to the material used (titanium)

Warnings and precautions

IMPORTANT NOTICE: THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING.

- For the safe and effective use of dental implants it is strongly suggested that specialised training be undertaken, including hands-on training to learn proper technique, biomechanical requirements and radiographic evaluations.
- Products must be secured against aspiration when handled intraorally. Aspiration of products may lead to infection or unplanned physical injury.

Responsibility for proper patient selection, adequate training, experience in the placement of implants and providing appropriate information for informed consent rests with the practitioner. Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone. Implant failure increases when implants are placed in irradiated bone as radiotherapy can result in progressive fibrosis of vessels and soft tissue, leading to diminished healing capacity.

Additionally, use of Zygomatic implants in bone tissue which has been irradiated as part of cancer therapy may result in the following:

- delayed or failed osseointegration of implants due to reduced bone vascularity, clinically expressed as osteoradionecrosis.
- tissue dehiscence and osteoradionecrosis.
- implant failure and loss.
- implant treatment of irradiated patients is dependent upon issues like the timing of implant placement in relation to the radiation therapy, anatomic site chosen for implant placement and radiation dosage at that site and consequent risk of osteoradionecrosis.

It is important to be aware and avoid damage to vital structures such as nerves, veins and arteries. Injuries to vital anatomical structures may cause serious complications such as injury to the eye, nerve damage and excessive bleeding. It is essential to protect the infraorbital nerve. Failing to identify actual measurements relative to the radiographic data could lead to complications.

New and experienced implant users should do training before using a new system or attempting to do a new treatment method. Take special care when treating patients who have local or systemic factors that could affect the healing of the bone and soft tissue (i.e., poor oral hygiene, uncontrolled diabetes, are on steroid therapy, smokers, infection in the nearby bone and patients who had orofacial radiotherapy).

Thorough screening of prospective implant candidates must be performed including:

- a comprehensive medical and dental history.
- visual and radiological inspection to determine adequate bone dimensions, anatomical landmarks, occlusal conditions and periodontal health.
- bruxism and unfavourable jaw relations must be taken into account.
- proper preoperative planning with a good team approach between well trained surgeons, restorative dentists and lab technicians is essential for successful implant treatment.
- minimising the trauma to the host tissue increases the potential for successful osseointegration.
- electrosurgery should not be attempted around metal implants as they are conductive.

Should the device not operate as intended, it must be reported to the manufacturer of the device. The contact information for the manufacturer of this device to report a change in performance is: sicomplaints@southernimplants.com.

Side effects

The clinical outcome of treatment is influenced by multiple variables and there are potential side effects and residual risks that are associated with the device group and may lead to the requirement of additional treatment. Furthermore, these side effects and residual risks can occur with varying possible severities and frequencies.

The common side effects and residual risks that are associated with the device include:

- pain or tenderness
- localized inflammation
- discomfort
- facial or periorbital hematoma or bruising
- mild to moderate nasal bleeding for 1 - 3 days

The less common and potentially more severe side effects and residual risks that are associated with the device include:

- peri-implantitis
- allergic reaction(s) to implant and/or abutment material
- intraoral soft tissue infection
- wound dehiscence or improper healing
- gingival hyperplasia
- subcutaneous malar emphysema
- infection in the maxillary sinus
- sinusitis
- anatomical landmark damage
- orbital cavity penetration
- perforation of the maxillary sinus
- perforation of the labial and lingual plates
- nerve damage resulting in permanent weakness, numbness, or pain
- nerve damage resulting in transient weakness, numbness, or pain
- formation of an oroantral fistula
- fat emboli formation
- breakage of the implant and/or abutment
- loosening of the abutment screw and/or retaining screw
- loosening of the implant
- insufficient levels of osseointegration or excessive bone loss
- additional surgical intervention (e.g., revision surgery or implant removal) surgery
- implant failure
- prosthetic failure

The combination of the two lists represents the side effects and residual risks associated with the Southern Implants® Zygomatic Implants.

Precaution: maintaining sterility protocol

1. Care must be taken to maintain the sterility of the implant by proper opening of the packaging and handling of the implant.
2. The outer rigid plastic box and the outside of the inner plastic blister tray-lid are not sterile; do not touch the outside with sterile gloves and do not place the plastic box or blister tray-lid onto the sterile field.
3. Inside the plastic box, the sealed inner plastic blister and peel back TYVEK lid is sterile only on the inside. The sealed blister is to be opened by an assistant (with non-sterile gloves). Remove the TYVEK lid and do not touch the sterile implant.
4. Follow the instructions illustrated in (Figure 1 to Figure 4) to remove the sterile implant, maintaining sterility and to attach the fixture mount and implant to the handpiece.

5. Maintain the sterility of the implant, after opening the tray and removing the implant, until placement in the surgical site.

NOTE: white gloves and background represent non-sterile items. Blue gloves and background represent sterile items.

Figure 1: to open implant package in the non-sterile field, with non-sterile gloves, tear the tamper-proof label to open the box. With nonsterile gloves

Figure 2: the sealed blister must be opened by an assistant (non-sterile gloves). Peel back the TYVEK lid and present the open blister to the surgeon.

Figure 3: without touching the outside of the blister, the surgeon removes the implant holder. Take care to not touch the implant.

Figure 4: engage the handpiece placement tool (I-CON-X) onto the fixture mount and separate the implant from the titanium clip on the implant holder

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established.

The contact information for the manufacturer of this device to report a serious incident is as follows:

sicomplaints@southernimplants.com.

Materials

Material type Commercially pure titanium (grade 4) according to ASTM F67 and ISO 5832-2
(UTS $\geq$ 900MPa)

Chemical Components	Residuals (Fe, O, C, N, H) in total	Titanium (Ti)
Composition, % (mass/mass)	<1.05	Balance

Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

Magnetic Resonance (MR) Safety

Nonclinical testing has demonstrated that the Southern Implants® dental implants, metallic abutments and prosthetic screws are MR conditional.

MRI Safety Information	
 MR Conditional	
A patient with the dental implants, metallic abutments and prosthetic screws from Southern Implants® may be safely scanned under the following conditions. Failure to follow these conditions may result in injury to the patient.	
Device Name	dental implants, metallic abutments and prosthetic screws
Static Magnetic Field Strength (B₀) [T]	1.5T or 3T
Maximum Spatial Field Gradient [T/m] [gauss/cm]	30 T/m (3000 gauss/cm)
Type of RF Excitation	Circularly polarized (CP) (i.e., quadrature driven)
RF Transmit Coil Type	Integrated Whole Body Transmit Coil
RF Receive Coil	Any receive-only RF coil may be used
Operating Mode	Normal Operating Mode
RF Conditions	Normal Operating mode: (Head SAR of 3.2 W/kg, 2 W/kg whole body)
Scan Duration	If the closest part of implant <25cm away from the isocenter in feet-head direction: a continuous scan time of 15 minutes will require a cooling delay of at least 5 minutes. If the closest part of implant $\geq$ 25cm away from the isocenter in feet-head direction: This is limited to 1 hour of scan duration.
MR Image Artifact	The presence of this implant may produce an image artifact.

Additional information:

1. Local transmit coil: Local transmit coil safety was not evaluated. Local transmit coil is only allowed if the implant is $\geq 25\text{cm}$ away from isocenter and for a wbSAR of up to 2W/kg .
2. Artifacts: MR Image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the implant. Tests following guidelines in ASTM F2119-24 showed artifact widths of $\leq 20\text{ mm}$ from the surface of the implant for the gradient echo sequence and $\leq 10\text{ mm}$ from the surface of the implant for the spin echo sequence.

Removable restorations should be taken out prior to scanning, as is done for watches, jewellery etc.

Should there be no MR symbol on the product label, please note that this device has not been evaluated for safety and compatibility in the MR environment.

Summary of Safety and Clinical Performance (SSCP)

As required by the European Medical Device Regulation (MDR; EU2017/745), a Summary of Safety and Clinical Performance (SSCP) is available for perusal with regard to Southern Implants® product ranges.

The relevant SSCP can be accessed at <https://ec.europa.eu/tools/eudamed>.

NOTE: the above website will be available upon the launch of the European Database on Medical Devices (EUDAMED).

Disclaimer of liability


This product is part of the Southern Implants® product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants® product range and take full responsibility for the correct indications and use of this product. Southern Implants® does not assume liability for damage due to incorrect use. Please note that some Southern Implants® products may not be cleared or released for sale in all markets.

Basic UDI

Product	Basic-UDI Number
Basic-UDI for Zygomatic Implants	60095440387194

Related literature and catalogues

CAT-2070 - Zygomatic Implants Product Catalogue

Symbols and warnings

* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

All rights reserved. Southern Implants®, the Southern Implants® logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants®. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale. It is the responsibility of the clinician to inspect the symbols that appear on the packaging of the product in use.

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