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
South Africa - Headquarters: 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

Subsidiaries

Australia

Southern Implants Australia
T: +61-(0)-8-9466-2627
E: info@southernimplants.com.au

Spain and Portugal

Southern Implants Iberica
T: +34 935 053 507
E: info@southernimplants.es

United Kingdom and Ireland

Southern Implants UK
T: +44-20-8899-6845 / 6 / 7
E: info@southernimplants.co.uk

USA and Canada

Southern Implants North America Inc.
T: +1-561-472-0990
E: customercare@southernimplants.com

Description

Southern Implants insertion tools are reusable instruments which are used in conjunction with Southern Implant's implant connections. These handpiece and wrench tools are available in various lengths. All implant insertion tools are provided non-sterile.

The handpiece versions have a latch dimension compatible to ISO 1797. This is in order to connect the insertion tools to the handpiece of an implant motor unit. Handpiece tools are available in:

- the standard version, which has a W&H hex on the shaft, allowing the handpiece to drive on the W&H hex as well as the latch.
- a universal version, which does not feature a W&H hex on the shaft, the handpiece will drive on the latch only.
- The wrench versions have a square connection, to be used with the I-WI-SS wrench converter connected to the Southern Implants torque wrench (I-TWS-B45/B100).

Intended use

The implant insertion tools are intended to aid in the treatment of partially or fully edentulous patients eligible for placement of one or more dental implants as a means of fixing a permanent or removable single crown, partial or full-arch dental prosthesis in the upper or lower jaw.

The implant insertion tools are intended to hold a dental implant (e.g., orthodontic implant screw) at its distal end during the application of the implant to the patient's oral cavity. It is typically a small metal shaft designed to be mounted onto the distal end of a manual screwdriver or a handpiece latch (ISO 1797 or W&H) whereby rotation is applied to introduce the implant into a surgical cavity.

The insertion tools are classified as medical devices and are multiple use devices.

Indications for use

The indications for use for the device are not dissimilar to that of dental implants.

Southern Implants dental implants are intended for both one and two-stage surgical procedures in the following situations and with the following clinical protocols:

- replacing single and multiple missing teeth in the mandible and maxilla.
- immediate placement in extraction sites and in situations with a partially or completely healed alveolar ridge.
- immediate loading in all indications, except in single tooth situations on implants shorter than 8 mm or in soft bone (type IV) where implant stability may be difficult to obtain, and immediate loading may not be appropriate.

Southern Implants (Pty) Ltd's manufactures and sells various reusable instruments that are intended to facilitate the insertion of restorations of its and dental implants by clinicians.

Intended user

Maxillo-facial Surgeons, General Dentists, Orthodontists, Periodontists, Prosthodontists, and other appropriately trained and experienced implant users.

Intended environment

The devices are intended to be used in a clinical environment such as an operating theatre or a dentist consultation room.

Intended patient population

The intended patient population for the device is not dissimilar to that of dental implant therapy.

The intended patient population for implant therapy is partially or fully edentulous patients requiring prosthetic dental restoration in the upper or lower jaw. Restorations may comprise single teeth, partial or full bridges and are fixed restorations.

Additionally, the intended population for some device variants includes those who have previously undergone dental implant therapy.

Compatibility information

Southern Implants' implants should be restored with Southern Implants' components. In the Southern Implants' range there are 8 implant connections. The product code and connection type can be identified by specific abbreviations in the product codes. Range identifiers are summarised in Table A. Refer to the individual product catalogues for details on implant connection sizes and replacement parts compatible with each listed insertion tool.

Table A – Compatible insertion tools for various implant ranges

Implant connection type	Insertion tool type	Compatible device
External Hex (EX)	W&H Hexed latch handpiece	Parts labelled I-CON-X and I-CON-XS (uses replacement part I-X-RING)
	ISO 1797 latch handpiece	Parts labelled I-CONU-X and I-CONU-XS (uses replacement part I-X-RING)
	Handheld hexed connection	Parts labelled I-IMP-INS-2 and I-ZYG-INS-2 (uses replacement part I-X-RING)
	Handheld square connection	Parts labelled I-IMP-INS-1 and I-ZYG-INS-1
External Hex Internal Drive (EX i)	W&H Hexed latch handpiece	Parts labelled I-HID-S/ M/ L (uses replacement part I-PBIT-H)
	ISO 1797 latch handpiece	Parts labelled I-HIDU-S/ M/ L (uses replacement part I-PBIT-H)
	Wrench bit	Parts labelled I-WIID-S/ L (uses replacement part I-PBIT-H)
Deep Conical (DC)	W&H Hexed latch handpiece	Parts labelled I-HDC3-GS/ M, I-HDC4-GS/ M and I-HDC5-GS/ M
	Wrench bit	Parts labelled I-WI-DC3-S/ M, I-WI-DC4-S/ M and I-WI-DC5-S/ M
	W&H Hexed latch handpiece connected to fixture mount	Parts labelled I-CON-X and I-CON-XS (uses replacement part I-X-RING)
	ISO 1797 latch handpiece connected to fixture mount	Parts labelled I-CONU-X and I-CONU-XS (uses replacement part I-X-RING)
TRI-NEX® (EL) (Lobe)	W&H Hexed latch handpiece	Parts labelled I-HLH-35GS/ 35M, I-HLH-43GS/ 43M, I-HLH-50GS/ 50M, I-HLH-60GS/ 60M, I-L43-12D/ -GS and I-L50-12D/ -GS (uses replacement part I-PBIT-L18 or I-PBIT-L2)
	ISO 1797 latch handpiece	Parts labelled I-HLHU-35S/ 35M, I-HLHU-43S/ 43M, I-HLHU-50S/ 50M and I-HLHU-60S/ 60M (uses replacement part I-PBIT-L18 or I-PBIT-L2)
	Wrench bit	Parts labelled I-WI-LH35S/ LH35M, I-WI-LH43S/ LH43M, I-WI-LH50S/ LH50M, I-WI-LH60S, I-WIL-43-12D and I-WIL50-12D (uses replacement part I-PBIT-L18 or I-PBIT-L2)
Internal Hex M-Series and PROVATA® (3M/ M/ Z)	W&H Hexed latch handpiece	Parts labelled I-H3M-M/ L, I-H-PRO-GS, I-HM-S/ M/ L and I-HZ-M/ L (uses replacement part I-PBIT-H16 or I-PBIT-L18)
	Wrench bit	Parts labelled I-WI-3M-S/ M/ L, I-WI-M-S/ M/ L, I-H-PRO12D-GS and I-WI-Z-S/ M/ L (uses replacement part I-PBIT-H16 or I-PBIT-L18)
Internal Octagon IT (ITS/ ITS6)	W&H Hexed latch handpiece	Parts labelled I-HITS and I-HITS-S (uses replacement part I-PBIT-H)
	ISO 1797 latch handpiece	Parts labelled I-HITSU and I-HITSU-S (uses replacement part I-PBIT-H)
	Wrench bit	Parts labelled I-WI-ITS and I-WI-ITS-S (uses replacement part I-PBIT-H)
	W&H Hexed latch handpiece connected to fixture mount	Parts labelled I-CON-IT

Single Platform (SP)	W&H Hexed latch handpiece	Parts labelled I-HSP-GSN/ GS and I-HSP-S/ M
	Wrench bit	Parts labelled I-WI-SP-S/ M
ILZ Mini Implants	W&H Hexed latch handpiece	Parts labelled I-HILZ-S/ M/ L
	Wrench bit	Parts labelled I-WI-ILZ-S
	Handheld PEEK cap	Part supplied assembled with the implant

Clinical benefits

The clinical benefits of Drills and Handpiece Devices are not dissimilar to that of dental implant therapy.

The clinical benefits of dental implant therapy include improved chewing function, speech, aesthetics and patient psychological wellbeing. Through this procedure patients can expect to have their missing teeth replaced and/or crowns restored.

The clinical benefits of implant insertion tools include the successful achievement of implant placement

Storage, cleaning and sterilisation

Limitations on reusable items

A direct value for reusable instruments cannot be given. Frequent processing may have minor effects on the instruments. The product life is normally determined by wear and damage during use, thus instruments if properly cared for and inspected after each use, can be reused many times. Maintain a checklist for these instruments recording the number of uses.

Prior to reprocessing the device it should be thoroughly inspected and tested to determine its suitability for reuse.

NOTE: during use, handle drills and instruments with sterile tweezers to minimise contamination of the instrument tray and risk of damage to sterile surgical gloves.

Containment

As soon as practically possible, remove all visible residue after use (bone, blood or tissue), by immersing the instrument in cold water (dried soil may be difficult to remove).

Precleaning

Disassemble instruments from hand pieces and all connecting parts from instruments in order to clean soil from obstructed areas. Remove the PEEK bits from placement tools. Rinse with lukewarm water for 3 minutes and remove hardened debris with a soft nylon brush. Avoid mechanical damage during cleaning.

Manual cleaning or automated cleaning

Prepare an ultrasonic bath with suitable detergent (i.e. Steritech instrument cleaner - 5% dilution), sonicate for 20 minutes (alternative methods can be used if proven by the end user). Rinse with purified/sterile water.

NOTE: always follow the instructions for use of the manufacturers for cleaning agents and disinfectants.

Load devices into a thermidisinfector. Run the cleaning and disinfection cycle, followed by the drying cycle.

Drying

Dry the instruments both inside and outside with filtered compressed air or single use, lint free wipes. Pack the instruments as quickly as possible after removal into the storage container. If additional drying is necessary, dry in a clean location.

Inspection

Do a visual inspection of the items to check for damages on instruments.

Packaging

Use the correct packaging material as indicated for steam sterilisation to ensure sterility is maintained. Double packaging is recommended. Where appropriate, the cleaned, disinfected and checked devices can be assembled and placed in instrument trays as applicable. Instrument trays can be double wrapped or placed in sterilisation bags.

Sterilisation

Southern Implants® recommends one of the following procedures to sterilise the restoration prior to use:

1. prevacuum sterilisation method: steam sterilise the abutments at 132°C (270°F) at 180 - 220 kPa for 4 minutes. Dry for at least 20 minutes in the chamber. Only an approved wrap or pouch for steam sterilisation must be used.
2. for users in the USA: prevacuum sterilisation method: wrapped, steam sterilise at 135°C (275°F) at 180 - 220 kPa for 3 minutes. Dry for 20 minutes in the chamber. Use a wrap or pouch that is cleared for the indicated steam sterilisation cycle.

NOTE: users in the USA must ensure that the steriliser, wrap or pouch, and all steriliser accessories are cleared by the FDA, for the intended sterilisation cycle.

Storage

Maintain packaging integrity to ensure sterility in storage. Packaging should be completely dry before storage to avoid corrosion and degradation of cutting edges.

Contraindications

Do not use in patients:

- who are medically unfit for dental implant procedures.
- where adequate numbers of implants could not be placed to achieve full functional support of the prosthesis.
- who are allergic or have hypersensitivity to pure titanium or titanium alloy (Ti-6Al-4V), gold, palladium, platinum or iridium.
- who are under the age of 18, have poor bone quality, blood disorders, infected implant site, vascular impairment, uncontrolled diabetes, drug or alcohol abuse, chronic high dose steroid therapy, anti-coagulant therapy, metabolic bone disease, radiotherapy treatment and sinus pathology.

Warnings and precautions

THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING.

- For the safe and effective use of dental implants it is strongly suggested that specialised training be undertaken, including hands-on training to learn proper technique, biomechanical requirements and radiographic evaluations.
- Products must be secured against aspiration when handled intraorally. Aspiration of products may lead to infection or unplanned physical injury.

Responsibility for proper patient selection, adequate training, experience in the placement of implants and providing appropriate information for informed consent rests with the practitioner. Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone. Implant failure increases when implants are placed in irradiated bone as radiotherapy can result in progressive fibrosis of vessels and soft tissue, leading to diminished healing capacity.

It is important to be aware and avoid damage to vital structures such as nerves, veins and arteries. Injuries to vital anatomical structures may cause serious complications such as injury to the eye, nerve damage and excessive bleeding. It is essential to protect the infraorbital nerve. Failing to identify actual measurements relative to the radiographic data could lead to complications.

New and experienced implant users should do training before using a new system or attempting to do a new treatment method. Take special care when treating patients who have local or systemic factors that could affect the healing of

the bone and soft tissue (i.e. poor oral hygiene, uncontrolled diabetes, are on steroid therapy, smokers, infection in the nearby bone and patients who had orofacial radiotherapy).

Thorough screening of prospective implant candidates must be performed including:

- a comprehensive medical and dental history.
- visual and radiological inspection to determine adequate bone dimensions, anatomical landmarks, occlusal conditions and periodontal health.
- bruxism and unfavourable jaw relations must be taken into account.
- proper preoperative planning with a good team approach between well trained surgeons, restorative dentists and lab technicians is essential for successful implant treatment.
- minimising the trauma to the host tissue increases the potential for successful osseointegration.
- electrosurgery should not be attempted around metal implants as they are conductive.

Should the device not operate as intended, it must be reported to the manufacturer of the device. The contact information for the manufacturer of this device to report a change in performance is: sicomplaints@southernimplants.com.

Side effects

The side effects of the use of the system are not dissimilar to those of dental implant therapy. Possible side effects to implant therapy include:

- pain
- swelling
- phonetic difficulties
- gingival inflammation

Less common but more persistent symptoms include, but are not limited to:

- allergic reaction(s) to implant and/or abutment material
- breakage of the implant and/or abutment
- loosening of the abutment screw and/or retaining screw
- infection requiring revision of the dental implant
- nerve damage resulting in permanent weakness, numbness, or pain
- histologic responses with possible macrophage and/or fibroblast involvement
- fat emboli formation
- loosening of the implant requiring revision surgery
- perforation of the maxillary sinus
- perforation of the labial and lingual plates
- bone loss possibly resulting in revision or removal of the implant.

Precaution: maintaining sterility protocol

Implants are packaged as follows:

1. An outer package consisting of a rigid, clear box which acts as protection for the inner package.
2. The inner package consisting of a blister pack (clear plastic-formed blister base with a TYVEK "peel-back" lid).
3. Within the inner package, there is a hollow tube which contains one implant suspended from a titanium ring, this ensures the implant never touches the inside of the plastic tube.
4. Labelling information is located on the surface of the peel-back lid and on the outside of the rigid box.

Care must be taken to maintain the sterility of the implant by proper opening of the packaging and handling of the implant.

1. Open the implant package in the non-sterile field, with non-sterile gloves, tear the address label to open the box.
2. With non-sterile gloves, remove the inner blister pack. Do not place the plastic box or blister pack-lid onto the sterile field. The contents of this inner package are sterile.

- The sealed blister is to be opened by an assistant (with nonsterile gloves): remove the TYVEK lid and drop or place the sterile tube onto the sterile field, open the tube cap and attach the implant placement tool onto the implant and carefully remove from the sterile tube. Do not touch the sterile implant.

Other sterile components are packed in a peel pouch or blister base with a “peel-back” lid. Labelling information is located on the bottom half of the pouch, inside the packet or on the surface of the peel-back lid. Sterility is assured unless the pouch is damaged or opened. Non-sterile components are supplied clean but not sterile in a peel pouch or blister base with peelback lid. Labelling information is located on the bottom half of the pouch or on the surface of the peel-back lid.

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established.

The contact information for the manufacturer of this device to report a serious incident is as follows: sicomplaints@southernimplants.com.

Materials

Material type Titanium Alloy (Ti-6Al-4V), Stainless Steel (DIN 1.4197 & SAE 316), Aluminium, Medical Grade PEEK

Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

Summary of Safety and Clinical Performance (SSCP)

As required by the European Medical Device Regulation (MDR; EU2017/745), a Summary of Safety and Clinical Performance (SSCP) is available for perusal with regard to Southern Implants® product ranges.

The relevant SSCP can be accessed at <https://ec.europa.eu/tools/eudamed>.

NOTE: the above website will be available upon the launch of the European Database on Medical Devices (EUDAMED).

Disclaimer of liability

This product is part of the Southern Implants® product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants® product range and take full responsibility for the correct indications and use of this product. Southern Implants® does not assume liability for damage due to incorrect use. Please note that some Southern Implants® products may not be cleared or released for sale in all markets.

Basic UDI

Product	Basic-UDI Number
Basic-UDI for Reusable Instruments	6009544038769E
Basic-UDI for Drills and Handpiece Devices	6009544038759C

Related literature and catalogues

CAT-2004 - Tri-Nex Implants Product Catalogue

CAT-2005 - IT Implants Product Catalogue

CAT-2010 - Osseointegrated Fixtures Catalogue

CAT-2020 - External Hex Implants Product Catalogue
 CAT-2042 - Deep Conical Implants Product Catalogue
 CAT-2043 - Internal Hex Implants Product Catalogue
 CAT-2060 - PROVATA® Implants Product Catalogue
 CAT-2069 - INVERTA® Implants Product Catalogue
 CAT-2070 - Zygomatic Implants Product Catalogue
 CAT-2087M - ILZ Mini Implants Product Catalogue
 CAT-2093 - Single Platform (SP1) Implants Product Catalogue

Symbols and warnings

											
Manufacturer: Southern Implants® 1 Albert Rd, P.O Box 605 IRENE, 0062, South Africa. Tel: +27 12 667 1046	CE mark 2797	Prescription device*	Sterilised using irradiation	Non-sterile	Use by date (mm-yy)	Do not reuse	Do not resterilise	Catalogue number	Batch code	Medical device	Authorised representative in the European Community
											
Authorised representative for Switzerland	Date of manufacture	Magnetic Resonance conditional	Magnetic Resonance safe	Single sterile barrier system with protective packaging inside	Single sterile barrier system	Consult instruction for use	Caution	Keep away from sunlight	Do not use if package is damaged		

* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

All rights reserved. Southern Implants®, the Southern Implants® logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants®. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale. It is the responsibility of the clinician to inspect the symbols that appear on the packaging of the product in use.

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