
<b>English</b>	<b>INSTRUCTIONS FOR USE: Southern Implants® Prosthetic Screws</b>
<b>Español</b>	<b>INSTRUCCIONES DE USO: Tornillos protéticos Southern Implants®</b>
<b>Italiano</b>	<b>ISTRUZIONI PER L'USO: Southern Implants® Viti protesiche</b>
<b>Français</b>	<b>MODE D'EMPLOI : Vis prothétiques Southern Implants®</b>
<b>Deutsch</b>	<b>GEBRAUCHSANWEISUNGEN: Southern Implants® Prothetiksrauben</b>
<b>Português</b>	<b>INSTRUÇÕES DE UTILIZAÇÃO: Southern Implants® Parafusos protéticos</b>

**South Africa - Headquarters:** 1 Albert Road, Irene, 0062, RSA  
T: +27-12-667-1046 | E: [info@southernimplants.com](mailto:info@southernimplants.com)

**Southern Implants Europe AB:** Holmgatan 30, S-791 71 Falun, Sweden  
T: +46 23 13300 | E: [ecrep@southernimplants.com](mailto:ecrep@southernimplants.com)

#### Subsidiaries

##### Australia

Southern Implants Australia  
T: +61-(0)-8-9466-2627  
E: [info@southernimplants.com.au](mailto:info@southernimplants.com.au)

##### Spain and Portugal

Southern Implants Iberica  
T: +34 935 053 507  
E: [info@southernimplants.es](mailto:info@southernimplants.es)

##### United Kingdom and Ireland

Southern Implants UK  
T: +44-20-8899-6845 / 6 / 7  
E: [info@southernimplants.co.uk](mailto:info@southernimplants.co.uk)

##### USA and Canada

Southern Implants North America Inc.  
T: +1-561-472-0990  
E: [customer care@southernimplants.com](mailto:customer care@southernimplants.com)

**Intended use**

Southern Implants® prosthetic retaining screws are intended to be used in the maxilla or mandible to connect an abutment or framework to an endosseous implant or abutment to restore chewing function for the patient.

**Description**

These are pre-manufactured screws that can connect an abutment direct to an endosseous implant, or a restoration on compact conical level. The Southern Implants prosthetic screws consist of gold and titanium screws. These components are supplied sterile and for single patient use. Screws are available in hexed, slotted, unigrip and a torx connection. Ensure the correct screw and screwdriver is used with one of Southern Implants 5 implant connections.

**Indications for use**

Southern Implants prosthetic screws are intended for both one and two-stage surgical procedures in the following situations and with the following clinical protocols:

- replacing single and multiple missing teeth in the mandible and maxilla.

- immediate placement in extraction sites and in situations with a partially or completely healed alveolar ridge.
- immediate loading in all indications, except in single tooth situations on implants shorter than 8mm or soft bone (type IV) where implant stability may be difficult to obtain and immediate loading may not be appropriate.

**Intended user**

Maxillo-facial Surgeons, General Dentists, Orthodontists, Periodontists, Prosthodontists and other appropriately trained and experienced implant users.

**Intended environment**

The devices are intended to be used in any standard autoclave as well as in the clinical environment such as an operating theatre or a dental consultation room.

**Intended patient population**

Patients that have lost one tooth or multiple teeth.

Table 1 - Torque Table for Southern Screws

External Hex	1.22 Hex (M2)		QUAD (M2)		SLOTTED (M2)					UNIGRIP (M2)				TORQUE: 30-40 Ncm Head Diameter: 2.70 mm		
	2 Series															
		TSH2	TSHZ2	BSH2	GSQ2	GSQZ2	TSS2	TSSZ2	GSS2	GSSZ2	BSS2*	TSU2	TSUZ2		GSU2	GSUZ2
3 Series																
	TSH3	TSHZ3	BSH3*	GSQ3	GSQZ3	TSS3	TSSZ3	GSS3	GSSZ3	BSS3*	TSU3	TSUZ3	GSU3	GSUZ3		
NOTE: screw TORQUE with PEEK Prosthetics for <math>\lt; \varnothing 4.0 \text{ mm}</math> implant interfaces: 15 Ncm >= <math>\varnothing 4.0 \text{ mm}</math> implant interfaces: 20 Ncm																
INPI	1.22 Hex (M1.6)		TORQUE: 32 Ncm Head Diameter: 2.40 mm		Piccolo	1.22 Hex (M1.6)		TORQUE: 25-32 Ncm Head Diameter: 2.15 mm								
			TSH16	BSH16*				TS-P-16	BS-P-16*							
NOTE: screw TORQUE with PEEK Prosthetics: 15 Ncm																

DC (Deep Conical)	1.22 Hex			1.22 Hex Passive Abutment Screws								
	M1.4	M1.6	M2	M1.4	M1.4	PA-DC4-16T	PA-DC4-16B*	PA-DC4-16T	PA-DC4-16B*	PA-DC5-20T	PA-DC5-20B*	
	TS-DC3-14	BS-DC3-14*	TS-DC4-16	BS-DC4-16*	TS-DC5-20	BS-DC5-20*	PA-DC3-14T	PA-DC3-14B* (anodized yellow)	PA-DC4-16T	PA-DC4-16B* (anodized blue)	PA-DC5-20T	PA-DC5-20B* (anodized purple)
	DC3 TORQUE: 15 Ncm Head Diameter: 2.10 mm		DC4 TORQUE: 20 Ncm Head Diameter: 2.35 mm		DC5 TORQUE: 25-35 Ncm Head Diameter: 2.35 mm		TORQUE: 15 Ncm Head Diameter: 2.20 mm		TORQUE: 20 Ncm Head Diameter: 2.40 mm		TORQUE: 25-32 Ncm Head Diameter: 2.80 mm	
NOTE: screw TORQUE with PEEK Prosthetics: DC3 interface: 15 Ncm DC4 interface: 20 Ncm DC5 interface: 20 Ncm												




Internal Hex (M-series and PROVATA®)	1.27 Hex (M1.6)	1.27 Hex (M1.8)	Quad (M1.8)	TORQUE: 32 Ncm Head Diameter: 2.25 mm	1.27 Hex (M1.6)	1.27 Hex Passive Abutment Screws (M1.8)			TORQUE: 32 Ncm Head Diameter: 2.60 mm
	TS-Z-16 BS-Z-16*	TS-Z-18 BS-Z-18*	GS-Q-18		PA-M-16T PA-M-16B	PA-M-18T PA-Q-18G PA-M-18B*			
NOTE: screw TORQUE with PEEK Prosthetics: 15 Ncm									

IT (Internal Octagon)	Torx (M2)			Torx Passive Abutment Screws (M2)			TORQUE: 32-40 Ncm Head Diameter: 2.70 mm
	TSIT2 GSIT2 BSIT2*	TS-IT-PA BS-IT-PA*					

TRI-NEX	Unigrip (M1.8)				TORQUE: 32 Ncm Head Diameter: 2.50 mm	Unigrip (M2)				TORQUE: 32-40 Ncm Head Diameter: 2.50 mm	
	TS-L-18 TS-L-18C (coated) GS-L-18 BS-L-18*	TS-L-20 TS-L-20C (coated) GS-L-20 BS-L-20*									
	NOTE: screw TORQUE with PEEK Prosthetics: 15 Ncm										
	Unigrip Passive Abutment Screws (M1.8)					TORQUE: 32 Ncm Head Diameter: 2.70 mm	Unigrip Passive Abutment Screws (M2)				
PA-L-18 PA-L-18G PA-L-18B*	PA-L-20 PA-L-20G PA-L-20B*										

Compact Conical Abutment Level Screws	Series 1 Screws (M1.4)			Slotted			Unigrip		TORQUE: 10-15 Ncm Head Diameter: 2.25 mm
	1.22 Hex	TSH1 GSH1 BSH1*	TSS1 GSS1 BSS1*	TSU1 GSU1					
	NOTE: screw TORQUE with PEEK Prosthetics: 10 - 15 Ncm								

Angled Compact Conical Abutment Screws	External Hex	Series 9 and 9Z Screws (M2)					TORQUE: 20 Ncm Head Diameter: 2.50 mm	Internal Hex (M-series and PROVATA®)	1.22 Hex		TORQUE: 32 Ncm Head Diameter: 2.20 mm
		Slotted	Unigrip	GSS9 TSU9 GSU9 TSUZ9 GSUZ9	TS-Z-MC-16 TS-Z-MC-18						
DC (Deep Conical)	TRI-NEX	1.22 Hex		Unigrip		TORQUE: 20 Ncm Head Diameter: 2.25 mm	Unigrip		TORQUE: 20 Ncm Head Diameter: 2.45 mm		
		TS-DC3-MC TS-DC4-MC	TU-MCL-18C (coated) GU-MCL-18	TU-MCL-20C (coated) GU-MCL-20							

<b>Procera System</b>	<b>Unigrip (M2)</b>		TORQUE: 32-40 Ncm Head Diameter: 2.50 mm	<b>Screw Connections</b>	<b>Hexed</b>	<b>Unigrip</b>	<b>Quad</b>	<b>Slotted</b>	<b>Torx</b>	G.. = Gold Alloy T.. = Gr. 5 Titanium Alloy B.. = Brass
	 TSP2	 TSP2C (coated)								

**NOTE:** • always ensure that the correct screw is used for the relevant implant and component.  
 • due to design revisions and changes, screw tips may be flat or rounded.  
 \*Blacked Brass Screws for Laboratory Use Only - Torque values do not apply.

The 1.22/1.27 Universal Hex drivers are available in hand-held, hand-piece and wrench insert versions. The advantage is that the same driver can be used on both 1.22 and 1.27 Hex driving interfaces (this does not include 0.9 Hex interfaces). The Universal Hex driver will allow users to easily swap between External Hex / DC (Deep Conical) and PROVATA® (Internal Hex) ranges. All laser markings indicate "22U", which stands for 1.22 Hex and the tool being Universal.

Product codes:

- hand-held: I-HD-22U-S/M/L
- hand-piece: I-HHD-22U-S/M/L
- wrench insert: I-WI-22U-S/M/L

Colour coding:

- hand-held: blue idler.
- hand-piece: blue band on tool shaft
- wrench insert: blue band on square feature of the wrench insert.

These tools are also laser marked with the tool size (i.e. S/M/L) - on the idler of the hand-held drivers and on the shaft of the hand-piece and wrench insert drivers to indicate shaft length.

Below provides information on which driver is used with Southern Implants' variety of implant systems.

**Table 2 - Southern Drivers**

		EXTERNAL HEX						DC (Deep Conical)		
		Hexed		Unigrip	Quad	Slotted	Hexed			
		COVER SCREW DRIVERS	HEALING ABUTMENT AND PROSTHETIC SCREW DRIVERS			PROSTHETIC SCREW DRIVERS			HEALING ABUTMENT AND PROSTHETIC SCREW DRIVERS	
<b>Handheld</b>		I-CS-HD/L 	I-HD-S/M/L 	I-HD-22U-S/M/L 	I-UGI-S/M/L 	I-QDI-S/M/L 	I-BD-S/M/L 	I-HD-S/M/L 	I-HD-22U-S/M/L 	
		I-HHD-09 	I-HHD-22S/M/L 	I-HHD-22U-S/M/L 	I-HUG-S/M/L 	I-HQD-S/M/L 	I-HBDS/M/L 	I-HHD-22S/M/L 	I-HHD-22U-S/M/L 	
		I-WI-09 	I-WI-22S/M/L 	I-WI-22U-S/M/L 	I-WI-UG-S/M/L 	I-WI-QS/M/L 	I-WI-BS/M/L 	I-WI-22S/M/L 	I-WI-22U-S/M/L 	

	Internal Hex (PROVATA® / M-Series)			TRI-NEX	IT (Internal Octagon)
	Hexed		Quad	Unigrip	Torx
	COVER SCREW, HEALING ABUTMENT AND PROSTHETIC SCREW DRIVERS		GOLD PROSTHETIC SCREW DRIVERS	COVER SCREW, HEALING ABUTMENT AND PROSTHETIC SCREW DRIVERS	COVER SCREW, HEALING ABUTMENT AND PROSTHETIC SCREW DRIVERS
Handheld					
Handpiece					
Wrench insert					

**Clinical procedure**

A proper clinical and radiological evaluation must be done to determine the bone dimensions and bone quality.

**Clinical benefits**

Patients can expect to have their missing teeth replaced and/or crowns restored. Screwdrivers are used in dental procedures or in dental implant crowns and bridges.

**Surgical procedure**

1. Select the appropriate screw for the implant or abutment based on the product catalogue.
2. Using the appropriate driver (Table 2), insert the retaining screw into the abutment or framework and connect the assembly to the implant or compact conical abutment.
3. Tighten the screw with an appropriate driver to the torque specified in Table 1.

**CAUTION:** do not exceed the recommended insertion torques for retaining screws. Overtightening may result in screw fracture.

**Before surgery**

All components, instruments and tooling used during the clinical or laboratory procedure must be maintained in good condition and care must be taken that instrumentation does not damage implants or other components.

**During surgery**

Care must be taken that parts are not swallowed or aspirated during any of the procedures, a rubber-dam application is recommended when appropriate. Care must be taken to apply the correct tightening torque of abutments and abutment screws.

**Post-surgery**

Regular patient follow-up and proper oral hygiene must be achieved to ensure favourable long-term results.

**Storage, cleaning and sterilisation**

The implants, cover screws, prosthetic screws and healing abutments are supplied sterile (sterilised by gamma irradiation) and intended for

single-use prior to the expiration date (see packaging label). Sterility is assured unless the container or seal is damaged or opened. If packaging is damaged do not use the product and contact your Southern representative or return to Southern Implants. Do not reuse implants, cover screws, temporary abutments, abutments and prosthetic screws. Re-using these components may result in:

- damage on the surface or critical dimensions, which may result in performance and compatibility degradation.
- adds the risk of cross-patient infection and contamination if single use items are reused.

Southern Implants does not accept any responsibility for complications associated with reused components

**Contraindications**

Do not use in patients:

- who are medically unfit for dental implant procedures.
- where adequate numbers of implants could not be placed to achieve full functional support of the prosthesis.
- who are allergic or have hypersensitivity to pure titanium or titanium alloy (Ti-6Al-4V), gold, palladium, platinum or iridium.
- who are under the age of 18, have poor bone quality, blood disorders, infected implant site, vascular impairment, uncontrolled diabetes, drug or alcohol abuse, chronic high dose steroid therapy, anti-coagulant therapy, metabolic bone disease, radiotherapy treatment and sinus pathology.

**Warnings**

THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING.

- For the safe and effective use of dental implants it is suggested that specialised training be undertaken, including hands-on training to learn proper technique, biomechanical requirements and radiographic evaluations.
- Responsibility for proper patient selection, adequate training, experience in the placement of implants and providing appropriate information for informed consent rests with the practitioner. Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone.

**Cautions**

New and experienced implant users should do training before using a new system or attempting to do a new treatment method. Take special care when treating patients who have local or systemic factors that could affect the healing of the bone and soft tissue (i.e. poor oral hygiene, uncontrolled diabetes, are on steroid therapy, smokers, infection in the nearby bone and patients who had oro-facial radiotherapy). Thorough screening of prospective implant candidates must be performed including:

- a comprehensive medical and dental history.
- visual and radiological inspection to determine adequate bone dimensions, anatomical landmarks, occlusal conditions and periodontal health.
- bruxism and unfavourable jaw relations must be taken into account.
- proper pre-operative planning with a good team approach between well trained surgeons, restorative dentists and lab technicians is essential for successful implant treatment.
- minimizing the trauma to the host tissue increases the potential for successful osseointegration.
- electro-surgery should not be attempted around metal implants, as they are conductive.

**Single use devices**

Do not reuse devices indicated for single use (use the device prior to the expiration date).

**Side effects**

Potential side effects and temporary symptoms: pain, swelling, phonetic difficulties, gingival inflammation. More persistent symptoms: the risks and complications with implants include, but are not limited to: (1) allergic reaction(s) to implant and/or abutment material; (2) breakage of the implant and/or abutment; (3) loosening of the abutment screw and/or retaining screw; (4) infection requiring revision of the dental implant; (5) nerve damage that could cause permanent weakness, numbness, or pain; (6) histologic responses possibly involving macrophages and/or fibroblasts; (7) formation of fat emboli; (8) loosening of the implant requiring revision surgery; (9) perforation of the maxillary sinus; (10) perforation of the labial and lingual plates and (11) bone loss possibly resulting in revision or removal.

**Breakage**

Implant and abutment fractures can occur when applied loads exceed the normal functional torque strength of the material. Potential overloading conditions may result from: deficiencies in implant numbers, lengths and/or diameters to adequately support a restoration, excessive cantilever length, incomplete abutment seating, abutment angles greater than 30 degrees, occlusal interferences causing excessive lateral forces, patient parafunction (e.g. bruxing, clenching), loss or changes in dentition or functionality, inadequate prosthesis fit and physical trauma. Additional treatment may be necessary when any of the above conditions are present to reduce the possibility of hardware complications or failure.

**Performance requirements**

It is the responsibility of the clinician to instruct the patient on all appropriate contraindications, side effects and precautions as well as the need to seek the services of a trained dental professional if there are any changes in the performance of the implant (e.g. looseness of the prosthesis, infection or exudate around the implant, pain or any other unusual symptoms that the patient has not been told to expect).

**Notice regarding serious incidents**

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established. The contact information for the manufacturer of this device to report a serious incident is as follows: [sicomplaints@southernimplants.com](mailto:sicomplaints@southernimplants.com)

**Materials**

Titanium retaining screws: Titanium alloy Ti-90%, Al-6%, V-4%  
 Gold retaining screws: Gold Alloy Au-61%, Ag-16.5%, Pt-13.5%, Cu -9%

**Disposal**

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

**Disclaimer of liability**

This product is part of the Southern Implants product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants product range and take full responsibility for the correct indications and use of this product. Southern Implants does not assume liability for damage due to incorrect use. Please note that some Southern Implants products may not be cleared or released for sale in all markets.

**Healing**

The healing time required for osseointegration depends on the individual and treatment protocol. It is the responsibility of the practitioner to decide when the implant can be restored. Good primary stability will govern if immediate loading can be done.

**Implant care and maintenance**

Potential implant patients should establish an adequate oral hygiene regime prior to implant therapy. Proper post-operative, oral hygiene and implant maintenance instructions must be discussed with the patient, as this will determine the longevity and health of the implants. The patient should maintain regular prophylaxis and evaluation appointments. It is the responsibility of the practitioner to decide when the implant can be restored. Good primary stability will govern if immediate loading can be done.

**MR conditional**

Non-clinical testing has demonstrated that the Southern Implants dental implants, metallic abutments and prosthetic screws are MR conditional. A patient with these devices can be safely scanned in a MR system meeting the following conditions:

- Static magnetic field of 1.5 Tesla and 3.0 Tesla only.
- Maximum spatial gradient magnetic field of 4000 Gauss/cm (40 T/m).
- Maximum MR system reported, head specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode) or whole body averaged specific absorption rate (wbSAR) of 1 W/kg.

Under the scan conditions defined above, the Southern Implants dental implants, abutments and prosthetic screws are expected to produce a maximum temperature rise of 5.8°C after 15 minutes of continuous scanning. In the non-clinical testing, the image artifact caused by the device extends approximately 20 mm from the Southern Implants dental implants, abutments and prosthetic screws, when imaged with a gradient echo pulse sequence and a 3.0 Tesla MRI system. Removable restorations should be taken out prior to scanning, as is done for watches, jewellery etc. Should there be no MR symbol on the product label, please note that this device has not been evaluated for safety and compatibility in the MR environment. This device has not been tested for heating or migration in the MR environment.


**Basic UDI**

Product	Basic-UDI Number
Basic-UDI for Metal Abutments	600954403872

**Related literature and catalogues**

- CAT-2004 - Tri-Nex Implants Product Catalogue
- CAT-2005 - IT Implants Product Catalogue
- CAT-2010 - Osseointegrated Fixtures Catalogue
- CAT-2020 - External Hex Implants Product Catalogue
- CAT-2042 - Deep Conical Implants Product Catalogue
- CAT-2043 - Internal Hex Implants Product Catalogue
- CAT-2060 - PROVATA® Implants Product Catalogue
- CAT-2069 - INVERTA® Implants Product Catalogue
- CAT-2070 - Zygomatic Implants Product Catalogue

### Symbols and Warnings

 Manufacturer: Southern Implants 1 Albert Rd, P.O Box 605 IRENE, 0062, South Africa. Tel: +27 12 667 1046		 Prescription device*	 Sterilization using Irradiation	 Consult instruction for use	 Use by date (mm-yy)	 Do not reuse	 Batch code	 Do not use if package is damaged	 Medical Device	 Authorize Representative in the European Community	 Catalog number	 Date of Manufacture
<p>* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.</p>												
						<p>Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.</p>						
<p>All rights reserved. Southern Implants®, the Southern Implants logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale.</p>												

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