English	INSTRUCTIONS FOR USE: Southern Implants® Prosthetic Screws
Español	INSTRUCCIONES DE USO: Tornillos Protésicos Southern Implants®
Italiano	ISTRUZIONI PER L'USO: Southern Implants® Viti Protesiche
Français	NOTICE D'UTILISATION : Vis prothétiques de Southern Implants®
Deutsch	GEBRAUCHSANWEISUNG: Southern Implants® Prothetische Schrauben
Português	INSTRUÇÕES DE UTILIZAÇÃO: Parafusos protéticos Southern Implants®
Ελληνικά	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ: Προσθετικές βίδες της Southern Implants®
Svenska	BRUKSANVISNING: Southern Implants® Protesskruvar
Română	INSTRUCȚIUNI DE UTILIZARE: Suruburi protetice Southern Implants®
Türkçe	KULLANMA TALİMATI: Southern Implants® Protez Vidaları

Southern Implants (Pty.) Ltd. 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

Subsidiaries

Australia
Southern Implants Australia
T: +61-(0)-8-9466-2627
E: info@southernimplants.com.au

Spain and Portugal
Southern Implants Iberica
T: +34 935 053 507
E: info@southernimplants.es

United Kingdom and Ireland
Southern Implants UK
T: +44-20-805-94490
E: info@southernimplants.co.uk

USA and Canada
Southern Implants North America Inc.
T: +1-561-472-0990
E: customer care@southernimplants.com

Description

The Southern Implants® Prosthetic Screws are pre-manufactured screws that can connect an abutment direct to an endosseous implant, or a restoration on compact conical level.

The Southern Implants® Prosthetic Screws consist of gold, titanium, and brass screws. The gold and titanium prosthetic screws are supplied sterile and are intended for single patient use. The brass prosthetic screws are for Laboratory use only, which are used to seat the prosthesis to the laboratory abutment analogue during prosthetic fabrication, and are not intended to come into contact with the patient. The screws are available in hexed, slotted, quad, unigrip and a torx connection. Ensure the correct screw and screwdriver is used with one of Southern Implants 6 implant connections.

Intended use

The Southern Implants® Metal Abutments are prosthetic components intended for connection to endosseous implant systems to aid in prosthetic rehabilitation and implant procedures.

Specifically, the Southern Implants® Prosthetic Screws are components of the endosseous implant system, intended to secure the system's components together to aid in prosthetic rehabilitation and implant procedures.

Indications for use

The Southern Implants® Metal Abutments are indicated for use as components of an endosseous implant system to support prosthetic rehabilitation in implant treatment procedures.

Specifically, the Southern Implants® Prosthetic Screws are indicated for use in the functional and aesthetic treatment of patients eligible for implant-retained prostheses.

Intended user

Dental technicians, Maxillofacial Surgeons, General Dentists, Orthodontists, Periodontists, Prosthodontists, and other appropriately trained and/or experienced medical professionals.

Intended environment

The devices are intended to be used in a dental laboratory as part of the restoration design and manufacture as well as in a clinical environment such as an operating theatre or a dentist consultation room

Intended patient population

Southern Implants® Metal Abutments are intended to be used in patients indicated and eligible, or otherwise not contraindicated, for implant-retained prosthetic restorations.

Compatibility information

Southern Implants' implants should be restored with Southern Implants' components. In the Southern Implants' range there are 6 implant connections. The implant code and connection type can be identified by specific abbreviations in the product codes. Table A summarises the screws available for each Southern range and the recommended torque values of the Southern screws.

Table A – Torque Tables for Southern Screws

Screw Type	External Hex Screws						Head Diameter
	1.22 Hex	Quad	Slotted	Unigrip	Torx	Torque	
2 Series (M2)	TSH2 TSHZ2 BSH2*	GSQ2 GSQZ2	TSS2 TSSZ2 GSS2 GSSZ2 BSS2*	TSU2 TSUZ2 GSU2 GSUZ2	-	32-40 Ncm (NOTE: screw TORQUE with PEEK Prosthetics for < Ø4.0 mm implant interfaces: 15 Ncm ≥ Ø4.0 mm implant interfaces: 20 Ncm)	2.70 mm
3 Series (M2)	TSH3 TSHZ3 BSH3*	GSQ3 GSQZ3	TSS3 TSSZ3 GSS3 GSSZ3 BSS3*	TSU3 TSUZ3 GSU3 GSUZ3	TSTZ5-ASC		2.40 mm

Screw Type	External Hex Screws						
	1.22 Hex	Quad	Slotted	Unigrip	Torx	Torque	Head Diameter
INPI (M1.6)	TSH16 BSH16*	-	-	-	-	32 Ncm	2.40 mm
Piccolo (M1.6)	TS-P-16 BS-P-16*	-	-	-	TS-P-16-ASC (Head Diameter – 2.20 mm)	25-32 Ncm (NOTE: screw TORQUE with PEEK Prosthetics for < Ø4.0 mm implant interfaces: 15 Ncm)	2.15 mm

Screw Type	DC (Deep Conical) Screws			
	1.22 Hex	Torx	Torque	Head Diameter
Standard Screws (M2)	TS-DC5-20 BS-DC5-20*	TS-DC5-ASC	32 Ncm (NOTE: screw TORQUE with PEEK Prosthetics DC5 interface: 20 Ncm)	2.35 mm
Standard Screws (M1.6)	TS-DC4-16 BS-DC4-16*	TS-DC4-ASC (Head Diameter – 2.30 mm)	30 Ncm (NOTE: screw TORQUE with PEEK Prosthetics DC4 interface: 20 Ncm)	2.35 mm
Standard Screws (M1.4)	TS-DC3-14 BS-DC3-14*	TS-DC3-ASC	20 Ncm (NOTE: screw TORQUE with PEEK Prosthetics DC5 interface: 15 Ncm)	2.10 mm
Passive Abutment Screws (M2)	PA-DC5-20T PA-DC5-20B* (anodised purple)	-	32 Ncm	2.80 mm
Passive Abutment Screws (M1.6)	PA-DC4-16T PA-DC4-16B* (anodised blue)	-	30 Ncm	2.40 mm
Passive Abutment Screws (M1.4)	PA-DC3-14T PA-DC3-14B* (anodised yellow)	-	20 Ncm	2.20 mm

Screw Type	Internal Hex (M-series and PROVATA®) Screws				
	1.27 Hex	Quad	Torx	Torque	Head Diameter
Standard Screws (M1.8)	TS-Z-18 BS-Z-18*	-	TS-Z-18-ASC (Head Diameter – 2.30mm)	32-40 Ncm (NOTE: screw TORQUE with PEEK Prosthetics: 15 Ncm)	2.25 mm
	-	GS-Q-18	-	32 Ncm (NOTE: screw TORQUE with PEEK Prosthetics: 15 Ncm)	
Passive Abutment Screws (M1.8)	PA-M-18T PA-M-18B*	PA-Q-18G	-	32 Ncm (NOTE: screw TORQUE with PEEK Prosthetics: 15 Ncm)	2.60 mm
Standard Screws (M1.6)	TS-Z-16 BS-Z-16*	-	TS-Z-16-ASC (Head Diameter – 2.20mm)	32 Ncm (NOTE: screw TORQUE with PEEK Prosthetics: 15 Ncm)	2.25 mm
Passive Abutment Screws (M1.6)	PA-M-16T PA-M-16B*	-	-		2.60 mm

Screw Type	SP1 (Single Platform) Screws			
	1.22 Hex	Torx	Torque	Head Diameter
Standard Screws (M1.6)	TS-SP-16 BS-SP-16*	TS-SP-ASC (Head Diameter – 2.20 mm)	32 Ncm	2.00 mm
Passive Abutment Screws (M1.6)	PA-SP-16T PA-SP-16B*	-		2.40 mm

Screw Type	IT (Internal Octagon) Screws		
	Torx	Torque	Head Diameter
Standard Screws (M2)	TSIT2 GSIT2 BSIT2*	32-40 Ncm	2.70 mm
	TS-IT-TIB TSIT2-ASC	32-40 Ncm	2.40 mm

Screw Type	IT (Internal Octagon) Screws		
	Torx	Torque	Head Diameter
Passive Abutment Screws (M2)	TS-IT-PA BS-IT-PA*	32-40 Ncm	2.70 mm

Screw Type	TRI-NEX Screws			
	Unigrip	Torx	Torque	Head Diameter
Standard Screws (M2)	TS-L-20 TS-L-20C (coated) GS-L-20 BS-L-20*	TS-L-20-ASC	32-40 Ncm (NOTE: screw TORQUE with PEEK Prosthetics: 20 Ncm)	2.50 mm
Standard Screws (M1.8)	TS-L-18 TS-L-18C (coated) GS-L-18 BS-L-18*	TS-L-18-ASC	32 Ncm (NOTE: screw TORQUE with PEEK Prosthetics: 15 Ncm)	2.50 mm
Passive Abutment Screws (M2)	PA-L-20 PA-L-20G PA-L-20B*	-	32-40 Ncm	2.70 mm
Passive Abutment Screws (M1.8)	PA-L-18 PA-L-18G PA-L-18B*	-	32 Ncm	2.70 mm

Screw Type	Compact Conical Abutment Screws					
	1.22 Hex	Slotted	Unigrip	Torx	Torque	Head Diameter
Abutment Level Screws	TSH1-NLT	-	-	-	10-15 Ncm (NOTE: screw TORQUE with PEEK Prosthetics: 15 Ncm)	2.0 mm
	TSH1 GSH1 BSH1*	TSS1 GSS1 BSS1*	TSU1 GSU1	TST1-ASC (Head Diameter – 2.30 mm)		2.25 mm
	TSH1-27	-	-	-		2.70 mm

Screw Type	Angled Compact Conical Abutment Screws				
	1.22 Hex	Slotted	Unigrip	Torque	Head Diameter
Series 9 and 9Z (M2)	-	GSS9	TSU9 GSU9 TSUZ9 GSUZ9	20 Ncm	2.50 mm
TRI-NEX (M2)	-	-	TU-MCL-20C (coated) GU-MCL-20		2.45 mm
External Hex, ZYGIN (M2)	-	-	TS-EX-MC TS-ZYG-MC		2.38 mm
Internal Hex (M-series and PROVATA®) (M1.8)	TS-Z-MC-18	-	-	32 Ncm	2.20 mm
TRI-NEX (M1.8)	-	-	TU-MCL-18C (coated) GU-MCL-18	20 Ncm	2.25 mm
Internal Hex (M-series and PROVATA®) (M1.6)	TS-Z-MC-16	-	-	32 Ncm	2.20 mm
DC (Deep Conical) (M1.6)	TS-DC4-MC	-	-	30 Ncm	2.00 mm
DC (Deep Conical) (M1.4)	TS-DC3-MC	-	-	20 Ncm	2.00 mm
SP1 (Single Platform) (M1.6)	TS-SP-MC	-	-	20 Ncm	2.00 mm

Screw Type	Procera System		
	Unigrip	Torque	Head Diameter
Standard Screws (M2)	TSP2 TSP2C (coated)	32-40 Ncm	2.50 mm

NOTE:

- always ensure that the correct screw is used for the relevant implant and component.
- due to design revisions and changes, screw tips may be flat or rounded.

**Blackened Brass Screws for Laboratory Use Only - Torque values do not apply.*

See CAT-8057 for information on which driver is used with Southern Implants' variety of implant systems. For further information see the relevant Southern Implants® Product Catalogues and the Prosthetic Manual (CAT-4100).

Clinical performance

The Southern Implants® Prosthetic Screws are designed to secure and fix the components of the implant system together. The clinical performance of these screws is primarily characterized by the degree of fixation achieved between the prosthetic device and the supporting structures.

Clinical benefits

The clinical benefits associated with the Southern Implants® Prosthetic Screws are intrinsically linked to the overall success of implant-supported prosthetic treatments and encompass both direct and indirect clinical benefits. The direct clinical benefits of the Southern Implants® Prosthetic Screws arise specifically from their use in the treatment process. These benefits include the enhancement of peri-implant soft tissue health and the aesthetic outcomes of the final restoration. The clinical benefits experienced by patients as a result of the implant-supported prosthetic treatments in which these abutments are utilized represent the indirect clinical benefits of the Southern Implants® Prosthetic Screws. While these benefits are substantial, they are not directly attributable to the devices themselves; rather, they reflect the contributions of the abutments to the overall treatment system. Additionally, successful treatment outcomes are anticipated to result in significant enhancements in patients' psychosocial well-being, self-esteem, and overall quality of life. Overall, the Southern Implants® Prosthetic Screws not only facilitate the mechanical aspects of dental restorations but also contribute to the holistic improvement of patient outcomes and experiences.

Storage, cleaning and sterilisation

The component is supplied sterile (sterilised by gamma irradiation) and intended for single-use prior to the expiration date (see packaging label). Sterility is assured unless the container or seal is damaged or opened. If packaging is damaged do not use the product and contact your Southern representative or return to Southern Implants®. The devices must be stored in a dry place at room temperature and not exposed to direct sunlight. Incorrect storage may influence device characteristics. Do not reuse components indicated for single-use only. Reusing these components may result in:

- damage to the surface or critical dimensions, which may result in performance and compatibility degradation.
- adds the risk of cross-patient infection and contamination if single-use items are reused.

Southern Implants® does not accept any responsibility for complications associated with reused single-use components.

Contraindications

The contraindications of all device groups used as part of the specific treatment or procedure apply. Therefore, the contraindications of the systems/medical devices utilized as part of implant surgery/therapy should be noted and the relevant documents consulted.

The contraindications specific to this device group include:

- patients who are medically unfit for implant or prosthetic related procedures
- patients that present with an allergy or hypersensitivity to the following materials: Titanium, Aluminium, Vanadium, Cobalt Chrome, Cobalt Chrome Molybdenum, Gold, Platinum, Palladium, Iridium, Silver and/or Copper

Warnings and precautions

IMPORTANT NOTICE: THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING

- To ensure the safe and effective use of dental implants, new technologies/systems and the Metal Abutment devices, it is strongly recommended that specialised training be undertaken. This training should include hands-on methods to gain competency on the proper technique, the system's biomechanical requirements and radiographic evaluations required for the specific system.
- Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone.
- Use of the device with incompatible or non-corresponding devices can result in the poor performance or failure of the device.
- When handling devices intraorally, it is imperative that they are adequately secured to prevent aspiration, as aspiration of products may lead to infection or physical injury.
- The use of non-sterile items can lead to secondary infections of the tissue or transfer infectious diseases.
- Failure to adhere to appropriate cleaning, re-sterilization, and storage procedures as outlined in the Instructions for Use can result in device damage, secondary infections, or patient harm.
- Exceeding the number of recommended uses for reusable devices can result in device damage, secondary infection or patient harm.
- The use of blunt drills may cause damage to the bone, potentially compromising osseointegration.

It is crucial to emphasize that training should be undertaken by both new and experienced implant users prior to utilizing a new system or attempting a new treatment method.

Patient Selection and Preoperative Planning

A comprehensive patient selection process and meticulous preoperative planning are essential for successful implant treatment. This process should involve consultation among a multidisciplinary team, including well-trained surgeons, restorative dentists, and laboratory technicians.

Patient screening should include, at a minimum, a thorough medical and dental history, as well as visual and radiological inspections to assess the presence of adequate bone dimensions, the positioning of anatomical landmarks, the presence of unfavourable occlusal conditions, and the periodontal health status of the patient.

For successful implant treatment, it is important to:

1. Minimize trauma to the host tissue, thereby increasing the potential for successful osseointegration.
2. Accurately identify measurements relative to radiographic data, as failure to do so may lead to complications.
3. Be vigilant in avoiding damage to vital anatomical structures, such as nerves, veins, and arteries. Injury to these structures may result in serious complications, including ocular injury, nerve damage, and excessive bleeding.

The responsibility for proper patient selection, adequate training and experience in implant placement, and the provision of appropriate information required for informed consent rests with the practitioner. By combining thorough screening of prospective implant candidates with a practitioner possessing a high level of competence in the use of the system, the potential for complications and severe side effects can be significantly reduced.

High Risk Patients

Special care should be exercised when treating patients with local or systemic risk factors that may adversely affect the healing of bone and soft tissue or otherwise increase the severity of side effects, the risk of complications, and/or the likelihood of implant failure. Such factors include:

- poor oral hygiene
- history of smoking/vaping/tobacco use
- history of periodontal disease
- history of orofacial radiotherapy**
- bruxism and unfavourable jaw relations

- use of chronic medications that may delay healing or increase the risk of complications including, but not limited to, chronic steroid therapy, anti-coagulant therapy, TNF- α blockers, bisphosphonate and cyclosporin

*** The potential for implant failure and other complications increases when implants are placed in irradiated bone, as radiotherapy can lead to progressive fibrosis of blood vessels and soft tissue (i.e., osteoradionecrosis), resulting in diminished healing capacity. Contributing factors to this increased risk include the timing of implant placement in relation to radiation therapy, the proximity of radiation exposure to the implant site, and the radiation dosage at that site.*

Should the device not operate as intended, it must be reported to the manufacturer of the device. The contact information or the manufacturer of this device to report a change in performance is: sicomplaints@southernimplants.com.

Side effects

The clinical outcome of treatment is influenced by multiple factors. The following side effects and residual risks are associated with the device group or the surgical procedures involved in its use. They may necessitate further treatment, revision surgery or additional visits to the relevant medical professional. Furthermore, these side effects and residual risks can vary in both severity and frequency.

Therefore, the full list of known potential undesirable side-effects and residual risks that have been identified in relation to the prosthetic screw device group include:

- Abutment fracture
- Abutment screw fracture or distortion
- Allergic reaction(s) to the abutment material
- Anaesthesia, paraesthesia, hyperesthesia, and hypoesthesia (transient or permanent)
- Bleeding
- Bleeding on probing
- Bruising
- Complications requiring revision surgery
- Difficulty in retrieving the abutment
- Gingival inflammation
- Gingival recession
- Hyperplastic tissue response
- Implant failure due to insufficient levels of osseointegration
- Improper implant positioning resulting in prosthetic compromise
- Infection (acute and/or chronic)
- Localized inflammation
- Loosening of the abutment screw and/or retaining screw
- Loss or damage to adjacent teeth
- Marginal bone loss
- Micromovements and implant instability
- Misfit or improper connection at the implant-abutment interface
- Overloading of the abutment/implant
- Pain, tenderness or discomfort
- Peri-implantitis, peri-mucositis or otherwise poor peri-implant soft tissue health
- Periodontal inflammation
- Phonetic difficulties
- Prosthetic failure
- Soft tissue irritation
- Sub-optimal aesthetic result
- Wound dehiscence or poor healing

Precaution: maintaining sterility protocol

The Southern Implants® Prosthetic Screws are packed in a peel pouch or blister base with a “peel-back” lid. Labelling information is located on the bottom half of the pouch, inside the packet or on the surface of the peel-back lid. Sterility is assured unless the pouch or blister is damaged or opened.

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established.

The contact information for the manufacturer of this device to report a serious incident is as follows:

sicomplaints@southernimplants.com.

Materials

Material type Titanium Alloy (Ti-6AL-4V and ASTM F136), gold alloy, brass (for laboratory screws only)

Titanium alloy (Ti-6Al-4V) according to ASTM F136 and ISO 5832-3

Chemical Components	Aluminium (Al)	Vanadium (V)	Residuals* (Fe, O, C, N, H) in total	Titanium (Ti)
Composition, % (mass/mass)	5.50 – 6.75	3.50 – 4.50	<0.55	Balance

*Where, Fe = Iron, O = Oxygen, C= Carbon, N = Nitrogen, H = Hydrogen

Gold Alloy Screws

Chemical Components	Gold (Au)	Silver (Ag)	Platinum (Pt)	Copper (Cu)	Residuals, combined total (individual)
Composition, % (mass/mass)	60.5 – 61.5	16.0 – 17.0	13.0 – 14.0	8.50 – 9.50	0.4 (≤0.1)

Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

MR safety

Nonclinical testing has demonstrated that the Southern Implants® dental implants, metallic abutments and prosthetic screws are MR conditional.

MRI Safety Information	
 MR Conditional	
A patient with the dental implants, metallic abutments and prosthetic screws from Southern Implants® may be safely scanned under the following conditions. Failure to follow these conditions may result in injury to the patient.	
Device Name	dental implants, metallic abutments and prosthetic screws
Static Magnetic Field Strength (B₀) [T]	1.5T or 3T
Maximum Spatial Field Gradient [T/m] [gauss/cm]	30 T/m (3000 gauss/cm)
Type of RF Excitation	Circularly polarized (CP) (i.e., quadrature driven)
RF Transmit Coil Type	Integrated Whole Body Transmit Coil

RF Receive Coil	Any receive-only RF coil may be used
Operating Mode	Normal Operating Mode
RF Conditions	Normal Operating mode: (Head SAR of 3.2 W/kg, 2 W/kg whole body)
Scan Duration	If the closest part of implant <25cm away from the isocenter in feet-head direction: a continuous scan time of 15 minutes will require a cooling delay of at least 5 minutes. If the closest part of implant ≥ 25cm away from the isocenter in feet-head direction: This is limited to 1 hour of scan duration.
MR Image Artifact	The presence of this implant may produce an image artifact.

Additional information:

1. Local transmit coil: Local transmit coil safety was not evaluated. Local transmit coil is only allowed if the implant is ≥ 25cm away from isocenter and for a wbSAR of up to 2W/kg.
2. Artifacts: MR Image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the implant. Tests following guidelines in ASTM F2119-24 showed artifact widths of ≤20 mm from the surface of the implant for the gradient echo sequence and ≤10 mm from the surface of the implant for the spin echo sequence.

Removable restorations should be taken out prior to scanning, as is done for watches, jewellery etc.

Should there be no MR symbol on the product label, please note that this device has not been evaluated for safety and compatibility in the MR environment.

Summary of Safety and Clinical Performance (SSCP)

As required by the European Medical Device Regulation (MDR; EU2017/745), a Summary of Safety and Clinical Performance (SSCP) is available for perusal with regard to Southern Implants® product ranges.

The relevant SSCP can be accessed at <https://ec.europa.eu/tools/eudamed>.

NOTE: the above website will be available upon the launch of the European Database on Medical Devices (EUDAMED). In the interim, email clinical@southernimplants.com to request a Summary of Safety and Clinical Performance (SSCP) for these devices.

Disclaimer of liability

This product is part of the Southern Implants® product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants® product range and take full responsibility for the correct indications and use of this product. Southern Implants® does not assume liability for damage due to incorrect use. Please note that some Southern Implants® products may not be cleared or released for sale in all markets.


Basic UDI

Product	Basic-UDI Number
Basic-UDI for Metal Abutments	60095440387296
Basic-UDI for Prosthetic Screws (Titanium)	60095440501886
Basic-UDI for Prosthetic Screws (Gold)	60095440501784

Related literature and catalogues

CAT-2004 - Tri-Nex Implants Product Catalogue
 CAT-2005 – Internal Octagon Implants Product Catalogue
 CAT-2010 - Osseointegrated Fixtures Catalogue
 CAT-2020 - External Hex Implants Product Catalogue
 CAT-2042 - Deep Conical Implants Product Catalogue
 CAT-2043 - Internal Hex Implants Product Catalogue
 CAT-2060 - PROVATA® Implants Product Catalogue
 CAT-2069 – Deep Conical INVERTA® Implants Product Catalogue
 CAT-2070 - Zygomax Implants Product Catalogue
 CAT-2093 - Single Platform Implants Product Catalogue
 CAT-2095 – External Hex INVERTA® Implants Product Catalogue
 CAT-4100 - Southern Implants® Prosthetic Manual

Symbols and warnings

* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

All rights reserved. Southern Implants®, the Southern Implants® logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants®. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale. It is the responsibility of the clinician to inspect the symbols that appear on the packaging of the product in use.

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