|                  |                                                                                 |
|------------------|---------------------------------------------------------------------------------|
| <b>English</b>   | <b>INSTRUCTIONS FOR USE: Southern Implants® Reusable Instruments</b>            |
| <b>Español</b>   | <b>INSTRUCCIONES DE USO: Instrumentos reutilizables Southern Implants®</b>      |
| <b>Italiano</b>  | <b>ISTRUZIONI PER L'USO: Southern Implants® Strumenti Riutilizzabili</b>        |
| <b>Français</b>  | <b>MODE D'EMPLOI : Instruments réutilisables Southern Implants®</b>             |
| <b>Deutsch</b>   | <b>GEBRAUCHSANWEISUNGEN: Southern Implants® Wiederverwendbare Instrumente</b>   |
| <b>Português</b> | <b>INSTRUÇÕES DE UTILIZAÇÃO: Southern Implants® Instrumentos reutilizáveis</b>  |
| <b>Ελληνικά</b>  | <b>ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ: Southern Implants® Επαναχρησιμοποιούμενα εργαλεία</b>        |
| <b>Svenska</b>   | <b>ANVISNINGAR FÖR ANVÄNDNING: Southern Implants® Återanvändbara Instrument</b> |
| <b>Română</b>    | <b>INSTRUCȚIUNI DE UTILIZARE: Instrumente reutilizabile Southern Implants®</b>  |

**South Africa - Headquarters:** 1 Albert Road, Irene, 0062, RSA  
T: +27-12-667-1046 | E: [info@southernimplants.com](mailto:info@southernimplants.com)

**Southern Implants Europe AB:** Holmgatan 30, S-791 71 Falun, Sweden  
T: +46 23 13300 | E: [ecrep@southernimplants.com](mailto:ecrep@southernimplants.com)

#### Subsidiaries

##### Australia

Southern Implants Australia  
T: +61-(0)-8-9466-2627  
E: [info@southernimplants.com.au](mailto:info@southernimplants.com.au)

##### Spain and Portugal

Southern Implants Iberica  
T: +34 935 053 507  
E: [info@southernimplants.es](mailto:info@southernimplants.es)

##### United Kingdom and Ireland

Southern Implants UK  
T: +44-20-8899-6845 / 6 / 7  
E: [info@southernimplants.co.uk](mailto:info@southernimplants.co.uk)

##### USA and Canada

Southern Implants North America Inc.  
T: +1-561-472-0990  
E: [customercare@southernimplants.com](mailto:customercare@southernimplants.com)

## Description

### Osteotomy Guide

An osteotomy guide is a non-sterile surgical instrument intended to be attached to a portion of bone (e.g. pelvis, femur, tibia, ulna) during an osteotomy procedure to outline the desired cut-line and mechanically support/guide the cutting device (e.g. oscillating cutting saw).

### Reusable Dental Implant Space/Depth/Angle Gauge

A reusable dental implant space/depth/angle gauge is used during the planning and post-osteotomy stages of a dental implantation procedure, prior to surgical implant placement/fixation, to measure various lengths (e.g. drill depth of the osteotomy cavity, mucosal height, distance between implants/teeth for determining the mesiodistal position of the implant) and/or angles (e.g. direction of the perforation in relation to the antagonist arch), to aid in selection of an appropriate implant and/or final prosthetic abutment to be installed. It is typically graduated and may be a hand-held instrument or a stand-alone device placed intraorally; some types are intended for use during x-ray imaging.

### Dental Osteotomes

A hand-held surgical instrument used for shaping and condensing bone in a dental osteotomy, typically when performing a dental implant surgical procedure. It is held by the user who will use a surgical mallet or hammer to manually impart an impacting force to the proximal end of the instrument. The device has a straight or bayonet-shaped handle leading to a round or concave tip at the distal end, which may be marked with graduations to indicate depth differences. It is used for ridge expanding, bone condensing and sinus floor elevation procedures. This device is made of high-grade stainless steel.

### Basic Dental Retractors

A hand-held, non-self-retaining dental instrument intended to be used to manually displace soft tissues of the oral cavity to improve their visualization and access, and to afford them protection during oral surgical procedures. It is a metallic device and includes a shaped/hooked blade with a handle/holding portion.

### Surgical Drill Guide, Reusable

A hand-held tubular surgical device designed to simultaneously place, angle, and guide a rotating surgical drill on hard tissue (e.g., bone), and to protect the surrounding soft tissue and dentition while the drill is in operation. It can have a handle at its proximal end for holding and manipulation, or it may be held by another device (e.g. orthopaedic jig, 3D-printed surgical guide). The device is dedicated to drill bit size, made of metal, and may have graduations along its length; it is not custom-made for a specific patient.

### Dental Torque Wrenches

A hand-held dental tool intended for the manual application of rotation to tighten/loosen/remove a device (e.g. dental implant, spatial register (localizer/tracker/calibrator)) fitted into or close to the oral cavity, and/or to a surgical instrument, in association with a dental surgical procedure (e.g. dental guided surgery registration). It may have a built-in ratchet function, with forward and reverse direction, and a calibrated scale to which a recommended torque can be set; it is not a screwdriver. It is typically made of metallic materials and may have a standard socket connection to accommodate different tool bits/adaptors.

### Wrench Adaptors

A component of a surgical wrench (typically ratchet type) intended to enable connection of the main body of the wrench to a variety of distal end pieces (bits). It is intended to be used during an invasive surgical procedure, and is not dedicated to dental use only. It may be referred to as a converter.

### Guide Stabilisers

A device used to retain a surgical guide to a placed implant, via an integral screw, in order to secure the surgical guide from shifting position while additional implant osteotomies are prepared and implants placed.

### Plastic Access Maintainers

A part used to maintain the screw chimney and access hole of a prosthetic component during various restoration manufacturing processes.

### Retention Components

A part intended to be assembled with an implant insertion tool or dental implant holding endpiece, to retain a dental implant or implant fixture mount at its distal end during the application of the implant to the patient's oral cavity.

### **Intended use**

Southern Implants' Reusable Surgical Instruments are intended to aid in the treatment of partially or fully edentulous patients eligible for placement of one or more dental implants as a means of fixing a permanent or removable single crown, partial or full-arch dental prosthesis in the upper or lower jaw.

The specific manner in which this intended use is fulfilled is discussed below:

### Osteotomy Guide

The Osteotomy Guide allows for improved preparation of the procedural site by outlining the desired cutline and mechanically supporting or guiding the cutting device.

### Reusable Dental Implant Space/Depth/Angle Gauge

Depth gauges are used to confirm that the depth of the hole drilled in bone for implant insertion is the correct depth.

Direction indicators are used to show the angle of the implant / implant prosthetic interface after a hole has been drilled in bone.

Profile gauges are used to during planning and post-osteotomy stages of implantation, prior implantation, to select the appropriate implants / abutments to use.

### Dental Osteotomes

The dental osteotome is a hand-held surgical instrument used to shape and condense the bone during implant site preparation allowing for the successful achievement of ridge expansion, bone condensing and sinus floor lift procedures.

### Basic Dental Retractors

The Basic Dental Retractors are used to manually displace the soft tissue of the oral cavity to improve visualisation and access to the implantation site.

### Surgical Drill Guide, Reusable

Drill guides are used to stabilise the drill shafts to position the drill accurately in the desired or planned position.

### Dental Torque Wrenches

The Dental torque wrench enables the manual application of torque to instruments.

### Wrench Adaptors

The wrench adapter enables the compatibility of the Dental Torque Wrench with various distal endpieces.

### Guide Stabilisers

The guide stabilisers aid in securing a surgical guide to prevent unintentional misalignment or repositioning during dental implant placement and surgical procedures.

### Plastic Access Maintainers

Plastic Access maintainers aid in maintaining the screw chimney and access hole of a prosthetic component during various restoration manufacturing processes.

### Retention Components

Retention Components aid in retaining a dental implant or implant fixture mount at its distal end during the application of the implant to the patient's oral cavity.

### Indications for use

The Reusable Instruments are indicated for use in cases of partial or full edentulousness, congenital tooth loss, or failing/problematic natural dentition, when a patient is dissatisfied with existing removable dentures, or as part of the revision/replacement of previous dental restorations.

### Intended user

The intended user for this system includes Maxillofacial Surgeons, General Dentists, Orthodontists, Periodontist, Prosthodontists and other appropriately trained and experienced implant users.

### Intended environment

This device is intended to be used in a clinical environment such as an operating theatre or a dentist consultation room.

### Intended patient population

The intended patient population for dental implant therapy includes partially or fully edentulous patients eligible or, otherwise, not contraindicated for implant placement that require prosthetic restoration or revision of existing restorations in the upper or lower jaw; where the planned restorations are to be fixed and comprise of single teeth, partial or full bridges

### Compatibility information

The reusable instrument can be identified by specific abbreviations in the product codes. Device identifiers are summarised in Table A.

**Table A – Reusable instrument device codes**

| Reusable Instrument Type                         | Device Product Code                                                                                           |
|--------------------------------------------------|---------------------------------------------------------------------------------------------------------------|
| Osteotomy Guides                                 | Parts labelled CG-F-RT-30/ 40/ 50 and I-FH-RT are osteotomy guides for Rohner Technique systems               |
| Dental Implant Space/Depth/Angle Gauge, Reusable | Part labelled I-DG are depth gauges for general dental implants                                               |
|                                                  | Part labelled I-ADG and I-DI are angle gauges for general dental implants                                     |
|                                                  | Parts labelled CH-I-DG and I-ZYG-DG-1 are depth gauges for zygomatic implants                                 |
|                                                  | Parts labelled I-ZYG-DI55 and ZYG-TR-55-35/ 45/ 52.5/ 52 are angle gauges for zygomatic implants              |
|                                                  | Parts labelled MAX-(Ø)-PG are angle gauges for MAX dental implants                                            |
|                                                  | Parts labelled LISA are space analyser gauges to measure the available space between teeth                    |
|                                                  | Parts labelled I-OST-(Ø)C/ P/ SU, I-FSD and I-HAMMER are for general dental implants                          |
| Dental Osteotomes                                | Parts labelled I-NIC/ -L are for general dental implants                                                      |
|                                                  | Parts labelled I-ZYG-RET are for zygomatic implants                                                           |
|                                                  | Parts labelled I-PT-RET are for pterygoid implants                                                            |
| Surgical Drill Guide, Reusable                   | Parts labelled DGS-L-(Ø) and DGS-S-(Ø) are drill guide spoons for general dental implants                     |
|                                                  | Parts labelled GSDS-(Ø) and I-TDS-(Ø) are drill stabilisers for general dental implants                       |
|                                                  | Parts labelled TDF-RT-(Ø) and ICF-RT are surgical guides for Rohner Technique systems                         |
| Dental Torque Wrenches                           | Parts labelled I-RATCHET-2, I-TW/ TWS/ TWA are torque wrenches for general dental implants                    |
| Wrench Adaptors                                  | Parts labelled I-WI-C-L/ S, I-WI-CNB/ CST and I-WI-SH/ SL/ SS are wrench adaptors for general dental implants |
| Guide Stabilisers                                | Parts labelled I-GSN-SP and I-GS-SP for Single Platform (SP1) implants                                        |
| Plastic Access Maintainers                       | Parts labelled PM22 are for general dental implants                                                           |
| Retention Components                             | Parts labelled I-X-RING are for use with the I-CON-X                                                          |
|                                                  | Parts labelled I-PBIT-L18/L2/2/H16/H are for implant insertion tools and dental implant holding endpieces     |

## Clinical benefits

The Reusable Instruments are used as part of dental implant therapy, a procedure in which a patient's missing teeth and/or crowns are restored. As a result, the benefits of the Reusable Instruments mirror those of dental implant therapy in general, including improved chewing function, speech, aesthetics, and patient psychological wellbeing.

### Osteotomy Guide

The Osteotomy Guide allows for improved preparation of the procedural site by outlining the desired cutline and mechanically supporting or guiding the cutting device.

### Reusable Dental Implant Space /Depth/Angle Gauge

The Reusable Dental Implant Space/Depth/Angle Gauge allows for improved selection of implants and/or the final prosthetic abutment-based geometry for the specific osteotomy.

### Dental Osteotomes

The Dental Osteotomes allow for the successful achievement of ridge expansion, bone condensing and sinus floor lift procedures.

### Basic Dental Retractors

The Basic Dental Retractors provide improved protection of the soft tissue of the mouth and surrounding areas during oral surgical procedures.

### Surgical Drill Guide, Reusable

Surgical Drill Guides provide improved protection of the surrounding soft tissue and dentition while a surgical drill is in operation.

### Dental Torque Wrenches

The Dental Torque Wrench allows for the successful achievement of tightening, loosening, or removal of a surgical instrument and/or a device (e.g., dental implant) fitted into or close to the oral cavity.

### Wrench Adaptors

The wrench adapter enables the Dental Torque Wrench to be used with various distal end pieces.

### Guide Stabilisers

The Guide Stabilisers provide improved alignment of the surgical guide and improved protection of the surrounding soft tissue and dentition while a surgical drill is in operation.

### Plastic Access Maintainers

Plastic Access maintainers improve prosthetic outcomes during various restoration manufacturing processes.

### Retention Components

Retention Components improve implant placement into the patient's oral cavity.

## Storage, cleaning and sterilisation

### Limitations on reusable items

A direct value for reusable instruments cannot be given. Frequent processing may have minor effects on the instruments. The product life is normally determined by wear and damage during use, thus instruments if properly cared for and inspected after each use, can be reused many times. Maintain a checklist for these instruments recording the number of uses.

Prior to reprocessing the device it should be thoroughly inspected and tested to determine its suitability for reuse.

**NOTE:** during use, handle drills and instruments with sterile tweezers to minimise contamination of the instrument tray and risk of damage to sterile surgical gloves.

**Containment**

As soon as practically possible, remove all visible residue after use (bone, blood or tissue), by immersing the instrument in cold water (dried soil may be difficult to remove).

**Precleaning**

Disassemble instruments from hand pieces and all connecting parts from instruments in order to clean soil from obstructed areas. Remove the PEEK bits from placement tools. Rinse with lukewarm water for 3 minutes and remove hardened debris with a soft nylon brush. Avoid mechanical damage during cleaning.

**Manual cleaning or automated cleaning**

Prepare an ultrasonic bath with suitable detergent (i.e. Steritech instrument cleaner - 5% dilution), sonicate for 20 minutes (alternative methods can be used if proven by the end user). Rinse with purified/sterile water.

**NOTE:** always follow the instructions for use of the manufacturers for cleaning agents and disinfectants.

Load devices into a thermoisinfectant. Run the cleaning and disinfection cycle, followed by the drying cycle.

**Drying**

Dry the instruments both inside and outside with filtered compressed air or single use, lint free wipes. Pack the instruments as quickly as possible after removal into the storage container. If additional drying is necessary, dry in a clean location.

**Inspection**

Do a visual inspection of the items to check for damages on instruments.

**Packaging**

Use the correct packaging material as indicated for steam sterilisation to ensure sterility is maintained. Double packaging is recommended. Where appropriate, the cleaned, disinfected and checked devices can be assembled and placed in instrument trays as applicable. Instrument trays can be double wrapped or placed in sterilisation bags.

**Sterilisation**

Southern Implants® recommends one of the following procedures to sterilise the device prior to use:

1. prevacuum sterilisation method: steam sterilise the device at 132°C (270°F) at 180 - 220 kPa for 4 minutes. Dry for at least 20 minutes in the chamber. Only an approved wrap or pouch for steam sterilisation must be used.
2. for users in the USA: prevacuum sterilisation method: wrapped, steam sterilise at 135°C (275°F) at 180 - 220 kPa for 3 minutes. Dry for 20 minutes in the chamber. Use a wrap or pouch that is cleared for the indicated steam sterilisation cycle.

**NOTE:** users in the USA must ensure that the steriliser, wrap or pouch, and all steriliser accessories are cleared by the FDA, for the intended sterilisation cycle.

**Storage**

Maintain packaging integrity to ensure sterility in storage. Packaging should be completely dry before storage to avoid corrosion and degradation of cutting edges.

**Contraindications**

Since the Reusable Instruments are required or utilized as part of dental implant surgery, the contraindications are not dissimilar to those of dental implant therapy as a whole. These contraindications include:

- patients medically unfit for oral surgical procedures
- cases where an adequate number of implants cannot be placed to provide full functional support of the prosthesis.
- patients under the age of 18 years.

- patients with poor bone quality.
- patients with blood disorders.
- presence of infection at the implant site.
- patients with vascular impairment.
- patients with uncontrolled diabetes.
- patients undergoing chronic high dose steroid therapy.
- patients with metabolic bone disease.
- patients with pure titanium, titanium alloy (Ti6Al4V), gold, palladium, polyether ether ketone (PEEK) or iridium allergies or hypersensitivities.

Other than the above, there are no contraindications unique to this system.

### Warnings and precautions

THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING

- For the safe and effective use of dental implants, new technologies/systems and the Reusable Instruments, it is strongly suggested that specialised training be undertaken. This training should include hands-on methods to gain competency on the proper technique, biomechanical requirements and radiographic evaluations required for the specific system.
- Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone.
- Use of the device with incompatible or non-corresponding devices can result in device failure or poor performance.
- When handled intraorally, devices must be adequately secured against aspiration. Aspiration of products may lead to infection or physical injury.
- Electrosurgery should not be attempted around metal implants as they are conductive.
- The use of non-sterile items can lead to secondary infections of the tissue or transfer infectious diseases.
- Failure to perform appropriate cleaning, re-sterilisation and storage procedures as per the Instructions for Use document can result in device damage or secondary infection.
- Exceeding the number of recommended uses for reusable devices can result in device damage, secondary infection or patient harm,
- Blunt drills may cause damage to the bone which could compromise osseointegration.
- For short implants, clinicians should closely monitor patients for any of the following conditions: peri implant bone loss, changes to implant's response to percussion, or radiographic changes in bone to implant contact along the implant's length. If the implant shows mobility or greater than 50% bone loss, the implant should be evaluated for possible removal. If the clinicians choose a short implant, then clinicians should consider a two-stage surgical approach, splinting a short implant to an additional implant, and placement of the widest possible fixture. Allow longer periods for osseointegration and avoid immediate loading.

**Note:** Training should be undertaken by both new and experienced implant users before using a new system or attempting to do a new treatment method.

It is important that thorough screening of prospective implant candidates is performed. This screening should include at least:

- a comprehensive medical and dental history.
- visual and radiological inspection to determine adequate bone dimensions, anatomical landmarks, occlusal conditions and periodontal health. Failure to accurately estimate lengths of drills relative to radiographic measurements can result in permanent injury to nerves or other vital structures.
- the evaluation and consideration of bruxism and unfavourable jaw relations.

Proper preoperative planning with a good team approach between well trained surgeons, restorative dentists and lab technicians is essential for successful implant treatment.

For successful implant treatment, it is important to:

- minimise the trauma to the host tissue as this increases the potential for successful osseointegration.

- be aware and avoid damage to vital anatomical structures such as nerves, veins and arteries. Injury to these structures may cause serious complications, including, injury to the eye, nerve damage and excessive bleeding. It is essential to protect the infraorbital nerve.
- identify actual measurements relative to the radiographic data as failing to do so could lead to complications.

Special care should be taken when treating patients with local or systemic risk factors that could affect the healing of the bone and soft tissue or otherwise increase the severity of the side effects, risk of complications and/or implant failure. These factors include:

- smoking.
- history of periodontal disease.
- poor oral hygiene.
- history of orofacial radiotherapy.\*\*
- bruxism and unfavourable jaw relations.

*\*\* The potential for implant failure and other complications increases when implants are placed in irradiated bone as radiotherapy can result in progressive fibrosis of vessels and soft tissue (i.e., osteoradionecrosis), leading to diminished healing capacity. Implant treatment of irradiated patients is dependent upon issues like the timing of implant placement in relation to the radiation therapy, anatomic site chosen for implant placement and radiation dosage at that site and consequent risk of osteoradionecrosis.*

The responsibility for proper patient selection, adequate training and experience in implant placement and the provision of the appropriate information required for informed consent rests with the practitioner.

Should the device not operate as intended, it must be reported to the manufacturer of the device. The contact information or the manufacturer of this device to report a change in performance is: [sicomplaints@southernimplants.com](mailto:sicomplaints@southernimplants.com)

### Side effects

Since the device is required or utilized as part of dental implant surgery, the expected side effects are not dissimilar to those of dental implant therapy as a whole. Common side effects of dental implant therapy include pain, local inflammation, phonetic difficulties and gingival inflammation. Other, less common side effects or complications to dental implant therapy include, but are not limited to: (1) wound dehiscence; (2) peri-implantitis; (3) transient weakness, numbness, and/or pain associated with mild nerve injury; (4) fat emboli formation; (5) marginal bone loss within acceptable limits; (6) allergic reaction(s) to implant and/or abutment material; (7) unspecified infection; (8) implant failure due to insufficient levels of osseointegration; (9) breakage of the implant and/or abutment; (10) loosening of the abutment screw and/or retaining screw; (11) complications requiring revision of the dental implant; (12) nerve damage resulting in permanent weakness, numbness, or pain; (13) loosening of the implant requiring revision surgery; (14) perforation of the maxillary sinus; (15) perforation of the labial and/or lingual plates and (16) bone loss possibly resulting in revision or removal of the implant.

### Precaution: maintaining sterility protocol

Implants are packaged as follows:

1. An outer package consisting of a rigid, clear box which acts as protection for the inner package.
2. The inner package consisting of a blister pack (clear plastic-formed blister base with a TYVEK "peel-back" lid).
3. Within the inner package, there is a hollow tube which contains one implant suspended from a titanium ring, this ensures the implant never touches the inside of the plastic tube.
4. Labelling information is located on the surface of the peel-back lid and on the outside of the rigid box.

Care must be taken to maintain the sterility of the implant by proper opening of the packaging and handling of the implant.

1. Open the implant package in the non-sterile field, with non-sterile gloves, tear the address label to open the box.
2. With non-sterile gloves, remove the inner blister pack. Do not place the plastic box or blister pack-lid onto the sterile field. The contents of this inner package are sterile.

- The sealed blister is to be opened by an assistant (with nonsterile gloves): remove the TYVEK lid and drop or place the sterile tube onto the sterile field, open the tube cap and attach the implant placement tool onto the implant and carefully remove from the sterile tube. Do not touch the sterile implant.

Other sterile components are packed in a peel pouch or blister base with a “peel-back” lid. Labelling information is located on the bottom half of the pouch, inside the packet or on the surface of the peel-back lid. Sterility is assured unless the pouch is damaged or opened. Non-sterile components are supplied clean but not sterile in a peel pouch or blister base with peel-back lid. Labelling information is located on the bottom half of the pouch or on the surface of the peel-back lid.

### Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established.

The contact information for the manufacturer of this device to report a serious incident is as follows: [sicomplaints@southernimplants.com](mailto:sicomplaints@southernimplants.com).

### Materials

Material type                      Stainless Steel (DIN 1.4197, SAE 316/ 316L/ 303/ 431/ 304), Stainless Steel Spring Wire (AISI 302 or DIN 1.4320), Commercially pure titanium (grade 2, 3, 4), Titanium Alloy (Ti-6Al-4V), Celeron/Tufnol, Delrin, Silicon, Radel PPSU, Medical Grade PEEK, Aluminium,

### Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

### Summary of Safety and Clinical Performance (SSCP)

As required by the European Medical Device Regulation (MDR; EU2017/745), a Summary of Safety and Clinical Performance (SSCP) is available for perusal with regard to Southern Implants® product ranges.

The relevant SSCP can be accessed at <https://ec.europa.eu/tools/eudamed>.

**NOTE:** the above website will be available upon the launch of the European Database on Medical Devices (EUDAMED).

### Disclaimer of liability

This product is part of the Southern Implants® product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants® product range and take full responsibility for the correct indications and use of this product. Southern Implants® does not assume liability for damage due to incorrect use. Please note that some Southern Implants® products may not be cleared or released for sale in all markets.

### Basic UDI

| Product                            | Basic-UDI Number |
|------------------------------------|------------------|
| Basic-UDI for Reusable Instruments | 6009544038769E   |

### Related literature and catalogues

CAT-2004 - Tri-Nex® Implants Product Catalogue  
 CAT-2005 - IT Implants Product Catalogue  
 CAT-2010 - Osseointegrated Fixtures Catalogue  
 CAT-2020 - External Hex Implants Product Catalogue  
 CAT-2042 - Deep Conical Implants Product Catalogue  
 CAT-2043 - Internal Hex Implants Product Catalogue  
 CAT-2060 - PROVATA® Implants Product Catalogue  
 CAT-2069 - INVERTA® Implants Product Catalogue  
 CAT-2070 - Zygomatic Implants Product Catalogue  
 CAT-2093 - Single Platform (SP1) Implants Product Catalogue

### Symbols and warnings

|                                                                                                                                |                                                                                     |                                                                                     |                                                                                     |                                                                                     |                                                                                     |                                                                                       |                                                                                       |                                                                                       |                                                                                       |                                                                                       |                                                                                       |
|--------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
|                                             |  |  |  |  |  |    |   |  |  |  |  |
| Manufacturer:<br>Southern Implants®<br>1 Albert Rd, P.O<br>Box 605 IRENE,<br>0062,<br>South Africa.<br>Tel: +27 12 667<br>1046 | CE mark<br>2797                                                                     | Prescription<br>device*                                                             | Sterilised<br>using<br>irradiation                                                  | Non-sterile                                                                         | Use by date<br>(mm-yy)                                                              | Do not reuse                                                                          | Do not<br>resterilise                                                                 | Catalogue<br>number                                                                   | Batch code                                                                            | Medical device                                                                        | Authorised<br>representative in the<br>European<br>Community                          |
|                                             |  |  |  |  |  |  |  |  |  |                                                                                       |                                                                                       |
| Authorised<br>representative<br>for Switzerland                                                                                | Date of<br>manufacture                                                              | Magnetic<br>Resonance<br>conditional                                                | Magnetic<br>Resonance<br>safe                                                       | Single sterile barrier system<br>with protective packaging<br>inside                | Single sterile<br>barrier system                                                    | Consult<br>instruction<br>for use                                                     | Caution                                                                               | Keep away from<br>sunlight                                                            | Do not use if<br>package is<br>damaged                                                |                                                                                       |                                                                                       |

\* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

All rights reserved. Southern Implants®, the Southern Implants® logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants®. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale. It is the responsibility of the clinician to inspect the symbols that appear on the packaging of the product in use.

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