English	INSTRUCTIONS FOR USE: Southern Implants® Reusable Instruments
Español	INSTRUCCIONES DE USO: Southern Implants® Reusable Instruments
Italiano	ISTRUZIONI PER L'USO: Southern Implants® Strumenti Riutilizzabili
Français	MODE D'EMPLOI : Instruments réutilisables de Southern Implants®
Deutsch	GEBRAUCHSANWEISUNG: Southern Implants® Wiederverwendbare Instrumente
Português	INSTRUÇÕES DE UTILIZAÇÃO: Southern Implants® Instrumentos reutilizáveis
Ελληνικά	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ: Southern Implants® Επαναχρησιμοποιούμενα εργαλεία
Svenska	BRUKSANVISNING: Southern Implants® Återanvändbara Instrument
Română	INSTRUCȚIUNI DE UTILIZARE: Instrumente reutilizabile Southern Implants®
Türkçe	KULLANMA TALİMATI: Southern Implants® Tekrar Kullanılabilir Aletler

South Africa - Headquarters: 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

EC REP

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

Subsidiaries

Australia

Southern Implants Australia
T: +61-(0)-8-9466-2627
E: info@southernimplants.com.au

Spain and Portugal

Southern Implants Iberica
T: +34 935 053 507
E: info@southernimplants.es

United Kingdom and Ireland

Southern Implants UK
T: +44-20 8059 4490
E: info@southernimplants.co.uk

USA and Canada

Southern Implants North America Inc.
T: +1-561-472-0990
E: customercare@southernimplants.com

Description

Cutting Guide

A cutting guide is a non-sterile surgical instrument intended to be attached to a portion of bone (e.g. pelvis, femur, tibia, ulna) during an osteotomy procedure to outline the desired cut-line and mechanically support/guide the cutting device (e.g. oscillating cutting saw).

Fibula Holder

The Fibula Holder has pointed clamps so as to securely hold the fibula for the surgeon to cut the bone. The fibula shall be elevated between 80mm to 100mm from the surface.

Implant Space Analyser Gauge

A reusable dental implant space/depth/angle gauge is used during the planning and post-osteotomy stages of a dental implantation procedure, prior to surgical implant placement/fixation, to measure various lengths (e.g. drill depth of the osteotomy cavity, mucosal height, distance between implants/teeth for determining the mesiodistal position of the implant) and/or angles (e.g. direction of the perforation in relation to the antagonist arch), to aid in selection of an appropriate implant and/or final prosthetic abutment to be installed. It is typically graduated and may be a hand-held instrument or a stand-alone device placed intraorally; some types are intended for use during x-ray imaging.

Angle Gauge

A reusable dental implant space/depth/angle gauge is used during the planning and post-osteotomy stages of a dental implantation procedure, prior to surgical implant placement/fixation, to measure various lengths (e.g. drill depth of the osteotomy cavity, mucosal height, distance between implants/teeth for determining the mesiodistal position of the implant) and/or angles (e.g. direction of the perforation in relation to the antagonist arch), to aid in selection of an appropriate implant and/or final prosthetic abutment to be installed. It is typically graduated and may be a hand-held instrument or a stand-alone device placed intraorally; some types are intended for use during x-ray imaging.

Depth Gauge

A reusable dental implant space/depth/angle gauge is used during the planning and post-osteotomy stages of a dental implantation procedure, prior to surgical implant placement/fixation, to measure various lengths (e.g. drill depth of the osteotomy cavity, mucosal height, distance between implants/teeth for determining the mesiodistal position of the implant) and/or angles (e.g. direction of the perforation in relation to the antagonist arch), to aid in selection of an appropriate implant and/or final prosthetic abutment to be installed. It is typically graduated and may be a hand-held instrument or a stand-alone device placed intraorally; some types are intended for use during x-ray imaging.

Reusable Bone Taps

A metal surgical instrument designed for cutting threads into bone so that the threads facilitate the insertion of bone screws without excessive torque. The screws anchor bone fragments, fixtures and/or other devices to the bone. It is made of titanium and has a helical thread intersected with longitudinal cutting grooves forming the cutting edges at its distal end and may be cannulated. The proximal end typically terminates in a square to which a tapping instrument or extension can be attached. It will typically be manually rotated in order to create the best possible thread.

Surgical Screwdrivers

The screwdrivers are intended to be used to apply torque to a prosthetic screw or abutment to affix the screw or abutment to the dental implant, or loosen and remove the component from the dental implant or abutment. The screwdrivers have a driving interface on one end that is compatible with the corresponding screw or abutment. This interface allows for pickup of the corresponding component and prevents relative rotation between the screw or abutment and the screwdriver when torque is applied. The screwdrivers are available in various lengths (short 5mm, medium 14mm, long 21mm) to allow for use in both the anterior and posterior where available space may differ.

Dental Implant Holding Endpiece

A device intended to hold a dental implant (e.g., orthodontic implant screw) at its distal end during the implantation of the implant to the patient's oral cavity. It is typically a small metal shaft designed to be mounted onto the distal end of a manual screwdriver whereby rotation is applied to introduce the implant into a surgical cavity.

Dental Osteotomes

A hand-held surgical instrument used for shaping and condensing bone in a dental osteotomy, typically when performing a dental implant surgical procedure. It is held by the user who will use a surgical mallet or hammer to manually impart an impacting force to the proximal end of the instrument. The device has a straight or bayonet-shaped handle leading to a round or concave tip at the distal end, which may be marked with graduations to indicate depth differences. It is used for ridge expanding, bone condensing and sinus floor elevation procedures. This device is made of high-grade stainless steel.

Bone Mills

A rotary surgical instrument designed to enlarge the diameter of the proximal portion of a hole drilled in bone. It is a shaft of metal with a cutting head, available in various configurations. The device is typically used to recess the head of a dental implant within bone (enable it to lie flush with or below the bone surface) or to help accommodate the implantation of a device in bone. Accessories (such as caps) may be required to connect to the implants, which align the countersinks and protect the implant's internal connection against damage from the rotary motion of the countersink.

Basic Dental Retractors

A hand-held, non-self-retaining dental instrument intended to be used to manually displace soft tissues of the oral cavity to improve their visualization and access, and to afford them protection during oral surgical procedures. It is a metallic device and includes a shaped/hooked blade with a handle/holding portion.

Screw Extractor Guides

A hand-held manual surgical tool designed to remove implanted screws that have a damaged screw recess by self-tapping into the defective recess/head and providing sufficient grip to rotate and remove the screw in the normal way. It typically has a conical tip that is placed into the implant connection to provide a channel through which the screw remover is guided during use. The device is also intended to protect the implant internal connection during the screw removal procedure.

Surgical Drill Guides

A surgical device designed to assist during surgical procedures to act as a physical guide for the position, depth and/or direction of the rotating surgical drill. It may also serve to protect the surrounding soft tissue and dentition while the drill is in operation. It can have a handle at its proximal end for holding and manipulation, or it may be held by another device (e.g. orthopaedic jig, 3D-printed surgical guide). The devices are dedicated to drill bit size, made of metal, and may have graduations along its length; it is not custom-made for a specific patient.

Dental Torque Wrenches

A hand-held dental tool intended for the manual application of rotation to tighten/loosen/remove a device (e.g. dental implant, spatial register (localizer/tracker/calibrator)) fitted into or close to the oral cavity, and/or to a surgical instrument, in association with a dental surgical procedure (e.g. dental guided surgery registration). It may have a built-in ratchet function, with forward and reverse direction, and a calibrated scale to which a recommended torque can be set; it is not a screwdriver. It is typically made of metallic materials and may have a standard socket connection to accommodate different tool bits/adaptors.

Fixture Mount Spanners

A Fixture Mount Spanner is intended to hold an implant in position while the fixture mount screw is loosened to prevent reverse torque from being transmitted to the implant.

Wrench Adaptors

A component of a surgical wrench (typically ratchet type) intended to enable connection of the main body of the wrench to a variety of distal end pieces (bits). It is intended to be used during an invasive surgical procedure, and is not dedicated to dental use only. It may be referred to as a converter.

Retention Components

A device intended to be assembled with an implant insertion tool, dental implant holding endpiece or drill extension, to retain a dental implant, implant fixture mount or drill during implant placement procedures.

Guide Stabilisers

A device used to retain a surgical guide to a placed implant, via an integral screw, in order to secure the surgical guide from shifting position while additional implant osteotomies are prepared and implants placed.

Space Maintainers

A device generally used on top of a compact conical abutment after implant placement to maintain soft tissue opening during preparation of the prosthesis.

Intended use

The Southern Implants® Reusable Instruments are reusable devices intended for use in the planning, osteotomy preparation and placement of implants/prosthetics where they aid in the surgical procedures associated with implant therapy.

Specifically, Southern Implants® Cutting Guide is intended to act as a physical guide during the bone cutting procedures associated with grafting.

The Southern Implants® Fibula Holder is intended to clamp a fibula graft and enable it to be cut to size for use in the grafting procedures.

The Southern Implants® Implant Space Analyser Gauge is intended to be used as a visual guide for the space that is available between dentition.

The Southern Implants® Angle Gauge is intended to be used as a visual guide to confirm or determine the direction or orientation of the prepared osteotomy.

The Southern Implants® Depth Gauge is intended to be used as a visual guide to confirm or determine the depth of the prepared osteotomy.

The Southern Implants® Reusable Bone Tap is intended to cut threads into the bone, facilitate the insertion of the target implant(s) and/or reduce the torque required for implant placement.

The Southern Implants® Surgical Screw Drivers are intended to apply or transmit torque to abutments and/or screws.

The Southern Implants® Dental Implant Holding Endpiece is intended to apply or transmit torque to an implant during implantation.

The Southern Implants® Dental Osteotomes are intended to shape and condense bone as part of ridge expansion and/or sinus floor elevation techniques.

The Southern Implants® Bone Mills are intended to remove or mill bone from around an implant platform.

The Southern Implants® Basic Dental Retractors are intended to displace the soft tissue of the oral cavity and improve visualization/access during implant bed preparation or implant placement.

The Southern Implants® Screw Extractor Guides are intended to be used as a physical guide for the Screw Remover during screw removal procedures.

The Southern Implants® Surgical Drill Guide is intended to act as a physical guide for the position, depth and direction of the rotating surgical drill.

The Southern Implants® Dental Torque Wrench is intended to apply or transmit torque to screws, instruments, implants or prosthetic devices to facilitate their tightening/loosening/removal.

The Southern Implants® Fixture Mount Spanner is intended to prevent reverse torque from being transmitted to the implant when the fixture mount screw is loosened.

The Southern Implants® Wrench Adaptor is intended to enable the connection of a torque wrench to a variety of instruments.

The Southern Implants® Retention Components are intended to be assembled with an implant insertion tool, dental implant holding endpiece or drill extension to retain a dental implant, implant fixture mount or drill during implant placement procedures.

The Southern Implants® Guide Stabilizer is intended to retain and secure the Surgical Drill Guide.

The Southern Implants® Space Maintainers are intended to maintain soft tissue opening after implant placement during preparation of the prosthesis.

Indications for use

The Southern Implants® Reusable Instruments are indicated for the procedures associated with the functional and aesthetic restoration of individuals indicated for implant-retained prostheses and the placement/use of Southern Implants® devices.

Specifically, the Southern Implants® Screw Extractor Guides are indicated for the removal of damaged or defective screws.

The Southern Implants® Surgical Drill Guides are indicated for guided surgery protocols.

The Southern Implants® Guide Stabilizer is indicated for guided surgery protocols.

Intended user

Dental Technicians, Maxillofacial Surgeons, General Dentists, Orthodontists, Periodontists, Prosthodontists, and other appropriately trained and/or experienced medical professionals.

Intended environment

The devices are intended to be used in a dental laboratory as part of the restoration design and manufacture as well as in a clinical environment such as an operating theatre or a dentist consultation room.

Intended patient population

The Southern Implants® Reusable Instruments are intended to be used in patients subject to implant treatment/revision.

Compatibility information

Table A – Reusable Instrument Device Codes

Reusable Instrument Type	Device Product Code
Cutting Guides and Fibula Holder	Parts labelled CG-F-RT-30/ 40/ 50 and I-FH-RT are osteotomy guides for Rohner Technique systems
Dental Implant Space/ Depth/ Angle Gauge	Parts labelled I-DG are depth gauges for general dental implants
	Parts labelled I-MD-SP are mucosal depth gauges for the Single Platform implants
	Parts labelled I-ADG and I-DI are angle gauges for general dental implants
	Parts labelled CH-I-DG and I-ZYG-DG-1 are depth gauges for zygomatic implants

	Parts labelled I-ZYG-DI55 and ZYG-TR-55-35/ 45/ 52.5/ 52 are angle gauges for zygomatic implants	
	Parts labelled LISA are space analyser gauges to measure the available space between teeth	
	Part labelled MAX-7/8/9/10-PG-xx and MAX-ROD are profile gauges for MAX implants	
Reusable Bone Taps (handheld)	Parts labelled IV-TR45xx/ TR50xx/ TR60xx are for use with the INVERTA implants (devices intended for up to 10 uses)	
Surgical Screwdrivers	Handheld Hexed	Parts labelled I-CS-HD/L, I-HD-S/M/L, I-HD-22U-S/M/L, I-HD-27-S/M/L
	Wrench Insert Hexed	Parts labelled I-WI-09, I-WI-22S/M/L, I-WI-22U-S/M/L, I-WI-27S/M/L
	Handheld Unigrip	Parts labelled I-UGI-S/M/L, I-UD-LH
	Wrench Insert Unigrip	Parts labelled I-WI-UG-S/M/L
	Handheld Quad	Parts labelled I-QDI-S/M/L, I-QD-LH
	Wrench Insert Quad	Parts labelled I-WI-QS/M/L
	Handheld Slotted/Blade	Parts labelled I-BD-S/M/L, I-BB-LH
	Wrench Insert Slotted/Blade	Parts labelled I-WI-BS/M/L
	Handheld Torx	Parts labelled I-SCS-S/M/L, I-HD-ASC-M/L
	Wrench Insert Torx	Parts labelled I-WI-SCS-S/M/L, I-WI-ASC-M/L
	Handheld Abutment Driver	Parts labelled I-AD, I-OAB-H
	Wrench Insert Abutment Driver	Parts labelled I-WI-A, I-WI-SA
Dental Implant Holding Endpieces	Handheld hexed connection	Parts labelled I-IMP-INS-2 and I-ZYG-INS-2 (uses replacement part I-X-RING) used with implants pre-mounted with fixture mounts
	Handheld square connection	Parts labelled I-IMP-INS-1 and I-ZYG-INS-1 used with implants pre-mounted with fixture mounts
	Wrench bit	Parts labelled I-WIID-S/ L (uses replacement part I-PBIT-H) used with External Hex Internal Drive implants
		Parts labelled I-WI-DC3-S/ M, I-WI-DC4-S/ M and I-WI-DC5-S/ M used with Deep Conical implants
		Parts labelled I-WI-LH35S/ LH35M and I-WIL-43-12D (uses replacement part I-PBIT-L18) used with TRI-NEX implants
		Parts labelled I-WI-LH43S/ LH43M, I-WI-LH50S/ LH50M, I-WI-LH60 and I-WIL50-12D (uses replacement part I-PBIT-L2) used with TRI-NEX implants
		Parts labelled I-WI-3M-S/ M/ L (uses replacement part I-PBIT-H16) used with PROVATA implants
		Parts labelled I-WI-M-S/ M/ L and I-WI-Z-S/ M/ L (uses replacement part I-PBIT-L18) used with Internal Hex M-Series and PROVATA implants
		Parts labelled I-WI-ITS and I-WI-ITS-S (uses replacement part I-PBIT-2) used with Internal Octagon IT implants
		Parts labelled I-WI-SP-S/ M/ L and SP-FMU /-L used with Single Platform implants
		Parts labelled I-WI-ILZ-S used with ILZ Mini implants
	Guided surgery fixture mounts	Parts labelled BA-GS-FMU, BAR12D-GS-FMU, IB-GS-FMU, IBN-GS-FMU, IBNT12D-GS-FMU, IBR12D-GS-FMU, PRO12D-GS-FMU are guided surgery fixture mounts for general dental implants
	Fixture mount extensions	Parts labelled I-FM-E, I-FME-XL, I-FME-XS, I-FME-XM, I-FM-H are fixture mounts extensions for general dental implants
	Fixture mount converters	Parts labelled I-ITS-FMU are for Internal Octagon IT implant fixture mounts
Dental Osteotomes	Parts labelled I-OST-(Ø)C/ P/ SU, I-FSD and I-HAMMER are for general dental implants	
Bone Mills	Handheld Bone Mills	Parts labelled I-HBM-DC30/35/40/50 are for Deep Conical implants
		Parts labelled I-HBML-35/43/50/60 are for Tri-Nex implants
		Parts labelled I-HBM-M-46/ M-56/ Z66 are for Internal Hex M-Series and PROVATA implants
	Bone Mill Caps	Parts labelled I-BNM/BPM/BM/BAM/BBBM-CAP are for External Hex implants
		Parts labelled I-BMZYG-CAP are for Zygomatic (ZYGIN/ZYGON) implants

Basic Dental Retractors	Parts labelled I-NIC/ -L are for general dental implants	
	Parts labelled I-ZYG-RET-1 are for zygomatic implants	
	Parts labelled I-PT-RET are for pterygoid implants	
Screw Extractor Guides	Parts labelled I-SRG-DC3/4/5 are for Deep Conical implants	
	Parts labelled I-SRG-EXT-HEX/ IBN or I-SRG-MSC-IP and I-SR4-GH are for External Hex implants	
	Parts labelled I-SRG-IT are for Internal Octagon IT implants	
	Parts labelled I-SRG-L-35/43/50/60 are for Tri-Nex implants	
	Parts labelled I-SRG-3M/M/Z are for Internal Hex M-Series and PROVATA implants	
	Parts labelled I-SRG-SP are for Single Platform implants	
Surgical Drill Guides	Drill Guide Spoons	Parts labelled DGS-L-(Ø) and DGS-S-(Ø) are drill guide spoons for general dental implants
	Drill Stabilizers	Parts labelled GSDS-(Ø) and I-TDS-(Ø) are drill stabilisers for general dental implants
	Rohner Technique Drill Guide	Parts labelled TDF-RT-(Ø) and ICF-RT are surgical guides for Rohner Technique systems
	Drill Stops	Parts labelled I-GS-TDS-20 are drill stops for general dental implants
	Offset Sleeves	Parts labelled I-GSSO4/5-(L) are offset sleeves for general dental implants
Dental Torque Wrenches	Parts labelled I-RATCHET-2, I-TW-45/ TWS-B45/ TWS-B100/ TWA-45 are torque wrenches for general dental implants	
Fixture Mount Spanners	Parts labelled I-SP-X and I-SP-IT are fixture mount spanners for general dental implants	
Wrench Adaptors	Parts labelled I-WI-C-L/ S, I-WI-CNB/ CST and I-WI-SH/ SL/ SS are wrench adaptors for general dental implants	
Retention Components	Parts labelled I-X-RING are for use with the I-CON-X	
	Parts labelled I-PBIT-L18/L2/2/H16/H are for implant insertion tools and dental implant holding endpieces	
Guide Stabilisers	Parts labelled I-GSN-SP and I-GS-SP for Single Platform (SP1) implants	
Space Maintainers	Parts labelled I-SMC6 are for use with Ø4.8 compact conical abutments	

Clinical performance

The Reusable instruments are intended to aid in the surgical procedures associated with implant therapy. Consequently, the clinical performance of the Reusable instruments is primarily defined by the success of the implant therapy. This performance can be quantitatively assessed through implant survival/success, which serves as an important indicator of the instruments' effectiveness and success in clinical practice.

Clinical benefits

This device aids in the successful outcome of surgical procedures and has no direct clinical benefit.

Storage, cleaning and sterilisation

For components supplied sterile (sterilised by gamma irradiation) and intended for multi-use prior to the expiration date (see packaging label), sterility is assured unless the container or seal is damaged or opened.

For all components, if the packaging is damaged, do not use the product and contact your Southern representative or return to Southern Implants®. The devices must be stored in a dry place at room temperature and not exposed to direct sunlight. Incorrect storage may influence device characteristics.

Limitations on reusable items

A direct value for reusable instruments cannot be given. Frequent processing may have minor effects on the instruments. The product life is normally determined by wear and damage during use, thus instruments if properly cared for and inspected after each use, can be reused many times. Maintain a checklist for these instruments recording the number of uses.

Prior to reprocessing the device, it should be thoroughly inspected and tested to determine its suitability for reuse.

The Bone Taps are the only instruments that have a life span of 10 uses per instrument.

NOTE: during use, handle drills and instruments with sterile tweezers to minimise contamination of the instrument tray and risk of damage to sterile surgical gloves.

Containment

As soon as practically possible, remove all visible residue after use (bone, blood or tissue), by immersing the instrument in cold water (dried soil may be difficult to remove).

Precleaning

Disassemble instruments from hand pieces and all connecting parts from instruments in order to clean soil from obstructed areas. Remove the PEEK bits from placement tools. Rinse with lukewarm water for 3 minutes and remove hardened debris with a soft nylon brush. Avoid mechanical damage during cleaning.

Manual cleaning or automated cleaning

Prepare an ultrasonic bath with suitable detergent (i.e. Steritech instrument cleaner - 5% dilution), sonicate for 20 minutes (alternative methods can be used if proven by the end user). Rinse with purified/sterile water.

NOTE: always follow the instructions for use of the manufacturers for cleaning agents and disinfectants.

Load devices into a thermodisinfectant. Run the cleaning and disinfection cycle, followed by the drying cycle.

Drying

Dry the instruments both inside and outside with filtered compressed air or single use, lint free wipes. Pack the instruments as quickly as possible after removal into the storage container. If additional drying is necessary, dry in a clean location.

Inspection

Do a visual inspection of the items to check for damages on instruments.

Packaging

Use the correct packaging material as indicated for steam sterilisation to ensure sterility is maintained. Double packaging is recommended. Where appropriate, the cleaned, disinfected and checked devices can be assembled and placed in instrument trays as applicable. Instrument trays can be double wrapped or placed in sterilisation bags.

Sterilisation

Southern Implants® recommends one of the following procedures to sterilise the device prior to use:

1. prevacuum sterilisation method: steam sterilise the device at 132°C (270°F) at 180 - 220 kPa for 4 minutes. Dry for at least 20 minutes in the chamber. Only an approved wrap or pouch for steam sterilisation must be used.
2. for users in the USA: prevacuum sterilisation method: wrapped, steam sterilise at 135°C (275°F) at 180 - 220 kPa for 3 minutes. Dry for 20 minutes in the chamber. Use a wrap or pouch that is cleared for the indicated steam sterilisation cycle.

NOTE: users in the USA must ensure that the steriliser, wrap or pouch, and all steriliser accessories are cleared by the FDA, for the intended sterilisation cycle.

Storage

Maintain packaging integrity to ensure sterility in storage. Packaging should be completely dry before storage to avoid corrosion and degradation of cutting edges.

Contraindications

The contraindications of the specific target implant and prosthetic devices apply. Therefore, the contraindications of the systems/medical devices utilized as part of implant surgery/therapy should be noted and the relevant documents consulted.

Contraindications specific to the Southern Implants® Reusable Instruments include:

- Allergies or hypersensitivity to chemical ingredients in the following materials: Titanium, Aluminium, Vanadium, Stainless Steel (DIN 1.4197, SAE 316, AISI 303, AISI 302 or DIN 1.4320), Silicon, Polyether Ether Ketone (PEEK), Acetal Copolymer (POM), Acetal homopolymer (POM-H) and Polyphenylsulfone (PPSU).

Warnings and precautions

IMPORTANT NOTICE: THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING

- To ensure the safe and effective use of dental implants, new technologies/systems and the Metal Abutment devices, it is strongly recommended that specialised training be undertaken. This training should include hands-on methods to gain competency on the proper technique, the system's biomechanical requirements and radiographic evaluations required for the specific system.
- Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone.
- Use of the device with incompatible or non-corresponding devices can result in the poor performance or failure of the device.
- When handling devices intraorally, it is imperative that they are adequately secured to prevent aspiration, as aspiration of products may lead to infection or physical injury.
- The use of non-sterile items can lead to secondary infections of the tissue or transfer infectious diseases.
- Failure to adhere to appropriate cleaning, re-sterilization, and storage procedures as outlined in the Instructions for Use can result in device damage, secondary infections, or patient harm.
- Exceeding the number of recommended uses for reusable devices can result in device damage, secondary infection or patient harm.
- The use of blunt drills may cause damage to the bone, potentially compromising osseointegration.

It is crucial to emphasize that training should be undertaken by both new and experienced implant users prior to utilizing a new system or attempting a new treatment method.

Patient Selection and Preoperative Planning

A comprehensive patient selection process and meticulous preoperative planning are essential for successful implant treatment. This process should involve consultation among a multidisciplinary team, including well-trained surgeons, restorative dentists, and laboratory technicians.

Patient screening should include, at a minimum, a thorough medical and dental history, as well as visual and radiological inspections to assess the presence of adequate bone dimensions, the positioning of anatomical landmarks, the presence of unfavourable occlusal conditions, and the periodontal health status of the patient.

For successful implant treatment, it is important to:

1. Minimize trauma to the host tissue, thereby increasing the potential for successful osseointegration.
2. Accurately identify measurements relative to radiographic data, as failure to do so may lead to complications.
3. Be vigilant in avoiding damage to vital anatomical structures, such as nerves, veins, and arteries. Injury to these structures may result in serious complications, including ocular injury, nerve damage, and excessive bleeding.

The responsibility for proper patient selection, adequate training and experience in implant placement, and the provision of appropriate information required for informed consent rests with the practitioner. By combining thorough screening of prospective implant candidates with a practitioner possessing a high level of competence in the use of the system, the potential for complications and severe side effects can be significantly reduced.

High Risk Patients

Special care should be exercised when treating patients with local or systemic risk factors that may adversely affect the healing of bone and soft tissue or otherwise increase the severity of side effects, the risk of complications, and/or the likelihood of implant failure. Such factors include:

- poor oral hygiene
- history of smoking/vaping/tobacco use
- history of periodontal disease
- history of orofacial radiotherapy**
- bruxism and unfavourable jaw relations
- use of chronic medications that may delay healing or increase the risk of complications including, but not limited to, chronic steroid therapy, ant-coagulant therapy, TNF- α blockers, bisphosphonate and cyclosporin

*** The potential for implant failure and other complications increases when implants are placed in irradiated bone, as radiotherapy can lead to progressive fibrosis of blood vessels and soft tissue (i.e., osteoradionecrosis), resulting in diminished healing capacity. Contributing factors to this increased risk include the timing of implant placement in relation to radiation therapy, the proximity of radiation exposure to the implant site, and the radiation dosage at that site.*

Should the device not operate as intended, it must be reported to the manufacturer of the device. The contact information or the manufacturer of this device to report a change in performance is: sicomplaints@southernimplants.com.

Side effects

The clinical outcome of treatment is influenced by various factors. The following side effects and residual risks are associated with the device group and may necessitate further treatment, revision surgery or additional visits to the relevant medical professional's office. Furthermore, these side effects and residual risks can occur with varying possible severities and frequencies.

- Pain, discomfort or tenderness
- Localized inflammation
- Soft tissue irritation
- Hematoma or bruising
- Minor bleeding
- Allergic or hypersensitivity reaction(s)
- Bruising
- Dental injury during surgery
- Fistula formation
- Gingival injury
- Infection (acute and/or chronic)
- Localized inflammation
- Loss or damage to adjacent dentition
- Pain or discomfort
- Peri-implantitis, peri-mucositis or otherwise poor peri-implant soft tissue health
- Soft tissue irritation
- Surgical side effects such as pain, inflammation, bruising and mild bleeding
- Wound dehiscence or poor healing

Precaution: maintaining sterility protocol

Devices supplied sterile are packed in a peel pouch or blister base with a "peel-back" lid. Labelling information is located on the bottom half of the pouch, inside the packet or on the surface of the peel-back lid. Sterility is assured unless the pouch is damaged or opened.

Devices supplied non-sterile are supplied clean but not sterile in a peel pouch, blister base with peel-back lid, zipper bag, or tin box with sponge insert. Labelling information is located on the bottom half of the pouch or on the surface of the peel-back lid.

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established.

The contact information for the manufacturer of this device to report a serious incident is as follows:

sicomplaints@southernimplants.com.

Materials

Device	Device Material
Osteotomy Guides (Cutting Guides & Fibula Holders)	Stainless Steel 431 (hardened to 45-55 HRC), Stainless Steel 304, Titanium Alloy Ti-6Al-4V (90% titanium, 6% aluminium, 4% vanadium) according to ASTM F136 and ISO 5832-3, and Medical Grade PEEK
Reusable Dental Implant Space/Depth/Angle Gauge	Commercially pure titanium (Grade 2, 3, 4) according to ASTM F67 and ISO 5832-2, Titanium Alloy Ti-6Al-4V (90% titanium, 6% aluminium, 4% vanadium) according to ASTM F136 and ISO 5832-3, Stainless Steel 316/316L/303
Bone Tap, Reusable	Commercially pure titanium (Grade 4) according to ASTM F67 and ISO 5832-2
Surgical screwdrivers, reusable	Stainless Steel (DIN 1.4197) (Harden to 48-52 HRC), Stainless Steel (SAE 316), Commercially pure titanium (Grade 4) according to ASTM F67 and ISO 5832-2, RADEL PPSU, Stainless Steel (SAE 316)
Dental Implant holding endpieces	Stainless Steel (DIN 1.4197), Stainless Steel 316, Commercially pure titanium (Grade 4) according to ASTM F67 and ISO 5832-2
Dental Osteotomes	Stainless Steel 316, Celeron/Tufnol
Bone Mills	Commercially pure titanium (Grade 4) according to ASTM F67 and ISO 5832-2, Titanium Alloy Ti-6Al-4V (90% titanium, 6% aluminium, 4% vanadium) according to ASTM F136 and ISO 5832-3 Handle - Stainless Steel (SAE 316) Idler - Commercially pure titanium (Grade 4) according to ASTM F67 and ISO 5832-2 Rivet - Stainless Steel (SAE 316)
Basic Dental Retractors	Aluminium, Stainless Steel
Screw extractor guides, reusable	Commercially pure titanium (Grade 4) according to ASTM F67 and ISO 5832-2, Titanium Alloy Ti-6Al-4V (90% titanium, 6% aluminium, 4% vanadium) according to ASTM F136 and ISO 5832-3
Surgical drill guides, reusable	Stainless Steel (AISI 303), Stainless Steel (SAE 316), Commercially pure titanium (Grade 4) according to ASTM F67 and ISO 5832-2, Titanium Alloy Ti-6Al-4V (90% titanium, 6% aluminium, 4% vanadium) according to ASTM F136 and ISO 5832-3
Dental Torque Wrenches	Stainless Steel (SAE 316)
Fixture Mount Spanners	Stainless Steel (SAE 316)
Wrench Adaptors	Stainless Steel DIN 1.4197, Titanium Alloy Ti-6Al-4V (90% titanium, 6% aluminium, 4% vanadium) according to ASTM F136 and ISO 5832-3, Silicon
Retention Components	Stainless Steel Spring Wire (AISI 302 or DIN 1.4320), Medical Grade PEEK, Silicone
Guide Stabilisers	Titanium Alloy Ti-6Al-4V (90% titanium, 6% aluminium, 4% vanadium) according to ASTM F136 and ISO 5832-3
Space Maintainers	Titanium Alloy Ti-6Al-4V (90% titanium, 6% aluminium, 4% vanadium) according to ASTM F136 and ISO 5832-3

Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

Summary of Safety and Clinical Performance (SSCP)

As required by the European Medical Device Regulation (MDR; EU2017/745), a Summary of Safety and Clinical Performance (SSCP) is available for perusal with regard to Southern Implants® product ranges.

The relevant SSCP can be accessed at <https://ec.europa.eu/tools/eudamed>.

NOTE: the above website will be available upon the launch of the European Database on Medical Devices (EUDAMED).

Disclaimer of liability

This product is part of the Southern Implants® product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants® product range and take full responsibility for the correct indications and use of this product. Southern Implants® does not assume liability for damage due to incorrect use. Please note that some Southern Implants® products may not be cleared or released for sale in all markets.

Basic UDI

Product	Basic-UDI Number
Basic-UDI for Reusable Instruments	6009544038769E

Related literature and catalogues

CAT-2004 - Tri-Nex® Implants Product Catalogue
 CAT-2005 - IT Implants Product Catalogue
 CAT-2010 - Osseointegrated Fixtures Catalogue
 CAT-2020 - External Hex Implants Product Catalogue
 CAT-2042 - Deep Conical Implants Product Catalogue
 CAT-2043 - Internal Hex Implants Product Catalogue
 CAT-2060 - PROVATA® Implants Product Catalogue
 CAT-2069 - Deep Conical INVERTA® Implants Product Catalogue
 CAT-2070 - Zygomatic Implants Product Catalogue
 CAT-2087M - ILZ Mini Implants Product Catalogue
 CAT-2088M - Nazalus Implants Product Catalogue and Surgical Manual
 CAT-2092 - Deep Conical Pterygoid Implants Product Catalogue
 CAT-2093 - Single Platform (SP1) Implants Product Catalogue
 CAT-2095 - External Hex INVERTA® Implants Product Catalogue
 CAT-2096 - External Hex Pterygoid Implants Product Catalogue

Symbols and warnings

											
Manufacturer: Southern Implants® 1 Albert Rd, P.O Box 605 IRENE, 0062, South Africa. Tel: +27 12 667 1046	CE mark 2797	Prescription device*	Sterilised using irradiation	Non-sterile	Use by date (mm-yy)	Do not reuse	Do not resterilise	Catalogue number	Batch code	Medical device	Authorised representative in the European Community
											
Authorised representative for Switzerland	Date of manufacture	Magnetic Resonance conditional	Magnetic Resonance safe	Single sterile barrier system with protective packaging inside	Single sterile barrier system	Consult instruction for use	Caution	Keep away from sunlight	Do not use if package is damaged		

* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

All rights reserved. Southern Implants®, the Southern Implants® logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants®. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale. It is the responsibility of the clinician to inspect the symbols that appear on the packaging of the product in use.

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