
English	INSTRUCTIONS FOR USE: Southern Implants® Osseointegrated Fixtures
Español	INSTRUCCIONES DE USO: Implantes Osteointegrados Southern Implants®
Italiano	ISTRUZIONI PER L'USO: Southern Implants® Fixture Osteointegrate
Français	MODE D'EMPLOI : Implants ostéo-intégrés de Southern Implants®
Deutsch	GEBRAUCHSANWEISUNGEN: Southern Implants® Osseointegrierte Implantate
Português	INSTRUÇÕES DE UTILIZAÇÃO: Fixações osseointegradas Southern Implants®
Ελληνικά	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ: Οστεοενσωματούμενα εξαρτήματα της Southern Implants®
Svenska	ANVISNINGAR FÖR ANVÄNDNING: Southern Implants® Osseointegrerad fixturer

South Africa - Headquarters: 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

Subsidiaries

Australia

Southern Implants Australia
T: +61-(0)-8-9466-2627
E: info@southernimplants.com.au

Spain and Portugal

Southern Implants Iberica
T: +34 935 053 507
E: info@southernimplants.es

United Kingdom and Ireland

Southern Implants UK
T: +44-20-8899-6845 / 6 / 7
E: info@southernimplants.co.uk

USA and Canada

Southern Implants North America Inc.
T: +1-561-472-0990
E: customercare@southernimplants.com

Description

The Osseointegrated Fixture is a self-tapping tapered implant made of commercially pure Grade 4 Titanium (UTS ≥ 900 MPa). All implants are surface-roughened by grit-blasting with aluminium oxide. The surface has an average Sa value of 1.4 microns.

Intended use

The Osseointegrated Fixtures are intended for surgical placement in the skull, mainly the mastoid process, periorbital and perinasal regions, as a means to provide rigid retention for cranio-maxillofacial epistheses.

These devices are intended for single use on a single patient.

Indications for use

The Osseointegrated Fixtures are indicated for use in the restorative treatment of patients with congenital or acquired cranio-maxillofacial defects to enable the attachment of external aesthetic restorations or epistheses in cases where other means of attachment are considered to be inadequate.

The devices are indicated for two-stage procedures only.

Intended user

The intended user for the Osseointegrated fixtures includes maxillofacial surgeons, plastic surgeons, anaplastologists (as part of the surgical planning and epistheses design) and other appropriately trained medical professionals.

Intended environment

This system is intended to be used in a clinical environment such as an operating theatre.

Intended patient population

The intended patient population for the Osseointegrated Fixtures includes those with acquired or congenital defects of the craniofacial or maxillofacial region that are eligible, or otherwise not contraindicated, for an extra-oral implant-retained episthesis.

Compatibility information

Southern Implants' implants should be restored with Southern Implants' components or compatible Steco Titanmagnetics Components. Use components that correspond to the connection type and prosthetic platform when restoring the Osseointegrated Fixtures. For further information see the Osseointegrated Fixtures Product Catalogue (CAT-2010).

Table A – Osseointegrated Fixtures implant range

Item Code	Implant Length Code	Implant Length (mm)	Implant Diameter (mm)
IE*	IE3	3	Ø3.75
	IE4	4	Ø3.75
	IE6	6	Ø3.75

**Osseointegrated Fixtures are packed with their corresponding cover screw*

Clinical benefits

The Osseointegrated Fixtures are used to provide a rigid retention method for a removable aesthetic restoration as part of the restorative treatment of maxillofacial and/or craniofacial defects. As such, the clinical benefits of these devices are closely related to aesthetics and the associated rigid retention of the episthesis (i.e., the prosthetic device used to restore facial defects). These clinical benefits include increased stability and resistance to episthesis movement, the provision of an alternative treatment in cases where traditional retention methods are insufficient or unsuitable, improved aesthetics and overall comfort to the patient, a decrease in the skin irritation associated with skin adhesives, and increased patient satisfaction and psychological well-being.

Storage, cleaning and sterilisation

The component is supplied sterile (sterilised by gamma irradiation) and intended for single-use prior to the expiration date (see packaging label). Sterility is assured unless the container or seal is damaged or opened. If packaging is damaged do not use the product and contact your Southern representative or return to Southern Implants®. The devices must be stored in a dry place at room temperature and not exposed to direct sunlight. Incorrect storage may influence device characteristics. Do not reuse components indicated for single-use only. Reusing these components may result in:

- damage to the surface or critical dimensions, which may result in performance and compatibility degradation.
- adds the risk of cross-patient infection and contamination if single-use items are reused.

Southern Implants® does not accept any responsibility for complications associated with reused single-use components.

Contraindications

Contraindications are circumstances, conditions, symptoms or factors that increase the risks associated with a medical procedure, drug or treatment where these can be considered to be either absolute or relative in nature. An absolute contraindication is one that results in the progression of treatment to be inadvisable without exception. However, in the case of relative contraindications, the user may proceed with the treatment if the necessary precautions have been taken or it is determined that the benefit to the patient would outweigh the increased risk. These cases represent high risk patients.

Absolute contraindications for Osseointegrated Fixtures include:

- cases where full functional support for the episthesis is not attainable (i.e., the number of implants that could be placed would provide insufficient support even if combined with other retention methods)
- a lack of oncological radicality (i.e., persistent cancer presence or high risk of cancer recurrence due to potentially incomplete eradication)
- conditions that prevent proper implant hygiene including, but not limited to, severe mental disorders (e.g., dementia), cachexia, drug or alcohol addiction or disability
- cases where there is an inability to adhere to a follow-up schedule
- patients who are medically unfit for implant surgery
- patients with pure titanium, titanium alloy (Ti6Al4V), gold, palladium, polyether ether ketone (PEEK) or iridium allergies or hypersensitivities
- presence of infection or pathology at the implant site

The relative contraindications for Osseointegrated Fixtures include:

- mild psychiatric disorders
- poor craniofacial bone condition (inadequate bone quality and/or quantity)
- compromised vasculature on the surgical bed
- presence of extensive soft tissue scarring
- history of orofacial radiation therapy following ablative oncologic surgery
- minors (parent/guardian to ensure adequate cleaning)
- lack of sufficient tissue covering
- smoking/vaping/tobacco use
- patients on chronic medications that may delay healing or increase the risk of complications including, but not limited to, chronic steroid therapy, ant-coagulant therapy, TNF- α blockers, Bisphosphonate and Cyclosporin
- decompensated/uncontrolled diabetes
- pathological skin conditions including psoriasis

Warnings and precautions

THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING.

- For the safe and effective use of the Osseointegrated Fixtures it is strongly suggested that specialised training be undertaken. This training should include hands-on methods to gain competency on the proper technique, biomechanical requirements and radiographic evaluations required for the specific system.
- Improper technique can result in an increased potential for implant failure, damage to nerves/vessels, loss of supporting bone or other more severe complications and side effects.
- Use of the device with incompatible or non-corresponding devices can result in device failure or poor performance.
- The use of non-sterile items can lead to secondary infections of the tissue or transfer infectious diseases.
- Failure to perform appropriate cleaning, re-sterilisation and storage procedures as per the Instructions for Use document can result in device damage or secondary infection.
- Blunt drills may cause damage to the bone which could compromise osseointegration.

It should be noted that training should be undertaken by *both* new and experienced implant users before using a new system or attempting a new treatment method.

It is important that thorough screening of prospective implant candidates be performed. This screening should include, at minimum, a comprehensive medical history taking in which special consideration is given to the aetiology of the defect and a visual/radiological inspection to determine the condition of the bone, bone dimensions and anatomical landmarks. Failure to accurately estimate lengths of drills relative to radiographic measurements can result in permanent injury to nerves or other vital structures.

Patient co-operation is a critical factor in the success of an extraoral implant and its importance should be thoroughly communicated to the patient during the initial screening process and all subsequent appointments. This co-operation includes the adherence to routine follow-ups and home care procedures throughout its lifecycle in order to maintain excellent hygiene around the implant site.

Once the patient has been screened and considered suitable for an implant-retained prosthesis, it is vital to ensure proper preoperative planning is performed. This preoperative planning should include craniofacial MRI and CT scan analysis, 3D reconstructions of the craniofacial bones, and clinical examination of the deformity. Pre-establishment of implant sites is crucial, during which consideration is given to bone parameters (thickness and shape), implant stability, minimization of the impact of unwanted forces on the implants, planning locations for replacement implants in case of implant failure at the originally selected sites, and maintaining the angles at which the implants are inserted so that other prosthetic components can subsequently be fixed to them.

Furthermore, this preoperative planning should include consultation between a multi-disciplinary team, including a prosthetist or anaplastologist. These factors are essential to ensure the most optimal location is selected for implant placement to maximise implant stability in the bone and promote a well-fitting prosthesis and are important considerations in the promotion of successful treatment outcomes.

For successful implant treatment, it is important to (1) minimise the trauma to the host tissue as this increases the potential for successful osseointegration, (2) identify actual measurements relative to the radiographic data as failing to do so could lead to complications, and (3) be aware and avoid damage to vital anatomical structures such as nerves, veins and arteries. Injury to these structures may cause serious complications, including injury to the eye, nerve damage and excessive bleeding.

It should be noted that the prognosis of the implant is dependent on implant site where the temporal region generally presents with better prognosis of survival and fewer complications and the orbital areas generally presents with a greater risk.

Special care should be taken when treating patients with local or systemic risk factors that could affect the healing of the bone and soft tissue or otherwise increase the severity of the side effects, risk of complications and/or implant failure.

These factors include:

- mild psychiatric disorders
- poor condition of the craniofacial bone (inadequate bone quality and/or quantity)
- compromised vascularity on the surgical bed
- presence of extensive soft tissue scarring
- history of orofacial radiation therapy following ablative oncologic surgery¹
- minors²
- lack of sufficient tissue covering
- smoking
- patients on chronic medications that may delay healing or increase the risk of complications including, but not limited to, chronic steroid therapy, anti-coagulant therapy, TNF- α blockers, Bisphosphates and Cyclosporin
- decompensated/uncontrolled diabetes
- pathological skin conditions including psoriasis

The age of the patient should also be considered as age-associated natural decline in vascularity, healing capacity and bone density/quality can affect the ability of the bone to fuse with the implant. However, factors such as overall bone health and the specific condition of the patient's bones will play a significant role.

Bacterial and/or recurrent meningitis is an increased risk in patients with cranio-facial anatomical abnormalities, therefore, this increased risk should be considered during the planning phase of any surgery or treatment modality provided to these patients where extra precautionary measures should be taken.

With the combination of the thorough screening of prospective implant candidates, a user with a high level of competence in the use of the system and strict adherence to aftercare instructions and hygiene maintenance measures, the potential for complications and severe side effects can be reduced.

The responsibility for proper patient selection, adequate training and experience in implant placement and the provision of the appropriate information required for informed consent rests with the practitioner.

¹ *The potential for implant failure increases when implants are placed in irradiated bone. This is due to the diminished healing capacity of the region which is secondary to the progressive fibrosis of vessels and soft tissue caused by radiotherapy. Additionally, the use of implants in irradiated bone may result in:*

- *delayed or failed osseointegration of implants due to reduced bone vascularity, clinically expressed as osteoradionecrosis.*
- *tissue dehiscence and osteoradionecrosis.*
- *implant failure and loss.*

The success of implant placement in irradiated patients is dependent on various factors, including the timing of implant placement in relation to the radiation therapy, the anatomic site of the implant placement as well as the radiation dosage at that site and consequent risk of osteoradionecrosis. It is recommended that these patients undergo hyperbaric oxygen treatment to ameliorate the reduction in the osteogenic potential and micro-vascularisation caused by irradiation.

² *Minors represent a relative contraindication which can be controlled by ensuring the parents/guardians are adequately educated in the hygiene required for implant maintenance.*

Should the device not operate as intended, it must be reported to the manufacturer of the device. The contact information for the manufacturer of this device to report a change in performance is: sicomplaints@southernimplants.com.

Side effects

The side effects associated with the implantation of the Osseointegrated Fixtures are comparable with the side effects associated with other forms of craniofacial surgery.

These include:

- Pain or tenderness
- Localized inflammation
- Soft tissue irritation
- Discomfort
- Hematoma or bruising
- Bleeding

Additionally, the normal side effects associated with anaesthesia should also be expected.

In addition to the side effects expected as a result of the implantation of the Osseointegrated Fixtures, there are also additional potential or foreseeable complications. These complications represent the relative risks associated with proceeding with the treatment of these patients where further intervention, medical treatment or removal of the implant may be required.

These include:

- Infection of the surrounding soft tissue
- Implant contact with underlying sigmoid sinus or dura when inserted in the mastoid region
- Anatomical landmark damage
- Nerve damage resulting in transient or permanent weakness, numbness, or pain
- Perforation of the dura
- Soft tissue complications including erythema, moisture secretion and granulation
- Skin overgrowth around/over abutment

Precaution: maintaining sterility protocol

Implants are packaged as follows:

1. An outer package consisting of a rigid, clear box which acts as protection for the inner package.
2. The inner package consisting of a blister pack (clear plastic-formed blister base with a TYVEK “peel-back” lid).
3. Within the inner package, there is a hollow tube which contains one implant suspended from a titanium ring, this ensures the implant never touches the inside of the plastic tube.
4. Labelling information is located on the surface of the peel-back lid and on the outside of the rigid box.

Care must be taken to maintain the sterility of the implant by proper opening of the packaging and handling of the implant.

1. Open the implant package in the non-sterile field, with non-sterile gloves, tear the address label to open the box.
2. With non-sterile gloves, remove the inner blister pack. Do not place the plastic box or blister pack-lid onto the sterile field. The contents of this inner package are sterile.
3. The sealed blister is to be opened by an assistant (with nonsterile gloves): remove the TYVEK lid and drop or place the sterile tube onto the sterile field, open the tube cap and attach the implant placement tool onto the implant and carefully remove from the sterile tube. Do not touch the sterile implant.

Other sterile components are packed in a peel pouch or blister base with a “peel-back” lid. Labelling information is located on the bottom half of the pouch, inside the packet or on the surface of the peel-back lid. Sterility is assured

unless the pouch is damaged or opened. Non-sterile components are supplied clean but not sterile in a peel pouch or blister base with peel-back lid. Labelling information is located on the bottom half of the pouch or on the surface of the peel-back lid.

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established.

The contact information for the manufacturer of this device to report a serious incident is as follows: sicomplaints@southernimplants.com.

Materials

Material type Commercially pure titanium (grade 4, ASTM F67 and ISO5832-2, UTS ≥ 900 MPa)

Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

MR safety

Nonclinical testing has demonstrated that the Southern Implants® osseointegrated fixtures, metallic abutments and prosthetic screws are MR conditional.

A patient with these devices can be safely scanned in a MR system meeting the following conditions:

- static magnetic field of 1.5 Tesla and 3.0 Tesla only.
- maximum spatial gradient magnetic field of 3000 Gauss/cm (30 T/m).
- maximum MR system reported SAR corresponding to Normal Operating mode for all landmarks (Head SAR of 3.2 W/kg for head landmark, 2 W/kg whole body, and appropriate partial body SAR for other landmarks). For imaging landmarks above the thorax, a continuous scan time of 15 minutes will require a cooling delay of at least 5 minutes.
- in the non-clinical testing, the image artifact caused by the device extends approximately 20 mm from the Southern Implants' osseointegrated fixtures, abutments and prosthetic screws, when imaged with a gradient echo pulse sequence and a 3.0 Tesla MRI system.

Removable restorations should be taken out prior to scanning, as is done for watches, jewellery etc.

Should there be no MR symbol on the product label, please note that this device has not been evaluated for safety and compatibility in the MR environment.

Summary of Safety and Clinical Performance (SSCP)

As required by the European Medical Device Regulation (MDR; EU2017/745), a Summary of Safety and Clinical Performance (SSCP) is available for perusal with regard to Southern Implants® product ranges.

The relevant SSCP can be accessed at <https://ec.europa.eu/tools/eudamed>.

NOTE: the above website will be available upon the launch of the European Database on Medical Devices (EUDAMED).

Disclaimer of liability

This product is part of the Southern Implants® product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants® product range and take full responsibility for the correct

indications and use of this product. Southern Implants® does not assume liability for damage due to incorrect use. Please note that some Southern Implants® products may not be cleared or released for sale in all markets.

Basic UDI

Product	Basic-UDI Number
Basic-UDI for Osseointegrated Fixtures	60095440422484

Related literature and catalogues

CAT-2010- Osseointegrated Fixtures Product Catalogue

Symbols and warnings

											
Manufacturer: Southern Implants® 1 Albert Rd, P.O Box 605 IRENE, 0062, South Africa. Tel: +27 12 667 1046	CE mark 2797	Prescription device*	Sterilised using irradiation	Non-sterile	Use by date (mm-yy)	Do not reuse	Do not resterilise	Catalogue number	Batch code	Medical device	Authorised representative in the European Community
											
Authorised representative for Switzerland	Date of manufacture	Magnetic Resonance conditional	Magnetic Resonance safe	Single sterile barrier system with protective packaging inside	Single sterile barrier system	Consult instruction for use	Caution	Keep away from sunlight	Do not use if package is damaged		

* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

All rights reserved. Southern Implants®, the Southern Implants® logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants®. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale. It is the responsibility of the clinician to inspect the symbols that appear on the packaging of the product in use.

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