
English	INSTRUCTIONS FOR USE: Southern Implants® Laboratory Components
Español	INSTRUCCIONES DE USO: Componentes de Laboratorio Southern Implants®
Italiano	ISTRUZIONI PER L'USO: Southern Implants® Componenti di Laboratorio
Français	MODE D'EMPLOI : Composants de laboratoire de Southern Implants®
Deutsch	GEBRAUCHSANWEISUNGEN: Southern Implants® Laborkomponenten
Português	INSTRUÇÕES DE UTILIZAÇÃO: Componentes de laboratório Southern Implants®
Ελληνικά	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ: Εργαστηριακά εξαρτήματα της Southern Implants®
Svenska	ANVISNINGAR FÖR ANVÄNDNING: Southern Implants® Laboratoriekomponenter
Română	INSTRUCȚIUNI DE UTILIZARE: Southern Implants® Componente laborator

South Africa - Headquarters: 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

Subsidiaries

Australia

Southern Implants Australia
T: +61-(0)-8-9466-2627
E: info@southernimplants.com.au

Spain and Portugal

Southern Implants Iberica
T: +34 935 053 507
E: info@southernimplants.es

United Kingdom and Ireland

Southern Implants UK
T: +44-20-8899-6845 / 6 / 7
E: info@southernimplants.co.uk

USA and Canada

Southern Implants North America Inc.
T: +1-561-472-0990
E: customercare@southernimplants.com

Description

Brass Screws

The brass screws are used in the dental laboratory as a cost-effective alternative to the final prosthetic screw. The brass screw is intended to allow the laboratory technician to replicate the seating of the restoration onto a laboratory analogue and ensure the fit of the final prosthetic screw through the screw access hole of the restoration.

Cylinders

The plastic cylinders are supplied with Gold abutments and Passive abutments and intended for use in the laboratory to maintain the screw access hole and, screw seat (passive abutments only) during waxing up procedures of a castable restoration.

Plastic Sleeves

The plastic sleeves are intended for use in the laboratory as an aid in waxing up a restoration, during conventional casting procedures.

Gold Bars

The gold bar is connected to two implants placed in the anterior to provide stability for dentures. Retention clips are placed into the denture that will clip over the gold bar to aid in stabilizing the denture.

Laboratory Analogues

The laboratory analogues are used to replicate the location, orientation and connection of the prosthetic interface of an implanted dental device, namely the implant or abutment, in a stone or printed model during fabrication processes.

Digital Laboratory Analogues

The digital laboratory analogues are used as dental implant replicas to be inserted into a 3D printed or milled model in order to duplicate the location, orientation and restorative platform of the implant placed in the mouth during CAD/CAM procedures.

Luting Screws

The luting screws are used to clamp the interfacial component of a passive abutment onto the laboratory analogue during the process of luting the casting onto the interfacial component. The luting screw maintains the screw access hole and creates a chimney in the resulting restoration.

Pin Gauges

The pin gauges are devices used to keep the screw channel of a restoration or crown open when packing porcelain around the pin, or when waxing up a restoration around the pin.

Polishing Caps

The polishing caps are used to hold a restoration after cementing has taken place, allowing the dental professional to polish off any excess cement.

Post Retention Components

The retention posts are used during the fabrication of removable dentures and are soldered into a bar for retention of a removable overdenture.

Shoulder Analogues

The shoulder analogues are used to replicate the location, orientation of a shouldered abutment and aid in impression taking.

Titanium Bars

The titanium bars are used to fabricate dentures by providing stability to dentures as they connect multiple implants together and provide retention points where retention posts are soldered on.

Analogue Screws

The analogue screws are used to secure the digital laboratory analogues in place within a 3D printed or milled model.

Analogue Try-in Abutments

The try-in abutments are used with implants and laboratory analogues to confirm the required orientation of abutments, prior to selecting the desired abutment for the restoration.

Reamers and Lapping Tools

The reamers are used to ensure the passive abutment screw access hole diameter is correct and that the screw seat is flat. The lapping tools are used with the castable UCLA hexed sleeve and UCLA non-hexed sleeve to improve the fitting surface of the cast abutment.

O-Rings

The O-rings are used as a processing aid during laboratory procedures to ensure that acrylic does not enter the retentive cap of the denture during curing by sealing the interface. The O-ring is then removed and discarded.

Plastic clips and spacers

The plastic spacers are used with cementable abutments to allow for a screw retained restoration to be made. The plastic clips are laboratory clips which interface with a ball abutment.

Anaplastology Components

The Anaplastology components are used in the laboratory as an aid in prosthetic restoration by anaplastologists.

Intended use

The Southern Implants' Laboratory Components are intended to aid in the treatment of partially or fully edentulous patients eligible for placement of one or more dental implants as a means of fixing a permanent or removable single crown, partial or full-arch dental prosthesis in the upper or lower jaw.

Indications for use

The Laboratory Components are indicated for use in cases of partial or full edentulousness, congenital tooth loss, or failing/problematic natural dentition, when a patient is dissatisfied with existing removable dentures, or as part of the revision/replacement of previous dental restorations.

Intended user

The intended user for this system includes Dental Technicians, Maxillo-facial Surgeons, General Dentists, Orthodontists, Periodontists, Prosthodontists, and other appropriately trained and experienced implant users.

Intended environment

This device is intended to be used in a dental laboratory as part of the restoration design and manufacture as well as in a clinical environment such as an operating theatre or a dentist consultation room.

Intended patient population

The intended patient population for dental implant therapy includes partially or fully edentulous patients eligible or, otherwise, not contraindicated for implant placement that require prosthetic restoration or revision of existing restorations in the upper or lower jaw; where the planned restorations are to be fixed and comprise of single teeth, partial or full bridges.

Compatibility information

The Laboratory Components can be identified by specific abbreviations in the product codes. Device identifiers are summarised in Table A. Table B summarises the brass screws available for each Southern range (corresponding prosthetic screws indicated in square brackets []) and the recommended torque values of the Southern screws.

Table A – Laboratory Component device codes

Implant connection type	Compatible device
Plastic Cylinders	Parts labelled ITS6-PC1 and ITS-PC1 for Internal Octagon (IT) engaging cylinders.
	Parts labelled ITS6-PC1ne and ITS-PC1ne for Internal Octagon (IT) non-engaging cylinders.
	Parts labelled IT-PE-x and IT6-PE-x where xx denotes the post height for Internal Octagon (IT) engaging cylinders.
	Parts labelled IT-PE-NE and IT6-PE-NE where xx denotes the post height for Internal Octagon (IT) non-engaging cylinders.
	Parts labelled PC-EL-xx and PC-NL-xx where xx denotes the corresponding platform diameter for TRI-NEX engaging and non-engaging cylinders respectively.
	Parts labelled PC-M-E and PC-M-NE for Internal Hex engaging and non-engaging cylinders respectively.
	Parts labelled SB1, SB5, SBA1, SBA5, SBBB1, SBBB5, SBN1, SBN5, SSD4 for External Hex cylinders
Plastic Sleeves	Parts labelled ITS6-PA-P and ITS-PA-P for Internal Octagon (IT) Passive Abutment engaging plastic sleeves.
	Parts labelled ITS6-PA-NEP and ITS-PA-NEP for Internal Octagon (IT) Passive Abutment non-engaging plastic sleeves.
	Parts labelled PA-DCx-P and PA-NDCx-P where x denotes the corresponding platform size for the Deep Conical Passive Abutment engaging and non-engaging plastic sleeves respectively.
	Parts labelled PA-EL-xP and PA-NL-xP where x denotes the corresponding platform size for the TRI-NEX Passive Abutment engaging and non-engaging plastic sleeves respectively.
	Parts labelled PA-3EM-SP, PA-EM-SP, PA-EM-SP45-P and PA-EZ-P for the Internal Hex Passive Abutment engaging plastic sleeves.
	Parts labelled PA-3NM-SP, PA-NM-SP, PA-NM-WP and PA-NZ-P for the Internal Hex Passive Abutment non-engaging plastic sleeves.
	Parts labelled PA-SP-P and PA-SP-PM-P for the Single Platform (SP1) Passive Abutment engaging plastic sleeves.
	Parts labelled PA-NSP-P for the Single Platform (SP1) Passive Abutment non-engaging plastic sleeves.
	Parts labelled PA-MC-48P and PMC-1 for the Ø4.8 Compact Conical Abutment interface.
	Parts labelled PA-MC-60P and PMCW-1 for the Ø6.0 Compact Conical Abutment interface.
	Parts labelled SP16PH, SB16P, SBA16PH, SBBB16PH, SBN16PH, SBNV16PH, SBV16PH and SMAX9PH for the External Hex Passive Abutment engaging plastic sleeves.
	Parts labelled SP-17P, SB17P, SBA17P, SBBB17P, SBN17P, SBV-17-9 and SMAX9P for the External Hex Passive Abutment non-engaging plastic sleeves.
	Parts labelled SB-DBAN-1 and SB-DBN-1 for the DBAN and DBN abutments.
	Parts labelled SLZ for ILZ implants.
	Parts labelled SI30-SLEEVE for SI30 abutment.
	Parts labelled SMP-1-S for the SMP1 abutment.
	Parts labelled GC3-P, GC-30P, GCBA-P, GCBBB-P, GCBN-P, GCB-P and GCMAX-P for the External Hex Gold Abutment plastic sleeves.
	Parts labelled GC-DCx-P where x denotes the corresponding platform size for the for the Deep Conical Gold Abutment plastic sleeves.
	Parts labelled GC-M-P and GC-Z-P for the Internal Hex Gold Abutment plastic sleeves.
	Parts labelled GC-PL-45 for the TRI-NEX Gold Abutment plastic sleeves.
	Parts labelled ITS6-GC1-P and ITS-GC1-P for the Internal Octagon (IT) Gold Abutment plastic sleeves.
	Parts labelled SB-2-6-P and SBBB-2-6-P for the Compact Conical Gold Abutment plastic sleeves.
	Parts labelled CC-L-35PC and CC-L-43PC for TRI-NEX Chrome Cobalt cylinder plastic sleeves.
	Parts labelled CC-PS-40 for External Hex Chrome Cobalt cylinder plastic sleeves.
	Parts labelled ITS6-CC1-PC and ITS-CC1-PC for Internal Octagon (IT) Chrome Cobalt cylinder plastic sleeves.
	Parts labelled ITS6-TA and ITS-TA for TSA6/ TSAF plastic sleeves.
Gold Bars	Parts labelled GD1 and GDC1.
Laboratory Analogues	Parts labelled LA-DCx where x denotes the corresponding platform size for the Deep Conical lab analogues.
	Parts labelled LA-L-x where x denotes the corresponding platform size for the TRI-NEX lab analogues.
	Parts labelled LA-3M, LA-M, LA-M-P45 and LA-Z for the Internal Hex lab analogues.
	Parts labelled LA-ILZ and LSLZ1 for the ILZ lab analogues.
	Parts labelled LITS4 and LITS6 for the Internal Octagon (IT) lab analogues.
	Parts labelled LP12, LSN12, LS12, LSA12, LBBBS12 and LMAX9 for the External Hex lab analogues.
	Parts labelled LSMC1, LSMW1, LSA7 and LS7 for the Compact Conical Abutment lab analogues.
	Parts labelled LAFx and LAF6-x where x denotes the corresponding post height for the TSAF lab analogues.
	Parts labelled LOAB-1, LOAB-H, AO2, PP2 and BP2 for Ball/ Overdenture abutment lab analogues.

	Parts labelled LSDBAN4 and LSDBN-4 for the DBAN and DBN abutment lab analogues.
	Parts labelled LSD4 and LSI3 for the Calcitek integral implant lab analogues.
	Parts labelled LSI1 for the Shoulder abutment lab analogues.
	Parts labelled LS1 for the Standard abutment lab analogues.
	Parts labelled LS6BBB1 for the A6BBB abutment lab analogues
	Parts labelled BCP1 for the CCP1 abutment lab analogues.
Digital Laboratory Analogues	Parts labelled LAD-DCx for where x denotes the corresponding platform size for the Deep Conical digital lab analogues.
	Parts labelled LAD-DCRx for where x denotes the corresponding platform size for the Deep Conical digital lab analogues and R indicates a reduced Deep Conical connection.
	Parts labelled LAD-L-x for where x denotes the corresponding platform size for the for the TRI-NEX digital lab analogues.
	Parts labelled LAD-3M, LAD-M, LAD-M-P45 and LAD-Z for the Internal Hex digital lab analogues.
	Parts labelled LAD-MC and LAD-MCW for the Compact Conical Abutment digital lab analogues.
	Parts labelled LAD-ITS6 and LAD-ITS for the Internal Octagon (IT) digital lab analogues.
	Parts labelled LAD-IP, LAD-IBN, LAD-IB, LAD-BA and LAD-BBB for the External Hex digital lab analogues.
	Parts labelled LAD-SP and LAD-SP-PM for the Single Platform (SP1) implant digital lab analogues.
Luting Screws	Parts labelled SP16LS, SB-16-S, and PA-xx-LS where xx denotes the corresponding implant range or passive abutment screw.
Pin Gauges	Parts labelled CP-H, CP-L-18, CP-L-20 and WP-GC-S3.
Polishing Caps	Parts labelled PPC-xx where xx denotes the corresponding implant range.
Post Retention Components	Parts labelled CCP1 and CCP1-1.4.
Shoulder Analogues	Parts labelled LT6-7/10, LT-7, LT-10 and LT-20.
Titanium Bars	Parts labelled MW-BAR195, BAS1, BAS1L, BAS2, BAS3 and BAS1P.
Analogue Screws	Parts labelled LAD-S.
Try-in Abutments	Parts labelled TRY-SYN-xx and TRY-SYN6-xx where xx denotes the corresponding abutment angulation aligned to either the octagon flat or corner.
Reamers	Parts labelled LT18-xx where xx denotes the reamer diameter.
Lapping Tools	Parts labelled I-LTNP.
O-Rings	Parts labelled OAB-R.
Plastic Spacers and Clips	Parts labelled PS1, PS2, PC1, PC2 and PC4.

Table B – Torque and Compatibility Tables for Southern Brass Screws

Screw Type	External Hex Screws			
	1.22 Hex	Slotted	Torque	Head Diameter
2 Series (M2)	BSH2* [TSH2; TSHZ2]	BSS2* [TSS2; TSSZ2; GSS2; GSSZ2]	5 Ncm	2.70 mm
3 Series (M2)	BSH3* [TSH3; TSHZ3]	BSS3* [TSS3; TSSZ3; GSS3; GSSZ3]		2.40 mm
INPI (M1.6)	BSH16* [TSH16]	-		2.40 mm
Piccolo (M1.6)	BS-P-16* [TS-P-16]	-		2.15 mm

Screw Type	DC (Deep Conical) Screws		
	1.22 Hex	Torque	Head Diameter
Standard Screws (M2)	BS-DC5-20* [TS-DC5-20]	5 Ncm	2.35 mm
Standard Screws (M1.6)	BS-DC4-16* [TS-DC4-16]		2.35 mm
Standard Screws (M1.4)	BS-DC3-14* [TS-DC3-14]		2.10 mm
Passive Abutment Screws (M2)	PA-DC5-20B* (anodised purple) [PA-DC5-20T]		2.80 mm

Passive Abutment Screws (M1.6)	PA-DC4-16B* (anodised blue) [PA-DC4-16T]		2.40 mm
Passive Abutment Screws (M1.4)	PA-DC3-14B* (anodised yellow) [PA-DC3-14T]		2.20 mm

Screw Type	Internal Hex (M-series and PROVATA®) Screws		
	1.27 Hex	Torque	Head Diameter
Standard Screws (M1.8)	BS-Z-18* [TS-Z-18]	5 Ncm	2.25 mm
Standard Screws (M1.6)	BS-Z-16* [TS-Z-16]		2.25 mm
Passive Abutment Screws (M1.8)	PA-M-18B* [PA-M-18T]		2.60 mm
Passive Abutment Screws (M1.6)	PA-M-16B* [PA-M-16T]		2.60 mm

Screw Type	SP1 (Single Platform) Screws		
	1.22 Hex	Torque	Head Diameter
Standard Screws (M1.6)	BS-SP-16* [TS-SP-16]	5 Ncm	2.00 mm
Passive Abutment Screws (M1.6)	PA-SP-16B* [PA-SP-16T]		2.40 mm

Screw Type	IT (Internal Octagon) Screws		
	Torx	Torque	Head Diameter
Standard Screws (M2)	BSIT2* [TSIT2 GSIT2]	5 Ncm	2.70 mm
Passive Abutment Screws (M2)	BS-IT-PA* [TS-IT-PA]		

Screw Type	TRI-NEX Screws		
	Unigrip	Torque	Head Diameter
Standard Screws (M2)	BS-L-20* [TS-L-20; TS-L-20C (coated); GS-L-20]	5 Ncm	2.50 mm
Standard Screws (M1.8)	BS-L-18* [TS-L-18; TS-L-18C (coated); GS-L-18]		2.50 mm
Passive Abutment Screws (M2)	PA-L-20B* [PA-L-20 ;PA-L-20G]		2.70mm
Passive Abutment Screws (M1.8)	PA-L-18B* [PA-L-18; PA-L-18G]		2.70mm

Screw Type	Compact Conical Abutment Screws			
	1.22 Hex	Slotted	Torque	Head Diameter
Abutment Level Screws	BSH1* [TSH1; GSH1]	BSS1* [TSS1; GSS1]	5 Ncm	2.25 mm

***NOTE:** The brass screws are blackened (unless otherwise indicated) and intended for Laboratory Use Only – the torque values indicated in Table B are applicable to the brass screws only. For torque values of the corresponding prosthetic screws, refer to CAT-8068 – INSTRUCTIONS FOR USE: Southern Implants Prosthetic Screws

Clinical benefits

The Laboratory Components are used as part of dental implant therapy, a procedure in which a patient's missing teeth and/or crowns are restored. As a result, the benefits of the Laboratory Components mirror those of dental implant therapy in general, including improved chewing function, speech, aesthetics, and patient psychological wellbeing.

Storage, cleaning and sterilisation

This component is supplied non-sterile and is indicated for single use. If packaging is damaged do not use the product and contact your Southern representative or return to Southern Implants®. The devices must be stored in a dry place at room temperature and not exposed to direct sunlight. Incorrect storage may influence device characteristics.

Do not reuse components indicated for single-use only. Re-using these components may:

- damage on the surface or critical dimensions, which may result in performance and compatibility degradation.
- adds the risk of cross-patient infection and contamination if single-use items are reused.

Southern Implants® does not accept any responsibility for complications associated with reused single-use components.

Southern Implants® recommends one of the following procedures to sterilise the restorations and non-sterile single-use components prior to use:

1. prevacuum sterilisation method: steam sterilise the abutments at 132°C (270°F) at 180 - 220 kPa for 4 minutes. Dry for at least 20 minutes in the chamber. Only an approved wrap or pouch for steam sterilisation must be used.
2. for users in the USA: prevacuum sterilisation method: wrapped, steam sterilise at 135°C (275°F) at 180 - 220 kPa for 3 minutes. Dry for 20 minutes in the chamber. Use a wrap or pouch that is cleared for the indicated steam sterilisation cycle.

NOTE: users in the USA must ensure that the steriliser, wrap or pouch, and all steriliser accessories are cleared by the FDA, for the intended sterilisation cycle.

Contraindications

Since the Laboratory Components are required or utilized as part of dental implant surgery, the contraindications are not dissimilar to those of dental implant therapy as a whole. These contraindications include:

- patients medically unfit for oral surgical procedures
- cases where an adequate number of implants cannot be placed to provide full functional support of the prosthesis.
- patients under the age of 18 years.
- patients with poor bone quality.
- patients with blood disorders.
- presence of infection at the implant site.
- patients with vascular impairment.
- patients with uncontrolled diabetes.
- patients undergoing chronic high dose steroid therapy.
- patients with metabolic bone disease.
- patients with pure titanium, titanium alloy (Ti6Al4V), gold, palladium, polyether ether ketone (PEEK) or iridium allergies or hypersensitivities.

Other than the above, there are no contraindications unique to this system.

Warnings and precautions

THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING

- For the safe and effective use of dental implants, new technologies/systems and the Laboratory Components, it is strongly suggested that specialised training be undertaken. This training should include hands-on

methods to gain competency on the proper technique, biomechanical requirements and radiographic evaluations required for the specific system.

- Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone.
- Use of the device with incompatible or non-corresponding devices can result in device failure or poor performance.
- When handled intraorally, devices must be adequately secured against aspiration. Aspiration of products may lead to infection or physical injury.
- Electrosurgery should not be attempted around metal implants as they are conductive.
- The use of non-sterile items can lead to secondary infections of the tissue or transfer infectious diseases.
- Failure to perform appropriate cleaning, re-sterilisation and storage procedures as per the Instructions for Use document can result in device damage or secondary infection.
- Exceeding the number of recommended uses for reusable devices can result in device damage, secondary infection or patient harm,
- Blunt drills may cause damage to the bone which could compromise osseointegration.
- For short implants, clinicians should closely monitor patients for any of the following conditions: peri implant bone loss, changes to implant's response to percussion, or radiographic changes in bone to implant contact along the implant's length. If the implant shows mobility or greater than 50% bone loss, the implant should be evaluated for possible removal. If the clinicians choose a short implant, then clinicians should consider a two-stage surgical approach, splinting a short implant to an additional implant, and placement of the widest possible fixture. Allow longer periods for osseointegration and avoid immediate loading.

Note: Training should be undertaken by both new and experienced implant users before using a new system or attempting to do a new treatment method.

It is important that thorough screening of prospective implant candidates is performed. This screening should include at least:

- a comprehensive medical and dental history.
- visual and radiological inspection to determine adequate bone dimensions, anatomical landmarks, occlusal conditions and periodontal health. Failure to accurately estimate lengths of drills relative to radiographic measurements can result in permanent injury to nerves or other vital structures.
- the evaluation and consideration of bruxism and unfavourable jaw relations.

Proper preoperative planning with a good team approach between well trained surgeons, restorative dentists and lab technicians is essential for successful implant treatment.

For successful implant treatment, it is important to:

- minimise the trauma to the host tissue as this increases the potential for successful osseointegration.
- be aware and avoid damage to vital anatomical structures such as nerves, veins and arteries. Injury to these structures may cause serious complications, including, injury to the eye, nerve damage and excessive bleeding. It is essential to protect the infraorbital nerve.
- identify actual measurements relative to the radiographic data as failing to do so could lead to complications.

Special care should be taken when treating patients with local or systemic risk factors that could affect the healing of the bone and soft tissue or otherwise increase the severity of the side effects, risk of complications and/or implant failure. These factors include:

- smoking.
- history of periodontal disease.
- poor oral hygiene.
- history of orofacial radiotherapy.**
- bruxism and unfavourable jaw relations.

**** The potential for implant failure and other complications increases when implants are placed in irradiated bone as radiotherapy can result in progressive fibrosis of vessels and soft tissue (i.e., osteoradionecrosis), leading to diminished healing capacity. Implant treatment of irradiated patients is dependent upon issues like the timing of implant placement in**

relation to the radiation therapy, anatomic site chosen for implant placement and radiation dosage at that site and consequent risk of osteoradionecrosis.

The responsibility for proper patient selection, adequate training and experience in implant placement and the provision of the appropriate information required for informed consent rests with the practitioner.

Should the device not operate as intended, it must be reported to the manufacturer of the device. The contact information or the manufacturer of this device to report a change in performance is: sicomplaints@southernimplants.com.

Side effects

Since the device is required or utilized as part of dental implant surgery, the expected side effects are not dissimilar to those of dental implant therapy as a whole. Common side effects of dental implant therapy include pain, local inflammation, phonetic difficulties and gingival inflammation. Other, less common side effects or complications to dental implant therapy include, but are not limited to: (1) wound dehiscence; (2) peri-implantitis; (3) transient weakness, numbness, and/or pain associated with mild nerve injury; (4) fat emboli formation; (5) marginal bone loss within acceptable limits; (6) allergic reaction(s) to implant and/or abutment material; (7) unspecified infection; (8) implant failure due to insufficient levels of osseointegration; (9) breakage of the implant and/or abutment; (10) loosening of the abutment screw and/or retaining screw; (11) complications requiring revision of the dental implant; (12) nerve damage resulting in permanent weakness, numbness, or pain; (13) loosening of the implant requiring revision surgery; (14) perforation of the maxillary sinus; (15) perforation of the labial and/or lingual plates and (16) bone loss possibly resulting in revision or removal of the implant.

Precaution: maintaining sterility protocol

Implants are packaged as follows:

1. An outer package consisting of a rigid, clear box which acts as protection for the inner package.
2. The inner package consisting of a blister pack (clear plastic-formed blister base with a TYVEK "peel-back" lid).
3. Within the inner package, there is a hollow tube which contains one implant suspended from a titanium ring, this ensures the implant never touches the inside of the plastic tube.
4. Labelling information is located on the surface of the peel-back lid and on the outside of the rigid box.

Care must be taken to maintain the sterility of the implant by proper opening of the packaging and handling of the implant.

1. Open the implant package in the non-sterile field, with non-sterile gloves, tear the address label to open the box.
2. With non-sterile gloves, remove the inner blister pack. Do not place the plastic box or blister pack-lid onto the sterile field. The contents of this inner package are sterile.
3. The sealed blister is to be opened by an assistant (with nonsterile gloves): remove the TYVEK lid and drop or place the sterile tube onto the sterile field, open the tube cap and attach the implant placement tool onto the implant and carefully remove from the sterile tube. Do not touch the sterile implant.

Other sterile components are packed in a peel pouch or blister base with a "peel-back" lid. Labelling information is located on the bottom half of the pouch, inside the packet or on the surface of the peel-back lid. Sterility is assured unless the pouch is damaged or opened. Non-sterile components are supplied clean but not sterile in a peel pouch or blister base with peel-back lid. Labelling information is located on the bottom half of the pouch or on the surface of the peel-back lid.

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established.

The contact information for the manufacturer of this device to report a serious incident is as follows: sicomplaints@southernimplants.com.

Materials

Material type Brass, Delrin (Acetal) POM, Gold, Commercially Pure Titanium (Grade 2, 3 and 4), ARCAP, Stainless Steel (SAE 316/ 303, DIN 1.4197), Aluminium, Medical Grade PEEK, Standard PEEK, HOSTAFORM, Chrome Cobalt, Silicone

Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

Summary of Safety and Clinical Performance (SSCP)

As required by the European Medical Device Regulation (MDR; EU2017/745), a Summary of Safety and Clinical Performance (SSCP) is available for perusal with regard to Southern Implants® product ranges.

The relevant SSCP can be accessed at <https://ec.europa.eu/tools/eudamed>.

NOTE: the above website will be available upon the launch of the European Database on Medical Devices (EUDAMED).

Disclaimer of liability

This product is part of the Southern Implants® product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants® product range and take full responsibility for the correct indications and use of this product. Southern Implants® does not assume liability for damage due to incorrect use. Please note that some Southern Implants® products may not be cleared or released for sale in all markets.

Basic UDI

Product	Basic-UDI Number
Basic-UDI for Laboratory Components	6009544038799L

Related literature and catalogues

CAT-2001 – Prosthetic and Laboratory Instructions Manual
 CAT-2004 – Tri-Nex® Implants Product Catalogue
 CAT-2005 – IT Implants Product Catalogue
 CAT-2010 – Osseointegrated Fixtures Catalogue
 CAT-2020 – External Hex Implants Product Catalogue
 CAT-2042 – Deep Conical Implants Product Catalogue
 CAT-2043 – Internal Hex Implants Product Catalogue
 CAT-2060 – PROVATA® Implants Product Catalogue
 CAT-2069 – INVERTA® Implants Product Catalogue
 CAT-2070 – Zygomatic Implants Product Catalogue
 CAT-2087M – ILZ Mini Implants Product Catalogue
 CAT-2093 – Single Platform (SP1) Implants Product Catalogue

Symbols and warnings

											
Manufacturer: Southern Implants® 1 Albert Rd, P.O Box 605 IRENE, 0062, South Africa. Tel: +27 12 667 1046	CE mark 2797	Prescription device*	Sterilised using irradiation	Non-sterile	Use by date (mm-yy)	Do not reuse	Do not resterilise	Catalogue number	Batch code	Medical device	Authorised representative in the European Community
											
	Authorised representative for Switzerland	Date of manufacture	Magnetic Resonance conditional	Magnetic Resonance safe	Single sterile barrier system with protective packaging inside	Single sterile barrier system	Consult instruction for use	Caution	Keep away from sunlight	Do not use if package is damaged	

* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

All rights reserved. Southern Implants®, the Southern Implants® logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants®. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale. It is the responsibility of the clinician to inspect the symbols that appear on the packaging of the product in use.

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