English	INSTRUCTIONS FOR USE: Southern Implants® Laboratory Components
Español	INSTRUCCIONES DE USO: Southern Implants® Laboratory Components
Italiano	ISTRUZIONI PER L'USO: Southern Implants® Componenti di Laboratorio
Français	MODE D'EMPLOI : Composants de laboratoire de Southern Implants®
Deutsch	GEBRAUCHSANWEISUNG: Southern Implants® Labor-Komponenten
Português	INSTRUÇÕES DE UTILIZAÇÃO: Componentes de laboratório Southern Implants®
Ελληνικά	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ: Εξαρτήματα εργαστηρίου της Southern Implants®
Svenska	BRUKSANVISNING: Southern Implants® Laboratoriekomponenter
Română	INSTRUCȚIUNI DE UTILIZARE: Componente de laborator Southern Implants®
Türkçe	KULLANMA TALİMATI: Southern Implants® Laboratuvar Bileşenleri

South Africa - Headquarters: 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

Subsidiaries

Australia
Southern Implants Australia
T: +61-(0)-8-9466-2627
E: info@southernimplants.com.au

Spain and Portugal
Southern Implants Iberica
T: +34 935 053 507
E: info@southernimplants.es

United Kingdom and Ireland
Southern Implants UK
T: +44-20 8059 4490
E: info@southernimplants.co.uk

USA and Canada
Southern Implants North America Inc.
T: +1-561-472-0990
E: customercare@southernimplants.com

EC REP

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

CH REP

MedEnvoy Switzerland: Gotthardstrasse 28, 6302 Zug, Switzerland

Description

Brass Screws

The brass screws are used in the dental laboratory as a cost-effective alternative to the final prosthetic screw. The brass screw is intended to allow the laboratory technician to replicate the seating of the restoration onto a laboratory analogue and ensure the fit of the final prosthetic screw through the screw access hole of the restoration.

Cylinders

The plastic cylinders are supplied with Gold abutments and Passive abutments and intended for use in the laboratory to maintain the screw access hole and, screw seat (passive abutments only) during waxing up procedures of a castable restoration.

Plastic Sleeves

The plastic sleeves are intended for use in the laboratory as an aid in waxing up a restoration, during conventional casting procedures.

Gold Bars

The gold bar is connected to two implants placed in the anterior to provide stability for dentures. Retention clips are placed into the denture that will clip over the gold bar to aid in stabilizing the denture.

Laboratory Analogues

The laboratory analogues are used to replicate the location, orientation and connection of the prosthetic interface of an implanted dental device, namely the implant or abutment, in a stone or printed model during fabrication processes.

Digital Laboratory Analogues

The digital laboratory analogues are used as dental implant replicas to be inserted into a 3D printed or milled model in order to duplicate the location, orientation and restorative platform of the implant placed in the mouth during CAD/CAM procedures.

Luting Screws

The luting screws are used to clamp the interfacial component of a passive abutment onto the laboratory analogue during the process of luting the casting onto the interfacial component. The luting screw maintains the screw access hole and creates a chimney in the resulting restoration.

Pin Gauges

The pin gauges are devices used to keep the screw channel of a restoration or crown open when packing porcelain around the pin, or when waxing up a restoration around the pin.

Polishing Caps

The polishing caps are used to hold a restoration after cementing has taken place, allowing the dental professional to polish off any excess cement.

Post Retention Components

The retention posts are used during the fabrication of removable dentures and are soldered into a bar for retention of a removable overdenture.

Shoulder Analogues

The shoulder analogues are used to replicate the location, orientation of a shouldered abutment and aid in impression taking.

Titanium Bars

The titanium bars are used to fabricate dentures by providing stability to dentures as they connect multiple implants together and provide retention points where retention posts are soldered on.

Analogue Screws

The analogue screws are used to secure the digital laboratory analogues in place within a 3D printed or milled model.

Analogue Try-in Abutments

The try-in abutments are used with implants and laboratory analogues to confirm the required orientation of abutments, prior to selecting the desired abutment for the restoration.

Reamers and Lapping Tools

The reamers are used to ensure the passive abutment screw access hole diameter is correct and that the screw seat is flat. The lapping tools are used with the castable UCLA hexed sleeve and UCLA non-hexed sleeve to improve the fitting surface of the cast abutment.

O-Rings

The O-rings are used as a processing aid during laboratory procedures to ensure that acrylic does not enter the retentive cap of the denture during curing by sealing the interface. The O-ring is then removed and discarded.

Plastic clips and spacers

The plastic spacers are used with cementable abutments to allow for a screw retained restoration to be made. The plastic clips are laboratory clips which interface with a ball abutment.

Anaplastology Components

The Anaplastology components are used in the laboratory as an aid in prosthetic restoration by anaplastologists.

Plastic Access Maintainers

The plastic access maintainers are used to maintain the screw chimney and screw access hole of a prosthetic component during various restoration manufacturing processes.

Intended use

The Southern Implants® Laboratory Components are intended to aid in the laboratory procedures associated with implant therapy.

Indications for use

The Southern Implants® Lab Components are indicated for the laboratory procedures associated with the functional and aesthetic restoration of individuals indicated for implant-retained prostheses.

Intended user

Dental Technicians and other appropriately trained and/or experienced medical professionals.

Intended environment

The devices are intended to be used in a dental laboratory as part of the restoration design and manufacture.

Intended patient population

The Southern Implants® Lab Components are not intended to come into contact with the patient and, therefore, do not have an intended patient population.

Compatibility information

The Laboratory Components can be identified by specific abbreviations in the product codes. Device identifiers are summarised in Table A. Table B summarises the brass screws available for each Southern range (corresponding prosthetic screws indicated in square brackets []) and the recommended torque values of the Southern screws.

Table A – Laboratory Component device codes

Implant connection type	Compatible device
Plastic Cylinders	Parts labelled ITS6-PC1 and ITS-PC1 for Internal Octagon (IT) engaging cylinders.
	Parts labelled ITS6-PC1ne and ITS-PC1ne for Internal Octagon (IT) non-engaging cylinders.
	Parts labelled IT-PE-x and IT6-PE-x where xx denotes the post height for Internal Octagon (IT) engaging cylinders.
	Parts labelled IT-PE-NE and IT6-PE-NE where xx denotes the post height for Internal Octagon (IT) non-engaging cylinders.
	Parts labelled PC-EL-xx and PC-NL-xx where xx denotes the corresponding platform diameter for TRI-NEX engaging and non-engaging cylinders respectively.
	Parts labelled PC-M-E and PC-M-NE for Internal Hex engaging and non-engaging cylinders respectively.
	Parts labelled SB1, SB5, SBA1, SBA5, SBBB1, SBBB5, SBN1, SBN5, SSD4 for External Hex cylinders
Plastic Sleeves	Parts labelled ITS6-PA-P and ITS-PA-P for Internal Octagon (IT) Passive Abutment engaging plastic sleeves.
	Parts labelled ITS6-PA-NEP and ITS-PA-NEP for Internal Octagon (IT) Passive Abutment non-engaging plastic sleeves.
	Parts labelled PA-DCx-P and PA-NDCx-P where x denotes the corresponding platform size for the Deep Conical Passive Abutment engaging and non-engaging plastic sleeves respectively.
	Parts labelled PA-EL-xP and PA-NL-xP where x denotes the corresponding platform size for the TRI-NEX Passive Abutment engaging and non-engaging plastic sleeves respectively.
	Parts labelled PA-3EM-SP, PA-EM-SP, PA-EM-WP, PA-EM-SP45-P and PA-EZ-P for the Internal Hex Passive Abutment engaging plastic sleeves.
	Parts labelled PA-3NM-SP, PA-NM-SP, PA-NM-WP and PA-NZ-P for the Internal Hex Passive Abutment non-engaging plastic sleeves.
	Parts labelled PA-SP-P and PA-SP-PM-P for the Single Platform (SP1) Passive Abutment engaging plastic sleeves.
	Parts labelled PA-NSP-P for the Single Platform (SP1) Passive Abutment non-engaging plastic sleeves.
	Parts labelled PA-MC-48P and PMC-1 for the Ø4.8 Compact Conical Abutment interface.
	Parts labelled PA-MC-60P and PMCW-1 for the Ø6.0 Compact Conical Abutment interface.
	Parts labelled SP16PH, SB16P, SBA16PH, SBBB16PH, SBN16PH, SBNV16PH, SBV16PH and SMAX9PH for the External Hex Passive Abutment engaging plastic sleeves.
	Parts labelled SP-17P, SB17P, SBA17P, SBBB17P, SBN17P, SBV-17-9 and SMAX9P for the External Hex Passive Abutment non-engaging plastic sleeves.
	Parts labelled SB-DBAN-1 and SB-DBN-1 for the DBAN and DBN abutments.
	Parts labelled SLZ for ILZ implants.
	Parts labelled SI30-SLEEVE for SI30 abutment.
	Parts labelled SMP-1-S for the SMP1 abutment.
	Parts labelled GC3-P, GC-30P, GCBA-P, GCBBB-P, GCBN-P, GCB-P and GCMAx-P for the External Hex Gold Abutment plastic sleeves.
	Parts labelled GC-DCx-P where x denotes the corresponding platform size for the for the Deep Conical Gold Abutment plastic sleeves.
	Parts labelled GC-M-P and GC-Z-P for the Internal Hex Gold Abutment plastic sleeves.
	Parts labelled GC-PL-45 for the TRI-NEX Gold Abutment plastic sleeves.
	Parts labelled ITS6-GC1-P and ITS-GC1-P for the Internal Octagon (IT) Gold Abutment plastic sleeves.
	Parts labelled SB-2-6-P and SBBB-2-6-P for the Compact Conical Gold Abutment plastic sleeves.
	Parts labelled CC-L-35PC and CC-L-43PC for TRI-NEX Chrome Cobalt cylinder plastic sleeves.
	Parts labelled CC-PS-40 for External Hex Chrome Cobalt cylinder plastic sleeves.
	Parts labelled ITS6-CC1-PC and ITS-CC1-PC for Internal Octagon (IT) Chrome Cobalt cylinder plastic sleeves.
	Parts labelled ITS6-TA and ITS-TA for TSA6/ TSAF plastic sleeves.
Gold Bars	Parts labelled GD1 and GDC1.
Laboratory Analogues	Parts labelled LA-DCx where x denotes the corresponding platform size for the Deep Conical lab analogues.
	Parts labelled LA-L-x where x denotes the corresponding platform size for the TRI-NEX lab analogues.
	Parts labelled LA-3M, LA-M, LA-M-P45 and LA-Z for the Internal Hex lab analogues.
	Parts labelled LA-ILZ and LSLZ1 for the ILZ lab analogues.
	Parts labelled LITS4 and LITS6 for the Internal Octagon (IT) lab analogues.

	Parts labelled LP12, LSN12, LS12, LSA12, LBBBS12 and LMAX9 for the External Hex lab analogues.
	Parts labelled LSMC1, LSMW1, LSA7 and LS7 for the Compact Conical Abutment lab analogues.
	Parts labelled LAFx and LAF6-x where x denotes the corresponding post height for the TSAF lab analogues.
	Parts labelled LOAB-1, LOAB-H, AO2, PP2 and BP2 for Ball/ Overdenture abutment lab analogues.
	Parts labelled LSDBAN4 and LSDBN-4 for the DBAN and DBN abutment lab analogues.
	Parts labelled LSD4 and LSI3 for the Calcitek integral implant lab analogues.
	Parts labelled LSI1 for the Shoulder abutment lab analogues.
	Parts labelled LS1 for the Standard abutment lab analogues.
	Parts labelled LS6BBB1 for the A6BBB abutment lab analogues
	Parts labelled BCP1 for the CCP1 abutment lab analogues.
Digital Laboratory Analogues	Parts labelled LAD-DCx for where x denotes the corresponding platform size for the Deep Conical digital lab analogues. Parts labelled LAD-DCRx for where x denotes the corresponding platform size for the Deep Conical digital lab analogues and R indicates a reduced Deep Conical connection.
	Parts labelled LAD-L-x for where x denotes the corresponding platform size for the for the TRI-NEX digital lab analogues.
	Parts labelled LAD-3M, LAD-M, LAD-M-P45 and LAD-Z for the Internal Hex digital lab analogues.
	Parts labelled LAD-MC and LAD-MCW for the Compact Conical Abutment digital lab analogues.
	Parts labelled LAD-ITS6 and LAD-ITS for the Internal Octagon (IT) digital lab analogues.
	Parts labelled LAD-IP, LAD-IBN, LAD-IB, LAD-BA and LAD-BBB for the External Hex digital lab analogues.
	Parts labelled LAD-SP and LAD-SP-PM for the Single Platform (SP1) implant digital lab analogues.
Luting Screws	Parts labelled SP16LS, SB-16-S, and PA-xx-LS where xx denotes the corresponding implant range or passive abutment screw.
Pin Gauges	Parts labelled CP-H, CP-L-18, CP-L-20 and WP-GC-S3.
Polishing Caps	Parts labelled PPC-xx where xx denotes the corresponding implant range.
Post Retention Components	Parts labelled CCP1 and CCP1-1.4.
Shoulder Analogues	Parts labelled LT6-7/10, LT-7, LT-10 and LT-20.
Titanium Bars	Parts labelled MW-BAR195, BAS1, BAS1L, BAS2, BAS3 and BAS1P.
Analogue Screws	Parts labelled LAD-S.
Try-in Abutments	Parts labelled TRY-SYN-xx and TRY-SYN6-xx where xx denotes the corresponding abutment angulation aligned to either the octagon flat or corner.
Reamers	Parts labelled LT18-xx where xx denotes the reamer diameter.
Lapping Tools	Parts labelled I-LTNP, I-LT6.
O-Rings	Parts labelled OAB-R.
Plastic Spacers and Clips	Parts labelled PS1, PS2, PC1, PC2 and PC4.
Anaplastology Components	Parts labelled V3-GC-93184, V2-GC-90772 and V3-PDA-6.5/14.
Plastic Access Maintainers	Parts labelled PM22.

Table B – Torque and Compatibility Tables for Southern Brass Screws

Screw Type	External Hex Screws			
	1.22 Hex	Slotted	Torque	Head Diameter
2 Series (M2)	BSH2* [TSH2; TSHZ2]	BSS2* [TSS2; TSSZ2; GSS2; GSSZ2]	5 Ncm	2.70 mm
3 Series (M2)	BSH3* [TSH3; TSHZ3]	BSS3* [TSS3; TSSZ3; GSS3; GSSZ3]		2.40 mm
INPI (M1.6)	BSH16* [TSH16]	-		2.40 mm
Piccolo (M1.6)	BS-P-16* [TS-P-16]	-		2.15 mm

Screw Type	Deep Conical (DC) Screws		
	1.22 Hex	Torque	Head Diameter
Standard Screws (M2)	BS-DC5-20* [TS-DC5-20]	5 Ncm	2.35 mm
Standard Screws (M1.6)	BS-DC4-16* [TS-DC4-16]		2.35 mm
Standard Screws (M1.4)	BS-DC3-14* [TS-DC3-14]		2.10 mm
Passive Abutment Screws (M2)	PA-DC5-20B* (anodised purple) [PA-DC5-20T]		2.80 mm
Passive Abutment Screws (M1.6)	PA-DC4-16B* (anodised blue) [PA-DC4-16T]		2.40 mm
Passive Abutment Screws (M1.4)	PA-DC3-14B* (anodised yellow) [PA-DC3-14T]		2.20 mm

Screw Type	Internal Hex (M-Series and PROVATA®) Screws		
	1.27 Hex	Torque	Head Diameter
Standard Screws (M1.8)	BS-Z-18* [TS-Z-18]	5 Ncm	2.25 mm
Standard Screws (M1.6)	BS-Z-16* [TS-Z-16]		2.25 mm
Passive Abutment Screws (M1.8)	PA-M-18B* [PA-M-18T]		2.60 mm
Passive Abutment Screws (M1.6)	PA-M-16B* [PA-M-16T]		2.60 mm

Screw Type	Single Platform (SP1) Screws		
	1.22 Hex	Torque	Head Diameter
Standard Screws (M1.6)	BS-SP-16* [TS-SP-16]	5 Ncm	2.00 mm
Passive Abutment Screws (M1.6)	PA-SP-16B* [PA-SP-16T]		2.40 mm

Screw Type	Internal Octagon (IT) Screws		
	Torx	Torque	Head Diameter
Standard Screws (M2)	BSIT2* [TSIT2 GSIT2]	5 Ncm	2.70 mm
Passive Abutment Screws (M2)	BS-IT-PA* [TS-IT-PA]		

Screw Type	TRI-NEX Screws		
	Unigrip	Torque	Head Diameter
Standard Screws (M2)	BS-L-20* [TS-L-20; TS-L-20C (coated); GS-L-20]	5 Ncm	2.50 mm
Standard Screws (M1.8)	BS-L-18* [TS-L-18; TS-L-18C (coated); GS-L-18]		2.50 mm
Passive Abutment Screws (M2)	PA-L-20B* [PA-L-20; PA-L-20G]		2.70mm
Passive Abutment Screws (M1.8)	PA-L-18B* [PA-L-18; PA-L-18G]		2.70mm

Screw Type	Compact Conical Abutment Screws			
	1.22 Hex	Slotted	Torque	Head Diameter
Abutment Level Screws	BSH1* [TSH1; GSH1]	BSS1* [TSS1; GSS1]	5 Ncm	2.25 mm

***NOTE:** The brass screws are blackened (unless otherwise indicated) and intended for Laboratory Use Only – the torque values indicated in Table B are applicable to the brass screws only. For torque values of the corresponding prosthetic screws, refer to CAT-8068 – INSTRUCTIONS FOR USE: Southern Implants Prosthetic Screws.

Clinical performance

The Southern Implants® Lab Components are intended to aid in the laboratory procedures associated with implant therapy. Consequently, the clinical performance of these devices is primarily defined by the success of the treatment. This performance can be quantitatively assessed through prosthetic survival/success, which serves as an important indicator of their effectiveness and success in clinical practice.

Clinical benefits

The Southern Implants® Lab Components play a crucial role in ensuring the successful outcome of laboratory-associated prosthetic procedures, though they do not provide a direct clinical benefit. However, the clinical benefits of the associated treatment, target implant, and/or prosthetic device(s) should be acknowledged, with relevant documentation reviewed as necessary.

Storage, cleaning and sterilisation

This component is supplied non-sterile and is indicated for single use. If packaging is damaged do not use the product and contact your Southern representative or return to Southern Implants®. The devices must be stored in a dry place at room temperature and not exposed to direct sunlight. Incorrect storage may influence device characteristics.

Do not reuse components indicated for single-use only. Re-using these components may:

- damage on the surface or critical dimensions, which may result in performance and compatibility degradation.
- adds the risk of cross-patient infection and contamination if single-use items are reused.

Southern Implants® does not accept any responsibility for complications associated with reused single-use components.

Contraindications

The Southern Implants® Lab Components do not come into contact with the patient and, therefore, no specific patient-related contraindications apply.

However, the contraindications of all device groups used as part of the specific treatment or procedure apply. Therefore, the contraindications of the systems/medical devices utilized as part of implant surgery/therapy should be noted and the relevant documents consulted.

Warnings and precautions

IMPORTANT NOTICE: THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING

- To ensure the safe and effective use of dental implants, new technologies/systems and the Lab Component devices, it is strongly recommended that specialised training be undertaken. This training should include hands-on methods to gain competency on the proper technique, the system's biomechanical requirements and radiographic evaluations required for the specific system.
- Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone.
- Use of the device with incompatible or non-corresponding devices can result in the poor performance or failure of the device.
- When handling devices intraorally, it is imperative that they are adequately secured to prevent aspiration, as aspiration of products may lead to infection or physical injury.

- The use of non-sterile items can lead to secondary infections of the tissue or transfer infectious diseases.
- Failure to adhere to appropriate cleaning, re-sterilization, and storage procedures as outlined in the Instructions for Use can result in device damage, secondary infections, or patient harm.
- Exceeding the number of recommended uses for reusable devices can result in device damage, secondary infection or patient harm.
- The use of blunt drills may cause damage to the bone, potentially compromising osseointegration.

It is crucial to emphasize that training should be undertaken by both new and experienced implant users prior to utilizing a new system or attempting a new treatment method.

Patient Selection and Preoperative Planning

A comprehensive patient selection process and meticulous preoperative planning are essential for successful implant treatment. This process should involve consultation among a multidisciplinary team, including well-trained surgeons, restorative dentists, and laboratory technicians.

Patient screening should include, at a minimum, a thorough medical and dental history, as well as visual and radiological inspections to assess the presence of adequate bone dimensions, the positioning of anatomical landmarks, the presence of unfavourable occlusal conditions, and the periodontal health status of the patient.

For successful implant treatment, it is important to:

1. Minimize trauma to the host tissue, thereby increasing the potential for successful osseointegration.
2. Accurately identify measurements relative to radiographic data, as failure to do so may lead to complications.
3. Be vigilant in avoiding damage to vital anatomical structures, such as nerves, veins, and arteries. Injury to these structures may result in serious complications, including ocular injury, nerve damage, and excessive bleeding.

The responsibility for proper patient selection, adequate training and experience in implant placement, and the provision of appropriate information required for informed consent rests with the practitioner. By combining thorough screening of prospective implant candidates with a practitioner possessing a high level of competence in the use of the system, the potential for complications and severe side effects can be significantly reduced.

High Risk Patients

Special care should be exercised when treating patients with local or systemic risk factors that may adversely affect the healing of bone and soft tissue or otherwise increase the severity of side effects, the risk of complications, and/or the likelihood of implant failure. Such factors include:

- poor oral hygiene
- history of smoking/vaping/tobacco use
- history of periodontal disease
- history of orofacial radiotherapy**
- bruxism and unfavourable jaw relations
- use of chronic medications that may delay healing or increase the risk of complications including, but not limited to, chronic steroid therapy, anti-coagulant therapy, TNF- α blockers, bisphosphonate and cyclosporin

*** The potential for implant failure and other complications increases when implants are placed in irradiated bone, as radiotherapy can lead to progressive fibrosis of blood vessels and soft tissue (i.e., osteoradionecrosis), resulting in diminished healing capacity. Contributing factors to this increased risk include the timing of implant placement in relation to radiation therapy, the proximity of radiation exposure to the implant site, and the radiation dosage at that site.*

Should the device not operate as intended, it must be reported to the manufacturer of the device. The contact information or the manufacturer of this device to report a change in performance is: sicomplaints@southernimplants.com.

Side effects

The Southern Implants® Lab Components do not come into contact with the patient and, therefore, no specific patient-related side effects apply.

However, the side effects of all device groups used as part of the specific treatment or procedure apply. Therefore, the side effects of the systems/medical devices utilized as part of implant surgery/therapy should be noted and the relevant documents consulted.

Precaution: maintaining sterility protocol

The devices are supplied clean but not sterile in a peel pouch or blister base with peel-back lid. Labelling information is located on the bottom half of the pouch or on the surface of the peel-back lid.

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established.

The contact information for the manufacturer of this device to report a serious incident is as follows:

sicomplaints@southernimplants.com.

Materials

Material type: Brass, Delrin (Acetal) POM, Gold Alloy, Commercially Pure Titanium Grade 2, 3 or 4 according to ASTM F67 and ISO 5832-2, Titanium Alloy Ti-6Al-4V (90% titanium, 6% aluminium, 4% vanadium) according to ASTM F136 and ISO 5832-2, ARCAP, Stainless Steel (SAE 316/ 303, DIN 1.4197), Aluminium, Medical Grade PEEK, Standard PEEK, HOSTAFORM, Chrome Cobalt, Silicone

Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

Summary of Safety and Clinical Performance (SSCP)

As required by the European Medical Device Regulation (MDR; EU2017/745), a Summary of Safety and Clinical Performance (SSCP) is available for perusal with regard to Southern Implants® product ranges.

The relevant SSCP can be accessed at <https://ec.europa.eu/tools/eudamed>.

NOTE: the above website will be available upon the launch of the European Database on Medical Devices (EUDAMED).

Disclaimer of liability

This product is part of the Southern Implants® product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants® product range and take full responsibility for the correct indications and use of this product. Southern Implants® does not assume liability for damage due to incorrect use. Please note that some Southern Implants® products may not be cleared or released for sale in all markets.


Basic UDI

Product	Basic-UDI Number
Basic-UDI for Laboratory Components	6009544038799L

Related literature and catalogues

CAT-2001 – Prosthetic and Laboratory Instructions Manual
 CAT-2004 – Tri-Nex® Implants Product Catalogue
 CAT-2005 – IT Implants Product Catalogue
 CAT-2010 – Osseointegrated Fixtures Catalogue
 CAT-2020 – External Hex Implants Product Catalogue
 CAT-2042 – Deep Conical Implants Product Catalogue
 CAT-2043 – Internal Hex Implants Product Catalogue
 CAT-2060 – PROVATA® Implants Product Catalogue
 CAT-2069 – Deep Conical INVERTA® Implants Product Catalogue
 CAT-2070 – Zygomatic Implants Product Catalogue
 CAT-2087M – ILZ Mini Implants Product Catalogue
 CAT-2088M – Nazalus Implants Product Catalogue and Surgical Manual
 CAT-2092 – Deep Conical Pterygoid Implants Product Catalogue
 CAT-2093 – Single Platform (SP1) Implants Product Catalogue
 CAT-2095 – External Hex INVERTA® Implants Product Catalogue
 CAT-2096 – External Hex Pterygoid Implants Product Catalogue

Symbols and warnings

* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

All rights reserved. Southern Implants®, the Southern Implants® logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants®. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale. It is the responsibility of the clinician to inspect the symbols that appear on the packaging of the product in use.

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