
English	INSTRUCTIONS FOR USE: Southern Implants® General Use Instruments
Español	INSTRUCCIONES DE USO: Instrumentos de Uso General Southern Implants®
Italiano	ISTRUZIONI PER L'USO: Southern Implants® Strumenti per Uso Generale
Français	MODE D'EMPLOI : Instruments à usage général Southern Implants®
Deutsch	GEBRAUCHSANWEISUNGEN: Southern Implants® Allgemeine Gebrauchsinstrumente
Português	INSTRUÇÕES DE UTILIZAÇÃO: Southern Implants® Instrumentos de uso geral
Ελληνικά	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ: Southern Implants® Εργαλεία γενικής χρήσης
Svenska	ANVISNINGAR FÖR ANVÄNDNING: Southern Implants® Instrument för Allmän Bruk

ARCHIVED

South Africa - Headquarters: 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

Subsidiaries

Australia

Southern Implants Australia
T: +61-(0)-8-9466-2627
E: info@southernimplants.com.au

Spain and Portugal

Southern Implants Iberica
T: +34 935 053 507
E: info@southernimplants.es

United Kingdom and Ireland

Southern Implants UK
T: +44-20-8899-6845 / 6 / 7
E: info@southernimplants.co.uk

USA and Canada

Southern Implants North America Inc.
T: +1-561-472-0990
E: customercare@southernimplants.com

Description

Dental implantation drilling template guide sleeves

A small tubular device intended to be inserted within the guide holes of a dental implantation drilling template (a device used to guide drilling in the jawbone during dental implant placement/fixation) to guide the drill and / or implant into position. Dental implantation drilling template guide sleeves are packaged in a peel pouch. These devices are supplied non-sterile and for single use only.

Dental implant internal thread shaping tools

A cylindrical rod with a spiralled cutting distal end, similar to a drill bit or a screw, intended to be used to shape and rectify the internal thread of a dental implant that has been implanted into bone. It is made of metal and is designed to be attached to a manual screwdriver. These devices are supplied non-sterile and for single use only.

Dental implant extractor

A manually-operated, dental device used to retrieve a dental implant from the oral cavity, typically due to damage (e.g., a broken collar) or malfunction of the dental implant. It is made of toughened, high-grade steel and has a cylindrical design with a coarse, left-handed flute. The flute has a steep pitch that spirals up the tapered working end and there is a hexagonal head at the proximal end that fits a torque wrench socket. The device is inserted into the implant which is to be removed and turned anticlockwise (counter clockwise). The anticlockwise torque exerted makes this device grip into the implant, which is then unscrewed as the device continues to turn anti-clockwise. This device are supplied sterile and for single use only.

Abutment retrieval tools

A screwdriver driven dental device used to retrieve a dental abutment, typically one with a short abutment post (e.g., a passive abutment) and conical connection (e.g., Single Platform SP1 connection), from a dental implant after prolonged loading when it is found that the abutment is jammed into the implant's connection. It is made of Titanium alloy (ASTM F136) and has a screw design with an elongated nib section intended to bottom out in the implant thread pilot hole and lift the abutment out of the implant without damaging the internal thread or connection of the implant. The abutment retrieval tool has external thread that corresponds with internal thread either manufactured into the abutment's screw access hole or tapped by the user. The device is inserted into the abutment and rotated clockwise until the seal between the abutment and implant's conical connection is broken and the abutment is dislodged and may be removed from the implant. These devices are supplied non-sterile and for single use only.

Intended use

The General Use Instruments are intended to aid in the treatment of partially or fully edentulous patients eligible for placement of one or more dental implants as a means of fixing a permanent or removable single crown, partial or full-arch dental prosthesis in the upper or lower jaw.

Indications for use

The General Use Instruments are indicated for use in the treatment of partial or full edentulousness, congenital tooth loss, or failing/problematic natural dentition, when a patient is dissatisfied with existing removable dentures, or as part of the revision/replacement of previous dental restorations.

Intended user

The intended user for this system includes Maxillofacial Surgeons, General Dentists, Orthodontists, Periodontist, Prosthodontists and other appropriately trained and experienced implant users.

Intended environment

This device is intended to be used in a clinical environment such as an operating theatre or a dentist consultation room.

Intended patient population

The intended patient population for dental implant therapy includes partially or fully edentulous patients eligible or, otherwise, not contraindicated for implant placement that require prosthetic restoration or revision of existing

restorations in the upper or lower jaw; where the planned restorations are to be fixed and comprise of single teeth, partial or full bridges.

Compatibility information

Southern Implants' implants should be restored with Southern Implants' components. In the Southern Implants' General Use Instruments range there are 4 product types.

Table A - Compatible

General Use Instrument type	Product Code	Description	Sterility Status
Dental implantation drilling template guide sleeves	I-GS-2513-L6	Guided surgery pin sleeve Ø2.80 outer diameter, Ø1.30 inner diameter, length 6mm	Non-sterile
	I-GS-3020-L4	Guided surgery sleeve Ø3.30 outer diameter, Ø2.0 inner diameter, length 4mm	Non-sterile
	I-GS-3020-L6	Guided surgery sleeve Ø3.30 outer diameter, Ø2.0 inner diameter, length 6mm	Non-sterile
	I-GS-5142-L4	Guided surgery sleeve Ø5.40 outer diameter, Ø4.15 inner diameter, length 4mm	Non-sterile
	I-GS-5142-L9	Guided surgery sleeve Ø5.40 outer diameter, Ø4.15 inner diameter, length 9mm	Non-sterile
	I-GS-6252-L4	Guided surgery sleeve Ø6.50 outer diameter, Ø5.15 inner diameter, length 4mm	Non-sterile
	I-GS-6252-L9	Guided surgery sleeve Ø6.50 outer diameter, Ø5.15 inner diameter, length 9mm	Non-sterile
	I-GSC-5142-L4	Guided surgery C-type sleeve Ø5.40 outer diameter, Ø4.20 inner diameter, length 4mm	Non-sterile
	I-GSC-5142-L9	Guided surgery C-type sleeve Ø5.40 outer diameter, Ø4.20 inner diameter, length 9mm	Non-sterile
	I-GSC-6252-L4	Guided surgery C-type sleeve Ø6.50 outer diameter, Ø5.20 inner diameter, length 4mm	Non-sterile
	I-GSC-6252-L9	Guided surgery C-type sleeve Ø6.50 outer diameter, Ø5.20 inner diameter, length 9mm	Non-sterile
	I-GSZ-4128-L8	Guided surgery sleeve for zygomatic implants Ø4.40 outer diameter, Ø2.80 inner diameter, length 8mm	Non-sterile
	I-GS-4542-CS	Guided surgery sleeve for CEREC Ø4.75 outer diameter, Ø4.15 inner diameter, length 4mm	Non-sterile
	I-GS-5342-CM	Guided surgery sleeve for CEREC Ø5.55 outer diameter, Ø4.15 inner diameter, length 4mm	Non-sterile
	I-GS-6052-CL	Guided surgery sleeve for CEREC Ø6.25 outer diameter, Ø5.15 inner diameter, length 4mm	Non-sterile
Dental implant internal thread shaping tools	I-M1.4-S	M1.4 thread tap, length short	Non-sterile
	I-M1.4	M1.4 thread tap, length long	Non-sterile
	I-M1.6-S	M1.6 thread tap, length short	Non-sterile
	I-M1.6	M1.6 thread tap, length long	Non-sterile
	I-M1.8-S	M1.8 thread tap, length short	Non-sterile
	I-M1.8	M1.8 thread tap, length long	Non-sterile
	I-M2-S	M2 thread tap, length short	Non-sterile
	I-M2	M2 thread tap, length long	Non-sterile
Dental implant extractor	I-IMP-REM	Implant extractor, wrench insert	Sterile
Abutment retrieval tools	I-ART-SP	Abutment retrieval tool/ screw to be used with the Single Platform SP1 passive abutments PA-SP & PA-NSP	Non-sterile

Clinical benefits

The General Use Instruments are used as part of dental implant therapy, a procedure in which a patient's missing teeth and/or crowns are restored. As a result, the benefits of the General Use Instruments mirror those of dental implant therapy in general, including improved chewing function, speech, aesthetics, and patient psychological wellbeing.

Storage, cleaning and sterilisation

APPLICABLE FOR STERILE COMPONENTS:

Some General Use Instruments are supplied sterile (sterilised by gamma irradiation) and intended for single-use prior to the expiration date (see packaging label). Sterility is assured unless the container or seal is damaged or opened. If packaging is damaged do not use the product and contact your Southern representative or return to Southern Implants®. The devices must be stored in a dry place at room temperature and not exposed to direct sunlight. Incorrect storage may influence device characteristics. Do not reuse components indicated for single-use only. Reusing these components may result in:

- damage to the surface or critical dimensions, which may result in performance and compatibility degradation.
- adds the risk of cross-patient infection and contamination if single-use items are reused.

Southern Implants® does not accept any responsibility for complications associated with reused single-use components.

APPLICABLE FOR NON-STERILE COMPONENTS:

Some General Use Instruments are supplied non-sterile and is indicated for single use. If packaging is damaged do not use the product and contact your Southern representative or return to Southern Implants®. The devices must be stored in a dry place at room temperature and not exposed to direct sunlight. Incorrect storage may influence device characteristics.

Do not reuse components indicated for single-use only. Re-using these components may:

- damage on the surface or critical dimensions, which may result in performance and compatibility degradation.
- adds the risk of cross-patient infection and contamination if single-use items are reused.

Southern Implants® does not accept any responsibility for complications associated with reused single-use components.

Southern Implants® recommends one of the following procedures to sterilise the restorations and non-sterile single-use components prior to use:

1. prevacuum sterilisation method: steam sterilise the abutments at 132°C (270°F) at 180 - 220 kPa for 4 minutes. Dry for at least 20 minutes in the chamber. Only an approved wrap or pouch for steam sterilisation must be used.
2. for users in the USA: prevacuum sterilisation method: wrapped, steam sterilise at 135°C (275°F) at 180 - 220 kPa for 3 minutes. Dry for 20 minutes in the chamber. Use a wrap or pouch that is cleared for the indicated steam sterilisation cycle.

NOTE: users in the USA must ensure that the steriliser, wrap or pouch, and all steriliser accessories are cleared by the FDA, for the intended sterilisation cycle.

Contraindications

Since the General Use Instruments are required or utilized as part of dental implant surgery, the complications for this device are not dissimilar to those of dental implant therapy as a whole. These contraindications include:

- patients medically unfit for oral surgical procedures
- cases where an adequate number of implants cannot be placed to provide full functional support of the prosthesis.
- patients under the age of 18 years.
- patients with poor bone quality.
- patients with blood disorders.
- presence of infection at the implant site.
- patients with vascular impairment.
- patients with uncontrolled diabetes.
- patients undergoing chronic high dose steroid therapy.
- patients with metabolic bone disease.
- patients with pure titanium, titanium alloy (Ti6Al4V), gold, palladium, polyether ether ketone (PEEK) or iridium allergies or hypersensitivities.

Other than the above, there are no contraindications unique to this system.

Warnings and precautions

THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING

- For the safe and effective use of dental implants, new technologies/systems and the <device name>, it is strongly suggested that specialised training be undertaken. This training should include hands-on methods to

gain competency on the proper technique, biomechanical requirements and radiographic evaluations required for the specific system.

- Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone.
- Use of the device with incompatible or non-corresponding devices can result in device failure or poor performance.
- When handled intraorally, devices must be adequately secured against aspiration. Aspiration of products may lead to infection or physical injury.
- Electrosurgery should not be attempted around metal implants as they are conductive.
- The use of non-sterile items can lead to secondary infections of the tissue or transfer infectious diseases.
- Failure to perform appropriate cleaning, resterilisation and storage procedures as per the Instructions for Use document can result in device damage or secondary infection.
- Blunt drills may cause damage to the bone which could compromise osseointegration.
- For short implants, clinicians should closely monitor patients for any of the following conditions: peri implant bone loss, changes to implant's response to percussion, or radiographic changes in bone to implant contact along the implant's length. If the implant shows mobility or greater than 50% bone loss, the implant should be evaluated for possible removal. If the clinicians choose a short implant, then clinicians should consider a two-stage surgical approach, splinting a short implant to an additional implant, and placement of the widest possible fixture. Allow longer periods for osseointegration and avoid immediate loading.

Note: Training should be undertaken by both new and experienced implant users before using a new system or attempting to do a new treatment method.

It is important that thorough screening of prospective implant candidates is performed. This screening should include at least:

- a comprehensive medical and dental history.
- visual and radiological inspection to determine adequate bone dimensions, anatomical landmarks, occlusal conditions and periodontal health. Failure to accurately estimate lengths of drills relative to radiographic measurements can result in permanent injury to nerves or other vital structures.
- the evaluation and consideration of bruxism and unfavourable jaw relations.

Proper preoperative planning with a good team approach between well trained surgeons, restorative dentists and lab technicians is essential for successful implant treatment.

For successful implant treatment, it is important to:

- minimise the trauma to the host tissue as this increases the potential for successful osseointegration.
- be aware and avoid damage to vital anatomical structures such as nerves, veins and arteries. Injury to these structures may cause serious complications, including, injury to the eye, nerve damage and excessive bleeding. It is essential to protect the infraorbital nerve.
- identify actual measurements relative to the radiographic data as failing to do so could lead to complications.

Special care should be taken when treating patients with local or systemic risk factors that could affect the healing of the bone and soft tissue or otherwise increase the severity of the side effects, risk of complications and/or implant failure. These factors include:

- smoking.
- history of periodontal disease.
- poor oral hygiene.
- history of orofacial radiotherapy.**
- bruxism and unfavourable jaw relations.

**** The potential for implant failure and other complications increases when implants are placed in irradiated bone as radiotherapy can result in progressive fibrosis of vessels and soft tissue (i.e., osteoradionecrosis), leading to diminished healing capacity. Implant treatment of irradiated patients is dependent upon issues like the timing of**

implant placement in relation to the radiation therapy, anatomic site chosen for implant placement and radiation dosage at that site and consequent risk of osteoradionecrosis.

The responsibility for proper patient selection, adequate training and experience in implant placement and the provision of the appropriate information required for informed consent rests with the practitioner.

Should the device not operate as intended, it must be reported to the manufacturer of the device. The contact information or the manufacturer of this device to report a change in performance is:

sicomplaints@southernimplants.com

Side effects

Since the device is required or utilized as part of dental implant surgery, the side effects of the device are not dissimilar to those of dental implant therapy as a whole. Common side effects of dental implant therapy include pain, inflammation, phonetic difficulties and gingival inflammation. Other, less common side effects or complications to dental implant therapy include, but are not limited to: (1) wound dehiscence; (2) peri-implantitis; (3) transient weakness, numbness, and/or pain associated with mild nerve injury; (4) fat emboli formation; (5) marginal bone loss within acceptable limits; (6) allergic reaction(s) to implant and/or abutment material; (7) unspecified infection; (8) implant failure due to insufficient levels of osseointegration; (9) breakage of the implant and/or abutment; (10) loosening of the abutment screw and/or retaining screw; (11) complications requiring revision of the dental implant; (12) nerve damage resulting in permanent weakness, numbness, or pain; (13) loosening of the implant requiring revision surgery; (14) perforation of the maxillary sinus; (15) perforation of the labial and/or lingual plates and (16) bone loss possibly resulting in revision or removal of the implant.

Precaution: maintaining sterility protocol

Implants are packaged as follows:

1. An outer package consisting of a rigid, clear box which acts as protection for the inner package.
2. The inner package consisting of a blister pack (clear plastic-formed blister base with a TYVEK "peel-back" lid).
3. Within the inner package, there is a hollow tube which contains one implant suspended from a titanium ring, this ensures the implant never touches the inside of the plastic tube.
4. Labelling information is located on the surface of the peel-back lid and on the outside of the rigid box.

Care must be taken to maintain the sterility of the implant by proper opening of the packaging and handling of the implant.

1. Open the implant package in the non-sterile field, with non-sterile gloves, tear the address label to open the box.
2. With non-sterile gloves, remove the inner blister pack. Do not place the plastic box or blister pack-lid onto the sterile field. The contents of this inner package are sterile.
3. The sealed blister is to be opened by an assistant (with nonsterile gloves): remove the TYVEK lid and drop or place the sterile tube onto the sterile field, open the tube cap and attach the implant placement tool onto the implant and carefully remove from the sterile tube. Do not touch the sterile implant.

Other sterile components are packed in a peel pouch or blister base with a "peel-back" lid. Labelling information is located on the bottom half of the pouch, inside the packet or on the surface of the peel-back lid. Sterility is assured unless the pouch is damaged or opened. Non-sterile components are supplied clean but not sterile in a peel pouch or blister base with peel-back lid. Labelling information is located on the bottom half of the pouch or on the surface of the peel-back lid.

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established.

The contact information for the manufacturer of this device to report a serious incident is as follows:

sicomplaints@southernimplants.com

Materials

Dental Implantation Drilling Template Guide Sleeves and
Abutment Retrieval Tools

Titanium alloy (Ti-6AL-4V, ASTM 136)

Dental Implant Internal Thread Shaping Tools and Dental
Implant Extractor

Stainless Steel 316 and Stainless Steel 431

Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

Summary of Safety and Clinical Performance (SSCP)

As required by the European Medical Device Regulation (MDR; EU2017/745), a Summary of Safety and Clinical Performance (SSCP) is available for perusal with regard to Southern Implants® product ranges.

The relevant SSCP can be accessed at <https://ec.europa.eu/tools/eudamed>.

NOTE: the above website will be available upon the launch of the European Database on Medical Devices (EUDAMED).

Disclaimer of liability

This product is part of the Southern Implants® product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants® product range and take full responsibility for the correct indications and use of this product. Southern Implants® does not assume liability for damage due to incorrect use. Please note that some Southern Implants® products may not be cleared or released for sale in all markets.

Basic UDI

Product	Basic-UDI Number
Basic-UDI for General Use Instruments	6009544039108M

Related literature and catalogues

CAT-2068 – SiGuided Product Catalogue
CAT-2089M – SIREAL Product Catalogue
CAT-2090M – Pilot Drill Guided Surgery Product Catalogue
CAT-2006 – Instrument Product Catalogue
CAT-1076 – Implant Removal Data Sheet
CAT-1224 – Abutment Retrieval Tool Data Sheet

Symbols and warnings

											
Manufacturer: Southern Implants® 1 Albert Rd, P.O Box 605 IRENE, 0062, South Africa. Tel: +27 12 667 1046	CE mark 2797	Prescription device*	Sterilised using irradiation	Non-sterile	Use by date (mm-yy)	Do not reuse	Do not resterilize	Catalogue number	Batch code	Medical device	Authorised representative in the European Community
											
Authorised representative for Switzerland	Date of manufacture	Magnetic Resonance conditional	Magnetic Resonance safe	Single sterile barrier system with protective packaging inside	Single sterile barrier system	Consult instruction for use	Caution	Keep away from sunlight	Do not use if package is damaged		

* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

All rights reserved. Southern Implants®, the Southern Implants® logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants®. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale. It is the responsibility of the clinician to inspect the symbols that appear on the packaging of the product in use.

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