English	INSTRUCTIONS FOR USE: Southern Implants® General Use Instruments
Español	INSTRUCCIONES DE USO: Instrumentos de uso general Southern Implants®
Italiano	ISTRUZIONI PER L'USO: Southern Implants® Strumenti per Uso Generale
Français	NOTICE D'UTILISATION : Instruments à usage général de Southern Implants®
Deutsch	GEBRAUCHSANWEISUNG: Southern Implants® Instrumente für den allgemeinen Gebrauch
Português	INSTRUÇÕES DE UTILIZAÇÃO: Southern Implants® Instrumentos de uso geral
Ελληνικά	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ: Εργαλεία γενικής χρήσης της Southern Implants®
Svenska	BRUKSANVISNING: Southern Implants® Instrument för allmän användning
Română	INSTRUCȚIUNI DE UTILIZARE: Southern Implants® Instrumente de utilizare generală
Türkçe	KULLANMA TALİMATI: Southern Implants® Genel Kullanıma Yönelik Aletler

Southern Implants (Pty.) Ltd. 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

Subsidiaries

Australia

Southern Implants Australia
T: +61-(0)-8-9466-2627
E: info@southernimplants.com.au

Spain and Portugal

Southern Implants Iberica
T: +34 935 053 507
E: info@southernimplants.es

United Kingdom and Ireland

Southern Implants UK
T: +44-20-805-94490
E: info@southernimplants.co.uk

USA and Canada

Southern Implants North America Inc.
T: +1-561-472-0990
E: customercare@southernimplants.com

Description

Dental implantation drilling template guide sleeves

A small tubular device intended to be inserted within the guide holes of a dental implantation drilling template (a device used to guide drilling in the jawbone during dental implant placement/fixation) to guide the drill and / or implant into position. These devices are supplied non-sterile and for single use only.

Dental implant internal thread shaping tools

A cylindrical rod with a spiralled cutting distal end, similar to a drill bit or a screw, intended to be used to shape and rectify the internal thread of a dental implant that has been implanted into bone. It is made of metal and is designed to be attached to a manual screwdriver. These devices are supplied non-sterile and for single use only.

Dental implant extractor

A manually-operated, dental device used to retrieve a dental implant from the oral cavity, typically due to damage (e.g., a broken collar) or malfunction of the dental implant. It is made of toughened, high-grade steel and has a cylindrical design with a coarse left-handed flute. The flute has a steep pitch that spirals up the tapered working end and there is a hexagonal head at the proximal end that fits a torque wrench socket. The device is inserted into the implant which is to be removed and turned anticlockwise (counter clockwise). The anticlockwise torque exerted makes this device grip into the implant, which is then unscrewed as the device continues to turn anti-clockwise. This device is supplied sterile and for single use only.

Abutment retrieval tools

A screwdriver driven dental device used to retrieve a dental abutment, typically one with a short abutment post (e.g., a passive abutment) and conical connection (e.g., Single Platform SP1 connection), from a dental implant after prolonged loading when it is found that the abutment is jammed into the implant's connection. It is made of Titanium alloy (ASTM F136) and has a screw design with an elongated nib section intended to bottom out in the implant thread pilot hole and lift the abutment out of the implant without damaging the internal thread or connection of the implant. The abutment retrieval tool has external thread that corresponds with internal thread either manufactured into the abutment's screw access hole or tapped by the user. The device is inserted into the abutment and rotated clockwise until the seal between the abutment and implant's conical connection is broken and the abutment is dislodged and may be removed from the implant. These devices are supplied non-sterile and for single use only.

Guide Pins

A device used to stabilize a drilling template guide during implant bed preparation or implant placement. They are used in conjunction with the drilling template guide sleeves. These devices are supplied non-sterile and for single use only.

Intended use

The Southern Implants® General Use Instruments are single-use devices intended to aid in the surgical procedures associated with implant therapy.

Specifically, the Southern Implants® Dental Implantation Drilling Template Guide Sleeve is intended for the visual control or physical guidance during implant bed preparation or implant placement.

The Southern Implants® Dental Implant Internal Thread Shaping Tool is intended to manually reshape the internal thread of an implant.

The Southern Implants® Dental Implant Extractor is intended to remove an endosseous implant through the application of manual reverse torque.

The Southern Implants® Abutment Retrieval Tool is intended to remove an abutment from an implant through the application of manual torque.

The Southern Implants® Guide Pins are intended for stabilizing a drilling template guide during implant bed preparation or implant placement.

Indications for use

The General Use Instruments are indicated for use in the treatment of partial or full edentulousness, congenital tooth loss, or failing/problematic natural dentition, when a patient is dissatisfied with existing removable dentures, or as part of the revision/replacement of previous dental restorations.

Intended user

Dental Technicians, Maxillofacial Surgeons, General Dentists, Orthodontists, Periodontists, Prosthodontists, and other appropriately trained and/or experienced medical professionals.

Intended environment

The devices are intended to be used in a dental laboratory as part of the restoration design and manufacture as well as in a clinical environment such as an operating theatre or a dentist consultation room.

Intended patient population

The Southern Implants® General Use Instruments are intended to be used in patients subject to implant treatment/revision.

Compatibility information

The Southern Implants' General Use Instruments and their identifiers are given in Table A.

Table A - Compatibility

General Use Instrument type	Product Code	Description	Sterility Status
Dental implantation drilling template guide sleeves	I-GS-2513-L6	Guided surgery pin sleeve Ø2.50 outer diameter, Ø1.30 inner diameter, length 6mm	Non-sterile
	I-GS-3020-L4	Guided surgery sleeve Ø3.00 outer diameter, Ø2.0 inner diameter, length 4mm	Non-sterile
	I-GS-3020-L6	Guided surgery sleeve Ø3.00 outer diameter, Ø2.0 inner diameter, length 6mm	Non-sterile
	I-GS-5142-L4	Guided surgery sleeve Ø5.10 outer diameter, Ø4.15 inner diameter, length 4mm	Non-sterile
	I-GS-5142-L9	Guided surgery sleeve Ø5.10 outer diameter, Ø4.15 inner diameter, length 9mm	Non-sterile
	I-GS-6252-L4	Guided surgery sleeve Ø6.20 outer diameter, Ø5.15 inner diameter, length 4mm	Non-sterile
	I-GS-6252-L9	Guided surgery sleeve Ø6.20 outer diameter, Ø5.15 inner diameter, length 9mm	Non-sterile
	I-GSC-5142-L4	Guided surgery C-type sleeve Ø5.10 outer diameter, Ø4.20 inner diameter, length 4mm	Non-sterile
	I-GSC-5142-L9	Guided surgery C-type sleeve Ø5.10 outer diameter, Ø4.20 inner diameter, length 9mm	Non-sterile
	I-GSC-6252-L4	Guided surgery C-type sleeve Ø6.20 outer diameter, Ø5.20 inner diameter, length 4mm	Non-sterile
	I-GSC-6252-L9	Guided surgery C-type sleeve Ø6.20 outer diameter, Ø5.20 inner diameter, length 9mm	Non-sterile
	I-GSZ-4128-L8	Guided surgery sleeve for zygomatic implants Ø4.10 outer diameter, Ø2.80 inner diameter, length 8mm	Non-sterile
	I-GS-4542-CS	Guided surgery sleeve for CEREC Ø4.45 outer diameter, Ø4.15 inner diameter, length 4mm	Non-sterile
	I-GS-5342-CM	Guided surgery sleeve for CEREC Ø5.25 outer diameter, Ø4.15 inner diameter, length 4mm	Non-sterile
	I-GS-6052-CL	Guided surgery sleeve for CEREC Ø5.95 outer diameter, Ø5.15 inner diameter, length 4mm	Non-sterile
	I-GSF-5142-L4	Guided surgery sleeve with anti-rotate flats, Ø5.10 outer diameter, Ø4.15 inner diameter, length 4mm	Non-sterile
	I-GSF-5142-L6	Guided surgery sleeve with anti-rotate flats, Ø5.10 outer diameter, Ø4.15 inner diameter, length 6mm	Non-sterile
	I-GSF-6252-L4	Guided surgery sleeve with anti-rotate flats, Ø6.20 outer diameter, Ø5.15 inner diameter, length 4mm	Non-sterile
	I-GSF-6252-L6	Guided surgery sleeve with anti-rotate flats, Ø6.20 outer diameter, Ø5.15 inner diameter, length 6mm	Non-sterile

	I-GSF-6252-L4-P5	Guided surgery sleeve with anti-rotate flats, Ø6.20 outer diameter, Ø5.15 inner diameter, length 4mm, 5 pack	Non-sterile
	I-GSCF-6252-L6	Guided surgery C-type sleeve with anti-rotate flats, Ø6.20 outer diameter, Ø5.15 inner diameter, length 6mm	Non-sterile
Dental implant internal thread shaping tools	I-M1.4-S	M1.4 thread tap, length short	Non-sterile
	I-M1.4	M1.4 thread tap, length long	Non-sterile
	I-M1.6-S	M1.6 thread tap, length short	Non-sterile
	I-M1.6	M1.6 thread tap, length long	Non-sterile
	I-M1.8-S	M1.8 thread tap, length short	Non-sterile
	I-M1.8	M1.8 thread tap, length long	Non-sterile
	I-M2-S	M2 thread tap, length short	Non-sterile
	I-M2	M2 thread tap, length long	Non-sterile
Dental implant extractor	I-IMP-REM	Implant extractor, wrench insert	Sterile
Abutment retrieval tools	I-ART-SP	Abutment retrieval tool/ screw to be used with the Single Platform SP1 passive abutments PA-SP & PA-NSP	Non-sterile
Guide Pins	I-D12T-GP	Guide Pin for Ø1.2 drill	Non-sterile

Clinical performance

The General Use Instruments are intended to aid in the surgical procedures associated with implant therapy. Consequently, the clinical performance of the General Use Instruments is primarily defined by the success of the implant therapy. This performance can be quantitatively assessed through implant survival/success, which serves as an important indicator of the instruments' effectiveness and success in clinical practice.

Clinical benefits

The Southern Implants® General Use Instruments play a crucial role in ensuring the successful outcome of surgical procedures, though they do not provide a direct clinical benefit. However, the clinical benefits of the associated treatment, target implant, and/or prosthetic device(s) should be acknowledged, with relevant documentation reviewed as necessary.

Storage, cleaning and sterilisation

Applicable for Sterile Components:

The component is supplied sterile (sterilised by gamma irradiation) and is intended for single use prior to the expiration date (see packaging label). Sterility is assured unless the container or seal is damaged or opened. If packaging is damaged do not use the product and contact your Southern representative or return to Southern Implants®. The devices must be stored in a dry place at room temperature and not exposed to direct sunlight. Incorrect storage may influence device characteristics.

Do not reuse components indicated for single-use only. Reusing these components may result in:

- damage to the surface or critical dimensions, which may result in performance and compatibility degradation.
- adds the risk of cross-patient infection and contamination if single-use items are reused.

Southern Implants® does not accept any responsibility for complications associated with reused single-use components.

Applicable for Non-Sterile Components:

This component is supplied non-sterile and is indicated for single use. If packaging is damaged do not use the product and contact your Southern representative or return to Southern Implants®. The devices must be stored in a dry place at room temperature and not exposed to direct sunlight. Incorrect storage may influence device characteristics.

Do not reuse components indicated for single-use only. Re-using these components may:

- damage on the surface or critical dimensions, which may result in performance and compatibility degradation.
- adds the risk of cross-patient infection and contamination if single-use items are reused.

Southern Implants® does not accept any responsibility for complications associated with reused single-use components.

Southern Implants® recommends one of the following procedures to sterilise the restorations and non-sterile single-use components prior to use:

1. pre-vacuum sterilisation method: steam sterilise the device at 132°C (270°F) at 180 - 220 kPa for 4 minutes. Dry for at least 20 minutes in the chamber. Only an approved wrap or pouch for steam sterilisation must be used.
2. for users in the USA: pre-vacuum sterilisation method: wrapped, steam sterilise at 135°C (275°F) at 180 - 220 kPa for 3 minutes. Dry for 20 minutes in the chamber. Use a wrap or pouch that is cleared for the indicated steam sterilisation cycle.

NOTE: users in the USA must ensure that the steriliser, wrap or pouch, and all steriliser accessories are cleared by the FDA, for the intended sterilisation cycle.

Contraindications

The contraindications of all device groups used as part of the specific treatment or procedure apply. Therefore, the contraindications of the systems/medical devices utilized as part of implant surgery/therapy should be noted and the relevant documents consulted.

The contraindications specific to this device group include:

- Allergies or hypersensitivity to chemical ingredients in the following materials: Titanium, Aluminium, Vanadium, Stainless Steel.

Warnings and precautions

IMPORTANT NOTICE: THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING

- To ensure the safe and effective use of dental implants, new technologies/systems and the General Use Instrument devices, it is strongly recommended that specialised training be undertaken. This training should include hands-on methods to gain competency on the proper technique, the system's biomechanical requirements and radiographic evaluations required for the specific system.
- Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone.
- Use of the device with incompatible or non-corresponding devices can result in the poor performance or failure of the device.
- When handling devices intraorally, it is imperative that they are adequately secured to prevent aspiration, as aspiration of products may lead to infection or physical injury.
- The use of non-sterile items can lead to secondary infections of the tissue or transfer infectious diseases.
- Failure to adhere to appropriate cleaning, re-sterilization, and storage procedures as outlined in the Instructions for Use can result in device damage, secondary infections, or patient harm.
- Exceeding the number of recommended uses for reusable devices can result in device damage, secondary infection or patient harm.
- The use of blunt drills may cause damage to the bone, potentially compromising osseointegration.

It is crucial to emphasize that training should be undertaken by both new and experienced implant users prior to utilizing a new system or attempting a new treatment method.

Patient Selection and Preoperative Planning

A comprehensive patient selection process and meticulous preoperative planning are essential for successful implant treatment. This process should involve consultation among a multidisciplinary team, including well-trained surgeons, restorative dentists, and laboratory technicians.

Patient screening should include, at a minimum, a thorough medical and dental history, as well as visual and radiological inspections to assess the presence of adequate bone dimensions, the positioning of anatomical landmarks, the presence of unfavourable occlusal conditions, and the periodontal health status of the patient.

For successful implant treatment, it is important to:

1. Minimize trauma to the host tissue, thereby increasing the potential for successful osseointegration.
2. Accurately identify measurements relative to radiographic data, as failure to do so may lead to complications.
3. Be vigilant in avoiding damage to vital anatomical structures, such as nerves, veins, and arteries. Injury to these structures may result in serious complications, including ocular injury, nerve damage, and excessive bleeding.

The responsibility for proper patient selection, adequate training and experience in implant placement, and the provision of appropriate information required for informed consent rests with the practitioner. By combining thorough screening of prospective implant candidates with a practitioner possessing a high level of competence in the use of the system, the potential for complications and severe side effects can be significantly reduced.

High Risk Patients

Special care should be exercised when treating patients with local or systemic risk factors that may adversely affect the healing of bone and soft tissue or otherwise increase the severity of side effects, the risk of complications, and/or the likelihood of implant failure. Such factors include:

- poor oral hygiene
- history of smoking/vaping/tobacco use
- history of periodontal disease
- history of orofacial radiotherapy**
- bruxism and unfavourable jaw relations
- use of chronic medications that may delay healing or increase the risk of complications including, but not limited to, chronic steroid therapy, anti-coagulant therapy, TNF- α blockers, bisphosphonate and cyclosporin

*** The potential for implant failure and other complications increases when implants are placed in irradiated bone, as radiotherapy can lead to progressive fibrosis of blood vessels and soft tissue (i.e., osteoradionecrosis), resulting in diminished healing capacity. Contributing factors to this increased risk include the timing of implant placement in relation to radiation therapy, the proximity of radiation exposure to the implant site, and the radiation dosage at that site.*

Should the device not operate as intended, it must be reported to the manufacturer of the device. The contact information or the manufacturer of this device to report a change in performance is: sicomplaints@southernimplants.com.

Side effects

The clinical outcome of treatment is influenced by multiple factors. The following side effects and residual risks are associated with the General Use Instruments device group or the surgical procedures involved in its use. They may necessitate further treatment, revision surgery or additional visits to the relevant medical professional. Furthermore, these side effects and residual risks can vary in both severity and frequency.

Therefore, the full list of known potential undesirable side-effects and residual risks that have been identified in relation to the General Use Instruments device group include:

- Allergic or hypersensitivity reaction(s)
- Bruising
- Dental injury during surgery
- Gingival injury or recession
- Infection (acute and/or chronic)
- Localized inflammation
- Pain or discomfort
- Soft tissue irritation
- Surgical side effects such as pain, inflammation, bruising and mild bleeding
- Wound dehiscence or poor healing
- Sub-optimal aesthetic result

Precaution: maintaining sterility protocol

Applicable to Sterile Components:

The devices supplied sterile are packed in a peel pouch or blister base with a “peel-back” lid. Labelling information is located on the bottom half of the pouch, inside the packet or on the surface of the peel-back lid. Sterility is assured unless the pouch or blister is damaged or opened.

Applicable to Non-Sterile Components:

The devices supplied non-sterile are supplied clean but not sterile in a peel pouch or blister base with peel-back lid. Labelling information is located on the bottom half of the pouch or on the surface of the peel-back lid.

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established.

The contact information for the manufacturer of this device to report a serious incident is as follows:

sicomplaints@southernimplants.com.

Materials

Dental Implantation Drilling Template Guide Sleeves:	Titanium alloy Ti-6AL-4V (90% titanium, 6% aluminium, 4% vanadium) according to ASTM 136 and ISO 5832-3
Abutment Retrieval Tools:	Titanium alloy Ti-6AL-4V (90% titanium, 6% aluminium, 4% vanadium) according to ASTM 136 and ISO 5832-3
Dental Implant Internal Thread Shaping Tools:	Stainless Steel 316
Dental Implant Extractor:	Stainless Steel 431
Guide Pins:	Stainless Steel 303 or 316

Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

Summary of Safety and Clinical Performance (SSCP)

As required by the European Medical Device Regulation (MDR; EU2017/745), a Summary of Safety and Clinical Performance (SSCP) is available for perusal with regard to Southern Implants® product ranges.

The relevant SSCP can be accessed at <https://ec.europa.eu/tools/eudamed>.

NOTE: the above website will be available upon the launch of the European Database on Medical Devices (EUDAMED). In the interim, email clinical@southernimplants.com to request a Summary of Safety and Clinical Performance (SSCP) for these devices.

Disclaimer of liability

This product is part of the Southern Implants® product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants® product range and take full responsibility for the correct indications and use of this product. Southern Implants® does not assume liability for damage due to incorrect use. Please note that some Southern Implants® products may not be cleared or released for sale in all markets.

Basic UDI

Product	Basic-UDI Number
Basic-UDI for General Use Instruments	6009544039108M

Related literature and catalogues

CAT-2068 – SIGuided Product Catalogue
 CAT-2089M – SIREAL Product Catalogue
 CAT-2090M – Pilot Drill Guided Surgery Product Catalogue
 CAT-1076 – Implant Removal Data Sheet
 CAT-1224 – Abutment Retrieval Tool Data Sheet

Symbols and warnings

Manufacturer

CE mark

Prescription device*

Sterilised using irradiation

Non-sterile

Use by date (mm-yy)

Do not reuse

Do not resterilise

Catalogue number

Batch code

Medical device

Authorised representative in the European Community

Authorised representative for Switzerland

Date of manufacture

Magnetic Resonance conditional

Magnetic Resonance safe

Single sterile barrier system with protective packaging inside

Single sterile barrier system

Consult instruction for use

Caution

Keep away from sunlight

Do not use if package is damaged

* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

All rights reserved. Southern Implants®, the Southern Implants® logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants®. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale. It is the responsibility of the clinician to inspect the symbols that appear on the packaging of the product in use.

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