English	INSTRUCTIONS FOR USE: Southern Implants® Titanium Blanks
Español	INSTRUCCIONES DE USO: Bloques de Titanio Southern Implants®
Italiano	ISTRUZIONI PER L'USO: Monconi vergini in titanio Southern Implants®
Français	MODE D'EMPLOI : Blocs titane de Southern Implants®
Deutsch	GEBRAUCHSANWEISUNG: Southern Implants®-Titanrohlinge
Português	INSTRUÇÕES DE UTILIZAÇÃO: Blocos de titânio da Southern Implants®
Ελληνικά	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ: Κενά Τιτανίου της Southern Implants®
Svenska	BRUKSANVISNING: Southern Implants® Titanblank
Română	INSTRUCȚIUNI DE UTILIZARE: Blocuri din titan Southern Implants®
Türkçe	KULLANMA TALİMATI: Southern Implants® Titanyum Bloklar

Southern Implants (Pty.) Ltd. 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

Subsidiaries

Australia

Southern Implants Australia
T: +61-(0)-8-9466-2627
E: info@southernimplants.com.au

Spain and Portugal

Southern Implants Iberica
T: +34 935 053 507
E: info@southernimplants.es

United Kingdom and Ireland

Southern Implants UK
T: +44-20-805-94490
E: info@southernimplants.co.uk

USA and Canada

Southern Implants North America Inc.
T: +1-561-472-0990
E: customercare@southernimplants.com

Description

Titanium Blanks are premanufactured prosthetic components intended to be connected to endosseous dental implants and act as a supporting base for permanent prosthetic restorations. The Titanium Blanks are indicated in Southern Implants' digital workflow scan files from desktop/intraoral scanners, CAD software, CAM software, ceramic material, milling machine and associated tooling and accessories. These abutments are supplied non-sterile and are intended for single-use only.

Intended use

The Southern Implants® Metal Abutments are prosthetic components intended for connection to endosseous implant systems to aid in prosthetic rehabilitation and implant procedures.

The Southern Implants® Direct Abutments are intended to provide a supportive base for connected prosthetic devices on endosseous implants.

Indications for use

The Southern Implants® Metal Abutments are indicated for use as components of an endosseous implant system to support prosthetic rehabilitation in implant treatment procedures.

The Southern Implants® Direct Abutments are indicated for connection to Southern Implants® fixtures as part of the functional and aesthetic restoration of the upper or lower jaw using a fixed, *permanent* prosthesis in cases where the endosseous implant system is eligible for functional loading.

Intended user

Dental technicians, Maxillofacial Surgeons, General Dentists, Orthodontists, Periodontists, Prosthodontists, and other appropriately trained and/or experienced medical professionals.

Intended environment

The devices are intended to be used in a dental laboratory as part of the restoration design and manufacture as well as in a clinical environment such as an operating theatre or a dentist consultation room.

Intended patient population

Southern Implants® Metal Abutments are intended to be used in patients indicated and eligible, or otherwise not contraindicated, for implant-retained prosthetic restorations.

Compatibility information

Southern Implants' implants should be restored with Southern Implants' components. In the Southern Implants' Titanium Blank range there are 7 implant connection types. The compatible implant and connection type can be identified by specific abbreviations in the product codes. Range identifiers are summarised in Table A.

For further information see the SI Digital Product Catalogue (CAT-2063).

Table A - Compatibility

Implant connection type	Compatible devices
External Hex (EX)	Parts labelled BL-EX(Ø)-11 for engaging items
TRI-NEX (EL) (Lobe)	Parts labelled BL-EL(Ø)-11 for engaging items
Deep Conical (DC)	Parts labelled BL-DC(Ø)-11 for engaging items
Internal Hex (M)	Parts labelled BL-M-11 for engaging items (used with Ø3.75, 4.20 and 5.00mm platforms)
	Parts labelled BL-M-P45-11 for engaging items (used with Ø5.0mm platform)
Internal Hex PROVATA® (3M/M/Z)	Parts labelled BL-3M-11 for engaging items (used with Ø3.3 mm platform)
	Parts labelled BL-M-11 for engaging items (used with Ø4.0, 5.0 and 6.0mm platforms)
	Parts labelled BL-M-P45-11 for engaging items (used with Ø5.0 and 6.0mm platforms)
	Parts labelled BL-Z-11 for engaging items (used with Ø6.0, 7.0, 8.0 and 9.0mm platforms)

Internal Octagon IT (ITS/ITS6)	Parts labelled BL-ITS-11 for engaging items (used with Ø4.8mm platforms)
	Parts labelled BL-IT6-11 for engaging items (used with Ø6.5mm platforms)
Single Platform (SP1)	Parts labelled BL-SP-11 for engaging items (used with all SP1 platforms)
	Parts labelled BL-SP-PM-11 for engaging items (used with Ø4.5, Ø5.0mm and Ø6.0 platforms)

Additionally, the Titanium Blanks are suitable for use with the Ø11.5mm Medentika® blank holders. Refer to the Medentika® website for information on compatible blank holders and other related equipment.

Clinical performance

The clinical performance of Southern Implants® Titanium Blanks is primarily defined by their ability to maintain stability and retention of the prosthetic restoration over time. This stability is achieved through the secure fixation of the restoration to the abutment, which can be accomplished via cementation or casting, as well as the mechanical fastening of the abutment to the implant through screws.

Clinical benefits

The clinical benefits associated with the Southern Implants® Titanium Blanks are intrinsically linked to the overall success of implant-supported prosthetic treatments and encompass both direct and indirect clinical benefits.

The direct clinical benefits of the Southern Implants® Titanium Blanks arise specifically from their use in the treatment process. These benefits include the enhancement of peri-implant soft tissue health and aesthetic outcomes of the final restoration. The clinical benefits experienced by patients as a result of implant-supported prosthetic treatments in which these abutments are utilized represent the indirect clinical benefits of the Southern Implants® Titanium Blanks. While these benefits are substantial, they are not directly attributable to the devices themselves; rather, they reflect the contributions of the abutments to the overall treatment system. Additionally, successful treatment outcomes are anticipated to result in significant enhancements in patients' psychosocial well-being, self-esteem, and overall quality of life. Overall, the Southern Implants® Titanium Blanks not only facilitate the mechanical aspects of dental restorations but also contribute to the holistic improvement of patient outcomes and experiences.

Storage, cleaning and sterilisation

This component is supplied non-sterile and is indicated for single-use. If packaging is damaged do not use the product and contact your Southern representative or return to Southern Implants®. The devices must be stored in a dry place at room temperature and not exposed to direct sunlight. Incorrect storage may influence device characteristics.

Do not reuse components indicated for single-use only. Re-using these components may:

- damage on the surface or critical dimensions, which may result in performance and compatibility degradation.
- adds the risk of cross-patient infection and contamination if single-use items are reused.

Southern Implants® does not accept any responsibility for complications associated with reused single-use components.

Southern Implants® recommends one of the following procedures to sterilise the restorations and non-sterile single-use components prior to use:

1. pre-vacuum sterilisation method: steam sterilise the abutments at 132°C (270°F) at 180 - 220 kPa for 4 minutes. Dry for at least 20 minutes in the chamber. Only an approved wrap or pouch for steam sterilisation must be used.
2. for users in the USA: pre-vacuum sterilisation method: wrapped, steam sterilise at 135°C (275°F) at 180 - 220 kPa for 3 minutes. Dry for 20 minutes in the chamber. Use a wrap or pouch that is cleared for the indicated steam sterilisation cycle.

NOTE: users in the USA must ensure that the steriliser, wrap or pouch, and all steriliser accessories are cleared by the FDA, for the intended sterilisation cycle.

Contraindications

The contraindications of all device groups used as part of the specific treatment or procedure apply. Therefore, the contraindications of the systems/medical devices utilized as part of the implant surgery/therapy should be noted and the relevant documents consulted.

The contraindications specific to this device group include:

- patients who are medically unfit for implant or prosthetic related procedures
- patients that present with an allergy or hypersensitivity to the following materials: Titanium, Aluminium, Vanadium, Cobalt Chrome, Cobalt Chrome Molybdenum, Gold, Platinum, Palladium, Iridium, Silver and/or Copper

Warnings and precautions

IMPORTANT NOTICE: THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING

- To ensure the safe and effective use of dental implants, new technologies/systems and the Metal Abutment devices, it is strongly recommended that specialised training be undertaken. This training should include hands-on methods to gain competency on the proper technique, the system's biomechanical requirements and radiographic evaluations required for the specific system.
- Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone.
- Use of the device with incompatible or non-corresponding devices can result in the poor performance or failure of the device.
- When handling devices intraorally, it is imperative that they are adequately secured to prevent aspiration, as aspiration of products may lead to infection or physical injury.
- The use of non-sterile items can lead to secondary infections of the tissue or transfer infectious diseases.
- Failure to adhere to appropriate cleaning, re-sterilization, and storage procedures as outlined in the Instructions for Use can result in device damage, secondary infections, or patient harm.
- Exceeding the number of recommended uses for reusable devices can result in device damage, secondary infection or patient harm.
- The use of blunt drills may cause damage to the bone, potentially compromising osseointegration.

It is crucial to emphasize that training should be undertaken by both new and experienced implant users prior to utilizing a new system or attempting a new treatment method.

Patient Selection and Preoperative Planning

A comprehensive patient selection process and meticulous preoperative planning are essential for successful implant treatment. This process should involve consultation among a multidisciplinary team, including well-trained surgeons, restorative dentists, and laboratory technicians.

Patient screening should include, at a minimum, a thorough medical and dental history, as well as visual and radiological inspections to assess the presence of adequate bone dimensions, the positioning of anatomical landmarks, the presence of unfavourable occlusal conditions, and the periodontal health status of the patient.

For successful implant treatment, it is important to:

1. Minimize trauma to the host tissue, thereby increasing the potential for successful osseointegration.
2. Accurately identify measurements relative to radiographic data, as failure to do so may lead to complications.
3. Be vigilant in avoiding damage to vital anatomical structures, such as nerves, veins, and arteries. Injury to these structures may result in serious complications, including ocular injury, nerve damage, and excessive bleeding.

The responsibility for proper patient selection, adequate training and experience in implant placement, and the provision of appropriate information required for informed consent rests with the practitioner. By combining thorough screening of prospective implant candidates with a practitioner possessing a high level of competence in the use of the system, the potential for complications and severe side effects can be significantly reduced.

High Risk Patients

Special care should be exercised when treating patients with local or systemic risk factors that may adversely affect the healing of bone and soft tissue or otherwise increase the severity of side effects, the risk of complications, and/or the likelihood of implant failure. Such factors include:

- poor oral hygiene
- history of smoking/vaping/tobacco use
- history of periodontal disease
- history of orofacial radiotherapy**
- bruxism and unfavourable jaw relations
- use of chronic medications that may delay healing or increase the risk of complications including, but not limited to, chronic steroid therapy, anti-coagulant therapy, TNF- α blockers, bisphosphonate and cyclosporin

*** The potential for implant failure and other complications increases when implants are placed in irradiated bone, as radiotherapy can lead to progressive fibrosis of blood vessels and soft tissue (i.e., osteoradionecrosis), resulting in diminished healing capacity. Contributing factors to this increased risk include the timing of implant placement in relation to radiation therapy, the proximity of radiation exposure to the implant site, and the radiation dosage at that site.*

Should the device not operate as intended, it must be reported to the manufacturer of the device. The contact information or the manufacturer of this device to report a change in performance is: sicomplaints@southernimplants.com.

Side effects

The clinical outcome of treatment is influenced by various factors and there are additional potential side effects and residual risks associated with the device group. These may necessitate further treatment, revision surgery or additional visits to the relevant medical professional's office. Furthermore, these side effects and residual risks can vary in both severity and frequency.

Therefore, the full list of known potential undesirable side-effects and residual risks that have been identified in relation to the Metal Abutments include:

- Abutment fracture
- Abutment screw fracture or distortion
- Allergic reaction(s) to the abutment material
- Anaesthesia, paraesthesia, hyperesthesia, and hypoesthesia (transient or permanent)
- Bleeding
- Bleeding on probing
- Bruising
- Complications requiring revision surgery
- Difficulty in retrieving the abutment
- Gingival inflammation
- Gingival recession
- Hyperplastic tissue response
- Implant failure due to insufficient levels of osseointegration
- Improper implant positioning resulting in prosthetic compromise
- Infection (acute and/or chronic)
- Localized inflammation
- Loosening of the abutment screw and/or retaining screw
- Loss or damage to adjacent teeth
- Marginal bone loss
- Micromovements and implant instability
- Misfit or improper connection at the implant-abutment interface
- Overloading of the abutment/implant
- Pain, tenderness or discomfort

- Peri-implantitis, peri-mucositis or otherwise poor peri-implant soft tissue health
- Periodontal inflammation
- Phonetic difficulties
- Prosthetic failure
- Soft tissue irritation
- Sub-optimal aesthetic result
- Wound dehiscence or poor healing

Precaution: maintaining sterility protocol

The Titanium Blanks are supplied clean but not sterile in a peel pouch or blister base with peel-back lid. Labelling information is located on the bottom half of the pouch or on the surface of the peel-back lid.

Restoration design restrictions

The Titanium Blanks can be altered to the desired anatomical shape provided the pre-milled implant connection is not altered. Always ensure that any sharp edges are rounded.

For the Titanium Blanks CAD software designs, the following restoration design restrictions are applicable per connection interface type (See Figure A and Table B):

Figure A – Illustration of Design Restrictions


Table B – Design Restrictions

	Connection Interface	Indication	Maximum Abutment Angulation (A)	Chimney Height (B)		Collar Height (C)		Platform Diameter (D)		Minimum Chimney Wall Thickness (E) (i.e. screw hole to outer surface)
				Max.	Min.	Max.	Min.	Max.	Min.	
External Hex	IP	1 mandible, 2	20°	8	4,5	5	1,5	4,5	3,85	0,4
	IBN	All*	20°	8 for anterior use 4.5 for all regions			0,6	4,8	4,3	
	IB	All		8				5,5	4,8	
	BA	All						6,5	5,5	
	BBB	All						7,5	6,5	

Deep Conical	DC3	1 mandible, 2	20°
	DC4	All	20°
	DCR50	All	
	DC5	All	
Internal Hex (PROVATA/ M-Series)	3M	1, 2	20°
	M	All	
	M-P45	All	
	Z	All	
Tri-Nex	EL35	All	20°
	EL43	All	
	EL50	All	
	EL60	All	
Internal Octagon (IT)	ITS	All	20°
	IT6	All	
Single Platform (SP1)	SP	All	20°
	SP-PM	All	

	5	1,5	4,5	3,85	
	5	0,7	4,8	3,85	
	3	0,85	5,5	5,5	
	5	0,75	6,5	5	
	5	0,6	1,5	4,5	3,85
			4,8	4,3	
			5,5	4,8	
			7,5	6	
	5	0,6	4,8	4,3	
	5		5,5	4,5	
	5		6,5	5,5	
	3		7,5	6,5	
	3	0,25	6	4,84	
	1,5		7,5	6,54	
	5	1,5	4,8	4,3	
			6,5	5,5	

*dependent on the selected maximum chimney height

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established.

The contact information for the manufacturer of this device to report a serious incident is as follows:

sicomplaints@southernimplants.com.

Materials

Titanium Blanks: Titanium alloy (Ti-6Al-4V) according to ASTM F136 and ISO 5832-3

Chemical Components	Aluminium (Al)	Vanadium (V)	Residuals* (Fe, O, C, N, H) in total	Titanium (Ti)
Composition, % (mass/mass)	5.50 – 6.75	3.50 – 4.50	<0.55	Balance

*Where, Fe = Iron, O = Oxygen, C = Carbon, N = Nitrogen, H = Hydrogen

Prosthetic Screws: Titanium alloy (Ti-6Al-4V) according to ASTM F136 and ISO 5832-3

Chemical Components	Aluminium (Al)	Vanadium (V)	Residuals* (Fe, O, C, N, H) in total	Titanium (Ti)
Composition, % (mass/mass)	5.50 – 6.75	3.50 – 4.50	<0.55	Balance

*Where, Fe = Iron, O = Oxygen, C = Carbon, N = Nitrogen, H = Hydrogen

Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

Magnetic Resonance (MR) Safety

Nonclinical testing has demonstrated that the Southern Implants® dental implants, metallic abutments and prosthetic screws are MR conditional.

MRI Safety Information	
 MR Conditional	
A patient with the dental implants, metallic abutments and prosthetic screws from Southern Implants® may be safely scanned under the following conditions. Failure to follow these conditions may result in injury to the patient.	
Device Name	dental implants, metallic abutments and prosthetic screws
Static Magnetic Field Strength (B₀) [T]	1.5T or 3T
Maximum Spatial Field Gradient [T/m] [gauss/cm]	30 T/m (3000 gauss/cm)
Type of RF Excitation	Circularly polarized (CP) (i.e., quadrature driven)
RF Transmit Coil Type	Integrated Whole Body Transmit Coil
RF Receive Coil	Any receive-only RF coil may be used
Operating Mode	Normal Operating Mode
RF Conditions	Normal Operating mode: (Head SAR of 3.2 W/kg, 2 W/kg whole body)
Scan Duration	If the closest part of implant <25cm away from the isocenter in feet-head direction: a continuous scan time of 15 minutes will require a cooling delay of at least 5 minutes. If the closest part of implant ≥ 25cm away from the isocenter in feet-head direction: This is limited to 1 hour of scan duration.
MR Image Artifact	The presence of this implant may produce an image artifact.

Additional information:

1. Local transmit coil: Local transmit coil safety was not evaluated. Local transmit coil is only allowed if the implant is ≥ 25cm away from isocenter and for a wbSAR of up to 2W/kg.
2. Artifacts: MR Image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the implant. Tests following guidelines in ASTM F2119-24 showed artifact widths of ≤20 mm from the surface of the implant for the gradient echo sequence and ≤10 mm from the surface of the implant for the spin echo sequence.

Removable restorations should be taken out prior to scanning, as is done for watches, jewellery etc.

Should there be no MR symbol on the product label, please note that this device has not been evaluated for safety and compatibility in the MR environment.

Summary of Safety and Clinical Performance (SSCP)

As required by the European Medical Device Regulation (MDR; EU2017/745), a Summary of Safety and Clinical Performance (SSCP) is available for perusal with regard to Southern Implants® product ranges.

The relevant SSCP can be accessed at <https://ec.europa.eu/tools/eudamed>.

NOTE: the above website will be available upon the launch of the European Database on Medical Devices (EUDAMED). In the interim, email clinical@southernimplants.com to request a Summary of Safety and Clinical Performance (SSCP) for these devices.

Disclaimer of liability

This product is part of the Southern Implants® product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants® product range and take full responsibility for the correct indications and use of this product. Southern Implants® does not assume liability for damage due to incorrect use. Please note that some Southern Implants® products may not be cleared or released for sale in all markets.


Basic UDI

Product	Basic-UDI Number
Basic-UDI for Metal Abutments	60095440387296
Basic-UDI for Direct Abut. Titanium	6009544050117Q

Related literature and catalogues

CAT-2004 – Tri-Nex® Implants Product Catalogue
 CAT-2005 – IT Implants Product Catalogue
 CAT-2020 – External Hex Implants Product Catalogue
 CAT-2042 – Deep Conical Implants Product Catalogue
 CAT-2043 – Internal Hex Implants Product Catalogue
 CAT-2060 – PROVATA® Implants Product Catalogue
 CAT-2069 – Deep Conical INVERTA® Implants Product Catalogue
 CAT-2095 – External Hex INVERTA® Implants Product Catalogue
 CAT-2093 – Single Platform (SP1) Implants Product Catalogue
 CAT-2063 – SI Digital Product Catalogue

Symbols and warnings

* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

All rights reserved. Southern Implants®, the Southern Implants® logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants®. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale. It is the responsibility of the clinician to inspect the symbols that appear on the packaging of the product in use.

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