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South Africa - Headquarters: 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

Subsidiaries

United Kingdom and Ireland
Southern Implants UK
T: +44-20 8059 4490
E: info@southernimplants.co.uk

USA and Canada
Southern Implants North America Inc.
T: +1-561-472-0990
E: customercare@southernimplants.com

EC REP

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

CH REP

MedEnvoy Switzerland: Gotthardstrasse 28, 6302 Zug, Switzerland

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Description

Angle Gauge

A reusable dental implant space/depth/angle gauge is used during the planning and post-osteotomy stages of a dental implantation procedure, prior to surgical implant placement/fixation, to measure various lengths (e.g. drill depth of the osteotomy cavity, mucosal height, distance between implants/teeth for determining the mesiodistal position of the implant) and/or angles (e.g. direction of the perforation in relation to the antagonist arch), to aid in selection of an appropriate implant and/or final prosthetic abutment to be installed. It is typically graduated and may be a hand-held instrument or a stand-alone device placed intraorally; some types are intended for use during x-ray imaging.

Depth Gauge

A reusable dental implant space/depth/angle gauge is used during the planning and post-osteotomy stages of a dental implantation procedure, prior to surgical implant placement/fixation, to measure various lengths (e.g. drill depth of the osteotomy cavity, mucosal height, distance between implants/teeth for determining the mesiodistal position of the implant) and/or angles (e.g. direction of the perforation in relation to the antagonist arch), to aid in selection of an appropriate implant and/or final prosthetic abutment to be installed. It is typically graduated and may be a hand-held instrument or a stand-alone device placed intraorally; some types are intended for use during x-ray imaging.

Surgical Screwdrivers

The screwdrivers are intended to be used to apply torque to a prosthetic screw or abutment to affix the screw or abutment to the dental implant, or loosen and remove the component from the dental implant or abutment. The screwdrivers have a driving interface on one end that is compatible with the corresponding screw or abutment. This interface allows for pickup of the corresponding component and prevents relative rotation between the screw or abutment and the screwdriver when torque is applied. The screwdrivers are available in various lengths (short 5 mm, medium 14 mm, long 21 mm) to allow for use in both the anterior and posterior where available space may differ.

Dental Implant Holding Endpiece

A device intended to hold a dental implant (e.g., orthodontic implant screw) at its distal end during the implantation of the implant to the patient's oral cavity. It is typically a small metal shaft designed to be mounted onto the distal end of a manual screwdriver whereby rotation is applied to introduce the implant into a surgical cavity.

Basic Dental Retractors

A hand-held, non-self-retaining dental instrument intended to be used to manually displace soft tissues of the oral cavity to improve their visualization and access, and to afford them protection during oral surgical procedures. It is a metallic device and includes a shaped/hooked blade with a handle/holding portion.

Screw Extractor Guides

A hand-held manual surgical tool designed to remove implanted screws that have a damaged screw recess by self-tapping into the defective recess/head and providing sufficient grip to rotate and remove the screw in the normal way. It typically has a conical tip that is placed into the implant connection to provide a channel through which the screw remover is guided during use. The device is also intended to protect the implant internal connection during the screw removal procedure.

Dental Torque Wrenches

A hand-held dental tool intended for the manual application of rotation to tighten/loosen/remove a device (e.g. dental implant, spatial register (localizer/tracker/calibrator)) fitted into or close to the oral cavity, and/or to a surgical instrument, in association with a dental surgical procedure (e.g. dental guided surgery registration). It may have a built-in ratchet function, with forward and reverse direction, and a calibrated scale to which a recommended torque can be set; it is not a screwdriver. It is typically made of metallic materials and may have a standard socket connection to accommodate different tool bits/adaptors.

Wrench Adaptors

A component of a surgical wrench (typically ratchet type) intended to enable connection of the main body of the wrench to a variety of distal end pieces (bits). It is intended to be used during an invasive surgical procedure, and is not dedicated to dental use only. It may be referred to as a converter.

Retention Components

A device intended to be assembled with an implant insertion tool, dental implant holding endpiece or drill extension, to retain a dental implant, implant fixture mount or drill during implant placement procedures.

Intended use

The Southern Implants® Reusable Instruments are reusable devices intended for use in the planning, osteotomy preparation and placement of implants/prosthetics where they aid in the surgical procedures associated with implant therapy.

Specifically, the Southern Implants® Angle Gauge is intended to be used as a visual guide to confirm or determine the direction or orientation of the prepared osteotomy.

The Southern Implants® Depth Gauge is intended to be used as a visual guide to confirm or determine the depth of the prepared osteotomy.

The Southern Implants® Surgical Screw Drivers are intended to apply or transmit torque to abutments and/or screws.

The Southern Implants® Dental Implant Holding Endpiece is intended to apply or transmit torque to an implant during implantation.

The Southern Implants® Basic Dental Retractors are intended to displace the soft tissue of the oral cavity and improve visualization/access during implant bed preparation or implant placement.

The Southern Implants® Screw Extractor Guides are intended to be used as a physical guide for the Screw Remover during screw removal procedures.

The Southern Implants® Dental Torque Wrench is intended to apply or transmit torque to screws, instruments, implants or prosthetic devices to facilitate their tightening/loosening/removal.

The Southern Implants® Wrench Adaptor is intended to enable the connection of a torque wrench to a variety of instruments.

The Southern Implants® Retention Components are intended to be assembled with an implant insertion tool, dental implant holding endpiece or drill extension to retain a dental implant, implant fixture mount or drill during implant placement procedures.

Indications for use

The Southern Implants® Reusable Instruments are indicated for the procedures associated with the functional and aesthetic restoration of individuals indicated for implant-retained prostheses and the placement/use of Southern Implants® devices.

Specifically, the Southern Implants® Screw Extractor Guides are indicated for the removal of damaged or defective screws.

Intended user

Dental Technicians, Maxillofacial Surgeons, General Dentists, Orthodontists, Periodontists, Prosthodontists, and other appropriately trained and/or experienced medical professionals.

Intended environment

The devices are intended to be used in a dental laboratory as part of the restoration design and manufacture as well as in a clinical environment such as an operating theatre or a dentist consultation room.

Intended patient population

The Southern Implants® Reusable Instruments are intended to be used in patients subject to implant treatment/revision.

Compatibility information

Table A – Reusable Instrument Device Codes

Reusable Instrument Type	Device Product Code	
Dental Implant Depth/ Angle Gauge	Parts labelled CH-I-DG and I-ZYG-DG-1	
	Parts labelled ZYG-TR-55-35/ 45/ 52	
Surgical Screwdrivers	Handheld Hexed	Parts labelled I-CS-HD, I-HD-M
Dental Implant Holding Endpieces	Handheld Hexed connection	Parts labelled I-IMP-INS-2 and I-ZYG-INS-2 (uses replacement part I-X-RING)
Basic Dental Retractors	Parts labelled I-ZYG-RET-1	
Screw Extractor Guides	Parts labelled I-SRG-EXT-HEX	
Dental Torque Wrenches	Parts labelled I-RATCHET-2	
Wrench Adaptors	Parts labelled I-WI-SH	
Retention Components	Parts labelled I-X-RING are replacement parts for use with the I-IMP-INS-2, I-ZYG-INS-2, I-CON-X	

Clinical performance

The Reusable instruments are intended to aid in the surgical procedures associated with implant therapy. Consequently, the clinical performance of the Reusable instruments is primarily defined by the success of the implant therapy. This performance can be quantitatively assessed through implant survival/success, which serves as an important indicator of the instruments' effectiveness and success in clinical practice.

Clinical benefits

This device aids in the successful outcome of surgical procedures and has no direct clinical benefit.

Storage, cleaning and sterilisation

For components supplied sterile (sterilised by gamma irradiation), sterility is assured unless the container or seal is damaged or opened. If packaging is damaged do not use the product and contact your Southern representative or return to Southern Implants®. The devices must be stored in a dry place at room temperature and not exposed to direct sunlight. Incorrect storage may influence device characteristics.

Limitations on reusable items

A direct value for reusable instruments cannot be given. Frequent processing may have minor effects on the instruments. The product life is normally determined by wear and damage during use, thus instruments if properly cared for and inspected after each use, can be reused many times. Maintain a checklist for these instruments recording the number of uses.

Prior to reprocessing the device, it should be thoroughly inspected and tested to determine its suitability for reuse.

NOTE: during use, handle drills and instruments with sterile tweezers to minimise contamination of the instrument tray and risk of damage to sterile surgical gloves.

Containment

As soon as practically possible, remove all visible residue after use (bone, blood or tissue), by immersing the instrument in cold water (dried soil may be difficult to remove).

Precleaning

Disassemble instruments from hand pieces and all connecting parts from instruments in order to clean soil from obstructed areas. Remove the PEEK bits from placement tools. Rinse with lukewarm water for 3 minutes and remove hardened debris with a soft nylon brush. Avoid mechanical damage during cleaning.

Manual cleaning or automated cleaning

Prepare an ultrasonic bath with suitable detergent (i.e. Steritech instrument cleaner - 5% dilution), sonicate for 20 minutes (alternative methods can be used if proven by the end user). Rinse with purified/sterile water.

NOTE: always follow the instructions for use of the manufacturers for cleaning agents and disinfectants.

Load devices into a thermodisinfectant. Run the cleaning and disinfection cycle, followed by the drying cycle.

Drying

Dry the instruments both inside and outside with filtered compressed air or single use, lint free wipes. Pack the instruments as quickly as possible after removal into the storage container. If additional drying is necessary, dry in a clean location.

Inspection

Do a visual inspection of the items to check for damages on instruments.

Packaging

Use the correct packaging material as indicated for steam sterilisation to ensure sterility is maintained. Double packaging is recommended. Where appropriate, the cleaned, disinfected and checked devices can be assembled and placed in instrument trays as applicable. Instrument trays can be double wrapped or placed in sterilisation bags.

Sterilisation

Southern Implants® recommends one of the following procedures to sterilise the device prior to use:

1. prevacuum sterilisation method: steam sterilise the device at 132°C (270°F) at 180 - 220 kPa for 4 minutes. Dry for at least 20 minutes in the chamber. Only an approved wrap or pouch for steam sterilisation must be used.
2. for users in the USA: prevacuum sterilisation method: wrapped, steam sterilise at 135°C (275°F) at 180 - 220 kPa for 3 minutes. Dry for 20 minutes in the chamber. Use a wrap or pouch that is cleared for the indicated steam sterilisation cycle.

NOTE: users in the USA must ensure that the steriliser, wrap or pouch, and all steriliser accessories are cleared by the FDA, for the intended sterilisation cycle.

Storage

Maintain packaging integrity to ensure sterility in storage. Packaging should be completely dry before storage to avoid corrosion and degradation of cutting edges.

Contraindications

The contraindications of the specific target implant and prosthetic devices apply. Therefore, the contraindications of the systems/medical devices utilized as part of implant surgery/therapy should be noted and the relevant documents consulted.

Contraindications specific to the Southern Implants® Reusable Instruments include:

- Allergies or hypersensitivity to chemical ingredients in the following materials: Titanium, Aluminium, Vanadium, Stainless Steel (DIN 1.4197, SAE 316, AISI 303, AISI 302 or DIN 1.4320), Silicon, Polyether Ether Ketone (PEEK), Acetal Copolymer (POM), Acetal homopolymer (POM-H) and Polyphenylsulfone (PPSU).

Warnings and precautions

IMPORTANT NOTICE: THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING

- To ensure the safe and effective use of dental implants, new technologies/systems and the Metal Abutment devices, it is strongly recommended that specialised training be undertaken. This training should include hands-on methods to gain competency on the proper technique, the system's biomechanical requirements and radiographic evaluations required for the specific system.
- Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone.
- Use of the device with incompatible or non-corresponding devices can result in the poor performance or failure of the device.
- When handling devices intraorally, it is imperative that they are adequately secured to prevent aspiration, as aspiration of products may lead to infection or physical injury.
- The use of non-sterile items can lead to secondary infections of the tissue or transfer infectious diseases.
- Failure to adhere to appropriate cleaning, re-sterilization, and storage procedures as outlined in the Instructions for Use can result in device damage, secondary infections, or patient harm.
- Exceeding the number of recommended uses for reusable devices can result in device damage, secondary infection or patient harm.
- The use of blunt drills may cause damage to the bone, potentially compromising osseointegration.

It is crucial to emphasize that training should be undertaken by both new and experienced implant users prior to utilizing a new system or attempting a new treatment method.

Patient Selection and Preoperative Planning

A comprehensive patient selection process and meticulous preoperative planning are essential for successful implant treatment. This process should involve consultation among a multidisciplinary team, including well-trained surgeons, restorative dentists, and laboratory technicians.

Patient screening should include, at a minimum, a thorough medical and dental history, as well as visual and radiological inspections to assess the presence of adequate bone dimensions, the positioning of anatomical landmarks, the presence of unfavourable occlusal conditions, and the periodontal health status of the patient.

For successful implant treatment, it is important to:

1. Minimize trauma to the host tissue, thereby increasing the potential for successful osseointegration.
2. Accurately identify measurements relative to radiographic data, as failure to do so may lead to complications.
3. Be vigilant in avoiding damage to vital anatomical structures, such as nerves, veins, and arteries. Injury to these structures may result in serious complications, including ocular injury, nerve damage, and excessive bleeding.

The responsibility for proper patient selection, adequate training and experience in implant placement, and the provision of appropriate information required for informed consent rests with the practitioner. By combining thorough screening of prospective implant candidates with a practitioner possessing a high level of competence in the use of the system, the potential for complications and severe side effects can be significantly reduced.

High Risk Patients

Special care should be exercised when treating patients with local or systemic risk factors that may adversely affect the healing of bone and soft tissue or otherwise increase the severity of side effects, the risk of complications, and/or the likelihood of implant failure. Such factors include:

- poor oral hygiene
- history of smoking/vaping/tobacco use
- history of periodontal disease
- history of orofacial radiotherapy**
- bruxism and unfavourable jaw relations

- use of chronic medications that may delay healing or increase the risk of complications including, but not limited to, chronic steroid therapy, ant-coagulant therapy, TNF- α blockers, bisphosphonate and cyclosporin

*** The potential for implant failure and other complications increases when implants are placed in irradiated bone, as radiotherapy can lead to progressive fibrosis of blood vessels and soft tissue (i.e., osteoradionecrosis), resulting in diminished healing capacity. Contributing factors to this increased risk include the timing of implant placement in relation to radiation therapy, the proximity of radiation exposure to the implant site, and the radiation dosage at that site.*

Should the device not operate as intended, it must be reported to the manufacturer of the device. The contact information or the manufacturer of this device to report a change in performance is: sicomplaints@southernimplants.com.

Side effects

The clinical outcome of treatment is influenced by various factors. The following side effects and residual risks are associated with the device group and may necessitate further treatment, revision surgery or additional visits to the relevant medical professional's office. Furthermore, these side effects and residual risks can occur with varying possible severities and frequencies.

- Pain, discomfort or tenderness
- Localized inflammation
- Soft tissue irritation
- Hematoma or bruising
- Minor bleeding
- Allergic or hypersensitivity reaction(s)
- Bruising
- Dental injury during surgery
- Fistula formation
- Gingival injury
- Infection (acute and/or chronic)
- Localized inflammation
- Loss or damage to adjacent dentition
- Pain or discomfort
- Peri-implantitis, peri-mucositis or otherwise poor peri-implant soft tissue health
- Soft tissue irritation
- Surgical side effects such as pain, inflammation, bruising and mild bleeding
- Wound dehiscence or poor healing

Precaution: maintaining sterility protocol

Devices supplied sterile are packed in a peel pouch or blister base with a "peel-back" lid. Labelling information is located on the bottom half of the pouch, inside the packet or on the surface of the peel-back lid. Sterility is assured unless the pouch is damaged or opened.

Devices supplied non-sterile are supplied clean but not sterile in a peel pouch, blister base with peel-back lid, zipper bag, or tin box with sponge insert. Labelling information is located on the bottom half of the pouch or on the surface of the peel-back lid.

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established.

The contact information for the manufacturer of this device to report a serious incident is as follows:
sicomplaints@southernimplants.com.

Materials

Material type Stainless Steel (DIN 1.4197, SAE 316/ 316L/ 303/ 431/ 304), Stainless Steel Spring Wire (AISI 302 or DIN 1.4320), Commercially pure titanium (Grade 2, 3, 4) according to ASTM F67 and ISO 5832-2, Titanium Alloy Ti-6Al-4V (90% titanium, 6% aluminium, 4% vanadium) according to ASTM F136 and ISO 5832-2, Aluminium, Polyphenylsulfone (PPSU)

Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

Summary of Safety and Clinical Performance (SSCP)

As required by the European Medical Device Regulation (MDR; EU2017/745), a Summary of Safety and Clinical Performance (SSCP) is available for perusal with regard to Southern Implants® product ranges.

The relevant SSCP can be accessed at <https://ec.europa.eu/tools/eudamed>.

NOTE: the above website will be available upon the launch of the European Database on Medical Devices (EUDAMED).

Disclaimer of liability

This product is part of the Southern Implants® product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants® product range and take full responsibility for the correct indications and use of this product. Southern Implants® does not assume liability for damage due to incorrect use. Please note that some Southern Implants® products may not be cleared or released for sale in all markets.


Basic UDI

Product	Basic-UDI Number
Basic-UDI for Reusable Instruments	6009544038769E

Related literature and catalogues

CAT-8049-STR - IFU - Straumann® ZAGA™ Zygomatic Implants
 CAT-8047-STR - IFU - Straumann® ZAGA™ Screw-Retained Abutments
 CAT-8051-STR - IFU - Straumann® ZAGA™ Cover Screws
 CAT-8048-STR - IFU - Straumann® ZAGA™ Drills & Handpiece Instruments
 CAT-8082-STR - IFU - Straumann® ZAGA™ Instrument Trays

Symbols and warnings

* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

All rights reserved. Southern Implants®, the Southern Implants® logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants®. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale. It is the responsibility of the clinician to inspect the symbols that appear on the packaging of the product in use.

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