|                 |                                                                    |
|-----------------|--------------------------------------------------------------------|
| English         | INSTRUCTIONS FOR USE: Straumann® ZAGA™ Instrument Trays            |
| Español         | INSTRUCCIONES DE USO: Bandejas para Instrumental ZAGA™ Straumann®  |
| Italiano        | ISTRUZIONI PER L'USO: Straumann® ZAGA™ Vassoi per Strumenti        |
| Français        | MODE D'EMPLOI : Plateaux d'instruments Straumann® ZAGA™            |
| Deutsch         | GEBRAUCHSANWEISUNGEN: Straumann® ZAGA™ Instrumententrays           |
| Português       | INSTRUÇÕES DE UTILIZAÇÃO: Implantes zigomáticos Straumann® ZAGA™   |
| Ελληνικά        | ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ: Δίσκοι οργάνων ZAGA™ της Straumann®                |
| Svenska         | ANVISNINGAR FÖR ANVÄNDNING: Straumann® ZAGA™ Instrumentbrickor     |
| Brazilian       | INSTRUÇÕES DE USO: Bandejas de instrumentos Straumann® ZAGA        |
| Hrvatski        | UPUTE ZA UPOTREBU: Straumann® ZAGA™ posude za instrumente          |
| Dansk           | BRUGSANVISNING: Straumann® ZAGA™ Instrumentbakker                  |
| Suomi           | KÄYTTÖOHJEET: Straumann® ZAGA™ -instrumenttialustat                |
| Nederlands (NL) | GEbruiksaanwijzing: Straumann® ZAGA™ instrumentenbakken            |
| Magyar          | HASZNÁLATI ÚTMUTATÓ: Straumann® ZAGA™ Műszertálcák                 |
| Lietuviškai     | NAUDOJIMO NURODYMAI: Straumann® ZAGA™ instrumentų dėklai           |
| Polski          | INSTRUKCJA UŻYTKOWANIA: Straumann® ZAGA™ tace na instrumenty       |
| Română          | INSTRUCȚIUNI DE UTILIZARE: Straumann® ZAGA™ Tăvi de instrumente    |
| Slovenščina     | NAVODILA ZA UPORABO: Pladnji za instrumente Straumann® ZAGA™       |
| Українська      | ІНСТРУКЦІЇ З ЕКСПЛУАТАЦІЇ: Лотки для інструментів Straumann® ZAGA™ |
| Русский         | ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ: Лотки для инструментов Straumann® ZAGA™  |
| Български       | ИНСТРУКЦИИ ЗА УПОТРЕБА: Тави за инструменти Straumann® ZAGA™       |
| Česky           | NÁVOD K POUŽITÍ: Podnosy na nástroje Straumann® ZAGA™              |
| Eesti           | KASUTUSJUHEND: Straumann® ZAGA™ instrumendialused                  |
| Íslenska        | NOTKUNARLEIÐBEININGAR: Straumann® ZAGA™ hljóðfærabakkar            |
| Latviski        | LIETOŠANAS INSTRUKCIJA: Straumann® ZAGA™ instrumentu paplātes      |
| Norsk           | BRUKSANVISNINGER: Straumann® ZAGA™ Instrumentbrett                 |
| Slovenčina      | NÁVOD NA POUŽITIE: Zásobníky nástrojov Straumann® ZAGA™            |
| Türkçe          | KULLANMA TALİMATI: Straumann® ZAGA™ Alet Tepsileri                 |
| Српски          | УПУТСТВО ЗА УПОТРЕБУ: Straumann® ZAGA™ Подметачи за инструмент     |

**South Africa - Headquarters:** 1 Albert Road, Irene, 0062, RSA  
T: +27-12-667-1046 | E: info@southernimplants.com

#### Subsidiaries

**United Kingdom and Ireland**  
Southern Implants UK  
T: +44-20 8059 4490  
E: info@southernimplants.co.uk

**USA and Canada**  
Southern Implants North America Inc.  
T: +1-561-472-0990  
E: customercare@southernimplants.com

EC REP

**Southern Implants Europe AB:** Holmgatan 30, S-791 71 Falun, Sweden  
T: +46 23 13300 | E: ecrep@southernimplants.com

CH REP

**MedEnvoy Switzerland:** Gotthardstrasse 28, 6302 Zug, Switzerland

## Description

The Southern Implants® Instrument Trays are reusable, multi-piece rigid containers or organizing trays designed for use in healthcare facilities to hold various surgical drills and instruments. Their primary purpose is to organize, contain, and transport non-porous, reusable medical devices for sterilization between uses. Each tray can consist of multiple components, including a base tray, a locking lid and an internal individualized insert tray, which are designed to integrate into a single unit that securely contains and protects the interior components during steam sterilization.

The components of the Southern Implants® Instrument Trays are perforated and/or have sufficient air gaps to facilitate steam sterilization. Additionally, select insert trays features retention sockets (grommets) made from silicone to allow steam penetration. The Instrument Trays are available in various dimensions enabling them to fit in any standard autoclave.

The lid is constructed from Radel R-5800 (a polymer resin). The base and insert trays are made from Radel R-5000.

The instruments to be sterilized in the proposed tray, are all non-porous devices and include dental surgical drills and tools. The instrument trays are intended to be reusable, however they cannot be repaired or serviced if damaged. Since the trays are perforated, an FDA cleared sterilization wrap must be used for sterilization purposes to maintain the sterility of the contents. The trays are reusable and the tray material allows repeated sterilization cycles.

## Intended use

The Southern Implants® Instrument Tray is a reusable medical device designed to hold various surgical drills and tools in order to organise, steam sterilise, and transport them between uses.

## Indications for use

The Southern Implants® Instrument Tray is indicated for the transport or sterilization of various reusable devices where it is covered in a steam sterilisable wrap and sterilized according to one of the following cycles:

- Pre-vacuum steam – At 132°C for 4 minutes with a 20 minute dry time.
- Pre-vacuum steam – At 135°C for 3 minutes with a 20 minute dry time.

## Intended user

Maxillofacial Surgeons, General Dentists, Orthodontists, Periodontists, Prosthodontists, and other appropriately trained and/or experienced medical professionals.

## Intended environment

The devices are intended to be used in a clinical environment such as an operating theatre or a dentist consultation room.

## Intended patient population

The Southern Implants® Instrument Trays do not have direct patient interaction.

## Compatibility information

The Southern Implants® Instrument Trays can be identified by specific abbreviations in the product codes. Range identifiers are summarised in Table A.

**Table A – Instrument Trays**

| Product Code | Corresponding System                       |
|--------------|--------------------------------------------|
| CH-I-ZYG     | Straumann Zygomatic Implant System         |
| I-ZYG-1      | Southern Implants Zygomatic Implant System |
| I-PROS-EG    | Applicable to all Southern Implant systems |

## Clinical performance

The Southern Implants® Instrument Trays aid in the successful organisation, sterilization, and transportation of reusable devices by having an appropriate geometry and material.

### Clinical benefits

The Southern Implants® Instrument Trays provide no direct clinical benefit for a patient; however, the instrument trays are able to provide effective sterilization of reusable instruments and improve air circulation which reduces moisture buildup on instruments and inhibit bacterial growth. The tray also protects the instruments from external damage which could affect their function and performance.

### Storage, cleaning and sterilisation

#### Limitations on reusable items

A direct value for reusable instruments cannot be given. Frequent processing may have minor effects on the instruments. The product life is normally determined by wear and damage during use, thus instruments if properly cared for and inspected after each use, can be reused many times. Maintain a checklist for these instruments recording the number of uses.

Prior to reprocessing the device, it should be thoroughly inspected and tested to determine its suitability for reuse.

**NOTE:** during use, handle drills and instruments with sterile tweezers to minimise contamination of the instrument tray and risk of damage to sterile surgical gloves.

#### Containment

As soon as practically possible, remove all visible residue after use (bone, blood or tissue), by immersing the instrument in cold water (dried soil may be difficult to remove).

#### Precleaning

Disassemble instruments from hand pieces and all connecting parts from instruments in order to clean soil from obstructed areas. Remove the PEEK bits from placement tools. Rinse with lukewarm water for 3 minutes and remove hardened debris with a soft nylon brush. Avoid mechanical damage during cleaning.

#### Manual cleaning or automated cleaning

Prepare an ultrasonic bath with suitable detergent (i.e. Steritech instrument cleaner - 5% dilution), sonicate for 20 minutes (alternative methods can be used if proven by the end user). Rinse with purified/sterile water.

**NOTE:** always follow the instructions for use of the manufacturers for cleaning agents and disinfectants.

Load devices into a thermodisinfectant. Run the cleaning and disinfection cycle, followed by the drying cycle.

#### Drying

Dry the instruments both inside and outside with filtered compressed air or single use, lint free wipes. Pack the instruments as quickly as possible after removal into the storage container. If additional drying is necessary, dry in a clean location.

#### Inspection

Do a visual inspection of the items to check for damages on instruments.

#### Packaging

Use the correct packaging material as indicated for steam sterilisation to ensure sterility is maintained. Double packaging is recommended. Where appropriate, the cleaned, disinfected and checked devices can be assembled and placed in instrument trays as applicable. Instrument trays can be double wrapped or placed in sterilisation bags.

#### Sterilisation

Southern Implants® recommends one of the following procedures to sterilise the device prior to use:

1. prevacuum sterilisation method: steam sterilise the device at 132°C (270°F) at 180 - 220 kPa for 4 minutes. Dry for at least 20 minutes in the chamber. Only an approved wrap or pouch for steam sterilisation must be used.

- for users in the USA: prevacuum sterilisation method: wrapped, steam sterilise at 135°C (275°F) at 180 - 220 kPa for 3 minutes. Dry for 20 minutes in the chamber. Use a wrap or pouch that is cleared for the indicated steam sterilisation cycle.

**NOTE:** users in the USA must ensure that the steriliser, wrap or pouch, and all steriliser accessories are cleared by the FDA, for the intended sterilisation cycle.

### Storage

Maintain packaging integrity to ensure sterility in storage. Packaging should be completely dry before storage to avoid corrosion and degradation of cutting edges.

### Contraindications

The Southern Implants® Instrument Trays are contraindicated for the sterilisation of non-porous loads. Furthermore, it is not recommended that the Instrument Trays be stacked during sterilization.

The Southern Implants® Instrument Trays do not come into contact with the patient and, therefore, no specific patient-related contraindications apply. However, since these devices are intended to be used as part of implant therapy/surgery, the contraindications of the associated systems/medical devices should be noted and the relevant documents consulted.

### Warnings and precautions

**IMPORTANT NOTICE:** THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING

- To ensure the safe and effective use of dental implants, new technologies/systems and the Metal Abutment devices, it is strongly recommended that specialised training be undertaken. This training should include hands-on methods to gain competency on the proper technique, the system's biomechanical requirements and radiographic evaluations required for the specific system.
- Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone.
- Use of the device with incompatible or non-corresponding devices can result in the poor performance or failure of the device.
- When handling devices intraorally, it is imperative that they are adequately secured to prevent aspiration, as aspiration of products may lead to infection or physical injury.
- The use of non-sterile items can lead to secondary infections of the tissue or transfer infectious diseases.
- Failure to adhere to appropriate cleaning, re-sterilization, and storage procedures as outlined in the Instructions for Use can result in device damage, secondary infections, or patient harm.
- Exceeding the number of recommended uses for reusable devices can result in device damage, secondary infection or patient harm.
- The use of blunt drills may cause damage to the bone, potentially compromising osseointegration.

It is crucial to emphasize that training should be undertaken by both new and experienced implant users prior to utilizing a new system or attempting a new treatment method.

### Patient Selection and Preoperative Planning

A comprehensive patient selection process and meticulous preoperative planning are essential for successful implant treatment. This process should involve consultation among a multidisciplinary team, including well-trained surgeons, restorative dentists, and laboratory technicians.

Patient screening should include, at a minimum, a thorough medical and dental history, as well as visual and radiological inspections to assess the presence of adequate bone dimensions, the positioning of anatomical landmarks, the presence of unfavourable occlusal conditions, and the periodontal health status of the patient.

For successful implant treatment, it is important to:

1. Minimize trauma to the host tissue, thereby increasing the potential for successful osseointegration.
2. Accurately identify measurements relative to radiographic data, as failure to do so may lead to complications.
3. Be vigilant in avoiding damage to vital anatomical structures, such as nerves, veins, and arteries. Injury to these structures may result in serious complications, including ocular injury, nerve damage, and excessive bleeding.

The responsibility for proper patient selection, adequate training and experience in implant placement, and the provision of appropriate information required for informed consent rests with the practitioner. By combining thorough screening of prospective implant candidates with a practitioner possessing a high level of competence in the use of the system, the potential for complications and severe side effects can be significantly reduced.

### High Risk Patients

Special care should be exercised when treating patients with local or systemic risk factors that may adversely affect the healing of bone and soft tissue or otherwise increase the severity of side effects, the risk of complications, and/or the likelihood of implant failure. Such factors include:

- poor oral hygiene
- history of smoking/vaping/tobacco use
- history of periodontal disease
- history of orofacial radiotherapy\*\*
- bruxism and unfavourable jaw relations
- use of chronic medications that may delay healing or increase the risk of complications including, but not limited to, chronic steroid therapy, anti-coagulant therapy, TNF- $\alpha$  blockers, bisphosphonate and cyclosporin

*\*\* The potential for implant failure and other complications increases when implants are placed in irradiated bone, as radiotherapy can lead to progressive fibrosis of blood vessels and soft tissue (i.e., osteoradionecrosis), resulting in diminished healing capacity. Contributing factors to this increased risk include the timing of implant placement in relation to radiation therapy, the proximity of radiation exposure to the implant site, and the radiation dosage at that site.*

Should the device not operate as intended, it must be reported to the manufacturer of the device. The contact information or the manufacturer of this device to report a change in performance is: [sicomplaints@southernimplants.com](mailto:sicomplaints@southernimplants.com).

### Side effects

Since the Instrument Trays do not come into direct contact with the patient, they do not pose inherent side effects or residual risks. However, as these devices are integral to surgical procedures, the typical side effects associated with such procedures may still be expected, including:

1. Pain, discomfort or tenderness
2. Localized inflammation
3. Soft tissue irritation
4. Hematoma or bruising
5. Minor bleeding

Additionally, due to the intended use of the Instrument Trays for the re-sterilization of reusable surgical instruments, a limited set of contamination-related complications may occur. These include persistent pain, inflammation, infection, peri-implantitis, peri-implant mucositis, and fistula formation

### Precaution: maintaining sterility protocol

The instrument trays are supplied clean but not sterile in a cardboard box. Labelling information is located on the box.

### Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established.

The contact information for the manufacturer of this device to report a serious incident is as follows:  
[sicomplaints@southernimplants.com](mailto:sicomplaints@southernimplants.com).

## Materials

|                       |                        |
|-----------------------|------------------------|
| Lid:                  | Radel R-5800           |
| Base and Insert Tray: | Radel R-5000           |
| Grommets:             | Medical Grade Silicone |

## Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

## Summary of Safety and Clinical Performance (SSCP)

As required by the European Medical Device Regulation (MDR; EU2017/745), a Summary of Safety and Clinical Performance (SSCP) is available for perusal with regard to Southern Implants® product ranges.

The relevant SSCP can be accessed at <https://ec.europa.eu/tools/eudamed>.

**NOTE:** the above website will be available upon the launch of the European Database on Medical Devices (EUDAMED).

## Disclaimer of liability

This product is part of the Southern Implants® product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants® product range and take full responsibility for the correct indications and use of this product. Southern Implants® does not assume liability for damage due to incorrect use. Please note that some Southern Implants® products may not be cleared or released for sale in all markets.


## Basic UDI

| Product                        | Basic-UDI Number |
|--------------------------------|------------------|
| Basic-UDI for Instrument Trays | 6009544039118P   |

## Related literature and catalogues

CAT-8049-STR - IFU - Straumann® ZAGA™ Zygomatic Implants  
 CAT-8047-STR - IFU - Straumann® ZAGA™ Screw-Retained Abutments  
 CAT-8051-STR - IFU - Straumann® ZAGA™ Cover Screws  
 CAT-8048-STR - IFU - Straumann® ZAGA™ Drills & Handpiece Instruments  
 CAT-8080-STR - IFU - Straumann® ZAGA™ Reusable Instruments  
 CAT-8083-STR - IFU - Straumann® ZAGA™ General Use Instruments

## Symbols and warnings

\* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licensed physician or dentist.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

All rights reserved. Southern Implants®, the Southern Implants® logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants®. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale. It is the responsibility of the clinician to inspect the symbols that appear on the packaging of the product in use.

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