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
South Africa - Headquarters: 1 Albert Road, Irene, 0062, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

Southern Implants Europe AB: Holmgatan 30, S-791 71 Falun, Sweden
T: +46 23 13300 | E: ecrep@southernimplants.com

Subsidiaries

Australia

Southern Implants Australia
T: +61-(0)-8-9466-2627
E: info@southernimplants.com.au

Spain and Portugal

Southern Implants Iberica
T: +34 935 053 507
E: info@southernimplants.es

United Kingdom and Ireland

Southern Implants UK
T: +44-20-8899-6845 / 6 / 7
E: info@southernimplants.co.uk

USA and Canada

Southern Implants North America Inc.
T: +1-561-472-0990
E: customercare@southernimplants.com

Description

Osteotomy Guide

An osteotomy guide is a non-sterile surgical instrument intended to be attached to a portion of bone (e.g. pelvis, femur, tibia, ulna) during an osteotomy procedure to outline the desired cut-line and mechanically support/guide the cutting device (e.g. oscillating cutting saw).

Reusable Dental Implant Depth/Angle Gauge

A reusable dental implant depth/angle gauge is used during the planning and post-osteotomy stages of a dental implantation procedure, prior to surgical implant placement/fixation, to measure various lengths (e.g. drill depth of the osteotomy cavity, mucosal height, distance between implants/teeth for determining the mesiodistal position of the implant) and/or angles (e.g. direction of the perforation in relation to the antagonist arch), to aid in selection of an appropriate implant and/or final prosthetic abutment to be installed. It is typically graduated and may be a hand-held instrument or a stand-alone device placed intraorally; some types are intended for use during x-ray imaging.

Dental Osteotomes

A hand-held surgical instrument used for shaping and condensing bone in a dental osteotomy, typically when performing a dental implant surgical procedure. It is held by the user who will use a surgical mallet or hammer to manually impart an impacting force to the proximal end of the instrument. The device has a straight or bayonet-shaped handle leading to a round or concave tip at the distal end, which may be marked with graduations to indicate depth differences. It is used for ridge expanding, bone condensing and sinus floor elevation procedures. This device is made of high-grade stainless steel.

Basic Dental Retractors

A hand-held, non-self-retaining dental instrument intended to be used to manually displace soft tissues of the oral cavity to improve their visualization and access, and to afford them protection during oral surgical procedures. It is a metallic device and includes a shaped/hooked blade with a handle/holding portion.

Surgical Drill Guide, Reusable

A hand-held tubular surgical device designed to simultaneously place, angle, and guide a rotating surgical drill on hard tissue (e.g., bone), and to protect the surrounding soft tissue and dentition while the drill is in operation. It can have a handle at its proximal end for holding and manipulation, or it may be held by another device (e.g. orthopaedic jig, 3D-printed surgical guide). The device is dedicated to drill bit size, made of metal, and may have graduations along its length; it is not custom-made for a specific patient.

Dental Torque Wrenches

A hand-held dental tool intended for the manual application of rotation to tighten/loosen/remove a device (e.g. dental implant, spatial register (localizer/tracker/calibrator)) fitted into or close to the oral cavity, and/or to a surgical instrument, in association with a dental surgical procedure (e.g. dental guided surgery registration). It may have a built-in ratchet function, with forward and reverse direction, and a calibrated scale to which a recommended torque can be set; it is not a screwdriver. It is typically made of metallic materials and may have a standard socket connection to accommodate different tool bits/adaptors.

Wrench Adaptors

A component of a surgical wrench (typically ratchet type) intended to enable connection of the main body of the wrench to a variety of distal end pieces (bits). It is intended to be used during an invasive surgical procedure, and is not dedicated to dental use only. It may be referred to as a converter.

Guide Stabilisers

A device used to retain a surgical guide to a placed implant, via an integral screw, in order to secure the surgical guide from shifting position while additional implant osteotomies are prepared and implants placed.

Intended use

Southern Implants' Reusable Surgical Instruments are intended to be used as part of surgical procedures and protocols to aid in the treatment of partially or fully edentulous patients eligible for placement of one or more dental implants as a means of fixing a permanent or removable single crown, partial or full-arch dental prosthesis in the upper or lower jaw.

The system constituents are classified as medical devices and are multiple use devices.

Osteotomy Guide

The Osteotomy Guide allows for improved preparation of the procedural site by outlining the desired cutline and mechanically supporting or guiding the cutting device.

Reusable Dental Implant Depth/Angle Gauge

The Reusable Dental Implant Depth/Angle Gauge allows for improved selection of implants and/or the final prosthetic abutment-based geometry for the specific osteotomy.

Dental Osteotomes

The dental osteotome is a hand-held surgical instrument used to shape and condense the bone during implant site preparation allowing for the successful achievement of ridge expansion, bone condensing and sinus floor lift procedures.

Basic Dental Retractors

The Basic Dental Retractors are used to manually displace the soft tissue of the oral cavity to improve visualisation and access to the implantation site as well as afford them protection during the surgical procedure.

Surgical Drill Guide, Reusable

The Surgical Drill Guide designed to simultaneously place, angle, and guide a rotating surgical drill on hard tissue (e.g., bone), and to protect the surrounding soft tissue and dentition while the drill is in operation.

Dental Torque Wrenches

The Dental torque wrench enables the manual application of rotation to tighten, loosen, or remove a surgical instrument and/or a device (e.g., dental implant) fitted into or close to the oral cavity in association with a dental surgical procedure.

Wrench Adaptors

The wrench adapter enables the compatibility of the Dental Torque Wrench with various distal endpieces.

Guide Stabilisers

The guide stabilisers aid in securing a surgical guide to prevent unintentional misalignment or repositioning during dental implant placement and surgical procedures.

Indications for use

Southern Implants' Reusable Surgical Instruments are reusable devices intended to assist with surgical procedures and facilitate the insertion and placement of dental implants and restorations by clinicians. Southern Implants dental implants are intended for both one and two-stage surgical procedures in the following situations and with the following clinical protocols:

- replacing single and multiple missing teeth in the mandible and maxilla.
- immediate placement in extraction sites and in situations with a partially or completely healed alveolar ridge.
- immediate loading in all indications, except in single tooth situations on implants shorter than 8 mm or in soft bone (type IV) where implant stability may be difficult to obtain, and immediate loading may not be appropriate.

Intended user

The intended user for the Reusable Surgical Instruments is not dissimilar to that of dental implant therapy.

The intended user for the Reusable Surgical Instruments, with respect to dental implant therapy, includes Dental Technicians, Maxillofacial Surgeons, General Dentists, Orthodontists, Periodontist, Prosthodontists and other appropriately trained and experienced implant users.

Intended environment

The intended environment for the Reusable Surgical Instruments device family is not dissimilar to that of dental implant therapy.

The intended environment of the Reusable Surgical Instruments system, with respect to dental implant therapy, includes a dental laboratory for the manufacture of the restoration and a clinical environment such as an operating theatre or a dentist consultation room.

Intended patient population

The intended patient population for the device is not dissimilar to that of dental implant therapy.

The intended patient population for implant therapy is partially or fully edentulous patients requiring prosthetic dental restoration in the upper or lower jaw. Restorations may comprise single teeth, partial or full bridges and are fixed restorations.

Additionally, the intended population for some device variants includes those who have previously undergone dental implant therapy.

Device code information

The reusable instrument can be identified by specific abbreviations in the product codes. Device identifiers are summarised in Table A.

Table A – Reusable instrument device codes

Reusable Instrument Type	Device Product Code
Osteotomy Guides	Parts labelled CG-F-RT-30/ 40/ 50 and I-FH-RT are osteotomy guides for Rohner Technique systems
Dental Implant Depth/Angle Gauge, Reusable	Part labelled I-DG are depth gauges for general dental implants
	Part labelled I-ADG and I-DI are angle gauges for general dental implants
	Parts labelled CH-I-DG and I-ZYG-DG-1 are depth gauges for zygomatic implants
	Parts labelled I-ZYG-DI55 and ZYG-TR-55-35/ 45/ 52.5/ 52 are angle gauges for zygomatic implants
	Parts labelled MAX-(Ø)-PG are angle gauges for MAX dental implants
Dental Osteotomes	Parts labelled I-OST-(Ø)C/ P/ SU, I-FSD and I-HAMMER are for general dental implants
Basic Dental Retractors	Parts labelled I-NIC/ -L are for general dental implants
	Parts labelled I-ZYG-RET are for zygomatic implants
	Parts labelled I-PT-RET are for pterygoid implants
Surgical Drill Guide, Reusable	Parts labelled DGS-L-(Ø) and DGS-S-(Ø) are drill guide spoons for general dental implants
	Parts labelled GSDS-(Ø) and I-TDS-(Ø) are drill stabilisers for general dental implants
	Parts labelled TDF-RT-(Ø) and ICF-RT are surgical guides for Rohner Technique systems
Dental Torque Wrenches	Parts labelled I-RATCHET-2, I-TW/ TWS/ TWA are torque wrenches for general dental implants
Wrench Adaptors	Parts labelled I-WI-C-L/ S, I-WI-CNB/ CST and I-WI-SH/ SL/ SS are wrench adaptors for general dental implants
Guide Stabilisers	Parts labelled I-GSN-SP and I-GS-SP for Single Platform (SP1) implants

Clinical benefits

The clinical benefits of the Reusable Surgical Instruments are similar to those of the General Dental Implants and Zygomatic Implants (i.e., dental implant therapy).

The clinical benefits of dental implant therapy include improved chewing function, speech, aesthetics and patient psychological wellbeing. Through this procedure patients can expect to have their missing teeth replaced and/or crowns restored.

Osteotomy Guide

The Osteotomy Guide allows for improved preparation of the procedural site by outlining the desired cutline and mechanically supporting or guiding the cutting device.

Reusable Dental Implant Depth/Angle Gauge

The Reusable Dental Implant Depth/Angle Gauge allows for improved selection of implants and/or the final prosthetic abutment-based geometry for the specific osteotomy.

Dental Osteotomes

The Dental Osteotomes allow for the successful achievement of ridge expansion, bone condensing and sinus floor lift procedures.

Basic Dental Retractors

The Basic Dental Retractors provide improved protection of the soft tissue of the mouth and surrounding areas during oral surgical procedures.

Surgical Drill Guide, Reusable

Surgical Drill Guides provide improved protection of the surrounding soft tissue and dentition while a surgical drill is in operation.

Dental Torque Wrenches

The Dental Torque Wrench allows for the successful achievement of tightening, loosening, or removal of a surgical instrument and/or a device (e.g., dental implant) fitted into or close to the oral cavity.

Wrench Adaptors

The wrench adapter enables the Dental Torque Wrench to be used with various distal end pieces.

Guide Stabilisers

The Guide Stabilisers provide improved alignment of the surgical guide and improved protection of the surrounding soft tissue and dentition while a surgical drill is in operation.

Storage, cleaning and sterilisation

Limitations on reusable items

A direct value for reusable instruments cannot be given. Frequent processing may have minor effects on the instruments. The product life is normally determined by wear and damage during use, thus instruments if properly cared for and inspected after each use, can be reused many times. Maintain a checklist for these instruments recording the number of uses.

Prior to reprocessing the device it should be thoroughly inspected and tested to determine its suitability for reuse.

NOTE: during use, handle drills and instruments with sterile tweezers to minimise contamination of the instrument tray and risk of damage to sterile surgical gloves.

Containment

As soon as practically possible, remove all visible residue after use (bone, blood or tissue), by immersing the instrument in cold water (dried soil may be difficult to remove).

Precleaning

Disassemble instruments from hand pieces and all connecting parts from instruments in order to clean soil from obstructed areas. Remove the PEEK bits from placement tools. Rinse with lukewarm water for 3 minutes and remove hardened debris with a soft nylon brush. Avoid mechanical damage during cleaning.

Manual cleaning or automated cleaning

Prepare an ultrasonic bath with suitable detergent (i.e. Steritech instrument cleaner - 5% dilution), sonicate for 20 minutes (alternative methods can be used if proven by the end user). Rinse with purified/sterile water.

NOTE: always follow the instructions for use of the manufacturers for cleaning agents and disinfectants.

Load devices into a thermidisinfector. Run the cleaning and disinfection cycle, followed by the drying cycle.

Drying

Dry the instruments both inside and outside with filtered compressed air or single use, lint free wipes. Pack the instruments as quickly as possible after removal into the storage container. If additional drying is necessary, dry in a clean location.

Inspection

Do a visual inspection of the items to check for damages on instruments.

Packaging

Use the correct packaging material as indicated for steam sterilisation to ensure sterility is maintained. Double packaging is recommended. Where appropriate, the cleaned, disinfected and checked devices can be assembled and placed in instrument trays as applicable. Instrument trays can be double wrapped or placed in sterilisation bags.

Sterilisation

Southern Implants® recommends one of the following procedures to sterilise the restoration prior to use:

1. prevacuum sterilisation method: steam sterilise the abutments at 132°C (270°F) at 180 - 220 kPa for 4 minutes. Dry for at least 20 minutes in the chamber. Only an approved wrap or pouch for steam sterilisation must be used.
2. for users in the USA: prevacuum sterilisation method: wrapped, steam sterilise at 135°C (275°F) at 180 - 220 kPa for 3 minutes. Dry for 20 minutes in the chamber. Use a wrap or pouch that is cleared for the indicated steam sterilisation cycle.

NOTE: users in the USA must ensure that the steriliser, wrap or pouch, and all steriliser accessories are cleared by the FDA, for the intended sterilisation cycle.

Storage

Maintain packaging integrity to ensure sterility in storage. Packaging should be completely dry before storage to avoid corrosion and degradation of cutting edges.

Contraindications

The contraindications for this system are not dissimilar to those of dental implant therapy.

Contraindications to implant therapy include:

- patients medically unfit for oral surgical procedures
- where an inadequate number of implants can be placed limiting the functional support of the prosthesis
- patients under the age of 18 years
- patients with poor bone quality
- patients with blood disorders
- presence of infection at the implant site
- patients with vascular impairment
- patients with uncontrolled diabetes
- patients with drug or alcohol abuse dependencies
- patients undergoing chronic high dose steroid therapy
- patients undergoing anti-coagulant therapy

- patients with metabolic bone disease
- patients undergoing radiotherapy treatment
- patients with pure titanium, titanium alloy (Ti6Al4V), gold, palladium, polyether ether ketone (PEEK) or iridium allergies or hypersensitivities

Other than the above, there are no contraindications unique to this system.

Warnings and precautions

THESE INSTRUCTIONS ARE NOT INTENDED AS A SUBSTITUTE FOR ADEQUATE TRAINING.

- For the safe and effective use of dental implants it is strongly suggested that specialised training be undertaken, including hands-on training to learn proper technique, biomechanical requirements and radiographic evaluations.
- Products must be secured against aspiration when handled intraorally. Aspiration of products may lead to infection or unplanned physical injury.

Responsibility for proper patient selection, adequate training, experience in the placement of implants and providing appropriate information for informed consent rests with the practitioner. Improper technique can result in implant failure, damage to nerves/vessels and/or loss of supporting bone. Implant failure increases when implants are placed in irradiated bone as radiotherapy can result in progressive fibrosis of vessels and soft tissue, leading to diminished healing capacity.

It is important to be aware and avoid damage to vital structures such as nerves, veins and arteries. Injuries to vital anatomical structures may cause serious complications such as injury to the eye, nerve damage and excessive bleeding. It is essential to protect the infraorbital nerve. Failing to identify actual measurements relative to the radiographic data could lead to complications.

New and experienced implant users should do training before using a new system or attempting to do a new treatment method. Take special care when treating patients who have local or systemic factors that could affect the healing of the bone and soft tissue (i.e. poor oral hygiene, uncontrolled diabetes, are on steroid therapy, smokers, infection in the nearby bone and patients who had orofacial radiotherapy).

Thorough screening of prospective implant candidates must be performed including:

- a comprehensive medical and dental history.
- visual and radiological inspection to determine adequate bone dimensions, anatomical landmarks, occlusal conditions and periodontal health.
- bruxism and unfavourable jaw relations must be taken into account.
- proper preoperative planning with a good team approach between well trained surgeons, restorative dentists and lab technicians is essential for successful implant treatment.
- minimising the trauma to the host tissue increases the potential for successful osseointegration.
- electrosurgery should not be attempted around metal implants as they are conductive.

Should the device not operate as intended, it must be reported to the manufacturer of the device. The contact information for the manufacturer of this device to report a change in performance is: sicomplaints@southernimplants.com.

Side effects

The side effects of the use of the system are not dissimilar to those of dental implant therapy. Possible side effects to implant therapy include:

- pain
- swelling
- phonetic difficulties
- gingival inflammation

Less common but more persistent symptoms include, but are not limited to:

- allergic reaction(s) to implant and/or abutment material
- breakage of the implant and/or abutment
- loosening of the abutment screw and/or retaining screw
- infection requiring revision of the dental implant
- nerve damage resulting in permanent weakness, numbness, or pain
- histologic responses with possible macrophage and/or fibroblast involvement
- fat emboli formation
- loosening of the implant requiring revision surgery
- perforation of the maxillary sinus
- perforation of the labial and lingual plates
- bone loss possibly resulting in revision or removal of the implant.

Precaution: maintaining sterility protocol

Implants are packaged as follows:

1. An outer package consisting of a rigid, clear box which acts as protection for the inner package.
2. The inner package consisting of a blister pack (clear plastic-formed blister base with a TYVEK “peel-back” lid).
3. Within the inner package, there is a hollow tube which contains one implant suspended from a titanium ring, this ensures the implant never touches the inside of the plastic tube.
4. Labelling information is located on the surface of the peel-back lid and on the outside of the rigid box.

Care must be taken to maintain the sterility of the implant by proper opening of the packaging and handling of the implant.

1. Open the implant package in the non-sterile field, with non-sterile gloves, tear the address label to open the box.
2. With non-sterile gloves, remove the inner blister pack. Do not place the plastic box or blister pack-lid onto the sterile field. The contents of this inner package are sterile.
3. The sealed blister is to be opened by an assistant (with nonsterile gloves): remove the TYVEK lid and drop or place the sterile tube onto the sterile field, open the tube cap and attach the implant placement tool onto the implant and carefully remove from the sterile tube. Do not touch the sterile implant.

Other sterile components are packed in a peel pouch or blister base with a “peel-back” lid. Labelling information is located on the bottom half of the pouch, inside the packet or on the surface of the peel-back lid. Sterility is assured unless the pouch is damaged or opened. Non-sterile components are supplied clean but not sterile in a peel pouch or blister base with peel-back lid. Labelling information is located on the bottom half of the pouch or on the surface of the peel-back lid.

Notice regarding serious incidents

Any serious incident that has occurred in relation with the device must be reported to the manufacturer of the device and the competent authority in the member state in which the user and/or patient is established.

The contact information for the manufacturer of this device to report a serious incident is as follows: sicomplaints@southernimplants.com.

Materials

Material type	Stainless Steel (DIN 1.4197, SAE 316/ 316L/ 303/ 431/ 304), Commercially pure titanium (grade 2, 3, 4), Titanium Alloy (Ti-6Al-4V), Celeron/Tufnol, Medical Grade PEEK, Aluminium,
---------------	--

Disposal

Disposal of the device and its packaging: follow local regulations and environmental requirements, taking different contamination levels into account. When disposing of spent items, take care of sharp drills and instruments. Sufficient PPE must be used at all times.

Summary of Safety and Clinical Performance (SSCP)

As required by the European Medical Device Regulation (MDR; EU2017/745), a Summary of Safety and Clinical Performance (SSCP) is available for perusal with regard to Southern Implants® product ranges.

The relevant SSCP can be accessed at <https://ec.europa.eu/tools/eudamed>.

NOTE: the above website will be available upon the launch of the European Database on Medical Devices (EUDAMED).

Disclaimer of liability

This product is part of the Southern Implants® product range and should only be used with the associated original products and according to the recommendations as in the individual product catalogues. The user of this product has to study the development of the Southern Implants® product range and take full responsibility for the correct indications and use of this product. Southern Implants® does not assume liability for damage due to incorrect use. Please note that some Southern Implants® products may not be cleared or released for sale in all markets.

Basic UDI

Product	Basic-UDI Number
Basic-UDI for Reusable Instruments	6009544038769E

Related literature and catalogues

CAT-2004 - Tri-Nex® Implants Product Catalogue
 CAT-2005 - IT Implants Product Catalogue
 CAT-2010 - Osseointegrated Fixtures Catalogue
 CAT-2020 - External Hex Implants Product Catalogue
 CAT-2042 - Deep Conical Implants Product Catalogue
 CAT-2043 - Internal Hex Implants Product Catalogue
 CAT-2060 - PROVATA® Implants Product Catalogue
 CAT-2069 - INVERTA® Implants Product Catalogue
 CAT-2070 - Zygomatic Implants Product Catalogue
 CAT-2093 - Single Platform (SP1) Implants Product Catalogue

Symbols and warnings

											
Manufacturer: Southern Implants® 1 Albert Rd, P.O Box 605 IRENE, 0062, South Africa. Tel: +27 12 667 1046	CE mark 2797	Prescription device*	Sterilised using irradiation	Non-sterile	Use by date (mm-yy)	Do not reuse	Do not resterilise	Catalogue number	Batch code	Medical device	Authorised representative in the European Community
											
Authorised representative for Switzerland	Date of manufacture	Magnetic Resonance conditional	Magnetic Resonance safe	Single sterile barrier system with protective packaging inside	Single sterile barrier system	Consult instruction for use	Caution	Keep away from sunlight	Do not use if package is damaged		

* Prescription device: Rx only. Caution: Federal Law restricts this device to sale by or on the order of a licenced physician or dentist.

Canada licence exemption: Please note that not all products may have been licensed in accordance with Canadian law.

All rights reserved. Southern Implants®, the Southern Implants® logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Southern Implants®. Product images in this document are for illustration purposes only and do not necessarily represent the product accurately to scale. It is the responsibility of the clinician to inspect the symbols that appear on the packaging of the product in use.

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