



The Deep Conical range surgery kit (I-DC-EG) includes all instruments to perform implant surgery for the DC range of implants. For cleaning and sterilization procedure guidelines refer to CAT-1039.



**Layout Information**

1	D-RB-MS	Round Burr	25	D-DCT5013	Tapered Drill Ø5.0 x 13mm
2	D-DC20	Twist Drill Ø2.0mm	26	D-DCT5015	Tapered Drill Ø5.0 x 15mm
3	D-DC25	Twist Drill Ø2.5mm	27	D-DC29 / 34 / 39 / 49	Twist Drills - Optional
4	D-DC27	Twist Drill Ø2.7mm	28	D-CSS-M	Countersink - Optional
5	D-DC32	Twist Drill Ø3.2mm	29	I-HDC3-S / M	Handpiece Insert Ø3.0 Implant, Short / Medium
6	D-DC37	Twist Drill Ø3.7mm	30	I-HDC4-S / M	Handpiece Insert Ø3.5 / Ø4.0 Implant, Short / Medium
7	D-DC42	Twist Drill Ø4.2mm	31	I-HDC5-S / M	Handpiece Insert Ø5.0 Implant, Short / Medium
8	D-DC47	Twist Drill Ø4.7mm	32	I-WI-DC3-S / M	Wrench Insert Ø3.0 Implant, Short / Medium
9	D-DCT3009	Tapered Drill Ø3.0 x 9mm	33	I-WI-DC4-S / M	Wrench Insert Ø3.5 / Ø4.0 Implant, Short / Medium
10	D-DCT3011	Tapered Drill Ø3.0 x 11mm	34	I-WI-DC5-S / M	Wrench Insert Ø5.0 Implant, Short / Medium
11	D-DCT3013	Tapered Drill Ø3.0 x 13mm	35	I-WI-C-S / L	Converter to Wrench, Short / Long
12	D-DCT3508	Tapered Drill Ø3.5 x 8mm	36	I-DE-MN	Drill Extension
13	D-DCT3509	Tapered Drill Ø3.5 x 9mm	37	I-FME-XS / M / L	Fixture Mount Extension, Short / Medium / Long
14	D-DCT3511	Tapered Drill Ø3.5 x 11mm	38	I-CON-X / XS	Connector to Handpiece / Short
15	D-DCT3513	Tapered Drill Ø3.5 x 13mm	39	I-DI	Direction Indicator
16	D-DCT3515	Tapered Drill Ø3.5 x 15mm	40	I-DI12d	Direction Indicator, 12° Angled
17	D-DCT4006	Tapered Drill Ø4.0 x 6mm	41	I-DG-DC	Depth Gauge
18	D-DCT4008	Tapered Drill Ø4.0 x 8mm	42	I-RATCHET-2	Ratchet Wrench
19	D-DCT4009	Tapered Drill Ø4.0 x 9mm	43	I-SP-X	Fixture Mount Spanner
20	D-DCT4011	Tapered Drill Ø4.0 x 11mm	44	I-HHD-22S / M / L	1.22 Hex Driver, Handpiece Insert, Short / Medium / Long
21	D-DCT4013	Tapered Drill Ø4.0 x 13mm	45	I-WI-22S / M / L	1.22 Hex Driver, Wrench Insert, Short / Medium / Long
22	D-DCT4015	Tapered Drill Ø4.0 x 15mm	46	I-WI-A	Abutment Driver, Wrench Insert
23	D-DCT5009	Tapered Drill Ø5.0 x 9mm	47	I-HD-S / M / L	1.22 Hex Driver, Short / Medium / Long
24	D-DCT5011	Tapered Drill Ø5.0 x 11mm	48	I-AD	Abutment Driver

Images are for illustration purposes only and do not necessarily accurately represent the product.

For contact information on your nearest distributor, visit

[www.southernimplants.com](http://www.southernimplants.com)