


External Hex Range



The above CIA abutments, for Externally Hexed Implants, all use 3 Series prosthetic screws.


Tri-NEX Range



Uses TS-L-18 / GS-L-18 /
TS-L-18C prosthetic screws


Uses TS-L-20 / GS-L-20 / TS-L-20C prosthetic screws

Internal Octagon (IT) Range



Uses TSIT2 / GSIT2 prosthetic screws

Internal Hex Range (M-Series / PROVATA)



Uses TS-Z-18 / GS-Q-18 prosthetic screws

 **Note: Not to be used with Cerec 1 and Cerec 2**

Images are for illustration purposes only and do not necessarily accurately represent the product. Not all products are cleared for sale in all countries.

www.southernimplants.com