


I-TRI-NEX-EG TRI-NEX Instrument Tray



Pilot Drills
D-3Spade-1.8M
D-RB-MS
D-12T-M15
D-16-T

Ø2.0mm Twist Drills
D-20T-M10
D-20T-M15

Dedicated Twist Drills
Ø2.5mm D-25T-M15 D-25T-M20
Ø2.9mm D-29T-M15 D-29T-M20
Ø3.0mm D-30T-M10 D-30T-M15
Ø3.3mm D-33T-M10 D-33T-M15
Ø4.0mm D-40T-M10 D-40T-M15

Universal 1.22/1.27 Hex Drivers
I-UGI-S/ML Handheld
I-HUG-S/ML Handpiece Insert
I-WI-UG-S/ML Wrench Insert

Counter Bore
D-CB

Counter Sink
D-CSS-M (Ø4mm)
D-CSS-5 (Ø5mm)

Compact Conical Abutment Drivers
I-AD Handheld
I-HAD Handpiece Insert
I-WI-A Wrench Insert

Primary Drills Ø3.5 Ø4.3 Ø5.0
Additional Drills Ø2.5 Ø3.5 Ø4.3
Drivers
Placement Tools
Drill / Implant length measure

INTERMEDIATE TAPERED DRILLS
Ø3.5 D-L-35-8 D-L-35-10 D-L-35-11.5 D-L-35-13 D-L-35-16
Ø4.3 D-L-43-8 D-L-43-10 D-L-43-11.5 D-L-43-13 D-L-43-16
Ø5.0 D-L-50-8 D-L-50-10 D-L-50-11.5 D-L-50-13 D-L-50-16
Ø6.0 D-L-60-8 D-L-60-10 D-L-60-11.5 D-L-60-13 D-L-60-16


PLACEMENT TOOLS
Ø3.5 I-HLH-35S / M (Handpiece insert) I-HLHU-35S / M (Universal) I-WI-LH-35S / M (Wrench insert)
Ø4.3 I-HLH-43S / M (Handpiece insert) I-HLHU-43S / M (Universal) I-WI-LH-43S / M (Wrench insert)
Ø5.0 I-HLH-50S / M (Handpiece insert) I-HLHU-S / M (Universal) I-WI-LH-50S / M (Wrench insert)
Ø6.0 I-HLH-60S / M (Handpiece insert) I-HLHU-S / M (Universal) I-WI-LH-60 (Wrench insert)

PLACEMENT TOOLS Co-Axis 12°
Ø4.3 I-L43-12d (Handpiece insert) I-WI-L43-12d (Wrench insert)
Ø5.0 I-L50-12d (Handpiece insert) I-WI-L50-12d (Wrench insert)

Images are for illustration purposes only and do not necessarily accurately represent the product. Not all products are cleared for sale in all countries.

SOUTHERNIMPLANTS.COM

I-TRI-NEX-EG TRI-NEX Instrument Tray



NOTE:

- The surgical kit has an intuitive layout to guide the surgeon through the drill sequence.
- Most Instruments available in various lengths.
- All instruments and tooling used during the procedure must be maintained in good condition, and cleaned and sterilized prior to use. Please consult the Southern Implants Cleaning and Sterilization Procedure Guidelines (CAT-1039) for guidance concerning the maintenance of drills, instruments, and surgical trays.