

Ferretti Healing Abutments

The Southern Implants special Ferretti Healing Abutments are available in the following variations:


TRI-NEX



Implant Diameter	Code	Major Diameter (mm)	Collar Height (mm)	Bone Mill	Bone Mill Cap
Ø4.3	S-HA-L-43-F4	6.5	4	I-BM-67	S-I-BM-L43CAP
	S-HA-L-43-F6	6.5	6	I-BM-67	S-I-BM-L43CAP
	S-HA-L-43W-F4	8	4	I-BAM-77	S-I-BAM-L43CAP
	S-HA-L-43W-F6	8	6	I-BAM-77	S-I-BAM-L43CAP



Implant Diameter	Code	Major Diameter (mm)	Collar Height (mm)	Bone Mill	Bone Mill Cap
Ø5.0	S-HA-L-50-F4	6.5	4	I-BAM-62/77	S-I-BAM-L50CAP
	S-HA-L-50-F6	6.5	6	I-BAM-62/77	S-I-BAM-L50CAP
	S-HA-L-50W-F4	8	4	I-BAM-77	S-I-BAM-L50CAP
	S-HA-L-50W-F6	8	6	I-BAM-77	S-I-BAM-L50CAP



Implant Diameter	Code	Major Diameter (mm)	Collar Height (mm)	Bone Mill	Bone Mill Cap
Ø6.0	S-HA-L-60-F4	6.5	4	I-BBBM-77	S-I-BBBM-L60CAP
	S-HA-L-60-F6	6.5	6	I-BBBM-77	S-I-BBBM-L60CAP
	S-HA-L-60W-F4	8	4	I-BBBM-77	S-I-BBBM-L60CAP
	S-HA-L-60W-F6	8	6	I-BBBM-77	S-I-BBBM-L60CAP

Images are for illustration purposes only and do not necessarily accurately represent the product. Not all products are cleared for sale in all countries.