






External Hex Titanium Abutments - used with 2 Series screws

These are premanufactured dental implant abutments that can either be connected direct to an endosseous implant or connect the prosthesis to a compact conical abutment for use as an aid in temporary prosthetic rehabilitation.

The Titanium abutments are provided sterile, however, it will no longer be sterile after modification.

Refer to the Instructions For Use for further information (CAT-8007 - Temporary Titanium Abutments).

<p>Ø3.0 mm interface</p> <p>TCP1H (engaging)</p> <p>TCP1NH (non-engaging)</p> 	<p>Ø3.25 mm interface</p> <p>TCBN1H (engaging)</p> <p>TCBN5H (engaging)</p> <p>TCBN1NH (non-engaging)</p> <p>TCBN5NH (non-engaging)</p> 	<p>Ø4.0 mm interface</p> <p>TCB1H (engaging)</p> <p>TCB5H (engaging)</p> <p>TCB1NH (non-engaging)</p> <p>TCB5NH (non-engaging)</p> 	
<p>Ø5.0 mm interface</p> <p>TCBA1H (engaging)</p> <p>TCBA5H (engaging)</p> <p>TCBA1NH (non-engaging)</p> <p>TCBA5NH (non-engaging)</p> 		<p>Ø6.0 mm interface</p> <p>TCBBB1H (engaging)</p> <p>TCBBB5H (engaging)</p> <p>TCBBB1NH (non-engaging)</p> <p>TCBBB5NH (non-engaging)</p> 	

2 Series screws (M2)

1.22 Hex



TSHZ2



BSH2*

QUAD



GSQZ2

SLOTTED



TSSZ2



GSSZ2



BSS2*

UNIGRIP



TSUZ2



GSUZ2

TORQUE:
32 - 40 Ncm
Head diameter:
2.70 mm

Screw Head Connections

Hex



Slotted



Unigrip



Quad



NOTE:

- due to design revisions screw tips may be flat or rounded.
- always ensure that the correct screw is used for the relevant implant and component.
- refer to CAT-8068 for alternative slotted 1 Series screws.
- * blackened and for laboratory use only.
- universal drivers are compatible with both 1.22 and 1.27 Hex screws:
 - I-HD22U-S/M/L
 - I-HHD-22U-S/M/L
 - I-WI-22U-S-/M/L