



**SOUTHERNIMPLANTS®**

Innovative Treatment Solutions

# OSSEOINTEGRATED FIXTURES

Restoring lives through precision anchorage.


**For Ears**



**For Eyes**



**For Noses**






# RESTORING MORE THAN FUNCTION

## RESTORING IDENTITY

Southern Implants connects  
precision engineering with  
life-changing outcomes.



Images courtesy of Brom Epithetik



### PROSTHETIC VERSATILITY

Compatible with  
mechanical or magnetic  
prosthetics.



### CONTROLLED PLACEMENT

Flanged collar prevents  
over-insertion in critical  
anatomical regions.




### HIGH STRENGTH TITANIUM

High Strength Grade 4  
Pure Titanium ( $\geq 920$  MPa)  
providing exceptional  
fatigue strength.


Available in two shapes




Available in 3 lengths




3 mm



4 mm



6 mm






IMPLANT DIMENSIONS AND INFORMATION

RANGE		MAJOR DIAMETER	PLATFORM DIAMETER (flange)	Prosthetic Colour Code	HEX WIDTH	HEX HEIGHT	COLLAR HEIGHT	THREAD PITCH	APEX DIAMETER	CYLINDRICAL OR TAPERED	IMPLANT LENGTHS		
											3	4	6
IE	Ø3.75 mm	3.75	4.80	●	2.70	0.7	0.40	0.6	2.80	C	✓	✓	✓
IET	Ø4.5 mm	4.5	4.07	●	2.70	0.7	0.25	0.5	2.45	T		✓	
IETi	Ø4.5 mm	4.5	4.07	●	2.70	0.7	0.25	0.5	2.45	T		✓	

**NOTE:** Not all products are cleared for sale in all regions. Please contact your local Southern Implants distributor for availability.



For more information scan below



to download the Osseointegrated catalogue  
or visit southernimplants.com



**SOUTHERNIMPLANTS®**  
Innovative Treatment Solutions

**South Africa – Headquarters** 1 Albert Road, Irene, RSA  
T: +27-12-667-1046 | E: info@southernimplants.com

**EC REP**  
Southern Implants Europe AB Holmgatan 30, S-791 71, Falun, Sweden  
E: ecrep@southernimplants.com

**Subsidiaries**

**Australia**  
Southern Implants Australia  
T: +61-8 -9466-2627  
E: info@southernimplants.com.au

**Spain and Portugal**  
Southern Implants Ibérica  
T: +34-935-053-507  
E: info@southernimplants.es

**United Kingdom and Ireland**  
Southern Implants UK  
T: +44-20-805-94490  
E: info@southernimplants.co.uk

**USA and Canada**  
Southern Implants North America Inc.  
T: +1-561-472-0990  
E: customercare@southernimplants.com

Southern Implants® are distributed world-wide, please visit southernimplants.com for a list of Distributors.  
For the latest and full assortment of Southern Implants' products, visit our current resources library at southernimplants.com.